



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____


PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 3 3S-SS \$9,047

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. DRIVE MOTORS

REQUIRED.

- | | |
|--|--|
| ☐ 6 mph Drive Motors Standard
Part: MOT147433 | ☐ 7 mph, 4 Pole Drive Motors..... \$1,686
Part: MOT1902455
Requires Q-Logic 3 electronics. |
| ☐ 6 mph, Accu-Trac Motor Technology \$1,387
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. | ☐ 7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,387
Part: MOT2008501. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

- | | | | |
|---|--|--|--|
| ☐  Red Corvette
Part: DGN1801101-PRD1 | ☐  Deep Purple
Part: DGN1801101-PPR5 | ☐  Midnight Blue
Part: DGN1801101-PBU10 | ☐  Green Machine
Part: DGN1801101-PGN11 |
| ☐  Back in Black
Part: DGN1801101-PBK1 | ☐  Raspberry Beret
Part: DGN1801101-PPK4 | ☐  Electric Blue
Part: DGN1801101-PBU9 | ☐  Military Green
Part: DGN1801101-PGN8 |
| ☐  Night Train
Part: DGN1801101-PBK3 | ☐  Pretty in Pink
Part: DGN1801101-PPK3 | ☐  Blue Jean
Part: DGN1801101-PBU11 | |
| ☐  Grey Street
Part: DGN1801101-PSL3 | ☐  Evolution Orange
Part: DGN1801101-POR7 | ☐  Lemon Crush
Part: DGN1801101-PYL4 | |


3.1b. Hydrographic Shrouds

- | | | |
|--|---|---|
| ☐  Military Digital Camo..... \$268
Part: DGN1801101-WCM7 | ☐  Faded Glory \$268
Part: DGN1801101-WFL2 | ☐  Carbon Fiber..... \$268
Part: DGN1801101-WCB1 |
|--|---|---|

4. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

- ☐ Switch Accent Lights No Charge
Part: DGN154117
Accent lights operate by pressing included single switch button.
- ☐ Base Mounted, LED, Full Light Package \$1,066
Part: XRF000129
Requires QLNE5 10 Key, Q-Logic 3 EX Joystick for Full Lights, or Q-Logic 3e Joystick.



5. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

- ☐ Quantum Backup Camera \$375
Part: MEC1910365
- 1. Optional Mounting Bracket**
- | | |
|---|---|
| ☐ Backup Camera Headrest Bracket..... \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | ☐ Backup Camera Track Mount Bracket \$28
Part: MEC2104227 |
|---|---|
- 2. Optional Display Mounting Bracket**
- ☐ Backup Camera Q-Logic Joystick Mount Bracket \$107
Part: MEC2106123



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

- | | |
|--|--|
| ☐ 3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair) | ☐ 3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2 |
| ☐ 4"x14" Knobby, Flat Free Black Tires \$592
Part: XRF000368 | ☐ 4"x14" Knobby, Pneumatic Black Tires..... \$592
Part: XRF000367 |

Drive wheels section continued on following page.

6. DRIVE WHEELS CONTINUED**6.2. Drive Wheel Accent Color**

Not available with knobby tires.

- | | | |
|---|--|---|
| ☐ Back in Black.....Standard
Part: DGN164245-PBK1 x2 | ☐ Pretty in PinkNo Charge
Part: DGN164245-PPK3 x2 | ☐ Green MachineNo Charge
Part: DGN164245-PGN11 x2 |
| ☐ Red Corvette.....No Charge
Part: DGN164245-PRD1 x2 | ☐ Evolution Orange.....No Charge
Part: DGN164245-POR7 x2 | ☐ Military GreenNo Charge
Part: DGN164245-PGN8 x2 |
| ☐ Night Train.....No Charge
Part: DGN164245-PBK3 x2 | ☐ Midnight BlueNo Charge
Part: DGN164245-PBU10 x2 | ☐ Military Digital Camo..... \$59 ea.
Part: DGN164245-WCM7 x2 |
| ☐ Grey Street.....No Charge
Part: DGN164245-PSL3 x2 | ☐ Electric BlueNo Charge
Part: DGN164245-PBU9 x2 | ☐ Faded Glory \$59 ea.
Part: DGN164245-WFL2 x2 |
| ☐ Deep Purple.....No Charge
Part: DGN164245-PPR5 x2 | ☐ Blue JeanNo Charge
Part: DGN164245-PBU11 x2 | ☐ Carbon Fiber \$59 ea.
Part: DGN164245-WCB1 x2 |
| ☐ Raspberry Beret.....No Charge
Part: DGN164245-PPK4 x2 | ☐ Lemon Crush.....No Charge
Part: DGN164245-PYL4 x2 | |

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

- | | |
|---|--|
| ☐ NF-22, Gel Cell, Interceptor 2 Battery \$362 ea.
Part: BAT2004376. HCPCS: E2361 | ☐ NF-22, Gel Cell, MK Battery \$441 ea.
Part: BATGEL1008. HCPCS: E2361 |
| ☐ NF-22, AGM, MK Battery..... \$441 ea.
Part: BAT121346. HCPCS: E2361 | |

8. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.

- | | | |
|--|---|---|
| ☐ Up to 200 lb User Weight
Part: PTOINDV3837 | ☐ 201 lb-250 lb User Weight
Part: PTOINDV3861 | ☐ 251 lb-300 lb User Weight
Part: PTOINDV3984 |
|--|---|---|

9. ELECTRONICS**REQUIRED.** Select a type of drive control from the following sections.**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- | | |
|--|---|
| ☐ QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators. | ☐ QLNES 10 Key Joystick\$429
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package. |
|--|---|

**9.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- | | |
|--|---|
| ☐ Q-Logic 3 EX Controller Module \$745
Part: XRF000002. HCPCS: E2377 | ☐ Q-Logic 3 Expandable Controller Harness.....\$467
Part: XRF000264. HCPCS: E2313 |
|--|---|

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

- | | |
|--|---|
| ☐ Q-Logic 3 EX Joystick\$857
Part: CTL167724 | ☐ Q-Logic 3 EX Joystick for Full Lights\$857
Part: CTL167725 |
| ☐ Q-Logic 3e Joystick\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls. | ☐ No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display. |

**9.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

- | |
|--|
| ☐ Magnetic to XLR Charger Adapter with Wall Holder \$61
Part: ELE2103635 |
|--|

Electronics section continued on following page.

9. ELECTRONICS CONTINUED**9.3. Q-Logic 3 EX Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.**

- | | | |
|---|---|--|
| ☐ Black, Soft TouchStandard
Part: DGN124307
Also matches Night Train shrouds. | ☐ Pretty in Pink \$64
Part: DGN124336-PPK3 | ☐ Green Machine \$64
Part: DGN124336-PGN11 |
| ☐ Red Corvette \$64
Part: DGN124336-PRD1 | ☐ Evolution Orange \$64
Part: DGN124336-POR7 | ☐ Military Green \$64
Part: DGN124336-PGN8 |
| ☐ Back in Black \$64
Part: DGN124336-PBK1 | ☐ Midnight Blue \$64
Part: DGN124336-PBU10 | ☐ Military Digital Camo \$64
Part: DGN124336-WCM7 |
| ☐ Grey Street \$64
Part: DGN124336-PSL3 | ☐ Electric Blue \$64
Part: DGN124336-PBU9 | ☐ Faded Glory \$64
Part: DGN124336-WFL2 |
| ☐ Deep Purple \$64
Part: DGN124336-PPR5 | ☐ Blue Jean \$64
Part: DGN124336-PBU11 | ☐ Carbon Fiber \$64
Part: DGN124336-WCB1 |
| ☐ Raspberry Beret \$64
Part: DGN124336-PPK4 | ☐ Lemon Crush \$64
Part: DGN124336-PYL4 | |

10. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|---|--|
| ☐ Mushroom Handle (B) \$155
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$174
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$155
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$179
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$195
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$179
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$155
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$174
Part: ACCASMB1856. HCPCS: E2323 | |

**11. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**

- | |
|--|
| ☐ Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617 |
|--|

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- | |
|---|
| ☐ Swing-away, Q-Logic 3 \$325
Part: ACC1908434. HCPCS: E1032 |
|---|

11.1c. Basic Swing Away Joystick Mount

- | | |
|--|---|
| ☐ Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032 | ☐ Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032 |
|--|---|



11.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

- | | |
|--|--|
| ☐ Retract4 Mount \$375
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. | ☐ Multi-Axis Mount \$129
Part: XRF000004
Allows height adjustment of joystick. |
|--|--|

11.2. Joystick Mount Location

- | | |
|--|---|
| ☐ Left Side Mount
Part: PTOINDV1225 | ☐ Right Side Mount
Part: PTOINDV1226 |
|--|---|



12. TRU BALANCE FLEX SEATING SYSTEM

Select one option from each section. Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.

12.1. TRU Balance Flex Frame Selection☐ TRU Balance Flex Seating System No Charge
Part: PTOINDV3853**12.2. TRU Balance Flex Seat Dimension - Width**

Maximum seat width is 22". Part number is SETWIDTH plus number following # below.

☐ 15" W (#1005) ☐ 16" W (#1006) ☐ 17" W (#1007) ☐ 18" W (#1008) ☐ 19" W (#1009) ☐ 20" W (#1010) ☐ 21" W (#1011) ☐ 22" W (#1012)**12.3. TRU Balance Flex Seat Dimension - Depth**

Maximum seat depth is 20". Part number is SETDEPTH plus number following # below.

☐ 15" D (#1005) ☐ 16" D (#1006) ☐ 17" D (#1007) ☐ 18" D (#1008) ☐ 19" D (#1009) ☐ 20" D (#1010)**12.4. TRU Balance Flex Seat to Floor**☐ 16.5" at front and rear
Part: PTOINDV3459 ☐ 17.5" at front with 2° posterior angle
Part: PTOINDV3460 ☐ 18.25" at front with 4° posterior angle
Part: PTOINDV3805**12.5. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

Cushion Model	Width	Depth
Cushion Top Foam V=Visco Top Foam: 1" thick, XXSoft viscoelastic foam topper. F=FRG Top Foam: Under 22"W, blue XXSoft FRG foam topper. 22"W and over, green MedSoft FRG foam topper. N=No Top Foam	Cushion Liner U=Urethane Liner: 2Way Stretch. P=Polyurethane Liner: 4Way Stretch. N=No Liner	Cushion Cover C=CoolCore: Chemical free, 4Way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. M=Mesh: Superior air circulation, non-slip. T=Stealth Tek: Waterproof with stretch capability.

1. Cushion Type

- | | |
|---|---|
| ☐ TRU-Comfort 2 SPP
Reference part: XRF000116. See cushion order form.
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608 | ☐ Premiere P
Reference part: XRF000122. See cushion order form.
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606 |
| ☐ Essence SPP
Reference part: XRF000118. See cushion order form.
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608 | ☐ Glacial SP
Reference part: XRF000114. See cushion order form.
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604 |
| ☐ Solution SPP
Reference part: XRF000113. See cushion order form.
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608 | ☐ Zen SP
Reference part: XRF000115. See cushion order form.
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604 |
| ☐ Spectrum Foam SPP
Reference part: XRF000112. See cushion order form.
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608 | ☐ Simplicity G
Reference part: XRF000111. See cushion order form.
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602 |
| ☐ Spectrum Gel SPP
Reference part: XRF000117. See cushion order form.
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608 | |

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- | | | |
|---|--|---|
| ☐ Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller | ☐ Medium Wedge Positioning Kit\$72
Part: POP1909006-M
14"W and smaller | ☐ Large Wedge Positioning Kit \$105
Part: POP1909006-L
14"W and smaller |
| ☐ Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger | ☐ Medium Wedge Positioning Kit\$72
Part: POP1712340-M
15"W and larger | ☐ Large Wedge Positioning Kit \$105
Part: POP1712340-L
15"W and larger |

TRU Balance Flex Seating System section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**3. Cushion Top Foam**

- ☐ Visco Top Foam
Reference part: XRF000380
- ☐ No Top Foam
Reference part: XRF000382
- ☐ FRG Top Foam
Reference part: XRF000381

4. Cushion Liner

- ☐ Stealth Urethane Liner
Reference part: XRF000383
- ☐ No Liner
Reference part: XRF000385
- ☐ Stealth Polyurethane Liner
Reference part: XRF000384

5. Cushion Cover

- ☐ CoolCore Cover
Reference part: XRF000386
- ☐ Stealth Tek Cover
Reference part: XRF000388
- ☐ Mesh Cover
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ☐ ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

- ☐ Additional Stealth Urethane Liner
Reference Part: XRF000390
- ☐ Additional Stealth Polyurethane Liner
Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

- ☐ Additional CoolCore Cover
Reference Part: XRF000392
- ☐ Additional Stealth Tek Cover
Reference Part: XRF000394
- ☐ Additional Mesh Cover
Reference Part: XRF000393

12.6. TRU Balance Flex Back Cane Selection

- ☐ Angle Adjustable Back Cane No Charge
Part: MEC140266
- ☐ 4" Endomorph (2" per side) Back Cane \$209
Part: XRF000272
- ☐ Flip Forward Back Cane \$188
Part: XRF000274
Not available with oxygen holder or cane mounted and cantilever armrests.
- ☐ 2" Mesomorph (1" per side) Back Cane \$205
Part: MEC1804736
- ☐ 2" Endomorph (1" per side) Back Cane \$209
Part: XRF000270

12.7. TRU Balance Flex Back Cane Height Selection

16" and 17" back canes are not available with cane mounted, flip up armrests.

- ☐ 14" High Canes
Part: FRMASMB15832
- ☐ 17" High Canes
Part: FRMASMB15835
- ☐ 20" High Canes
Part: FRMASMB15838
- ☐ 23" High Canes
Part: FRMASMB15841
- ☐ 15" High Canes
Part: FRMASMB15833
- ☐ 18" High Canes
Part: FRMASMB15836
- ☐ 21" High Canes
Part: FRMASMB15839
- ☐ 24" High Canes
Part: FRMASMB15842
- ☐ 16" High Canes
Part: FRMASMB15834
- ☐ 19" High Canes
Part: FRMASMB15837
- ☐ 22" High Canes
Part: FRMASMB15840
- ☐ 25" High Canes
Part: FRMASMB15843

12.8. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Sport Back Selection

- ☐ Sport Back Cushion, 21" Wide or Less \$402
Part: XRF000361
- ☐ Sport Back Cushion, 22" Wide or More \$616
Part: XRF000362

1b. TRU-Comfort Back

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

- ☐ TRU-Comfort Back, Back Cane Option, Less than 22" Wide \$759
Part: TC4V10W21WBC. HCPCS: E2620
- ☐ TRU-Comfort Back, Back Cane Option, 22" Wide or Greater ... \$872
Part: TC4V22W32WBC. HCPCS: E2621

TRU Balance Flex Seating System section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1c. Q-Back 2.0, AL Series, Active Contour**

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$568
Part: QB2AL14AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$568
Part: QB2AL18AT. HCPCS: E2613 |
| ☐ Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$568
Part: QB2AL14AT. HCPCS: E2613 | ☐ Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$568
Part: QB2AL18AX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$568
Part: QB2AL16AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703
Part: QB2AL18AXX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$568
Part: QB2AL16AT. HCPCS: E2613 | ☐ Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$568
Part: QB2AL20AS. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$568
Part: QB2AL16AX. HCPCS: E2613 | ☐ Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$568
Part: QB2AL20AT. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
Part: QB2AL16AXX. HCPCS: E2613 | ☐ Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$568
Part: QB2AL20AX. HCPCS: E2613 |
| ☐ Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$568
Part: QB2AL18AS. HCPCS: E2613 | ☐ Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703
Part: QB2AL20AXX. HCPCS: E2613 |

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL14DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL18DT. HCPCS: E2620 |
| ☐ Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL14DT. HCPCS: E2620 | ☐ Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL18DX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL16DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18DXX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL16DT. HCPCS: E2620 | ☐ Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL20DS. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL16DX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL20DT. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16DXX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL20DX. HCPCS: E2620 |
| ☐ Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL18DS. HCPCS: E2620 | ☐ Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20DXX. HCPCS: E2620 |

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|--|
| ☐ Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
Part: QB2AL16ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
Part: QB2AL16ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620 |

1f. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613 | ☐ AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613 |
|---|---|

1g. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620 | ☐ AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT16DT. HCPCS: E2620 | ☐ AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620 | |

1h. Q-Back 2.0, AFT Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|--|
| ☐ AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
Part: QB2AFT16EXX. HCPCS: E2620 | |
|---|--|

TRU Balance Flex Seating System section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1i. Omit Back Selection**

☐ Omit Back..... No Charge
Part: PTOINDV1287

12.9. Back Width Package

REQUIRED with Sport and TRU-Comfort back selections.

☐ 16"-20" Back Width Package ☐ 18"-22" Back Width Package
Part: XRF000038 Part: XRF000039

12.10. Overall Back Height

☐ 14" Overall Back Height Part: SETHGHT1027	☐ 22" Overall Back Height Part: SETHGHT1035
☐ 15" Overall Back Height Part: SETHGHT1028	☐ 23" Overall Back Height Part: SETHGHT1036
☐ 16" Overall Back Height Part: SETHGHT1029	☐ 24" Overall Back Height Part: SETHGHT1037
☐ 17" Overall Back Height Part: SETHGHT1030	☐ 25" Overall Back Height Part: SETHGHT1038
☐ 18" Overall Back Height Part: SETHGHT1031	☐ 26" Overall Back Height Part: SETHGHT1039
☐ 19" Overall Back Height Part: SETHGHT1032	☐ 27" Overall Back Height Part: SETHGHT1040
☐ 20" Overall Back Height Part: SETHGHT1033	☐ 28" Overall Back Height Part: SETHGHT1041
☐ 21" Overall Back Height Part: SETHGHT1034	☐ 29" Overall Back Height Part: SETHGHT1042

		Back Width										
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	14"	P	T	P	T	-	-	-	-	-	-	-
	15"	P	T	P	T	-	-	-	-	-	-	-
	16"	P	T	P	T	P	T	P	T	P	T	P
	17"	P	T	P	T	P	T	P	T	P	T	P
	18"	P	T	P	T	P	T	P	T	P	T	P
	19"	P	T	P	T	P	T	P	T	P	T	P
	20"	-	-	P	T	P	T	P	T	P	T	P
	21"	-	-	P	T	P	T	P	T	P	T	P
	22"	-	-	P	T	P	T	P	T	P	T	P
	23"	-	-	P	T	P	T	P	T	P	T	P
	24"	-	-	P	T	P	T	P	T	P	T	P
	25"	-	-	P	T	P	T	P	T	P	T	P
	26"	-	-	P	T	P	T	P	T	P	T	P
	27"	-	-	P	T	P	T	P	T	P	T	P
	28"	-	-	P	T	P	T	P	T	P	T	P
	29"	-	-	P	T	P	T	P	T	P	T	P

Legend: P Sport Back T TRU-Comfort Back

Seat Back Availability Chart

12.11. TRU Balance Flex Armrests

Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests..... No Charge L R Left part: ARM1803253. Right part: ARM156742	☐ ☐ Desk Length, Drop-In Armrests No Charge L R Left part: ARM1803254. Right part: ARM156743
--	---

1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests.....\$64 L R Left part: XRF000371. Right part: XRF000372	☐ ☐ Desk Length, Flip Back Armrests\$64 L R Left part: XRF000373. Right part: XRF000374
---	--

1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$161 ea. L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	☐ ☐ Desk Length, Two-Post Armrests\$161 ea. L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
---	--

1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests\$81 ea. L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973	☐ ☐ Desk Length, Single Post Armrests\$81 ea. L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973
--	--

12.12. TRU Balance Flex Armpads

Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads

☐ Straight, Full Length Left and Right No Charge Part: ARM134275	☐ Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
☐ Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	☐ Straight, Desk Length Left and Right..... No Charge Part: ARM134254

1b. Waterfall Armpads

☐ Waterfall, Full Length Left and Right..... No Charge Part: ARM134285	☐ Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
☐ Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	☐ Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads

☐ Falcon Gel, Full Length Left and Right.....\$317 Part: ARM136438	☐ Falcon Gel, Desk Length Left, Full Length Right..... \$294 Part: ARM136439
☐ Falcon Gel, Full Length Left, Desk Length Right..... \$294 Part: ARM136442	☐ Falcon Gel, Desk Length Left and Right \$270 Part: ARM136437

13. HEADRESTS**Optional.** Select a headrest pad and bracket.**13.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads, Single Point Adjustable**☐ Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)\$179
Part: CP08CCSPM. HCPCS: E0955☐ Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)\$179
Part: CP10CCSPM. HCPCS: E0955☐ Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)\$179
Part: CP14CCSPM. HCPCS: E0955☐ Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)\$179
Part: CP06CCSPM. HCPCS: E0955**2. Headrest Brackets**☐ Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033☐ Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**13.2. Flip-down Headrest Bracket**☐ Stealth Flip-down Headrest Bracket (FDM350).....\$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

13.3. Q Cam Lock**Optional.**☐ Q Cam Lock.....\$51
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

14. LAP BELTS**REQUIRED.** Select one lap belt below.
☐ 38" Lap Belt\$100
Part: ACCBELT1006

 ☐ 50" Lap Belt\$100
Part: ACCBELT1000

 ☐ 60" Lap Belt\$100
Part: ACCBELT1002

 ☐ 70" Lap Belt\$100
Part: ACCBELT1003

 ☐ 80" Lap Belt\$100
Part: ACCBELT1004

 ☐ 90" Lap Belt\$100
Part: ACCBELT1001
15. LEGRESTS**Optional.** Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**15.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

☐ Base Mounted Foot Platform Standard
Part: FRMASMB16077☐ Omit base mounted foot platform No Charge
Part: PTOINDV3674☐ Bell Shaped, Base Mounted Foot Platform \$252
Part: RIGASMB7120041**15.1b. Basic Center Mount Foot Platform****1. Basic Center Mount Foot Platform Assembly**☐ Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates☐ Short Extension, 5"x9" Footplates \$184
Part: XRF000312. HCPCS: K0040☐ Long Extension, 5"x9" Footplates \$184
Part: XRF000314. HCPCS: K0040☐ Medium Extension, 5"x9" Footplates \$184
Part: XRF000313. HCPCS: K0040**3. Basic Center Mount Calf Pads**☐ Small Calf Pads, Basic Center Mount.....\$161
Part: RIG2104369☐ Large Calf Pads, Basic Center Mount\$161
Part: RIG2202997
Only available with Long Extension.**4. Basic Center Mount Additional Accessory**☐ Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following pages.

15. LEGRESTS CONTINUED**15.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in later section. Select legrest length in following section.

1. Foot Platform Assembly

- ☐ Quick Adjustable Center Mount Foot Platform \$268
Part: XRF000311. HCPCS: K0040
- ☐ Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

2. Lower Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

- ☐ Pediatric Extension, 5"x8" Footplates
Part: RIG2404015
- ☐ Long Extension, 5"x11" Footplates
Part: RIG2404019
- ☐ Medium Extension, 5"x8" Footplates
Part: RIG2404016
- ☐ Medium Extension, 5"x11" Footplates
Part: RIG2404018
- ☐ Long Extension, 5"x8" Footplates
Part: RIG2404017

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

- ☐ Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018
- ☐ Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019
- ☐ Small Calf Pads, Quick Adjustable Center Mount Only, \$161
6"Wx6"H (pair)
Part: XRF000156
- ☐ Large Calf Pads, Quick Adjustable Center Mount Only, \$161
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

- ☐ Small Calf Panel, 18"Wx8"H \$188
Part: FRMASMB7110012
- ☐ Large Calf Panel, 24"Wx8"H \$242
Part: RIG123825
Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Assembly required. May require seating configuration adjustments to avoid front caster interference.

- ☐ Width Extension for 8" \$70
Part: RIG2404919
- ☐ Width Extension for 11" \$70
Part: RIG2404920
- ☐ Outside Edge \$70
Part: RIG2404921
Not available with width extensions.
- ☐ Outside Edge for Width Extension \$70
Part: RIG2404922
Available with footplates and width extension.
- ☐ Front and Rear Cap \$16
Part: RIG2404923

5. Center Mount Additional Accessories

- ☐ Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q
- ☐ Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q
- ☐ Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q
- ☐ Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q
- ☐ Footplate Connector \$16
Part: RIGASMB7120023
- ☐ Footplate Riser \$77
Part: XRF000348
Not available with Long Extension. Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments.
- ☐ Footplate Riser Calf Panel \$51
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

- ☐ 5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108
- ☐ 5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

15.1d. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$88 ea.
P L R Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.
- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$88 ea.
P L R Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.

Legrests section continued on following pages.

15. LEGRESTS CONTINUED**2. Lower Extension**☐ ☐ ☐

11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ ☐ ☐

13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ ☐ ☐

15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐

5"Wx8"D Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐

6"Wx11"D¹ Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐

6"Wx8"D¹ Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐

8"Wx11"D² Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

☐

Therafin 31515 Multi-Axis Amputee Support \$259 ea.
 Part: INDPART2822. HCPCS: E1020

☐

Therafin 32387 Amputee Support Gel Cover \$183 ea.
 Part: INDPART2823. HCPCS: K0108

15.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper☐ ☐ ☐

Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension☐ ☐ ☐

11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐

14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

☐

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

☐

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)☐ ☐ ☐

7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐

8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐

5"Wx8"D Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐

6"Wx11"D¹ Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐

6"Wx8"D¹ Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐

8"Wx11"D² Footplates with Heel Loops \$92 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories☐

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$134
 Part: INDPART2242. HCPCS: E0953

15.2. Legrest Accessories☐

Heel Loops (pair) \$54
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

☐

Wheel Bumper \$38
 Part: XRF000100

15.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length☐

2" L
 (#3989)

☐

5" L
 (#3992)

☐

8" L
 (#3680)

☐

9.5" L
 (#3683)

☐

11" L
 (#3686)

☐

12.5" L
 (#3689)

☐

14" L
 (#3692)

☐

15.5" L
 (#3695)

☐

17" L
 (#3698)

☐

18.5" L
 (#3701)

☐

3" L
 (#3990)

☐

6" L
 (#3993)

☐

8.5" L
 (#3681)

☐

10" L
 (#3684)

☐

11.5" L
 (#3687)

☐

13" L
 (#3690)

☐

14.5" L
 (#3693)

☐

16" L
 (#3696)

☐

17.5" L
 (#3699)

☐

19" L
 (#3702)

☐

4" L
 (#3991)

☐

7" L
 (#3994)

☐

9" L
 (#3682)

☐

10.5" L
 (#3685)

☐

12" L
 (#3688)

☐

13.5" L
 (#3691)

☐

15" L
 (#3694)

☐

16.5" L
 (#3697)

☐

18" L
 (#3700)

Legrests section continued on following pages.

15. LEGRESTS CONTINUED**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

- | | | | | | | | | | | |
|--|--|---|---|---|---|---|---|---|---|---|
| ☐ 7" L
(#3470) | ☐ 8.5" L
(#3473) | ☐ 10" L
(#3476) | ☐ 11.5" L
(#3479) | ☐ 13" L
(#3482) | ☐ 14.5" L
(#3485) | ☐ 16" L
(#3488) | ☐ 17.5" L
(#3491) | ☐ 19" L
(#3494) | ☐ 20.5" L
(#3497) | ☐ 22" L
(#3669) |
| ☐ 7.5" L
(#3471) | ☐ 9" L
(#3474) | ☐ 10.5" L
(#3477) | ☐ 12" L
(#3480) | ☐ 13.5" L
(#3483) | ☐ 15" L
(#3486) | ☐ 16.5" L
(#3489) | ☐ 18" L
(#3492) | ☐ 19.5" L
(#3495) | ☐ 21" L
(#3498) | ☐ 22.5" L
(#3670) |
| ☐ 8" L
(#3472) | ☐ 9.5" L
(#3475) | ☐ 11" L
(#3478) | ☐ 12.5" L
(#3481) | ☐ 14" L
(#3484) | ☐ 15.5" L
(#3487) | ☐ 17" L
(#3490) | ☐ 18.5" L
(#3493) | ☐ 20" L
(#3496) | ☐ 21.5" L
(#3600) | |

16. POWER OPTION ELECTRONICS**REQUIRED for articulating foot platform.**

- | | |
|---|---|
| ☐ Single Actuator Through Controller Interface\$2,206
Part: XRF000335. HCPCS: E2310 | ☐ Single Actuator Through Controller and Interface (AAM)..... \$3,101
Part: XRF000068. HCPCS: E2310 |
|---|---|

17. SEAT MOUNTED LED LIGHTS**Optional.** Requires QLNE5 10 key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

- | |
|--|
| ☐ Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130
Not available with Base Mounted, LED, Full Light Package. |
|--|

18. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

- | | |
|--|---|
| ☐ Flex Seat Cup Holder No Charge
Part: ACCASMB2634
Not available on TRU Balance Flex seating without two-post armrests. | ☐ Quantum Backpack\$171
Part: ACC1905853 |
| ☐ Oxygen Holder for Flex Seat \$214
Part: ACC125006. HCPCS: E2208
Not available with Manual Recline. | ☐ Grip Drink Holder (20085) \$15.50
Part: ACC152993 |
| ☐ Walker Holder \$111
Flex part: FRMASMB7031 | ☐ Accessory Bar\$135
Part: XRF000128 |
| ☐ Cane/Crutch Holder \$107
Flex part: ACC152515. HCPCS: E2207 | ☐ Large, Black, Weather Cover \$202
Part: ACCCOVR1021 |
| ☐ Flag Holder \$55
Flex part: ACCASMB2415 | ☐ Privacy Bag\$43
Part: ACC2005723 |
| ☐ Saddle Bag\$40
Part: DWR1020D010 | ☐ Transfer Handles (pair) \$214
Part: XRF000396. HCPCS: E0705 |
| ☐ Stealth Cell Phone/Media Device Holder \$255
Part: ACC7120006 | ☐ Personal Item Pouch\$86
Part: ACC2101017 |
| ☐ Stealth Large Tablet Holder \$625.28
Part: ACC141818 | ☐ Remote Stop Switch \$1,321
Part: XRF000379 |

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

- | | |
|---|--|
| ☐ Aftermarket Accessory Mounting Bracket, Slotted \$54 ea.
Part: ACC2407046 | ☐ Aftermarket Accessory Mounting Bracket, Non-slotted .. \$54 ea.
Part: ACC2405496 |
|---|--|

19. BRANCH PATCHES**Optional.**

- | | | |
|---|---|--|
| ☐ Air Force Patch\$21
Part: ACC140240 | ☐ Coast Guard Patch\$21
Part: ACC140239 | ☐ Navy Patch\$21
Part: ACC140244 |
| ☐ Army Patch\$21
Part: ACC140241 | ☐ Marines Patch\$21
Part: ACC140243 | ☐ Purple Heart Patch\$21
Part: ACC140238 |



20. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**20.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

- ☐ Swing-away, UniLink, Back Cane, pair (ULPK404104B00)\$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034

2. Pads

- ☐ Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250
Part: INDPART2663. HCPCS: E0956
- ☐ Flat Pads, 3"x5", pair (ULP846050320)..... \$250
Part: ST-ULP846050320-Q. HCPCS: E0956
- ☐ Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
Part: INDPART2664. HCPCS: E0956
- ☐ Flat Pads, 4"x6", pair (ULP835050320)..... \$250
Part: ST-ULP835050320-Q. HCPCS: E0956

3. Installation

- ☐ Install Laterals
Part: PTOINDV3499
- ☐ Do Not Install
Part: PTOINDV3500

20.2. Stealth Unilink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

- ☐ UniLink, 5" Height Adjustable, \$265 ea.
Removable Mounting Hardware
Part: ST-UTTP5-Q. HCPCS: HCPCS: E1028/E1034^
- ☐ Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: HCPCS: E1028/E1034^
- ☐ UniLink, 7" Height Adjustable, \$276 ea.
Removable Mounting Hardware
Part: ST-UTTP7-Q. HCPCS: HCPCS: E1028/E1034^
- ☐ Flip Down with RIGHT Full Surface Contact..... \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: HCPCS: E1028/E1034^
- ☐ UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: HCPCS: E1028/E1034^
- ☐ UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

- ☐ UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^
- ☐ UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
- ☐ UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

- ☐ UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q
- ☐ UniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

- ☐ Install on Left
Part: PTOINDV3585
- ☐ Install on Right
Part: PTOINDV3586
- ☐ Do Not Install
Part: PTOINDV3719

20.3. Stealth Positioning Chest Harness

- ☐ Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960
- ☐ Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960
- ☐ Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960
- ☐ Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960

20.4. Stealth Positioning Chest Straps

- ☐ Chest Strap, Medium w/ Belt Anchor (IFCS4113) \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960
- ☐ Chest Strap, Large w/ Belt Anchor (IFCS4114)..... \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960

20.5. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

- ☐ UniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028
- ☐ UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028




Positioning components section continued on following page.

20. POSITIONING COMPONENTS CONTINUED

2. Facial Lateral Pads and Accessories

- ☐ Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815
- ☐ Facial Pad, 4"x5x2", Cool Core (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108
- ☐ Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea.
Part: ST-920M-Q. HCPCS: K0108
- ☐ Swing-Away Joystick Mount,..... \$276 ea.
18" Rod (C)
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest



20.6. Stealth Elbow Support

1. Elbow Support Bracket

- ☐ Left, Swing-Away, Elbow Support Bracket\$199 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028
- ☐ Right, Swing-Away, Elbow Support Bracket.....\$199 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

- ☐ 3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108
- ☐ 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108

20.7. Stealth Foot Boxes

- ☐ Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954
- ☐ Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954
- ☐ Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954



20.8. Stealth Positioning Shoe Holder

- ☐ Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952
- ☐ Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952
- ☐ Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

21. LAP TRAY

21.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- ☐ Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950
- ☐ Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950
- ☐ Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

- ☐ Swing Away Mounting Hardware.....\$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

Compatible

Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.