



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked For: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

☐ Quantum Q6 Edge 3 Stretto 3S-SS..... \$9,047

Part: Q6 EDGE 3 STRETTO 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

☐ Accu-Trac Motor Technology \$1,387

Part: XRF000268. HCPCS: K0108

3. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**3.1a. Standard Shrouds**

- | | | | |
|--|---|---|---|
| ☐  Red Corvette
Part: DGN1900199-PRD1 | ☐  Deep Purple
Part: DGN1900199-PPR5 | ☐  Midnight Blue
Part: DGN1900199-PBU10 | ☐  Green Machine
Part: DGN1900199-PGN11 |
| ☐  Back in Black
Part: DGN1900199-PBK1 | ☐  Raspberry Beret
Part: DGN1900199-PPK4 | ☐  Electric Blue
Part: DGN1900199-PBU9 | ☐  Military Green
Part: DGN1900199-PGN8 |
| ☐  Night Train
Part: DGN1900199-PBK3 | ☐  Pretty in Pink
Part: DGN1900199-PPK3 | ☐  Blue Jean
Part: DGN1900199-PBU11 | |
| ☐  Grey Street
Part: DGN1900199-PSL3 | ☐  Evolution Orange
Part: DGN1900199-POR7 | ☐  Lemon Crush
Part: DGN1900199-PYL4 | |

3.1b. Hydrographic Shrouds

- | | |
|--|--|
| ☐  Military Digital Camo \$268
Part: DGN1801101-WCM7 | ☐  Faded Glory \$268
Part: DGN1801101-WFL2 |
|--|--|

4. BASE MOUNTED LIGHTS**Optional.** Select one option.

- ☐ **Switch Accent Lights** No Charge
Part: DGN1810041
Accent lights operate by pressing included single switch button. Accent lighting is not available with Single Actuator Control through Toggle Box.

**5. QUANTUM BACKUP CAMERA****Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

- ☐ **Quantum Backup Camera** \$375
Part: MEC1910365

1. Optional Camera Mounting Bracket

- | | |
|---|--|
| ☐ Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | ☐ Backup Camera Track Mount Bracket \$28
Part: MEC2104227 |
|---|--|

2. Optional Display Mounting Bracket

- ☐ **Backup Camera Q-Logic Joystick Mount Bracket** \$107
Part: MEC2106123

**6. DRIVE WHEELS****REQUIRED.** Select a type of drive wheel. Accent colors are available for 14" tires.**6.1. Drive Wheels**

- | | |
|---|---|
| ☐ 3"x14" Pneumatic Tires No Charge
Part: WHL1713925 | ☐ 2"x12.5" Pneumatic Tires No Charge
Part: WHL1810666 |
| ☐ 3"x14" Flat Free Tires No Charge
Part: WHL1704058 | ☐ 2"x12.5" Flat Free Tires No Charge
Part: WHL1810665 |
| ☐ 4"x14" Knobby, Flat Free Black Tires \$592
Part: XRF000368 | ☐ 4"x14" Knobby, Pneumatic Black Tires \$592
Part: XRF000367 |

6.2. Drive Wheel Accent Color for 14" Tires

Not available with 12.5" tires.

- | | | | |
|--|---|---|--|
| ☐ Back in Black No Charge
Part: DGN164245-PBK1 x2 | ☐ Deep Purple No Charge
Part: DGN164245-PPR5 x2 | ☐ Midnight Blue No Charge
Part: DGN164245-PBU10 x2 | ☐ Green Machine No Charge
Part: DGN164245-PGN11 x2 |
| ☐ Red Corvette No Charge
Part: DGN164245-PRD1 x2 | ☐ Raspberry Beret No Charge
Part: DGN164245-PPK4 x2 | ☐ Electric Blue No Charge
Part: DGN164245-PBU9 x2 | ☐ Military Green No Charge
Part: DGN164245-PGN8 x2 |
| ☐ Night Train No Charge
Part: DGN164245-PBK3 x2 | ☐ Pretty in Pink No Charge
Part: DGN164245-PPK3 x2 | ☐ Blue Jean No Charge
Part: DGN164245-PBU11 x2 | ☐ Military Digital Camo \$59 ea.
Part: DGN164245-WCM7 x2 |
| ☐ Grey Street No Charge
Part: DGN164245-PSL3 x2 | ☐ Evolution Orange No Charge
Part: DGN164245-POR7 x2 | ☐ Lemon Crush No Charge
Part: DGN164245-PYL4 x2 | ☐ Faded Glory \$59 ea.
Part: DGN164245-WFL2 x2 |

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

- ☐ NF-22, Gel Cell, Interceptor 2 Battery \$362 ea.
Part: BAT2004376. HCPCS: E2361
Uses 8 AMP charger
- ☐ NF-22, AGM, MK Flat Top Battery \$441 ea.
Part: BAT2003982. HCPCS: E2361
Uses 8 AMP charger
- ☐ Gel Cell, 40AH, Interceptor 2 Battery.....\$195 ea.
Part: BAT2100316. HCPCS: K0108
Uses 5 AMP charger

8. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.

- ☐ Up to 150 lb User Weight
Part: PTOINDV3977
- ☐ 151-250 lb User Weight
Part: PTOINDV3985
- ☐ 250-300 lb User Weight
Part: PTOINDV3984

9. ELECTRONICS**REQUIRED.** Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).**9.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- ☐ QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators.
- ☐ QLNES 10 Key Joystick\$429
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**9.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

- ☐ Q-Logic 3 EX Controller Module \$745
Part: XRF000002. HCPCS: E2377
- ☐ Q-Logic 3 Expandable Controller Harness.....\$467
Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

- ☐ Q-Logic 3 EX Joystick\$857
Part: CTL167724
- ☐ Q-Logic 3 EX Joystick for Full Lights\$857
Part: CTL167725
- ☐ Q-Logic 3e Joystick\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.
- ☐ No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display

**9.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.


- ☐ Magnetic to XLR Charger Adapter with Wall Holder \$61
Part: ELE2103635

9.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

- | | | |
|--|---|--|
| ☐ Black, Soft TouchStandard
Part: DGN124307
Also matches Night Train shrouds. | ☐ Raspberry Beret..... \$64
Part: DGN124336-PPK4 | ☐ Lemon Crush..... \$64
Part: DGN124336-PYL4 |
| ☐ Red Corvette..... \$64
Part: DGN124336-PRD1 | ☐ Pretty in Pink \$64
Part: DGN124336-PPK3 | ☐ Green Machine \$64
Part: DGN124336-PGN11 |
| ☐ Back in Black..... \$64
Part: DGN124336-PBK1 | ☐ Evolution Orange..... \$64
Part: DGN124336-POR7 | ☐ Military Green \$64
Part: DGN124336-PGN8 |
| ☐ Grey Street \$64
Part: DGN124336-PSL3 | ☐ Midnight Blue \$64
Part: DGN124336-PBU10 | ☐ Military Digital Camo..... \$64
Part: DGN124336-WCM7 |
| ☐ Deep Purple \$64
Part: DGN124336-PPR5 | ☐ Electric Blue \$64
Part: DGN124336-PBU9 | ☐ Faded Glory \$64
Part: DGN124336-WFL2 |
| | ☐ Blue Jean \$64
Part: DGN124336-PBU11 | |

10. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|--|---|
| ☐ Mushroom Handle (B) \$155
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$174
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$155
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$179
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$195
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$179
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$155
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$174
Part: ACCASMB1856. HCPCS: E2323 | |

**11. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**

- ☐
- Fixed, In-line, Bracket for Q-Logic 3 Standard
-
- Part: FRMASMB10617

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- ☐
- Swing-away, Q-Logic 3 \$325
-
- Part: ACC1908434. HCPCS: E1032

11.1c. Basic Swing Away Joystick Mount

- ☐
- Basic Swing Away, Q-Logic, Left Side \$245
-
- Part: ACC2201041. HCPCS: E1032

- ☐
- Basic Swing Away, Q-Logic, Right Side \$245
-
- Part: ACC2201042. HCPCS: E1032



11.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

- | | |
|---|---|
| ☐ Retract4 Mount \$375
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. | ☐ Multi-Axis Mount \$129
Part: XRF000004
Allows height adjustment of joystick. |
|---|---|

11.2. Joystick Mount Location

- | | |
|---|--|
| ☐ Left Side Mount
Part: PTOINDV1225 | ☐ Right Side Mount
Part: PTOINDV1226 |
|---|--|

**12. TRU BALANCE FLEX SEATING SYSTEM**

Seat width and depth can be adjusted in 1" increments from 15"W-20"W and 16"D-20"D. Maximum seat size 20"W x 20"D.

12.1. TRU Balance Flex Frame Selection**REQUIRED.**

- ☐
- TRU Balance Flex Seating System No Charge
-
- Part: PTOINDV3853

12.2. Occupied Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Omitting a backrest selection requires a provider to facilitate an aftermarket backrest that best fits the consumer's needs. The third-party backrest if applicable must be occupied transit approved per ISO/RESNA standards. When modifying a power chair, including the addition of third-party componentry, to accommodate the medical needs of the power chair occupant, a risk assessment in conformance with ISO 14971, as outlined in ISO 7176, should be performed.

- ☐
- WC-19 Occupied Transit Option \$295
-
- Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**12.3. Seat Width****REQUIRED.** Maximum seat width is 20". Part number is SETWDTH plus number following # below.

- | | | | | | | | | | | |
|---|---|---|---|---|---|---|---|---|---|---|
| ☐ 10" W
(#1000) | ☐ 11" W
(#1001) | ☐ 12" W
(#1002) | ☐ 13" W
(#1003) | ☐ 14" W
(#1004) | ☐ 15" W
(#1005) | ☐ 16" W
(#1006) | ☐ 17" W
(#1007) | ☐ 18" W
(#1008) | ☐ 19" W
(#1009) | ☐ 20" W
(#1010) |
|---|---|---|---|---|---|---|---|---|---|---|

12.4. Seat Depth**REQUIRED.** Maximum seat depth is 20". Part number is SETDPH plus number following # below.

- | | | | | | |
|---|---|---|---|---|---|
| ☐ 10" D
(#1000) | ☐ 12" D
(#1002) | ☐ 14" D
(#1004) | ☐ 16" D
(#1006) | ☐ 18" D
(#1008) | ☐ 20" D
(#1010) |
| ☐ 11" D
(#1001) | ☐ 13" D
(#1003) | ☐ 15" D
(#1005) | ☐ 17" D
(#1007) | ☐ 19" D
(#1009) | |

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**12.5. Seat to Floor Height**☐ **Lowest Seat to Floor**
Part: PTOINDV3558☐ **Middle Seat to Floor**
Part: PTOINDV3559☐ **Highest Seat to Floor**
Part: PTOINDV3560

High	18.5" at front with 4" posterior angle
Med	17.5" at front with 2" posterior angle
TRU Balance Flex Seat to Floor Heights	

12.6. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

Cushion Model	Width	Depth
Cushion Top Foam V=Visco Top Foam: 1" thick, XXSoft viscoelastic foam topper. F=FRG Top Foam: Under 22"W, blue XXSoft FRG foam topper. 22"W and over, green MedSoft FRG foam topper. N=No Top Foam	Cushion Liner U=Urethane Liner: 2Way Stretch. P=Polyurethane Liner: 4Way Stretch. N=No Liner	Cushion Cover C=CoolCore: Chemical free, 4Way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. M=Mesh: Superior air circulation, non-slip. T=Stealth Tek: Waterproof with stretch capability.

1. Cushion Type
☐ **TRU-Comfort 2 SPP**
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

☐ **Essence SPP**
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

☐ **Solution SPP**
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

☐ **Spectrum Foam SPP**
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

☐ **Spectrum Gel SPP**
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

☐ **Premiere P**
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

☐ **Glacial SP**
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

☐ **Zen SP**
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

☐ **Simplicity G**
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602
2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

☐ **Small Wedge Positioning Kit**\$39
Part: POP1909006-S
14"W and smaller

☐ **Medium Wedge Positioning Kit**\$72
Part: POP1909006-M
14"W and smaller

☐ **Large Wedge Positioning Kit**\$105
Part: POP1909006-L
14"W and smaller

☐ **Small Wedge Positioning Kit**\$39
Part: POP1712340-S
15"W and larger

☐ **Medium Wedge Positioning Kit**\$72
Part: POP1712340-M
15"W and larger

☐ **Large Wedge Positioning Kit**\$105
Part: POP1712340-L
15"W and larger
3. Cushion Top Foam
☐ **Visco Top Foam**
Reference part: XRF000380

☐ **No Top Foam**
Reference part: XRF000382

☐ **FRG Top Foam**
Reference part: XRF000381
4. Cushion Liner
☐ **Stealth Urethane Liner**
Reference part: XRF000383

☐ **No Liner**
Reference part: XRF000385

☐ **Stealth Polyurethane Liner**
Reference part: XRF000384
5. Cushion Cover
☐ **CoolCore Cover**
Reference part: XRF000386

☐ **Stealth Tek Cover**
Reference part: XRF000387

☐ **Mesh Cover**
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	—
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**6. Optional Base Insert**

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ☐ **ABS Flat Rigidizer**
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

- ☐ **Additional Stealth Urethane Liner** ☐ **Additional Stealth Polyurethane Liner**
Reference Part: XRF000390 Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

- ☐ **Additional CoolCore Cover** ☐ **Additional Stealth Tek Cover**
Reference Part: XRF000392 Reference Part: XRF000394
- ☐ **Additional Mesh Cover**
Reference Part: XRF000393

12.7. TRU Balance Flex Back Cane Selection**REQUIRED.**

- ☐ **Angle Adjustable Back Cane** No Charge
Part: MEC140266
- ☐ **Flip Forward Back Cane** \$188
Part: XRF000274
Not available with oxygen holder, or the flip back, cane mount armrest.
- ☐ **2" Mesomorph (1" per side) Back Cane** \$205
Part: MEC1804736
- ☐ **2" Endomorph (1" per side) Back Cane** \$209
Part: XRF000270
- ☐ **4" Endomorph (2" per side) Back Cane** \$209
Part: XRF000272

12.8. TRU Balance Flex Back Cane Height Selection**REQUIRED.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

- ☐ **14" High Canes** Part: FRMASMB15832
- ☐ **15" High Canes** Part: FRMASMB15833
- ☐ **16" High Canes** Part: FRMASMB15834
- ☐ **17" High Canes** Part: FRMASMB15835
- ☐ **18" High Canes** Part: FRMASMB15836
- ☐ **19" High Canes** Part: FRMASMB15837
- ☐ **20" High Canes** Part: FRMASMB15838
- ☐ **21" High Canes** Part: FRMASMB15839
- ☐ **22" High Canes** Part: FRMASMB15840
- ☐ **23" High Canes** Part: FRMASMB15841
- ☐ **24" High Canes** Part: FRMASMB15842
- ☐ **25" High Canes** Part: FRMASMB15843

12.9. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Sport Back Selection

- ☐ **Sport Back Cushion, 21" Wide or Less** \$402
Part: XRF000361
- ☐ **Sport Back Cushion, 22" Wide or More** \$616
Part: XRF000362

1b. TRU-Comfort Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

- ☐ **TRU-Comfort Back, Back Cane Option, Less than 22" Wide** \$759
Part: TC4V10W21WBC. HCPCS: E2620
- ☐ **TRU-Comfort Back, Back Cane Option, 22" Wide or Greater** ... \$872
Part: TC4V22W32WBC. HCPCS: E2621

		Back Width										
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	14"	P	T	P	T	P	T	-	-	-	-	-
	15"	P	T	P	T	P	T	-	-	-	-	-
	16"	P	T	P	T	P	T	P	T	P	T	P
	17"	P	T	P	T	P	T	P	T	P	T	P
	18"	P	P	P	T	P	T	P	T	P	T	P
	19"	P	P	P	T	P	T	P	T	P	T	P
	20"	-	-	P	T	P	T	P	T	P	T	P
	21"	-	-	P	T	P	T	P	T	P	T	P
	22"	-	-	P	T	P	T	P	T	P	T	P
	23"	-	-	P	T	P	T	P	T	P	T	P
	24"	-	-	P	T	P	T	P	T	P	T	P
	25"	-	-	P	T	P	T	P	T	P	T	P
	26"	-	-	P	T	P	T	P	T	P	T	P
	27"	-	-	P	T	P	T	P	T	P	T	P
	28"	-	-	P	T	P	T	P	T	P	T	P
	29"	-	-	P	T	P	T	P	T	P	T	P

☐ Sport Back ☐ TRU-Comfort Back

Seat Back Availability Chart

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1c. Q-Back 2.0, AL Series, Active Contour**

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL14AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL18AT. HCPCS: E2613 |
| ☐ Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL14AT. HCPCS: E2613 | ☐ Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL18AX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL16AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18AXX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL16AT. HCPCS: E2613 | ☐ Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL20AS. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL16AX. HCPCS: E2613 | ☐ Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL20AT. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16AXX. HCPCS: E2613 | ☐ Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL20AX. HCPCS: E2613 |
| ☐ Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL18AS. HCPCS: E2613 | ☐ Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20AXX. HCPCS: E2613 |

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL14DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL18DT. HCPCS: E2620 |
| ☐ Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL14DT. HCPCS: E2620 | ☐ Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL18DX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL16DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18DXX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL16DT. HCPCS: E2620 | ☐ Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL20DS. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL16DX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL20DT. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16DXX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL20DX. HCPCS: E2620 |
| ☐ Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL18DS. HCPCS: E2620 | ☐ Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20DXX. HCPCS: E2620 |

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|--|
| ☐ Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL16ES. HCPCS: E2620 | ☐ Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL16ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620 |
| ☐ Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620 |
| ☐ Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620 | ☐ Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620 |

1f. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613 | ☐ AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613 |
|---|---|

1g. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| ☐ AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT16DS. HCPCS: E2620 | ☐ AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT18DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT16DT. HCPCS: E2620 | ☐ AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620 |
| ☐ AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620 | |

TRU Balance Flex seating section continued on following pages.

12. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1h. Q-Back 2.0, AFT Series, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section.

☐ AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline\$703
 Part: QB2AFT16EXX. HCPCS: E2620
1i. Omit Back Selection
☐ Omit Back.....No Charge
 Part: PTOINDV1287
12.10. Back Width Package**REQUIRED with Sport and TRU-Comfort back selections.**
☐ 16"-20" Back Width Package
 Part: XRF000038

☐ 18"-22" Back Width Package
 Part: XRF000039
12.11. Overall Back Height**REQUIRED with back selection.** Refer to the back availability chart for stock order back size ranges.

☐ 14" Height Part: SETHGHT1027	☐ 18" Height Part: SETHGHT1031	☐ 22" Height Part: SETHGHT1035	☐ 26" Height Part: SETHGHT1039
☐ 15" Height Part: SETHGHT1028	☐ 19" Height Part: SETHGHT1032	☐ 23" Height Part: SETHGHT1036	☐ 27" Height Part: SETHGHT1040
☐ 16" Height Part: SETHGHT1029	☐ 20" Height Part: SETHGHT1033	☐ 24" Height Part: SETHGHT1037	☐ 28" Height Part: SETHGHT1041
☐ 17" Height Part: SETHGHT1030	☐ 21" Height Part: SETHGHT1034	☐ 25" Height Part: SETHGHT1038	☐ 29" Height Part: SETHGHT1042

12.12. Armrests**REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU Balance Flex Seat, Drop-In Armrests**

Includes full length, flat armpads. Available with 14" deep seats and greater.

☐ Full Length, Drop-In Armrests..... No Charge L R Left part: ARM1803253. Right part: ARM156742	☐ Desk Length, Drop-In Armrests No Charge L R Left part: ARM1803254. Right part: ARM156743
---	--

1b. Pediatric, Height Adjustable, Removable Armrests

Available with straight armpads only.

☐ Desk Length, Pediatric Armrests \$81
 L R Left part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973
1c. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ Full Length, Flip Back Armrests.....\$64 L R Left part: XRF000371. Right part: XRF000372	☐ Desk Length, Flip Back Armrests\$64 L R Left part: XRF000373. Right part: XRF000374
--	---

1d. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ Full Length, Two-Post Armrests.....\$161 ea. L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	☐ Desk Length, Two-Post Armrests\$161 ea. L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
--	---

1e. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ Full Length, Single Post Armrests \$81 ea. L R Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973	☐ Desk Length, Single Post Armrests \$81 ea. L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973
--	--

12.13. Armpads**REQUIRED with armrest selection except Drop-In Armrests.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**

☐ Straight, Full Length Left and Right No Charge Part: ARM134275	☐ Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
☐ Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	☐ Straight, Desk Length Left and Right..... No Charge Part: ARM134254

1b. Waterfall Armpads

☐ Waterfall, Full Length Left and Right..... No Charge Part: ARM134285	☐ Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
☐ Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	☐ Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads

☐ Falcon Gel, Full Length Left and Right.....\$317 Part: ARM136438	☐ Falcon Gel, Desk Length Left, Full Length Right..... \$294 Part: ARM136439
☐ Falcon Gel, Full Length Left, Desk Length Right..... \$294 Part: ARM136442	☐ Falcon Gel, Desk Length Left and Right \$270 Part: ARM136437

13. HEADRESTS**Optional.** Select a headrest pad and bracket.**13.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads, Single Point Adjustable**☐ Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)\$179
Part: CP08CCSPM. HCPCS: E0955☐ Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)\$179
Part: CP10CCSPM. HCPCS: E0955☐ Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)\$179
Part: CP14CCSPM. HCPCS: E0955☐ Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)\$179
Part: CP06CCSPM. HCPCS: E0955**2. Headrest Brackets**☐ Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033☐ Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**13.2. Q Cam Lock****Optional.**☐ Q Cam Lock.....\$51
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

14. LAP BELTS**REQUIRED.** Select one lap belt below.☐ 38" Lap Belt
Part: ACCBELT1006☐ 50" Lap Belt
Part: ACCBELT1000☐ 60" Lap Belt
Part: ACCBELT1002☐ 70" Lap Belt
Part: ACCBELT1003☐ 80" Lap Belt
Part: ACCBELT1004☐ 90" Lap Belt
Part: ACCBELT1001**15. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**15.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform Type☐ Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium
8"-10"	10"-14"
Basic Center Mount Extension Range	

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates☐ Short Extension, 5"x9" FootplatesNo Charge
Part: MEC1906675☐ Medium Extension, 5"x9" FootplatesNo Charge
Part: MEC1906680**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

☐ Small Calf Pads, Basic Center Mount.....\$161
Part: RIG2104369**4. Basic Center Mount Additional Accessories**☐ Basic Center Mount Footplate Connector\$16
Part: RIG2104240**15.1b. Center Mount Foot Platform**

Select legrest length in following section.

1. Foot Platform Assembly☐ Quick Adjustable Center Mount Foot Platform\$268
Part: XRF000311. HCPCS: K0040☐ Power Articulating Foot Platform\$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040**2. Extension and Angle Adjustable Footplates**

Pediatric extensions are not compatible with Power Articulating Foot Platform.

☐ Pediatric Extension, 5"x8" Footplates
Part: RIG2404015☐ Long Extension, 5"x8" Footplates
Part: RIG2404017☐ Medium Extension, 5"x8" Footplates
Part: RIG2404016☐ Long Extension, 5"x11" Footplates
Part: RIG2404019☐ Medium Extension, 5"x11" Footplates
Part: RIG2404018

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

☐ Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018☐ Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019☐ Small Calf Pads, Quick Adjustable Center Mount Only,\$161
6"Wx6"H (pair)
Part: XRF000156☐ Large Calf Pads, Quick Adjustable Center Mount Only,\$161
8" taper to 6"Wx8"H (pair)
Part: XRF000157

Legrests seating section continued on following pages.

15. LEGRESTS CONTINUED**3b. Center Mount Calf Panel**

- ☐ Small Calf Panel, 18"Wx8"H..... \$188
Part: FRMASMB7110012
- ☐ Large Calf Panel, 24"Wx8"H \$242
Part: RIG123825
Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

- ☐ Width Extension for 8" \$70
Part: RIG2404919
- ☐ Width Extension for 11" \$70
Part: RIG2404920
- ☐ Outside Edge \$70
Part: RIG2404921
Not available with width extensions.
- ☐ Outside Edge for Width Extension \$70
Part: RIG2404922
Available with footplates and width extension.
- ☐ Front and Rear Cap \$16
Part: RIG2404923

5. Center Mount Additional Accessories

- ☐ Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q
- ☐ Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q
- ☐ Footplate Connector \$16
Part: RIGASMB7120023
- ☐ Footplate Riser \$77
Part: XRF000348
Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments. Not available with long extensions.
- ☐ Footplate Riser Calf Panel \$51
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

- ☐ 5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108
- ☐ 5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108
- ☐ 5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

15.1c. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$88 ea.
☐ ☐ ☐ Pair part: XRF000082. Left: XRF000083. Right: XRF000084
☐ ☐ ☐ Not available with 11"-13" extensions.
- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$88 ea.
☐ ☐ ☐ Pair part: XRF000085. Left: XRF000086. Right: XRF000087
☐ ☐ ☐ For use only with 11"-13" lower extensions.

2. Lower Extension

- ☐ ☐ ☐ 11"-13" Length
☐ ☐ ☐ Pair part: XRF000088.
☐ ☐ ☐ Left: XRF000089.
☐ ☐ ☐ Right: XRF000090.
- ☐ ☐ ☐ 13"-16.5" Length
☐ ☐ ☐ Pair part: XRF000094.
☐ ☐ ☐ Left: XRF000095.
☐ ☐ ☐ Right: XRF000096.
- ☐ ☐ ☐ 15.5"-19.5" Length
☐ ☐ ☐ Pair part: XRF000097.
☐ ☐ ☐ Left: XRF000098.
☐ ☐ ☐ Right: XRF000099.

3. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- ☐ ☐ ☐ 5"Wx8"D with Heel Loops \$92 ea.
☐ ☐ ☐ Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
☐ ☐ ☐ HCPCS: K0040+E095
- ☐ ☐ ☐ 6"Wx11"D¹ with Heel Loops \$92 ea.
☐ ☐ ☐ Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
☐ ☐ ☐ HCPCS: K0040+E0951
- ☐ ☐ ☐ 6"Wx8"D¹ with Heel Loops \$92 ea.
☐ ☐ ☐ Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
☐ ☐ ☐ HCPCS: K0040+E0951
- ☐ ☐ ☐ 8"Wx11"D² with Heel Loops \$92 ea.
☐ ☐ ☐ Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
☐ ☐ ☐ HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- ☐ ☐ ☐ Therafin 31515 Multi-Axis Amputee Support \$259 ea.
☐ ☐ ☐ Part: INDPART2822. HCPCS: E1020
- ☐ ☐ ☐ Therafin 32387 Amputee Support Gel Cover \$183 ea.
☐ ☐ ☐ Part: INDPART2823. HCPCS: K0108

15.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

- ☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
☐ ☐ ☐ Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
☐ ☐ ☐ Select lower extension and footplates

Legrests section continued on following pages.

15. LEGRESTS CONTINUED

2a. Lower Extension

☐ ☐ ☐
P L R

11"-14" Length

Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
For reinforced only

☐ ☐ ☐
P L R

14"-17" Length

Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
For reinforced only

☐ ☐ ☐
P L R

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length

Pair part: XRF000101. Left: XRF000102. Right: XRF000103
For angle adjustable only

☐ ☐ ☐
P L R

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length

Pair part: XRF000104. Left: XRF000105. Right: XRF000106
For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

☐
Left Residual Limb Support with Calf Pad..... \$342
Part: FRMASMB5301. HCPCS: E1020

☐
Right Residual Limb Support with Calf Pad \$342
Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐
P L R

7.25"Wx6"D Reinforced Footplate..... Standard

Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐
P L R

8.5"Wx6"D Reinforced Footplate Standard

Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐
P L R

5"Wx8"D with Heel Loops..... \$92 ea.

Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

☐ ☐ ☐
P L R

6"Wx11"D¹ with Heel Loops..... \$92 ea.

Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

☐ ☐ ☐
P L R

6"Wx8"D¹ with Heel Loops..... \$92 ea.

Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

☐ ☐ ☐
P L R

8"Wx11"D² with Heel Loops..... \$92 ea.

Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

4. Elevating Legrest Accessories

☐
Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$134
Part: INDPART2242. HCPCS: E0953

15.2. Legrest Accessories

☐
Heel Loops (pair)..... \$54
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi axis footplates when noted.

☐
Wheel Bumper..... \$38
Part: XRF000100

15.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

☐ **2" L**
(#3989)

☐ **4" L**
(#3991)

☐ **6" L**
(#3993)

☐ **8" L**
(#3680)

☐ **9" L**
(#3682)

☐ **10" L**
(#3684)

☐ **11" L**
(#3686)

☐ **12" L**
(#3688)

☐ **13" L**
(#3690)

☐ **14" L**
(#3692)

☐ **15" L**
(#3694)

☐ **3" L**
(#3990)

☐ **5" L**
(#3992)

☐ **7" L**
(#3994)

☐ **8.5" L**
(#3681)

☐ **9.5" L**
(#3683)

☐ **10.5" L**
(#3685)

☐ **11.5" L**
(#3687)

☐ **12.5" L**
(#3689)

☐ **13.5" L**
(#3691)

☐ **14.5" L**
(#3693)

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

☐ **7" L**
(#3470)

☐ **8" L**
(#3472)

☐ **9" L**
(#3474)

☐ **10" L**
(#3476)

☐ **11" L**
(#3478)

☐ **12" L**
(#3480)

☐ **13" L**
(#3482)

☐ **14" L**
(#3484)

☐ **15" L**
(#3486)

☐ **7.5" L**
(#3471)

☐ **8.5" L**
(#3473)

☐ **9.5" L**
(#3475)

☐ **10.5" L**
(#3477)

☐ **11.5" L**
(#3479)

☐ **12.5" L**
(#3481)

☐ **13.5" L**
(#3483)

☐ **14.5" L**
(#3485)

16. POWER OPTION ELECTRONICS

REQUIRED for articulating foot platform.
☐
Single Actuator Through Controller Interface \$2,206
Part: XRF000335. HCPCS: E2310

☐
Single Actuator Through Controller and Interface (AAM)..... \$3,101
Part: XRF000068. HCPCS: E2310

17. SEAT MOUNTED LED LIGHTS

Optional. Requires QLNE5 10 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

☐
Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130
Not available with Base Mounted, LED, Full Light Package.

18. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

- | | |
|--|--|
| ☐ Flex Seat Cup Holder No Charge
Part: ACCASMB2634
Not available on TRU Balance Flex seating without two-post armrests. | ☐ Grip Drink Holder (20085) \$15.50
Part: ACC152993 |
| ☐ Oxygen Holder for Flex Seat \$214
Part: ACC125006. HCPCS: E2208 | ☐ Accessory Bar \$135
Part: XRF000128 |
| ☐ Cane/Crutch Holder \$107
Flex part: ACC152515. HCPCS: E2207 | ☐ Large, Black, Weather Cover \$202
Part: ACCCOVR1021 |
| ☐ Cell Phone Holder for 12" Lift Captain Seat \$70
Part: ACC1611300 | ☐ Transfer Handles (pair) \$214
Part: XRF000396. HCPCS: E0705
Not available with Heavy Duty armrests. |
| ☐ Flag Holder \$55
Part: ACCASMB2415 | ☐ Privacy Bag \$43
Part: ACC2005723 |
| ☐ Saddle Bag \$40
Part: DWR1020D010 | ☐ Personal Item Pouch \$86
Part: ACC2101017 |
| ☐ Stealth Cell Phone/Media Device Holder \$255
Part: ACC7120006 | ☐ Remote Stop Switch \$1,321
Part: XRF000379 |
| ☐ Quantum Backpack \$171
Part: ACC1905853 | |

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

- | | |
|---|--|
| ☐ Aftermarket Accessory Mounting Bracket, Slotted \$54 ea.
QTY Part: ACC2407046 | ☐ Aftermarket Accessory Mounting Bracket, Non-slotted .. \$54 ea.
QTY Part: ACC2405496 |
|---|--|

19. BRANCH PATCHES**Optional.**

- | | | |
|--|--|---|
| ☐ Air Force Patch \$21
Part: ACC140240 | ☐ Coast Guard Patch \$21
Part: ACC140239 | ☐ Navy Patch \$21
Part: ACC140244 |
| ☐ Army Patch \$21
Part: ACC140241 | ☐ Marines Patch \$21
Part: ACC140243 | ☐ Purple Heart Patch \$21
Part: ACC140238 |

**20. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**20.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

- | |
|---|
| ☐ Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034 |
|---|

2. Pads

- | | |
|---|---|
| ☐ Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956 | ☐ Flat Pads, 3"x5", pair (ULP846050320) \$250
Part: ST-ULP846050320-Q. HCPCS: E0956 |
| ☐ Curved Pads, 4"x6", pair (TWBL-4X6) \$250
Part: INDPART2664. HCPCS: E0956 | ☐ Flat Pads, 4"x6", pair (ULP835050320) \$250
Part: ST-ULP835050320-Q. HCPCS: E0956 |

3. Installation

- | | |
|--|--|
| ☐ Install Laterals
Part: PTOINDV3499 | ☐ Do Not Install
Part: PTOINDV3500 |
|--|--|

20.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

- | | |
|---|--|
| ☐ UniLink, 5" Height Adjustable, \$265 ea.
QTY Removable Mounting Hardware
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^ | ☐ Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^ |
| ☐ UniLink, 7" Height Adjustable, \$276 ea.
QTY Removable Mounting Hardware
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^ | ☐ Flip Down with RIGHT Full Surface Contact \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^ |
| ☐ UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^ | ☐ UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^ |

Positioning components section continued on following page.

20. POSITIONING COMPONENTS CONTINUED**2. UniLink Pelvic/Thigh Pads**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

- ☐ **UniLink, 4" x 4" Pad** \$106 ea.
 QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^
- ☐ **UniLink, 4" x 12" Pad** \$106 ea.
 QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
- ☐ **UniLink, 4" x 8" Pad** \$106 ea.
 QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

- ☐ **UniLink, 1.5" Link** \$81 ea.
 QTY Part: ST-UIL001150A000-Q
- ☐ **UniLink, 2" Link** \$81 ea.
 QTY Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

- ☐ **Install on Left**
 Part: PTOINDV3585
- ☐ **Install on Right**
 Part: PTOINDV3586
- ☐ **Do Not Install**
 Part: PTOINDV3719

20.3. Stealth Position Chest Harness

- ☐ **Upper Classic, Medium (IFPS1103)** \$154 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960
- ☐ **Upper Contour, Medium (IFPS2103)** \$154 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960
- ☐ **Upper Classic, Large (IFPS1104)** \$154 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960
- ☐ **Upper Contour, Large (IFPS2104)** \$154 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

20.4. Stealth Position Chest Straps

- ☐ **Chest Strap, Medium** \$97 ea.
 w/ Belt Anchor, (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960
- ☐ **Chest Strap, Large** \$97 ea.
 w/ Belt Anchor, (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

20.5. Stealth Swing-Away Headrest Mounts


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

- ☐ **UniLink Facial Lateral, Swing Away, Right** \$318 ea.
 Part: ST-SAR2A-R-Q. HCPCS: E1028
- ☐ **UniLink Facial Lateral, Swing Away, Left** \$318 ea.
 Part: ST-SAR2A-L-Q. HCPCS: E1028
- ☐ **UniLink Facial Lateral, Dual Swing Away** \$492 ea.
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028
- ☐ **UniLink Facial Lateral, Dual Swing Away** \$492 ea.
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
- ☐ **UniLink Facial Lateral, Dual Swing Away** \$434 ea.
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

2. Facial Lateral Pads and Accessories

- ☐ **Egg Switch Mount (ESM) (A)** \$57 ea.
 QTY Part: INDPART2815
- ☐ **Facial Pad, 4".5x2", Cool Core (900M-Q)** \$61 ea.
 QTY Part: ST-900M-Q. HCPCS: K0108
- ☐ **Facial Pad, 6"x2", Cool Core (920M-Q)** \$61 ea.
 QTY Part: ST-920M-Q. HCPCS: K0108
- ☐ **Swing-Away Joystick Mount, 18" Rod (C)** \$276 ea.
 QTY Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest

**20.6. Stealth Elbow Support****1. Elbow Support Bracket**

- ☐ **Left, Swing-Away, Elbow Support Bracket, Tube** \$199 ea.
 Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028
- ☐ **Right, Swing-Away, Elbow Support Bracket, Tube** \$199 ea.
 Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

- ☐ **3"Wx3"H Elbow Pad (ARMS-ES-3X3)** \$57 ea.
 QTY Part: POS124065. HCPCS: K0108
- ☐ **4"Wx4"H Elbow Pad (ARMS-ES-4X4)** \$57 ea.
 QTY Part: POS124067. HCPCS: K0108

20.7. Stealth Foot Boxes

- ☐ **Stealth Small Foot Boxes, 5"Wx8"Dx7"H** \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954
- ☐ **Stealth Large Foot Boxes, 7"Wx12"Dx11"H** \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954
- ☐ **Stealth Medium Foot Boxes, 6"Wx10"Dx9"H** \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

20.8. Stealth Positioning Shoe Holder

- ☐ **Small Shoe Holder, 3/4" Webbing (IFSH4211)** \$87 ea.
 QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952
- ☐ **Large Shoe Holder, 3/4" Webbing (IFSH4213)** \$87 ea.
 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952
- ☐ **Medium Shoe Holder, 3/4" Webbing (IFSH4212)** \$87 ea.
 QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Positioning components section continued on following page.

21. LAP TRAY

21.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- ☐ Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950
- ☐ Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950
- ☐ Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

- ☐ Swing Away Mounting Hardware..... \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.