



## Q6 Edge® Z Order Form with Contoured Seating

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
 Provider Name: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_  
 PO Number: \_\_\_\_\_  
 Quote Name/Marked for: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Q6 edge<sup>®</sup> Z



### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.


HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

- A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_
- B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_
- C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

### 1. BASE MODEL











**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors, 8 amp off board charger, and switch operated fender lights. Dimensions of base without legrests: 26.25"W x 35.5"L.

Quantum Q6 Edge Z 3S-SS.....\$8,290  
 Part: Q6 EDGE Z 3S-SS. HCPCS: K0848  
 Solid seat pan standard with no power option.

Quantum Q6 Edge Z 3S-C.....\$8,290  
 Part: Q6 EDGE Z 3S-C. HCPCS: K0849  
 Captain seat standard.

**2. BASE COLOR**

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- |   |  |  |  |
|---|--|--|--|
|  <b>Red Corvette</b><br>Part: DGN1610108-PRD1  |  <b>Night Train</b><br>Part: DGN1610108-PBK3    |  <b>Evolution Orange</b><br>Part: DGN1610108-POR7 |  <b>Green Machine</b><br>Part: DGN1610108-PGN11 |
|  <b>Electric Blue</b><br>Part: DGN1610108-PBU9 |  <b>Pretty in Pink</b><br>Part: DGN1610108-PPK3 |  <b>Grey Street</b><br>Part: DGN1610108-PSL3      |  |
|  <b>Back in Black</b><br>Part: DGN1610108-PBK1 |  <b>Deep Purple</b><br>Part: DGN1610108-PPR5    |  <b>Midnight Blue</b><br>Part: DGN1610108-PBU10   |  |

**3. DRIVE MOTORS**

**REQUIRED.**

- 6 mph Drive Motors ..... Standard  
Part: MOT150186

**4. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera ..... \$368  
Part: MEC1910365  
Available for Q-Logic 3 electronics only. NE/NE+ coming soon.



**1. Optional Camera Mounting Bracket**

- Backup Camera Headrest Bracket..... \$21      Backup Camera Track Mount Bracket ..... \$27  
Part: MEC1901231      Part: MEC2104227

**2. Optional Display Mounting Bracket**

- Backup Camera Q-Logic Joystick Mount Bracket ....\$105      Backup Camera NE/NE+ Joystick Mount Bracket ....\$105  
Part: MEC2106123      Part: MEC2106132

**5. DRIVE WHEELS**

**REQUIRED.** Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

**5.1. Drive Wheels**

- |  |  |
|--|--|
| 14" Run Flat Tires ..... Standard<br>Part: PTOINDV3510                           | 4"W Flat Free, Knobby, Extended Traction Tires ..... \$290 ea.<br>Part: WHL121698 x2 |
| 14" Run Flat Tires for Drive Wheel Accents ..... Standard<br>Part: WHL1703421 x2 | 4"W Pneumatic, Knobby, Extended Traction Tires ..... \$290 ea.<br>Part: WHL153737 x2 |

**5.2. Drive Wheel Accent Color**

**Optional.** Only available WHL1703421 above. Sold in Pairs (parts below will be quantity of 2).

- |   |   |  |
|---|---|--|
| Back in Black Accent Color..... \$63<br>Part: XRF000179 | Deep Purple Accent Color..... \$63<br>Part: XRF000184       | Electric Blue Accent Color ..... \$63<br>Part: XRF000188 |
| Red Corvette Accent Color ..... \$63<br>Part: XRF000180 | Pretty in Pink Accent Color ..... \$63<br>Part: XRF000185   | Green Machine Accent Color ..... \$63<br>Part: XRF000190 |
| Night Train Accent Color ..... \$63<br>Part: XRF000181  | Evolution Orange Accent Color ..... \$63<br>Part: XRF000186 |  |
| Grey Street Accent Color..... \$63<br>Part: XRF000182   | Midnight Blue Accent Color..... \$63<br>Part: XRF000187     |  |

**6. BATTERIES**

**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

- NF-22, Gel Cell, Interceptor 2 Battery ..... \$341 ea.      Group 24, 70 Amp Gel Cell, MK Battery ..... \$520 ea.  
Part: BAT2004376. HCPCS: E2361      Part: BATGEL1009. HCPCS: E2363

**7. ELECTRONICS**

**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

**7.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

- |  |  |
|--|--|
| NE+ 6-Key Joystick .....Standard<br>Part: CTLDC1560<br>Control up to 2 actuators through joystick. | NE+ 8-Key Joystick .....\$420<br>Part: CTLDC1561<br>Control up to 2 actuators and Seat Mounted LED Light Package through joystick. |
| No Joystick, Setup for NE+ Joystick ..... No Charge<br>Part: XRF000009                             |  |



Electronics continued on following page.

**7. ELECTRONICS CONTINUED**

**7.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

**1. Q-Logic 3 Required Components**

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

- |  |  |
|--|--|
| Q-Logic 3 EX Controller Module .....\$730<br>Part: XRF000002. HCPCS: E2377 | Q-Logic 3 Expandable Controller Harness.....\$447<br>Part: XRF000264. HCPCS: E2313 |
|--|--|



**2. Q-Logic 3 Joystick**

- |   |   |
|---|---|
| Q-Logic 3 EX Joystick .....\$840<br>Part: CTL167724   | Q-Logic 3 EX Joystick for Full Lights .....\$840<br>Part: CTL167725 |
| Q-Logic 3e Joystick .....\$735<br>Part: CTL1910600<br>Can operate full light package. Not available with alternative or specialty controls. | No Joystick, Setup for Q-Logic<br>Part: XRF000010                   |

**7.2. Q-Logic 3 Joystick Armor**

**REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2.

- |  |   |   |
|--|---|---|
| Black, Soft Touch ..... Standard<br>Part: DGN124307<br>Also matches Night Train shrouds. | Back in Black ..... \$63<br>Part: DGN124336-PBK1    | Grey Street ..... \$63<br>Part: DGN124336-PSL3    |
| Red Corvette ..... \$63<br>Part: DGN124336-PRD1  | Pretty in Pink ..... \$63<br>Part: DGN124336-PPK3   | Midnight Blue ..... \$63<br>Part: DGN124336-PBU10 |
| Electric Blue ..... \$63<br>Part: DGN124336-PBU9   | Evolution Orange ..... \$63<br>Part: DGN124336-POR7 | Green Machine ..... \$63<br>Part: DGN124336-PGN11 |
|  | Deep Purple ..... \$63<br>Part: DGN124336-PPR5      |   |

**8. JOYSTICK HANDLES**

**REQUIRED.** The NE Plus are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

- |  |  |
|--|--|
| Mushroom Handle (B) ..... \$90<br>Part: ACCASMB1860. HCPCS: E2323                    | Extended I-Handle (G) ..... \$90<br>Part: ACCASMB1857. HCPCS: E2323          |
| Chin Cup (C) ..... \$90<br>Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | Bodypoint 3" Goal Post Handle ..... \$104<br>Part: PLSKNCP1057. HCPCS: E2323 |
| T-Bar (D) ..... \$126<br>Part: ACCASMB1858. HCPCS: E2323                             | Bodypoint 4" Goal Post Handle ..... \$104<br>Part: PLSKNCP1058. HCPCS: E2323 |
| Medium Ball Handle (E) ..... \$69<br>Part: ACCASMB1859. HCPCS: E2323                 | Q-Logic 3 Carrot Handle (A) ..... No Charge<br>Part: DGN124308               |
| Soft Ball Handle (F) ..... \$80<br>Part: ACCASMB1856. HCPCS: E2323                   |  |




**9. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**9.1a. Fixed, In-line Joystick Bracket**

- |   |  |
|---|--|
| Fixed, In-line, Bracket for Q-Logic 3 .....Standard<br>Part: FRMASMB10617 | Fixed, In-line, Bracket for NE/NE+ .....Standard<br>Part: FRMASMB12744 |
|---|--|



**9.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- |   |  |
|---|--|
| Swing-away, Q-Logic 3 ..... \$294<br>Part: ACC1908434. HCPCS: E1028 | Swing-away, NE/NE+ ..... \$294<br>Part: ACC1908437. HCPCS: E1028 |
|---|--|

**9.1c. Basic Swing Away Joystick Mount**

- |  |   |
|--|---|
| Basic Swing Away, NE/NE+, Left Side ..... \$245<br>Part: ACC2201049. HCPCS: E1028  | Basic Swing Away, NE/NE+, Right Side ..... \$245<br>Part: ACC2201050. HCPCS: E1028  |
| Basic Swing Away, Q-Logic, Left Side ..... \$245<br>Part: ACC2201043. HCPCS: E1028 | Basic Swing Away, Q-Logic, Right Side ..... \$245<br>Part: ACC2201044. HCPCS: E1028 |

Joystick mounting brackets section continued on following page.

**9. JOYSTICK MOUNTING BRACKETS CONTINUED**

**9.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

**Retract4 Mount**.....\$368  
 Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.


**Multi-Axis Mount** ..... \$126

Part: XRF000004

Allows height adjustment of joystick.

**Stealth Swing-Away/Flip Down Bracket** .....\$250  
 Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.



Retract4



Multi-Axis Joystick Mount

**9.2. Joystick Mount Location**

**Left Side Mount**

Part: PTOINDV1225

**Right Side Mount**

Part: PTOINDV1226

**10. CONTOURED SEATING**

**Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

**10.1. Occupied Transit, WC-19**

**Optional.** Not available with Black Vinyl, High Back Seat

**WC-19 Occupied Transit Option for Contoured Seat** ..... \$368

Part: KIT122609

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

**10.2a. Q Captain Seat**

**16"W x 16"D** ..... No Charge  
 Part: SET2005145

**18"W x 20"D** ..... No Charge  
 Part: SET2005151

**20"W x 22"D** ..... No Charge  
 Part: SET2005909

**16"W x 18"D** ..... No Charge  
 Part: SET2005146

**18"W x 22"D** ..... No Charge  
 Part: SET2005908

**22"W x 18"D** ..... No Charge  
 Part: SET2005155

**16"W x 20"D** ..... No Charge  
 Part: SET2005147

**20"W x 16"D** ..... No Charge  
 Part: SET2005152

**22"W x 20"D** ..... No Charge  
 Part: SET2005156

**18"W x 16"D** ..... No Charge  
 Part: SET2005149

**20"W x 18"D** ..... No Charge  
 Part: SET2005153

**22"W x 22"D** ..... No Charge  
 Part: SET2005157

**18"W x 18"D** ..... No Charge  
 Part: SET2005150

**20"W x 20"D** ..... No Charge  
 Part: SET2005154

**10.2b. Q Captain Seat with Solid Seat Pan**

**16"W x 16"D** ..... No Charge  
 Part: SET2005394

**18"W x 20"D** ..... No Charge  
 Part: SET2005399

**20"W x 22"D** ..... No Charge  
 Part: SET2005911

**16"W x 18"D** ..... No Charge  
 Part: SET2005395

**18"W x 22"D** ..... No Charge  
 Part: SET2005910

**22"W x 18"D** ..... No Charge  
 Part: SET2005403

**16"W x 20"D** ..... No Charge  
 Part: SET2005396

**20"W x 16"D** ..... No Charge  
 Part: SET2005400

**22"W x 20"D** ..... No Charge  
 Part: SET2005404

**18"W x 16"D** ..... No Charge  
 Part: SET2005397

**20"W x 18"D** ..... No Charge  
 Part: SET2005401

**22"W x 22"D** ..... No Charge  
 Part: SET2005405

**18"W x 18"D** ..... No Charge  
 Part: SET2005398

**20"W x 20"D** ..... No Charge  
 Part: SET2005402

**10.1c. Black Vinyl, High Back Seat**

**24"W x 20"D** ..... No Charge  
 Part: SETASMB2387

**24"W x 22"D** ..... No Charge  
 Part: SETASMB2388

**24"W x 24"D** ..... No Charge  
 Part: SETASMB2389

**10.2. Seat to Floor Height**

**REQUIRED.**

**Lowest Seat to Floor**  
 Part: PTOINDV3558

**Middle Seat to Floor**  
 Part: PTOINDV3559

**Highest Seat to Floor**  
 Part: PTOINDV3560

**10.3. Matching Seat Cushion**

**Optional with Q Captain Seat with Solid Seat Pan.** Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
 Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608

TRU-Comfort 2 SPP, Matching Seat Size  
 Reference part: XRF000116. MSRP: \$420/\$525



**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
 Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size  
 Reference part: XRF000117. MSRP: \$483/\$557



Contoured Seating continued on following pages.


**10. CONTOURED SEATING CONTINUED**

**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: *ST-SOLSPPXXYY*. HCPCS: *E2607/E2608*

Solution SPP, Matching Seat Size  
Reference part: *XRF000113*. MSRP: *\$394/\$452*




**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: *ST-SPCFSPXXYY*. HCPCS: *E2607/E2608*

Spectrum Foam SPP, Matching Seat Size  
Reference part: *XRF000112*. MSRP: *\$483/\$557*



**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-PREMPXXYY*. HCPCS: *E2605/E2606*

Premiere P, Matching Seat Size  
Reference part: *XRF000122*. MSRP: *\$341/\$499*



**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-ZENSPXXYY*. HCPCS: *E2603/E2604*

Zen SP, Matching Seat Size  
Reference part: *XRF000115*. MSRP: *\$331/\$446*



**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Part series: *ST-ESSSPXXYY*. HCPCS: *E2607/E2608*

Essence SPP, No Wedge, Matching Seat Size  
Reference part: *XRF000118*. MSRP: *\$420/\$525*

Essence SPP, Small Wedge, Matching Seat Size  
Reference part: *XRF000119*. MSRP: *\$457/\$562*

Essence SPP, Medium Wedge, Matching Seat Size  
Reference part: *XRF000120*. MSRP: *\$488/\$593*

Essence SPP, Large Wedge, Matching Seat Size  
Reference part: *XRF000121*. MSRP: *\$520/\$625*




**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: *ST-GLSPXXYY*. HCPCS: *E2603/E2604*

Glacial SP, Matching Seat Size  
Reference part: *XRF000114*. MSRP: *\$347/\$462*




**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.

Part series: *ST-SIMGXXYY*. HCPCS: *E2601/E2602*

Simplicity G, Matching Seat Size  
Reference part: *XRF000111*. MSRP: *\$116/\$173*



**10.4a. Armrests for Q Captain Seat**

Optional. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Q-Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat armpads.

Full Length, Left Armrest ..... \$90 ea.  
Part: *MEC2201304*. HCPCS: *K0020*

Full Length, Right Armrest..... \$90 ea.  
Part: *MEC2201303*. HCPCS: *K0020*

**1b. Drop-In, Height Adjustable Armrests (T-Arm)**

Includes flat armpads.

Full Length, Drop-In Armrests..... \$79 ea.  
 L  R Left part: *XRF000298*. Right part: *XRF000299*. HCPCS: *E0973*

Desk Length, Drop-In Armrests ..... \$79 ea.  
 L  R Left part: *XRF000300*. Right part: *XRF000301*. HCPCS: *E0973*

**1c. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests ..... \$79 ea.  
 L  R Left part: *ARM2000906*. Right part: *ARM2000907*. HCPCS: *E0973*

Desk Length, Single Post Armrests ..... \$79 ea.  
 L  R Left part: *ARM2000908*. Right part: *ARM2000909*. HCPCS: *E0973*

**1d. Bariatric, Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Bariatric Armrests, 8"-12" H ..... \$90 ea.  
 L  R Left part: *FRMASMB6728*. Right part: *FRMASMB6729*. HCPCS: *E0973*

Full Length, Bariatric Armrests, 10"-14"H ..... \$90 ea.  
 L  R Left part: *FRMASMB6730*. Right part: *FRMASMB6731*. HCPCS: *E0973*

**1e. Heavy Duty, Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 8"-12" H ..... \$90 ea.  
 L  R Left part: *XRF000253*. Right part: *XRF000254*. HCPCS: *E0973*

Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.  
 L  R Left Part: *XRF000257*. Right part: *XRF000258*. HCPCS: *E0973*

**1f. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

Pediatric, Removable Armrests ..... \$79 ea.  
 L  R Left Part: *ARM142876*. Right part: *ARM142877*. HCPCS: *E0973*

**10.4b. Armrests for Black Vinyl, High Back Seat**

Optional. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1. Flip-Up, Height Adjustable Armrests**

Includes full length, flat armpads.

Full Length, Left Armrest ..... \$90  
Part: *FRMASMB8614*

Full Length, Right Armrest..... \$90  
Part: *FRMASMB8615*

Desk Length, Left Armrest..... \$90  
Part: *FRMASMB8612*

Desk Length, Right Armrest ..... \$90  
Part: *FRMASMB8613*

Contoured Seating continued on following pages.

**9. CONTOURED SEATING CONTINUED**

**10.5. Armrest Height**

**REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.**

**2. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

**3. Right Armrest Height**

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

**10.6. Armpads for Q Captain Seat**

**REQUIRED except with Flip-Up and Pediatric armrests.**

**1a. Straight Armpads for Single Post Armrests**

Straight, Full Length Left and Right ..... No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right ..... No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right ..... No Charge Part: ARM134281	Straight, Desk Length Left and Right ..... No Charge Part: ARM134254

**1b. Waterfall Armpads for Single Post Armrests**

Waterfall, Full Length Left and Right ..... No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right ..... No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right ..... No Charge Part: ARM134289	Waterfall, Desk Length Left and Right ..... No Charge Part: ARM134282

**1c. Falcon Gel Armpads for Single Post Armrests**

Falcon Gel, Full Length Left and Right .....\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right .....\$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right ..... \$289.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right ..... \$265.02 Part: ARM136437

**1d. Armpads for Bariatric or Heavy Duty Armrests**

Bariatric Armpads ..... No Charge Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads ..... No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
Straight Armpads ..... No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests	

**10.7. Quad Release Lever**

**Optional. Available with Single Post and Pediatric armrests.**

Quad Release Lever, Pair.....No Charge Part: FRMASMB10517
--

**11. LAP BELTS**

**REQUIRED. Select one lap belt below.**

**11.1a. Standard Lap Belts**

<input type="checkbox"/> 38" Lap Belt Part: ACCBELT1006	<input type="checkbox"/> 50" Lap Belt Part: ACCBELT1000	<input type="checkbox"/> 60" Lap Belt Part: ACCBELT1002	<input type="checkbox"/> 70" Lap Belt Part: ACCBELT1003	<input type="checkbox"/> 80" Lap Belt Part: ACCBELT1004	<input type="checkbox"/> 90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

**11.1b. Optional Lap Belts**

Stealth Positioning, 60" Single Pull, Push Button, Medium.... \$47.38 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning, 60" Single Pull, Push Button, ..... \$47.38 Non-padded 2" Belt (IFPB5310) Part: ST-IFPB5310-Q
Stealth Positioning, 60" Single Pull, Push Button, Large ..... \$47.38 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Part: ACCBELT1016 Required with ISO 7176-19 Compliant Transit Option.
Stealth Positioning, 60" Single Pull, Push Button, ..... \$47.38 Extra Large, Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q	

**12. BASE MOUNTED FOOT PLATFORM**

**REQUIRED. Select a foot platform option. May interfere with seat mounted legrests.**

Base Mounted Foot Platform ..... Standard Part: FRMASMB16077 Not available with center mount foot platforms.	Bell Shaped, Base Mounted Foot Platform ..... \$247 Part: RIGASMB7120041
Omit Base Mounted Foot Platform ..... No Charge Part: PTOINDV3674	

**13. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**13.1a. Basic Center Mount Foot Platform**

**Available with Q-Captain seat only.** Select legrest length in following section.

**1. Basic Center Mount Foot Platform Assembly**

Basic Center Mount.....No Charge  
 Part: XRF000323

**2. Basic Center Mount Extension and Angle Adjustable Footplates**

Short Extension, 5"x9" Footplates ..... \$180 Part: XRF000312. HCPCS: K0040	Long Extension, 5"x9" Footplates ..... \$180 Part: XRF000314. HCPCS: K0040
Medium Extension, 5"x9" Footplates ..... \$180 Part: XRF000313. HCPCS: K0040	

**3. Basic Center Mount Calf Pads**

Small Calf Pads, Basic Center Mount.....\$158  
 Part: RIG2104369

**4. Basic Center Mount Additional Accessory**

Basic Center Mount Footplate Connector ..... \$16  
 Part: RIG2104240

**13.1b. Center Mount Foot Platform**

**Available with Q-Captain seat only.** Select legrest length in following section.

**1. Center Mount Foot Platform Assembly**

Quick Adjustable Center Mount Foot Platform.....\$263  
 Part: XRF000311. HCPCS: K0040

**2. Quick Adjustable Center Mount Lower Extension and Angle Adjustable Footplates**

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x8" Footplates Part: RIG150203	Long Extension, 5"x8" Footplates Part: RIG151506
Pediatric Extension, 5"x8" Footplates Part: RIG152578	Medium Extension, 5"x11" Footplates Part: RIG150212	Long Extension, 5"x11" Footplates Part: RIG151508
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x6" Footplates Part: RIG151505	

**3a. Center Mount Calf Pads**

Small Calf Pads, Manually Adjustable Center Mount Only, ..... \$158 6"Wx6"H (pair) Part: XRF000156	Large Calf Pads, Manually Adjustable Center Mount Only, ..... \$158 8" taper to 6"Wx8"H (pair) Part: XRF000157
--	--

**3b. Center Mount Calf Panel**

Small Calf Panel for Manually Adjustable, 18"Wx8"H ..... \$184 Part: FRMASMB7110012	Large Calf Panel for Manually Adjustable, 24"Wx8"H ..... \$237 Part: RIG123825 Only available on widths of 20" or greater.
--	--

**4. Center Mount Modular Footplate Extension Components**

**Optional.** Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates ..... \$69 Part: RIG152517	4" Outside Edge (pair) ..... \$69 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
Width Extension for 11"D Footplates ..... \$69 Part: RIG152518	6" Outside Edge (pair) ..... \$69 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
2" Bell-Shaped Length Extension ..... \$69 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

**5. Center Mount Additional Accessories**

Small Reversed Dartex Calf Pad Cover (pair) ..... \$158 Part: RIG133717	Large Foam Panel Cover (pair) ..... \$268.77 Part: ST-CSCP0927-Q
Large Reversed Dartex Calf Pad Cover (pair) ..... \$158 Part: RIG133719	Footplate Connector ..... \$16 Part: RIGASMB7120023
Small Reversed Dartex Panel Cover (pair) ..... \$210 Part: RIG134673	Footplate Riser ..... \$75 Part: XRF000348 Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.
Large Reversed Dartex Panel Cover (pair) ..... \$210 Part: RIG134744	Footplate Riser Calf Panel ..... \$50 Part: FRMASMB7110011 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".
Small Foam Calf Pad Cover (pair) ..... \$63.74 Part: ST-CSCS0775-Q	
Large Foam Calf Pad Cover (pair) ..... \$78.72 Part: ST-CSCS0985-Q	
Small Foam Panel Cover (pair) ..... \$167.57 Part: ST-CSCP0920-Q	

Legrests section continued on following page.

**13. LEGRESTS CONTINUED**

**6. Stealth Gel Footplate Covers**

- 5"x8" Tapered Gel Footplate Covers ..... \$211.15  
Part: ST-FPG5008T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Covers ..... \$231.75  
Part: ST-FPG5011T-Q. HCPCS: K0108

- 5"x8" Tapered Gel Footplate Covers for Width Extension ... \$231.75  
Part: ST-FPG6508T-Q. HCPCS: K0108
- 5"x11" Tapered Gel Footplate Cover for Width Extension ... \$257.50  
Part: ST-FPG7011T-Q. HCPCS: K0108

**13.1c. 60°-90° Swing-Away Legrests**

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1a. Swing-Away Legrest Upper**

- P  L  R 60°-90° Swing-Away Legrest ..... No Charge  
Pair part: RIG120062. Left: FRMASMB17558.  
Right: FRMASMB17559
- P  L  R 60°-90° Swing-Away Legrest, Short ..... No Charge  
Pair part: RIG136481. Left: RIG132607.  
Right: RIG132609.

**1b. 70° Swing-Away Legrest Uppers**

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

- 70° Swing-Away Legrest (pair) ..... \$172  
Part: FRMASMB7841  
Not available with 11"-13" extensions.

**2. Lower Extension**

- P  L  R 11"-13" Length  
Pair part: XRF000088.  
Left: XRF000089.  
Right: XRF000090.  
N/A with Multi-axis.
- P  L  R 13"-16.5" Length  
Pair part: XRF000094.  
Left: XRF000095.  
Right: XRF000096.  
N/A with Multi-Axis.
- P  L  R 11"-16.5" Length  
Pair part: XRF000091.  
Left: XRF000092.  
Right: XRF000093.  
Multi-Axis Only.
- P  L  R 15.5"-19.5" Length  
Pair part: XRF000097.  
Left: XRF000098.  
Right: XRF000099.

60°-90° Swing-Away Legrest Compatibility Matrix	Contoured Seat			Synergy Seat		
	Low	Medium	High	Low	Medium	High
Extruded, 13" - 16.5"	✓	✓	✓	!	!	✓
Extruded, 15.5" - 19.5"	!	✓	✓	X	!	!
Angle Adjustable, 13" - 16.5"	✓	✓	✓	!	!	✓
Angle Adjustable, 15.5" - 19.5"	!	✓	✓	X	!	!
Multi-Axis, 11" - 16.5"	✓	✓	✓	!	!	✓
Multi-Axis, 15.5" - 19.5"	!	✓	✓	X	!	!

11"-13" Extensions cause no interference at all seat-to-floor heights.

- This rigging configuration causes no interference issues at this seat-to-floor height.
- This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
- This rigging configuration does not work with this seat-to-floor height.

**3a. Extruded Footplate(s)**

- P  L  R 7.25"Wx6"D Extruded Footplate ..... Standard  
Part pair: FRMASMB7847. Left: FRMASMB7861  
Right: FRMASMB7862

- P  L  R 8.5"Wx6"D Extruded Footplate ..... Standard  
Part pair: FRMASMB7848. Left: FRMASMB7863  
Right: FRMASMB7864  
May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s) with Heel Loop(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- P  L  R 4"Wx6"D ..... \$90 ea.  
Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560  
HCPCS: K0040+E0951
- P  L  R 5"Wx6"D ..... \$90 ea.  
Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561  
HCPCS: K0040+E0951
- P  L  R 5"Wx8"D ..... \$90 ea.  
Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562  
HCPCS: K0040+E0951

- P  L  R 6"Wx8"D<sup>1</sup> ..... \$90 ea.  
Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563  
HCPCS: K0040+E0951
- P  L  R 6"Wx11"D<sup>1</sup> ..... \$90 ea.  
Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564  
HCPCS: K0040+E0951
- P  L  R 8"Wx11"D<sup>2</sup> ..... \$90 ea.  
Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565  
HCPCS: K0040+E0951

**3c. Multi-Axis Footplate(s) with Heel Loops(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

- P  L  R 4"Wx6"D ..... \$260 ea.  
Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552  
HCPCS: K0040+E0951
- P  L  R 5"Wx6"D ..... \$260 ea.  
Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553  
HCPCS: K0040+E0951
- P  L  R 5"Wx8"D ..... \$260 ea.  
Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554  
HCPCS: K0040+E0951

- P  L  R 6"Wx8"D<sup>1</sup> ..... \$260 ea.  
Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555  
HCPCS: K0040+E0951
- P  L  R 6"Wx11"D<sup>1</sup> ..... \$260 ea.  
Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558  
HCPCS: K0040+E0951
- P  L  R 8"Wx11"D<sup>2</sup> ..... \$260 ea.  
Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559  
HCPCS: K0040+E0951

**4. 60°-90° Swing-Away Legrest Accessories**

- Gel Neoprene Footrest Sleeve (pair) ..... \$84  
Part: INDPART2243. HCPCS: K0108

**5. Amputee Supports with Swing-Away Upper**

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- Therafin 31515 Multi-Axis Amputee Support ..... \$230 ea.  
Part: INDPART2822. HCPCS: E1020
- Therafin 32387 Amputee Support Gel Cover ..... \$155 ea.  
Part: INDPART2823. HCPCS: K0108

Legrests section continued on following page.



**13. LEGRESTS CONTINUED**

**13.1d. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

P  L  R Elevating Legrest Assembly ..... \$171 ea.  
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990  
 Select lower extension and footplates

**2. Lower Extension**

P  L  R **11"-14" Length**  
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030  
 For reinforced only

P  L  R **8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**  
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103  
 For angle adjustable only

P  L  R **14"-17" Length**  
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035  
 For reinforced only

P  L  R **13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**  
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106  
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

This rigging configuration causes no interference issues at this seat-to-floor height.

This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

**3a. Reinforced Cast Footplate(s)**

P  L  R **7.25"Wx6"D Reinforced Footplate**..... Standard  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P  L  R **8.5"Wx6"D Reinforced Footplate** ..... Standard  
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

P  L  R **4"Wx6"D** ..... \$90 ea.  
 Pair parts: FRMASMB8812 x2.  
 HCPCS: K0040

P  L  R **6"Wx8"D<sup>1</sup> with Heel Loops**..... \$90 ea.  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795  
 HCPCS: K0040+E0951

P  L  R **5"Wx6"D**..... \$90 ea.  
 Pair parts: FRMASMB8811 x2.  
 HCPCS: K0040

P  L  R **6"Wx11"D<sup>1</sup> with Heel Loops** ..... \$90 ea.  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797  
 HCPCS: K0040+E0951

P  L  R **5"Wx8"D with Heel Loops**..... \$90 ea.  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792  
 HCPCS: K0040+E0951

P  L  R **8"Wx11"D<sup>2</sup> with Heel Loops** ..... \$90 ea.  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799  
 HCPCS: K0040+E0951

**4. Elevating Legrest Accessories**

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) ..... \$131.25  
 Part: INDPART2242. HCPCS: E0953

**13.2. Legrest Accessories**

Heel Loops (pair)..... \$53  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$37  
 Part: XRF000100

**13.3. Set Legrest Extension Length**

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

**1a. Centermount Extension Length**

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

**14. SEAT MOUNTED LED LIGHTS**

**Optional.** Requires Q Captain Seat and Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package .....	\$1,045
<i>Part: XRF000130</i>	
Not available with Base Mounted, LED, Full Light Package.	

**15. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

**15.1. Q Captain Seat or Black Vinyl, High Back Seat Accessories**

Saddle Bag .....	\$40	Grip Drink Holder (20085) .....	\$15.50
<i>Part: DWR1020D010</i>		<i>Part: ACC152993</i>	
Quantum Backpack .....	\$168	Q-Logic EX/NE Programmer .....	\$573
<i>Part: ACC1905853</i>		<i>Part: ELEASMB6976</i>	
Large, Black, Weather Cover .....	\$202	Q-Logic PC Programming Disk with Cable .....	\$445
<i>Part: ACCCOVR1021</i>		<i>Part: ELEASMB5216</i>	
Medium, Black, Weather Cover .....	\$168	Personal Item Pouch .....	\$84
<i>Part: ACCCOVR1022</i>		<i>Part: ACC2101017</i>	
Grip Lap Board (20005).....	\$38.95		
<i>Part: ACC152992</i>			

**15.2. Q Captain Seat Only Accessories**

Cup Holder, C Style .....	\$32	Reflector Kit .....	\$53
<i>Part: ACCASMB2288</i>		<i>Part: ACC124670</i>	
For use with Flip-Up armrests.			
Cup Holder, Bariatric Arms or Gel Armpads .....	\$32	Transfer Handles .....	\$105 ea.
<i>Part: ACCASMB2632</i>		<i>Part: MEC2205680. HCPCS: E0705</i>	
Available with Flip-Up armrests only. Not available with Elevating Legrests.			
Cup Holder, Tube Style Armrests.....	\$32	Clothing Guard (Pair).....	\$237
<i>Part: XRF000221</i>		<i>Part: ACC125007</i>	
For HD Removable, Single Post, or Pediatric Arms.			
Oxygen Holder.....	\$118	Privacy Bag.....	\$42
<i>Part: ACC2006353. HCPCS: E2208</i>		<i>Part: ACC1907570</i>	
Cane/Crutch Holder .....	\$110		
<i>Part: ACC2006352. HCPCS: E2207</i>			
Rear Basket .....	\$118		
<i>Part: ACC2006351</i>			

**15.3. Black Vinyl, High Back Seat Only Accessories**

Contoured Seat Cup Holder.....	\$32	Rear Basket for High Back Seat .....	\$118
<i>Part: ACCASMB2287</i>		<i>Part: FRMASMB1280</i>	
Oxygen Holder for High Back Seat .....	\$118	Flag Holder for High Back Seat.....	\$40
<i>Part: FRMASMB1269. HCPCS: E2208</i>		<i>Part: ACCASMB1907</i>	
Walker Holder for High Back Seat .....	\$103	Stealth Cell Phone/Media Device Holder .....	\$205
<i>Part: FRMASMB1271</i>		<i>Part: ACC7120006</i>	
Cane/Crutch Holder for High Back Seat.....	\$110	Stealth Large Tablet Holder .....	\$400
<i>Part: FRMASMB1270. HCPCS: E2207</i>		<i>Part: ACC7120005</i>	

**16. BRANCH PATCHES**

**Optional.**

Air Force Patch .....	\$21	Coast Guard Patch .....	\$21	Navy Patch .....	\$21
<i>Part: ACC140240</i>		<i>Part: ACC140239</i>		<i>Part: ACC140244</i>	
Army Patch .....	\$21	Marines Patch.....	\$21	Purple Heart Patch .....	\$21
<i>Part: ACC140241</i>		<i>Part: ACC140243</i>		<i>Part: ACC140238</i>	



Military Patches

**17. POSITIONING COMPONENTS**

Optional, Q Captain Seat only. More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.

**17.1. Stealth Unilink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, .....\$252.21 ea.  
Removable Mounting Hardware  
Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable, .....\$262.50 ea.  
Removable Mounting Hardware  
Part: ST-UTTP7-Q. HCPCS: E1028



**2. UniLink Pelvic/Thigh Pads**

UniLink, 4" x 4" Pad .....\$101.25 ea.  
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad .....\$101.25 ea.  
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad.....\$101.25 ea.  
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link .....\$77.20 ea.  
Part: ST-UIL001150A000-Q

UniLink, 2" Link.....\$77.20 ea.  
Part: ST-UIL001200A000-Q

**4. UniLink Pelvic/Thigh Guide Installation**

Install on Left  
Part: PTOINDV3585

Install on Right  
Part: PTOINDV3586

Do Not Install  
Part: PTOINDV3719

**17.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) .....\$267.75 ea.  
Part: INDPART2804. HCPCS: E1028

Flip Down with LEFT Full Surface Contact .....\$267.75 ea.  
Feature (TWBADD-L & LPHW407)  
Part: POS140428. HCPCS: E1028

Removable Bracket (LP430).....\$270.38 ea.  
Part: INDPART2805. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact.....\$267.75 ea.  
Feature (TWBADD-R & LPHW407)  
Part: POS140429. HCPCS: E1028

Removable Bracket with 1" Offset .....\$363.61 ea.  
(LP430 & LHW-121)  
Part: INDPART2854. HCPCS: E1028

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122).....\$27.04  
Part: INDPART2821. HCPCS: K0108  
Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) .....\$54.08  
Part: POS140430. HCPCS: K0108  
Only available with the flip down full surface bracket.

**3. Pelvic/Thigh Guide Pad**

Short Pad, 4"x4" (TWBADD-4X4).....\$77.25 ea.  
Part: INDPART2801. HCPCS: E0953/E0956

Long Pad, 4"x12" (TWBADD-4X12).....\$77.25 ea.  
Part: INDPART2803. HCPCS: E0953/E0956

Medium Pad, 4"x8" (TWBADD-4X8).....\$77.25 ea.  
Part: INDPART2802. HCPCS: E0953/E0956

**4. Pelvic/Thigh Guide Installation**

Install on Left  
Part: PTOINDV3585

Install on Right  
Part: PTOINDV3586

Do Not Install  
Part: PTOINDV3719

**17.3. Stealth Elbow Support**

**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube .....\$189.26 ea.  
Style Arms (ARMS-ES-L)  
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube .....\$189.26 ea.  
Style Arms (ARMS-ES-R)  
Part: POS124066. HCPCS: E1028

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) .....\$54.08 ea.  
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$54.08 ea.  
Part: POS124067. HCPCS: K0108

**17.4. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H.....\$513.71  
(SPLE107) (pair)  
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H.....\$513.71  
(SPLE109) (pair)  
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H .....\$513.71  
(SPLE108) (pair)  
Part: INDPART2787. HCPCS: E0954

Positioning components section continued on following page.

**17. POSITIONING COMPONENTS CONTINUED**

**17.5. Stealth Positioning Shoe Holder**

Small Shoe Holder, 3/4" Webbing (IFSH4211).....\$82.40 ea.  
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952  
QTY \_\_\_\_\_

Medium Shoe Holder, 3/4" Webbing (IFSH4212) .....\$82.40 ea.  
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952  
QTY \_\_\_\_\_

Large Shoe Holder, 3/4" Webbing (IFSH4213) .....\$82.40 ea.  
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952  
QTY \_\_\_\_\_