



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Jazzy® 600 ES Order Form with Synergy® Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



FDA Class II Medical Device*

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 2 power base. 300 lb. weight capacity with Dynamic joystick, 5-amp off-board charger, and 14" flat-free drive wheels. Dimensions of base without legrests: 27.75"Wx35.5"L. *Pride FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

☐ Jazzy 600ES 2S-SS, Cardinal Red..... \$5,359
Part: JAZZY600 ES 2S-SS. HCPCS: K0822

Solid seat pan standard with no power option. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

Account Number: _____ Date: _____ Marked for: _____

2. BASE MOUNTED FOOT PLATFORM

Optional. May interfere with seat mounted legrests.

☐ Base Mounted Foot Platform No Charge
Part: FRMASMB16077

3. UNOCCUPIED SECUREMENT

Optional.

☐ Unoccupied Securement \$105
Part: KITASMB1582

4. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

☐ NF-22, 55 Ah Battery \$347 ea.
Part: BATLIQ1018. HCPCS: E2361

5. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

5.1a. Fixed, In-line Joystick Bracket

☐ Fixed, In-line, Bracket Standard
Part: MEC2002220

1. In-line Joystick Bracket Location

☐ Left Side Mount ☐ Right Side Mount
Part: PTOINDV1225 Part: PTOINDV1226

5.1b. Swing-Away (Retractable) Joystick Bracket

☐ Left, Swing-away, Joystick Bracket \$258 ☐ Right, Swing-away, Joystick Bracket \$258
Part: KIT2002912. HCPCS: E1032 Part: KIT2003310. HCPCS: E1032

6. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

☐ Mushroom Handle \$90 Part: ACCASMB1860. HCPCS: E2323	☐ Soft Ball Handle \$80 Part: ACCASMB1856. HCPCS: E2323
☐ Chin Cup \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	☐ Extended I-Handle \$90 Part: ACCASMB1857. HCPCS: E2323
☐ T-Bar \$126 Part: ACCASMB1858. HCPCS: E2323	☐ Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
☐ Medium Ball Handle \$69 Part: ACCASMB1859. HCPCS: E2323	☐ Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323

7. SYNERGY SEATING SYSTEM

Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 20"W x 20"D.

7.1. Synergy Seating Frame Selection

REQUIRED.

☐ Synergy Seating System \$456
Part: PTOINDV2995

7.2. Seat Width

REQUIRED. Part number is SETWDTH plus number following # below.

☐ 10" W ☐ 11" W ☐ 12" W ☐ 13" W ☐ 14" W ☐ 15" W ☐ 16" W ☐ 17" W ☐ 18" W ☐ 19" W ☐ 20" W
(#1000) (#1001) (#1002) (#1003) (#1004) (#1005) (#1006) (#1007) (#1008) (#1009) (#1010)

7.3. Seat Width Package

REQUIRED.

☐ 10"-14" Width Package ☐ 14"-18" Width Package ☐ 16"-20" Width Package
Part: FRMASMB18011 Part: FRMASMB18012 Part: FRMASMB18013

7.4. Seat Depth

REQUIRED. Part number is SETDPTH plus number following # below.

☐ 10" D ☐ 11" D ☐ 12" D ☐ 13" D ☐ 14" D ☐ 15" D ☐ 16" D ☐ 17" D ☐ 18" D ☐ 19" D ☐ 20" D
(#1000) (#1001) (#1002) (#1003) (#1004) (#1005) (#1006) (#1007) (#1008) (#1009) (#1010)

7.5. Seat Depth Package

REQUIRED.

☐ 10"-14" Depth Package ☐ 14"-18" Depth Package ☐ 15"-20" Depth Package
Part: FRMASMB6981 Part: FRMASMB17823 Part: FRMASMB7732

7.6. Accessory Bar

Optional.

☐ Accessory Bar \$132
Part: XRF000128

Synergy seating systems section continued on following pages.

7. SYNERGY SEATING SYSTEM CONTINUED**7.7. Synergy Seating Back Cane Selection****REQUIRED.**

- | | |
|--|---|
| ☐ Synergy Back Cane No Charge
Part: FRMASMB15845 | ☐ 2" Mesomorph (1" per side) Back Cane \$205
Part: FRMASMB12067 |
| ☐ Flip Forward Back Cane \$184
Part: FRMASMB15846
Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests. | ☐ 2" Endomorph (1" per side) Back Cane \$205
Part: FRMASMB12115 |
| | ☐ 4" Endomorph (2" per side) Back Cane \$205
Part: FRMASMB12139 |

7.8. Synergy Seating Back Cushion**REQUIRED.** Select one type of back. See availability chart for sizes.**1a. Synergy Back Selection**

- | |
|--|
| ☐ Synergy Back \$289
Part: XRF000363 |
|--|

1b. Sport Back Selection

- | | |
|---|---|
| ☐ Sport Back Cushion, 21" Wide or Less \$394
Coming soon. Part: XRF000361 | ☐ Sport Back Cushion, 22" Wide or More \$604
Coming soon. Part: XRF000362 |
|---|---|

1c. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

- | | |
|---|--|
| ☐ TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620 | ☐ TRU-Comfort 4 Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621 |
|---|--|

1d. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL14AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL18AT. HCPCS: E2613 |
| ☐ Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL14AT. HCPCS: E2613 | ☐ Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL18AX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL16AS. HCPCS: E2613 | ☐ Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18AXX. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL16AT. HCPCS: E2613 | ☐ Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL20AS. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL16AX. HCPCS: E2613 | ☐ Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL20AT. HCPCS: E2613 |
| ☐ Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16AXX. HCPCS: E2613 | ☐ Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL20AX. HCPCS: E2613 |
| ☐ Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL18AS. HCPCS: E2613 | ☐ Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20AXX. HCPCS: E2613 |

1e. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| ☐ Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL14DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL18DT. HCPCS: E2620 |
| ☐ Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL14DT. HCPCS: E2620 | ☐ Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL18DX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL16DS. HCPCS: E2620 | ☐ Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18DXX. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL16DT. HCPCS: E2620 | ☐ Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL20DS. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL16DX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL20DT. HCPCS: E2620 |
| ☐ Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16DXX. HCPCS: E2620 | ☐ Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL20DX. HCPCS: E2620 |
| ☐ Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL18DS. HCPCS: E2620 | ☐ Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20DXX. HCPCS: E2620 |

Synergy seating section continued on following pages.

7. SYNERGY SEATING SYSTEM CONTINUED**1f. Q-Back 2.0, AL Series, Extra Deep Contour**

Select a back type and the specific overall back height in subsequent section.

- ☐ Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620
- ☐ Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
Part: QB2AL16ES. HCPCS: E2620
- ☐ Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
Part: QB2AL16ET. HCPCS: E2620
- ☐ Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620
- ☐ Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620
- ☐ Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620
- ☐ Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620
- ☐ Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620
- ☐ Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620
- ☐ Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620
- ☐ Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620

1g. Q-Back 2.0, AFT Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

- ☐ AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613
- ☐ AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613

1h. Q-Back 2.0, AFT Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

- ☐ AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620
- ☐ AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT16DT. HCPCS: E2620
- ☐ AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620
- ☐ AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620
- ☐ AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620

1i. Q-Back 2.0, AFT Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

- ☐ AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AFT16EXX. HCPCS: E2620

1j. Omit Back Selection

- ☐ Omit Back.....No Charge
Part: PTOINDV1287

7.9. Synergy Seating Back Cane Height**REQUIRED except if a TRU-Comfort back is desired.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

- ☐ 14" High Canes
Part: FRMASMB15832
- ☐ 15" High Canes
Part: FRMASMB15833
- ☐ 16" High Canes
Part: FRMASMB15834
- ☐ 17" High Canes
Part: FRMASMB15835
- ☐ 18" High Canes
Part: FRMASMB15836
- ☐ 19" High Canes
Part: FRMASMB15837
- ☐ 20" High Canes
Part: FRMASMB15838
- ☐ 21" High Canes
Part: FRMASMB15839
- ☐ 22" High Canes
Part: FRMASMB15840
- ☐ 23" High Canes
Part: FRMASMB15841
- ☐ 24" High Canes
Part: FRMASMB15842
- ☐ 25" High Canes
Part: FRMASMB15843

7.10. Synergy Seating Overall Back Height**REQUIRED.**

- ☐ 14" Overall Back Height
Part: XRF000136
- ☐ 15" Overall Back Height
Part: XRF000137
- ☐ 16" Overall Back Height
Part: XRF000138
- ☐ 17" Overall Back Height
Part: XRF000139
- ☐ 18" Overall Back Height
Part: XRF000140
- ☐ 19" Overall Back Height
Part: XRF000141
- ☐ 20" Overall Back Height
Part: XRF000142
- ☐ 21" Overall Back Height
Part: XRF000143
- ☐ 22" Overall Back Height
Part: XRF000144
- ☐ 23" Overall Back Height
Part: XRF000145
- ☐ 24" Overall Back Height
Part: XRF000146
- ☐ 25" Overall Back Height
Part: XRF000147
- ☐ 26" Overall Back Height
Part: SETHGHT1039
- ☐ 27" Overall Back Height
Part: SETHGHT1040
- ☐ 28" Overall Back Height
Part: SETHGHT1041
- ☐ 29" Overall Back Height
Part: SETHGHT1042

		Back Width											
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
Overall Height	14"	P	T	P	T	P	T	-	-	-	-	-	-
	15"	P	T	P	T	P	T	-	-	-	-	-	-
	16"	P	T	P	T	P	T	P	T	P	T	P	T
	17"	P	T	P	T	P	T	P	T	P	T	P	T
	18"	P	P	P	T	P	T	P	T	P	T	P	T
	19"	P	P	P	T	P	T	P	T	P	T	P	T
	20"	-	-	P	T	P	T	P	T	P	T	P	T
	21"	-	-	P	T	P	T	P	T	P	T	P	T
	22"	-	-	P	T	P	T	P	T	P	T	P	T
	23"	-	-	P	T	P	T	P	T	P	T	P	T
	24"	-	-	P	T	P	T	P	T	P	T	P	T
	25"	-	-	P	T	P	T	P	T	P	T	P	T
	26"	-	-	P	T	P	T	P	T	P	T	P	T
	27"	-	-	P	T	P	T	P	T	P	T	P	T
	28"	-	-	P	T	P	T	P	T	P	T	P	T
	29"	-	-	P	T	P	T	P	T	P	T	P	T

☐ Sport Back ☐ TRU-Comfort Back

Seat Back Availability Chart

Synergy seating section continued on following pages.

7. SYNERGY SEATING SYSTEM CONTINUED**7.11. Synergy Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

- | | |
|---|---|
| ☐ ☐ Full Length, Flip Back Armrests..... Standard | ☐ ☐ Desk Length, Flip Back Armrests Standard |
| ☐ ☐ Left Part: FRMASMB7038. Right part: FRMASMB7039 | ☐ ☐ Left part: FRMASMB7040. Right part: FRMASMB7041 |

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- | | |
|---|---|
| ☐ ☐ Full Length, Two-Post Armrests.....\$158 ea. | ☐ ☐ Desk Length, Two-Post Armrests\$158 ea. |
| ☐ ☐ Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 | ☐ ☐ Left part: XRF000053. Right part: XRF000054. HCPCS: E0973 |

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

- | | |
|---|---|
| ☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea. | ☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea. |
| ☐ ☐ Left Part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973 | ☐ ☐ Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973 |

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

- | | |
|---|---|
| ☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea. | ☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$90 ea. |
| ☐ ☐ Left Part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973 | ☐ ☐ Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973 |

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

- | |
|---|
| ☐ ☐ Pediatric, Removable Armrests \$79 ea. |
| ☐ ☐ Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973 |

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

- | | |
|---|---|
| ☐ ☐ Full Length, Cantilever Armrests\$103 ea. | ☐ ☐ Desk Length, Cantilever Armrests.....\$103 ea. |
| ☐ ☐ Left Part: FRMASMB16188. Right part: FRMASMB16189 | ☐ ☐ Left part: FRMASMB16190. Right part: FRMASMB16191 |

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- | | |
|---|---|
| ☐ ☐ Full Length, Single Post Armrests \$79 ea. | ☐ ☐ Desk Length, Single Post Armrests \$79 ea. |
| ☐ ☐ Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973 | ☐ ☐ Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973 |

7.12. Synergy Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

- | | |
|--|---|
| ☐ Straight Armpads..... Standard | ☐ Waterfall Armpads..... No Charge |
| ☐ Part: XRF000148 | ☐ Part: XRF000149 |
| ☐ No Armpad Selection No Charge | |
| ☐ Part: PTOINDV1373 | |

7.13. Joystick Bracket Receiver**REQUIRED.**

- | | |
|---|--|
| ☐ Standard Fixed Receiver Standard | ☐ Stealth Height Adjustable Receiver..... \$236 |
| ☐ Part: FRMASMB16062 | ☐ Part: XRF000064. HCPCS: K0108 |

7.14. Lap Belts**REQUIRED.****1. Standard Lap Belts**

- | | | | | | |
|--|--|--|--|--|--|
| ☐ 38" Lap Belt | ☐ 50" Lap Belt | ☐ 60" Lap Belt | ☐ 70" Lap Belt | ☐ 80" Lap Belt | ☐ 90" Lap Belt |
| ☐ Part: ACCBELT1006 | ☐ Part: ACCBELT1000 | ☐ Part: ACCBELT1002 | ☐ Part: ACCBELT1003 | ☐ Part: ACCBELT1004 | ☐ Part: ACCBELT1001 |

2. Optional Lap Belts

- | | |
|---|---|
| ☐ Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q | ☐ Stealth Positioning 60" Single Pull, Push Button, \$50
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q |
| ☐ Stealth Positioning 60" Single Pull, Push Button, Large \$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q | ☐ 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with WC-19 Transit Option.
Part: ACCBELT1016 |

7.15. Seat to Floor Height**REQUIRED.**

- | | |
|--|--|
| ☐ 16.875" Seat to Floor | ☐ 17.875" Seat to Floor |
| ☐ Part: PTOINDV3653 | ☐ Part: PTOINDV3654 |

Synergy seating section continued on following page.

7. SYNERGY SEATING SYSTEM CONTINUED**7.17. Headrests, Comfort Plus with Cool Core****Optional.****1. Headrest Pads, Single Point Adjustable**

- ☐ Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)\$179
Part: CP08CCSPM. HCPCS: E0955
- ☐ Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)\$179
Part: CP14CCSPM. HCPCS: E0955
- ☐ Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)\$179
Part: CP10CCSPM. HCPCS: E0955
- ☐ Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)\$179
Part: CP06CCSPM. HCPCS: E0955

2. Headrest Brackets

- ☐ Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033
- ☐ Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033

7.18. Flip-down Headrest Bracket**Optional.**

- ☐ Stealth Flip-down Headrest Bracket (FDM350)..... \$170
Part: INDPART2219
The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

8. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**8.1a. Basic Center Mount Foot Platform****1. Basic Center Mount Foot Platform Assembly**

- ☐ Basic Center Mount.....No Charge
Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

- ☐ Short Extension, 5"x9" Footplates \$180
Part: XRF000332
- ☐ Medium Extension, 5"x9" Footplates \$180
Part: XRF000333
- ☐ Long Extension, 5"x9" Footplates \$180
Part: XRF000334

Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High
Basic			
Short Extension	8"-10"	8"-10"	8"-10"
Medium Extension	10"-13"	10"-14"	10"-14"
Long Extension	N/A	14"	14"-15"

3. Basic Center Mount Calf Pads

- ☐ Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369
- ☐ Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Extension.

4. Basic Center Mount Additional Accessory

- ☐ Basic Center Mount Footplate Connector \$16
Part: RIG2104240

8.1b. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat\$86 ea.
P L R Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.
- ☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat\$86 ea.
P L R Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.

2. Lower Extension

- ☐ ☐ ☐ 11"-13" Length
P L R Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.
- ☐ ☐ ☐ 13"-16.5" Length
P L R Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.
- ☐ ☐ ☐ 15.5"-19.5" Length
P L R Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

Legrests section continued on following page.

8. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | | | |
|--|---|--|--|
| ☐ ☐ ☐
P L R | 5"Wx8"D Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 | ☐ ☐ ☐
P L R | 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 |
| ☐ ☐ ☐
P L R | 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 | ☐ ☐ ☐
P L R | 8"Wx11"D² Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 |

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- | | |
|--|---|
| ☐ Therafin 31515 Multi-Axis Amputee Support \$230 ea.
Part: INDPART2822. HCPCS: E1020 | ☐ Therafin 32387 Amputee Support Gel Cover \$155 ea.
Part: INDPART2823. HCPCS: K0108 |
|--|---|

8.1c. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

- | | |
|--|--|
| ☐ ☐ ☐
P L R | Elevating Legrest Assembly \$171 ea.
Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: K0195
<i>Select lower extension and footplates</i> |
|--|--|

2. Lower Extension

- | | | | |
|--|--|--|---|
| ☐ ☐ ☐
P L R | 11"-14" Length
Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
<i>For reinforced only</i> | ☐ ☐ ☐
P L R | 14"-17" Length
Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
<i>For reinforced only</i> |
| ☐ ☐ ☐
P L R | 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
Pair part: XRF000101. Left: XRF000102. Right: XRF000103
<i>For angle adjustable only</i> | ☐ ☐ ☐
P L R | 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
Pair part: XRF000104. Left: XRF000105. Right: XRF000106
<i>For angle adjustable only</i> |

3a. Reinforced Cast Footplate(s)

- | | | | |
|--|--|--|--|
| ☐ ☐ ☐
P L R | 7.25"Wx6"D Reinforced Footplate Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056 | ☐ ☐ ☐
P L R | 8.5"Wx6"D Reinforced Footplate Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
<i>May cause interference on seat widths less than 18."</i> |
|--|--|--|--|

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

- | | | | |
|--|---|--|--|
| ☐ ☐ ☐
P L R | 5"Wx8"D Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 | ☐ ☐ ☐
P L R | 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 |
| ☐ ☐ ☐
P L R | 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 | ☐ ☐ ☐
P L R | 8"Wx11"D² Footplates with Heel Loops \$90 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 |

4. Elevating Legrest Accessories

- | |
|--|
| ☐ Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
Part: INDPART2242. HCPCS: E0953 |
|--|

8.2. Legrest Accessories

- | | |
|--|--|
| ☐ Heel Loops (pair) \$53
Part: FRMASMB7873. HCPCS: E0951
<i>Standard on angle adjustable and multi axis footplates when noted.</i> | ☐ Wheel Bumper \$37
Part: XRF000100 |
|--|--|

8.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

- | | | | | | | | | | | |
|---|---|--|--|--|--|--|--|--|--|--|
| ☐ 7" L
(#3470) | ☐ 8.5" L
(#3473) | ☐ 10" L
(#3476) | ☐ 11.5" L
(#3479) | ☐ 13" L
(#3482) | ☐ 14.5" L
(#3485) | ☐ 16" L
(#3488) | ☐ 17.5" L
(#3491) | ☐ 19" L
(#3494) | ☐ 20.5" L
(#3497) | ☐ 22" L
(#3669) |
| ☐ 7.5" L
(#3471) | ☐ 9" L
(#3474) | ☐ 10.5" L
(#3477) | ☐ 12" L
(#3480) | ☐ 13.5" L
(#3483) | ☐ 15" L
(#3486) | ☐ 16.5" L
(#3489) | ☐ 18" L
(#3492) | ☐ 19.5" L
(#3495) | ☐ 21" L
(#3498) | ☐ 22.5" L
(#3670) |
| ☐ 8" L
(#3472) | ☐ 9.5" L
(#3475) | ☐ 11" L
(#3478) | ☐ 12.5" L
(#3481) | ☐ 14" L
(#3484) | ☐ 15.5" L
(#3487) | ☐ 17" L
(#3490) | ☐ 18.5" L
(#3493) | ☐ 20" L
(#3496) | ☐ 21.5" L
(#3600) | |

9. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

- | | |
|--|---|
| ☐ Synergy Seat Cup Holder No Charge
Part: ACCASMB2634 | ☐ Stealth Cell Phone/Media Device Holder \$215
Part: ACC7120006 |
| ☐ Oxygen Holder \$158
Part: ACC125006. HCPCS: E2208 | ☐ Stealth Large Tablet Holder \$625.28
Part: ST-URPKTAB0000-Q |
| ☐ Walker Holder \$109
Part: FRMASMB7031 | ☐ Grip Drink Holder (20085) \$15.50
Part: ACC152993 |
| ☐ Cane/Crutch Holder \$116
Part: FRMASMB4822. HCPCS: E2207 | ☐ Large, Black, Weather Cover \$202
Part: ACCCOVR1021 |
| ☐ Rear Basket \$124
Part: FRMASMB4821 | ☐ Transfer Handles (pair) \$210
Part: XRF000396. HCPCS: E0705 |
| ☐ Flag Holder \$55
Part: ACCASMB2415 | |
| ☐ Saddle Bag \$40
Part: DWR1020D010 | |

10. BRANCH PATCHES**Optional.**

- | | | |
|--|--|---|
| ☐ Air Force Patch \$21
Part: ACC140240 | ☐ Coast Guard Patch \$21
Part: ACC140239 | ☐ Navy Patch \$21
Part: ACC140244 |
| ☐ Army Patch \$21
Part: ACC140241 | ☐ Marines Patch \$21
Part: ACC140243 | ☐ Purple Heart Patch \$21
Part: ACC140238 |

11. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**11.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

- | |
|---|
| ☐ Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034 |
|---|

2. Pads

- | | |
|---|---|
| ☐ Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956 | ☐ Flat Pads, 3"x5", pair (ULP846050320) \$250
Part: ST-ULP846050320-Q. HCPCS: E0956 |
| ☐ Curved Pads, 4"x6", pair (TWBL-4X6) \$250
Part: INDPART2664. HCPCS: E0956 | ☐ Flat Pads, 4"x6", pair (ULP835050320) \$250
Part: ST-ULP835050320-Q. HCPCS: E0956 |

3. Installation

- | | |
|--|--|
| ☐ Install Laterals
Part: PTOINDV3499 | ☐ Do Not Install
Part: PTOINDV3500 |
|--|--|

11.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

- | | |
|---|--|
| ☐ UniLink Flip Down Hip/Thigh Support, Left \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^ | ☐ UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^ |
| ☐ Flip Down with LEFT Full Surface Contact \$281
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^ | ☐ Flip Down with RIGHT Full Surface Contact \$281
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^ |
| ☐ UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^ | ☐ UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^ |

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

- | | |
|---|--|
| ☐ UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^ | ☐ UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^ |
| ☐ UniLink, 4" x 8" Pad \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^ | |

3. UniLink Pelvic/Thigh Additional Links

- | | |
|--|--|
| ☐ UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q | ☐ UniLink, 2" Link \$81 ea.
Part: ST-UIL001200A000-Q |
|--|--|

4. UniLink Pelvic/Thigh Guide Installation

- | | | |
|---|--|--|
| ☐ Install on Left
Part: PTOINDV3585 | ☐ Install on Right
Part: PTOINDV3586 | ☐ Do Not Install
Part: PTOINDV3719 |
|---|--|--|

Positioning components section continued on following page.

11. POSITIONING COMPONENTS CONTINUED**11.3. Stealth Positioning Chest Harness**

- | | |
|--|---|
| ☐ Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960 | ☐ Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960 |
| ☐ Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960 | ☐ Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960 |

11.4. Stealth Positioning Chest Straps

- | | |
|--|---|
| ☐ Chest Strap, Medium w/ Belt Anchor (IFCS4113) \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960 | ☐ Chest Strap, Large w/ Belt Anchor (IFCS4114) \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960 |
|--|---|

11.5. Stealth Swing-Away Headrest Mounts


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

- | | |
|--|---|
| ☐ UniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028 | ☐ UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028 |
| ☐ UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028 | ☐ UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028 |
| ☐ UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028 | |

2. Facial Lateral Pads and Accessories

- | | |
|--|--|
| ☐ Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815 | ☐ Swing-Away Joystick Mount,..... \$276 ea.
18" Rod (C)
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest |
| ☐ Facial Pad, 4.5"x2" (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108 | |
| ☐ Facial Pad, 6"x2" (920M-Q).....\$61 ea.
Part: ST-920M-Q. HCPCS: K0108 | |

**11.6. Stealth Elbow Support****1. Elbow Support Bracket**

- | | |
|---|---|
| ☐ Left, Swing-Away, Elbow Support Bracket\$199 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028 | ☐ Right, Swing-Away, Elbow Support Bracket.....\$199 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028 |
|---|---|

2. Elbow Support Pad

- | | |
|--|---|
| ☐ 3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108 | ☐ 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108 |
|--|---|

11.7. Stealth Foot Boxes

- | | |
|---|--|
| ☐ Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954 | ☐ Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954 |
| ☐ Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954 | |

11.8. Stealth Positioning Shoe Holder

- | | |
|--|---|
| ☐ Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952 | ☐ Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952 |
| ☐ Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952 | |

12. STEALTH SEAT CUSHION

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

Cushion Model ####-###XXYY	Width	Depth
Cushion Top Foam V=Visco Top Foam: 1" thick, XXSoft viscoelastic foam topper. F=FRG Top Foam: Under 22"W, blue XXSoft FRG foam topper. 22"W and over, green MedSoft FRG foam topper. N=No Top Foam	Cushion Liner U=Urethane Liner: 2Way Stretch. P=Polyurethane Liner: 4Way Stretch. N=No Liner	Cushion Cover C=CoolCore: Chemical free, 4Way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. M=Mesh: Superior air circulation, non-slip. T=Stealth Tek: Waterproof with stretch capability.

1. Cushion Type

- | | |
|---|---|
| ☐ TRU-Comfort 2 SPP
Reference part: XRF000116. See cushion order form.
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608 | ☐ Premiere P
Reference part: XRF000122. See cushion order form.
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606 |
| ☐ Essence SPP
Reference part: XRF000118. See cushion order form.
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608 | ☐ Glacial SP
Reference part: XRF000114. See cushion order form.
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604 |
| ☐ Solution SPP
Reference part: XRF000113. See cushion order form.
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608 | ☐ Zen SP
Reference part: XRF000115. See cushion order form.
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604 |
| ☐ Spectrum Foam SPP
Reference part: XRF000112. See cushion order form.
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608 | ☐ Simplicity G
Reference part: XRF000111. See cushion order form.
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602 |
| ☐ Spectrum Gel SPP
Reference part: XRF000117. See cushion order form.
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608 | |

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

- | | | |
|---|--|---|
| ☐ Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller | ☐ Medium Wedge Positioning Kit\$72
Part: POP1909006-M
14"W and smaller | ☐ Large Wedge Positioning Kit \$105
Part: POP1909006-L
14"W and smaller |
| ☐ Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger | ☐ Medium Wedge Positioning Kit\$72
Part: POP1712340-M
15"W and larger | ☐ Large Wedge Positioning Kit \$105
Part: POP1712340-L
15"W and larger |

3. Cushion Top Foam

- | | |
|---|--|
| ☐ Visco Top Foam
Reference part: XRF000380 | ☐ No Top Foam
Reference part: XRF000382 |
| ☐ FRG Top Foam
Reference part: XRF000381 | |

4. Cushion Liner

- | | |
|---|---|
| ☐ Stealth Urethane Liner
Reference part: XRF000383 | ☐ No Liner
Reference part: XRF000385 |
| ☐ Stealth Polyurethane Liner
Reference part: XRF000384 | |

5. Cushion Cover

- | | |
|---|--|
| ☐ CoolCore Cover
Reference part: XRF000386 | ☐ Stealth Tek Cover
Reference part: XRF000388 |
| ☐ Mesh Cover
Reference part: XRF000387 | |

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

- ☐
- ABS Flat Rigidizer**
-
- Reference part: XRF000389
-
- MSRP: \$48-\$69

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

Stealth seat cushion section continued on following page.

Account Number: _____ Date: _____ Marked for: _____

12. STEALTH SEAT CUSHION CONTINUED

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

☐ _____ Additional Stealth Urethane Liner
QTY Reference Part: XRF000390

☐ _____ Additional Stealth Polyurethane Liner
QTY Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

☐ _____ Additional CoolCore Cover
QTY Reference Part: XRF000392

☐ _____ Additional Stealth Tek Cover
QTY Reference Part: XRF000394

☐ _____ Additional Mesh Cover
QTY Reference Part: XRF000393

©2026 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

