

Jazzy® 1122

Order Form / Price List
Revised: December 2005



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

BASE ONLY

- Jazzy 1122 K0011 \$5530
*Programmable speed up to 6 MPH,
300 lb. weight capacity.*
- High-torque Motors NC
Top speed of 4.5 MPH

COLOR

- Champagne Std.
- Candy Apple Red \$52
- Viper Blue \$52
- Onyx Black \$52

CONTROLLER

- PG Electronics**
- VSI 70 amp Std.
 - VSI 6-Key 70 amp¹¹ E2399 \$200
 - 3-Key Remote Plus 70 amp^{11,19} \$430
(Select Box for Profile) E2320
 - 1 Profile
 - Multiple Profile
 - 7-Key Accessory Remote Plus 70 amp^{11,19,21} E2320 \$713
 - Remote Plus 70 amp with Speed Dial^{11,15} E2320 \$627
 - Remote Plus 70 amp with Toggles^{11,15} E2320 \$660
 - Remote Plus 70 amp with Switch Jacks^{11,15} E2320 \$615
 - Microdrive, 70 amp^{9,15} E2320 \$995
- Dynamic Electronics^{11,15}**
- Europa G805 70 amp E2320 \$430
 - Europa G80T 70 amp with Toggles E2320 \$765
 - Europa G80A 70 amp with Switch Jacks E2320 \$765

SPECIALTY CONTROLS

- See Specialty Controls Order Form⁴

JOYSTICK HANDLE OPTIONS

- PG Controls**
- PG Handle Std.
 - Mushroom Handle E2323 \$71
 - T-Bar E2323 \$107
 - Chin Cup E2323 \$71
 - Extended I-Handle E2323 \$71
 - Soft Ball Handle E2323 \$50
 - Medium Ball Handle E2323 \$50

JOYSTICK HANDLES CONT.

- Dynamic Controls**
- Ball Knob Std.
 - Sponge Ball E2323 \$35
 - Chin Cup E2323 \$30
 - DL Knob E2323 \$20
 - Carrot Knob E2323 \$20
 - Straight Bar (Chrome) E2323 \$45
 - Finger Tube (Chrome) E2323 \$45
 - Round T-Bar (Chrome) E2323 \$45
 - Flat T-Bar (Chrome) E2323 \$45
 - Wrist Plate E2323 \$45

PROGRAMMER

- Remote Plus/Pilot Programmer \$423
- Dynamic Programmer \$423

CONTROLLER MOUNT

- Fixed, In-line, Joystick Mount Std.
 - Left
 - Right
- Quick Disconnect, Joystick Mount²⁰ NC
 - Left
 - Right
- Swing-away Joystick Mount E1028 \$196
 - Left
 - Right
- Multi-axis Joystick Mount K0108 \$100
 - Left
 - Right*(Fixed, does not swing-away)*

DRIVE WHEELS

- 14" Pneumatic (pair) Std.
- 14" Flat-Free (pair) K0093 \$213

BATTERIES (2 required)

- Group 24, 70 amp, Gel E2363 \$315
 - Group 34, 60 amp, Gel K0108 \$300
- Required w/ Super Low Tilt*
- NF-22, 50-amp Batteries E2361
 - Maintenance-free AGM \$225
 - Maintenance-free Gel Cell \$265
- Required with power tilt / recline elevating system*
- Ship unit with batteries installed \$200
 - Ship unit with batteries, but do not install NC

BATTERY CHARGER

- 8 amp Off-board Charger Std.
 - 5 amp On & 8 amp Off-board Charger \$220
- Off-board charger required with super low tilt*

FRONT RIGGINGS

- Foot Platform, 9.5" deep NC
- Foot Platform, 11.5" deep Std.
- No Front Rigging NC

FRONT RIGGINGS CONT.

- Rigging Type (pair)**
- Swing-away footrests \$162
 - Elevating Leg Rests E0990 \$325
- Select footplate size for above:
- 7 1/4 x 6
 - 8 1/2 x 6
- Adjustable Calf Pads w/ ELRS (1 hole std.)
- 1 hole
 - 3 holes*
 - 6 holes*
- *Potential for interference w/ shrouds*
- Heavy-Duty Drop-Ins, 70^{07,8,22} K0108 \$420
 - Heavy-Duty Drop-Ins, 80^{07,8,22} K0108 \$420

Footplates (pair)

- Reinforced Footplate Std.
 - 7 1/4 x 6
 - 8 1/2 x 6
- 5"x8" Angle-Adj. Footplate¹⁷ K0040 \$150
- 6"x8" Angle-Adj. Footplate³⁶ K0040 \$150
- 6"x11" Angle-Adj. Footplate³⁶ K0040 \$150
- 8"x11" Angle-Adj. Footplate²⁸ K0040 \$150
- Wheel Bumpers (Ang-adj. only) \$25

Extensions¹²

- 11"-14" Ext. NC
- 14"-17" Ext. Std.
- 17"-20" Ext. NC

Rigging Options

- Heel Loops (pair)¹⁴ E0951 \$20
- Gel Neoprene Footrest Sleeve (pair)¹⁶ K0108 \$65
- ELR/ALR Adductor Buttons (pr) K0108 \$100
- POSAlinc Knee Adductor (pair) for S/A legrests K0108 \$464
- Residual limb support left E1020 \$300
- Residual limb support right E1020 \$300

CUSHIONS⁵

- Simplicity Cushion
 - 16 x 16 - 20 x 18 E2601 \$98
- Solution Cushion
 - 10 x 10 - 20 x 18 E2607 \$380
 - 20 x 20 - 20 x 22 E2607 \$430
 - 22 x 16 - 24 x 24 E2608 \$430
- Spectrum Cushion
 - Gel
 - 10 x 10 - 20 x 18 E2605 \$480
 - 20 x 20 - 20 x 24 E2605 \$530
 - 22 x 16 - 24 x 24 E2606 \$530
 - Visco
 - 10 x 10 - 20 x 18 E2607 \$480
 - 20 x 20 - 20 x 24 E2607 \$530
 - 22 x 16 - 24 x 24 E2608 \$530
- Structure Cushion
 - Gel
 - Visco
 - 10 x 10 - 20 x 18 E2607 \$560
 - 20 x 20 - 20 x 22 E2607 \$610

Jazzy® 1122

Order Form / Price List
Revised: December 2005

POWER ELEVATING SEAT^{13,23,24,26}

- Power Elevating Seat Through Toggle \$975
- Power Elevating Seat Through Remote Controller²⁷ \$1375
- Lap Belt Std.

ACCESSORIES

- Rear Basket^{1,2} \$103
- Cane/Crutch Holder^{1,2} \$98
- Oxygen Holder^{1,2} K0104 \$103
- Walker Holder^{1,2,3} \$92
- Flag Holder \$17
- Dust Cover \$192
 - Gray
 - Burgundy
- Battery Spacers²⁹ \$75
- Saddle bag (contoured seating only) \$28
- Manual Hand Brake NC
- Cup holder \$17
 - Left
 - Right
- Lap Belt E0978 \$38
 - Pediatric 38" ²⁵
 - 50"
 - 60"
 - 70"
 - 80"
 - 90"

*SEAT DIMENSION CODING FOR CONTOURED SEATING 20" AND LARGER IS K0108

STANDARD SEATING

- High-Back with Headrest Std.
- High-Back with Headrest^{34,35} w/ Solid Seat Pan Std.
- Black Vinyl
- Gray Vinyl
- Beige Vinyl
- Gray Fabric
- Brown Fabric

High-Back Seating Dimensions (Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.
<input type="checkbox"/> 20"*	\$325	<input type="checkbox"/> 20"*	\$454
<input type="checkbox"/> 22"*	\$325	<input type="checkbox"/> 22"*	\$454

22" W AND ABOVE WILL NOT BE LTD RECLINE.

SEATING

- Medium-Back NC
- Medium-Back w/ Solid Seat Pan^{34,35} NC
- Brown Vinyl
- Beige Vinyl
- Gray Vinyl
- Black Vinyl
- Gray Fabric
- Brown Fabric

Medium-Back Seating Dimensions (Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.

ARMRESTS

For Medium and High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable removable³⁰ E0973 \$170

Armpads

- Desk Length, 10" NC
 - Left
 - Right
- Full Length, 14" NC
 - Left
 - Right
- Black Vinyl Armpads Std.

MANUAL RECLINER^{26,32,33,}

- Manual Recline High Back E1226 \$380
 - Manual Recline High Back³⁴ w/ Solid Seat Pan E1226 \$380
(180 degrees of Recline)
(Dimensions 18"W x 18"D)
 - Gray Vinyl
 - Beige Vinyl
- | Width: | | Depth: | |
|-------------------------------|-------|-------------------------------|-------|
| <input type="checkbox"/> 16" | NC | <input type="checkbox"/> 16" | NC |
| <input type="checkbox"/> 18" | Std. | <input type="checkbox"/> 18" | Std. |
| <input type="checkbox"/> 20"* | \$325 | <input type="checkbox"/> 20"* | \$454 |

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length.

REHAB SEATING

- Omit seat
- See Specialty Seat Order Form
- See Synergy Seating Order Form
- See TRU-Balance Power Tilt Order Form
- See TB Power Recline Order Form
- See TB Power Tilt & Recline Order Form
- See Motion Concepts Order Form

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility. Contact your Quantum Representative for details.

Footnotes:

1. Rear mounted accessories are N/A w/ manual recline high-back seating.
2. Includes mounting bracket and hardware.
3. Walker holder is available on contoured, semi-recline, and Synergy seat.
4. May require special order, call for lead-time, availability, and pricing.
5. See Synergy cushion order form to select sizes.
6. N/A w/ manual recline high-back.
7. Standard w/ heel loops and angle-adjustable footplates, check size below.
8. Available on Synergy Seating systems and Motion Concepts TRx only.
9. Requires multi-axis joystick mount; No charge for mount.
10. Group 34 low profile batteries are required for use with the TRU-Balance Super Low Power Tilt.
11. Needed for controlling power positioning seating systems.
12. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
13. High-torque motors have a half-speed inhibitor. High-speed motors have a quarter-speed inhibitor.
14. Std. when angle-adjustable footplates are selected.
15. Available w/ five individually programmable drive profiles.
16. Available w/ swing-away footrests, 70°-80° heavy-duty drop-ins.
17. Only avail. on 14" W or wider seats.
18. NF-22 batteries are necessary w/ power elevating seat in combination with TB Power Tilt/Recline
19. Std. w/ one programmable drive profile.
20. N/A w/ waterfall armpads; only available on Quantum Seating Systems.
21. Order only when controlling lights through controller.
22. May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Customer Service.
23. Will raise min. seat-to-floor height of unit 1.5"
24. N/A w/ any power tilt or tilt/ recline units.
25. N/A w/ limited recline seating.
26. Power elevating seat is not available with manual recliner.
27. Includes appropriate actuator control module.
28. Only avail. on 20" W or wider seats.
29. Available when NF-22 batteries are selected w/o power-elevating seat.
30. Not available on Manual Recliner.
31. Arms are upholstered to match the color and / or fabric of the seat and back.
32. Rear mount accessories are not N/A.
33. Special measurements are N/A.
34. Includes black vinyl seat pan, back available in all listed colors.
35. 18" x 18" std. Special order available. Call Customer Service.
36. Only avail. on 16" W or wider seats.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.