

Jazzy 600

Group 3 Single Power
Order Form
250lb Weight Capacity



Pride Mobility Products Corporation
182 Susquehanna Ave., Exeter, PA 18643
Phone: (866) 800-2002 • Fax: (866) 707-3422




Discount: _____
Terms: _____
Freight: _____
Taken By: _____

Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____
Fax: _____
Date: _____
PO#: _____
Marked For: _____
Ship To Address: _____
City: _____ State: _____ Zip: _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Items Designated with a Circular Checkbox (⊙) are not available through Quantum West

BASE MODEL			EXPANDABLE CONTROLS		
Description	Price	HCPCS Code	Description	Price	HCPCS Code
<input type="checkbox"/> Jazzy 600 3SP-SS	\$5,363	K0856	<i>All Joysticks or Omni Plus must be used in conjunction with PG Controller and Harness for Expandable Controller</i>		
COLOR			<input type="checkbox"/> PG Controller	\$700	E2377
<input type="checkbox"/> Candy Apple Red			<input type="checkbox"/> Harness for Expandable Controller	\$350	E2313
<input type="checkbox"/> Viper Blue			<input type="checkbox"/> Remote Plus 3-Key	\$700	
<input type="checkbox"/> Onyx Black			<input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile		
<input type="checkbox"/> Silver			<input type="checkbox"/> Remote Plus 3-Key ¹	\$700	
<input type="checkbox"/> Emerald Green			With Switch Jacks		
<input type="checkbox"/> Sapphire Blue			<input type="checkbox"/> Remote Plus 3-Key ¹	\$700	
BATTERY CHARGER			With Speed Dial		
Description	Price		<input type="checkbox"/> Remote Plus 3-Key ¹	\$700	
<input type="checkbox"/> 5 Amp Off-board	Standard		With Toggles		
DRIVE WHEELS			<input type="checkbox"/> Remote Plus 3-Key,	\$700	
Description	Price		With Environmental Controls		
<input type="checkbox"/> 14" Pneumatic (pair)	No Charge		<input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile		
<input type="checkbox"/> 14" Flat-Free (pair)	\$213		<input type="checkbox"/> Remote Plus 3-Key, With Toggles	\$700	
BASE MOUNTED FOOT PLATFORM			<input type="checkbox"/> Remote Plus 3-Key, With Toggles	\$700	
Description	Price		And Environmental Controls ¹		
<input type="checkbox"/> Base Mounted Foot Platform	Standard		<input type="checkbox"/> Remote Plus 3-Key, With Switch	\$700	
<input type="checkbox"/> Omit Base Mounted Foot Platform	No Charge		Jacks And Environmental Controls		
LED LIGHT PACKAGE			<input type="checkbox"/> Microdrive Joystick ¹	\$700	
Description	Price				
<input type="checkbox"/> Through Remote Plus or NE+ 8-Key w/ Lights	\$995				
NON-EXPANDABLE CONTROLS					
Description	Price				
<input type="checkbox"/> 75 Amp NE	Standard				
<input type="checkbox"/> 75 Amp NE+ 6-Key For Up To 2 Actuators Thru Joystick	\$200				
<input type="checkbox"/> 75 Amp NE+ 8-Key For Lights	\$400				
<input type="checkbox"/> P&G 6-Key VR2, 60 Amp	\$200				

1. Standard with 5 individual drive profiles.

PROGRAMMING		
Description	Price	HCPCS Code
<input type="checkbox"/> P&G Programmer 	\$423	

JOYSTICK MOUNT		
Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline ¹	Standard	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-away	\$210	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
Style: <input type="checkbox"/> Inline <input type="checkbox"/> Offset		
Length: <input type="checkbox"/> 4" Ext. <input type="checkbox"/> 7" Ext. ²		
<input type="checkbox"/> Multi-axis Fixed	\$100	K0108
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Not available with Q-logic EX Stand Alone Joystick.
2. In iQuantum choose in addition to Swing-Away 4" mount.

JOYSTICK HANDLES		
Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup	\$71	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$69	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

TRU-BALANCE® POWER POSITIONING		
Description	Price	HCPCS Code
<input type="checkbox"/> TRU-Balance® Power Tilt (55°) <i>55 Degree Center Of Gravity (CG) Power Tilt Seating System. 250lbs Weight Capacity</i>	\$4,995	E1002

SEATING SPECIFICATIONS		
Seat Width:		
<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"
<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"		
Seat Width Packages:		
<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 15"W-20"W	
Seat Depth:		
<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"
<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"		
Seat Depth Package:		
<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 15"D-20"D	
Back Cane Height:		
<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"
<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"
<input type="checkbox"/> 23"	<input type="checkbox"/> 24"	<input type="checkbox"/> 25"

BACK OPTIONS		
Description	Price	HCPCS Code
Tilt Back:		
<input type="checkbox"/> Synergy® Shape Back <i>Positioning Back</i>	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
<input type="checkbox"/> Omit Back	-\$100	
<input type="checkbox"/> Synergy® Back <i>with 1.5" Thick Back Cushion</i>	Standard	
Overall Back: <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
<i>Choose Lateral Extensions¹:</i>		
<input type="checkbox"/> 6" for 13"-16" H Back	\$25	
<input type="checkbox"/> 10" for 17"-18" H Back	\$25	
<input type="checkbox"/> 12" For 19"-22" H Back	\$25	
<input type="checkbox"/> 16" For 23"-25" H Back	\$25	
<input type="checkbox"/> TRU-Comfort™ Back <i>4-way Stretch Material Cover</i>		
Width: <input type="checkbox"/> 16" W <input type="checkbox"/> 18" W <input type="checkbox"/> 20" W	\$595	E2620
Overall Back: <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
Height ² : <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 17" W <input type="checkbox"/> 19" W	\$595	E2620
Overall Back: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Height ² : <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 22" W ³ <input type="checkbox"/> 24" W ³	\$695	E2621
Overall Back: <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Height ² :		

Tru Comfort Color Shroud Option:	
<input type="checkbox"/> Black <i>(Standard color option)</i>	
Metallic Colors: \$100	
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Silver
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Emerald Green
<input type="checkbox"/> Sapphire Blue	

Tilt Back Options:	
<input type="checkbox"/> Standard Back Canes	Standard
<input type="checkbox"/> Omit Push Handles	\$125
<input type="checkbox"/> Flip Forward Back Canes ⁴	\$150
<input type="checkbox"/> Omit Push Handles	\$125
<input type="checkbox"/> Offset (Meso) ⁵	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	
<input type="checkbox"/> Inset (Endo) ⁵	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	

1. Lateral extensions are required to mount lateral thoracic supports.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
3. Available with offset (meso) brackets.
4. Not available with Attendant Manual Recline, Vent Tray, or O2 Holder. May interfere with Lateral Supports and Waterfall Armrests.
5. Dual Post Flip-back Armrests interfere with offset canes. Not available with Manual Recline.

ARMRESTS	
Description	Price
Tilt Armrests:	
<input type="checkbox"/> Flip Back, Cane Mounted Armrests ¹	No Charge
<input type="checkbox"/> 2-Post, Flip-back, Height Adjustable, Removable ²	No Charge
<input type="checkbox"/> Single Post Height Adjustable Removable ³	No Charge
<input type="checkbox"/> Quad Release Lever	
<input type="checkbox"/> Removable, Single Post, Height Adjustable, Heavy Duty	No Charge
Height: <input type="checkbox"/> 8" H - 12" H <input type="checkbox"/> 10" H - 14" H	
<input type="checkbox"/> Pediatric, Removable ⁴	No Charge
<input type="checkbox"/> Quad Release Lever	
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁵	\$195
<i>Angle Adjustable: -30° through +30°</i>	

1. Interferes with Swing-Away laterals.
2. Height adjustment from 9"-13" in 1" increments. Interferes with swing-away laterals. Not available with a depth less than 16".
3. Height adjustable from 11" - 17" in 1/2" increments.
4. Height adjustment form 7.75" - 9.75" in 1/2" increments. Not available with full length armrests.
5. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward Back Canes, Cane Clamp Pushbar, or Shape Back. Call for details when selecting after market arm troughs.

ARMPADS		
Description	Price	
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Full Length (2"x14") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Desk Length (2"x10") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$160 pr.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Full Length Gel Armpad	\$56.69 ea.	
Mounting Plate (Required)*		
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$125 pr.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Desk Length Gel Armpad	\$53.52 ea.	

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

HEADRESTS		
Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
TRU-Comfort™ Headrest¹:		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350)	\$210	
Used in Addition to Removable Headrest Bracket		
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

1. Measurement is taken from the inside of the headrest.


NON-POWER LEGRESTS		
Description	Price	
<input type="checkbox"/> 60-90° Swing-away Legrests¹	Standard	
(Adjustable in 5° increments)		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Standard Extruded Footplates:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ²	\$150	
Angle Adjustable Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Angle-Adjustable Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		
<input type="checkbox"/> Multi-Axis Footplates	\$495	
Multi-Axis Extensions:		
<input type="checkbox"/> 15.5"-19.5"		
Multi-Axis Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		

1. Possible caster wheel interference when adjusted to less than 70 degrees.
 2. Includes Heel Loops
 3. Causes interference on seat widths less than 18".
 4. Causes interference on seat widths less than 20".

NON-POWER LEGRESTS (CONT'D)		
Description	Price	HCPCS Code
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests⁵	\$162	
<input type="checkbox"/> Angle Adjustable Footplates ⁶		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		
<input type="checkbox"/> Heavy Duty Drop-In Legrests⁵	\$420	
<input type="checkbox"/> Angle Adjustable Footplates ⁶	\$150	
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁷ <input type="checkbox"/> 6" x 11" ⁷ <input type="checkbox"/> 8" x 11" ⁸		
<input type="checkbox"/> Elevating Legrests⁵	\$325	E0990
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length:		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates ⁶	\$150	
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ⁷ <input type="checkbox"/> 6" x 11" ⁷ <input type="checkbox"/> 8" x 11" ⁸		
<input type="checkbox"/> Calf Pad Adjustment		
<input type="checkbox"/> 1 Hole <input type="checkbox"/> 3 Hole <input type="checkbox"/> 6 Hole		
<input type="checkbox"/> High Mount, Clamp-On Footrest*		
<input type="checkbox"/> Angle Adjustable		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

5. Extensions are measured from the seatpan to footplate.
 6. Includes Heel Loops
 7. Causes interference on seat widths less than 18".
 8. Causes interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM	
Description	Price
<input type="checkbox"/> Adjustable Center Mount Foot Platform	\$250
Includes 10"W x 11"D Split Foot Platform	
Extensions: ¹ <input type="checkbox"/> 10"-16" Ext. <input type="checkbox"/> 13"-19" Ext ²	
Calf Supports ³ :	
<input type="checkbox"/> Calf Pads (pair)	
<input type="checkbox"/> Small (6"H x 6"W)	\$150
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150
<input type="checkbox"/> Calf Panel (ea.)	
<input type="checkbox"/> Small (8"H x 18"W)	\$150
<input type="checkbox"/> Large (8"H x 24"W) ⁴	\$200
Optional Modular Footplate Extension Components/Kits ^{3,5}	
<input type="checkbox"/> 4" Width Extension Kit	\$50
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50



1. Extension length is measured from top of the seat pan to footplate.
 2. Not available on Recline only.
 3. Choose in Rigging Accessory Section in iQuantum.
 4. Only available on 20"W or greater.
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

ADJUSTABLE CENTER MOUNT FOOT PLATFORM (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/> 6-Piece Extension Stock Kit (Includes all above pieces)	\$95	
Additional Optional Accessories ^{1,2} :		
<input type="checkbox"/> Footplate Connector	\$49	

1. Choose in Rigging Accessory Section in iQuantum.
2. Assembly required. May require seating configuration adjustments to avoid front caster interference.

POWER POSITIONING ELECTRONICS

Description	Price	HCPCS Code
Single Actuator Electronics:		
<input type="checkbox"/> Tilt Only Through Toggle	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through VR2 6-Key (Tilt Only)	No Charge	E2310
<input type="checkbox"/> Through Single Switch (Remote+ and Tilt Only)	\$500	
<input type="checkbox"/> Through 2-Directional Control Box w/ Dual Switch Jacks	\$1,175	
<input type="checkbox"/> Through Remote Plus Controller	\$1,175	E2310
<input type="checkbox"/> Through NE+ Controller	\$1,175	E2310
Multiple Actuator Electronics:		
<i>One TRU Balance Power Positioning Function Plus Power Foot Rigging</i>		
<input type="checkbox"/> Tilt and Power AFP Through NE+	\$2,375	E2311
<input type="checkbox"/> Through Remote Plus Controller	\$2,375	E2311
<input type="checkbox"/> Remote Plus Through 4-Way Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Remote Plus Through 4-Way Button Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Remote Plus Through Quad Push-button	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		

CENTER MOUNT ARTICULATING FOOT PLATFORM

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Articulating Foot Platform	\$2895	K0108
<i>Includes 10"W x 11"D Split Foot Platform</i>		
<i>Extensions¹: <input type="checkbox"/> 10"-16" Ext.</i>		
<i>Calf Supports²:</i>		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ³	\$200	
<i>Optional Modular Footplate Extension Components/Kits^{2,4}</i>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/> 6-Piece Extension Stock Kit (Includes all above pieces)	\$95	
Additional Optional Accessories ⁴ :		
<input type="checkbox"/> Footplate Connector	\$49	

1. Extension length is measured from top of the seat pan to footplate.
2. Choose in Rigging Accessory Section in iQuantum.
3. Only available on 20"W or greater.
4. Assembly required. May require seating configuration adjustments to avoid front caster interference.

POWER ELEVATING LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 70° Articulating Swing-Away Legrests, 4" Articulation Range	\$2,395	E1010
<i>Quad Release Lever</i>		
<i>Calf Pads¹:</i>		
<input type="checkbox"/> 6" x 5.35"	<input type="checkbox"/> 8" x 5.86"	
<i>Angle Adjustable Extension Length²:</i>		
<input type="checkbox"/> 11"-13"	<input type="checkbox"/> 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"	<input type="checkbox"/> 5" x 6"	<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 8" ³	<input type="checkbox"/> 6" x 11" ³	<input type="checkbox"/> 8" x 11" ⁴
<i>Operational Method:</i>		
<input type="checkbox"/> Independent	<input type="checkbox"/> Combined	
<input type="checkbox"/> 70° Articulating Drop-In Legrests	\$2,395	E1010
<i>4" Range <input type="checkbox"/> 6" Range <input type="checkbox"/></i>		
<i>Calf Pads¹:</i>		
<input type="checkbox"/> 6" x 5.35"	<input type="checkbox"/> 8" x 5.86"	
<i>Angle Adjustable Extension Length²:</i>		
<input type="checkbox"/> 12"-14.5"	<input type="checkbox"/> 14.5"-18"	
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"	<input type="checkbox"/> 5" x 6"	<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 8" ³	<input type="checkbox"/> 6" x 11" ³	<input type="checkbox"/> 8" x 11" ⁴
<i>Operational Method:</i>		
<input type="checkbox"/> Independent	<input type="checkbox"/> Combined	

1. Choose in Rigging Accessory Section in iQuantum.
2. Extensions are measured from top of seat pan to footplate.
3. Only available with seat width 18" or greater.
4. Only available with seat width 20" or greater.

SYNERGY® CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion General Use		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
Solution@ 1 Cushion Skin Protection:		
<input type="checkbox"/> 14" x 14" - 20" x 20"	\$380	E2603
Solution@ Cushion Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$380	E2607
Spectrum@ Gel Cushion Positioning		
<input type="checkbox"/> 12" x 10" - 20" x 20"	\$480	E2605
Spectrum@ Visco Cushion Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$480	E2607
Spectrum@ Air Cushion Adjustable Skin Protection		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$480	E2622
Spectrum@ Air Contour Cushion Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 20"	\$499	E2624
TRU-Comfort™ Seat Cushion General Use		
<i>Includes Allure leather-like cover</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$380	E2601
TRU-Comfort™ Plus Seat Cushion Skin Protection		
<i>Includes 4-way stretch material cover</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$380	E2603

FOOT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) ¹	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$65	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$512	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair)*	\$574	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Bodypoint Ankle Huggers (pair)*	\$134	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium		

1. Standard with angle adjustable footplates.
2. Available only with swing-away legrests.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

LATERALS **LAP BELTS**

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)*		
<i>Bracket:</i>		
<input type="checkbox"/> 125° Swing-Away Bracket	\$450	E1028
<i>For Use on TRU-Comfort Backs w/ Recline</i>		
<input type="checkbox"/> Swing-Away Rail Mounted	\$450	E1028
<i>For Use on Synergy Back and Tru-Comfort (Synergy and Tilt Only)</i>		
<i>Pad:</i>		
<input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals*		
<i>Side(s):</i>		
<input type="checkbox"/> Left Side		<input type="checkbox"/> Right Side
<i>Choose Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$250	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$225	E1028
<i>Choose Pad:</i>		
<input type="checkbox"/> Short Pad, 4" x 4" (ea.)	\$50	E0956
<input type="checkbox"/> Medium Pad, 4" x 8" (ea.)	\$50	E0956
<input type="checkbox"/> Long Pad, 4" x 12" (ea.)	\$50	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ²	\$25	K0108

1. Swing-away Laterals Interfere with flip-up armrests.
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.

Description	Price
<input type="checkbox"/> 38" Lap Belt	Standard
<input type="checkbox"/> 50" Lap Belt	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Padded Lap Belt	\$115
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$87
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width	\$92
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width*	\$112
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width*	\$116

BATTERY OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> NF-22, 50 AGM	\$255	E2361
<input type="checkbox"/> NF-22, 50Gel Cell	\$295	E2361

MISCELLANEOUS POSITIONING COMPONENTS

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket*	\$172	E1028
<i>Pad:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$135	E0957
<input type="checkbox"/> Stealth Foot Boxes (pair)*	\$475	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)*	\$166	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/> Otto Bock Arm Channels (ea.)*	\$104	E2209
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Otto Bock Flat Hand Pad(ea.)*	\$52	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Medium <input type="checkbox"/> Large		

Please see positioning components order form for additional types, sizes and details.

SEAT TO FLOOR HEIGHT

Seating System	Seat-To-Floor Heights			
TRU-Balance® Tilt	18.125"	19.125"	20.125"	21.125"

Seat to floor height are +/- 1/8" with 18"W x 18"D with Flat-Free tires (add 3/8" for pneumatic tires). Increased widths may increase lowest seat to floor heights.

NOTES:

ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Quantum Backpack ¹	\$60	
<input type="checkbox"/> Oxygen Holder ²	\$150	E2208
<input type="checkbox"/> Weather Cover, Medium, Black	\$158	
<input type="checkbox"/> Cup Holder	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Tablet Holder	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder	\$140	
<input type="checkbox"/> Cell Phone/Media Holder	No Charge	

1. Backpack may interfere w/ other rear accessories.
 2. Oxygen holder on Tilt with Tru Comfort back requires SOF. N/A on Power Recline.