


EFFECTIVE NOVEMBER 15TH 2006.

Jazzy® 610 Series
 Order Form / Price List
 Revised: February 2007

Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636



This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:		Address:		Discount:
Phone Contact:		Marked For:		Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

<input type="checkbox"/> JAZZY 610 2S-C CONFIGURATION #2-1	
<input type="checkbox"/> Jazzy 610 2S-C Red Group 2 / Standard K0823 w/ captain's seating 300 lbs. weight capacity	\$5,526 -18" W x 18" D Gray vinyl medium back seat -Full length armrests -75 amp Flight controller -10" Solid tires -4 amp on-board charger -60" lapbelt -Foot platform
<input type="checkbox"/> U-1(pair) Installed	\$224

<input type="checkbox"/> JAZZY 610 2S-C CONFIGURATION #2-2	
<input type="checkbox"/> Jazzy 610 2S-C Red Group 2 / Standard K0823 w/ captain's seating 300 lbs. weight capacity	\$5,801 -18" W x 18" D Ltd. Recline Gray vinyl high-back -Full length armrests -75 amp Flight controller -10" Solid tires -4 amp on-board charger -60" lapbelt -Foot platform
<input type="checkbox"/> U-1(pair) Installed	\$224

<input type="checkbox"/> JAZZY 610 2S-SS CONFIGURATION #2-3	
<input type="checkbox"/> Jazzy 610 2S-SS Red Group 2 / Standard K0822 w/ sling/solid seating *Select cushion for the solid seat pan in appropriate section or review specific cushion order form. 300 lbs. weight capacity	\$5,801 -18" W x 18" D Ltd. Recline Gray vinyl high-back seat with solid seat pan -Full length armrests -75 amp Flight controller -10" Solid tires -4 amp on-board charger -60" lapbelt -Foot platform
<input type="checkbox"/> U-1(pair) Installed	\$224

<input type="checkbox"/> JAZZY 610 2SP-SS (single power) REFERENCE GUIDE #2-4	
<input type="checkbox"/> Jazzy 610 2SP-SS Group 2 / Single Power K0835 Seat Options: <input type="checkbox"/> TB Tilt E1002 250 lbs. weight capacity	\$5,726 -Flip-back, cane mounted, full length arms -50 amp, 6-Key PG VSI controller -10" Solid tires -4 amp on-board charger -60" lapbelt -Swing-away legrests -Reinforced footplates -Tilt through 6-key VSI <input type="checkbox"/> Omit Solid Tires (-\$213) <input type="checkbox"/> Omit 6-Key VSI (-\$200)

<input type="checkbox"/> JAZZY 610 2SP-C (single power) REFERENCE GUIDE #2-5	
<input type="checkbox"/> Jazzy 610 2SP-C Group 2 / Single Power K0836 250 lbs. weight capacity *For use when a single power option is added at a later date.	5,726 -18" W x 18" D Gray vinyl med-back seat -Full length armrests -50 amp 6-key PG VSI -10" Solid tires -4 amp on-board charger -60" lapbelt -Foot Platform <input type="checkbox"/> Omit Solid Tires (-\$213) <input type="checkbox"/> Omit 6-Key VSI (-\$200)

SPECIFIC CONFIGURATION

<input type="checkbox"/> Jazzy 610	\$5,313
<input type="checkbox"/> Jazzy 610 2S-C #2-1	\$5,526
<input type="checkbox"/> Jazzy 610 2S-C #2-2	\$5,801
<input type="checkbox"/> Jazzy 610 2S-SS	\$5,801
<input type="checkbox"/> Jazzy 610 2SP-SS	\$5,726
<input type="checkbox"/> Jazzy 610 2SP-C	\$5,726
<input type="checkbox"/> Omit contoured seat ²⁴	(-\$275)

When configuring power chair use reference guides to determine accessories included in the MSRP pricing package

Model number with suffix must be used for all chairs billed through Medicare.

COLOR

<input type="checkbox"/> Champagne	Std.
<input type="checkbox"/> Candy Apple Red	\$52
<input type="checkbox"/> Viper Blue	\$52
<input type="checkbox"/> Onyx Black	\$52

BATTERIES (2 REQUIRED)

<input type="checkbox"/> U-1, Installed	E2365	\$224pr
---	-------	---------

Jazzy® 610 SeriesOrder Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

CONTROLLER

- Pride Flight 75 amp¹ Std.
Cannot be used w/ power positioning systems.

PG Electronics

- VSI Integral 50 amp^{1,16,28} NC
 VSI Integral 6 Key 50 amp^{1,3} \$200
(Applicable for use w/ Tilt or PES)

Expandable Controllers:

(Applicable for use w/ Tilt or PES)

- 3-Key Remote Plus 70 amp^{1,3} E2377 \$1025
(Select Box for Profile)
 1 Profile Multiple Profile
- 7-Key Access. Remote Plus
70 amp^{1,3,4} E2377 \$1300
- Remote Plus 70 amp with
Speed Dial E2377 \$1300
- Remote Plus 70 amp
with Toggles^{1,3} E2377 \$1300
- Remote Plus 70 amp with
Switch Jacks^{1,3} E2377 \$1300
- Microdrive, 70 amp^{1,3,6} E2377 \$1300

PROGRAMMER

- Flight Programmer \$423
- PG Programmer \$423

CONTROLLER MOUNT

- Fixed, In-Line Joystick Mount Std.
 Left Right
- Swing-away Joystick Mount²⁷ E1028 \$196
 4" L²⁹ 7" L
 Left Right
- Quick Disconnect, Joystick Mount³³ NC
 Left Right
- Multi-axis Joystick Mount K0108 \$100
 Left Right
(Fixed, does not swing-away)
- Vertical Adj. Clamp⁷ K0108 \$95
Provides vertical joystick adjustment
w/o adjusting armrest height.

JOYSTICK HANDLE OPTIONS**Flight Controls**

- Flight Knob Std.
- Sponge Ball E2323 \$35
- Chin Cup E2323 \$30
- DL Knob E2323 \$20
- Carrot Knob E2323 \$20
- Straight Bar (Chrome) E2323 \$45
- Finger Tube (Chrome) E2323 \$45
- Round T-Bar (Chrome) E2323 \$45
- Flat T-Bar (Chrome) E2323 \$45
- Wrist Plate E2323 \$45

PG Controls

- VSI Knob Std.
- Mushroom Handle E2323 \$71
- T-Bar E2323 \$107
- Chin Cup E2323 \$71
- Extended I-Handle E2323 \$71
- Soft Ball Handle E2323 \$50
- Medium Ball Handle E2323 \$50
- Goal Post Handle E2323 \$66
 3" W 4" W

FRONT RIGGINGS

- Foot Platform Std.
- Rigging Type (pair)**
- Swing-away footrests²⁵ \$162
- Elevating Leg Rests E0990 \$325

Select footplate size for above:

- 7 1/4 x 6 8 1/2 x 6

Adjustable Calf Pads w/ ELRS (1 hole std.)

- 1 hole 3 holes* 6 holes*

*Potential for interference w/ shrouds

- HD Drop-Ins, 70^{07,13,22} K0108 \$420
(HD Drop-ins N/A w/ contoured or Specialty seating)

- Power-elevating, high-pivot^{30,22} \$1695
Legrests, (calf pads included)³² E1010
- Independent / Combined Combined³⁰
- Power-elevating, Articulating^{30,31,32} E1010 \$2395
Articulating Range: 4" L Std.
 6" L (Drop-in only) NC
- Set at:** 70°
- Independent / Combined³⁰ Combined³⁰
- Swing-away Power-elevating, E1010
Articulating^{25,30,31} \$2395
Articulating Range: 4" L Std.
 6" L (Drop-in only) NC
- Set at:** 70°
- Quad Release Lever NC
- Independent / Combined³⁰ Combined
- Calf Pads Left & Right (6 x 5.35) Std.
- Calf Pads Left & Right (8.00 x 5.86) NC

Footplates (pair)

- Reinforced Footplate Std.
 7 1/4 x 6 8 1/2 x 6
- 5"x8" Angle-Adj. Footplate¹⁰ \$150
- 6"x8" Angle-Adj. Footplate⁹ \$150
- 6"x11" Angle-Adj. Footplate⁹ \$150
- 8"x11" Angle-Adj. Footplate¹¹ \$150
- Wheel Bumpers (Ang-adj. only) NC

Extensions¹²

- 11"-14" Ext. NC
- 14"-17" Ext. Std.
- High mount Clamp-on; Ang. Adj. \$326
 Small Large
- High mount Clamp-on; Multi-axis \$412
 Small Large

- High mount footplatform²⁶ K0108 \$450
 Mid-wheel 6 foot platform
Lower Extensions
- 10"-14" L
- 14"-17" L

Rigging Options

- Heel Loops (pair)¹³ E0951 \$20
- Gel Neoprene Footrest
Sleeve (pair)⁸ K0108 \$65
- ELR / ALR Adductor Buttons (pr) K0108 \$100
- POSAlinc Knee Adductor (pair)
for S/A legrests K0108 \$464
- Residual limb support left E1020 \$300
- Residual limb support right E1020 \$300

Jazzy® 610 Series

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

DRIVE WHEELS

<input type="checkbox"/> 10" Pneumatic (pair) ²⁸	NC
<input type="checkbox"/> 10" Flat-Free (pair)	\$213

SEATING

Med-Back Std.

Med-Back & Solid Seat Pan¹⁴ Std.

18" W x 18" D

Gray Vinyl Beige Vinyl

High-Back w/ Headrest \$275

High-Back w/ Headrest & Solid Seat Pan¹⁴ \$275

18" W x 18" D

18" W x 20" D

Gray Vinyl Beige Vinyl

**Custom seat dimensions are available.
Please contact the Custom Group for
pricing and leadtime.**

ARMRESTS

For Medium and High-Back Seats

Flip-up, removable Std.

Flip-up, height-adjustable removable²⁰ \$170
E0973

Armpads

Desk Length, 10" NC

Left NC Right NC

Full Length, 14" NC

Left Std. Right Std.

POWER ELEVATING SEAT^{16,17,18,19}

Power Elevating Seat \$975

Through 6-key VSI NC

Lap Belt Std.

MANUAL RECLINER¹⁵

Manual Recline High Back E1226 \$600

Manual Recline High Back¹⁴
w/ Solid Seat Pan E1226 \$600
(180 degrees of Recline)

(Dimensions 18"W x 18"D)

Gray Vinyl Beige Vinyl

Width: **Depth:**

16" NC 16" NC

18" Std. 18" Std.

20" NC 20" NC

**When ordering a lap belt with this seat, an
additional 10" is added to accommodate**

CUSHIONS³

Simplicity Cushion
16 x 16 - 20 x 20 E2601 \$98

Solution Cushion
10 x 10 - 20 x 18 E2607 \$380

20 x 20 - 20 x 20 E2607 \$430

22 x 16 - 22 x 20 E2608 \$430

Spectrum Cushion
 Gel

10 x 10 - 20 x 18 E2605 \$480

20 x 20 - 20 x 20 E2605 \$530

22 x 16 - 22 x 20 E2606 \$530

Visco

10 x 10 - 20 x 18 E2607 \$480

20 x 20 - 20 x 20 E2607 \$530

22 x 16 - 22 x 20 E2608 \$530

Structure Cushion
 Gel Visco

10 x 10 - 20 x 18 E2607 \$560

20 x 20 - 20 x 20 E2607 \$610

22 x 16 - 22 x 20 E2608 \$610

TRU-Comfort Cushion
16 x 16 - 20 x 18 \$380

20 x 20 - 20 x 22 \$380

22 x 16 - 22 x 24 \$480

Footnotes for pages 1, 2 & 3

1. Std. w/ one programmable drive profile.
2. See Synergy Cushion order forms for detailed pricing, sizes & codes.
3. Needed for power positioning systems.
4. Needed for lighting packages.
5. Avail. w/ five individual drive profiles.
6. Requires multi-axis joystick mount; N/C for mount.
7. N/A on Specialty, contoured seating or waterfall armpads.
8. Avail. w/ S/A footrests & 70 degree HD drop-ins.
9. Only avail. on 16" W or greater seats.
10. Only avail. on 14" W or greater seats.
11. Only avail. on 20" W or greater seats.
12. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
13. Std. w/ ang-adj. footplates.
14. Includes black vinyl seat pan, back available in all listed colors.
15. 18" x 18" std. Special measurements are N/A.
16. N/A w/ power positioning systems
17. Includes 1/2 speed inhibit.
18. STF will increase.
19. Manual recliner is N/A w/ PES.
20. N/A on manual recliner.
21. Rear mount accessories are N/A.
22. May cause interference with front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
23. Applicable when ordering no seat, specialty, Synergy seat or any power positioning system.
24. Applicable when ordering no seat, specialty, Synergy seat or any power positioning system.
25. Please note S/A footrests are standard when ordered with Power Positioning.
26. Std. calf pad is not depth adjustable.
27. Can only mount one armrest accessory in combination w/ S/A joystick.
28. If ordered w/ reference guides or configuration packages discounts may apply.
29. N/A with waterfall armpads and dual post armrests.
30. Requires Remote Plus.
31. Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
32. Std. calf pad is not depth adjustable.
33. N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.

Jazzy® 610 Series

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) \$275

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" NC
- 18" Std.
- 20" NC

Specialty Seating Kit Depth

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | |

Specialty Seating Kit Back Cane Height

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 20" NC |
| <input type="checkbox"/> 22" NC | |

Back Selection

- See Synergy Shape Back Order Form^{1,2,3}
- Solid Plate Back up to 20"W¹⁹ \$275
- See Positioning Components Order Form for additional backs

Back Height Solid Pan⁴

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

Sling Back

- 16" wide \$98
- 18" wide \$98
- 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
- 20" height NC
- 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

- 2-Post, flip-back⁶ Std.
- Flip-back, cane mounted^{5,6} \$95ea
- Removable, adjustable, quick height^{3,20} E0973 \$275ea
- Removable, height adjustable²⁰ heavy duty armrests E0973 \$150ea
- 8"H-12"H 10"H-14"H
- Pediatric, removable,²⁰ adjustable, quick height E0973 \$250ea

Armrest Pads

- Full Length, straight, 2"x14" NC
 - Left Right
- Desk Length, straight, 2"x10" NC
 - Left Right
- Full Length, Waterfall, 2"x14"⁷ NC
 - Left Right
- Desk Length, Waterfall, 2"x10"^{7,8} NC
 - Left Right

HEADRESTS⁹

- Stealth^{10,11} E0955
 - Large (11"wide x 6" high) \$210
 - Small (9 3/4" wide x 3" high) \$210
 - Stealth, swing-away (Flip-Down)
 - Head Rest Bracket E1028 \$150
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications^{10,11} E0955
 - Large (11"wide x 6" high) \$275
 - Small (9 3/4" wide x 3" high) \$275
 - Fixed, Offset Mounting Bracket NC

Synergy Seating Systems¹²

- Synergy Seating™ System \$434
(Price with Base)
- Synergy Manual Recline Pkg.^{13,14,15,16} \$700
Maximum weight 300 lbs. E1226
- Synergy Manual Tilt Pkg.^{14,17} \$550
(Max Dimensions 18"W x 18"D) K0108
Maximum weight 160 lbs.
N/A w/ elevating seat
- Synergy Quick Release¹⁸ \$200
(Max Dimensions 20"W x 20"D)
Std. rear mount accessories are N/A

Synergy Quick Release is in addition to the MSRP of the Synergy Seat.

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 10" NC | <input type="checkbox"/> 11" NC |
| <input type="checkbox"/> 12" NC | <input type="checkbox"/> 13" NC |
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | |

Seat Width Package

- 10"W-14"W 14"W-18"W
- 16"W-20"W

Seat Width Package Manual Tilt

- 10"W-12"W 12"W-14"W
- 14"W-16"W 16"W-18"W

DEPTH

Set Seat Depth At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 10" NC | <input type="checkbox"/> 11" NC |
| <input type="checkbox"/> 12" NC | <input type="checkbox"/> 13" NC |
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | |

Seat Depth Package

- 10"D-14"D 14"D-18"D
- 16"D-20"D

Seat Depth Package Manual Tilt

- 10"D-14"D 14"D-16"D
- 16"D-18"D

Footnotes for pages 4:

1. 14" & 22" width are N/A.
2. Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
3. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
4. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
5. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
6. S/A laterals interfere w/ flip-back arms.
7. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
8. S/A joystick mount N/A.
9. Only available w/ solid backs.
10. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
11. Comes with swing-away / removable hardware; HCPCS code is E1028.
12. Manual systems do not include a headrest.
13. N/A w/ flip-back, cane mounted arms.
14. Solid back not included.
15. Includes removable, adjustable, quick-height armrests.
16. N/A w/ Shape Back.
17. Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
18. N/A w/ Synergy manual recline.
19. 20" D is N/A.
20. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Jazzy® 610 Series

Order Form / Price List
Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.

BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

BACK OPTIONS

See Synergy Shape Back Order Form^{2,3,4}

Solid Plate Back up to 20" W \$275

Back Height:¹

<input type="checkbox"/> 18"H	<input type="checkbox"/> 19"H	<input type="checkbox"/> 20"H	<input type="checkbox"/> 21"H
<input type="checkbox"/> 22"H	<input type="checkbox"/> 23"H	<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H

Sling Back^{1,5} \$98

Back height will be the same as back can height. Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back¹⁶

Back Height:	<input type="checkbox"/> 22"	<input type="checkbox"/> 25"
<input type="checkbox"/> Back up to 20" W		\$595
<input type="checkbox"/> 22" W		\$695

ARMRESTS

Flip-back, cane mounted^{6,7,8} Std.

2-Post, flip-back, height^{6,7,9} E0973 \$150ea adjustable, removable armrests

Removable, adjustable^{10, 19} E0973 \$75ea quick height

Quad Release Lever NC

(3" of height adjustment in 3/8" increments)

Removable, height adjustable E0973 \$75ea heavy duty armrests¹⁹

8"H-12"H 10"H-14"H

Pediatric removable¹⁹ E0973 \$75ea adjustable, quick height.

Match Shroud Color Black

Quad Release Lever NC

(7 3/4" - 10 3/4" height adjustment in 1/2"

increments from seat pan to top of armrest)

ARMPADS

Full length, Straight, 2" X 14" Std.

Left Right

Desk Length, Straight, 2" X 10" Std.

Left Right

Full Length, Waterfall, 2" X 14"¹¹ NC

Left Right

Desk Length, Waterfall, 2" X 10"^{11,12} NC

Left Right

HEAD REST¹⁵

Stealth^{13,14} E0955

Large (11" wide x 6" high) \$210

Stealth^{13,14} E0955

Small (9 3/4" wide x 3" high) \$210

Stealth, swing-away (Flip-Down)

Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications^{13,14} E0955

Large (11" wide x 6" high) \$275

Medical Applications^{13,14} E0955

Small (9 3/4" wide x 3" high) \$275

Fixed, Offset Mounting Bracket NC

TRU-Comfort^{13,14} E0955

8" Wide w/ Cobra hardware \$417

10" Wide w/ Cobra hardware \$417

Bracket for after market headrest \$180

TRU-BALANCE 2 TILT

TRU-Balance 2 Power Tilt E1002 \$4995

55 degree center of gravity (CG) power tilt seating system. 250lbs weight capacity.

SYSTEM SPECIFICATIONS

Set Seat Width At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Width Package

<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 16"W-20"W
------------------------------------	------------------------------------

Set Seat Depth At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Depth Package

<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 16"D-20"D
------------------------------------	------------------------------------

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

Set Back Cane Height At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

Seat Overall Back Height¹

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

Back Type

See Synergy Shape Back Order Form

Solid Plate Back Std.

Attendant manual recline^{16,17,18} \$750

180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.

TRU-Comfort Back¹⁶

Back Height: 22" 25"

Back up to 20" W \$595

22" W \$695

Omit Back -(\$100)

Footnotes for page 4:

1. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
2. 14" & 22" width are N/A.
3. Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
4. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
5. N/A w/ Synergy manual systems.
6. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
7. S/A laterals interfere w/ flip-back arms.
8. N/A w/ Synergy manual recline.
9. May cause interference with certain dimension configurations. Call Quantum Customer Service.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
11. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
12. N/A w/ S/A joystick mount.
13. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
14. Comes with swing-away / removable hardware; HCPCS code is E1028.
15. Available w/ solid backs only.
16. N/A w/ flip-back, cane mounted arms.
17. N/A w/ Shape Back.
18. Includes price for removable, adjustable, quick-height arms & 2-post arms.
19. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Jazzy® 610 Series

Order Form / Price List

Revised: February 2007

EFFECTIVE NOVEMBER 15TH 2006.**ARMRESTS**

<input type="checkbox"/> Flip-back, cane mounted ^{1,2,3}	Std.
<i>Mounted on back canes with full width spreader bar.</i>	
<input type="checkbox"/> 2-Post, flip-back, hgt. adj, removable ^{1,2}	NC
<input type="checkbox"/> Removable, adj. quick height ^{4,19}	NC
<input type="checkbox"/> Quad Release Lever	NC
<i>(3" of height adjustment in 3/8" increments)</i>	
<input type="checkbox"/> Removable, single-post, ¹⁹ height adjustable, heavy duty	NC
<input type="checkbox"/> 8"-12" Adjustment	
<input type="checkbox"/> 10"-14" Adjustment	
<input type="checkbox"/> Pediatric, removable, ¹⁹ adjustable, quick height	NC
<input type="checkbox"/> Match Shroud Color <input type="checkbox"/> Black	
<input type="checkbox"/> Quad Release Lever	NC

Armrest Pads

<input type="checkbox"/> Full Length, straight, 2"x14"	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Desk Length, straight, 2"x10"	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Full Length, Waterfall, 2"x14" ¹⁵	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Desk Length, Waterfall, 2"x10" ^{15,6}	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	

POWER POSITIONING ELECTRONICS**PG Electronics**

<input type="checkbox"/> Through Toggle Box ^{6,7}	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right (mounts to side of armrest)	
<input type="checkbox"/> Through Single Switch Tilt Module ^{8,9} E2310	\$500
<input type="checkbox"/> Through 6 Key VSI Controller ¹¹	NC

Single Motor Controller:¹¹

<input type="checkbox"/> Through Remote Plus Controller E2310*	\$1175
<input type="checkbox"/> Through ALM w/ 2 Jacks ¹⁰ E2310*	\$1175

*E2310 coding may vary due to configuration

HEAD REST

<input type="checkbox"/> Stealth ^{12,13} E0955	
Large (11" wide x 6" high)	\$210
<input type="checkbox"/> Stealth ^{12,13} E0955	
Small (9 3/4" wide x 3" high)	\$210
<input type="checkbox"/> Stealth, swing-away (Flip-Down)	
Head Rest Bracket E1028	\$150
<i>Interfaces w/ any headrest that uses 1/2" square down-tube.</i>	
<input type="checkbox"/> Medical Applications ^{12,13}	
Large (11" wide x 6" high) E0955	\$275
<input type="checkbox"/> Medical Applications ^{12,13}	
Small (9 3/4" wide x 3" high) E0955	\$275
<input type="checkbox"/> Omit Head Rest - To be supplied by provider	
<input type="checkbox"/> Fixed, Offset Mounting Bracket	NC
<input type="checkbox"/> TRU-Comfort ^{12,13} E0955	
<input type="checkbox"/> 8" Wide w/ Cobra hardware	\$417
<input type="checkbox"/> 10" Wide w/ Cobra hardware	\$417
<input type="checkbox"/> Bracket for after market headrest	\$180

ACCESSORIES

<input type="checkbox"/> Rear Basket ^{14,15,16}	\$103
<input type="checkbox"/> Cane/Crutch Holder ^{14,15,16} E2207	\$98
<input type="checkbox"/> Oxygen Holder ^{14,15,16} E2208	\$103
<input type="checkbox"/> Oxygen Holder ^{14,15,17,18} E2208	\$150
<i>(TRU-Balance Tilt & Synergy Seat only)</i>	
<input type="checkbox"/> Stroller handle	\$85
<i>(Synergy Manual Tilt only)</i>	
<input type="checkbox"/> Walker Holder ^{14,15,16} K0108	\$92
<input type="checkbox"/> Flag Holder (SOF on rehab seating)	\$17
<input type="checkbox"/> Weather Cover, Medium, Black	\$158
<input type="checkbox"/> Saddle bag(contoured seating only)	\$28
<input type="checkbox"/> Cup holder	\$17
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Lap Belt E0978	Std.
<input type="checkbox"/> 38" ¹⁶ <input type="checkbox"/> 50" <input type="checkbox"/> 60"	
<input type="checkbox"/> 70" <input type="checkbox"/> 80" <input type="checkbox"/> 90"	

Footnotes for page 6:

1. N/A w/ flip-back, cane mounted arms.
2. N/A w/ Shape Back.
3. Includes price for removable, adjustable, quick-height arms & 2-post arms.
4. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
5. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
6. N/A w/ S/A joystick mount.
7. N/A w/ waterfall armpads.
8. Switch not included.
9. PG electronics only.
10. Consult specialty control order form for a wide variety of switch types and sizes-- 2 required.
11. N/A w/ power elevating seat.
12. Comes with swing-away / removable hardware; HCPCS code is E1028.
13. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
14. Includes mounting bracket and hardware.
15. N/A w/ TB Tilt.
16. N/A w/ manual recline high back.
17. N/A w/ flip-bac arms.
18. Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
19. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Lowest seat-to-floor height	
Specialty Seat:	16.75"
Synergy Seat:	16.625"
TB 2 Tilt:	16.875"
High-back Seat:	20.5"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service
Includes: liftgate service and delivery of your product(s) inside your facility
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.