



A DIVISION OF
QUANTUM
LIFE BEYOND LIMITS™

May 2014

Kids ROCK™ Size 2 Order Form

Maximum Weight Capacity 115 lbs. 14"Wx16"Dx22"H Seating Capacity.

Customer Service: 800-800-8586 (toll free) Email: orders@pridemobility.com

Date: _____ <input type="checkbox"/> Quote <input type="checkbox"/> Order Provider Acct #: _____ Provider: _____ Provider Contact: _____ Provider Address: _____ Provider City: _____ ST: _____ ZIP: _____ Provider Phone: () _____ Fax: () _____ Confirmation Email: _____	SHIPPING INFORMATION PO#: _____ Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: () _____ Fax: () _____
--	--

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

MEASUREMENT WORKSHEET

Mark For: _____
 Age: _____
 Height: _____ ft. _____ in.
 Weight: _____ lbs.




Figure 1 - Client Measurements

- A. Sitting Depth _____
- B. Trunk Depth _____
- C. Trunk Width _____
- D. Hip Width _____
- E. Top of Shoulder to Seat _____
- F. Lower Leg Length _____

Kids ROCK Size 2 chair, Tilt-in-Space

Includes: Tilt in space, non-dynamic, mobility base, Silver metallic frame color, standard configuration.

<input checked="" type="checkbox"/>	E1233	Kids ROCKi Chair - Size 2	\$2,650.00
-------------------------------------	-------	---------------------------	------------

Transit Option

3222-5100 Seating System and/or 3222-5200 Back System Required for ISO 7176-19 Compliance. Transit Option must be factory installed and is not retrofitable.

<input type="checkbox"/>	Transit Option	\$200.00
--------------------------	----------------	----------

Upholstery Color

Upholstery Color applies to: back cushions, lateral pelvic/thigh supports, calf protector

<input type="checkbox"/>	Black	<input type="checkbox"/>	Pink/Black
<input type="checkbox"/>	Vibrant Red/Black	<input type="checkbox"/>	Blue/Black
<input type="checkbox"/>	Green/Black		

Reaction™ Dynamic Seating



Requires 3222-5800 Angle Adjustable Leg Rests

<input type="checkbox"/>	E2295	Kids Rock™ Active Component	\$695.00
--------------------------	-------	-----------------------------	----------

Upper Body Springs

<input type="checkbox"/>	Soft - Pair	STD
<input type="checkbox"/>	Medium - Pair	N/C
<input type="checkbox"/>	Hard - Pair	N/C
<input type="checkbox"/>	No Springs	N/C

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

Front Rigging

3222-5800 Required with REACTION™ System. Calf Protector included with all Front Rigging Option.

<input type="checkbox"/>	3222-5700		90° Fixed Front Rigging	STD	Lower Body Springs			
<input type="checkbox"/>	3222-5800	E0990	Angle Adj. Elevating Leg Rest	\$225.00	<input type="checkbox"/>	3222-7310C	Soft - Pair	STD
					<input type="checkbox"/>	3222-7350C	Medium - Pair	N/C
					<input type="checkbox"/>	3222-7380C	Hard - Pair	N/C
					<input type="checkbox"/>	3222-7300	No Springs	N/C

Foot Plate

Select Foot Plate Installation: Std(10.5-14".5 Leg Drop) Inverted(7.5-11" Leg Drop) Long(14"-16" Leg Drop)

<input type="checkbox"/>			90° Fixed Foot Plate	N/C	<input type="checkbox"/>	K0040	Angle Adjust Flip Up Foot Plate	\$165.00
<input type="checkbox"/>	K0040		Bolt Angle Adjustable Foot Plate	\$140.00	<input type="checkbox"/>	K0040	Quick Angle Adj. Foot Plate -18° to +45°Dorsal/Plantar Flexion; 8 Positions; Easy Pull Adj.	\$165.00
±8° Dorsal/Plantar Flexion; -8° - 0 - +8°								

Kids ROCK 2 Adjustable Depth Seating System

<input checked="" type="checkbox"/>	K0108	8-17" Adj. Depth Seat Pan	\$195.00
-------------------------------------	-------	---------------------------	----------

Kids UP LifeForm™ Contoured Seat Cushion Selections

<i>Seat Cushions Incl: 2 Positioning Wedges, Abductor Wedge; Black Mesh Cover</i>				Ischial Insert Selection Required			
<input type="checkbox"/>	E2294	8"-11" Depth	\$365.00	<input type="checkbox"/>		Open Cell Foam	STD
<input type="checkbox"/>	E2294	11"-14" Depth	\$365.00	<input type="checkbox"/>	K0108	Soft Memory Foam	\$50.00
<input type="checkbox"/>	E2294	14"-17" Depth	\$365.00	<input type="checkbox"/>	K0108	Medium Memory Foam	\$50.00

Kids UP LifeForm™ Seating System Accessories

<input type="checkbox"/>		Extra Positioning Wedges (4)	\$35.00
<input type="checkbox"/>		Pressure Relief Layer - Sunmate (14" x 16" x 1/2")	\$60.00
<input type="checkbox"/>		Wedge for Seat Cushion (14" x 8" x 1")	\$40.00

Synergy® Seat Cushion Options

<input type="checkbox"/>	E2607	Synergy® Solution Cushion 14"x12"	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution Cushion 14"x14"	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution Cushion 14"x16"	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution Cushion 14"x18"	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution 1 Seat Cushion 14"x14" W/Gel Layer	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution 1 Seat Cushion 14"x16" W/Gel Layer	\$380.00
<input type="checkbox"/>	E2607	Synergy® Solution 1 Seat Cushion 14"x18" W/Gel Layer	\$380.00

Back System Options

<input type="checkbox"/>	E2293	Configurable, Contoured, Adjustable Height Solid I-Back System	\$585.00
<input type="checkbox"/>		Back Canes Only. (Not ISO 7176-19 Compliant)	STD

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

Kids UP Lifeform™ Back Cushion Options

Select 1 if 3222-5200 is Selected. Back Cushions include contouring wedges and lumbar support

<input type="checkbox"/>	12"-17" Contoured Back Cushion. (I-Back)	STD
<input type="checkbox"/>	16"-22" Contoured Back Cushion. (I-Back)	STD
<input type="checkbox"/>	Omit Back Cushion	Less \$135.00

Back Cushion Embroidery

(Embroidery is applied to colored panel)

<input type="checkbox"/>	Block Letters	<input type="checkbox"/>	<i>Script Letters</i>
--------------------------	----------------------	--------------------------	-----------------------

Inscription, Please place one letter per box

--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

<input type="checkbox"/>	Blue Thread	<input type="checkbox"/>	Green Thread
<input type="checkbox"/>	Red Thread	<input type="checkbox"/>	Purple Thread
<input type="checkbox"/>	Black Thread	<input type="checkbox"/>	White Thread

Armrest Options

<input type="checkbox"/>	Basic Fixed Height Armrest Thin Profile - Pair	STD	<input type="checkbox"/>	E0973 Adjustable Height Armrest Full Length Pad - Pair *	\$195.00
<input type="checkbox"/>	Fixed Height Armrest Desk Length Pad - Pair *	\$170.00	<input type="checkbox"/>	E0973 Adjustable Height Armrest Thin Profile Pad - Pair	\$195.00
<input type="checkbox"/>	E0973 Adjustable Height Armrest Desk Length Pad - Pair *	\$195.00	<i>Armrests with * not for use with 22" or 24" Rear Wheels</i>		

Tray Options

<input type="checkbox"/>	E0950	Upper Extremity Positioning Platform ~ Color - Clear	\$128.00	<input type="checkbox"/>	E0950	Upper Extremity Positioning Platform ~ Color - Gray	\$128.00
<input type="checkbox"/>		Angle/Depth Adj. Support Hardware - Req. for all trays	\$144.00				

Rear Wheels

<input type="checkbox"/>	16" Composite Rear Wheels	STD	<input type="checkbox"/>	24" Composite Rear Wheels (only avail. w/ 8" Casters). <i>Requires 3222-5600</i>	N/C
<input type="checkbox"/>	22" Composite Rear Wheels. <i>Requires 3222-5600</i>	N/C			

Rear Tires

<input type="checkbox"/>	Full Polyurethane w/Treads	STD			
<input type="checkbox"/>	Pneumatic	STD			
<input type="checkbox"/>	E2213 Pneumatic w/Airless Inserts	\$88.00			

Hand Rims

Available with 22" and 24" rear wheels only.

<input type="checkbox"/>	Chrome Hand Rims (Not Avail. With 16" Wheels)	N/C	<input type="checkbox"/>	K0108 Plastic Coated Hand Rim Upgrade (Not Avail. w/16" Wheels)	\$92.00	
				<input type="checkbox"/>	Omit Hand Rims	N/C

Rear Anti-Tips

<input checked="" type="checkbox"/>	E0971	Rear Anti-Tips with Wheels	\$135.00
-------------------------------------	-------	----------------------------	----------

22"/24" Wheel Axle Options

<input type="checkbox"/>	Adj. Position Axle System	\$122.00
--------------------------	---------------------------	----------

Attendant Wheel Lock Option

Requires 16" rear wheels and 8" casters

<input type="checkbox"/>	Attendant Wheel Lock	\$100.00
--------------------------	----------------------	----------

Caster Options

Seat pan to floor height is 15.5" with 5" caster, and 17.5" with 8" caster.

<input type="checkbox"/>	5" x 1.5" Poly	N/C	<input type="checkbox"/>	8" x 2" Poly	N/C
--------------------------	----------------	-----	--------------------------	--------------	-----

Lateral Pelvic/Thigh Supports

<input type="checkbox"/>	E0956	Left Pelvic/Thigh Support Pkg.	\$120.00	<input type="checkbox"/>	E0956	Right Pelvic/Thigh Support Pkg.	\$120.00
--------------------------	-------	--------------------------------	----------	--------------------------	-------	---------------------------------	----------

Lifeform™ Anterior Pelvic Positioning Belts

Medial 4-point pelvic positioning belt recommended with REACTIONi seating. Positioning belts not for use as a occupant restraint during transit in a motor vehicle.

<input type="checkbox"/>	E0978	LifeFormi 2-point 1" Metal PB Buckle, 3.5" pads	\$88.00	<input type="checkbox"/>	E0978	LifeFormi 2-point 1.5" Metal PB Buckle, 6.5" pads	\$90.00
<input type="checkbox"/>	E0978	LifeFormi Medial 4-point 1" Metal PB Buckle, 3.5" pads	\$110.00	<input type="checkbox"/>	E0978	LifeFormi Medial 4-point 1.5" Metal PB Buckle, 6.5" pads	\$110.00
<input type="checkbox"/>	E0978	LifeFormi Medial 4-point 1" Metal PB Buckle, 4.5" pads	\$110.00				

Lifeform™ Anterior Trunk Support

Anterior Trunk Supports not for use as a occupant restraint during transit in a motor vehicle.

LifeForm™ Contoured Dynamic				LifeForm™ Zipper Opening Dynamic			
<input type="checkbox"/>	E0960	Contoured - Early Intervention	\$150.00	<input type="checkbox"/>	E0960	Zipper Opening - S	\$150.00
<input type="checkbox"/>	E0960	Contoured - XS	\$150.00	<input type="checkbox"/>	E0960	Zipper Opening - M	\$150.00
<input type="checkbox"/>	E0960	Contoured - S	\$150.00	<input type="checkbox"/>	E0960	Zipper Opening - L	\$150.00
LifeForm™ Contoured Structured				LifeForm™ Zipper Opening Structured			
<input type="checkbox"/>	E0960	Contoured Structured - S	\$150.00	<input type="checkbox"/>	E0960	Structured Zipper Opening - S	\$150.00
				<input type="checkbox"/>	E0960	Structured Zipper Opening - M	\$150.00
				<input type="checkbox"/>	E0960	Structured Zipper Opening - L	\$150.00

Lifeform™ Swing-Away Abductor Hardware and Pads

4-way stretch cover included with Lifeform Abductor pads.

<input type="checkbox"/>	E1028	ROCK 2 S/A Abduct. Hardware	\$220.00	<input type="checkbox"/>	E0957	Small Wedge(2"-3"W x 4"L x 3"D)	\$140.00
				<input type="checkbox"/>	E0957	Med. Wedge (2"-4.5" x 5"L x 3"D)	\$140.00
				<input type="checkbox"/>	E0957	Lg. Wedge (3"-5" x 5"L x 4"D)	\$140.00

Lifeform™ Headrest and Support Hardware

<input type="checkbox"/>	E0955	8" Curved Oval Pad	\$195.00	<input type="checkbox"/>	E1028	LifeFormi Removable Headrest Mounting Hardware (Required for all Headrest Options)	\$215.00
<input type="checkbox"/>	E0955	10" Curved Oval Pad	\$195.00				
<input type="checkbox"/>	E0955	14" Curved Oval Pad	\$195.00	<input type="checkbox"/>		Universal Head Rest Mount Plate Only*	\$50.00
<input type="checkbox"/>	E0955	Narrow Adjustable Arc Pad	\$280.00				
<input type="checkbox"/>	E0955	Standard Adjustable Arc Pad	\$280.00				

*Accommodates several other manufacturers headrests. Call for details.

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

Lifeform™ Ankle and Foot Positioning

Stealth Hard Side Foot Stabilizers				LifeForm™ Ankle Hugger Foot Stabilizers				
<input type="checkbox"/>	K0108	Small, Pair, Left and Right	\$166.00	<input type="checkbox"/>	3610-5230L	K0108	Left - Small	\$72.00
<input type="checkbox"/>	K0108	Medium, Pair, Left and Right	\$166.00	<input type="checkbox"/>	3610-5230R	K0108	Right - Small	\$72.00
<input type="checkbox"/>	K0108	Large, Pair, Left and Right	\$166.00	<input type="checkbox"/>	3610-5240L	K0108	Left - Medium	\$72.00
				<input type="checkbox"/>	3610-5240R	K0108	Right - Medium	\$72.00
				<input type="checkbox"/>	3610-5250L	K0108	Left - Large	\$72.00
				<input type="checkbox"/>	3610-5250R	K0108	Right - Large	\$72.00

Lifeform™ Swing Away Lateral Support Hardware

Choose your bracket and CURVED pad size				Choose your bracket and PLANAR pad size			
<input type="checkbox"/>	E0956	Left 1.5" pediatric bracket lateral thoracic support package, Curved Pad	\$110.00	<input type="checkbox"/>	E0956	Left 1.5" pediatric bracket lateral thoracic support package, Planar Pad	\$110.00
<input type="checkbox"/>	E0956	Right 1.5" pediatric bracket lateral thoracic support package, Curved Pad	\$110.00	<input type="checkbox"/>	E0956	Right 1.5" pediatric bracket lateral thoracic support package, Planar Pad	\$110.00
<input type="checkbox"/>	E0956	Left 2" bracket lateral thoracic support package, Curved Pad	\$135.00	<input type="checkbox"/>	E0956	Left 2" bracket lateral thoracic support package, Planar Pad	\$135.00
<input type="checkbox"/>	E0956	Right 2" bracket lateral thoracic support package, Curved Pad	\$135.00	<input type="checkbox"/>	E0956	Right 2" bracket lateral thoracic support package, Planar Pad	\$135.00
Choose a pad size: <input type="checkbox"/> 3"x4" <input type="checkbox"/> 4"x4" <input type="checkbox"/> 4"x5" <input type="checkbox"/> 4"x6"				Choose a pad size: <input type="checkbox"/> 3"x4" <input type="checkbox"/> 4"x4" <input type="checkbox"/> 4"x5" <input type="checkbox"/> 4"x6"			
Choose Swing-Away Hardware (required)				Choose Swing-Away Hardware (required)			
<i>Requires 3610-3110, 3120, 3210, 3220.</i>				<i>Requires 3610-3110, 3120, 3210, 3220.</i>			
<input type="checkbox"/>	E1028	Left 1.5" DoubleSwingi pediatric swing-away addition	\$175.00	<input type="checkbox"/>	E1028	Left 1.5" DoubleSwingi pediatric swing-away addition	\$175.00
<input type="checkbox"/>	E1028	Right 1.5" DoubleSwingi pediatric swing-away addition	\$175.00	<input type="checkbox"/>	E1028	Right 1.5" DoubleSwingi pediatric swing-away addition	\$175.00
<input type="checkbox"/>	E1028	Left 2" DoubleSwingi swing-away bracket	\$175.00	<input type="checkbox"/>	E1028	Left 2" DoubleSwingi swing-away bracket	\$175.00
<input type="checkbox"/>	E1028	Right 2" DoubleSwingi swing-away bracket	\$175.00	<input type="checkbox"/>	E1028	Right 2" DoubleSwingi swing-away bracket	\$175.00

LifeForm™ Swing-Away Width and Depth Extensions (optional)		
Options Available only with Swing-Away Hardware.		
<input type="checkbox"/>	Depth Ext. + 1", pediatric, Left	N/C
<input type="checkbox"/>	Depth Ext. + 1", pediatric, Right	N/C
<input type="checkbox"/>	Width Ext. + 1", pediatric, Left	N/C
<input type="checkbox"/>	Width Ext. + 1", pediatric, Right	N/C
<input type="checkbox"/>	Depth Ext. + 1", 2" Brkt., Left	N/C
<input type="checkbox"/>	Depth Ext. + 1", 2" Brkt., Right	N/C
<input type="checkbox"/>	Width Ext. + 1", 2" Brkt., Left	N/C
<input type="checkbox"/>	Width Ext. + 1", 2" Brkt., Right	N/C

LifeForm™ Fixed Lateral Support Hardware

Choose your bracket and CURVED pad size				Choose your bracket and PLANAR pad size			
<input type="checkbox"/>	E0956	Left 1.5" fixed pediatric bracket lateral thoracic support package, Curved Pad	\$110.00	<input type="checkbox"/>	E0956	Left 1.5" fixed pediatric bracket lateral thoracic support package, Planar Pad	\$110.00
<input type="checkbox"/>	E0956	Right 1.5" fixed pediatric bracket lateral thoracic support package, Curved Pad	\$110.00	<input type="checkbox"/>	E0956	Right 1.5" fixed pediatric bracket lateral thoracic support package, Planar Pad	\$110.00
<input type="checkbox"/>	E0956	Left 2" fixed bracket lateral thoracic support package, Curved Pad	\$135.00	<input type="checkbox"/>	E0956	Left 2" fixed bracket lateral thoracic support package, Planar Pad	\$135.00
<input type="checkbox"/>	E0956	Right 2" fixed bracket lateral thoracic support package, Curved Pad	\$135.00	<input type="checkbox"/>	E0956	Right 2" fixed bracket lateral thoracic support package, Planar Pad	\$135.00
Choose a pad size: <input type="checkbox"/> 3"x4" <input type="checkbox"/> 4"x4" <input type="checkbox"/> 4"x5" <input type="checkbox"/> 4"x6"				Choose a pad size: <input type="checkbox"/> 3"x4" <input type="checkbox"/> 4"x4" <input type="checkbox"/> 4"x5" <input type="checkbox"/> 4"x6"			

Accessories

<input type="checkbox"/>	ROCK Sunshade	\$265.00	<input type="checkbox"/>	Underchair Cargo Tray (carries up to 10 lbs)	\$150.00
Replacement Seat Cushion Covers					
<input type="checkbox"/>	E2610	Replacement Cover 8"-11" Seat Cushion	\$102.00		
<input type="checkbox"/>	E2610	Replacement Cover 11"-14" Seat Cushion	\$102.00		
<input type="checkbox"/>	E2610	Replacement Cover 14"-17" Seat Cushion	\$102.00		

Additional Support Accessories

Additional Support Accessories NOT tested with Transit. Must remove prior to tie down.							
<input type="checkbox"/>	K0105	IV Support Pole	\$150.00	<input type="checkbox"/>	E1029	Battery/Vent. Accessory Tray	\$295.00
<input type="checkbox"/>	E2208	Oxygen Tank Holder	\$150.00	Requires 16" Rear Wheels and 8" casters. (Does not include vent. mount for LTV)			
<input type="checkbox"/>	E1029	Ventilator Tray for LTV Ventilator Series Only Requires 16" Rear Wheels and 8" Casters	\$370.00				

Custom Notes:

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.