

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643
Sales: Phone: (866) 800-2002 • Fax: (866) 707-3422

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:	P.O. #:	Discount:	
Address:	Address:	Terms:		
Phone Contact:	Marked For:	Freight:		
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.			All pricing is MSRP.	Taken By:

☐ QUANTUM 1107 2SP-SS (single power) REFERENCE GUIDE #1

☐ Quantum 1107 2SP-SS Group 2 / Single Power **K0835** \$6,321
200 lbs. weight capacity

Seat Options:

☐ TB Tilt E1002
☐ TB Recline E1004

Options:

- Flip-back, cane mounted, full length arms
- 70 amp, 6-Key PG VSI controller
- 10" Solid tires
- 5 amp off-board charger
- 60" lapbelt
- Swing-away legrests
- Reinforced footplates
- Tilt through 6-key VSI
- ☐ Omit Solid Tires (-\$213)
- ☐ Omit 6-Key VSI (-\$200)
- ☐ Omit contoured seat (-\$275)

Select seating , front rigging options in appropriate section & base color.

☐ QUANTUM 1107 2SP-C (single power) REFERENCE GUIDE #2

☐ Quantum 1107 2SP-C Group 2 / Single Power **K0836** \$6,321
200 lbs. weight capacity

***The following single power options are available at a later date:**

1. Power Tilt
2. Power Recline

Options:

- 18" W x 18" D Ltd. Recline Gray vinyl high-back
- Full length armrests
- 70 amp 6-key PG VSI
- 10" Solid tires
- 5 amp off-board charger
- 60" lapbelt
- Foot Platform
- ☐ Omit Solid Tires (-\$213)
- ☐ Omit 6-Key VSI (-\$200)
- ☐ Omit contoured seat (-\$275)

Select seating , front rigging options in appropriate section & base color.

EXPANDABLE CONTROLS

All joysticks must be used in conjunction w/ either a PG or Q-Logic controller

CONTROLLER OPTIONS

- ☐ Q-Logic Controller E2377 \$695
- ☐ PG Controller E2377 \$695

All joysticks or enhanced display must be used in conjunction w/ either a PG or Q-Logic controller

JOYSTICK OPTIONS

Q-Logic Electronics

- ☐ 70 amp Q-Logic Hand Control^{3,5,29} E2399 \$700
- ☐ 70 amp Q-Logic Stand Alone Control^{3,29,30} E2399 \$700

- ☐ Q-logic Enhanced Display²⁷ \$1720

All joysticks or enhanced display must be used in conjunction w/ either a PG or Q-Logic controller

PG Electronics

- ☐ 3-Key Remote Plus 70 amp^{1,3,22} E2399 \$700
(Select Box for Profile)

- ☐ 1 Profile ☐ Multiple Profile

- ☐ 7-Key Access. Remote Plus 70 amp^{1,3,4,29} E2399 \$700

- ☐ Remote Plus 70 amp with Speed Dial^{3,5,29} E2399 \$700

- ☐ Remote Plus 70 amp with Toggles^{3,5,29} E2399 \$700

- ☐ Remote Plus 70 amp with Switch Jacks^{3,5,29} E2399 \$700

- ☐ Microdrive, 70 amp^{3,5,6,29} E2399 \$700
- ☐ Omni Plus²⁷ \$1721.63

(Requires a secondary control to operate the power chair)

Joystick includes electronics, cables, & junction boxes.

When configuring power chair use reference guides to determine accessories included in the MSRP pricing package

Model number with suffix must be used for all chairs billed through Medicare.

SPECIFIC CONFIGURATION

- ☐ Quantum 1107 \$5,908
- ☐ Quantum 1107 2SP-SS \$6,321
- ☐ Quantum 1107 2SP-C \$6,321
- ☐ Omit contoured seat²² (-\$275)

300 lbs. weight capacity w/o power positioning systems

COLOR

- ☐ Candy Apple Red NC
- ☐ Viper Blue NC
- ☐ Onyx Black NC
- ☐ Sunburst Orange NC

BATTERIES (2 required)

- ☐ U-1, 32 amp E2365 \$133ea

NON-EXPANDABLE CONTROLS

- ☐ Pride Flight 75 amp¹ NC
- ☐ VSI Integral 70 amp^{1,26} NC
- ☐ 6-Key VSI Integral 70 amp^{1,3} \$200

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007

PROGRAMMER

<input type="checkbox"/> Flight Programmer	\$423
<input type="checkbox"/> PG Programmer	\$423
<input type="checkbox"/> Q-Logic Programmer	\$423

CONTROLLER MOUNT

<input type="checkbox"/> Fixed, In-Line Joystick Mount	Std.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Quick Disconnect, Joystick Mount ^{24,32}	NC
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Swing-away Joystick Mount ²⁵	E1028 \$196
<input type="checkbox"/> 4" L ² <input type="checkbox"/> 7" L	
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Multi-axis Joystick Mount	K0108 \$100
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<i>(Fixed, does not swing-away)</i>	
<input type="checkbox"/> Vertical Adj. Clamp ⁸	K0108 \$95
Provides vertical joystick adjustment w/o adjusting armrest height.	

JOYSTICK HANDLE OPTIONS

Q-Logic & PG Controls

<input type="checkbox"/> Std. Knob		Std.
<input type="checkbox"/> Mushroom Handle	E2323	\$71
<input type="checkbox"/> T-Bar	E2323	\$107
<input type="checkbox"/> Chin Cup	E2323	\$71
<input type="checkbox"/> Extended I-Handle	E2323	\$71
<input type="checkbox"/> Soft Ball Handle	E2323	\$50
<input type="checkbox"/> Medium Ball Handle	E2323	\$50
<input type="checkbox"/> Goal Post Handle	E2323	\$66
<input type="checkbox"/> 3" W <input type="checkbox"/> 4" W		

Flight Controls

<input type="checkbox"/> Flight Knob		Std.
<input type="checkbox"/> Sponge Ball	E2323	\$35
<input type="checkbox"/> Chin Cup	E2323	\$30
<input type="checkbox"/> DL Knob	E2323	\$20
<input type="checkbox"/> Carrot Knob	E2323	\$20
<input type="checkbox"/> Straight Bar (Chrome)	E2323	\$45
<input type="checkbox"/> Finger Tube (Chrome)	E2323	\$45
<input type="checkbox"/> Round T-Bar (Chrome)	E2323	\$45
<input type="checkbox"/> Flat T-Bar (Chrome)	E2323	\$45
<input type="checkbox"/> Wrist Plate	E2323	\$45

DRIVE WHEELS

<input type="checkbox"/> 10" Pneumatic (pair) ²⁶	NC
<input type="checkbox"/> 10" Flat-Free (pair)	\$213

Quantum[®] 1107 Series

Order Form / Price List
Revised: November 2007

Legrest measurements have changed. Lengths are now measured according to industry standard from the seat pan to the back of the footplate

NEW SWING-AWAY LEGRESTS

<input type="checkbox"/> 70° Swing-away footrests Style #8 ³⁴	\$162
<input type="checkbox"/> 70° Tapered S/A footrests Style #8 ³⁴	\$162
<i>(N/C with power positioning systems)</i>	
<input type="checkbox"/> Quad Release Lever	NC
S/A Standard Extensions:	
<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
<i>(Lengths over 17" may cause interference)</i>	
<input type="checkbox"/> Aluminum Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 8.5 x 6" ⁹ (Tapered legs require 18"W or greater seat)	

Angle Adjustable Ext.:	\$150
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<i>(Lengths over 17" may cause interference)</i>	

Angle Adjustable Footplates (pair)

<input type="checkbox"/> 4"x6" ³³ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

Multi-Axis Ext. (ball & socket):	\$495
<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	

(Lengths over 17" may cause interference)

Multi-Axis Footplates (pair)

<input type="checkbox"/> 4"x6" ³³ (Tapered legs require 14"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 22"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

FRONT RIGGINGS

<input type="checkbox"/> Elevating Leg Rests E0990	\$325
Adj. Calf Pads w/ ELRS (1 hole std.)	
<input type="checkbox"/> 1 hole <input type="checkbox"/> 3 holes* <input type="checkbox"/> 6 holes*	
*Potential for interference w/ shrouds	

<input type="checkbox"/> HD Drop-Ins, 70° ^{7,8,15}	\$420
<i>N/A w/ contoured or Specialty seating</i>	
Angle Extensions:	
<input type="checkbox"/> 12" - 14.5" Ext.	NC
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	
<i>(Lengths over 17" may cause interference)</i>	

<input type="checkbox"/> Swing-away, Pediatric	\$162
<input type="checkbox"/> 4"x6" <input type="checkbox"/> 5"x6" Angle-Adj. (pair)	
<input type="checkbox"/> Drop-ins, Pediatric	\$162
<input type="checkbox"/> 4"x6" <input type="checkbox"/> 5"x6" Angle-Adj. (pair)	

FRONT RIGGINGS CONT.

<input type="checkbox"/> Foot Platform, 9.5" deep ^{20,31}	NC
<input type="checkbox"/> Platform Bracket Small	NC
<input type="checkbox"/> Platform Bracket Large ²¹	NC
<input type="checkbox"/> Foot Platform, 11.5" deep ^{20,31}	NC
<input type="checkbox"/> Platform Bracket Small	NC
<input type="checkbox"/> Platform Bracket Large ²¹	NC
<input type="checkbox"/> No Front Rigging	NC
Power-elevating, high-pivot^{16,17} E1010 \$1695	
Legrests, (calf pads included) ¹⁹	
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Power-elevating, Articulating ^{13,16,17} E1010 \$2395	
Articulating Range: <input type="checkbox"/> 4" L Std. <input type="checkbox"/> 6" L NC	
Set at: <input type="checkbox"/> 70° <input type="checkbox"/> 80°	
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Swing-away Power-elevating, E1010	
Articulating ^{13,16,17,28} \$2395	
Articulating Range: <input type="checkbox"/> 4" L Std.	
Set at: <input type="checkbox"/> 70° <input type="checkbox"/> 80°	
<input type="checkbox"/> Quad Release Lever	NC
<input type="checkbox"/> Independent <input type="checkbox"/> Combined ¹⁶ <input type="checkbox"/> Synchronized ¹⁸	
<input type="checkbox"/> Calf Pads Left & Right (6 x 5.35)	Std.
<input type="checkbox"/> Calf Pads Left & Right (8.00 x 5.86)	NC
Footplates (pair)	
<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁰ <input type="checkbox"/> 8.5 x 6" ⁹	
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ¹⁰	\$150
<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ⁹	\$150
<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ¹¹	\$150
<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$25
Extensions:¹²	
<input type="checkbox"/> 11"-14" Ext. ²⁸ (\$125 w/ power legs)	NC
<input type="checkbox"/> 14"-17" Ext	Std.
<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	\$326
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount Clamp-on; Multi-axis	\$412
<input type="checkbox"/> Small <input type="checkbox"/> Large	
<input type="checkbox"/> High mount footplatform ¹⁹ K0108 \$450	
<input type="checkbox"/> Mid-wheel 6 foot platform	
Lower Extensions	
<input type="checkbox"/> 10"-14" L	
<input type="checkbox"/> 14"-17" L	

Rigging Options

<input type="checkbox"/> Heel Loops (pair) ¹³ E0951	\$20
<input type="checkbox"/> Gel Neoprene Footrest K0108	
Sleeve (pair) ¹⁴ \$65	
<input type="checkbox"/> ELR / ALR Adductor Buttons (pair) K0108	\$100
<input type="checkbox"/> POSALinc Knee Adductor (pair)	
for S/A legrests K0108	\$464
<input type="checkbox"/> Residual limb support left E1020	\$300
<input type="checkbox"/> Residual limb support right E1020	\$300

Footnotes for pages 1, 2 & 3:

- Std. w/ one programmable drive profile.
- N/A with waterfall armpads and dual post armrests.
- Needed for controlling power positioning systems.
- Needed for lighting packages.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- Std. w/ heel loops & ang-adj. footplates, check size below.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Customer Service.
- Requires the addition of an ALM.
- Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
- Only avail. on TB Recline.
- Std. calf pad is not depth adjustable.
- No charge when substituting for S/A legrests.
- Allows for the depth adjustability range to be offset 1".
- Applicable when ordering no seat, specialty, Synergy seat or any power positioning system.
- Please note swing away footrests are standard when ordered with Power Positioning.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- If ordered w/ reference guides or configuration packages discounts may apply.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- 11" - 14" Ext. are N/A w/ power S/A legrests.
- Separately reimbursable through Medicare for multiple power units.
- Std. w/ three drive profiles.
- Fixed bracket, does not elevate.
- N/A w/ 2-post, desk length arms w/ a toggle.
- Only avail. on 12" W or greater seats.
- Please note swing away footrests are standard when ordered with Power Positioning.

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007

STANDARD SEATING

- Medium-Back** NC
- Medium-Back w/ Solid Seat Pan^{1,2}** NC
- Brown Vinyl
- Beige Vinyl
- Gray Vinyl
- Black Vinyl
- Gray Fabric
- Brown Fabric

Medium-Back Seating Dimensions (Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.

UPGRADED SEATING

- Ltd. Recline High-Back w/Headrest** NC
- Ltd. Recline High-Back with Headrest^{1,2} w/ Solid Seat Pan** NC
- Black Vinyl
- Gray Vinyl
- Beige Vinyl
- Gray Fabric
- Brown Fabric

High-Back Seating Dimensions (Dimensions 18"W x 18"D)

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 20"	NC

ARMRESTS

For Medium and High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable removable⁴ \$170 E0973

Armpads

- Desk Length, 10" NC
- Left NC
- Right NC
- Full Length, 14" NC
- Left Std.
- Right Std.

MANUAL RECLINER ^{2,5}

- Manual Recline High Back** \$600
- Manual Recline High Back¹ w/ Solid Seat Pan** E1226 \$600
- (180 degrees of Recline)^{2,6}**
- (Dimensions 18"W x 18"D)*

- Gray Vinyl
 - Beige Vinyl
- | | | | |
|------------------------------|------|------------------------------|------|
| Width: | | Depth: | |
| <input type="checkbox"/> 16" | NC | <input type="checkbox"/> 16" | NC |
| <input type="checkbox"/> 18" | Std. | <input type="checkbox"/> 18" | Std. |
| <input type="checkbox"/> 20" | NC | <input type="checkbox"/> 20" | NC |

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS⁶

- Simplicity Cushion**
- 16 x 16 - 20 x 20 E2601 \$98
- Solution 1 Cushion**
- 14 x 14 - 20 x 18 E2603 \$380
- 20 x 20 E2603 \$430
- Solution Cushion**
- 10 x 10 - 20 x 18 E2607 \$380
- 20 x 20 E2607 \$430
- Spectrum Cushion**
- Gel
- 10 x 10 - 20 x 18 E2605 \$480
- 20 x 20 E2605 \$530
- Visco
- 10 x 10 - 20 x 18 E2607 \$480
- 20 x 20 - 20 x 20 E2607 \$530
- Structure 2 Cushion**
- 14 x 16 - 20 x 18 K0734 \$480
- 20 x 20 K0734 \$580

SPECIALTY SEATING

- Specialty Seat Kit (limited adjustability)** NC

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" NC
- 18" Std.
- 20" NC

Specialty Seating Kit Depth

- 16" Std.
- 17" NC
- 18" Std.
- 19" NC
- 20" NC

Specialty Seating Kit Back Cane Height

- 18" Std.
- 20" NC
- 22" NC

Back Selection

- See Synergy Shape Back Order Form^{11,12,13}**
- Solid Plate Back up to 20" W ¹⁷ \$275
- See Positioning Components Order Form for additional backs

Back Height Pan ¹⁴

- 18" Std.
- 19" NC
- 20" NC
- 21" NC
- 22" NC
- 23" NC
- 24" NC
- 25" NC

Sling Back

- 16" wide \$98
- 18" wide \$98
- 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
- 20" height NC
- 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

- 2-Post, fixed flip-back^{13,16} Std.
- Flip-back, cane mounted^{13,16} \$95ea
- Removable, adjustable, quick height^{10,18} E0973 \$275ea
- 11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments*
- Removable, height adjustable¹⁸
- heavy duty armrests E0973 \$150ea
- 8"H-12"H NC
- 10"H-14"H NC

Armrest Pads

- Full Length, straight, 2"x14" NC
- Left NC
- Right NC
- Desk Length, straight, 2"x10" NC
- Left NC
- Right NC
- Full Length, Waterfall, 2"x14"¹⁹ NC
- Left NC
- Right NC
- Desk Length, Waterfall, 2"x10"^{18,9} NC
- Left NC
- Right NC

HEADRESTS ⁷

- Stealth ^{3,15} E0955**
- Large (11" wide x 6" high) \$210
- Stealth ^{3,15} E0955**
- Small (9 3/4" wide x 3" high) \$210
- Stealth, swing-away (Flip-Down)**
- Head Rest Bracket E1028 \$150
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications ^{3,15} E0955**
- Large (11" wide x 6" high) \$275
- Medical Applications ^{3,15} E0955**
- Small (9 3/4" wide x 3" high) \$275
- Fixed, Offset Mounting Bracket NC

Footnotes for pages 4:

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special measurements are not available.
3. Comes with swing-away / removable hardware; HCPCS code is E1028.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.
6. See Synergy Cushion order forms for detailed pricing, sizes & codes.
7. Only available w/ solid backs.
8. S/A joystick mount N/A.
9. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adj. range is 11" - 17".
11. 14" & 22" width are N/A.
12. Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
13. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
14. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
15. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
16. S/A laterals interfere w/ flip-back arms.
17. 20" D is N/A.
18. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007

Synergy Seating Systems¹¹

- Synergy Seating™ System \$434
(Price with Base)
- Synergy Manual Recline Pkg.^{12,13,14,15,21} \$700
Maximum weight 300 lbs. E1226
- Synergy Manual Tilt Pkg.^{13,16} \$550
(Max Dimensions 18"W x 18"D) K0108
Maximum weight 160 lbs.

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Width Package

- 10"W-14"W
- 14"W-18"W
- 16"W-20"W

Seat Width Package Manual Tilt

- 10"W-12"W
- 12"W-14"W
- 14"W-16"W
- 16"W-18"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 10" NC	<input type="checkbox"/> 11" NC
<input type="checkbox"/> 12" NC	<input type="checkbox"/> 13" NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Depth Package

- 10"D-14"D
- 14"D-18"D
- 16"D-20"D

Seat Depth Package Manual Tilt

- 10"D-14"D
- 14"D-16"D
- 16"D-18"D

BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

BACK OPTIONS

- See Synergy Shape Back Order Form^{1,2,3}
- Solid Plate Back up to 20" W²¹ \$275
- Back Height:**⁴
 - 18"H
 - 19"H
 - 20"H
 - 21"H
 - 22"H
 - 23"H
 - 24"H
 - 25"H
- Sling Back^{4,18} \$98

Back height will be the same as back cane height. Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back¹²

- Back Height:**
 - 22" 25"
 - Back up to 20" W E2620 \$595
 - 22" W E2621 \$695
- TRU-Comfort Cushion
 - 16 x 16 - 20 x 18 E2601 \$380
 - 20 x 20 - 20 x 22 E2601 \$380
 - 22 x 16 - 22 x 24 E2602 \$480

See TRU-Comfort Back & Cushion order forms for specific sizes

ARMRESTS

- Flip-back, cane mounted^{3,5,18} Std.
- 2-Post, flip-back, height^{3,5,19} E0973 \$150ea
adjustable, removable armrests
(9"-13" of height adjustment in 1" increments)
- Single Post Height adj.^{6,20} E0973 \$75ea
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments
- Quad Release Lever NC

- Removable, height adjustable E0973 \$75ea
heavy duty armrests²⁰

- 8"H-12"H 10"H-14"H
- Pediatric removable²⁰ E0973 \$75ea
adjustable, quick height.

- Match Shroud Color Black
- Quad Release Lever NC

(7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)

ARMPADS

- Full length, Straight, 2" X 14" Std.
 - Left Right
- Desk Length, Straight, 2" X 10" Std.
 - Left Right
- Full Length, Waterfall, 2" X 14"⁷ NC
 - Left Right
- Desk Length, Waterfall, 2" X 10"^{7,17} NC
 - Left Right

HEAD REST⁸

- Stealth^{9,10} E0955
 - Large (11" wide x 6" high) \$210
 - Small (9 3/4" wide x 3" high) \$210
- Stealth, swing-away (Flip-Down)
 - Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

- Medical Applications^{9,10} E0955
 - Large (11" wide x 6" high) \$275

- Medical Applications^{9,10} E0955
 - Small (9 3/4" wide x 3" high) \$275

- Fixed, Offset Mounting Bracket NC

- TRU-Comfort^{9,10} E0955
 - 8" Wide w/ Cobra hardware \$417
 - 10" Wide w/ Cobra hardware \$417
 - Bracket for after market headrest \$180

Footnotes for pages 5:

1. 14" & 22" width are N/A.
2. Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
3. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
4. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
5. S/A laterals interfere w/ flip-back arms.
6. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
7. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
8. Available w/ solid backs only.
9. Comes with swing-away / removable hardware; HPCPS code is E1028.
10. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
11. Manual systems do not include a headrest.
12. N/A w/ flip-back, cane mounted arms.
13. Solid back not included.
14. Includes removable, adjustable, quick-height armrests.
15. N/A w/ Shape Back.
16. Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
17. N/A w/ S/A joystick mount.
18. N/A w/ Synergy manual systems.
19. May cause interference with certain dimension configurations. Call Quantum Customer Service.
20. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
21. Solid plate is N/A w/manual recline.

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007

RECOMMENDED INFORMATION

Height: _____
Weight: _____
Leg Length: _____
Lower Leg Length: _____

TRU-BALANCE 2 TILT

TRU-Balance 2 Power Tilt **E1002** \$4995
55 degree center of gravity (CG) power
tilt seating system. 300lbs weight capacity.

TRU-BALANCE 2 RECLINE

TRU-Balance Power Recline **E1004** \$5595
Full Drive lock-out at 20 degrees
300 lbs. weight capacity

SYSTEM SPECIFICATIONS**Tower Selection (Tilt or Tilt/Recline only)**

- Adjustable-angle Std.
Will raise minimum STF by 7/8".
 Fixed NC

Set Seat Width At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | |

Seat Width Package

- 14"W-18"W 16"W-20"W

Set Seat Depth At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | |

Seat Depth Package

- 14"D-18"D 16"D-20"D

Push Bar selection must match back width;
center tube standard when selecting flip-up arms.

Set Back Cane Height At:

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" Std. | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" NC | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

Seat Overall Back Height³

- | | |
|-----------------------------------|---------------------------------|
| <input type="checkbox"/> 14" NC | <input type="checkbox"/> 15" NC |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 17" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 19" NC |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 21" NC |
| <input type="checkbox"/> 22" NC | <input type="checkbox"/> 23" NC |
| <input type="checkbox"/> 24" NC | <input type="checkbox"/> 25" NC |

BACKS

- See Synergy Shape Back Order Form⁶
 Solid Plate Back⁴ Std.
 Attendant manual recline^{1,4,11,23,24} \$750
180 degrees of recline with 2 1/2"
raised pivot point. Controlled through
dual gas struts. 300 lb. weight capacity.

(For use on TB Tilt Systems only)

- TRU-Balance Reclining Back Std.
(For use on Recline Systems only)

 TRU-Comfort Back¹³

- Back Height:** 22" 25"
 Back up to 20" W **E2620** \$595
 22" W **E2621** \$695
 TRU-Comfort Cushion
16 x 16 - 20 x 18 **E2601** \$380
20 x 20 - 20 x 22 **E2601** \$380
22 x 16 - 22 x 24 **E2602** \$480
 Omit Back -(\$100)

See TRU-Comfort Back & Cushion
order forms for specific sizes

ARMRESTS

- Flip-back, cane mounted^{5,9,14} Std.
Mounted on back canes with full width spreader bar.

- 2-Post, flip-back, hgt. adj, removable^{9,15} NC
(9"-13" of height adjustment in 1" increments)
 Single Post Height adj.^{15,16} NC
*11"-17" of height adj. 3" on top in 1/2"
increments & 3" on bottom in 1" increments*
 Quad Release Lever NC
(11"-17" of height adj in 3/8" increments)

- Removable, single-post,¹⁶
height adjustable, heavy duty NC
 8"-12" Adjustment
 10"-14" Adjustment

- Pediatric, removable,¹⁶ NC
adjustable, quick height
*(7 3/4" - 9 3/4" height adjustment in 1/2"
increments from seat pan to top of armrest)*
 Quad Release Lever NC
 Match Shroud Kid Color Black

- Reclining, Flip-back⁹ Std.
 8"-9" H 10"-11" H
(For use on Recline Systems only)

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
 Desk Length, straight, 2"x10" NC
 Left Right

- Full Length, Waterfall, 2"x14"¹⁷ NC
 Left Right
 Desk Length, Waterfall, 2"x10"^{2,17} NC
 Left Right

HEAD REST

- Stealth^{18,19} **E0955**
Large (11" wide x 6" high) \$210
 Stealth^{18,19} **E0955**
Small (9 3/4" wide x 3" high) \$210

- Stealth, swing-away (Flip-Down)
Head Rest Bracket **E1028** \$150
Interfaces w/ any headrest that uses 1/2" square down-tube.

- Medical Applications^{18,19} **E0955**
Large (11" wide x 6" high) \$275

- Medical Applications^{18,19} **E0955**
Small (9 3/4" wide x 3" high) \$275

- Omit Head Rest - To be supplied by provider

- Fixed, Offset Mounting Bracket NC

- Reclining Head Rest (6.5" x 7")^{18,19} **E0955** \$180
(For use on Recline Systems only)

- TRU-Comfort^{18,19} **E0955**
 8" Wide w/ Cobra hardware \$417
 10" Wide w/ Cobra hardware \$417

- Bracket for after market headrest \$180

POWER POSITIONING ELECTRONICS**Select Electronics to be used**

- Q-Logic
 PG Remote Plus
 PG 4-Key VSI
 PG 6-Key VSI

SINGLE POWER ELECTRONICS**Power Tilt or Power Recline**

- Through Toggle Box²⁰ (Tilt only) NC
 Left Right (mounts to side of armrest)

- Through 6-Key VSI (Tilt only) NC

- Through Single Switch (1/8" Jack)
Tilt Module^{21,22} (Tilt only) \$500

- Through 2-Directional Control Box
w/ dual Switch Jacks (1/8" Jack)^{21,22} \$1175

- Through Q-Logic Controller **E2310** \$1175

- Through Remote Plus Controller **E2310** \$1175

Multiple Motor/Actuator Electronics*:

* At least one TB2 seat function
plus power riggings

- Operation Q-Logic Controller **E2311** \$2375

- Operation Remote Plus Controller **E2311** \$2375

- Operation through toggle \$2375
(Select one)

- 4-way Toggle
 4-way Button Toggle
 Quad Push-Button Control

Toggle box mounts to side of armrest.

- Left Right

Quantum® 1107 Series

Order Form / Price List
Revised: November 2007

POSITIONING COMPONENTS

AES Laterals

<input type="checkbox"/> Thoracic Curved-S ^{18,19}	\$532pr
<input type="checkbox"/> Thoracic Curved-M ^{18,19}	\$532pr
<input type="checkbox"/> Thoracic Curved-L ^{18,19}	\$532pr
<input type="checkbox"/> Thoracic Tapered-S ^{18,19}	\$532pr
<input type="checkbox"/> Thoracic Tapered-M ^{18,19}	\$532pr
<input type="checkbox"/> Thoracic Tapered-L ^{18,19}	\$532pr
<input type="checkbox"/> Hip & Thigh Removable Short ^{18,19}	\$474pr
<input type="checkbox"/> Hip & Thigh Removable Long ^{18,19}	\$474pr
<input type="checkbox"/> Thoracic S/A Curved Rail-S ^{18,19,25}	\$532pr
<input type="checkbox"/> Thoracic S/A Curved Rail-M ^{18,19,25}	\$532pr
<input type="checkbox"/> Thoracic S/A Curved Rail-L ^{18,19,25}	\$532pr

Stealth

<input type="checkbox"/> Lateral Thoracic Curved-S/A-Pediatric	\$650pr
Flush to 1" Offset, Back cane mount ^{11,18,19,26}	
<input type="checkbox"/> Lateral Thoracic Curved-S/A-Medium	\$650pr
Flush to 1" Offset, Back cane mount ^{11,18,19,26}	
<input type="checkbox"/> Lateral Thoracic Curved-S/A-Large	\$650pr
Flush to 1" Offset, Back cane mount ^{11,18,19,26}	

ACCESSORIES

<input type="checkbox"/> Rear Basket ^{6,7,12}	\$103
<input type="checkbox"/> Cane/Crutch Holder ^{6,7,12}	\$98
<input type="checkbox"/> Oxygen Holder ^{6,7,12,13} E2208	\$103
<input type="checkbox"/> Oxygen Holder ^{6,7,10,11,13} E2208	\$150
(TRU-Balance 2 Tilt & Synergy Seat only)	
<input type="checkbox"/> Stroller handle	\$85
(Synergy Manual Tilt only)	
<input type="checkbox"/> Accessory Mounting Bracket	\$50
(1" adaptor tube for use w/ TB Power Recline)	
<input type="checkbox"/> Walker Holder ^{2,6,7,12}	\$92
<input type="checkbox"/> Flag Holder (SOF on rehab seating)	\$17
<input type="checkbox"/> Weather Cover, Large, Black	\$202
<input type="checkbox"/> Saddle bag(contoured seating only)	\$28
<input type="checkbox"/> Cup holder	\$17
<input type="checkbox"/> Left <input type="checkbox"/> Right	
N/A on desk length arm with S/A mount or toggle box.	
<input type="checkbox"/> Lap Belt	Std.
<input type="checkbox"/> Pediatric 38" N/A w/ limited recline seating	
<input type="checkbox"/> 50" ¹⁸ <input type="checkbox"/> 60" <input type="checkbox"/> 70" <input type="checkbox"/> 80" <input type="checkbox"/> 90"	

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.

Footnotes for pages 5,6 & 7:

1. Includes price for removable, adjustable, quick-height arms & 2-post arms.
2. S/A joystick mount N/A.
3. Measured from seat pan to top of back pan w/ 3-4" gap on tilt & 7-8" gap on recline to accommodate the cushion.
4. N/A w/ flip-back, cane mounted arms.
5. Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
6. Includes mounting bracket and hardware.
7. N/A on reclining seats.
8. N/A on manual recline high-back.
9. S/A laterals interfere w/ flip-back arms.
10. Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
11. S/A laterals interfere w/ flip back arms.
12. N/A w/ tilt systems.
13. N/A w/ flip-back arms.
14. N/A w/ Attendant manual recline.
15. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
16. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
17. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
18. Comes with swing-away / removable hardware; HCPCS code is E1028.
19. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
20. N/A w/ waterfall armpads.
21. Consult specialty control order form for a wide variety of switch types and sizes-- 2 required.
22. Does not include switch.
23. N/A w/ Shape Back.
24. Solid plate is N/A w/manual recline.
25. For use w/ recline systems only.
26. For use with Synergy Seat and TB Tilt only.

Lowest seat-to-floor height

Specialty Seat:	16.75"
Synergy Seat:	15.5"
TB 2 Tilt:	16.25"
TB Recline:	15.5"
High-back Seat:	19.125"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires).
Increased widths may increase seat-to-floor height. Contact customer service for details.
Listed measurements are in the lowest position.
It does not include elevating systems.