

Quantum® 500

Order Form / Price List
Revised: December 2005

*SEAT DIMENSION CODING FOR SPECIALTY & SYNERGY SEATING 20" AND LARGER IS K0108

SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) \$70

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" \$20
- 18" Std.
- 20" \$20

Specialty Seating Kit Depth

- 16" Std. 17" \$20
- 18" \$20 19" \$20
- 20" \$20

Specialty Seating Kit Back Cane Height

- 18" Std. 20" NC
- 22" NC

Back Selection

See Synergy Shape Back Order Form¹⁴

Solid Curved Back \$286

Full width, fully T-nutted, 1" foam with naugahyde sides and back and a rubatex face.

Back Height Solid Curved⁵

- 18" Std. 19" NC
- 20" NC 21" NC
- 22" NC 23" NC
- 24" NC 25" NC

Sling Back

- 16" wide \$98
- 18" wide \$98
- 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
- 20" height NC
- 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

- 2-Post, flip-back^{8,19} Std.
- Flip-back, cane mounted^{8,19} \$95ea
- Removable, adjustable, E0973 quick height¹⁰ \$275ea
- Removable, single-post, fixed E0973 height adjustable, heavy duty \$150ea
 - 8"-12" Adjustment
 - 10"-14" Adjustment
- Pediatric, removable, E0973 \$250ea adjustable, quick height

Armrest Pads

- Full Length, straight, 2"x14" NC
 - Left Right
- Desk Length, straight, 2"x10" NC
 - Left Right
- Full Length, Waterfall, 2"x14"^{11,15} NC
 - Left Right
- Desk Length, Waterfall, 2"x10"^{11,15} NC

*SEAT DIMENSION CODING FOR SPECIALTY & SYNERGY SEATING 20" AND LARGER IS K0108

Synergy Seating Systems³

- Synergy Seating™ System \$434 (Price with Base) K0108
- Synergy Manual Recline Pkg.^{1,2,3,18} \$700 Maximum weight 300 lbs. E1226
- Synergy Manual Tilt Pkg.² \$550 (Max Dimensions 18" W x 18" D)

Maximum weight 160 lbs. K0108
Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

- 10" \$20 11" \$20
- 12" \$20 13" \$20
- 14" \$20 15" \$20
- 16" \$20 17" \$20
- 18" Std. 19" \$20
- 20" \$20

Seat Width Package

- 10"W-14"W 14"W-18"W
- 16"W-20"W

Seat Width Package Manual Tilt

- 10"W-12"W 12"W-14"W
- 14"W-16"W 16"W-18"W

DEPTH

Set Seat Depth At:

- 10" \$20 11" \$20
- 12" \$20 13" \$20
- 14" \$20 15" \$20
- 16" Std. 17" \$20
- 18" \$20 19" \$20
- 20" \$20

Seat Depth Package

- 10"D-14"D 14"D-18"D
- 16"D-20"D

Seat Depth Package Manual Tilt

- 10"D-14"D 14"D-16"D
- 16"D-18"D

BACK CANES

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

- 14" NC 15" NC
- 16" NC 17" NC
- 18" Std. 19" NC
- 20" NC 21" NC
- 22" NC 23" NC
- 24" NC 25" NC

BACK OPTIONS

See Synergy Shape Back Order Form

Solid Curved Back up to 20" W \$281

Back Height:⁵

- 18"H 19"H 20"H 21"H
- 22"H 23"H 24"H 25"H
- Sling Back⁶ \$98

Back height will be the same as back can height.

Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

- Flip-back, cane mounted^{7,9,19} Std.
- 2-Post, flip-back, height^{8,9,19} E0973 \$150ea adjustable, removable armrests
- Removable, adjustable¹⁰ E0973 \$75ea quick height (3" of height adjustment in 1/2" increments)
- Removable, height adjustable E0973 \$75ea heavy duty armrests
 - 8"H-12"H 10"H-14"H
- Pediatric removable E0973 \$75ea adjustable, quick height. (7 3/4" - 10 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)

ARMPADS

- Full length, Straight, 2" X 14" Std.
 - Left Right
- Desk Length, Straight, 2" X 10" Std.
 - Left Right
- Full Length, Waterfall, 2" X 14"¹¹ NC
 - Left Right
- Desk Length, Waterfall, 2" X 10"^{11,12} NC
 - Left Right

HEADRESTS FOR SPECIALTY AND SYNERGY SEATING⁴

Stealth^{16,17} E0955 Large (11" wide x 6" high) \$210

Stealth^{16,17} E0955 Small (9 3/4" wide x 3" high) \$210

Stealth, swing-away (Flip-Down) Head Rest Bracket E1028 \$150

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications^{16,17} E0955 Large (11" wide x 6" high) \$275

Medical Applications^{16,17} E0955 Small (9 3/4" wide x 3" high) \$275

Fixed, Offset Mounting Bracket NC

Footnotes for page 3:

1. Cannot use flip-back, cane mounted arms.
2. Solid back is not included.
3. Includes removable, adjustable, quick height armrests.
4. N/A with sling back.
5. Measured from seat pan to top of back with a gap for the cushion.
6. N/A on Synergy manual tilt or recline.
7. N/A on Synergy Manual Recline.
8. May cause interference with certain dimension configurations. Call Quantum Customer Service.
9. Call for configuration avail. w/ Shape back.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17" from seat pan to top of armrest pad.
11. Only available on removable, adjustable, quick height arms and 2-post arms.
12. Swing-away joystick mount is N/A.
13. Manual systems do not include a headrest.
14. 14" & 22" widths are N/A.
15. Only available on removable, adjustable, quick height arms and pediatric arms.
16. Comes with removable hardware; HCPCS code is E1028.
17. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
18. N/A w/ Shape Back.
19. Stealth S/A laterals interfere w/ flip back arms.

Quantum® 500
Order Form / Price List
Revised: December 2005

TRU-BALANCE TILT

TRU-Balance Power Tilt E1002 \$4995
 55 degree center of gravity (CG) power tilt seating system.

SYSTEM SPECIFICATIONS

Tower Selection

Adjustable-angle Std.

Set Seat Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Width Package

14"W-18"W 16"W-20"W

Set Seat Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Depth Package

14"D-18"D 16"D-20"D

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

Set Back Cane Height At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Seat Overall Back Height²

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Back Type

- See Synergy Shape Back Order Form**
- Solid Curved Back Std.
 Full width, fully T-nutted, 1" foam with naugahyde sides and back and a rubatex face.
- Attendant manual recline^{3,4} \$750
 180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.
- Omit Back -(\$100)

ARMRESTS

- Flip-back, cane mounted^{5,9,10} Std.
Mounted on back canes with full width spreader bar.
- 2-Post, flip-back, hgt. adj, removable^{9,10} NC
- Removable, adj. quick height⁶ NC
- Removable, single-post, height adjustable, heavy duty NC
 - 8"-12" Adjustment
 - 10"-14" Adjustment
- Pediatric, removable, adjustable, quick height NC

Armrest Pads

- Full Length, straight, 2"x14" NC
 - Left Right
- Desk Length, straight, 2"x10" NC
 - Left Right
- Full Length, Waterfall, 2"x14" NC
 - Left Right
- Desk Length, Waterfall, 2"x10"⁸ NC
 - Left Right

ELECTRONICS FOR TILT

- PG 6-KEY VSI NC
- Toggle Box: mounts to side of armrest. NC
 - Left Right

HEAD REST

- Stealth**^{1,7} E0955
 Large (11"wide x 6" high) \$210
- Stealth**^{1,7} E0955
 Small (9 3/4" wide x 3" high) \$210
- Stealth**, swing-away (Flip-Down) Head Rest Bracket E1028 \$150
 Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications**^{1,7}
 Large (11"wide x 6" high) E0955 \$275
- Medical Applications**^{1,7}
 Small (9 3/4" wide x 3" high) E0955 \$275
- Omit Head Rest - To be supplied by provider
- Fixed, Offset Mounting Bracket NC
- Bracket for after market headrest \$180

Lowest seat-to-floor height	
Specialty Seat:	17"
Synergy Seat:	17"
TB Tilt:	19.375"
Contoured Seat:	21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Footnotes for page 4:

1. Comes with removable hardware; HCPCS code is E1028.
2. Measured from seat pan to top of back pan with a gap to accommodate the cushion.
3. N/A w/ flip-back, cane mounted arms.
4. Includes price for removable, adjustable, quick-height arms and 2-post arms.
5. N/A w/ Attendant manual recline.
6. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11"- 17" from seat pan to top of armrest pad.
7. These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
8. Swing-away joystick mount is N/A.
9. Call for configuration availability w/ Shape back.
10. Stealth S/A laterals interfere w/ flip back

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility
 Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.