



DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

**BASE MODEL** **EXPANDABLE CONTROLS**

Description	Price	HCPCS Code
<input type="checkbox"/> Quantum 610 2S-SS	\$5,971	K0822
<input type="checkbox"/> Quantum 610 2S-C	\$5,971	K0823
<input type="checkbox"/> Quantum 610 2SP-C	\$6,171	K0836
<input type="checkbox"/> Quantum 610 2SP-SS	\$6,171	K0835

**COLORS**

**Colors:**

Candy Apply Red     Sunburst Orange  
 Viper Blue         Pearl Pink  
 Onyx Black

**BASE MOUNTED FOOT PLATFORM**

Description	Price
<input type="checkbox"/> Base Mounted Foot Platform Installed	Standard
<input type="checkbox"/> No Base Mounted Foot Platform	No Charge
<input type="checkbox"/> Bell Shaped Foot Platform	\$235

**DRIVE WHEELS**

Description	Price
<input type="checkbox"/> 10" Flat-Free (pair)	Standard
<input type="checkbox"/> 10" Pneumatic (pair)	No Charge




**BATTERY CHARGER**





Description	Price
<input type="checkbox"/> 5 Amp Off-board Charger	Standard

**POWER ELEVATING SEAT**



Description	Price	HCPCS Code
<input type="checkbox"/> 5" Power Elevating Seat Thru VSI 6-Key	\$975	E2300

**NON-EXPANDABLE CONTROLS**

Description	Price
<input type="checkbox"/> PG 4-Key VR2 	Standard
<input type="checkbox"/> PG 6-Key VR2 	\$200
<input type="checkbox"/> PG 6-Key VSI <i>For Use With Power Elevate</i> 	\$200

Description	Price
<i>Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls</i>	
<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller 	\$695
<input type="checkbox"/> Expandable Controller Harness	\$425
<input type="checkbox"/> Q-Logic 2 EX Joystick 	\$700
<input type="checkbox"/> Q-Logic 2 EX Joystick for Lights 	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port 	\$700
<b>Armor Color Options:</b>	
<input type="checkbox"/> Black Soft Touch (Standard color option)	
<b>Colors: \$50</b>	
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Sunburst Orange
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Pearl Pink

**PROGRAMMING**

Description	Price
<input type="checkbox"/> Q-Logic Handheld Programmer 	\$545
<input type="checkbox"/> Q-Logic Programming Cable w/ Disk <sup>1</sup>	\$750
<input type="checkbox"/> PG Handheld Programmer 	\$423

1. Works with Q-Logic 2 EX only

**JOYSTICK HANDLES**

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$85	E2323
<input type="checkbox"/> T-Bar	\$120	E2323
<input type="checkbox"/> Chin Cup	\$85	E2324 <sup>1</sup>
<input type="checkbox"/> Extended I-Handle	\$85	E2323
<input type="checkbox"/> Soft Ball Handle	\$65	E2323
<input type="checkbox"/> Medium Ball Handle	\$65	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$99	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$99	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

**JOYSTICK MOUNTING BRACKETS**

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline <sup>1</sup>	Standard	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-away	\$245	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
Offset: <input type="checkbox"/> Inline <input type="checkbox"/> Offset		
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$350	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Multi-axis	\$100	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Not available with Q-Logic 2 EX Stand Alone Joystick

**CONTOURED SEATING**

**SYNERGY® BACK OPTIONS**

Description	Price	HCPCS Code
<b>Depth Adjustable High Back:</b>		
<input type="checkbox"/> Limited Recline Depth Adjustable High Back w/ Headrest <sup>1</sup>	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 18"W x 18-20"D <input type="checkbox"/> 20"W x 18-20"D		

<b>High Back with Solid Seat Pan:</b>		
<input type="checkbox"/> Limited Recline High-Back with Headrest & Solid Seat Pan	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 16"W x 20"D		
<input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 18"W x 18"D <input type="checkbox"/> 18"W x 20"D		
<input type="checkbox"/> 20"W x 16"D <input type="checkbox"/> 20"W x 18"D <input type="checkbox"/> 20"W x 20"D		

<b>Armrests:</b>		
<input type="checkbox"/> Flip-Up, Height-Adjustable, Removable	\$85 ea.	E0973

<b>Armpads:</b>		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length <sup>2</sup>		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length <sup>2</sup>		

**Special Order Seat**

\_\_\_\_\_ "W x \_\_\_\_\_ "D

Seat size can not exceed 20"W x 20"D

\*\*Contact your inside sales specialist for pricing\*\*

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part.
2. Selection of desk length armpads on a depth-adjustable high back requires them to be added as an additional part.

**MANUAL RECLINER**

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back with Headrest	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back with Headrest and Solid Seat Pan <sup>1</sup>	\$600	E1226

<b>Upholstery:</b>		
<input type="checkbox"/> Gray Vinyl		
<b>Seat Width:</b>		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
<b>Seat Depth:</b>		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
<b>Armpads:</b>		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		

1. Includes black vinyl seatpan. Back available in listed colors.

**SYNERGY® SEATING SYSTEMS**

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy Seating Frame	\$434	
<input type="checkbox"/> Synergy Manual Recline Package <sup>1</sup>	\$700	E1226

<b>Seat Width:</b>					
<input type="checkbox"/> 10" <sup>2</sup>	<input type="checkbox"/> 11" <sup>3</sup>	<input type="checkbox"/> 12" <sup>3</sup>	<input type="checkbox"/> 13" <sup>3</sup>	<input type="checkbox"/> 14"	
<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input checked="" type="checkbox"/> 19"	<input checked="" type="checkbox"/> 20"

<b>Seat Width Package:</b>		
<input type="checkbox"/> 10"W-14"W <input type="checkbox"/> 14"W-18"W <input type="checkbox"/> 16"W-20"W		

<b>Seat Depth:</b>					
<input type="checkbox"/> 10"	<input type="checkbox"/> 11"	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input checked="" type="checkbox"/> 14"	
<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input checked="" type="checkbox"/> 19"	<input checked="" type="checkbox"/> 20"

<b>Seat Depth Package:</b>		
<input type="checkbox"/> 10"D-14"D <input type="checkbox"/> 14"D-18"D <input type="checkbox"/> 15"D-20"D		

<b>Back Cane Height:</b>				
<input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <sup>3</sup> <input type="checkbox"/> 17" <sup>4</sup>				
<input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21"				
<input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24" <input type="checkbox"/> 25"				

1. Not available with Flip Back Arms, Cantilever Arms, Shape Back, AFP or Power Elevate. Requires 22"H canes. Pivot point 2.5" above seat pan.
2. SOF width is required when choosing a power elevating seat on widths 10"-13".
3. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

Description	Price	HCPCS Code
<input type="checkbox"/> Sling Back	\$98	
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Back: <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Choose Lateral Extensions <sup>1</sup> :		
<input type="checkbox"/> 6" for 13"-16" H Back	\$25	
<input type="checkbox"/> 10" for 17"-18" H Back	\$25	
<input type="checkbox"/> 12" For 19"-22" H Back	\$25	
<input type="checkbox"/> 16" For 23"-25" H Back	\$25	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width: <input type="checkbox"/> 16" W <input type="checkbox"/> 18" W <input type="checkbox"/> 20" W	\$675	E2620
Overall Back: <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
Height <sup>2</sup> : <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 17" W <input type="checkbox"/> 19" W	\$675	E2620
Overall Back: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Height <sup>2</sup> : <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 22"W <sup>3</sup>	\$775	E2621
Overall Back: <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Height <sup>2</sup> :		

**Synergy® Back Cane Options:**

<input type="checkbox"/> Standard Back Canes	Standard
<input type="checkbox"/> Omit Push Handles <sup>4</sup>	\$125
<input type="checkbox"/> Flip Forward Back Canes <sup>5</sup>	\$150
<input type="checkbox"/> Omit Push Handles	\$125
<input type="checkbox"/> Offset (Meso) <sup>6</sup>	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	
<input type="checkbox"/> Inset (Endo)	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	

1. Lateral extensions are required to mount lateral thoracic supports.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
3. Available with offset (meso) brackets.
4. Not available on 14"H canes.
5. Not available with Vent Tray, O2 Holder, cane mounted or cantilever armrests. Also not available w/ a Synergy manual recline or manual tilt. May interfere with Lateral Supports and Waterfall Armpads.
6. Dual Post Flip-back Armrests interfere with offset canes.

**SYNERGY® ARMRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> Flip-back, Cane Mounted <sup>1</sup>	No Charge	
<input type="checkbox"/> 2-post, Flip-back, Height-adjustable, Removable <sup>2,3</sup>	\$150 ea.	E0973
<input type="checkbox"/> Single-post, Height-adjustable, Removable <sup>3,4</sup> <input type="checkbox"/> Quad Release Lever	\$75 ea.	E0973
<input type="checkbox"/> Removable, Single-post, Height Adjustable, Heavy Duty <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H	\$75 ea.	E0973
<input type="checkbox"/> Pediatric, Removable <sup>5</sup> <input type="checkbox"/> Quad Release Lever	\$75 ea.	E0973
<input type="checkbox"/> Bariatric, Removable, Height Adjustable <sup>6</sup> <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H	\$75 ea.	E0973
<input type="checkbox"/> Height Adjustable Cantilever Armrests <sup>7</sup> Angle-Adjustable: -30° through +30°	\$195	

1. Interferes with swing-away laterals. Not available on flip forward canes or manual recline.
2. Interferes with swing-away laterals. Height adjustable from 9"-13" in 1" increments.
3. Not available with a depth less than 16".
4. Height adjustable from 11"-17" in 1/2" increments
5. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available
6. Bariatric Armrests require Bariatric Armpads
7. Height adjustment from 11-17" depends on cane height. Not available with Manual Recline, Flip Forward back canes & Shape Back.

**SYNERGY® ARMPADS**

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14") <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard
<input type="checkbox"/> Straight, Desk Length (2"x10") <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Full Length (2"x14") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Desk Length (2"x10") <sup>1</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")* <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$92.50 ea. \$56.69 ea.
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")* <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$70.00 ea. \$53.52 ea.
<input type="checkbox"/> Flat, Bariatric <sup>2</sup> <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.
2. Bariatric Armrests require Bariatric Armpads

**LAPBELTS**

Description	Price
<input type="checkbox"/> 38" Lap Belt <sup>1</sup>	Standard
<input type="checkbox"/> 50" Lap Belt <sup>1</sup>	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Non-Retractable (Sure-Lok) <sup>2</sup>	Standard
<input type="checkbox"/> 60" Padded Lap Belt <sup>3</sup>	\$120
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard

1. Includes mounting hardware. Not available with Manual Recliner.
2. Only needed with the occupied transit kit.
3. Only available with Synergy seating.

**I-FIT PELVIC BELTS**

Description	Price
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded 1.5" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Large (9.5") Padded 2" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Extra Large (11.5") Padded 2" Belt	\$86

**HEADRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
<b>Stealth:</b>		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$64.20	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$80.25	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$151.94	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$151.94	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$151.94	E0955
<input type="checkbox"/> Removable Headrest Bracket <sup>2</sup> (HMO475P)	\$160.50	E1028
<b>Aftermarket Bracket:</b>		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$220.50	E1028

**NON-POWERED LEGRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>60-90° Swing-away Legrests<sup>1</sup></b> (Adjustable in 5° increments)	<b>\$162</b>	
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
<i>Standard Extruded Extensions:</i>		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Standard Extruded Footplates:</i>		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates <sup>2</sup>	\$85 ea.	
<i>Angle Adjustable Extensions:</i>		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Angle-Adjustable Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		
<input type="checkbox"/> <b>High Mount, Clamp-On Footrest*</b>		
<input type="checkbox"/> Angle Adjustable		
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> <b>Heavy Duty Drop-In Legrests<sup>5</sup></b>	<b>\$420</b>	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$85 ea.	
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		

1. Possible caster wheel interference when adjusted to less than 70 degrees.
2. Includes Heel Loops
3. Causes interference on seat widths less than 18".
4. Causes interference on seat widths less than 20".
5. Only available with Synergy seating.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**NON-POWERED LEGRESTS (CONT'D)**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Elevating Legrests<sup>1</sup></b>	\$325	K0195
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length:		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$85 ea.	
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>		
<input type="checkbox"/> <b>90° Pediatric Swing-Away Legrests<sup>1</sup></b>	\$162	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		

1. Extensions are measured from the seatpan to footplate.  
 2. Includes heel loops.  
 3. Causes interference on seat widths less than 18".  
 4. Causes interference on seat widths less than 20".

**ADJUSTABLE CENTER MOUNT FOOT PLATFORM**

Description	Price
<input type="checkbox"/> <b>Adjustable Center Mount Foot Platform<sup>1</sup></b>	\$250
Footplates:	
<input type="checkbox"/> 10"W 11"D Split Footplates	
<input type="checkbox"/> 10"W 8"D Split Footplates	
<input type="checkbox"/> 10"W 6"D Split Footplates	
Extensions: <sup>2</sup>	
<input type="checkbox"/> 8"-12" Ext. <input type="checkbox"/> 10"-14" Ext.	
<input type="checkbox"/> 10"-16" Ext.	
Calf Supports <sup>3</sup> :	
<input type="checkbox"/> Calf Pads (pair)	
<input type="checkbox"/> Small (6"H x 6"W)	\$150
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150
<input type="checkbox"/> Calf Panel (ea.)	
<input type="checkbox"/> Small (8"H x 18"W)	\$150
<input type="checkbox"/> Large (8"H x 24"W) <sup>4</sup>	\$200
Optional Modular Footplate Extension Components/Kits <sup>3,5</sup>	
<input type="checkbox"/> 4" Width Extension Kit	\$50
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/> 4" Side Guard <sup>6</sup>	\$50
<input type="checkbox"/> 6" Side Guard <sup>7</sup>	\$50
Additional Optional Accessories <sup>5</sup> :	
<input type="checkbox"/> Footplate Connector	\$15
<input type="checkbox"/> Wheel Bumper	\$12.50

1. Only available with Synergy seating.  
 2. Extension length is measured from top of the seat pan to footplate.  
 3. Choose in Rigging Accessory Section in iQuantum.  
 4. Only available on 20"W or greater.  
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.  
 6. Only available with 10"W x 6"D Split Footplates  
 7. Only available with 10"W X 11"D and 10"W x 8"D Split Footplates

**ELECTRONICS FOR POWER FOOT RIGGINGS**

Description	Price	HCPCS Code
<input type="checkbox"/> Through VR2 6-Key	No Charge	
<input type="checkbox"/> Through Q-Logic 2 EX	\$1,895	E2310

**CENTER MOUNT POWER ARTICULATING FOOT PLATFORM**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Center Mount Power Articulating Foot Platform<sup>1</sup></b>	\$2995	See note <sup>2</sup>
Footplates:		
<input type="checkbox"/> 10"W 11"D Split Footplates		
<input type="checkbox"/> 10"W 8"D Split Footplates		
<input type="checkbox"/> 10"W 6"D Split Footplates		
Extensions <sup>3</sup> :		
<input type="checkbox"/> 10"-16" Ext.		
Calf Supports <sup>4</sup> :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) <sup>5</sup>	\$200	
Optional Modular Footplate Extension Components/Kits <sup>3,6</sup>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Side Guard <sup>7</sup>	\$50	
<input type="checkbox"/> 6" Side Guard <sup>8</sup>	\$50	
Additional Optional Accessories <sup>5</sup> :		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumpers	\$12.50 ea.	

1. Only available with Synergy seating.  
 2. Please verify the appropriate billing code to be used for this item with all third-party payors.  
 3. Extension length is measured from top of the seat pan to footplate.  
 4. Choose in Rigging Accessory Section in iQuantum.  
 5. Only available on 20"W or greater.  
 6. Assembly required. May require seating configuration adjustments to avoid front caster interference.  
 7. Only available with 10"W x 6"D Split Footplates  
 8. Only available with 10"W x 11"D and 10"W x 8"D split footplates

**FRONT RIGGING ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumper	\$12.50 ea.	
<input type="checkbox"/> Heel Loops (pair) <sup>1</sup>	\$50	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) <sup>2</sup>	\$80	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Standard with Angle Adjustable or Multi-Axis footplates.  
 2. Only available with Swing-Away legrests.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**SEAT CUSHIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Simplicity Cushion</b> <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<b>Solution® 1 Cushion</b> <i>Skin Protection:</i>		
<input type="checkbox"/> 14" x 14" - 20" x 20"	\$430	E2603
<b>Solution® Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$430	E2607
<b>Spectrum® Gel Cushion</b> <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 20"	\$530	E2605
<b>Spectrum® Visco Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 20"	\$530	E2607
<b>TRU-Comfort 2™ 4-Way Stretch Seat Cushion</b> <i>Skin Protection</i> <i>Includes 4-way stretch material cover</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$400	E2603

**ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> Rear Basket <sup>1</sup>	\$118	
<input type="checkbox"/> Cane/Crutch Holder <sup>1</sup>	\$110	E2207
<input type="checkbox"/> Weather Cover, Medium, Black	\$168	
<input type="checkbox"/> Cup Holder	\$20 <sup>2</sup>	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Saddle Bag	\$35	
<input type="checkbox"/> Flag Holder	\$25	
<input type="checkbox"/> Oxygen Holder, Contoured Seating <sup>1</sup>	\$118	E2208
<input type="checkbox"/> Oxygen Holder, Synergy Seat <sup>3</sup>	\$150	E2208
<input type="checkbox"/> Walker Holder <sup>4</sup>	\$103	
<input type="checkbox"/> Tablet Holder <sup>4</sup>	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder <sup>5</sup>	\$140	

1. Includes mounting hardware. Not available on Manual Recliner.  
 2. No charge on Synergy seating.  
 3. Not available with Synergy manual recline or manual tilt. Includes mounting hardware.  
 4. Only available with Synergy seating.

**ADDITIONAL ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Stealth Foot Boxes (pair)</b>	\$498.75	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> Large		
<input type="checkbox"/> <b>Stealth Padded Shoe Holders (pair)</b>	\$266	K0108
<input type="checkbox"/> Small		
<input type="checkbox"/> Medium		
<input type="checkbox"/> X-Large		
<input type="checkbox"/> <b>Stealth Elbow Support</b>	\$183.75 ea.	E1028
<i>Side:</i>		
<input type="checkbox"/> Bracket - Left		
<input type="checkbox"/> Bracket - Right		
<i>Pad:</i>		
<input type="checkbox"/> 3" x 3" (ARMS-ES-3x3)	\$52.50 ea.	K0108
<input type="checkbox"/> 4" x 4" (ARMS-ES-4x4)	\$52.50 ea.	K0108
<input type="checkbox"/> <b>Stealth i-Fit Dynamic Chest Straps</b>	\$142	E0960
<input type="checkbox"/> Upper Classic, Medium (IFPS1103)		
<input type="checkbox"/> Upper Classic, Large (IFPS1104)		
<input type="checkbox"/> Upper Contour, Medium (IFPS2103)		
<input type="checkbox"/> Upper Contour, Large (IFPS2104)		
<input type="checkbox"/> <b>Stealth i-Fit Chest Strap Belt Anchor</b>	\$90	
<input type="checkbox"/> Medium (IFPS4113)		
<input type="checkbox"/> Large (IFPS4114)		

**LATERALS**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Stealth Thoracic Laterals (pair)*</b>		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away Rail Mounts	\$472.50	E1028
<i>For Use on Synergy Back and Tru-Comfort</i>		
<i>Pad:</i>		
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$210	E0956
<input type="checkbox"/> <b>Stealth Pelvic/Thigh Laterals*</b>		
<i>Side:</i>		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$262.50	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$400	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$236.25	E1028
<input type="checkbox"/> Flip Down, Full Surface Contact Bracket (ea.)	\$236.26	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$75	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$75	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$75	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) <sup>2</sup>	\$26.25	K0108
<input type="checkbox"/> 2 1/4" Depth/Height Adapter Plate (ea.) <sup>3</sup>	\$26.25	K0108

1. Swing-away laterals interfere with flip-up armrests.  
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.  
 3. Only available on the flip down bracket.

**ADDITIONAL HEADREST ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Stealth Facial Laterals<sup>1</sup></b>		
<input type="checkbox"/> Right Side, 4.5" x 2" Pad (CPS9-R-HW & 900-Q)	\$347.75	E1028
<input type="checkbox"/> Left Side, 4.5" x 2" Pad (CPS9-L-HW & 900-Q)	\$347.75	E1028
<input type="checkbox"/> Right and Left 4.5" x 2" Pad (CPS9-B-HW & 900-Q)	\$453.75	E1028
<input type="checkbox"/> Egg Switch Mount (ESM)	\$52.50 ea.	
<input type="checkbox"/> Facial Pad, 6" x 2" (Q920-Q)	\$53.50 ea.	K0108
<input type="checkbox"/> Figure Eight Bracket with 13" Rod and Egg Switch Mount (F8678-13-Q)	\$220	K0108
<input type="checkbox"/> Swing-away Joystick Mount, 18" Rod (SUS18-Q)	\$250	K0108

1. Only available with Stealth Comfort Series Headrests.

**BATTERY OPTIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> U1 32 Amp, AGM	\$133 ea.	E2365
<i>Installation:</i>		
<input type="checkbox"/> Batteries Installed <input type="checkbox"/> Do Not Install		

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum