



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Alternative Controls and Electronics Order Form

Account Number: _____ Date: _____ PO Number: _____
Provider Name: _____ Marked for: _____
Contact: _____ Ship to Address: _____
Phone: _____ Fax: _____ City: _____ State: _____ Zip: _____
Email: _____

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com).


HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice.

1. PROPORTIONAL SECONDARY CONTROLS

1.1. Mo-Vis™ Joysticks

Requires enhanced display.

- | | |
|---|--|
| ☐ Mo-Vis™ All Round \$1,703
Joystick (IDM-ALLROUND-9)
Part: ST-IDMALLRND9-Q. HCPCS: E2373
May be used as chin, standard, or attendant joystick. | ☐ Mo-Vis™ Micro Joystick (IDM-MICRO-9) \$3,492
Part: STIDM-MICRO-9Q. HCPCS: E2312
Designed to be controlled by the finger, lip, or tongue.
Supplied with 2 tips for driving option: a ball or cup mount. |
| ☐ Mo-Vis™ Multi Joystick (IDM-MULTI-9) \$3,492
Part: STIDM-MULTI-9Q. HCPCS: E2312
Designed to be controlled by the chin, finger, or lip. | |



1.2. Stealth Joysticks

- | | |
|--|--|
| ☐ Stealth Precision Mini Proportional \$3,492
Joystick (IDPMP-9)
Part: ST-IDPMP-9-Q. HCPCS: E2312
Designed to work for users with minimal range of motion and strength. Additional tops available through SOF. | ☐ Stealth Mushroom Joystick \$2,158
(MJM250Q)
Part: ACCASMB2380 |
|--|--|



1.3. Brackets

1. Joystick to Gatlin Mounting Bracket

- | | |
|---|--|
| ☐ Stealth Gatlin Mount Bracket for \$199
Stand Alone Joystick (GAT452CJ)
Part: ST-GAT452CJ-Q | ☐ Stealth Gatlin Mount Bracket for \$199
Mo-Vis Multi Joystick (GAT461M)
Part: ST-GAT461M-Q |
| ☐ Stealth Gatlin Mount Bracket for \$199
NE, NE+, or Q-Logic (GAT456QL)
Part: ST-GAT456QL-Q | ☐ Stealth Gatlin Mount Bracket for \$203
Mini Proportional Joystick (GAT486)
Part: ST-GAT486-Q |
| ☐ Stealth Gatlin Mount Bracket for \$199
Mo-Vis Micro Joystick (GAT451M)
Part: ST-GAT451M-Q | ☐ Stealth Gatlin Mount Bracket for \$142
Mo-Vis All Round Joystick (GAT472P-M)
Part: ST-GAT472P-M-Q |
| ☐ Stealth Gatlin Mount Bracket for \$203
Enhanced Display (GAT454ED)
Part: ST-GAT454ED-Q | |



2. Stealth Adjustable Tray and Bracket

- | | |
|---|--|
| ☐ Tray Mount for Mo-Vis Micro Joystick \$511
Left Side Control (GAT440M-51M-L)
Part: ST-GAT44051ML-Q | ☐ Tray Mount for Mo-Vis Micro Joystick \$511
Right Side Control (GAT440M-51M-R)
Part: ST-GAT44051MR-Q |
|---|--|



3. Stealth i-Drive Control Harness

Available with Mo-Vis Multi and Micro Joystick.

- | | |
|---|--|
| ☐ Pediatric, i-Drive Control Harness (IDCH-P) \$937
Part: ST-IDCH-P-Q. HCPCS: K0108 | ☐ Medium, i-Drive Control Harness (IDCH-M) \$937
Part: ST-IDCH-M-Q. HCPCS: K0108 |
| ☐ Small, i-Drive Control Harness (IDCH-S) \$937
Part: ST-IDCH-S-Q. HCPCS: K0108 | ☐ Large, i-Drive Control Harness (IDCH-L) \$937
Part: ST-IDCH-L-Q. HCPCS: K0108 |



Proportional secondary controls section continued on following page.

1. PROPORTIONAL SECONDARY CONTROLS CONTINUED**4. Gatlin Mounts**


- ☐ Stealth Gatlin Mount, Cane Mount, \$653
for Chin Control (A) (GAT406C)
Part: FRMASMB8304
Includes a receiver block for tube-style mounting.

- ☐ Stealth 8 Position, Gatlin Midline Mount (B) \$653
(GAT400Q)
Part: ST-GAT400Q-Q
Requires a receiver block, not included.

5. Other Mounts/Brackets

- ☐ TRU-Balance 2 Power Recline..... \$50.57
Accessory Bracket
Part: FRMASMB7322
Allows cane mounted accessories to be mounted to a TRU-Balance 2 Recline system.

- ☐ Therafin Softline Bridge and Collar (32685) \$313.95
Part: INDPART2152
- ☐ Therafin Semi-Rigid Bridge and Collar (32669)..... \$227
Part: INDPART2895

**2. SWITCH CONTROLS****2.1. Head Array**

- i-Drive 4.0 Bluetooth Tri-array Series,..... \$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330

- i-Drive 4.0 Bluetooth Pro Series,..... \$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330

- Adult, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330

- Junior, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330

- Adult, i-Drive Pro Series Head Array \$4,825
(IDHBT300-1)
Part: ST-IDHBT300-1-Q. HCPCS: E2330
Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.

**2.2. Sip & Puff Module****1. Sip & Puff Module**

- ☐ Q-Logic 3 Sip & Puff Module \$1,990
Part: CTL168056. HCPCS: E2325
Select Sip & Puff component below.

- ☐ Q-Logic 2 Sip & Puff Module \$2,185
Part: KITASMB1697. HCPCS: E2325
Select Sip & Puff component below.

**2. Sip & Puff Straw and Mount**

- ☐ Stealth Sip and Puff for Comfort Series,..... \$318
Left Mount (SP625-CL)
Part: ST-SP625-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.


- ☐ Stealth Sip and Puff for i-Drive Head \$318
Array (SP625-U)
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive head arrays.

- ☐ Stealth Sip and Puff for Comfort Series,..... \$318
Right Mount (SP625-CR)
Part: ST-SP625-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

**3. Sip & Puff Accessories**

- ☐ Sip & Puff Straws, 12 pack (SP625-S12)..... \$35
Part: ST-SP625-S12-Q

- ☐ TRU-Balance 2 Power Recline..... \$50.57
Accessory Bracket
Part: FRMASMB7322
Allows cane mounted accessories to be mounted to a TRU-Balance 2 Recline system.

**2.3. Egg Switches**

- ☐ Red Egg Switch (ESRED) \$87
Part: SWTMCRO1077

- ☐ Yellow Egg Switch (ESYEL) \$87
Part: SWTMCRO1080

- ☐ Green Egg Switch (ESGRN) \$87
Part: SWTMCRO1078

- ☐ Black Egg Switch (ESBLK)..... \$87
Part: SWTMCRO1081

- ☐ Blue Egg Switch (ESBLU)..... \$87
Part: SWTMCRO1079

**2.4. Ablenet Switches**

- ☐ Round, Mini Cup Switch w/ Stickers (58080)..... \$93
Part: SWTMCRO1064
Single input, mechanical button switch. Includes 5 stickers.

- ☐ Micro Light Switch (58500) \$104
Part: SWTMCRO1044
Single input, feather touch, mechanical button switch.

- ☐ Black Spec Switch (100-SPBK)..... \$85
Part: SWTMCRO1116

- ☐ Pillow Switch (58750) \$125
Part: SWTMCRO1052

- ☐ Yellow Spec Switch (100-SPY)..... \$85
Part: SWTMCRO1115

- Single input, feather touch, mechanical button switch with a smooth, soft surface.



Switch controls continued on following page.

2. SWITCH CONTROLS CONTINUED**2.5. Mo-Vis Twister Satellite Buttons**

Single input, mechanical button switch with light activation force activated at any point on the switch. The Satellite Twister and its attached bracket only fits into the housing of the the All-Round Joystick.

- | | |
|--|--|
| ☐ Red Twister Satellite Button (IDM-17)\$449
Part: ST-IDM-17-Q | ☐ Blue Twister Satellite Button (IDM-20)\$449
Part: ST-IDM-20-Q |
| ☐ Black Twister Satellite Button (IDM-18)\$449
Part: ST-IDM-18-Q | ☐ Green Twister Satellite Button (IDM-21)\$449
Part: ST-IDM-21-Q |
| ☐ Yellow Twister Satellite Button (IDM-19).....\$449
Part: ST-IDM-19-Q | |

**2.6. Mo-Vis Twister Basic Buttons**

Single input, mechanical button switch with light activation force activated at any point on the switch.

- | | |
|--|---|
| ☐ Yellow Twister Basic Button (IDM-24) \$313
Part: ST-IDM-24-Q | ☐ Blue Twister Basic Button (IDM-27) \$313
Part: ST-IDM-27-Q |
| ☐ Black Twister Basic Button (IDM-25) \$313
Part: ST-IDM-25-Q | ☐ Green Twister Basic Button (IDM-28) \$313
Part: ST-IDM-28-Q |
| ☐ Red Twister Basic Button (IDM-26) \$313
Part: ST-IDM-26-Q | |

**2.7. Mo-Vis Twister Dome Buttons**

Includes fastener to mount to flat surfaces.

- | | |
|--|--|
| ☐ Red Twister Dome Button (IDM-53) \$195
Part: ST-IDM-53-Q | ☐ Green Twister Dome Button (IDM-56) \$195
Part: ST-IDM-56-Q |
| ☐ Black Twister Dome Button (IDM-54) \$195
Part: ST-IDM-54-Q | ☐ Yellow Twister Dome Button (IDM-57).....\$195
Part: ST-IDM-57-Q |
| ☐ Blue Twister Dome Button (IDM-55) \$195
Part: ST-IDM-55-Q | |

**2.8. Single Switch Accessories**

- | | |
|--|--|
| ☐ Single Switch Extension, 6'\$30
(Abelnet:43330)
Part: INDPART2824 | ☐ 5 to 1 Switch Adapter \$327
Part: HARUSHD2101
Connects up to 5 single switches to connect to 1 9-pin (DB9) connection to act as a multiple switch device. |
| ☐ 4-Switch to AAM Adapter\$69
Part: DWR1265H064
Allows 4 separate mechanical switches to be connected to AAM to control seat functions. | ☐ 2 to 1 Switch Adapter \$15
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch. |

3. ENVIRONMENTAL CONTROLS**3.1. ECU Module**

- | | |
|--|--|
| ☐ Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057 | ☐ Q-Logic 2 ECU Module Kit \$2,520
Part: KITASMB1696
Includes ECU Module and Bus Cable. |
|--|--|

3.2. Additional Accessories

- | | |
|---|---|
| ☐ ECU Bus Cable for Q-Logic 2.....\$447
Part: DWR1111H046 | ☐ Q-Logic 2 Multiplier Harness\$232
Part: DWR1111H053
Provides 2 additional bus connections. |
|---|---|


**4. STEALTH SWING-AWAY MOUNTS AND BRACKETS****4.1. Swing-Away Mount**

- | | |
|--|--|
| ☐ Swing-Away Joystick Mount with 18" Rod (A) \$276
Part: ST-SUS18-Q
Attaches to the Stealth Comfort Series headrests. | |
|--|--|

4.2. Brackets for Swing-Away Mount

These brackets mount a joystick to the Swing-Away Joystick Mount with 18" rod or similar hardware.

- | | |
|--|---|
| ☐ Swing-Away Joystick Mount for Micro \$170
Joystick (SM600-651) (B)
Part: ST-SM600-651-Q | ☐ Swing-Away Joystick Mount for\$145
Precision Mini Proportional Joystick (SM671) (D)
Part: ST-SM671-Q |
| ☐ Swing-Away Joystick Mount for Multi \$170
Joystick (SM600-661) (C)
Part: ST-SM600-661-Q | |



Account Number: _____ Date: _____ Marked for: _____

5. USB CHARGER WITH MOUNTING HARDWARE

The USB charger module can be connected anywhere on the system bus and provides a dedicated charging port for mobile devices according to USB Battery Charging Specifications Revision 1.2. It is intended to charge only one device at a time at a current of 1A or less. The USB charger module is not to be used with passive USB hubs or splitters.

- ☐ Q-Logic 2 USB Mobile Device Charger \$184
Part: ELE124122
- ☐ Q-Logic 3 USB Mobile Device Charger \$184
Part: ELE168069



6. POWERTECH VENT POWER MODULE

6.1. Powertech Vent Power Module Kit

- ☐ Powertech LTV Style Vent Power Module Kit.....\$2,220
Part: ELEASMB7110034. HCPCS: K0108
Powers a LTV style vent from the wheelchair's batteries. Includes compatible components below.
- ☐ Powertech Trilogy Vent Power Module Kit\$2,220
Part: ELE143268. HCPCS: K0108
Powers a Trilogy vent from the wheelchair's batteries. Includes compatible components below.

6.2. Powertech Vent Individual Items

- ☐ Powertech Vent Module (06076) \$1,950
Part: CTLDC8110005
- ☐ Powertech Vent Harness, 24V Input (060774)\$146.18
Part: HAR8110008
- ☐ Powertech Vent Harness, 12V Output (060770) \$219.28
Part: HAR8110009
- ☐ Powertech Trilogy Vent Harness, 12V Output \$219.28
(060771)
Part: ELE143265

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2024 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



[youtube.com/
QuantumRehab](https://www.youtube.com/QuantumRehab)



[facebook.com/
QuantumRehab](https://www.facebook.com/QuantumRehab)



[twitter.com/
quantumrehab](https://www.twitter.com/quantumrehab)



[instagram.com/
quantum_rehab](https://www.instagram.com/quantum_rehab)