



**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: [quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)

## Q6 Edge® 2.0 Order Form with TRU-Balance® 3 Positioning - ODJFS Only

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
Email: \_\_\_\_\_  
PO Number: \_\_\_\_\_  
Marked for: \_\_\_\_\_  
Ship to Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_



### INTRODUCTION

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Selections with HCPCS codes may have a part number noted with the selection. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.

If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

**Thank you for choosing Quantum®!**


### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

## BASE MODEL

**Required.** Select a power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum® Q6 Edge® 2.0 3S-SS ..... \$7,595  
Part: Q6 EDGE 2.0 3S-SS MFG Sugg HCPCS: K0848/K0011  
Solid seat pan standard with no power option.

Quantum® Q6 Edge® 2.0 3MP-SS ..... \$7,595  
Part: Q6 EDGE 2.0 3MP-SS MFG Sugg HCPCS: K0861/K0011  
Multiple actuator power option.

Quantum® Q6 Edge® 2.0 3SP-SS ..... \$7,595  
Part: Q6 EDGE 2.0 3SP-SS MFG Sugg HCPCS: K0856/K0011  
Single power option.

## BASE COLOR

**Required.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



## DRIVE MOTORS

**Required.**

6 mph Drive Motors ..... Standard

Accu-Trac Motor Technology ..... \$1,295  
HCPCS: K0108

Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

## BASE MOUNTED FOOT PLATFORM

**Required.** May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform ..... Standard

Omit base mounted foot platform ..... No Charge

Bell Shaped, Base Mounted Foot Platform ..... \$235

## DRIVE WHEELS

**Required.** Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

3"x14" Flat Free Black Tires ..... Standard

3"x14" Pneumatic Black Tires ..... No Charge

4"x14" Knobby, Flat Free Black Tires ..... \$551  
Not available with iLevel.

4"x14" Knobby, Pneumatic Black Tires ..... \$551  
Not available with iLevel.

### Rim Color

Rim colors other than the standard are available through a special order and extended lead time.

Standard Rims ..... Standard

Match Rims to Base Color ..... \$85

Non-matching Rims ..... \$85  
Select color below

Red Corvette

Back in Black

Pretty in Pink

Lime Light

Midnight Blue

Yellow Blaze

Electric Blue

White Satin

Deep Purple

Grey Street

Orange Crush

## CASTER WHEEL SIZE

**Required.**

6" Caster Wheels ..... Standard

8" Caster Wheels ..... \$150  
Allows smoother transitions over thresholds and obstacles.

## BASE LIGHTING PACKAGES

**Optional.** Select one type of optional lighting packages below. Lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights ..... \$250  
Fender lights that operate only by pressing included single switch button.

Automatic Fender Lights ..... \$275  
Fender lights that power on when the controller is activated and can be operated by the included single switch button.



Fender Lights

## ELECTRONICS

**Required.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Q-Logic Controls Order Form](#).

### Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick ..... Standard  
Available with Tilt Only through Toggle.

NE+ 6-Key Joystick (1 Actuator)..... No Charge  
Control of one actuator through Joystick on single power base model (3SP-SS) only.

No Joystick, Setup for NE+ Joystick... No Charge

NE+ 6-Key Joystick (2 Actuators) ..... \$200  
Control up to 2 actuators through joystick.

NE+ 8-Key Joystick (2 Actuators) ..... \$400  
Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

### Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

#### Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module ..... \$695  
HCPCS: E2377

Q-Logic 2 Expandable Controller..... \$425  
Harness  
Part: DWR1111H065 HCPCS: E2313



Q-Logic 2 EX Controller Module

#### Q-Logic 2 Joystick

Q-Logic 2 EX Joystick..... \$700

Q-Logic 2 Ex Joystick for Seat Mounted, LED Light Package ..... \$700

No Joystick, Setup for Q-Logic

Stand Alone Joystick with ..... \$700  
External Charger Port

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



Stand Alone Joystick

#### Q-Logic 2 Joystick Shroud

**Required with Q-Logic 2 EX Joystick selection above.** Reference colors shown on page 2.

Black, Soft Touch ..... Standard

Non-matching Joystick Shroud ..... \$50

Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze



Q-Logic 2 EX Joystick with Lime Light

## JOYSTICK HANDLES

**Optional.** The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$71  
Part: ACCASMB1860 HCPCS: E2323

Chin Cup (C) ..... \$71  
Part: ACCASMB1861 HCPCS: E2323  
(E2324 when chin control)

T-Bar (D) ..... \$107  
Part: ACCASMB1858 HCPCS: E2323

Medium Ball Handle (E)..... \$50  
Part: ACCASMB1859 HCPCS: E2323

Soft Ball Handle (F)..... \$50  
Part: ACCASMB1856 HCPCS: E2323

Extended I-Handle (G)..... \$71  
Part: ACCASMB1857 HCPCS: E2323

Bodypoint 3" Goal Post Handle..... \$69  
Part: PLSKNCP1057 HCPCS: E2323

Bodypoint 4" Goal Post Handle..... \$69  
Part: PLSKNCP1058 HCPCS: E2323



Joystick Handles

## JOYSTICK MOUNTING BRACKETS

**Required.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket or See Positioning Positioning Components Page or [Order Form](#)

### Fixed, Inline Joystick Bracket

Fixed, Inline, Left Side Mount ..... Standard

Fixed, Inline, Right Side Mount ..... Standard

### Swing-Away (Retractable) Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount..... \$245  
NE/NE+ Part: ACC142526 HCPCS: E1028  
Q-Logic: ACC142524

Swing-away, 1" Offset, Left Side Mount ..... \$245  
NE/NE+ Part: ACC142747 HCPCS: E1028  
Q-Logic: ACC142742

Swing-away, Inline, Right Side Mount ..... \$245  
NE/NE+ Part: ACC142527 HCPCS: E1028  
Q-Logic: ACC142525

Swing-away, 1" Offset, Right Side Mount.... \$245  
NE/NE+ Part: ACC142748 HCPCS: E1028  
Q-Logic: ACC142743

### Retract4 Swing-Away Arm Kit

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Retract4, Left Side Mount ..... \$250  
NE/NE+ Part: MEC138876 HCPCS: E1028  
Q-Logic: MEC138860

Retract4, Right Side Mount..... \$250  
NE/NE+ Part: MEC138876 HCPCS: E1028  
Q-Logic: MEC138860

### Multi-Axis

Allows height adjustment of joystick.

Multi-axis, Left Side Mount.....\$100  
NE/NE+ Part: ACC142528 HCPCS: E1028  
Q-Logic: ACC131216


Multi-axis, Right Side Mount .....\$100  
NE/NE+ Part: ACC142528 HCPCS: E1028  
Q-Logic: ACC131216

### Stealth Swing Away/Flip Down Bracket


Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket ..... \$275  
Left Side Mount  
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028  
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028


Stealth Swing-Away/Flip Down Bracket ..... \$275  
Right Side Mount  
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028  
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028



Swing Away Joystick  
Inline Mount (Right)



Swing Away Joystick  
Offset Mount (Right)



Retract4



Multi-Axis Joystick  
Mount



Swing-Away/Flip Down  
Bracket (Stand Alone  
Hardware Shown)

## ELECTRONICS ACCESSORIES

**Optional.**

USB Mobile Device Charger .....\$175  
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

Multiplier Harness .....\$221.21

Q-Logic EX/NE Programmer ..... \$545

Q-Logic PC Programming Disk with Cable... \$750



USB Media  
Device Charger

## TRU-BALANCE® 3 SEATING

**Required.** Select one seating system. Be sure to complete the Seat Dimension, Seat to Floor, and back selections.

### TRU-Balance® 3 Configuration

TRU-Balance® 3 Static Seat .....\$2,295  
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance® 3 Power Recline ..... \$7,595  
Part: PTOINDV3444 HCPCS: E1004  
175° of power recline.

TRU-Balance® 3 Power Tilt .....\$6,995  
Part: PTOINDV3443 HCPCS: E1002  
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance® 3 Power Tilt & Recline .....\$11,895  
Part: PTOINDV3446 HCPCS: E1007  
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

### 10" Power Lift or Manual Recline

**Optional.** Select one option below.

TRU-Balance® 3 10" Lift with iLevel® .....\$2,495  
Part: PTOINDV3773 HCPCS: E2300  
10" of scissor lift with iLevel® System. 22"W x 22"D maximum seat size.

Attendant Manual Recline ..... \$500  
Part: PTOINDV3501 HCPCS: E1226  
Not available with NE+ joystick with Tilt through toggle. Available with Static Seat or Power Tilt. Not available with Power Lift.



## TRU-BALANCE® 3 SEATING CONTINUED

### Transit Kit

**Optional.** The Q6 Edge® 2.0 is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option..... \$200

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://www.ercwts.org>.

### Seat Dimensions

**Required with seat selection.** Select width and depth package and enter desired set width in that package range, available in 1" increments. Configurations disallowed or not shown may be available through special order. Maximum seat size is 22"W x 22"D.

#### Seat Width

12"-16" Width Package

14"-16" Width Package

16"-20" Width Package

20"-24" Width Package  
(Maximum Width is 22")

Set Width at: \_\_\_\_\_ (Max seat width 22")

#### Seat Depth

12"-16" Depth Package

14"-16" Depth Package

16"-20" Depth Package

20"-24" Depth Package  
(Maximum Depth is 22")

Set Depth at: \_\_\_\_\_ (Max seat depth 22")

### Seat to Floor Height

**Required.** Select a seat to floor height compatible with the seat functions previously selected.

16.5" Seat to Floor

Not available with Power Tilt or  
Power Lift configurations.

17.5" Seat to Floor

18.5" Seat to Floor

19.5" Seat to Floor

Not available with Static Seat  
or Power Recline Only.

### Back Selection

**Required.** Match choice in this section with the options in the back cushion configuration section.

Static Back ..... Standard  
Select a back cushion configuration in the subsequent section.

Back Cane Option, Omit Back Pan/ Cushion .....(-\$125)  
Available only with Static Back and Power Tilt (1 actuator) selections.  
Select cane type and the back cane height below.

Flip Forward Back ..... \$150  
Available with Static Back and Power Tilt (1 actuator). Only available  
with 2-post, flip back or pediatric armrests. Not available with 12"-14"  
width or with attendant manual recline. Select a back cushion type in  
the subsequent section.

#### Back Cane

Select with Back Cane Option selected above.

Angle Adjustable..... Standard

1" Endo ..... \$35

Flip Forward ..... \$25

2" Endo ..... \$35

#### Back Cane Height

Select with Back Cane Option selected above.

16"

17"

20"

21"

24"

25"

### Back Cushion Configuration

**Required except when back cane option selected.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with Special Order and extended lead time.

Sport Back..... No Charge  
Complete following width package and size.

TRU-Comfort 4-way Stretch, 20" Width or less..... \$595  
Complete following width package and size. HCPCS: E2620

ADI Q-Back, Active Contour ..... \$525  
Complete following width package and size. HCPCS: E2613

TRU-Comfort 4-way Stretch, 22" Width or greater ..... \$695  
Complete following width package and size.. HCPCS: E2621

ADI Q-Back, Deep Contour..... \$595  
Complete following width package and size. HCPCS: E2615

Omit Back Pan/Cushion..... No Charge  
Select a following width package.

ADI Q-Back, Extra Deep Contour..... \$650  
Complete following width package and size. HCPCS: E2620  
20" Wide Extra Deep Contour is not compatible with Gel and Waterfall  
Armpads, Endo/Meso Downposts, 2 Post Arms, or Ottobock arm blocks  
unless wider arm bars are installed on the chair.

Aftermarket Back Interface ..... \$200  
Enter the preset back width below. Overall seat depth may need to be  
adjusted to accommodate certain backs. May interfere with gel/wa-  
terfall armpads and arm through. Does not include brackets to mount  
aftermarket backs. A Special Order is required on back widths less  
than 16".



## TRU-BALANCE® 3 SEATING CONTINUED

### Back Cushion Width

Available with +/-4" of seat width.

Match Seat Width Selections

10"W-14"W Package

16"W-20"W Package

20"W-24"W Package

Set Width at: \_\_\_\_\_

### Back Cushion Overall Height

**Static Seat and Power Tilt available height range:** 16"-25"H. *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.*

**Power Recline available height range:** 19"-28"H. *Odd heights include 6" gap and even heights include 7" gap.* See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

	Tilt or Static Width												Recline Width											
	12"	14"	16"	17"	18"	19"	20"	21"	22"	24"			12"	14"	16"	17"	18"	19"	20"	21"	22"	24"		
Overall Height	14"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	15"	ST	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
	16"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-	-	-	-	-
	17"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-	-	-	-	-
	18"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	-	-	-	-	-	-	-	-	-	-	-
	19"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	ST	ST	ST	S	ST	S	S
	20"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	S	ST	ST	ST	ST	ST	ST	ST	ST	S	ST	S	S
	21"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	ST	ST	ST	ST	ST	ST	ST	S	ST	S	S
	22"	-	-	SQT	ST	SQT	ST	SQT	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	ST	SQT	S	ST	S	S
	23"	-	-	ST	ST	ST	ST	ST	S	ST	ST	-	-	SQT	ST	SQT	ST	SQT	ST	SQT	S	ST	S	S
	24"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	ST	SQT	S	ST	S	S
	25"	-	-	S	S	ST	S	ST	S	S	S	-	-	SQT	ST	SQT	ST	SQT	ST	SQT	S	ST	S	S
	26"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	ST	ST	S	ST	S	S
	27"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	ST	S	S	S	S
	28"	-	-	-	-	-	-	-	-	-	-	-	-	S	S	ST	S	ST	S	ST	S	S	S	S

**S** Sport Back      **ST** Sport or TRU-Comfort Back

**SQT** Sport, Q-Back, or TRU-Comfort

Sport, Q-Back, and TRU-Comfort Back Availability Chart

Set Overall Height at: \_\_\_\_\_

### Aftermarket Back Interface Preset Width

Set Interface Width at: \_\_\_\_\_ (Range 16"-24" Wide)

### TRU-Balance® 3 Back Shroud

**Required.** The back shroud is a feature of the TRU-Balance® 3 seating system and is available in matching shroud colors. The back shroud is not available with the back cane configuration. Reference colors shown on page 2.

Back in Black

Electric Blue

Pretty in Pink

Lime Light

Silver Bullet

Orange Crush

Red Corvette

White Satin

Deep Purple

Grey Street

Midnight Blue

Yellow Blaze

### Additional TRU-Balance® 3 Back Options

**Optional.**

TRU-Comfort Lateral Positioning Wedge, pair .....\$30

Available with TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Privacy Flap, 9.75" .....\$35

Standard on recline seating systems.

Privacy Flap, 12.75" .....\$35

Standard on recline seating systems.

See Comfort Company Back order form

## ARMRESTS

**Required.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

### TRU-Balance® 3 Flip-Back Armrests

Choose a left and right TRU-Balance® 3 armrest and set height and optional downpost if desired. Different Left (L) and Right (R) armrest height packages are allowed. Not available with Back Cane Option.

**L** **R**

TRU-Balance® 3 Flip-Back, 8"-12"H..... Standard

**L** **R**

TRU-Balance® 3 RIGHT Flip-Back, 12"-16"H ..... Standard

Set LEFT height at: \_\_\_\_\_ Set RIGHT height at: \_\_\_\_\_

TRU-Balance® 3 Armrest Downpost (Pair) ..... \$175

Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available on 14" Width seat or back.

### Two-Post Armrests

Select the armrest and set the left and right heights. These armrests may interfere with offset backs/ back canes.

2-Post, Flip Back, Height Adjustable, ..... No Charge  
and Removable Armrests (9"-13" Height)

Not available with a depth less than 16". Required with flip forward back and laptrays.

Set LEFT height at: \_\_\_\_\_ Set RIGHT height at: \_\_\_\_\_

## ARMRESTS CONTINUED

### Heavy Duty, Removable Armrests

Select a type of HD removable armrests and the set height. Only one height package is selectable as a normal configuration.

8"-12"H, Single Post, Height Adjustable, .....\$75 ea. Removable, Heavy Duty (Sold as pair) HCPCS: E0973	10"-14"H, Single Post, Height Adjustable, .....\$75 ea. Removable, Heavy Duty (Sold as pair) HCPCS: E0973
---	--

Set LEFT height at: \_\_\_\_\_ Set RIGHT height at: \_\_\_\_\_

### ARMREST SIZE

**Required.** Select a Left (L) and Right (R) armrest size. This also affects the armpad selection below.

<input type="checkbox"/> L <input type="checkbox"/> R 14" Full Length	<input type="checkbox"/> L <input type="checkbox"/> R 10" Desk Length
---	---

### ARMPADS

**Required.** Select to omit an armpad selection or select a left and right armpad length and a type of left and right armpad. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or ..... No Charge  
See Positioning Components Page or Order Form

### Foam Armpads

Select a left and right standard type armpad (mix or match) OR choose the gel armpads that follow.

<input type="checkbox"/> L <input type="checkbox"/> R Straight, Armpad ..... Standard	<input type="checkbox"/> L <input type="checkbox"/> R Waterfall, Armpad..... No Charge
---	--

### Falcon Gel Armpads

Sold only in pairs. Select an armpad set below. Armpads are 3" wide. The mounting plate allows the armpad to be offset from the armrest tube.

Falcon Gel Left and Right Armpad: Full Length .....\$92.50 ea Desk Length .....\$70.00 ea For use with TB3 standard armrests.	Falcon Gel Left and Right Armpad with Mounting Plate: Full Length ..... \$148.02 ea Desk Length ..... \$126.20 ea For use with 2-Post, Heavy Duty Removable, or Pediatric Armrests.
--	--

### JOYSTICK BRACKET RECEIVER

**Required.**

Standard Fixed Receiver ..... Standard
Stealth Height Adjustable Receiver..... \$200 Tube Arm Part: ST-ARM250C-Q HCPCS: K0108 Track Arm Part: ST-ARM250T-Q HCPCS: K0108



### HEADRESTS


**Optional.** Select a headrest pad and bracket.

#### Headrest Pads

Stealth 8" Comfort Plus (CP450) .....\$142 Part: POS148484 HCPCS: E0955	Stealth Small Contoured (9.75"Wx3"H)..... \$60 (STL-OTSM250) Part: SETHEAD1052 HCPCS: E0955
Stealth 10" Comfort Plus (CP250)(A) .....\$142 Part: SETHEAD1057 HCPCS: E0955	Stealth Large Contoured (11"Wx6"H) (B) ..... \$75 (STL-OTLG150) Part: SETHEAD1053 HCPCS: E0955
Stealth 14" Comfort Plus (CP150) .....\$142 Part: INDPART2327 HCPCS: E0955	
Pediatric Comfort Plus (6"Wx4.5"H) .....\$142 (CPN350) Part: SETHEAD1069 HCPCS: E0955	

#### Headrest Bracket

Stealth's Worlds Best Removable.....\$210 Headrest Bracket (TWB485-1) (C) Part: FRMASMB15999 HCPCS: E1028	Stealth Multi-Axis Removable Headrest .....\$150 Bracket (HMO475P) (D) Part: HDWASMB1084 HCPCS: E1028
---	---



## LAP BELTS

**Required.** Select a lap belt below.

### Standard Lap Belts

38" Lap Belt	50" Lap Belt	60" Lap Belt	70" Lap Belt	80" Lap Belt	90" Lap Belt
--------------	--------------	--------------	--------------	--------------	--------------

### Optional Lap Belts

i-Fit 60" Single Pull, Push Button, Medium .....\$86 (7.5") Padded, 1.5" Belt (IFPB5213)	i-Fit 60" Single Pull, Push Button, Extra Large .....\$86 Padded (11.5"), 2" Belt (IFPB5315)
i-Fit 60" Single Pull, Push Button, Large .....\$86 Padded (9.5"), 2" Belt (IFPB5314)	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Required with ICO 7176-19 Compliant Transit Option.

## MANUFACTURER INSTALLED ACCESSORIES

**Optional.**

Trilogy Vent Tray ..... \$400  
Part: ACC136791 HCPCS: E1029

Compact Vent Tray ..... \$1,395  
Part: ACC133889 HCPCS: E1029

Push Handles (pair) .....\$175  
Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Accessory Bar .....\$100  
Standard with the SCIM module, Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.



Trilogy Vent Tray



Compact Vent Tray

## POWER LEGREST CONTROL

**Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)

Articulating Foot Platform (Requires One Actuator Control)

Combined Legrests (Requires 1 actuator control)

## POWER POSITIONING ELECTRONICS

**Required, except with Static Seat without power options.** Select the desired control type for the seating system. Please note the available actuator controls. The 4-way toggle and 4 function push button are mounted on the same side as the joystick.

Single Actuator Control, Through Toggle ..... No Charge  
Only available with Power Tilt only. Manual recline is not available with NE Tilt through toggle. Fender lights are not available with tilt through toggle. Select toggle side.

Left Side Toggle

Right Side Toggle

Through NE+ 6-Key Joystick with either:

Single Actuator Control.....\$1,895  
Quote will select this price/code with 1 actuator. HCPCS: E2310

Two Actuator Control .....\$2,895  
Quote will select this price/code with 2 actuators. HCPCS: E2311

Single Actuator Through Q-Logic 2 Controller .....\$1,895  
and Interface (AM1) HCPCS: E2310

Multiple Actuator Through Q-Logic 2 Controller.....\$2,895  
and Interface (AAM) HCPCS: E2311

Multiple Actuator Through 4-Way Toggle, .....\$2,895  
Q-Logic 2 Controller, and Interface (AAM) HCPCS: E2311  
Bracket mounts to outside of the armpad.

Multiple Actuators Through 4-Way Button Toggle, .....\$2,895  
Q-Logic 2 Controller, and Interface (AAM) HCPCS: E2311  
Bracket mounts to outside of the armpad.

Through 4 Function Push Button and Q-Logic 2 Controller and Interface with either:

Single Actuator Control (AAM) .....\$2,895  
Quote will select this price/code with 1 actuator. HCPCS: E2310

Multiple Actuator Control (AAM) .....\$2,895  
Quote will select this with multiple actuators. HCPCS: E2311  
Only available with TB3 armrests.

## Attendant Control

**Optional.** Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control ..... \$500  
Part: ELE125264 HCPCS: E2331

## Seat Mounted, LED Light Package

**Optional.** Available only with NE+ 8 Key Joystick or Q-Logic 2 Ex Joystick for LED Lights.

Seat Mounted, LED Light Package..... \$995



Seat Mounted LED Lights (One Shown)



## LEGRESTS

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

### Center Mount Foot Platforms

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section.

Power Articulating Foot Platform ..... \$2,995  
Please verify the appropriate billing code to be used for this item with all third-party payors.

Manually Adjustable Center Mount ..... \$250  
Foot Platform



### Center Mount Lower Extension

8"-12" Extension

N/A with Power Articulation.

10"-16" Extension

13"-19" Extension

Set Length at: \_\_\_\_\_

### Center Mount Footplates

5"x6" Footplates

5"x8" Footplates

5"x11" Footplates

Center Mount Foot Platform Compatibility Matrix	TRU-Balance® 3 Power Recline			TRU-Balance® 3 Power Tilt			TRU-Balance® 3 Power Lift			
	Low	Medium	High	Low	Medium	High	Low	Medium	High	
10"-16" Length Extension	!	!	✓	✓	✓	✓	✓	✓	✓	✓ No interference issues at this seat-to-floor height.
13"-19" Length Extension	!	!	!	!	!	!	!	!	!	! May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

### Center Mount Calf Support

**Optional.** Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pad (pair) 6"Wx6"H ..... \$150

Small Calf Panel 18"Wx8"H ..... \$150

Large Calf Pad (pair) 8" taper to 6"Wx8"H ..... \$150

Large Calf Panel 24"Wx8"H ..... \$200

Only available on widths of 20" or greater.

### Modular Footplate Extension Components

**Optional.** Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension, 8"D and 11"D Footplates ..... \$50

4" Outside Edge ..... \$50

Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension ..... \$50

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

6" Outside Edge ..... \$50

Available with 8"D and 11"D footplates and 11"D footplate extension.

### Additional Accessories

Reversed Dartex Calf Pad or Panel Cover (pair) ..... \$120

Footplate Connector ..... \$15

Wheel Bumpers (pair) ..... \$25

### 70° Power Articulating, Swing Away Legrests

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper, lower extension, and footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

### Configuration Verification

Left and Right Legrests will Match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

### Power Articulating, Swing Away Legrest

☐ L ☐ R 70° Articulating, Swing-Away Legrest ..... \$1,197.50 ea.  
HCPCS: E1010

### Lower Extension

☐ L ☐ R 11"-13" Length

☐ L ☐ R 14.5"-18" Length

Set Length at: \_\_\_\_\_

## LEGRESTS CONTINUED

### Angle Adjustable Footplate

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ 4"Wx6"D    ☐ ☐ 5"Wx6"D    ☐ ☐ 5"Wx8"D    ☐ ☐ 6"Wx8"D<sup>1</sup>    ☐ ☐ 6"Wx11"D<sup>1</sup>    ☐ ☐ 8"Wx11"D<sup>2</sup>  
 L R                      L R                      L R                      L R                      L R                      L R

### Accessories

Calf Pad (pair) 6"Wx6"H ..... No Charge

Gel Padded Knee Abductor Buttons (pair) ..... \$125  
 INDPART2242 HCPCS: K0108

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	TRU-Balance® 3 Power Recline			TRU-Balance® 3 Power Tilt			TRU-Balance® 3 Power Lift			
	Low	Medium	High	Low	Medium	High	Low	Medium	High	
Angle Adjustable, 11" - 13"	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓ This rigging configuration causes no interference issues at this seat-to-floor height.
Angle Adjustable, 14.5" - 18"	!	!	!	!	!	✓	!	!	✓	! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 17.5" - 21"	X	X	X	X	X	!	X	!	!	X This rigging configuration does not work with this seat-to-floor height.

### 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

### Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

### Swing-Away Legrest Upper

☐ ☐ 60°-90° Swing-Away Legrest ..... No Charge  
 L R Lower extensions with minimum 11" length utilize a short style legrest upper.

### Lower Extension

☐ ☐ 11"-13" Length    ☐ ☐ 11"-16.5" Length    ☐ ☐ 13"-16.5" Length    ☐ ☐ 15.5"-19.5" Length  
 L R N/A with Multi-axis    L R Only with Multi-Axis    L R N/A with Multi-Axis    L R

Set Length at: \_\_\_\_\_

### Extruded Footplate(s)

☐ ☐ 7.25"Wx6"D Extruded Footplate ..... Standard    ☐ ☐ 8.5"Wx6"D Extruded Footplate ..... Standard  
 L R May cause interference on seat widths less than 18."

### Angle Adjustable or Multi Axis Selection

☐ ☐ Angle Adjustable with Heel Loops ..... \$75 ea.  
 L R (Select footplate size) HCPCS: K0040    ☐ ☐ Multi-Axis (Select Footplate Size) ..... \$247.50 ea  
 L R HCPCS: K0040

### Angle Adjustable or Multi-Axis Footplate Size

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ 4"Wx6"D    ☐ ☐ 5"Wx6"D    ☐ ☐ 5"Wx8"D    ☐ ☐ 6"Wx8"D<sup>1</sup>    ☐ ☐ 6"Wx11"D<sup>1</sup>    ☐ ☐ 8"Wx11"D<sup>2</sup>  
 L R                      L R                      L R                      L R                      L R                      L R

### 60°-90° Swing-Away Legrest Accessories

Heel Loops (pair) ..... \$38  
 Part: FRMASMB7873 HCPCS: E0951  
 Standard on angle adjustable and multi-axis footplates.

Wheel Bumpers (pair) ..... \$25

Gel Neoprene Footrest Sleeve (pair) ..... \$80  
 Part: INDPART2243 HCPCS: K0108

Therafin 31515 Multi-Axis Amputee Support, LEFT ..... \$230  
 Part: INDPART2822 HCPCS: E1020

Therafin 31515 Multi-Axis Amputee Support, RIGHT ..... \$230  
 Part: INDPART2822 HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover, LEFT ..... \$150  
 Part: INDPART2823 HCPCS: K0108

Therafin 32387 Amputee Support Gel Cover, RIGHT .... \$150  
 Part: INDPART2823 HCPCS: K0108

## LEGRESTS CONTINUED

60°-90° Swing-Away Legrest Compatibility Matrix	TRU-Balance® 3 Power Recline			TRU-Balance® 3 Power Tilt			TRU-Balance® 3 Power Lift			
	Low	Medium	High	Low	Medium	High	Low	Medium	High	
Extruded, 13" - 16.5"	!	!	✓	!	✓	✓	✓	✓	✓	✓ This rigging configuration causes no interference issues at this seat-to-floor height.
Extruded, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!	! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 13" - 16.5"	!	!	✓	!	!	!	!	!	!	
Angle Adjustable, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!	
Multi-Axis, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!	X This rigging configuration does not work with this seat-to-floor height.

## Elevating Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

### Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

### Elevating Legrest

☐ Elevating Legrest Assembly ..... \$162.50 ea.  
☐ ☐ Select lower extension and footplates HCPCS: E0990

### Lower Extension

☐ ☐ 11"-14" Length ☐ ☐ 12"-14.5" Length ☐ ☐ 14"-17" Length ☐ ☐ 14.5"-18" Length  
☐ ☐ For reinforced only ☐ ☐ For angle adjustable only ☐ ☐ For reinforced only ☐ ☐ For angle adjustable only

Set Length at: \_\_\_\_\_

### Reinforced Cast Footplate(s)

☐ ☐ 7.25"Wx6"D Reinforced Footplate ..... Standard ☐ ☐ 8.5"Wx6"D Reinforced Footplate ..... Standard  
☐ ☐ May cause interference on seat widths less than 18."

### Angle Adjustable Footplate(s)

☐ ☐ Angle Adjustable with Heel Loops ..... \$75 ea.  
☐ ☐ (Select footplate size below) HCPCS: K0040

### Angle Adjustable Size

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ 4"Wx6"D ☐ ☐ 5"Wx6"D ☐ ☐ 5"Wx8"D ☐ ☐ 6"Wx8"D<sup>1</sup> ☐ ☐ 6"Wx11"D<sup>1</sup> ☐ ☐ 8"Wx11"D<sup>2</sup>

### Accessories

Heel Loops (pair) ..... \$38  
Part: FRMASMB7873 HCPCS: E0951  
Standard on angle adjustable footplates.

Wheel Bumpers (pair) ..... \$25  
Falcon LPW Gel Padded Adductor Buttons (pair) ..... \$125  
Part: INDPART2243 HCPCS: K0108

## CUSHION QUICK PICK

**Optional.** General size selections are shown below. Refer to the [specific cushion order form](#) for detailed pricing, sizes, codes, and part numbers. Check the non-matching box and attached a cushion order form for a size other than the previously selected seat width and depth, or for a cushion not shown here.

Match seat selection width and depth

Non-matching cushion size, [see cushion order form](#)

### Simplicity Cushion (General Use)

16"Wx16"D to 20"Wx20"D ..... \$98 HCPCS: E2601  
22"Wx16"D to 24"Wx24"D ..... \$98 HCPCS: E2602

### Solution® Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D ..... \$380 HCPCS: E2607  
22"Wx16"D to 24"Wx24"D ..... \$430 HCPCS: E2608

## CUSHION QUICK PICK CONTINUED

### Solution® 1 Cushion (Skin Protection)

14"Wx14"D to 20"Wx24"D.....	\$380	22"Wx16"D to 24"Wx24"D.....	\$430
	HCPCS: E2603		HCPCS: E2604

### Spectrum® Gel Cushion (Positioning)

12"Wx10"D to 20"Wx24"D.....	\$480	22"Wx16"D to 24"Wx24"D.....	\$530
	HCPCS: E2605		HCPCS: E2606

### Spectrum® Foam Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D.....	\$480	22"Wx16"D to 24"Wx24"D.....	\$530
	HCPCS: E2607		HCPCS: E2608

### TRU-Comfort™ 2 4-Way Stretch Cushion (Skin Protection)

16"Wx14"D to 20"Wx24"D.....	\$380	22"Wx16"D to 22"Wx24"D.....	\$480
	HCPCS: E2603		HCPCS: E2604

### TRU-Comfort™ Sport Seat Cover


Available as an additional cover option for both TRU-Comfort seat cushions. Only available in sizes below.

16"Wx16"D.....	\$100	20"Wx20"D.....	\$100
18"Wx18"D.....	\$100		

## PROVIDER INSTALLED ACCESSORIES

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU-Balance® 3 Cup Holder .....	No Charge	TRU-Balance® 3 Front Mounted.....	\$50
For use with TRU-Balance 3 (track style) armrests.		Cell Phone Holder (E)	
Cup Holder (Tube-style Armrest) .....	No Charge	Large, Black, Weather Cover .....	\$158
For use with tube-style armrests.		TRU-Balance® 3 Reflector Kit .....	\$25
TRU-Balance® 3 Backpack Holder (Pair) (A)...	\$25	TRU-Balance® 3 Personal Item Hook .....	\$25
Not available with back cane option.		Not available with back cane option.	
TRU-Balance® 3 Hydration System (B).....	\$175	TRU-Balance® 3 Rear View Mirror .....	\$25
TRU-Balance® 3 Glove Box (C) .....	\$100	TRU-Balance® 3 Oxygen Holder .....	\$150
TRU-Balance® 3 Transfer Bar, .....	\$175		HCPCS: E2208
Bolt-On (Pair) .....	HCPCS: E0705	TRU-Balance® 3 Tablet Holder .....	\$375
TRU-Balance® 3 Transfer Bar, .....	\$175	TRU-Balance® 3 Clothing Guard (Pair).....	\$175
Quick Release (Pair) (D) .....	HCPCS: E0705		



## BATTERIES

**Required.** Select a type of battery below, or to omit the batteries. These batteries are Group NF-22, maintenance free. The powerchair requires 2 batteries.

NF-22, 50 Amp AGM, UPG Battery .....	\$295 ea.	NF-22, 50 Amp AGM, MK Battery.....	\$295 ea.
Part: BATLIQ1018 .....	HCPCS: E2361	Part: BAT121346 .....	HCPCS: E2361
NF-22, 50 Amp Gel Cell, UPG Battery .....	\$295 ea.	NF-22, 50 Amp Gel Cell, MK Battery.....	\$295 ea.
Part: BATGEL1003 .....	HCPCS: E2361	Part: BATGEL1008 .....	HCPCS: E2361

No batteries. Provider will assume responsibility of batteries.

### Battery Installation

Please choose if batteries will be installed, if selected above.

Install Batteries .....	No Charge	Do NOT Install Batteries.....	No Charge
-------------------------	-----------	-------------------------------	-----------

## MANUFACTURER INSTALLED POSITIONING COMPONENTS

**Optional.** More options are available on the [Positioning Components Order Form](#).

### Stealth Thoracic Laterals

Swing-away laterals may interfere with TB3 and 2-post, flip up armrests. These laterals are not available with the back cane option.

#### Bracket

Swing-away, Rail Mounts, pair (TWBSTL)..... \$450  
Part: INDPART2896 HCPCS: E1028



#### Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$200  
Part: INDPART2663 HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$200  
Part: INDPART2664 HCPCS: E0956

#### Installation

Install Laterals

Do Not Install

### Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

#### Bracket

Flip Down Bracket (TWBTLADD-HW)..... \$225 ea.  
— Part: INDPART2804 HCPCS: E1028

Removable Bracket (LP430)..... \$250 ea.  
— Part: INDPART2805 HCPCS: E1028

Flip Down with LEFT Full Surface Contact ..... \$225 ea.  
Feature (TWBADD-L & LPHW407)  
Part: POS140428 HCPCS: E1028

Removable Bracket with 1" Offset ..... \$325 ea.  
— (LP430 & LHW-121)  
Part: INDPART2854 HCPCS: E1028

Flip Down with RIGHT Full Surface Contact ..... \$225 ea.  
Feature (TWBADD-R & LPHW407)  
Part: POS140429 HCPCS: E1028

#### Adapter

2" Height Adapter Plate (LHW-122)..... \$25  
Part: INDPART2821 HCPCS: K0108

2 1/4" Adapter Plate (LHW-129) ..... \$50  
Part: POS140430 HCPCS: K0108

Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

Only available with the flip down full surface bracket.

#### Pad

Short Pad, 4"x4" (TWBADD-4X4) ..... \$50 ea.  
— Part: INDPART2801 HCPCS: E0956

Long Pad, 4"x12" (TWBADD-4X12) ..... \$50 ea.  
— Part: INDPART2803 HCPCS: E0956

Medium Pad, 4"x8" (TWBADD-4X8) ..... \$50 ea.  
— Part: INDPART2802 HCPCS: E0956

#### Installation

Install on Left

Install on Right

Do Not Install

## NON-INSTALLED POSITIONING COMPONENTS

**Optional.** More options are available on the [Positioning Components Order Form](#).

### Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$475  
(SPLE107) (pair)  
Part: INDPART2786 HCPCS: K0108

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$475  
(SPLE109) (pair)  
Part: INDPART2788 HCPCS: K0108

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$475  
(SPLE108) (pair)  
Part: INDPART2787 HCPCS: K0108

### Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D ..... \$266  
(SPLE111) (pair)  
Part: INDPART2783 HCPCS: K0108

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D ..... \$266  
(SPLE114) (pair)  
Part: INDPART2785 HCPCS: K0108

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D ..... \$266  
(SPLE112) (pair)  
Part: INDPART2784 HCPCS: K0108



## NON-INSTALLED POSITIONING COMPONENTS CONTINUED

### Stealth Elbow Support

- Left, Swing-Away, Elbow Support Bracket..... \$175 ea.  
 — Tube Style Arms (ARMS-ES-L)  
*Part: POS124064 HCPCS: E1028*
- Left, Fixed Elbow Support Bracket..... \$175 ea.  
 — TRU-Balance 3 Arms (ARMS-ESFTB3-L)  
*Part: POS133465*

- Right, Swing-Away, Elbow Support Bracket..... \$175 ea.  
 — Tube Style Arms (ARMS-ES-R)  
*Part: POS124066 HCPCS: E1028*
- Right, Fixed Elbow Support Bracket ..... \$175 ea.  
 — TRU-Balance 3 Arms (ARMS-ESFTB3-R)  
*Part: POS133466*

### Pad

- 3"Wx3"H Elbow Pad (ARMS-ES-3X3) .....\$50 ea.  
 — *Part: POS124065 HCPCS: K0108*

- 4"Wx4"H Elbow Pad (ARMS-ES-4X4) .....\$50 ea.  
 — *Part: POS124067 HCPCS: K0108*

### Stealth i-Fit Dynamic Chest Straps

- Upper Classic, Medium (IFPS1103) ..... \$142  
*Part: ST-IFPS1103-Q HCPCS: E0960*
- Upper Classic, Large (IFPS1104) ..... \$142  
*Part: ST-IFPS1104-Q HCPCS: E0960*
- Upper Contour, Medium (IFPS2103) ..... \$142  
*Part: ST-IFPS2103-Q HCPCS: E0960*

- Upper Contour, Large (IFPS2104) ..... \$142  
*Part: ST-IFPS2104-Q HCPCS: E0960*
- Chest Strap, Medium, Belt Anchor (IFCS4113).....\$90
- Chest Strap, Large, Belt Anchor (IFCS4114).....\$90

### Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.


- Facial Lateral, 4.5"x2" Pad, Right Side, ..... \$325  
 for Comfort Series (CPS9-R-HW & 900-Q)  
*Part: INDPART2647 HCPCS: E1028*
- Facial Lateral, 4.5"x2" Pad, Left Side ..... \$325  
 for Comfort Series (CPS9-L-HW & 900-Q)  
*Part: INDPART2648 HCPCS: E1028*

- Facial Lateral, 4.5"x2" Pads, Left and Right Side ..... \$425  
 for Comfort Series (CPS9-B-HW & 900-Q)  
*Part: INDPART2649 HCPCS: E1028*

### Other Accessories

- Egg Switch Mount (ESM) (A) ..... \$50 ea.  
 — *Part: INDPART2815*
- Figure 8 Bracket with 13" Rod and..... \$250  
 Egg Switch Mount (B)  
*Part: ST-F8678-13-Q HCPCS: K0108*

- Facial Pad, 6"x2" (920-Q) ..... \$50  
*Part: ST-920-Q HCPCS: K0108*
- Swing-Away Joystick Mount, 18" Rod (C) .... \$250  
*Part: ST-SUS18-Q HCPCS: K0108*  
 Allows mounting of facial pads and attaches to the headrest



## SPECIALTY CONTROLS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

### Enhanced Display

- Q-Logic 2 Enhanced Display Kit ..... \$1,884.93  
 Select a color option below. Includes enhanced display, flexible goose-neck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



### Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

- Black, Soft Touch ..... Standard      Match Joystick Shroud to Base Color..... \$50  
*Part: See color below.*
- Non-matching Joystick Shroud ..... \$50  
 Select color below
- |               |                |               |               |
|---------------|----------------|---------------|---------------|
| Red Corvette  | White Satin    | Lime Light    | Midnight Blue |
| Electric Blue | Pretty in Pink | Grey Street   | Orange Crush  |
| Back in Black | Deep Purple    | Silver Bullet | Yellow Blaze  |

## SPECIALTY CONTROLS CONTINUED

### Head Array

Adult i-Drive Tri Array Series, ..... \$4,175  
3-Switch with Egg Switch for Mode (IDH200-1)  
Part: CTL135940 HCPCS: E2330  
Features 3 embedded proximity sensors for drive control and an egg switch for mode change.



### Other

See Specialty Controls Order Form

Multiplier Harness .....\$221.21

Red Egg Switch (ESRED) ..... \$77

Blue Egg Switch (ESBLU) ..... \$77


Black Egg Switch (ESBLK) ..... \$77

Green Egg Switch (ESGRN)..... \$77

Yellow Egg Switch (ESYEL) ..... \$77

SCIM Module Kit ..... \$575

Allows use of joystick display for specialty devices. Available with Q-Logic 2 EX Joystick selection above. Includes the TRU-Balance® 3 accessory bar. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



### Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.