



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® HD Order Form with TRU-Balance® 3 Positioning - Group 3, ODJFS Only

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Selections with HCPCS codes may have a part number noted with the selection. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.

If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



BASE MODEL

REQUIRED. Select a power base. 450 lb. weight capacity with standard 4.7 mph programmable motors, 8 amp off board charger, and unoccupied transit loops. Dimensions of base without legrests: 26"Wx36.75"L.

Quantum Q6 Edge HD 3HD-SS\$8,595

Part: Q6 EDGEHD 3HD-SS Mfg Sugg HCPCS: K0850/K0014

Solid seat pan standard with no power option.

Quantum Q6 Edge HD 3MPHD-SS\$8,595

Part: Q6 EDGEHD 3MPHD-SS Mfg Sugg HCPCS: K0862/K0014

Multiple power option.

Quantum Q6 Edge HD 3SPHD-SS\$8,595

Part: Q6 EDGEHD 3SPHD-SS Mfg Sugg HCPCS: K0858/K0014

Single power option.

BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



Red Corvette



Electric Blue



Back in Black



White Satin



Pretty in Pink



Deep Purple



Lime Light



Grey Street



Silver Bullet



Midnight Blue



Orange Crush



Yellow Blaze



To be
determined
at a later
date.

BASE MOUNTED FOOT PLATFORM

REQUIRED. Select a foot platform option. May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard

9.45"W x 10"D platform.

Omit base mounted foot platform No Charge

Base Mounted, Bell Shaped Foot Platform No Charge

10.12"W tapering out to 14"W x 13.06"D platform.

DRIVE WHEELS

REQUIRED. Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

3" x 14" Flat Free Tires Standard

4"W Flat Free, Knobby, Extended Traction Tires \$551

3" x 14" Pneumatic Tires No Charge

4"W Pneumatic, Knobby, Extended \$551
Traction Tires

Rim Color

Rim colors other than the standard are available through a special order and extended lead time.

Standard Rims Standard

Match Rims to Base Color\$50

Non-matching Rims.....\$50

Select color below

Red Corvette

Back in Black

Pretty in Pink

Lime Light

Midnight Blue

Yellow Blaze

Electric Blue

White Satin

Deep Purple

Grey Street

Orange Crush

BASE LIGHTING PACKAGES

Optional. Select one type of optional lighting packages below. Lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights \$250

Fender lights that operate only by pressing included single switch button.

Automatic Fender Lights..... \$275

Fender lights that power on when the controller is activated and can be operated by the included single switch button.



Fender Lights

ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Forms](#).

Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE 4-Key Joystick Standard
Available with Static Seat or Tilt Only through Toggle.

NE+ 6-Key Joystick (1 Actuator)..... No Charge
Control of one actuator through Joystick on single power base model (3SPHD-SS) only.

No Joystick, Setup for NE+ Joystick... No Charge

NE+ 6-Key Joystick (2 Actuators) \$200
Control up to 2 actuators through joystick.

NE+ 8-Key Joystick (2 Actuators) \$400
Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module \$695
HCPCS: E2377

Q-Logic 2 Expandable Controller..... \$425
Harness
Part: DWR1111H065 HCPCS: E2313



Q-Logic 2 EX
Controller Module

Q-Logic 2 Joystick

Q-Logic 2 EX Joystick..... \$700

Q-Logic 2 Ex Joystick for Seat..... \$700
Mounted, LED Light Package

No Joystick, Setup for Q-Logic

Stand Alone Joystick with \$700
External Charger Port

Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.



Stand Alone Joystick

Q-Logic 2 Joystick Shroud

REQUIRED with Q-Logic 2 EX Joystick selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Match Joystick Shroud to Base Color..... \$50

Non-matching Joystick Shroud \$50

Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze



Q-Logic 2 EX Joystick
with Lime Light

JOYSTICK HANDLES

Optional. The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$71
Part: ACCASMB1860 HCPCS: E2323

Soft Ball Handle (F)..... \$50
Part: ACCASMB1856 HCPCS: E2323

Chin Cup (C) \$71
Part: ACCASMB1861 HCPCS: E2323
(E2324 when chin control)

Extended I-Handle (G)..... \$71
Part: ACCASMB1857 HCPCS: E2323

T-Bar (D) \$107
Part: ACCASMB1858 HCPCS: E2323

Bodypoint 3" Goal Post Handle..... \$69
Part: PLSKNCP1057 HCPCS: E2323

Medium Ball Handle (E)..... \$50
Part: ACCASMB1859 HCPCS: E2323

Bodypoint 4" Goal Post Handle..... \$69
Part: PLSKNCP1058 HCPCS: E2323



Joystick Handles

JOYSTICK MOUNTING BRACKETS

REQUIRED. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket or See Positioning Positioning Components Page or [Order Form](#)

Fixed, Inline Joystick Bracket

Fixed, Inline, Left Side Mount Standard

Fixed, Inline, Right Side Mount Standard

Swing-Away (Retractable) Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount..... \$245
NE/NE+ Part: ACC142526 HCPCS: E1028
Q-Logic: ACC142524

Swing-away, 1" Offset, Left Side Mount \$245
NE/NE+ Part: ACC142747 HCPCS: E1028
Q-Logic: ACC142742

Swing-away, Inline, Right Side Mount \$245
NE/NE+ Part: ACC142527 HCPCS: E1028
Q-Logic: ACC142525

Swing-away, 1" Offset, Right Side Mount.... \$245
NE/NE+ Part: ACC142748 HCPCS: E1028
Q-Logic: ACC142743

Retract4 Swing-Away Arm Kit

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Retract4, Left Side Mount \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860

Retract4, Right Side Mount..... \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860

Multi-Axis

Allows height adjustment of joystick.

Multi-axis, Left Side Mount.....\$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216

Multi-axis, Right Side Mount\$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216

Stealth Swing Away/Flip Down Bracket

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket \$275
Left Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028

Stealth Swing-Away/Flip Down Bracket \$275
Right Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028

ELECTRONICS ACCESSORIES

Optional.

USB Mobile Device Charger\$175
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

Multiplier Harness\$221.21
Q-Logic EX/NE Programmer \$545
Q-Logic PC Programming Disk with Cable ... \$750

TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™

REQUIRED. Select one seating system. Be sure to complete the Seat Dimension, Seat to Floor, and back selections.

TRU-Balance® 3 Configuration

TRU-Balance 3 Static Seat\$2,295
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Heavy Duty, Power Tilt \$7,995
Part: PTOINDV3643 HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Heavy Duty, Power Recline\$8,595
Part: PTOINDV3645 HCPCS: E1004
175° of power recline.

TRU-Balance 3 Heavy Duty, Power Tilt & Recline \$12,895
Part: PTOINDV3644 HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.





10" Power Lift - iLevel System

Optional. Select one option below.

TRU-Balance 3 Heavy Duty, 10", Power Lift.....\$2,495
with iLevel®
Part: PTOINDV3797 HCPCS: E2300
400 lb weight capacity. 24"W x 24"D maximum seat size.

EXTRA
STABILITY TECHNOLOGY™

QUANTUM®
iLevel



TRU-BALANCE® 3 SEATING CONTINUED

Seat Dimensions

REQUIRED with seat selection. Select width and depth package and enter desired set width and depth in that available package range, available in 1" increments. Configurations disallowed or not shown may be available through special order. Maximum non-lift seat size is 32"W x 25"D. Maximum lift seat size is 24"Wx25"D. Lift seat depth of 25" D is attained by using the 20"-24" Depth Package.

Seat Width

20"-24" Width Package

24"-28" Width Package

28"-32" Width Package

Set Width at: _____ (From 20" to 32", Max Lift Seat Width = 24")

Seat Depth

16"-20" Depth Package
Minimum Seat Depth is 20"

20"-24" Depth Package

24"-28" Depth Package
Not available with Power Lift.
Maximum Seat Depth is 25"

Set Depth at: _____ (From 20" to 25")

Seat to Floor Height

REQUIRED. Select a seat to floor height compatible with the seat functions previously selected.

Static Seat Only or Recline Only

16.5" Seat to Floor

17.5" Seat to Floor

18.5" Seat to Floor

Tilt Only, Tilt and Recline, or Any Power Lift Configuration

18.5" Seat to Floor

19.5" Seat to Floor

20.5" Seat to Floor

Back Selection

REQUIRED. Match choice in this section with the options in the back cushion configuration section.

Static Back Standard
Select a back cushion configuration in the subsequent section.

Back Cane Option, Omit Back Pan/ Cushion (-\$125)
Available only with Static Back and Power Tilt (1 actuator) selections.
Select cane type and the back cane height below.

Back Cane

Select when Back Cane Option, Omit Back Pan/ Cushion is selected above.

Angle Adjustable Standard

1" Endo \$35

Flip Forward \$25

2" Endo \$35

Back Cane Height

Select when Back Cane Option, Omit Back Pan/ Cushion is selected above.

16"

17"

20"

21"

24"

25"

Back Cushion Configuration

REQUIRED except when back cane option selected. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with Special Order and extended lead time.

Sport Back..... No Charge
Complete following width package and size.

Omit Back Pan/Cushion..... No Charge
Select a following width package.

TRU-Comfort 4-way Stretch, 20" Width or less..... \$595
Complete following width package and size. HCPCS: E2620

Aftermarket Back Interface \$200
Enter the preset back width below. Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount aftermarket backs. A Special Order is required on back widths less than 16".

TRU-Comfort 4-way Stretch, 22" Width or greater \$695
Complete following width package and size.. HCPCS: E2621

TRU-BALANCE® 3 SEATING CONTINUED

Back Cushion Width

Available with +/-4" of seat width. Reference chart for available back width.

Match Seat Width Selections

16"W-20"W Package

20"W-24"W Package

24"W-28"W Package

28"W-32"W Package

Set Width at: _____

Back Cushion Overall Height

Static Seat and Power Tilt available height range: 16"-25"H. *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.*

Power Recline available height range: 19"-28"H. *Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.*

Set Overall Height at: _____

Aftermarket Back Interface Preset Width

Set Interface Width at: _____ (Range 16"-24" Wide)

TRU-Balance® 3 Back Shroud

REQUIRED. The back shroud is a feature of the TRU-Balance® 3 seating system and is available in matching shroud colors. The back shroud is not available with the back cane configuration. Reference colors shown on page 2.

Back in Black	Electric Blue	Pretty in Pink	Lime Light	Silver Bullet	Orange Crush
Red Corvette	White Satin	Deep Purple	Grey Street	Midnight Blue	Yellow Blaze

Additional TRU-Balance® 3 Back Options

Optional.

TRU-Comfort Lateral Wing, pair.....\$30
Adds lateral wing to Sport Backs. Standard with TRU-Comfort Backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Privacy Flap, 9.75"\$35
Standard on recline seating systems.

Privacy Flap, 12.75"\$35
Standard on recline seating systems.

See Comfort Company Back order form

ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

TRU-Balance® 3 Flip-Back Armrests

Choose a left and right TRU-Balance® 3 armrest and set height and optional downpost if desired. Different Left (L) and Right (R) armrest height packages are allowed. Not available with Back Cane Option.

<input type="checkbox"/> L <input type="checkbox"/> R	TRU-Balance 3 Flip-Back, 8"-12"H..... Standard	<input type="checkbox"/> L <input type="checkbox"/> R	TRU-Balance 3 Flip-Back, 12"-16"H Standard
---	--	---	--

Set LEFT height at: _____ Set RIGHT height at: _____

TRU-Balance 3 Armrest Downpost (Pair).....\$175 pr.
Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available on 14" Width seat or back.

	Tilt or Static Back Width												Recline Back Width											
	18"	19"	20"	21"	22"	24"	26"	28"	30"	32"			18"	19"	20"	21"	22"	24"	26"	28"	30"	32"		
Overall Back Height	16"	ST	ST	ST	S	ST	ST	S	S	S			-	-	-	-	-	-	-	-	-	-	-	-
	17"	ST	ST	ST	S	ST	ST	S	S	S			-	-	-	-	-	-	-	-	-	-	-	-
	18"	ST	ST	ST	S	ST	ST	ST	ST	S			-	-	-	-	-	-	-	-	-	-	-	-
	19"	ST	ST	ST	S	ST	ST	ST	ST	S			ST	ST	ST	S	ST	S	S	S	S	S	S	S
	20"	ST	ST	ST	S	ST	ST	ST	ST	S			ST	ST	ST	S	ST	S	S	S	S	S	S	S
	21"	ST	ST	ST	S	ST	ST	ST	ST	S			ST	ST	ST	S	ST	S	S	S	S	S	S	S
	22"	ST	ST	ST	S	ST	ST	ST	ST	S			ST	ST	ST	S	ST	S	S	S	S	S	S	S
	23"	ST	ST	ST	S	ST	ST	S	S	S			ST	ST	ST	S	ST	S	S	S	S	S	S	S
	24"	ST	S	ST	S	S	S	S	S	S			ST	ST	ST	S	ST	S	ST	ST	ST	ST	S	S
	25"	ST	S	ST	S	S	S	S	S	S			ST	ST	ST	S	ST	S	ST	ST	ST	ST	S	S
	26"	-	-	-	-	-	-	-	-	-			ST	ST	ST	S	ST	ST	S	S	S	S	S	S
	27"	-	-	-	-	-	-	-	-	-			ST	S	ST	S	S	S	S	S	S	S	S	S
	28"	-	-	-	-	-	-	-	-	-			ST	S	ST	S	S	S	S	S	S	S	S	S

☐ Sport Back ☐ Sport or TRU-Comfort Back

Sport and TRU-Comfort Back Availability Chart

ARMRESTS CONTINUED

Heavy Duty, Removable Armrests

Select a type of HD removable armrests and the set height. Only one height package is selectable as a normal configuration.

8"-12"H, Single Post, Height Adjustable,\$75 ea.
Removable, Heavy Duty (Sold as pair) HCPCS: E0973

10"-14"H, Single Post, Height Adjustable,\$75 ea.
Removable, Heavy Duty (Sold as pair) HCPCS: E0973

Set LEFT height at: _____ Set RIGHT height at: _____

ARMREST SIZE

REQUIRED. Select a Left (L) and Right (R) armrest size. This also affects the armpad selection below.

☐ L ☐ R 14" Full Length

☐ L ☐ R 10" Desk Length

ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad length and a type of left and right armpad. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or See Positioning Components Page or Order Form..... No Charge

Foam Armpads

Select a left and right standard type armpad (mix or match) OR choose the gel armpads that follow.

☐ L ☐ R Straight, Armpad Standard

☐ L ☐ R Waterfall, Armpad..... No Charge

Falcon Gel Armpads

Sold only in pairs. Select an armpad set below. Armpads are 3" wide. The mounting plate allows the armpad to be offset from the armrest tube.

Falcon Gel Left and Right Armpad:
Full Length\$92.50 ea
Desk Length\$70.00 ea
For use with TB3 standard armrests.

Falcon Gel Left and Right Armpad with Mounting Plate:
Full Length \$148.02 ea
Desk Length \$126.20 ea
For use with 2-Post, Heavy Duty Removable, or Pediatric Armrests.

JOYSTICK BRACKET RECEIVER

REQUIRED.

Standard Fixed Receiver Standard

Stealth Height Adjustable Receiver..... \$200
Tube Arm Part: ST-ARM250C-Q HCPCS: K0108
Track Arm Part: ST-ARM250T-Q HCPCS: K0108



HEADRESTS

Optional. Select a headrest pad and bracket.

Headrest Pads

Stealth 8" Comfort Plus\$142
Part: POS148484 HCPCS: E0955

Stealth 10" Comfort Plus (A)\$142
Part: SETHEAD1057 HCPCS: E0955

Stealth 14" Comfort Plus.....\$142
Part: INDPART2327 HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) \$142
Part: SETHEAD1069 HCPCS: E0955


Stealth Small Contoured (9.75"Wx3"H)..... \$60
Part: SETHEAD1052 HCPCS: E0955

Stealth Large Contoured (11"Wx6"H) (B)..... \$75
Part: SETHEAD1053 HCPCS: E0955

Headrest Bracket

Stealth's Worlds Best Removable.....\$210
Headrest Bracket (C)
Part: FRMASMB15999 HCPCS: E1028

Stealth Removable Headrest Bracket (D).....\$150
Part: HDWASMB1084 HCPCS: E1028



LAP BELTS

REQUIRED. Select a lap belt below.

Standard Lap Belts

38" Lap Belt	50" Lap Belt	60" Lap Belt	70" Lap Belt	80" Lap Belt	90" Lap Belt	No Lap Belt
--------------	--------------	--------------	--------------	--------------	--------------	-------------

Optional Lap Belts

i-Fit 60" Single Pull, Push Button,\$86 Medium (7.5") Padded, 1.5" Belt	i-Fit 60" Single Pull, Push Button,\$86 Extra Large Padded (11.5"), 2" Belt
i-Fit 60" Single Pull, Push Button,\$86 Large Padded (9.5"), 2" Belt	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Required with ICO 7176-19 Compliant Transit Option.

MANUFACTURER INSTALLED ACCESSORIES

Optional.

Push Handles (pair)\$175
No charge with Power Tilt or Power Recline systems. Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Accessory Bar\$100
Standard with the SCIM module, Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.

Trilogy Vent Tray \$400
Part: ACC136791 HCPCS: E1029

Compact Vent Tray \$1,395
Part: ACC133889 HCPCS: E1029



Trilogy Vent Tray



Compact Vent Tray

POWER LEGREST CONTROL

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)

Articulating Foot Platform (Requires One Actuator Control)

Combined Legrests (Requires 1 actuator control)

POWER POSITIONING ELECTRONICS

REQUIRED, except with Static Seat without power options. Select the desired control type for the seating system. Please note the available actuator controls. The 4-way toggle and 4 function push button are mounted on the same side as the joystick.

Power Tilt Only through Toggle

Single Actuator Control, Through Toggle (Power Tilt Only) No Charge
Only available with Power Tilt only. Manual recline is not available with NE Tilt through toggle. Fender lights are not available with tilt through toggle. Select toggle side.

Left Side Toggle

Right Side Toggle

NE+ 6-Key Controls

Single Actuator Control Through \$1,895
NE+ 6-Key Joystick HCPCS: E2310

Two Actuator Control Through \$2,895
NE+ 6-Key Joystick HCPCS: E2311

Q-Logic 3 Controls

Single Actuator Through Q-Logic 2 \$1,895
Controller and Interface (AM1) HCPCS: E2310

Multiple Actuator Through Q-Logic 2 \$2,895
Controller and Interface (AAM) HCPCS: E2311

Multiple Actuator Through 4-Way Toggle, \$2,895
Q-Logic 2 Controller and Interface (AAM)
Bracket mounts to outside of the armpad. HCPCS: E2311

Multiple Actuators Through 4-Way Button Toggle, \$2,895
Q-Logic 2 Controller and Interface (AAM)
Bracket mounts to outside of the armpad. HCPCS: E2311

Single Actuator Control Through 4 Function \$2,895
Push Button and Q-Logic 2 Controller and Interface (AAM)
Only available with TB3 armrests. HCPCS: E2310

Multiple Actuator Control Through 4 Function \$2,895
Push Button and Q-Logic 2 Controller and Interface (AAM)
Only available with TB3 armrests. HCPCS: E2311



Switch Controls

Attendant Control


Optional. Attendant control will be no charge when ordered with Power Tilt and Recline System.

Attendant Control \$500
Part: ELE125264 HCPCS: E2331

Seat Mounted, LED Light Package

Optional. Available only with NE+ 8 Key Joystick or Q-Logic 2 Ex Joystick for LED Lights.

Seat Mounted, LED Light Package \$995



Seat Mounted LED Lights (One Shown)

LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

Heavy Duty, 70°, Drop-In, Swing-Away Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

Heavy Duty, 70°, Swing-Away Legrest Upper

☐ ☐ Heavy Duty, 70° Drop-In, No Charge
☐ ☐ Swing-Away Legrests

Lower Extension

☐ ☐ 12"-14.5" Length ☐ ☐ 14.5"-18" Length
☐ ☐ ☐ ☐

Set Length at: _____

Angle Adjustable Footplates with Heel Loops

☐ ☐ 6"Wx8"D\$75 ea.
☐ ☐ HCPCS: K0040

☐ ☐ 6"Wx11"D\$75 ea.
☐ ☐ HCPCS: K0040

May cause interference on seat widths less than 18".

☐ ☐ 8"Wx11"D\$75 ea.
☐ ☐ HCPCS: K0040

May cause interference on seat widths less than 20".

Accessories

Wheel Bumpers (pair)\$25

Heavy Duty Elevating Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

Elevating Legrest

☐ ☐ Elevating Legrest Assembly \$162.50 ea.
☐ ☐ Select lower extension and footplates HCPCS: E0990

Lower Extension

☐ ☐ 11"-14" Length ☐ ☐ 12"-14.5" Length ☐ ☐ 14"-17" Length ☐ ☐ 14.5"-18" Length
☐ ☐ For reinforced only ☐ ☐ For angle adjustable only ☐ ☐ For reinforced only ☐ ☐ For angle adjustable only

Set Length at: _____

Reinforced Cast Footplate(s)

☐ ☐ 7.25"Wx6"D Reinforced Footplate Standard
☐ ☐

☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
☐ ☐ May cause interference on seat widths less than 18".

Angle Adjustable Footplate(s)

☐ ☐ Angle Adjustable with Heel Loops\$75 ea.
☐ ☐ (Select footplate size below) HCPCS: K0040

Angle Adjustable Size

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ 4"Wx6"D ☐ ☐ 5"Wx6"D ☐ ☐ 5"Wx8"D ☐ ☐ 6"Wx8"D¹ ☐ ☐ 6"Wx11"D¹ ☐ ☐ 8"Wx11"D²
☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐ ☐

Accessories

Heel Loops (pair)\$38
Part: FRMASMB7873 HCPCS: E0951
Standard on angle adjustable footplates.

Wheel Bumpers (pair)\$25

Falcon LPW Gel Padded Adductor Buttons (pair) \$125
Part: INDPART2242 HCPCS: K0108

LEGRESTS CONTINUED

Center Mount Foot Platforms

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section.

Power Articulating Foot Platform \$2,995
Please verify the appropriate billing code to be used for this item with all third-party payors.

Manually Adjustable Center Mount \$250
Foot Platform



Center Mount Lower Extension

10"-16" Extension 13"-19" Extension

Set Length at: _____

Center Mount Footplates

5"x6" Footplates 5"x8" Footplates 5"x11" Footplates

Center Mount Calf Support

Optional. Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pad (pair) 6"Wx6"H \$150	Small Calf Panel 18"Wx8"H \$150
Large Calf Pad (pair) 8" taper to 6"Wx8"H \$150	Large Calf Panel 24"Wx8"H \$200
	Only available on widths of 20" or greater.

Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension, 8"D and 11"D Footplates\$50	4" Outside Edge\$50
	Only available with 6"D footplates and 8"D footplate extension.
2" Bell-Shaped Length Extension\$50	6" Outside Edge.....\$50
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	Available with 8"D and 11"D footplates and 11"D footplate extension.

Additional Accessories

Reversed Dartex Calf Pad or Panel Cover (pair) \$120	Footplate Connector \$15
Wheel Bumpers (pair)\$25	

CUSHION QUICK PICK

Optional. General size selections are shown below. Refer to the [specific cushion order form](#) for detailed pricing, sizes, codes, weight capacity, and part numbers. Check the non-matching box and attached a cushion order form for a size other than the previously selected seat width and depth, or for a cushion not shown here.

Match seat selection width and depth

Non-matching cushion size, [see cushion order form](#) or cushion part number: _____

Cushion Selection

See cushion order form for detailed sizing and pricing.

Simplicity Cushion\$98	Spectrum® Gel Cushion.....\$480-\$580
General use. (Max 350 lbs) HCPCS: E2601	Positioning. (Max 350 lbs) HCPCS: E2607-E2608
Solution® Cushion.....\$380-\$430	Spectrum® Foam Cushion\$480-\$530
Skin protection and positioning. (Max 350 lbs) HCPCS: E2607-E2608	Skin protection and positioning. (Max 650 lbs) HCPCS: E2605-E2606
Solution® 1 Cushion.....\$380-\$430	TRU-Comfort 2 4-Way Stretch Cushion\$380-\$480
Skin protection. (Max 350 lbs) HCPCS: E2603-E2604	Skin protection. (Max 450 lbs) HCPCS: E2607-E2608

TRU-Comfort Sport Seat Cover

Available as an additional cover option for TRU-Comfort seat cushions. Only available in sizes below.

16"Wx16"D \$100	20"Wx20"D \$100
18"Wx18"D \$100	

PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU-Balance 3 Cup Holder No Charge
For use with TRU-Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest) No Charge
For use with tube-style armrests.

TRU-Balance 3 Backpack Holder (Pair) (A) \$25
Not available with back cane option.

TRU-Balance 3 Hydration System (B) \$175

TRU-Balance 3 Glove Box (C) \$100

TRU-Balance 3 Transfer Bar, \$175
Quick Release (Pair) (D) HCPCS: E0705

TRU-Balance 3 Front Mounted \$50
Cell Phone Holder (E)

Large, Black, Weather Cover \$158

TRU-Balance 3 Reflector Kit \$25


TRU-Balance 3 Personal Item Hook \$25
Not available with back cane option.

TRU-Balance 3 Rear View Mirror \$25

TRU-Balance 3 Oxygen Holder \$150
HCPCS: E2208

TRU-Balance 3 Tablet Holder \$375

TRU-Balance 3 Clothing Guard (Pair) \$175



BATTERIES

REQUIRED. Select a type of battery below, or to omit the batteries. The powerchair requires 2 batteries.

Group 24, 70 Amp Gel Cell, UPG Battery \$352 ea.
Part: BATGEL1004 HCPCS: E2363

Group 24, 70 Amp Gel Cell, MK Battery \$352 ea.
Part: BATGEL1009 HCPCS: E2363

No batteries. Provider will assume responsibility of batteries.

Battery Installation

Please choose if batteries will be installed, if selected above.

Install Batteries No Charge

Do NOT Install Batteries No Charge

BRANCH PATCHES

Optional.

Air Force Patch \$20

Marines Patch \$20

Arm Patch \$20

Navy Patch \$20

Coast Guard Patch \$20

Purple Heart Patch \$20



MANUFACTURER INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Thoracic Laterals

Swing-away laterals may interfere with TB3 and 2-post, flip up armrests. These laterals are not available with the back cane option.

Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$450
Part: INDPART2896 HCPCS: E1028



Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$200
Part: INDPART2663 HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$200
Part: INDPART2664 HCPCS: E0956

Installation

Install Laterals

Do Not Install

MANUFACTURER INSTALLED POSITIONING COMPONENTS CONTINUED

Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

Bracket

Flip Down Bracket (TWBTLADD-HW)..... \$225 ea.
QTY Part: INDPART2804 HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$225 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428 HCPCS: E1028

Flip Down with RIGHT Full Surface Contact \$225 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429 HCPCS: E1028

Removable Bracket (LP430)..... \$250 ea.
QTY Part: INDPART2805 HCPCS: E1028

Removable Bracket with 1" Offset \$325 ea.
(LP430 & LHW-121)
QTY Part: INDPART2854 HCPCS: E1028

Adapter

2" Height Adapter Plate (LHW-122)..... \$25
Part: INDPART2821 HCPCS: K0108

Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129)..... \$50
Part: POS140430 HCPCS: K0108

Only available with the flip down full surface bracket.

Pad

Short Pad, 4"x4" (TWBADD-4X4) \$50 ea.
QTY Part: INDPART2801 HCPCS: E0956

Medium Pad, 4"x8" (TWBADD-4X8) \$50 ea.
QTY Part: INDPART2802 HCPCS: E0956

Long Pad, 4"x12" (TWBADD-4X12) \$50 ea.
QTY Part: INDPART2803 HCPCS: E0956

Installation

Install on Left

Install on Right

Do Not Install

NON-INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$475
(SPLE107) (pair)
Part: INDPART2786 HCPCS: K0108

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$475
(SPLE108) (pair)
Part: INDPART2787 HCPCS: K0108

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$475
(SPLE109) (pair)
Part: INDPART2788 HCPCS: K0108

Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783 HCPCS: K0108

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
(SPLE112) (pair)
Part: INDPART2784 HCPCS: K0108

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
(SPLE114) (pair)
Part: INDPART2785 HCPCS: K0108

Stealth Elbow Support

Left, Swing-Away, Elbow Support Bracket..... \$175 ea.
QTY Tube Style Arms (ARMS-ES-L)
Part: POS124064 HCPCS: E1028

Left, Fixed Elbow Support Bracket..... \$175 ea.
QTY TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465

Right, Swing-Away, Elbow Support Bracket..... \$175 ea.
QTY Tube Style Arms (ARMS-ES-R)
Part: POS124066 HCPCS: E1028

Right, Fixed Elbow Support Bracket \$175 ea.
QTY TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466

Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$50 ea.
QTY Part: POS124065 HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$50 ea.
QTY Part: POS124067 HCPCS: K0108

NON-INSTALLED POSITIONING COMPONENTS CONTINUED

Stealth i-Fit Dynamic Chest Straps

Upper Classic, Medium (IFPS1103) \$142
Part: ST-IFPS1103-Q HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142
Part: ST-IFPS1104-Q HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$142
Part: ST-IFPS2103-Q HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142
Part: ST-IFPS2104-Q HCPCS: E0960

Chest Strap, Medium with Belt Anchor (IFCS4113) \$90
Part: ST-IFCS4113-Q HCPCS: E0960

Chest Strap, Large with Belt Anchor (IFCS4114) \$90
Part: ST-IFCS4114-Q HCPCS: E0960

Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.

Facial Lateral, 4.5"x2" Pad, Right Side, \$325
for comfort Series (CPS9-R-HW & 900-Q)
Part: INDPART2647

Facial Lateral, 4.5"x2" Pad, Left Side \$325
for comfort Series (CPS9-L-HW & 900-Q)
Part: INDPART2648

Facial Lateral, 4.5"x2" Pads, Left and Right Side \$425
for comfort Series (CPS9-B-HW & 900-Q)
Part: INDPART2649


Other Accessories

Egg Switch Mount (ESM) (A) \$50 ea.
Part: INDPART2815

Figure 8 Bracket with 13" Rod and..... \$250
Egg Switch Mount (B)
Part: ST-F8675-13-Q

Facial Pad, 6"x2" (920-Q) \$50
Part: ST-920-Q HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$250
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest



SPECIALTY CONTROLS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

Enhanced Display

Q-Logic 2 Enhanced Display Kit \$1,884.93
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Q-Logic 2 Enhanced Display Shroud

REQUIRED with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Match Joystick Shroud to Base Color \$50
Part: See color below.

Non-matching Joystick Shroud \$50
Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze

SPECIALTY CONTROLS CONTINUED

Head Array

Adult i-Drive Tri Array Series, 3-Switch with Egg Switch for Mode (IDH200-1)..... 4,175

Part: CTL135940

HCPCS: E2330

Features 3 embedded proximity sensors for drive control and an egg switch for mode change.



Adult Tri-Array

Other

Multiplier Harness.....\$221.21

QTY Part: DWR1111H053

Red Egg Switch (ESRED) \$77

QTY Part: SWTMCRO1077

Blue Egg Switch (ESBLU) \$77

QTY Part: SWTMCRO1079

Black Egg Switch (ESBLK) \$77

QTY Part: SWTMCRO1081

Green Egg Switch (ESGRN)..... \$77


QTY Part: SWTMCRO1078

Yellow Egg Switch (ESYEL)..... \$77

QTY Part: SWTMCRO1080

SCIM Module Kit \$575

Allows use of joystick display for specialty devices. Available with Q-Logic 2 EX Joystick selection above. Includes the TRU-Balance® 3 accessory bar. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This is also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



facebook.com/QuantumRehab



youtube.com/user/QuantumRehab