



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® Z Order Form with Pediatric TRU-Balance® 3 Positioning - Group 3, ODJFS Only

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Selections with HCPCS codes may have a part number noted with the selection. Other part numbers can be found on the quote acknowledgment after this order form is submitted. Please contact Quantum® Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.

If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



BASE MODEL

REQUIRED. Select a power base. 300 lb. weight capacity with standard 6mph programmable motors, 8 amp off board charger, and switch operated fender lights. Dimensions of base without legrests: 26.25"Wx35.5"L.

Quantum Q6 Edge Z 3S-SS \$7,895

Part: Q6 EDGE Z 3S-SS

MFG Sugg HCPCS: K0848/K0011

Solid seat pan standard with no power option.

Quantum Q6 Edge Z 3MP-SS \$7,895

Part: Q6 EDGE Z 3MP-SS

MFG Sugg HCPCS: K0861/K0011

Multiple power option.

Quantum Q6 Edge Z 3SP-SS \$7,895

Part: Q6 EDGE Z 3SP-SS

MFG Sugg HCPCS: K0856/K0011

Single power option.

BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.



Red Corvette



Electric Blue



Back in Black



White Satin



Pretty in Pink



Deep Purple



Lime Light



Grey Street



Silver Bullet



Midnight Blue



Orange Crush



Yellow Blaze



To be
determined
at a later
date.

DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard

Community Use Package \$2,125

300 lbs. weight capacity, 4-pole high speed motor package, 4" wide flat-free or pneumatic extended traction tires. Increases overall base width by 1.5".

High Speed Motor Package \$1,574

Programmable speed up to 8 mph. High speed motors decrease per-charge range. Not available with iLevel seating.

Accu-Trac Motor Technology \$1,295

HCPCS: K0108

Requires Q-Logic 2 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. Programmable speed up to 6 mph.

BASE MOUNTED FOOT PLATFORM

REQUIRED. Select a foot platform option. May interfere with seat mounted legrests.

Base Mounted Foot Platform Standard

Not available with center mount foot platforms.

Omit base mounted foot platform No Charge

DRIVE WHEELS

REQUIRED. Select a type of drive wheels. Knobby tires increase overall base width by 1.5".

14" Run Flat Tires Standard

4"W Flat Free, Knobby, Extended Traction Tires \$551

4"W Pneumatic, Knobby, Extended Traction Tires \$551

ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

NE+ 6-Key Joystick (1 Actuator) Standard

Control of one actuator through Joystick on single power base model (3SP-SS) only.

NE+ 6-Key Joystick (2 Actuators) No Charge

Control up to 2 actuators through joystick.

NE+ 8-Key Joystick (2 Actuators) \$400

Control up to 2 actuators and Seat Mounted LED Light Package through joystick.



NE+ 6 Key Joystick

ELECTRONICS CONTINUED

Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers. This sections allows a selection of common configurations.

Q-Logic 2 Required Components

Q-Logic 2 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected later in this form, Customer Service will quote the order with these items.

Q-Logic 2 EX Controller Module	\$695	Q-Logic 2 Expandable Controller.....	\$425
<i>HCPCS: E2377</i>		<i>Part: DWR1111H065</i>	<i>HCPCS: E2313</i>



Q-Logic 2 Joystick

Q-Logic 2 EX Joystick.....	\$700	Stand Alone Joystick with	\$700
Q-Logic 2 Ex Joystick for Seat.....	\$700	External Charger Port	
Mounted, LED Light Package		Up to 3 drive profiles and 1 seating profile when used alone. Use in conjunction with Enhanced Display for up to 5 profiles and auxiliary functions such as Bluetooth Mouse and IR. Includes charger port mounting bracket and multiplier box harness.	
No Joystick, Setup for Q-Logic			



Q-Logic 2 Joystick Shroud

REQUIRED with Q-Logic 2 EX Joystick selection above. Reference colors shown on page 2.

Black, Soft Touch	Standard	Match Joystick Shroud to Base Color.....	\$50
Non-matching Joystick Shroud	\$50		
Select color below			
Red Corvette	White Satin	Lime Light	Midnight Blue
Electric Blue	Pretty in Pink	Grey Street	Orange Crush
Back in Black	Deep Purple	Silver Bullet	Yellow Blaze



JOYSTICK HANDLES

Optional. The NE Plus and Q-Logic 2 EX Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)	\$71	Soft Ball Handle (F).....	\$50
<i>Part: ACCASMB1860</i>	<i>HCPCS: E2323</i>	<i>Part: ACCASMB1856</i>	<i>HCPCS: E2323</i>
Chin Cup (C)	\$71	Extended I-Handle (G).....	\$71
<i>Part: ACCASMB1861</i>	<i>HCPCS: E2323</i>	<i>Part: ACCASMB1857</i>	<i>HCPCS: E2323</i>
	<i>(E2324 when chin control)</i>		
T-Bar (D)	\$107	Bodypoint 3" Goal Post Handle.....	\$69
<i>Part: ACCASMB1858</i>	<i>HCPCS: E2323</i>	<i>Part: PLSKNCP1057</i>	<i>HCPCS: E2323</i>
Medium Ball Handle (E).....	\$50	Bodypoint 4" Goal Post Handle.....	\$69
<i>Part: ACCASMB1859</i>	<i>HCPCS: E2323</i>	<i>Part: PLSKNCP1058</i>	<i>HCPCS: E2323</i>



JOYSTICK MOUNTING BRACKETS

REQUIRED. Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket or See Positioning Positioning Components Page or [Order Form](#)



Fixed, Inline Joystick Bracket

Fixed, Inline, Left Side Mount	Standard	Fixed, Inline, Right Side Mount	Standard
--------------------------------------	----------	---------------------------------------	----------

Swing-Away (Retractable) Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" is used with TRU-Balance 3 arms when the joystick is inline with armpad. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

Swing-away, Inline, Left Side Mount.....	\$245	Swing-away, 1" Offset, Left Side Mount	\$245
<i>NE/NE+ Part: ACC142526</i>	<i>HCPCS: E1028</i>	<i>NE/NE+ Part: ACC142747</i>	<i>HCPCS: E1028</i>
<i>Q-Logic: ACC142524</i>		<i>Q-Logic: ACC142742</i>	
Swing-away, Inline, Right Side Mount	\$245	Swing-away, 1" Offset, Right Side Mount....	\$245
<i>NE/NE+ Part: ACC142527</i>	<i>HCPCS: E1028</i>	<i>NE/NE+ Part: ACC142748</i>	<i>HCPCS: E1028</i>
<i>Q-Logic: ACC142525</i>		<i>Q-Logic: ACC142743</i>	




JOYSTICK MOUNTING BRACKETS CONTINUED

Retract4 Swing-Away Arm Kit

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Retract4, Left Side Mount \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860

Retract4, Right Side Mount \$250
NE/NE+ Part: MEC138876 HCPCS: E1028
Q-Logic: MEC138860



Multi-Axis

Allows height adjustment of joystick.

Multi-axis, Left Side Mount \$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216

Multi-axis, Right Side Mount \$100
NE/NE+ Part: ACC142528 HCPCS: E1028
Q-Logic: ACC131216



Stealth Swing Away/Flip Down Bracket

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket \$275
Left Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028

Stealth Swing-Away/Flip Down Bracket \$275
Right Side Mount
Stand Alone Part: ST-TWBM480CJ-Q HCPCS: E1028
Q-Logic Part: ST-TWBM480QL-Q HCPCS: E1028



ELECTRONICS ACCESSORIES

Optional.

USB Mobile Device Charger \$175
Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

Multiplier Harness \$221.21

Q-Logic EX/NE Programmer \$545

Q-Logic PC Programming Disk with Cable ... \$750



TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™

REQUIRED. Select one seating system. Be sure to complete the Seat Dimension, Seat to Floor, and back selections.

TRU-Balance® 3 Configuration

TRU-Balance 3 Static Seat \$2,295
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Recline \$7,595
Part: PTOINDV3444 HCPCS: E1004
175° of power recline.

TRU-Balance 3 Power Tilt \$6,995
Part: PTOINDV3443 HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Tilt & Recline \$11,895
Part: PTOINDV3446 HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

10" Power Lift or Manual Recline

Optional. Select one option below.

TRU-Balance 3 10" Lift with iLevel® \$2,495
Part: PTOINDV3773 HCPCS: E2300
10" of scissor lift with iLevel® System. 22"W x 22"D maximum seat size.

Attendant Manual Recline* \$500
Part: PTOINDV3501 HCPCS: E1226

*Available through Special Order. Not available with NE+ joystick with Tilt through toggle. Available with Static Seat or Power Tilt. Not available with Power Lift.



Occupied Transit Kit

Optional. The Q6 Edge Z is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option \$200

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://www.ercwts.org>.

TRU-BALANCE® 3 SEATING CONTINUED

Seat Dimensions

REQUIRED with seat selection. Select width and depth package and enter desired set width and depth in that package range, available in 1" increments. Configurations disallowed or not shown may be available through special order.

Seat Width

12"-16" Width Package

14"-16" Width Package

Set Width at: _____

Seat Depth

12"-16" Depth Package

14"-16" Depth Package

Set Depth at: _____

Seat to Floor Height

REQUIRED. Select a seat to floor height compatible with the seat functions previously selected.

Static Seat or Recline Only

17.5625" Seat to Floor

18.5625" Seat to Floor

19.5625" Seat to Floor

Tilt Only or Tilt and Recline

18.125" Seat to Floor

19.125" Seat to Floor

20.125" Seat to Floor

Power Lift (with or without Tilt or Recline)

19.375" Seat to Floor

20.375" Seat to Floor

21.375" Seat to Floor

Back Selection

REQUIRED. Match choice in this section with the options in the back cushion configuration section.

Static Back Standard
Select a back cushion configuration in the subsequent section.

Back Cane Option, Omit Back Pan/ Cushion(-\$125)
Available only with Static Back and Power Tilt (1 actuator) selections.
Select cane type and the back cane height below.

Flip Forward Back \$150
Available with Static Back and Power Tilt (1 actuator). Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width or with attendant manual recline. Select a back cushion type in the subsequent section.

Back Cane

Select with Back Cane Option selected above.

Angle Adjustable Standard

1" Endo \$35

Flip Forward \$25

2" Endo \$35

Back Cane Height

Select with Back Cane Option selected above.

16"

17"

20"

21"

24"

25"

Back Cushion Configuration

REQUIRED except when back cane option selected. Select a configuration type and applicable dimensions. Other size backs may be available with Special Order and extended lead time. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

Sport Back No Charge
Complete following width package and size..

TRU-Comfort 4-way Stretch, 20" Width or less \$595
Complete following width package and size. HCPCS: E2620

Omit Back Pan/Cushion No Charge
Complete following width package.

Aftermarket Back Interface \$200
Enter the preset back width below. Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/water-fall arm pads and arm through. Does not include brackets to mount aftermarket backs. A Special Order is required on back widths less than 16".

Back Cushion Width

Available with +/-4" of seat width.

10"W-14"W Package

16"W-20"W Package

Set Width at: _____

	Tilt or Static Width			Recline Width		
	12"	14"	16"	12"	14"	16"
14"	S TC	S	-	-	-	-
15"	S TC	S	-	-	-	-
16"	S	S TC	S TC	S	S	S TC
17"	S	S TC	S TC	S	S	S TC
18"	-	-	S TC	-	-	S TC
19"	-	-	S TC	-	-	S TC
20"	-	-	S TC	-	-	S TC
21"	-	-	S TC	-	-	S TC
22"	-	-	S TC	-	-	S TC
23"	-	-	S TC	-	-	S TC
24"	-	-	S TC	-	-	S

S Sport Back S TC Sport or TRU-Comfort Back

Sport and TRU-Comfort Back Availability Chart

TRU-BALANCE® 3 SEATING CONTINUED

Back Cushion Overall Height

Static Seat and Power Tilt available height range: 14"-24". *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.* **Power Recline available height range:** 16"-24". *Odd heights include 6" gap and even heights include 7" gap.* Please see the Q6Edge Z TRU-Balance 3 order form for larger sizes. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

Set Overall Height at: _____

Aftermarket Back Interface Preset Width

Set Interface Width at: _____ (Range 16"-24" Wide)

TRU-Balance® 3 Back Shroud

REQUIRED. The back shroud is a feature of the TRU-Balance® 3 seating system and is available in matching shroud colors. The back shroud is not available with the back cane configuration. Reference colors shown on page 2.

Back in Black	Electric Blue	Pretty in Pink	Lime Light	Silver Bullet	Orange Crush
Red Corvette	White Satin	Deep Purple	Grey Street	Midnight Blue	Yellow Blaze

Additional TRU-Balance® 3 Back Options

Optional.

TRU-Comfort Lateral Wings, pair\$30
Adds lateral wing to Sport Backs. Standard with TRU-Comfort Backs.
Not available on power tilt backs 16"-17"H with limited adjustability
on 18" H. Not available on power recline back 19"-20"H with limited
adjustability on 21"H.

Privacy Flap, 9.75"\$35
Standard on recline seating systems.

Privacy Flap, 12.75"\$35
Standard on recline seating systems.

See Comfort Company Back order form

ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

TRU-Balance® 3 Flip-Back Armrests

Choose a left and right TRU-Balance® 3 armrest and set height and optional downpost if desired. Different Left (L) and Right (R) armrest height packages are allowed. Not available with Back Cane Option.

<input type="checkbox"/> L <input type="checkbox"/> R	TRU-Balance 3 Flip-Back, 8"-12"H.....Standard	<input type="checkbox"/> L <input type="checkbox"/> R	TRU-Balance 3 Flip-Back, 12"-16"HStandard
---	---	---	---

Set LEFT height at: _____ Set RIGHT height at: _____

Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

Pediatric, Armrests\$75 ea.
HCPCS: E0973

Set LEFT height at: _____ Set RIGHT height at: _____

ARMREST SIZE

REQUIRED. Select a Left (L) and Right (R) armrest size. This also affects the armpad selection below.

<input type="checkbox"/> L <input type="checkbox"/> R	14" Full Length	<input type="checkbox"/> L <input type="checkbox"/> R	10" Desk Length
---	-----------------	---	-----------------

ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad length and a type of left and right armpad. Please note that waterfall and gel armpads may cause interference with the seat back.

No Armpad Selection and/or No Charge
See Positioning Components Page or Order Form

Foam Armpads

Select a left and right standard type armpad (mix or match) OR choose the gel armpads that follow.

<input type="checkbox"/> L <input type="checkbox"/> R	Straight, Armpad Standard	<input type="checkbox"/> L <input type="checkbox"/> R	Waterfall, Armpad..... No Charge
---	---------------------------------	---	----------------------------------

ARMPADS CONTINUED

Falcon Gel Armpads

Sold only in pairs. Select an armpad set below. Armpads are 3" wide. The mounting plate allows the armpad to be offset from the armrest tube.

Falcon Gel Left and Right Armpad:

Full Length\$92.50 ea
Desk Length\$70.00 ea

For use with TB3 standard armrests.

Falcon Gel Left and Right Armpad with Mounting Plate:

Full Length \$148.02 ea
Desk Length \$126.20 ea

For use with 2-Post, Heavy Duty Removable, or Pediatric Armrests.

JOYSTICK BRACKET RECEIVER

REQUIRED.

Standard Fixed Receiver Standard

Stealth Height Adjustable Receiver \$200

Tube Arm Part: ST-ARM250C-Q

HCPCS: K0108

Track Arm Part: ST-ARM250T-Q

HCPCS: K0108



Height Adjustable
Receiver (Tube Arm)



Height Adjustable
Receiver (Track Arm)

HEADRESTS

Optional. Select a headrest pad and bracket.

Headrest Pads

Stealth 8" Comfort Plus\$142

Part: POS148484

HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H)\$142

Part: SETHEAD1069

HCPCS: E0955

Stealth 10" Comfort Plus (A)\$142

Part: SETHEAD1057

HCPCS: E0955

Stealth Small Contoured (9.75"Wx3"H) \$60

Part: SETHEAD1052

HCPCS: E0955

Headrest Bracket

Stealth's Worlds Best Removable.....\$210

Headrest Bracket (C)

Part: FRMASMB15999

HCPCS: E1028

Pediatric Offset (HMO475P-17) (E)\$150


Part: HDWASMB1166

HCPCS: E1028

Stealth Removable Headrest Bracket (D)\$150

Part: HDWASMB1084

HCPCS: E1028



Headrests

LAP BELTS

REQUIRED. Select a lap belt below.

Standard Lap Belts

38" Lap Belt

50" Lap Belt

60" Lap Belt

Optional Lap Belts

i-Fit 60" Single Pull, Push Button,\$86
Medium (7.5") Padded, 1.5" Belt

i-Fit 60" Single Pull, Push Button,\$86
Large Padded (9.5"), 2" Belt

i-Fit 60" Single Pull, Push Button,\$86
Extra Large Padded (11.5"), 2" Belt

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ICO 7176-19 Compliant Transit Option.

MANUFACTURER INSTALLED ACCESSORIES

Optional.

Push Handles (pair)\$175

No charge with Power Tilt or Power Recline systems. Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.

Accessory Bar\$100

Standard with the SCIM module, Trilogy vent tray, and oxygen holder. May be added to mount additional accessory to back of unit. Not available with back cane option.

Trilogy Vent Tray \$400

Part: ACC136791

HCPCS: E1029


Compact Vent Tray \$1,395

Part: ACC133889

HCPCS: E1029



Trilogy Vent Tray



Compact Vent Tray

POWER LEGREST CONTROL

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)

Articulating Foot Platform (Requires 1 actuator control)

Combined Legrests (Requires 1 actuator control)

POWER POSITIONING ELECTRONICS

REQUIRED, except with Static Seat and no power options. Select the desired control type for the seating system. Please note the available actuator controls. The 4-way toggle and 4 function push button are mounted on the same side as the joystick.

Single Actuator Control, Through No Charge Toggle (Power Tilt Only)*

*Only available with Power Tilt only though Special Order. Manual recline is not available with NE Tilt through toggle. Fender lights are not available with tilt through toggle. Select toggle side.

Left Side Toggle

Right Side Toggle

Through NE+ 6-Key Joystick with either:

Single Actuator Control..... \$1,895
Select this price/code with 1 actuator. HCPCS: E2310

Two Actuator Control..... \$2,895
Select this price/code with 2 actuators. HCPCS: E2311

Single Actuator Through Q-Logic 2 \$1,895
Controller and Interface (AM1) HCPCS: E2310

Multiple Actuator Through Q-Logic 2 \$2,895
Controller and Interface (AAM) HCPCS: E2311

Multiple Actuator Through 4-Way..... \$2,895
Toggle, Q-Logic 2 Controller and Interface (AAM) HCPCS: E2311
Bracket mounts to outside of the armpad.

Multiple Actuators Through 4-Way..... \$2,895
Button Toggle, Q-Logic 2 Controller and Interface (AAM) HCPCS: E2311
Bracket mounts to outside of the armpad.

Through 4 Function Push Button and Q-Logic 2 Controller and Interface with either:

Single Actuator Control (AAM) \$2,895
Select this price/code with 1 actuator. HCPCS: E2310

Multiple Actuator Control (AAM) \$2,895
Select this with multiple actuators. HCPCS: E2311
Only available with TB3 armrests.



Attendant Control

Optional. Attendant control will be no charge when ordered with Tilt and Recline System.

Attendant Control \$500
Part: ELE125264 HCPCS: E2331

Seat Mounted, LED Light Package

Optional. Available only with NE+ 8 Key Joystick or Q-Logic 2 Ex Joystick for LED Lights.

Seat Mounted, LED Light Package..... \$995



LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

Center Mount Foot Platforms

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section.

Power Articulating Foot Platform \$2,995
Please verify the appropriate billing code to be used for this item with all third-party payors.

Manually Adjustable Center Mount \$250
Foot Platform

Center Mount Lower Extension

8"-12" Extension 13"-19" Extension
N/A with Power Articulation.
10"-16" Extension

Set Length at: _____

Center Mount Footplates

5"x6" Footplates 5"x8" Footplates 5"x11" Footplates

Center Mount Calf Support

Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pad (pair) 6"Wx6"H.....\$150

Large Calf Pad (pair) 8" taper to 6"Wx8"H..\$150

Small Calf Panel 18"Wx8"H\$150

Large Calf Panel 24"Wx8"H..... \$200
Only available on widths of 20" or greater.

Center Mount Foot Platform Compatibility Matrix for TRU-Balance® 3	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
10"-16" Length Extension	!	✓	✓	✓	✓	✓	✓	✓	✓
13"-19" Length Extension	!	!	!	!	!	!	!	!	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.



LEGRESTS CONTINUED

Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension, 8"D and 11"D Footplates\$50

2" Bell-Shaped Length Extension\$50

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

4" Outside Edge\$50

Only available with 6"D footplates and 8"D footplate extension.

6" Outside Edge.....\$50

Available with 8"D and 11"D footplates and 11"D footplate extension.

Additional Accessories

Reversed Dartex Calf Pad or Panel Cover (pair) \$120

Footplate Connector \$15

Wheel Bumpers (pair)\$25

90° Pediatric Swing Away Legrests

Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will Match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

Pediatric, Swing Away Legrest

☐ ☐ 90° Pediatric, Swing-Away Legrest.....\$81 ea.
L R

Lower Extension

☐ ☐ 7"-10" Length
L R

Set Length at:_____

Angle Adjustable Footplate with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ 4"Wx6"D\$75 ea.
L R HCPCS: K0040

☐ ☐ 5"Wx6"D\$75 ea.
L R HCPCS: K0040

Accessories

Wheel Bumpers (pair)\$25

Elevating Legrests

Select a configuration, a legrest upper, a legrest extension, and type of footplate. Left (L) and right (R) configurations are selected by checking the appropriate boxes.

Configuration Verification

Left and Right Legrests will match

Left and Right Legrests will not match

One Legrest (Left or Right as Marked)

Elevating Legrest

☐ ☐ Elevating Legrest Assembly \$162.50 ea.
L R Select lower extension and footplates HCPCS: E0990

Lower Extension

☐ ☐ 11"-14" Length ☐ ☐ 14"-17" Length
L R For reinforced only L R For reinforced only

☐ ☐ 12"-14.5" Length ☐ ☐ 14.5"-18" Length
L R For angle adjustable only L R For angle adjustable only


Set Length at:_____


Reinforced Cast Footplate(s)

☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
L R

☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
L R May cause interference on seat widths less than 18."

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-14.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	✓	!	✓	✓

 This rigging configuration causes no interference issues at this seat-to-floor height.

 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

LEGRESTS CONTINUED

Angle Adjustable Footplate(s)

Angle Adjustable with Heel Loops.....\$75 ea.	
L R	(Select footplate size below) HCPCS: K0040

Angle Adjustable Size

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

L R	4"Wx6"D	L R	5"Wx6"D	L R	5"Wx8"D	L R	6"Wx8"D ¹	L R	6"Wx11"D ¹	L R	8"Wx11"D ²
-----	---------	-----	---------	-----	---------	-----	----------------------	-----	-----------------------	-----	-----------------------

Accessories

Heel Loops (pair)	\$38
Part: FRMASMB7873	HCPCS: E0951
Standard on angle adjustable footplates.	

Wheel Bumpers (pair)	\$25
Gel Padded Adductor Buttons (pair).....	\$125
Part: INDPART2242	HCPCS: K0108

CUSHION QUICK PICK

Optional. General size selections are shown below. Refer to the [specific cushion order form](#) for detailed pricing, sizes, codes, and part numbers. Check the non-matching box and attached a cushion order form for a size other than the previously selected seat width and depth, or for a cushion not shown here.

Match seat selection width and depth

Non-matching cushion size, [see cushion order form](#)

Solution® Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D.....	\$380	22"Wx16"D to 24"Wx24"D	\$430
	HCPCS: E2607		HCPCS: E2608

Solution® 1 Cushion (Skin Protection)

14"Wx14"D to 20"Wx24"D.....	\$380	22"Wx16"D to 24"Wx24"D	\$430
	HCPCS: E2603		HCPCS: E2604

Spectrum® Gel Cushion (Positioning)

12"Wx10"D to 20"Wx24"D.....	\$480	22"Wx16"D to 24"Wx24"D	\$530
	HCPCS: E2605		HCPCS: E2606

Spectrum® Foam Cushion (Skin Protection and Positioning)

10"Wx10"D to 20"Wx24"D.....	\$480	22"Wx16"D to 24"Wx24"D	\$530
	HCPCS: E2607		HCPCS: E2608

TRU-Comfort™ 2 4-Way Stretch Cushion (Skin Protection)

16"Wx14"D to 20"Wx24"D	\$380	22"Wx16"D to 22"Wx24"D	\$480
	HCPCS: E2603		HCPCS: E2604

TRU-Comfort™ Sport Seat Cover

Available as an additional cover option for both TRU-Comfort seat cushions. Only available in sizes below.

16"Wx16"D.....	\$100	20"Wx20"D	\$100
18"Wx18"D	\$100		

PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU-Balance 3 Cup Holder No Charge
For use with TRU-Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest) No Charge
For use with tube-style armrests.

TRU-Balance 3 Backpack Holder (Pair) (A) \$25
Not available with back cane option.

TRU-Balance 3 Hydration System (B)\$175

TRU-Balance 3 Glove Box (C).....\$100

TRU-Balance 3 Transfer Bar,\$175
Quick Release (Pair) (D) HCPCS: E0705

TRU-Balance 3 Front Mounted \$50
Cell Phone Holder (E)

Large, Black, Weather Cover.....\$158

TRU-Balance 3 Reflector Kit \$25


TRU-Balance 3 Personal Item Hook..... \$25
Not available with back cane option.

TRU-Balance 3 Rear View Mirror \$25

TRU-Balance 3 Oxygen Holder.....\$150
HCPCS: E2208

TRU-Balance 3 Tablet Holder \$375

TRU-Balance 3 Clothing Guard (Pair)\$175



BATTERIES

REQUIRED. Select a type of battery below, or to omit the batteries. The powerchair requires 2 batteries.

Group 24, 70 Amp Gel Cell, UPG Battery..... \$295 ea.
Part: BATGEL1004 HCPCS: E2363

Group 24, 70 Amp Gel Cell, MK Battery \$295 ea.
Part: BATGEL1009 HCPCS: E2363

No batteries. Provider will assume responsibility of batteries.

Battery Installation

Please choose if batteries will be installed, if selected above.

Install Batteries No charge

Do NOT install batteries..... No Charge

BRANCH PATCHES

Optional.

Air Force Patch..... \$20

Marines Patch \$20

Arm Patch \$20

Navy Patch \$20

Coast Guard Patch..... \$20

Purple Heart Patch..... \$20



MANUFACTURER INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Thoracic Laterals

Swing-away laterals may interfere with TB3 and 2-post, flip up armrests. These laterals are not available with the back cane option.

Bracket

Swing-away, Rail Mounts, pair (TWBSTL)..... \$450
Part: INDPART2896 HCPCS: E1028



Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$200
Part: INDPART2663 HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$200
Part: INDPART2664 HCPCS: E0956

Installation

Install Laterals

Do Not Install

Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

Bracket

Flip Down Bracket (TWBTLADD-HW)..... \$225 ea.
Part: INDPART2804 HCPCS: E1028

Removable Bracket (LP430)..... \$250 ea.
Part: INDPART2805 HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$225 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428 HCPCS: E1028

Removable Bracket with 1" Offset \$325 ea.
(LP430 & LHW-121)
Part: INDPART2854 HCPCS: E1028

Flip Down with RIGHT Full Surface Contact \$225 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429 HCPCS: E1028

Adapter

2" Height Adapter Plate (LHW-122)..... \$25
Part: INDPART2821 HCPCS: K0108

2 1/4" Adapter Plate (LHW-129) \$50
Part: POS140430 HCPCS: K0108

Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

Only available with the flip down full surface bracket.

Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4) \$50 ea.
Part: INDPART2801 HCPCS: E0956

Long Pad, 4"x12" (TWBADD-4X12) \$50 ea.
Part: INDPART2803 HCPCS: E0956

Medium Pad, 4"x8" (TWBADD-4X8) \$50 ea.
Part: INDPART2802 HCPCS: E0956

Installation

Install on Left

Install on Right

Do Not Install

NON-INSTALLED POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$475
(SPLE107) (pair)
Part: INDPART2786 HCPCS: K0108

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$475
(SPLE108) (pair)
Part: INDPART2787 HCPCS: K0108

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$475
(SPLE109) (pair)
Part: INDPART2788 HCPCS: K0108

Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783 HCPCS: K0108

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
(SPLE112) (pair)
Part: INDPART2784 HCPCS: K0108

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
(SPLE114) (pair)
Part: INDPART2785 HCPCS: K0108

Stealth Elbow Support

Left, Swing-Away, Elbow Support Bracket..... \$175 ea.
QTY Tube Style Arms (ARMS-ES-L)
Part: POS124064 HCPCS: E1028

Left, Fixed Elbow Support Bracket \$175 ea.
QTY TRU-Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465

Right, Swing-Away, Elbow Support Bracket..... \$175 ea.
QTY Tube Style Arms (ARMS-ES-R)
Part: POS124066 HCPCS: E1028

Right, Fixed Elbow Support Bracket \$175 ea.
QTY TRU-Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466

Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$50 ea.
QTY Part: POS124065 HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$50 ea.
QTY Part: POS124067 HCPCS: K0108

Stealth i-Fit Dynamic Chest Straps

Upper Classic, Medium (IFPS1103) \$142
Part: ST-IFPS1103-Q HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142
Part: ST-IFPS1104-Q HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$142
Part: ST-IFPS2103-Q HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142
Part: ST-IFPS2104-Q HCPCS: E0960

Chest Strap, Medium, Belt Anchor (IFCS4113) \$90

Chest Strap, Large, Belt Anchor (IFCS4114) \$90

Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.

Facial Lateral, 4.5"x2" Pad, Right Side, \$325
for comfort Series (CPS9-R-HW & 900-Q)
Part: INDPART2647

Facial Lateral, 4.5"x2" Pad, Left Side \$325
for comfort Series (CPS9-L-HW & 900-Q)
Part: INDPART2648

Facial Lateral, 4.5"x2" Pads, Left and Right Side \$425
for comfort Series (CPS9-B-HW & 900-Q)
Part: INDPART2649


Other Accessories

Egg Switch Mount (ESM) (A) \$50 ea.
QTY Part: INDPART2815

Figure 8 Bracket with 13" Rod and..... \$250
Egg Switch Mount (B)
Part: ST-F8675-13-Q HCPCS: K0108

Facial Pad, 6"x2" (920-Q) \$50
Part: ST-920-Q HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$250
Part: ST-SUS18-Q HCPCS: K0108
Allows mounting of facial pads and attaches to the headrest



SPECIALTY CONTROLS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

Enhanced Display

Q-Logic 2 Enhanced Display Kit\$1,884.93

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 2 Enhanced Display Shroud

Required with Q-Logic 2 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Match Joystick Shroud to Base Color..... \$50

Non-matching Joystick Shroud \$50

Select color below

Red Corvette

White Satin

Lime Light

Midnight Blue

Electric Blue

Pretty in Pink

Grey Street

Orange Crush

Back in Black

Deep Purple

Silver Bullet

Yellow Blaze

Head Array

Pediatric i-Drive Tri Array Series, 3-Switch with Egg Switch for Mode (IDH250-1)\$4,175

Part: CTL135956

HCPCS: E2330

Features 3 embedded proximity sensors for drive control and an egg switch for mode change.



Pediatric Tri-Array

Other

Multiplier Harness.....\$221.21

Part: DWR1111H053

Yellow Egg Switch (ESYEL)..... \$77

Part: SWTMCR01080

Red Egg Switch (ESRED) \$77

Part: SWTMCR01077

SCIM Module Kit \$575

Allows use of joystick display for specialty devices. Available with Q-Logic 2 EX Joystick selection above. Includes the TRU-Balance® 3 accessory bar. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

Blue Egg Switch (ESBLU) \$77

Part: SWTMCR01079

Black Egg Switch (ESBLK) \$77

Part: SWTMCR01081

Green Egg Switch (ESGRN)..... \$77

Part: SWTMCR01078



SCIM Module

REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.



facebook.com/QuantumRehab



youtube.com/user/QuantumRehab