

Jazzy Select

Order Form / Price List
ODJSF Only
Revised: February 2011



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

Account #:	Date:	Order No.:
Bill To:	Ship To:	P.O. #:
Address:	Address:	Discount:
Phone Contact:	Marked For:	Terms:
<p>This product is available for purchase online at www.mypride.com</p> <p>DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.</p>		<p>All pricing is MSRP.</p> <p>Freight:</p> <p>Taken By:</p>

BASE MODEL

- Jazzy Select -C \$5,496
Comes standard with Solid Tires
K0823/K0011 Mfg. Suggested
- Jazzy Select-SS \$5,496
Comes standard with Solid Tires
K0822/K0011 Mfg. Suggested

COLOR

- Jazzy Red NC
- Jazzy Blue NC

BATTERIES (2 required)

- U-1 Installed E2365 \$133ea

JOYSTICK MOUNT

- Fixed, In-Line Joystick Mount Std.
 Left Right
- Swing-away Joystick Mount E1028 \$210
 4" 7"
 Left Right

PROGRAMMER

- PG GC Programmer \$423

FRONT RIGGINGS

- Foot Platform (Depth-Adj.) Std.
- Rigging Type (pair)
Select footplate size for S/A & ELRs:
 71/4 x 6 81/2 x 6
- Swing-away footrests \$162
 11"-13" Ext. NC
 13"-16.5" Ext. Std.
 15.5"-19.5" Ext. NC
- Elevating Leg Rests E0990/K0195⁵ \$325
Adjustable Calf Pads w/ ELRS (1 hole std.)
 1 hole 3 holes* 6 holes*

*Potential for interference w/ shrouds

Extensions²

- 11"-14" Ext. NC
- 14"-17" Ext. Std.

Rigging Options

- Heel Loops (pair) E0951 \$38

SEATING

- Med-Back Gray Vinyl Std.
- Med-Back, Gray Vinyl & Solid Seat Pan³ NC
 18" W x 18" D
- Ltd. Recline, Gray Vinyl, High-Back w/ Headrest, **Depth Adjustable** \$275
 18" W x 18" - 20" D
 20" W x 18" - 20" D

SOLID SEAT PAN SELECTION

- Ltd. Recline, Gray Vinyl, High-Back w/ Headrest & Solid Seat Pan³ \$275
 18" W x 18" D
 18" W x 20" D
 20" W x 18" D
 20" W x 20" D

Custom seat dimensions are available.
Please contact the Custom Group for pricing and leadtime.

ARMRESTS

- For High-Back Seats**
- Flip-up, removable Std.
- Flip-up, height-adjustable E0973 \$170 removable
- Armpads**
- Desk Length, 10" NC
 Left NC Right NC
- Full Length, 14" ⁴ NC
 Left Std. Right Std.

Special Order Seat

_____ W x _____ D

Seat sizes can not exceed 20 W x 20 D
\$450 MSRP

Footnotes for page 1:

1. Can only mount one armrest accessory in combination w/ S/A joystick.
2. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
3. Includes black vinyl seat pan.
4. Full length arms required on Comfort & DepthAdjustable High-Back seating options.
5. E0990 Code used for purchase in CBA. K0195 code used for capped rental.

Jazzy Select

Order Form / Price List

ODJSF Only

Revised: February 2011

CUSHIONS⁴

<input type="checkbox"/> Simplicity Cushion		
16 x 16 - 20 x 20	E2601	\$98
<input type="checkbox"/> Solution 1 Cushion		
14 x 14 - 20 x 18	E2603	\$380
20 x 20	E2603	\$430
<input type="checkbox"/> Solution Cushion		
10 x 10 - 20 x 18	E2607	\$380
20 x 20	E2607	\$430
<input type="checkbox"/> Spectrum Cushion		
<input type="checkbox"/> Gel		
10 x 10 - 20 x 18	E2605	\$480
20 x 20	E2605	\$530
<input type="checkbox"/> Visco		
10 x 10 - 20 x 18	E2607	\$480
20 x 20 - 20 x 20	E2607	\$530
<input type="checkbox"/> Spectrum Air		
14 x 12 - 20 x 24	E2622	\$480
22 x 16 - 24 x 24	E2623	\$530
<input type="checkbox"/> Spectrum Air Contour		
14 x 12 - 20 x 24	E2624	\$499
22 x 16 - 24 x 24	E2625	\$549

ACCESSORIES

<input type="checkbox"/> Rear Basket ⁵		\$103
<input type="checkbox"/> Cane/Crutch Holder ⁵	E2207	\$98
<input type="checkbox"/> Oxygen Holder ⁵	E2208	\$103
<input type="checkbox"/> Walker Holder ⁵	K0108	\$92
<input type="checkbox"/> Weather Cover, Medium, Black		\$158
<input type="checkbox"/> Cup Holder		\$17
	<input type="checkbox"/> Right <input type="checkbox"/> Left	
<input type="checkbox"/> Lap Belt		Std.
<input type="checkbox"/> 38" <input type="checkbox"/> 50" <input type="checkbox"/> 60" <input type="checkbox"/> 70" <input type="checkbox"/> 80" <input type="checkbox"/> 90"		

SPECIALTY SEATINGSpecialty Seat Kit (limited adjustability) **\$275****SYSTEM SPECIFICATIONS**Specialty Seating Kit Width

<input type="checkbox"/> 16" NC	<input type="checkbox"/> 20" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 17" NC

Specialty Seating Kit Depth

<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	

Specialty Seating Kit Back Cane Height

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 20" NC
<input type="checkbox"/> 22" NC	

BACKS

- See Synergy Shape Back Order Form⁶
- Solid Plate Back up to 20"W⁷ **\$275**
- Omit Back
- See Positioning Components Order Form for additional backs

Back Height Solid Pan⁷

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Sling Back

<input type="checkbox"/> 16" wide	\$98
<input type="checkbox"/> 18" wide	\$98
<input type="checkbox"/> 20" wide	\$98

Sling Back Height

(Hgt must match cane hgt.)	
<input type="checkbox"/> 18" height	NC
<input type="checkbox"/> 20" height	NC
<input type="checkbox"/> 22" height	NC

ARMRESTS

- 2-Post, fixed flip-back^{8,9} Std.
- Flip-back, cane mounted^{8,9} \$95ea
- Removable, adjustable, quick height¹⁰ **E0973** \$275ea
- 11"-17" of height adj 3" on top in 1/2" increments & 3" on bottom in 1" increments
- Removable, height adjustable¹⁰ **E0973** \$150ea
- heavy duty armrests
- 8"H-12"H 10"H-14"H

Armrest Pads

- Full Length, straight, 2"x14" NC
- Left Right
- Desk Length, straight, 2"x10" NC
- Left Right
- Full Length, Waterfall, 2"x14"¹¹ NC
- Left Right
- Desk Length, Waterfall, 2"x10"^{3,11} NC
- Left Right

HEAD REST

- Stealth **E0955**
- Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
- Stealth Comfort Plus **E0955**
- Large (14"W) Small (10" W) \$205
- Removable Head Rest Bracket **E1028** \$210
- Stealth, swing-away (Flip-Down) Head Rest Bracket **E1028** \$210
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications **E0955**
- Large (11"W x 6" H)
- Small (9 3/4" W x 3" H) \$210
- Removable Head Rest Bracket **E1028** \$235

Footnotes for page 2:

- Can only mount one armrest accessory in combination w/ S/A joystick.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and/or cushion.
- S/A joystick mount N/A
- See Synergy cushion order form for sizes and pricing.
- Includes mounting bracket and hardware.
- Shape Back and w/ flip-back arms may interfere in the same configuration, call for details.
- 20" D is N/A.
- Call for availability with Shape back.
- S/A laterals interfere w/ flip-back arms.
- Desk length armrests w/ an accessory limit the accessibility to the height adjustment lever.
- N/A w/ flip back, cane mounted arms.

Jazzy Select

Order Form / Price List
ODJSF Only
Revised: February 2011

Synergy Seating Systems

Synergy Seating™ System \$434

WIDTH

Set Seat Width At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Width Package

14"W-18"W 16"W-20"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Depth Package

14"D-18"D 16"D-20"D

Please consult your Quantum Specialist or contact the Custom Group before ordering

BACK CANES¹

<input type="checkbox"/> Back Canes w/ push handles	Std.
<input type="checkbox"/> Back Canes w/out push handles	\$125
<input type="checkbox"/> Flip Forward ¹³	\$150
<input type="checkbox"/> Offset (Meso) \$195	<input type="checkbox"/> Inset (Endo) \$195
<input type="checkbox"/> 2" (1" per side) <input type="checkbox"/> 4" (2" per side)	

N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .
N/A w/ any reclining systems.

Actual Back Cane Height is reflected in measurement

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

BACK OPTIONS

- See Synergy Shape Back Order Form^{2,3}
- Solid Plate Back up to 20" W \$275
- Omit Back
- Back Height:⁴
 - 18"H 19"H 20"H 21"H
 - 22"H 23"H 24"H 25"H
 - Sling Back⁵ \$98

Back height will be the same as back cane height.
Includes gap for cushion.
Push Bar selection must match back width; center tube standard when selecting flip-up arms.

- TRU-Comfort Back ^{4,6}
 - Std. New Plus
 - Back Height: 22" 25"
 - Back up to 20" W **E2620** \$595
 - 22" W **E2621** \$695
 - TRU-Comfort Cushion
 - 16 x 16 - 20 x 18 **E2601** \$380
 - 20 x 20 - 20 x 22 **E2601** \$380
 - 22 x 16 - 22 x 24 **E2602** \$480
 - TRU-Comfort Plus Cushion
 - 16 x 16 - 20 x 18 **E2603** \$380
 - 20 x 20 - 20 x 22 **E2603** \$380
 - 22 x 16 - 22 x 24 **E2604** \$480
 - Omit Back -(\$100)

See TRU-Comfort Back & Cushion

order forms for specific sizes

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

ARMRESTS

- Flip-back, cane mounted ^{3,7} Std.
- 2-Post, flip-back, height ^{3,7,8} **E0973** \$150ea adjustable, removable armrests (9"-13" of height adjustment in 1" increments)
- Single Post height adj. ^{9,10} **E0973** \$75ea quick height, removable (11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.
 - Quad Release Lever NC (11"-17" of height adj in 3/8" increments)
- Removable, height adjustable **E0973** \$75ea heavy duty armrests⁹
 - 8"H-12"H 10"H-14"H
- Pediatric removable⁹ **E0973** \$75ea adjustable, quick height.
 - Match Shroud Color Black
 - Quad Release Lever NC (7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)
- Cantilever Armrests ¹⁴ SOF Angle-Adjustable: - 30° through + 30°

ARMPADS

- Full length, Straight, 2" X 14" Std.
 - Left Right
- Desk Length, Straight, 2" X 10" Std.
 - Left Right
- Full Length, Waterfall, 2" X 14" ¹¹ NC
 - Left Right
- Desk Length, Waterfall, 2" X 10" ^{11,12} NC
 - Left Right

HEAD REST⁵

- Stealth **E0955**
 - Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
 - Stealth Comfort Plus **E0955**
 - Large (14"W) Small (10" W) \$205
 - Removable Head Rest Bracket **E1028** \$210
 - Stealth, swing-away (Flip-Down)
 - Head Rest Bracket **E1028** \$210
- Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications **E0955**
 - Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230
 - Removable Head Rest Bracket **E1028** \$235
 - Fixed, Offset Mounting Bracket NC
 - TRU-Comfort **E0955**
 - 8" Wide 10" Wide \$205
 - Cobra Removable Head Rest Bracket **E1028** \$210
 - Removable Bracket for **E1028** \$210 after market headrest

Footnotes for page 3:

1. Height is measured from seat pan to top of cane.
2. 14" & 22" widths are N/A.
3. Shape Back & w/ flip-back arms may interfere in the same configuration call for details.
4. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
5. Headrests are N/A w/ sling backs.
6. N/A w/ flip back, cane mounted arms.
7. S/A laterals interfere w/ flip-back arms.
8. May cause interference with certain dimension configurations. Call Customer Service.
9. Desk length armrests w/ an accessory limit the accessibility to the height adjustment lever.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Available adj. range is 11" - 17".
11. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
12. S/A joystick mount N/A.
13. N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and waterfall armpads.
14. Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting aftermarket arm troughs or Universal Swivel w/ desk length arms.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.