

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Positioning Components Order Form - ODJFS Only

Account Number: _____ Date: _____ PO Number: _____
 Provider Name: _____ Marked for: _____
 Contact: _____ Ship to Address: _____
 Phone: _____ Fax: _____ City: _____ State: _____ Zip: _____
 Email: _____

This form is interactive when viewed as a PDF in most software. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. All prices are on-chair MSRP. Prices, specifications, part numbers, and availability are subject to change without notice.

STEALTH THORACIC LATERALS


Both a bracket and pad must be selected for a complete assembly.

Bracket

Select a bracket and a pad for complete assembly.

Swing Away, Rail Mount Bracket \$450 pr.
 QTY (TWBSTL)
 Part: INDPART2896 HCPCS: E1028
 For use with all TRU-Balance® 3 back options, Synergy Back, and TRU Comfort 2 on Synergy systems. May interfere with flip-up or cantilever armrests.

125° Swing Away Bracket..... \$450 pr.
 QTY (TWBTLTCL)
 Part: INDPART2677 HCPCS: E1028
 For use with Shape Backs or with TRU Comfort on TRU-Balance® 2 Recline systems. May interfere with flip-up or cantilever armrests.



Swing away, rail mount bracket

Lateral Extension

6" Lateral Extension..... \$25 ea.
 QTY Part: FRMASMB16274
 For back heights 13"-16" overall.

12" Lateral Extension \$25 ea.
 QTY Part: FRMASMB16298
 For back heights 19"-22" overall.

10" Lateral Extension \$25 ea.
 QTY Part: FRMASMB16280
 For back heights 17"-18" overall.

16" Lateral Extension \$25 ea.
 QTY Part: FRMASMB16281
 For back heights 23"-25" overall.



6" Lateral Extension

Pads

4"H x 6"D Flat Pads with \$200 pr.
 QTY Leatherette Cover (TWBLP4x6)
 Part: INDPART2675 HCPCS: E0956

4"Hx6"D Curved Pads with..... \$200 pr.
 QTY Dartex Cover (TWBL4x6)
 Part: INDPART2664 HCPCS: E0956

3"H x 5"D Flat Pads with \$200 pr.
 QTY Leatherette Cover (TWBLP3x5)
 Part: INDPART2676 HCPCS: E0956

3"Hx5"D Curved Pads with..... \$200 pr.
 QTY Dartex Cover (TWBL3x5)
 Part: INDPART2663 HCPCS: E0956



Lateral Pad

STEALTH PELVIC/THIGH LATERALS

Bracket

Brackets may interfere with some armrests configurations. Select a bracket and a pad for complete assembly.

Flip Down Bracket.....\$225 ea.
 (TWBSTLADD-HW)
 Part: INDPART2804 HCPCS: E1028
 For use with all TRU-Balance 3 back options, Synergy Back, and TRU Comfort 2 on Synergy systems. May interfere with flip-up or cantilever armrests.

Removable Bracket (LP430).....\$250 ea.
 Part: INDPART2805 HCPCS: E1028
 Height adapter plate is needed with TRU-Comfort seat cushions.


Removable Bracket with.....\$325 ea.
 1" Offset (LP430 & LHW-121)
 Part: INDPART2854 HCPCS: E1028

Left, Flip Down w/ Full Surface\$225 ea.
 Contact Feature (TWBADD-L & LPHW407)
 Part: POS140428 HCPCS: E1028

Right, Flip Down w/ Full Surface\$225 ea.
 Contact Feature (TWBADD-R & LPHW407)
 Part: POS140429 HCPCS: E1028




Left, Adjustable,
Flip-Down Pelvic/Thigh
Lateral



Removable Pelvic/Thigh
Bracket



Adapter Link



2" Depth Adapter

Adapters

Recommend with a 4" thick or greater seat cushion.

Adapter Link, 2 1/4" Height..... \$25 ea.
 (LHW-129)
 Part: POS140430 HCPCS: K0108
 For use with Adjustable, Flip Down brackets only.

2" Depth/ Height Adapter Plate..... \$25 ea.
 (LHW-122)
 Part: INDPART2821 HCPCS: K0108
 For use with brackets other than the Adjustable Flip-down.

Pads

4"Hx4"D Short (TWBADD-4x4) \$50 ea.
 Part: INDPART2801 HCPCS: E0956

4"Hx8"D Medium (TWBADD-4x8) \$50 ea.
 Part: INDPART2802 HCPCS: E0956

4"Hx12"D Long (TWBADD-4x12) \$50 ea.
 Part: INDPART2803 HCPCS: E0956

SHOULDER HARNESSES

Bodypoint

Shoulder Harness, Rear Pull, Medium, 13"-16" \$124 ea.
 (SH220M)
 Part: INDPART2744 HCPCS: E0960



Shoulder Harness

CHEST HARNESSES/ STRAPS

i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea.
 Part: ST-IFPS1103-Q HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
 Part: ST-IFPS1104-Q HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$142 ea.
 Part: ST-IFPS2103-Q HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
 Part: ST-IFPS2104-Q HCPCS: E0960



Contoured and Classic
Chest Harness

i-Fit Chest Straps

Chest Strap, Medium w/ Belt Anchor \$90 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q HCPCS: E0960



Chest Strap

PELVIC BELTS

i-Fit

60" Single Pull, Push Button, \$86 ea.
 Medium (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q HCPCS: E0978

60" Single Pull, Push Button, \$86 ea.
 Large Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q HCPCS: E0978

60" Single Pull, Push Button, \$86 ea.
 Extra Large Padded (11.5"), 2" Belt
 (IFPB5315)
 Part: ST-IFPB5315-Q HCPCS: E0978



i-Fit Padded Belt

PELVIC BELTS CONTINUED

Innovative Concepts and Bodypoint

Innovative Concepts: Lap Belt, Standard Pull,\$84 ea.
60", Padded, Push Button (PRI-11041-22)

Part: ACCBELT1032

HCPCS: E0978

Bodypoint: Pelvic Positioning Belt,\$87 ea.
Padded, Push Button, Medium, 9"-15" Hip Width
(HB205M-A1)

Part: ACCBELT1009

HCPCS: E0978

Bodypoint: Pelvic Positioning Belt,\$92 ea.
Padded, Push Button, Large, 13"-19" Hip Width
(HB205L-A1)

Part: ACCBELT1010

HCPCS: E0978

KNEE ABDUCTOR

Therafin

Hardware

Flip Down Abductor Hardware (50180) .. \$172 ea.

Part: INDPART2755

HCPCS: E1028

Pads

Small, 3.5"W to 3.25"W x 4"L (50285).... \$135 ea.

Part: INDPART2751

HCPCS: E0957

Medium, 4"W to 3.25"W x 6"L (50286) ... \$135 ea.

Part: INDPART2752

HCPCS: E0957

Small, Gel 2.25"W x 4.25"L (50288) \$145 ea.

Part: INDPART2754

HCPCS: E0957

Large, 4.5"W to 2.75"W x 6"L (50287) ... \$135 ea.

Part: INDPART2753

HCPCS: E0957



KNEE ADDUCTOR

Therafin Adductors

Hardware

Left, Flip Down Adductor, Long \$308.82 ea.
Arm with Clamp and Rail Mount (32944)

Part: INDPART2912

HCPCS: E1028

Right, Flip Down Adductor, Long..... \$308.82 ea.
Arm with Clamp and Rail Mount (32945)

Part: INDPART2913

HCPCS: E1028

Pads

3" x 4" Small Pad with Cover (50490) \$105 ea.

Part: INDPART2914

HCPCS: E0956

4" x 6" Med. Pad with Cover (50491) \$105 ea.

Part: INDPART2915

HCPCS: E0956



AEL Adductors

Posalinc, 1" Tubing Mount \$512 pr.
(13515LX/ 13515RX)

Part: INDPART2433

HCPCS: K0108

Available with swing away legrests only.

Posalinc, Knob Release, Rail Mount \$574 pr.
(13555LX/ 13555RX)

Part: INDPART2825

HCPCS: K0108

May interfere with some armrest configurations.



AMPUTEE SUPPORT

Therafin

Support

Multi-Axis Support, 7" x 10" (31515) \$230 ea.

Part: INDPART2822

HCPCS: E1020

Swing away upper legrest is required.

Optional Accessory

Amputee Support Gel Cover (32391).....\$135 ea.

Part: INDPART2823

HCPCS: K0108



BARIATRIC CALF SUPPORTS

Bodypoint

Bariatric Calf Panel, Small, 22"W\$135 ea.
Fits 20"W-24"W (CM165S)

Part: FRMASMB17538

HCPCS: E0995

Bariatric Calf Panel, Large, 30"W\$135 ea.
Fits 28"W-32"W (CM165L)

Part: FRMASMB17540

HCPCS: E0995

Bariatric Calf Panel, Medium, 26"W\$135 ea.
Fits 24"W-28"W (CM165M)

Part: FRMASMB17539

HCPCS: E0995

CENTER MOUNT CALF PAD/ FOOTPLATE COVERS

Innovative Concepts

Pad Covers

Calf Pad Covers, Air Mesh, Small \$120 pr.
(PRI-30005)
Part: RIG133716

Calf Pad Covers, Air Mesh, Large \$120 pr.
(PRI-30022)
Part: RIG133718

Calf Pad Covers, Reversed Dartex, \$120 pr.
Small (PRI-30006)
Part: RIG133717

Calf Pad Covers, Reversed Dartex, \$120 pr.
Large (PRI-30023)
Part: RIG133719



Calf Pad Covers

Calf Panel Covers

2 pairs of the side covers are needed with the large calf panel.

Calf Panel Cover, Air Mesh \$120 ea.
(PRI-30042)
Part: RIG133722

Calf Panel Side Covers, Air Mesh \$120 pr.
(PRI-30055)
Part: RIG133724

Calf Panel Cover, Reversed Dartex \$120 ea.
(PRI-30043)
Part: RIG133723

Calf Panel Side Covers, Reversed \$120 pr.
Dartex (PRI-30056)
Part: RIG133725



Calf Panel Covers

Footplate Covers

Center Mount Footplate Covers, \$120 pr.
Reversed Dartex (PRI-30036)
Part: RIG133721

LEG STRAPS

AES

18"W, Single Leg Strap \$48 ea.
(PLS1000-18)
Part: INDPART2289 HCPCS: K0038

20"W, Single Leg Strap \$48 ea.
(PLS1000-20)
Part: INDPART2288 HCPCS: K0038



Leg Strap

TOE/ ANKLE STRAPS

Innovative Concepts

Toe/Ankle Strap, Adult (PRI-11006) \$50 pr.
Part: INDPART2942 HCPCS: E0952

Toe/Ankle Strap, Adult, Padded \$60 pr.
(PRI-110003)
Part: INDPART2944 HCPCS: E0952



Toe/Ankle Strap

AES

Medium, Toe/Ankle Strap (TST1500) \$24 pr.
Part: INDPART2293 HCPCS: E0952

Large, Toe/Ankle Strap (TST2000) \$24 pr.
Part: INDPART2294 HCPCS: E0952

STEALTH FOOT BOXES

Small, Single Leg Foot Boxes, \$475 pr.
5"Wx8"Dx7"H (SPLE107)
Part: INDPART2786 HCPCS: K0108

Medium, Single Leg Foot Boxes, \$475 pr.
6"Wx10"Dx9"H (SPLE108)
Part: INDPART2787 HCPCS: K0108

Large, Single Leg Foot Boxes, \$475 pr.
7"Wx12"Dx11"H (SPLE109)
Part: INDPART2788 HCPCS: K0108



Stealth Foot Boxes

STEALTH SHOE HOLDER

Small Shoe Holders, Padded Buckle \$266 pr.
3.75"Wx6.25"D (SPLE111)
Part: INDPART2783 HCPCS: K0108

Medium Shoe Holders, Padded \$266 pr.
Buckle, 4"Wx8.25"D (SPLE112)
Part: INDPART2784 HCPCS: K0108

X-Large Shoe Holders, Padded \$266 pr.
Buckle, 5"Wx10.75"D (SPLE114)
Part: INDPART2785 HCPCS: K0108



Stealth Shoe Holders

ANKLE HUGGERS

Innovative Concepts

Medium Ankle Huggers, Side \$155 pr.
Release Strap (PRI-110452-M)
Part: INDPART2938 HCPCS: K0108

Large Ankle Huggers, Side Release \$155 pr.
Strap (PRI-110452-L)
Part: INDPART2939 HCPCS: K0108



Bodypoint

Medium Ankle Huggers, Side \$134 pr.
Release Buckle (FT220M)
Part: INDPART2416 HCPCS: K0108

Large Ankle Huggers, Side Release \$134 pr.
Buckle (FT220L)
Part: INDPART2748 HCPCS: K0108

HIGH MOUNT CLAMP ON FOOTRESTS

Bodypoint

Angle Adjustable Footplates, Small, \$356 pr.
4"W x 7.75"D, Fits 1" Tubing
(FT101/ HDW406-2)
Part: FRMASMB3257 HCPCS: K0040

Fulcrum Series, Multi-Axis, Large, \$460 pr.
4.75"W x 9.25"D, Fits 1" Tubing
(FT106/HDW406-2)
Part: FRMASMB3572 HCPCS: K0040



Angle Adjustable Footplates, Large, \$356 pr.
4.75"W x 9.25"D, Fits 1" Tubing
(FT102/ HDW406-2)
Part: FRMASMB3258 HCPCS: K0040

ARM PADS

Falcon

Gel Armpad, Full Length, \$80 ea.
3.5"W x 14"L (314G)
Part: ARM136408

Gel Armpad, Desk Length, \$75 ea.
3.5"W x 10"L (310G)
Part: ARM136409



Motion Concepts

Flat Armpad, Full Length, \$85 ea.
3.5"W x 14"L (TRM0606)
Part: SETLRAM1167

ARM CHANNELS

Otto Bock

Armpad mounting hardware is required. Refer to the armrest hardware section on following page.

Medium Forearm Channel, \$104 ea.
4.75"W x 11.5"L (436A2=2-7)
Part: INDPART2123 HCPCS: E2209

One Piece Arm Channel, \$104 ea.
4.75"W x 18"L (436A3=2-7)
Part: INDPART2125 HCPCS: E2209



Large Forearm Channel, \$104 ea.
4.75"W x 14"L (436A2=1-7)
Part: INDPART2124 HCPCS: E2209



Motion Concepts

Motion concepts armpad mounting hardware is required. Refer to the armrest hardware section on following page.

Left Arm Trough (EAT-L) \$100 ea.
Part: TRXARM1053 HCPCS: E2209

Right Arm Trough (EAT-R) \$100 ea.
Part: TRXARM1054 HCPCS: E2209

HAND PADS

Otto Bock Flat Pad

Medium Flat Hand Pad, 4.5"W x 7"L \$52 ea.
(436P5=2-7)
Part: INDPART2127

Large Flat Hand Pad, 5.5"W x 9"L \$52 ea.
(436P5=1-7)
Part: INDPART2126



HAND PADS CONTINUED

Otto Bock Palm Extensor

Left, Medium, Palm Extensor, \$52 ea.
 QTY 5.75"W x 8"L (436P2=L2-7)
 Part: INDPART2130

Right, Medium, Palm Extensor, \$52 ea.
 QTY 5.75"W x 8"L (436P2=R2-7)
 Part: INDPART2131

Left, Large, Palm Extensor, \$52 ea.
 QTY 5.75"W x 9.25"L (436P2=L1-7)
 Part: INDPART2128

Right, Large, Palm Extensor, \$52 ea.
 QTY 5.75"W x 9.25"L (436P2=R1-7)
 Part: INDPART2129



Motion Concepts

Flat Hand Pad Left or Right \$85 ea.
 QTY (TRM0325)
 Part: TRXARM1055



ARMREST HARDWARE

Otto Bock Accessories


Armrest and Handpad Strap \$30 ea.
 QTY Package (436S1=1-7)
 Part: INDPART2816 HCPCS: K0108

Elevating Swivel Unit, Black \$182 ea.
 QTY (476E52=SK035)
 Part: INDPART2633 HCPCS: K0108



Motion Concepts

Arm Trough Mounting Hardware \$38 ea.
 QTY (EMH)
 Part: HDWASMB1096



Arm Trough/ Gel Armpad Mounting Hardware

Desk Length Arm Mounting \$56.20 ea.
 QTY Hardware
 Part: HDWASMB1188

Full Length Arm Mounting \$56.52 ea.
 QTY Hardware
 Part: HDWASMB1189

Desk Length Recline or Single \$56.20 ea.
 QTY Post Armrest Mounting Hardware
 Part: HDWASMB1190

Full Length Recline or Single \$56.52 ea.
 QTY Post Armrest Mounting Hardware
 Part: HDWASMB1191

TRU-Balance 3 Hand Pad or Arm \$5.92 ea.
 QTY Channel Mounting Hardware
 Part: HDW131075

STEALTH ELBOW SUPPORT

Select a bracket and pad for complete assembly.


Bracket

Swing-away Elbow Support, Left \$175 ea.
 Tube-Style Bracket (A) (ARMS-ES-L)
 Part: POS124064 HCPCS: E1028

Fixed Elbow Support, Left \$175 ea.
 Track-Style Bracket (B) (ARMS-ESFTB3-L)
 Part: POS133465

Swing-away Elbow Support, Right \$175 ea.
 Tube-Style Bracket (ARMS-ES-R)
 Part: POS124066 HCPCS: E1028

Fixed Elbow Support, Right \$175 ea.
 Track-Style Bracket (B) (ARMS-ESFTB3-R)
 Part: POS133466



Pad

Elbow Support Pad, 3"x3" \$50 ea.
 QTY (ARMS-ES 3x3)
 Part: POS124065 HCPCS: K0108

Elbow Support Pad, 4"x4" \$50 ea.
 QTY (ARMS-ES 4x4)
 Part: POS124067 HCPCS: K0108

HEADRESTS AND MOUNTS

Stealth Headrests


Comfort Plus Headrest Pads

8" Comfort Plus (CP450) \$142 ea.
 Part: POS148484 HCPCS: E0955

10" Comfort Plus (CP250)(A) \$142 ea.
 Part: SETHEAD1057 HCPCS: E0955

14" Comfort Plus (CP150) \$142 ea.
 Part: INDPART2327 HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) \$142 ea.
 (CPN350)
 Part: SETHEAD1069 HCPCS: E0955



HEADRESTS AND MOUNTS CONTINUED

Contoured Headrest Pads


Small Contoured (9.75"Wx3"H)..... \$60 ea.
(STL-OTSM250)
Part: SETHEAD1052 HCPCS: E0955

Large Contoured (11"Wx6"H) \$75 ea.
(STL-OTLG150) (B)
Part: SETHEAD1053 HCPCS: E0955

Headrest Bracket

Stealth's Worlds Best Removable.....\$210 ea.
Headrest Bracket (TWB485)(C)
Part: FRMASMB15999 HCPCS: E1028

Removable Headrest Bracket\$150 ea.
(HMO475P)(D)
Part: HDWASMB1084 HCPCS: E1028



Stealth Swing Away Headrest Mounts

These Swing-Away headrest mounts are only available with the Stealth Comfort Series Headrests.

Facial Lateral, 4.5"x2" Pad,.....\$325 ea.
Right Side, for Comfort Series
(CPS9-R-HW & 900-Q)
Part: INDPART2647

Facial Laterals, 4.5"x2" Pad, Left \$425 ea.
and Right Side for Comfort Series
(CPS9-B-HW & 900-Q)
Part: INDPART2649

Facial Lateral, 4.5"x2" Pad,.....\$325 ea.
Left Side for Comfort Series
(CPS9-L-HW & 900-Q)
Part: INDPART2648



Other Accessories


Egg Switch Mount (ESM) (A) \$50 ea.
Part: INDPART2815

Swing-Away Joystick Mount,\$250 ea.
18" Rod (C)
Part: ST-SUS18-Q

Figure 8 Bracket with 13" Rod.....\$250 ea.
and Egg Switch Mount (B)
Part: ST-F8678-13-Q

Allows mounting of facial pads and attaches to the headrest

Facial Pad, 6"x2" (920-Q) \$50 ea.
Part: ST-920-Q HCPCS: K0108



LAP TRAYS

Therafin

These trays are compatible with rehab style armrests only. Reference the subsequent charts to help determine lap tray compatibility.

Medium, Clear Lap Tray, Top Drop.....\$352 ea.
Clamps, PVC Rim, Right Side Controller
Cutout, (T15R12A)
Part: FRMASMB11233 HCPCS: E0950
Fits seats measuring 18" to 20.5" armpad to armpad measurement.

X-Large, Clear Lap Tray, Top Drop\$399 ea.
Clamps, PVC Rim, Right Side Controller
Cutout, (T17R12A)
Part: INDPART2644 HCPCS: E0950
Fits seats measuring 22" to 24.5" armpad to armpad measurement.

Large, Clear Lap Tray, Top Drop.....\$378 ea.
Clamps, PVC Rim, Right Side Controller
Cutout, (T16R12A)
Part: FRMASMB11232 HCPCS: E0950
Fits seats measuring 20" to 22.5" armpad to armpad measurement.



LAP TRAYS CONTINUED

For Synergy and TRU-Balance 2, add the number below to the seat width to determine the overall armrest to armrest width. For TRU-Balance 3 armrests, add the measurement below to the back width.

			TRU-Balance 2 Tilt	TRU-Balance 2 with TRU-Comfort	TRU-Balance 2 Tilt and Recline	TRU-Balance 2 Tilt and Recline with TRU Comfort
Armrest and Arm Pad	Cane Mounted	Std/ Gel	5.5"	--	--	--
		Waterfall	5.25"	--	--	--
	Single Post	Std/ Gel	4.5"	4.5"	--	--
		Waterfall	4.25"	4.55"	--	--
	Dual Post	Std/ Gel	5.5"	5.5"	--	--
		Waterfall	5.25"	5.5"	--	--
	Cantilever	Std/ Gel	6.25"	7.25"	7.75"	7.75"
		Waterfall	6"	7"	7.5"	7.5"
	TRU-Balance 2 Recline	Std/ Gel	--	--	6.5"	6.5"
		Waterfall	--	--	6.25"	6.25"
TRU-Balance 2 Armpad to Armpad Measurements						

			Synergy	Synergy with TRU-Comfort
Armrest and Arm Pad	Cane Mounted	Std/ Gel	5.5”	--
		Waterfall	5.25”	--
	Single Post	Std/ Gel	4.5”	4.5”
		Waterfall	4.25”	4.55”
	Dual Post	Std/ Gel	5.5”	5.5”
		Waterfall	5.25”	5.5”
	Cantilever	Std/ Gel	6.25”	7.25”
		Waterfall	6”	7”
Synergy Seat Armpad to Armpad Measurements				

		TRU-Balance 3
Dual Post Armrest	Std/ Gel	5.5"
	Waterfall	5.25"
TRU-Balance 3 Armrests	Std/ Gel	4.5"
	Waterfall	4.25"
TRU-Balance 3 Armpad to Armpad Measurements		

ROHO SEAT CUSHION

Quadro Select, High Profile, \$454 ea.
16"X x 18"D (QS910C)
Part: SETBASE2569 HCPCS: E2624

Quadro Select, High Profile, \$454 ea.
20"X x 18"D (QS1110C)
Part: SETBASE2570 HCPCS: E2624

Quadro Select, High Profile, \$454 ea.
18"X x 18"D (QS1010C)
Part: SETBASE2555 HCPCS: E2624



PLANAR MOUNTING KIT

Planar Mounting Kit \$75 ea.
Part: SOF

Used to mount aftermarket seating to Synergy or TRU-Balance seating frames.

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2016 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

