

Quantum® 610 Series

Group 2 Standard
Order Form/ODJFS Only
300 lbs. weight capacity



ACCOUNT INFO

Account Number: _____
 Provider Name: _____
 Contact Name: _____
 Phone #: _____
 Fax #: _____
 E-mail Address: _____

Date: _____
 PO#: _____
 Marked For: _____
 3 Column Pricing on Quote/Order: Yes No
 Quote #: _____

Ship To: _____
 Address: _____
 City: _____ State: _____ Zip: _____
 Lift Gate Needed: Yes No

Taken By: _____

SET-TO-SPEC*

Please Supply Client Information:

Client Height: _____
 Client Weight: _____

Mount Hip/Thigh Supports: Yes No
 Mount Lateral Thoracic Supports: Yes No

Armrest Height from Pan: _____

Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

Seating System	Seat-To-Floor Height	
	Contoured Seat	20.5"
Synergy Seat	16.5"	17.5"

*Seat to floor height are +/- 1/8" with 18"W x 18"D with solid tires (add 3/8" for pneumatic tires).
 Increased widths may increase lowest seat to floor heights.*

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.

NOTES:

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASE MODEL		
Description	Price	HCPCS Code
<input type="checkbox"/> Quantum 610 2S-SS	\$5,971	K0822/K0011
<input type="checkbox"/> Quantum 610 2S-C	\$5,971	K0823K0011
<input type="checkbox"/> Quantum 610 2SP-C	\$6,171	K0836/K0011
<input type="checkbox"/> Quantum 610 2SP-SS	\$6,171	K0835/K0011

COLORS	
Colors:	
<input type="checkbox"/> Candy Apply Red	<input type="checkbox"/> Sunburst Orange
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Pearl Pink
<input type="checkbox"/> Onyx Black	



BASE MOUNTED FOOT PLATFORM	
Description	Price
<input type="checkbox"/> Base Mounted Foot Platform Installed	Standard
<input type="checkbox"/> No Base Mounted Foot Platform	No Charge
<input type="checkbox"/> Bell Shaped Foot Platform	\$235

DRIVE WHEELS	
Description	Price
<input type="checkbox"/> 10" Flat-Free (pair)	Standard
<input type="checkbox"/> 10" Pneumatic (pair)	No Charge

BATTERY CHARGER	
Description	Price
<input type="checkbox"/> 5 Amp Off-board Charger	Standard





POWER ELEVATING SEAT		
Description	Price	HCPCS Code
<input type="checkbox"/> 5" Power Elevating Seat Thru VSI 6-Key	\$975	E2300

NON-EXPANDABLE CONTROLS	
Description	Price

<input type="checkbox"/> PG 4-Key VR2	Standard	
<input type="checkbox"/> PG 6-Key VR2	\$200	
<input type="checkbox"/> PG 6-Key VSI <i>For Use With Power Elevate</i>	\$200	



EXPANDABLE CONTROLS	
Description	Price

Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls

<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller	\$695	
<input type="checkbox"/> Expandable Controller Harness	\$425	
<input type="checkbox"/> Q-Logic 2 EX Joystick	\$700	
<input type="checkbox"/> Q-Logic 2 EX Joystick for Lights	\$700	
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port	\$700	

Armor Color Options:	
<input type="checkbox"/> Black Soft Touch (Standard color option)	
Colors: \$50	
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Sunburst Orange
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Pearl Pink

PROGRAMMING	
Description	Price

<input type="checkbox"/> Q-Logic Handheld Programmer	\$545	
<input type="checkbox"/> Q-Logic Programming Cable w/ Disk ¹	\$750	
<input type="checkbox"/> PG Handheld Programmer	\$423	

1. Works with Q-Logic 2 EX only

JOYSTICK HANDLES		
Description	Price	HCPCS Code

<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup	\$71	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$69	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

JOYSTICK MOUNTING BRACKETS		
Description	Price	HCPCS Code

<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline ¹	Standard	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-away	\$245	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
Offset: <input type="checkbox"/> Inline <input type="checkbox"/> Offset		
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$250	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Multi-axis	\$100	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Not available with Q-Logic 2 EX Stand Alone Joystick

CONTOURED SEATING

SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
Depth Adjustable High Back:		
<input type="checkbox"/> Limited Recline Depth Adjustable High Back w/ Headrest ¹	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 18"W x 18-20"D <input type="checkbox"/> 20"W x 18-20"D		

High Back with Solid Seat Pan:		
<input type="checkbox"/> Limited Recline High-Back with Headrest & Solid Seat Pan	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 16"W x 20"D		
<input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 18"W x 18"D <input type="checkbox"/> 18"W x 20"D		
<input type="checkbox"/> 20"W x 16"D <input type="checkbox"/> 20"W x 18"D <input type="checkbox"/> 20"W x 20"D		

Armrests:		
<input type="checkbox"/> Flip-Up, Height-Adjustable, Removable	\$85 ea.	E0973

Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length ²		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length ²		

Special Order Seat

_____ "W x _____ "D

Seat size can not exceed 20"W x 20"D

Contact your inside sales specialist for pricing

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part.
2. Selection of desk length armpads on a depth-adjustable high back requires them to be added as an additional part.

MANUAL RECLINER

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back with Headrest	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back with Headrest and Solid Seat Pan ¹	\$600	E1226

Upholstery:		
<input type="checkbox"/> Gray Vinyl		
Seat Width:		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
Seat Depth:		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		

1. Includes black vinyl seatpan. Back available in listed colors.

SYNERGY® SEATING SYSTEMS

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy Seating Frame	\$434	
<input type="checkbox"/> Synergy Manual Recline Package ¹	\$700	E1226

Seat Width:					
<input type="checkbox"/> 10" ² <input type="checkbox"/> 11" ³ <input type="checkbox"/> 12" ³ <input type="checkbox"/> 13" ³ <input type="checkbox"/> 14"					
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input checked="" type="checkbox"/> 19" <input checked="" type="checkbox"/> 20"					

Seat Width Package:		
<input type="checkbox"/> 10"W-14"W <input type="checkbox"/> 14"W-18"W <input type="checkbox"/> 16"W-20"W		

Seat Depth:				
<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input checked="" type="checkbox"/> 14"				
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input checked="" type="checkbox"/> 19" <input checked="" type="checkbox"/> 20"				

Seat Depth Package:		
<input type="checkbox"/> 10"D-14"D <input type="checkbox"/> 14"D-18"D <input type="checkbox"/> 15"D-20"D		

Back Cane Height:				
<input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" ³ <input type="checkbox"/> 17" ⁴				
<input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21"				
<input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24" <input type="checkbox"/> 25"				

1. Not available with Flip Back Arms, Cantilever Arms, Shape Back, AFP or Power Elevate. Requires 22"H canes. Pivot point 2.5" above seat pan.
2. SOF width is required when choosing a power elevating seat on widths 10"-13".
3. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

Description	Price	HCPCS Code
<input type="checkbox"/> Sling Back	\$98	
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Back: <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Choose Lateral Extensions ¹ :		
<input type="checkbox"/> 6" for 13"-16" H Back	\$25	
<input type="checkbox"/> 10" for 17"-18" H Back	\$25	
<input type="checkbox"/> 12" For 19"-22" H Back	\$25	
<input type="checkbox"/> 16" For 23"-25" H Back	\$25	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width: <input type="checkbox"/> 16" W <input type="checkbox"/> 18" W <input type="checkbox"/> 20" W	\$595	E2620
Overall Back: <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
Height ² : <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 17" W <input type="checkbox"/> 19" W	\$595	E2620
Overall Back: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Height ² : <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 22"W ³	\$695	E2621
Overall Back: <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Height ² :		

Synergy® Back Cane Options:

<input type="checkbox"/> Standard Back Canes	Standard
<input type="checkbox"/> Omit Push Handles ⁴	\$125
<input type="checkbox"/> Flip Forward Back Canes ⁵	\$150
<input type="checkbox"/> Omit Push Handles	\$125
<input type="checkbox"/> Offset (Meso) ⁶	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	
<input type="checkbox"/> Inset (Endo)	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	

1. Lateral extensions are required to mount lateral thoracic supports.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
3. Available with offset (meso) brackets.
4. Not available on 14"H canes.
5. Not available with Vent Tray, O2 Holder, cane mounted or cantilever armrests. Also not available w/ a Synergy manual recline or manual tilt. May interfere with Lateral Supports and Waterfall Armpads.
6. Dual Post Flip-back Armrests interfere with offset canes.

SYNERGY® ARMRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> Flip-back, Cane Mounted ¹	No Charge	
<input type="checkbox"/> 2-post, Flip-back, Height-adjustable, Removable ^{2,3}	\$150 ea.	E0973
<input type="checkbox"/> Single-post, Height-adjustable, Removable ^{3,4} <input type="checkbox"/> Quad Release Lever	\$75 ea.	E0973
<input type="checkbox"/> Removable, Single-post, Height Adjustable, Heavy Duty <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H	\$75 ea.	E0973
<input type="checkbox"/> Pediatric, Removable ⁵ <input type="checkbox"/> Quad Release Lever	\$75 ea.	E0973
<input type="checkbox"/> Bariatric, Removable, Height Adjustable ⁶ <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H	\$75 ea.	E0973
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁷ Angle-Adjustable: -30° through +30°	\$195	

1. Interferes with swing-away laterals. Not available on flip forward canes or manual recline.
2. Interferes with swing-away laterals. Height adjustable from 9"-13" in 1" increments.
3. Not available with a depth less than 16".
4. Height adjustable from 11"-17" in 1/2" increments
5. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available
6. Bariatric Armrests require Bariatric Armpads
7. Height adjustment from 11-17" depends on cane height. Not available with Manual Recline, Flip Forward back canes & Shape Back.

SYNERGY® ARMPADS

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14") <input type="checkbox"/> Left <input type="checkbox"/> Right	Standard
<input type="checkbox"/> Straight, Desk Length (2"x10") <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Full Length (2"x14") ¹ <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Waterfall, Desk Length (2"x10") ¹ <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")* <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$92.50 ea. \$56.69 ea.
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")* <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$70.00 ea. \$53.52 ea.
<input type="checkbox"/> Flat, Bariatric ² <input type="checkbox"/> Left <input type="checkbox"/> Right	No Charge

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.
2. Bariatric Armrests require Bariatric Armpads

LAPBELTS

Description	Price
<input type="checkbox"/> 38" Lap Belt ¹	Standard
<input type="checkbox"/> 50" Lap Belt ¹	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Non-Retractable (Sure-Lok) ²	Standard
<input type="checkbox"/> 60" Padded Lap Belt ³	\$115
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$87
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width ³	\$92
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width*	\$112
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width*	\$116

1. Includes mounting hardware. Not available with Manual Recliner.
2. Only needed with the occupied transit kit.
3. Only available with Synergy seating.

I-FIT PELVIC BELTS

Description	Price
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded 1.5" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Large (9.5") Padded 2" Belt	\$86
<input type="checkbox"/> i-Fit 60" Single Pull, Push Button, Extra Large (11.5") Padded 2" Belt	\$86

HEADRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$60	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$75	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$150	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350) Used in Addition to Removable Headrest Bracket	\$210	
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

NON-POWERED LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 60-90° Swing-away Legrests¹ (Adjustable in 5° increments)	\$162	
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
<i>Standard Extruded Extensions:</i>		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Standard Extruded Footplates:</i>		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ²	\$75 ea.	
<i>Angle Adjustable Extensions:</i>		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
<i>Angle-Adjustable Footplates:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		
<input type="checkbox"/> High Mount, Clamp-On Footrest*		
<input type="checkbox"/> Angle Adjustable		
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> Heavy Duty Drop-In Legrests⁵	\$420	
<input type="checkbox"/> Angle Adjustable Footplates ²	\$75 ea.	
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		

1. Possible caster wheel interference when adjusted to less than 70 degrees.
2. Includes Heel Loops
3. Causes interference on seat widths less than 18".
4. Causes interference on seat widths less than 20".
5. Only available with Synergy seating.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

NON-POWERED LEGRESTS (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> Elevating Legrests¹	\$325	K0195
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length:		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates ²	\$75 ea.	
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests¹	\$162	
<input type="checkbox"/> Angle Adjustable Footplates ²		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		

1. Extensions are measured from the seatpan to footplate.
 2. Includes heel loops.
 3. Causes interference on seat widths less than 18".
 4. Causes interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM

Description	Price
<input type="checkbox"/> Adjustable Center Mount Foot Platform¹	\$250
Footplates:	
<input type="checkbox"/> 10"W 11"D Split Footplates	
<input type="checkbox"/> 10"W 8"D Split Footplates	
<input type="checkbox"/> 10"W 6"D Split Footplates	
Extensions: ²	
<input type="checkbox"/> 8"-12" Ext. <input type="checkbox"/> 10"-14" Ext.	
<input type="checkbox"/> 10"-16" Ext.	
Calf Supports ³ :	
<input type="checkbox"/> Calf Pads (pair)	
<input type="checkbox"/> Small (6"H x 6"W)	\$150
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150
<input type="checkbox"/> Calf Panel (ea.)	
<input type="checkbox"/> Small (8"H x 18"W)	\$150
<input type="checkbox"/> Large (8"H x 24"W) ⁴	\$200
Optional Modular Footplate Extension Components/Kits ^{3,5}	
<input type="checkbox"/> 4" Width Extension Kit	\$50
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/> 4" Side Guard ⁶	\$50
<input type="checkbox"/> 6" Side Guard ⁷	\$50
Additional Optional Accessories ³ :	
<input type="checkbox"/> Footplate Connector	\$15
<input type="checkbox"/> Wheel Bumpers, pair	\$25

1. Only available with Synergy seating.
 2. Extension length is measured from top of the seat pan to footplate.
 3. Choose in Rigging Accessory Section in iQuantum.
 4. Only available on 20"W or greater.
 5. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 6. Only available with 10"W x 6"D Split Footplates
 7. Only available with 10"W x 11"D and 10"W x 8"D Split Footplates

ELECTRONICS FOR POWER FOOT RIGGINGS

Description	Price	HCPCS Code
<input type="checkbox"/> Through VR2 6-Key	No Charge	
<input type="checkbox"/> Through Q-Logic 2 EX	\$1,895	E2310

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Power Articulating Foot Platform¹	\$2995	See note ²
Footplates:		
<input type="checkbox"/> 10"W 11"D Split Footplates		
<input type="checkbox"/> 10"W 8"D Split Footplates		
<input type="checkbox"/> 10"W 6"D Split Footplates		
Extensions ³ :		
<input type="checkbox"/> 10"-16" Ext.		
Calf Supports ⁴ :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ⁵	\$200	
Optional Modular Footplate Extension Components/Kits ^{3,6}		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Side Guard ⁷	\$50	
<input type="checkbox"/> 6" Side Guard ⁸	\$50	
Additional Optional Accessories ⁵ :		
<input type="checkbox"/> Footplate Connector	\$15	
<input type="checkbox"/> Wheel Bumpers, pair	\$25	

1. Only available with Synergy seating.
 2. Please verify the appropriate billing code to be used for this item with all third-party payors.
 3. Extension length is measured from top of the seat pan to footplate.
 4. Choose in Rigging Accessory Section in iQuantum.
 5. Only available on 20"W or greater.
 6. Assembly required. May require seating configuration adjustments to avoid front caster interference.
 7. Only available with 10"W x 6"D Split Footplates
 8. Only available with 10"W x 11"D and 10"W x 8"D split footplates

FRONT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) ¹	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$80	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Standard with Angle Adjustable or Multi-Axis footplates.
 2. Only available with Swing-Away legrests.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

SEAT CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required <i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion <i>General Use</i> <input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
Solution® 1 Cushion <i>Skin Protection:</i> <input type="checkbox"/> 14" x 14" - 20" x 20"	\$380	E2603
Solution® Cushion <i>Skin Protection & Positioning</i> <input type="checkbox"/> 10" x 10" - 20" x 20"	\$380	E2607
Spectrum® Gel Cushion <i>Positioning</i> <input type="checkbox"/> 12" x 10" - 20" x 20"	\$480	E2605
Spectrum® Visco Cushion <i>Skin Protection & Positioning</i> <input type="checkbox"/> 10" x 10" - 20" x 20"	\$480	E2607
Spectrum® Air <i>Adjustable Skin Protection</i> <input type="checkbox"/> 14" x 12" - 20" x 20"	\$480	E2622
Spectrum® Air Contour <i>Adjustable Skin Protection & Positioning</i> <input type="checkbox"/> 14" x 12" - 20" x 20"	\$499	E2624
TRU-Comfort 2™ 4-Way Stretch Seat Cushion <i>Skin Protection</i> <i>Includes 4-way stretch material cover</i> <input type="checkbox"/> 16" x 16" - 20" x 20"	\$380	E2603

ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Rear Basket ¹	\$103	
<input type="checkbox"/> Cane/Crutch Holder ¹	\$98	E2207
<input type="checkbox"/> Weather Cover, Medium, Black	\$158	
<input type="checkbox"/> Cup Holder <i>Side:</i> <input type="checkbox"/> Left <input type="checkbox"/> Right	\$17 ²	
<input type="checkbox"/> Saddle Bag	\$28	
<input type="checkbox"/> Flag Holder	\$17	
<input type="checkbox"/> Oxygen Holder, Contoured Seating ¹	\$103	E2208
<input type="checkbox"/> Oxygen Holder, Synergy Seat ³	\$150	E2208
<input type="checkbox"/> Walker Holder ⁴	\$92	
<input type="checkbox"/> Tablet Holder ⁴	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder ⁵	\$140	

1. Includes mounting hardware. Not available on Manual Recliner.
 2. No charge on Synergy seating.
 3. Not available with Synergy manual recline or manual tilt. Includes mounting hardware.
 4. Only available with Synergy seating.

ADDITIONAL ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Foot Boxes (pair) <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$475	K0108
<input type="checkbox"/> Stealth Padded Shoe Holders (pair) <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large	\$266	K0108
<input type="checkbox"/> Stealth Elbow Support <i>Side:</i> <input type="checkbox"/> Bracket - Left <input type="checkbox"/> Bracket - Right <i>Pad:</i> <input type="checkbox"/> 3" x 3" (ARMS-ES-3x3) <input type="checkbox"/> 4" x 4" (ARMS-ES-4x4)	\$175 ea.	E1028
<input type="checkbox"/> Stealth i-Fit Dynamic Chest Straps <input type="checkbox"/> Upper Classic, Medium (IFPS1103) <input type="checkbox"/> Upper Classic, Large (IFPS1104) <input type="checkbox"/> Upper Contour, Medium (IFPS2103) <input type="checkbox"/> Upper Contour, Large (IFPS2104)	\$142	E0960
<input type="checkbox"/> Stealth i-Fit Chest Strap Belt Anchor <input type="checkbox"/> Medium (IFPS4113) <input type="checkbox"/> Large (IFPS4114)	\$90	

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

LATERALS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)* <i>Bracket:</i> <input type="checkbox"/> Swing-Away Rail Mounts <i>For Use on Synergy Back and Tru-Comfort</i>	\$450	E1028
<i>Pad:</i> <input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover <input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover <input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover <input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200 \$200 \$200 \$200	E0956 E0956 E0956 E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals* <i>Side:</i> <input type="checkbox"/> Left Side <input type="checkbox"/> Right Side <i>Bracket:</i> <input type="checkbox"/> Removable Bracket (ea.) <input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.) <input type="checkbox"/> Flip Down Bracket (ea.) <input type="checkbox"/> Flip Down, Full Surface Contact Bracket (ea.)	\$250 \$325 \$225 \$225	E1028 E1028 E1028 E1028
<i>Pad:</i> <input type="checkbox"/> Short Pad, 4"x4" (ea.) <input type="checkbox"/> Medium Pad, 4"x8" (ea.) <input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$50 \$50 \$50	E0956 E0956 E0956
<i>Adapter:</i> <input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ² <input type="checkbox"/> 2 1/4" Depth/Height Adapter Plate (ea.) ³	\$25 \$25	K0108 K0108

1. Swing-away laterals interfere with flip-up armrests.
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.
 3. Only available on the flip down bracket.

ADDITIONAL HEADREST ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Facial Laterals¹ <input type="checkbox"/> Right Side, 4.5" x 2" Pad (CPS9-R-HW & 900-Q) <input type="checkbox"/> Left Side, 4.5" x 2" Pad (CPS9-L-HW & 900-Q) <input type="checkbox"/> Right and Left 4.5" x 2" Pad (CPS9-B-HW & 900-Q)	\$325 \$325 \$425	E1028 E1028 E1028
<input type="checkbox"/> Egg Switch Mount (ESM) <input type="checkbox"/> Facial Pad, 6" x 2" (Q920-Q) <input type="checkbox"/> Figure Eight Bracket with 13" Rod and Egg Switch Mount (F8678-13-Q) <input type="checkbox"/> Swing-away Joystick Mount, 18" Rod (SUS18-Q)	\$50 ea. \$50 ea. \$250 \$250	 K0108 K0108

1. Only available with Stealth Comfort Series Headrests.

BATTERY OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> U1 32 Amp, AGM <i>Installation:</i> <input type="checkbox"/> Batteries Installed <input type="checkbox"/> Do Not Install	\$133 ea.	E2365