

Quantum® Q6 Edge 2.0 with
TB3 Power Positioning
Group 3



QUANTUM REHAB®

Pride Mobility Products Corporation®
182 Susquehanna Ave., Exeter, PA 18643
Phone: (866)800-2002 • Fax: (866) 707-3422
Email: quantumorders@pridemobility.com

Discount: _____
Terms: _____
Freight: _____
Taken By: _____

Account Number: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
Date: _____
PO#: _____
Marked For: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

Note Section

For Quantum Use:

For Provider Use:

Set-to-Spec

This section must be filled out, otherwise the selection will default to the mid-range measurements.

Provide back height, armrest height, center mount height and legrest height.


Set Back Height: _____

Set Armrest Height Left: _____ Right: _____


Set Center Mount Length: _____

Set Legrest Length Left: _____ Right: _____

Height: _____ Weight: _____



A. Shoulder Width _____
B. Chest Width _____
C. Hip Width _____
D. External Knee Width _____



E. Maximum Sitting Height _____
F. Shoulder Height _____
G. Axilla Height _____
H. PSIS Height _____
I. Trunk Depth _____
J. Ischial Depth _____
K. Thigh Depth _____
L. Lower Leg Length _____
M. Elbow Height _____

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
Base Model		
<input type="checkbox"/> Quantum Q6 Edge 2.0 3SP-SS <i>Single power option. 300 lbs. wt capacity. Comes standard with 6 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24" W x 35.5" L.</i>	\$6,508	K0856
<input type="checkbox"/> Quantum Q6 Edge 2.0 3MP-SS <i>Multiple power option base. 300 lbs. wt capacity. Comes standard with 6 mph programmable motors, 8 amp off board charger. Dimensions of base without legrests are 24" W x 35.5" L.</i>	\$6,308	K0861

Colors

Metallic Colors:

<input type="checkbox"/> Red Corvette	<input type="checkbox"/> Black in Black	<input type="checkbox"/> Pretty in Pink	<input type="checkbox"/> Lime Light	<input type="checkbox"/> Silver Bullet	<input type="checkbox"/> Orange Crush
<input type="checkbox"/> Electric Blue	<input type="checkbox"/> White Satin	<input type="checkbox"/> Deep Purple	<input type="checkbox"/> Grey Street	<input type="checkbox"/> Midnight Blue	<input type="checkbox"/> Yellow Blaze

Base Mounted Foot Platform

<input type="checkbox"/> Base Mounted Foot Platform <i>(Not available with center mount foot platforms)</i>	Standard
<input type="checkbox"/> No Foot Platform	No Charge

Drive Wheels

<input type="checkbox"/> 3" x 14" Flat Free, Black Tires	Standard
<input type="checkbox"/> 4" x 14" Flat Free, Black Tires, Knobby, Extended Traction Tires <i>Increases overall base width by 1.5"</i>	\$551
<input type="checkbox"/> 3" x 14" Pneumatic Black Tires	No Charge
<input type="checkbox"/> 4" x 14" Pneumatic, Black Tires, Knobby, Extended Traction Tires <i>Increases overall base width by 1.5"</i>	\$551

Rim Color:

<input type="checkbox"/> Silver Bullet	<input type="checkbox"/> Matching	<input type="checkbox"/> Non-Matching
--	-----------------------------------	---------------------------------------

Rim Color Options: \$85 *(Only available through SOF)*

<input type="checkbox"/> Red Corvette	<input type="checkbox"/> Black in Black	<input type="checkbox"/> Pretty in Pink	<input type="checkbox"/> Lime Light	<input type="checkbox"/> Midnight Blue	<input type="checkbox"/> Yellow Blaze
<input type="checkbox"/> Electric Blue	<input type="checkbox"/> White Satin	<input type="checkbox"/> Deep Purple	<input type="checkbox"/> Grey Street	<input type="checkbox"/> Orange Crush	

Front Caster Forks

<input type="checkbox"/> Single Sided Forks <i>Coming Soon</i>	Standard
<input type="checkbox"/> Dual Sided Forks	No Charge

Batteries

<input type="checkbox"/> NF-22, 50 Amp Batteries		
<input type="checkbox"/> Maintenance Free AGM (MK) (each)	\$295	E2361
<input type="checkbox"/> Maintenance Free Gel Cell (MK) (each)	\$295	E2361
<input type="checkbox"/> No Batteries		





Installation:

<input type="checkbox"/> Batteries Installed	<input type="checkbox"/> Do Not Install
--	---

LED Light Package











Description	Price	HCPCS Code
<input type="checkbox"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995	

Non-Expandable Controls

Description	Price	HCPCS Code
<input type="checkbox"/>  75 Amp NE	Standard	
<input type="checkbox"/>  75 Amp NE+ 6-Key on Tilt Only <i>For up to 2 actuators through joystick on single power base only</i>	No Charge	
<input type="checkbox"/>  75 Amp NE+ 6-Key <i>For up to 2 Actuators through joystick</i>	\$200	
<input type="checkbox"/>  75 Amp NE+ 8-Key <i>For Lights</i>	\$400	








Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description		Price	HCPCS Code								
Expandable Controls											
<input type="checkbox"/>	 100 Amp Q-Logic 2 EX Controller <i>Required with Q-Logic Electronics</i>	\$695	E2377								
<input type="checkbox"/>	 Expandable Controller Harness <i>Required with Q-Logic Electronics</i>	\$350	E2313								
<input type="checkbox"/>	 Q-Logic 2 EX Joystick	\$700									
<input type="checkbox"/>	 Q-Logic 2 EX Joystick <i>For Lights</i>	\$700									
Q-Logic Color Shroud Options:											
<input type="checkbox"/>	<input checked="" type="radio"/> Black <i>Soft Touch Standard</i>	<input type="checkbox"/> Matching	<input type="checkbox"/> Non-Matching								
Joystick Color Shroud: \$50 <i>(Available through SOF, not available with 4"W tires)</i>											
<input type="checkbox"/>	<input checked="" type="radio"/> Red Corvette	<input type="checkbox"/>	<input checked="" type="radio"/> Black in Black	<input type="checkbox"/>	<input checked="" type="radio"/> Pretty in Pink	<input type="checkbox"/>	<input checked="" type="radio"/> Lime Light	<input type="checkbox"/>	<input checked="" type="radio"/> Silver Bullet	<input type="checkbox"/>	<input checked="" type="radio"/> Orange Crush
<input type="checkbox"/>	<input checked="" type="radio"/> Electric Blue	<input type="checkbox"/>	<input checked="" type="radio"/> White Satin	<input type="checkbox"/>	<input checked="" type="radio"/> Deep Purple	<input type="checkbox"/>	<input checked="" type="radio"/> Grey Street	<input type="checkbox"/>	<input checked="" type="radio"/> Midnight Blue	<input type="checkbox"/>	<input checked="" type="radio"/> Yellow Blaze
<input type="checkbox"/>	 Stand Alone Joystick with External Charger Port	\$700									
<i>Features:</i>											
<ul style="list-style-type: none">Up to 3 drive profiles and 1 seating profile when used alone.Use in conjunction w/ enhanced display for up to 5 profiles.Auxiliary functions including Bluetooth Mouse and IR.											
<input type="checkbox"/>	 Q-Logic 2 Enhanced Display Kit	\$1,884.93									
<i>Includes:</i>											
<ul style="list-style-type: none">Enhanced DisplayFlexible GooseneckGooseneck Mounting Bracket											
<i>Features:</i>											
<ul style="list-style-type: none">Specialty controls interface w/ 3.5" LCD displayUp to 5 profiles to control drive, seat and auxiliary functionsIntegrated single switch scanner, Bluetooth mouse emulator and IR controlsBluetooth connectivity for mobile devices											
Enhanced Display Color Options:											
<input type="checkbox"/>	<input checked="" type="radio"/> Black <i>Soft Touch Standard</i>	<input type="checkbox"/> Matching	<input type="checkbox"/> Non-Matching								
Display Color Shroud: \$50 <i>(Available through SOF, not available with 4"W tires)</i>											
<input type="checkbox"/>	<input checked="" type="radio"/> Red Corvette	<input type="checkbox"/>	<input checked="" type="radio"/> Black in Black	<input type="checkbox"/>	<input checked="" type="radio"/> Pretty in Pink	<input type="checkbox"/>	<input checked="" type="radio"/> Lime Light	<input type="checkbox"/>	<input checked="" type="radio"/> Silver Bullet	<input type="checkbox"/>	<input checked="" type="radio"/> Orange Crush
<input type="checkbox"/>	<input checked="" type="radio"/> Electric Blue	<input type="checkbox"/>	<input checked="" type="radio"/> White Satin	<input type="checkbox"/>	<input checked="" type="radio"/> Deep Purple	<input type="checkbox"/>	<input checked="" type="radio"/> Grey Street	<input type="checkbox"/>	<input checked="" type="radio"/> Midnight Blue	<input type="checkbox"/>	<input checked="" type="radio"/> Yellow Blaze
<input type="checkbox"/>	 Adult i-Drive Tri Array Series, 3-Switch w/ Egg Switch for Mode	\$4,175	E2330								
<i>Stealth: IDH200-1</i>											
<i>Features:</i>											
3 embedded proximity sensors for drive control, and an egg switch for mode change											
<input type="checkbox"/>	 Egg Switches	\$77									
<input type="checkbox"/>	<input checked="" type="radio"/> Red <i>Stealth: ESRED</i>	<input type="checkbox"/>	<input checked="" type="radio"/> Black <i>Stealth: ESBK</i>	<input type="checkbox"/>	<input checked="" type="radio"/> Yellow <i>Stealth: ESYEL</i>						
<input type="checkbox"/>	<input checked="" type="radio"/> Blue <i>Stealth: ESBLU</i>	<input type="checkbox"/>	<input checked="" type="radio"/> Green <i>Stealth: ESGRN</i>								
<input type="checkbox"/>	 USB Mobile Device Charger	\$175									
Available with NE, NE+ and Q-Logic 2 EX controllers; Can only charge 1 device at a time. A multiplier harness is required with all Static and Tilt only seating systems. The selection of more than 3 modules requires an additional multiplier.											
<input type="checkbox"/>	 SCIM Module Kit	\$575									
Allows use of joystick display for specialty controls; Up to 5 profiles to control drive, seat and auxiliary functions; Excludes sip-n-puff and single switch functions; A switch is required to operate the input device attached to the SCIM module; Includes accessory bar on TB3; <i>Only available with Q-Logic 2.</i>											
<input type="checkbox"/>	Click Here for Alternative Drive Devices										
Joystick Handles											
<input type="checkbox"/>	Standard Carrot Handle	Standard									
<input type="checkbox"/>	Mushroom Handle	\$71	E2323								
<input type="checkbox"/>	T-Bar	\$107	E2323								
<input type="checkbox"/>	Chin Cup <i>Use HCPCS Code E2323 when Chin Cup not used in conjunction with a chin control.</i>	\$71	E2324								
<input type="checkbox"/>	Extended I-Handle	\$71	E2323								
<input type="checkbox"/>	Soft Ball Handle	\$50	E2323								
<input type="checkbox"/>	Medium Ball Handle	\$50	E2323								
<input type="checkbox"/>	Bodypoint 3 Inch Goal Post Handle	\$69	E2323								
<input type="checkbox"/>	Bodypoint 4 Inch Goal Post Handle	\$69	E2323								



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
Joystick Mounting Brackets		
<input type="checkbox"/> Fixed, Inline <i>Not available with Q-Logic 2 EX Stand Alone Joystick.</i>	Standard	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-Away	\$210	E1028
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Inline <input type="checkbox"/> Offset		
<input type="checkbox"/> Retractable Swing-Away Joystick Arm Kit <i>Full 180° Movement</i>	\$250	E1028
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Multi-Axis	\$100	E1028
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> No Joystick Mount		
Programming		
<input type="checkbox"/>  Q-Logic EX/NE Programmer	\$545	
<input type="checkbox"/> Q-Logic PC Programming Disk w/ Cable <i>Works with Q-Logic 2 EX only.</i>	\$750	
TRU-Balance®3 Power Positioning		
<input type="checkbox"/> Tru-Balance®3 Power Tilt <i>50 degree center of gravity (CG) power tilt seating system</i>	\$5,395	E1002
<input type="checkbox"/> Tru-Balance®3 Power Recline <i>175 degrees of power recline</i>	\$5,995	E1004
<input type="checkbox"/> Tru-Balance®3 Power Tilt & Recline <i>50 Degree center of gravity (CG) power tilt seating system with 175 degrees of power recline; Full drive lock-out at 20 degrees</i>	\$10,195	E1007
NEW!!! iLevel Option:		
<input type="checkbox"/> Tru-Balance®3 10" Power Adjustable Seat Lift with iLevel <i>iLevel technology enables the unit to drive up to 3 Mph while elevated with 10" of lift. Requires Q-Logic 2 EX. Not available with manual recline back.</i>	\$1,895	E2300
<input type="checkbox"/> Omit iLevel technology		
TRU-Balance®3 Manual Positioning Option		
<input type="checkbox"/> Attendant Manual Recline <i>Available w/ Tilt only. N/A w/ NE Tilt through toggle, AFP or Back Cane Option. Push handles are required.</i>	\$500	E1226
Power Positioning Electronics		
Single Actuator Electronics:		
<input type="checkbox"/> Through Toggle Box <i>Tilt Only. Manual Recline back is not available through NE tilt through toggle.</i>	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through NE+ 6 Key	\$1175	E2310
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$1175	E2310
<input type="checkbox"/>  Through Q-Logic 2 EX Controller and Single Function with 4-Button Actuator <i>Includes AAM and 4-Button actuator module. Only available with TB3 armrests.</i>	\$2,375	E2310
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
Dual Actuator Electronics:		
<input type="checkbox"/> Through NE+ 6-Key Controller	\$1,175/\$2,375	E2310/E2311 ¹
Multiple Actuator Electronics:		
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2,375	E2311
<input type="checkbox"/>  Through Q-Logic 2 EX Controller and 4-Button Actuator <i>Mounts to same side as controller. Includes AAM and 4-Button actuator module. Only available with TB3 armrests.</i>	\$2,375	E2311
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/>  Q-Logic 2 EX Through 4-way Toggle ²	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/>  Q-Logic 2 EX Through 4-way Button Toggle ²	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Omit attendant control <i>(Recline Only)</i>	-\$250	

1. HCPCS codes apply to covered actuators only.

2. Bracket mounts to the outside of the armpad.



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning



Description	Price	HCPSC Code
TRU-Balance®3 Seat Specifications		
Seat Width:		
<input type="checkbox"/> 12"W ^{1,3}	<input type="checkbox"/> 15"W ^{2,3}	<input type="checkbox"/> 18"W
<input type="checkbox"/> 13"W ^{2,3}	<input type="checkbox"/> 16"W	<input type="checkbox"/> 19"W
<input type="checkbox"/> 14"W ³	<input type="checkbox"/> 17"W	<input type="checkbox"/> 20"W
<input type="checkbox"/> 21"W	<input type="checkbox"/> 22"W	
Seat Depth:		
<input type="checkbox"/> 12"D ⁴	<input type="checkbox"/> 15"D	<input type="checkbox"/> 18"D
<input type="checkbox"/> 13"D ⁴	<input type="checkbox"/> 16"D	<input type="checkbox"/> 19"D
<input type="checkbox"/> 14"D	<input type="checkbox"/> 17"D	<input type="checkbox"/> 20"D
<input type="checkbox"/> 21"D	<input type="checkbox"/> 22"D	
Seat Width Package:		
<input type="checkbox"/> 12"W-16"W ³	<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W
Seat Depth Package:		
<input type="checkbox"/> 12"D-16"D ³	<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"D-24"D

1. 16" - 17" depth is not available
2. Only available with non-matching back width

3. Not available with iLevel
4. Only available with back cane option

TRU-Balance®3 Seat to Floor Height		
Tilt or Tilt/Recline:	Recline:	Adjustable Seat Lift:
<input type="checkbox"/> 17.5" Seat to Floor	<input type="checkbox"/> 16.5" Seat to Floor	<input type="checkbox"/> 17.5" Seat to Floor
<input type="checkbox"/> 18.5" Seat to Floor	<input type="checkbox"/> 17.5" Seat to Floor	<input type="checkbox"/> 18.5" Seat to Floor
<input type="checkbox"/> 19.5" Seat to Floor	<input type="checkbox"/> 18.5" Seat to Floor	<input type="checkbox"/> 19.5" Seat to Floor

TRU-Balance®3 TILT Back Type		
<input type="checkbox"/> Static Back		Standard
<input type="checkbox"/> Flip Forward Back	Only available w/ 2-post, flip back or pediatric armrest. Not available on 12"-14"widths or with attendant manual recline.	\$150
<input type="checkbox"/> Back Cane Option	Only available with "Omit Back" and 2-post, flip back or pediatric armrests. Not available with recline.	(-\$125)
Back Cane Type:		
<input type="checkbox"/> Angle Adjustable		Standard
<input type="checkbox"/> Flip Forward		\$35
<input type="checkbox"/> 1" Endo		\$25
<input type="checkbox"/> 2" Endo		\$25
Back Cane Height:		
<input type="checkbox"/> 16"H	<input type="checkbox"/> 17"H	<input type="checkbox"/> 20"H
<input type="checkbox"/> 21"H	<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H



TRU-Balance®3 TILT Back Specifications		
<input type="checkbox"/>		TRU-Comfort Back (Positioning)
		<input type="checkbox"/> 4-way stretch
		<input type="checkbox"/> Air mesh
		<input type="checkbox"/> 20"W or less
		<input type="checkbox"/> 22"W or greater
<input type="checkbox"/>		Sport Back
		Size:
<input type="checkbox"/> 12"Wx14"H	<input type="checkbox"/> 17"Wx16"H	<input type="checkbox"/> 18"Wx20"H
<input type="checkbox"/> 12"Wx15"H	<input type="checkbox"/> 17"Wx17"H	<input type="checkbox"/> 18"Wx21"H
<input type="checkbox"/> 14"Wx16"H	<input type="checkbox"/> 17"Wx18"H	<input type="checkbox"/> 18"Wx22"H
<input type="checkbox"/> 14"Wx17"H	<input type="checkbox"/> 17"Wx19"H	<input type="checkbox"/> 18"Wx23"H
<input type="checkbox"/> 16"Wx16"H	<input type="checkbox"/> 17"Wx20"H	<input type="checkbox"/> 18"Wx24"H
<input type="checkbox"/> 16"Wx17"H	<input type="checkbox"/> 17"Wx21"H	<input type="checkbox"/> 18"Wx25"H
<input type="checkbox"/> 16"Wx18"H	<input type="checkbox"/> 17"Wx22"H	<input type="checkbox"/> 19"Wx16"H
<input type="checkbox"/> 16"Wx19"H	<input type="checkbox"/> 17"Wx23"H	<input type="checkbox"/> 19"Wx17"H
<input type="checkbox"/> 16"Wx20"H	<input type="checkbox"/> 18"Wx16"H	<input type="checkbox"/> 19"Wx18"H
<input type="checkbox"/> 16"Wx21"H	<input type="checkbox"/> 18"Wx17"H	<input type="checkbox"/> 19"Wx19"H
<input type="checkbox"/> 16"Wx22"H	<input type="checkbox"/> 18"Wx18"H	<input type="checkbox"/> 19"Wx20"H
<input type="checkbox"/> 16"Wx23"H	<input type="checkbox"/> 18"Wx19"H	<input type="checkbox"/> 19"Wx21"H
<input type="checkbox"/> 19"Wx22"H	<input type="checkbox"/> 19"Wx22"H	<input type="checkbox"/> 20"Wx16"H
<input type="checkbox"/> 19"Wx23"H	<input type="checkbox"/> 19"Wx23"H	<input type="checkbox"/> 20"Wx17"H
<input type="checkbox"/> 20"Wx16"H	<input type="checkbox"/> 20"Wx18"H	<input type="checkbox"/> 20"Wx19"H
<input type="checkbox"/> 20"Wx17"H	<input type="checkbox"/> 20"Wx19"H	<input type="checkbox"/> 20"Wx20"H
<input type="checkbox"/> 20"Wx18"H	<input type="checkbox"/> 20"Wx20"H	<input type="checkbox"/> 20"Wx21"H
<input type="checkbox"/> 20"Wx19"H	<input type="checkbox"/> 20"Wx21"H	<input type="checkbox"/> 20"Wx22"H
<input type="checkbox"/> 20"Wx20"H	<input type="checkbox"/> 20"Wx22"H	<input type="checkbox"/> 20"Wx23"H
<input type="checkbox"/> 20"Wx21"H	<input type="checkbox"/> 20"Wx23"H	<input type="checkbox"/> 20"Wx24"H
<input type="checkbox"/> 20"Wx22"H	<input type="checkbox"/> 20"Wx24"H	<input type="checkbox"/> 20"Wx25"H
<input type="checkbox"/> 20"Wx23"H	<input type="checkbox"/> 21"Wx16"H ¹	<input type="checkbox"/> 21"Wx17"H ¹
<input type="checkbox"/> 20"Wx24"H	<input type="checkbox"/> 21"Wx17"H ¹	<input type="checkbox"/> 21"Wx18"H ¹
<input type="checkbox"/> 20"Wx25"H	<input type="checkbox"/> 21"Wx18"H ¹	<input type="checkbox"/> 21"Wx19"H ¹
	<input type="checkbox"/> 21"Wx19"H ¹	<input type="checkbox"/> 21"Wx20"H ¹
	<input type="checkbox"/> 21"Wx20"H ¹	<input type="checkbox"/> 21"Wx21"H ¹
	<input type="checkbox"/> 21"Wx21"H ¹	<input type="checkbox"/> 21"Wx22"H ¹
	<input type="checkbox"/> 21"Wx22"H ¹	<input type="checkbox"/> 21"Wx23"H ¹
	<input type="checkbox"/> 21"Wx23"H ¹	<input type="checkbox"/> 21"Wx24"H ¹
	<input type="checkbox"/> 21"Wx24"H ¹	<input type="checkbox"/> 21"Wx25"H ¹
	<input type="checkbox"/> 22"Wx16"H	<input type="checkbox"/> 22"Wx17"H
	<input type="checkbox"/> 22"Wx17"H	<input type="checkbox"/> 22"Wx18"H
	<input type="checkbox"/> 22"Wx18"H	<input type="checkbox"/> 22"Wx19"H
	<input type="checkbox"/> 22"Wx19"H	<input type="checkbox"/> 22"Wx20"H
	<input type="checkbox"/> 22"Wx20"H	<input type="checkbox"/> 22"Wx21"H
	<input type="checkbox"/> 22"Wx21"H	<input type="checkbox"/> 22"Wx22"H
	<input type="checkbox"/> 22"Wx22"H	<input type="checkbox"/> 22"Wx23"H
	<input type="checkbox"/> 22"Wx23"H	<input type="checkbox"/> 22"Wx24"H
	<input type="checkbox"/> 22"Wx24"H	<input type="checkbox"/> 22"Wx25"H
	<input type="checkbox"/> 23"Wx16"H ¹	<input type="checkbox"/> 23"Wx17"H ¹
	<input type="checkbox"/> 23"Wx17"H ¹	<input type="checkbox"/> 23"Wx18"H ¹
	<input type="checkbox"/> 23"Wx18"H ¹	<input type="checkbox"/> 23"Wx19"H ¹
	<input type="checkbox"/> 23"Wx19"H ¹	<input type="checkbox"/> 23"Wx20"H ¹
	<input type="checkbox"/> 23"Wx20"H ¹	<input type="checkbox"/> 23"Wx21"H ¹
	<input type="checkbox"/> 23"Wx21"H ¹	<input type="checkbox"/> 23"Wx22"H ¹
	<input type="checkbox"/> 23"Wx22"H ¹	<input type="checkbox"/> 23"Wx23"H ¹
	<input type="checkbox"/> 23"Wx23"H ¹	<input type="checkbox"/> 23"Wx24"H ¹
	<input type="checkbox"/> 23"Wx24"H ¹	<input type="checkbox"/> 23"Wx25"H ¹
	<input type="checkbox"/> 24"Wx16"H	<input type="checkbox"/> 24"Wx17"H
	<input type="checkbox"/> 24"Wx17"H	<input type="checkbox"/> 24"Wx18"H
	<input type="checkbox"/> 24"Wx18"H	<input type="checkbox"/> 24"Wx19"H
	<input type="checkbox"/> 24"Wx19"H	<input type="checkbox"/> 24"Wx20"H
	<input type="checkbox"/> 24"Wx20"H	<input type="checkbox"/> 24"Wx21"H
	<input type="checkbox"/> 24"Wx21"H	<input type="checkbox"/> 24"Wx22"H
	<input type="checkbox"/> 24"Wx22"H	<input type="checkbox"/> 24"Wx23"H
	<input type="checkbox"/> 24"Wx23"H	

Back Width Package:			
<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W	<input type="checkbox"/> 25"W-28"W

1. Not available with TRU-Comfort Back



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPSC Code
TRU-Balance®3 RECLINE OR TILT/RECLINE Back Specifications		
+/-4" of Seat Width. Odd heights = 6" Gap, even heights = 7" gap. Measurements listed is overall back height.		
<input type="checkbox"/>  TRU-Comfort Back (Positioning) <input type="checkbox"/> 4-way stretch <input type="checkbox"/> Air mesh <input type="checkbox"/> 20"W or less <input type="checkbox"/> 22"W or greater	\$595	E2620
<input type="checkbox"/>  Sport Back	\$695	E2621
<input type="checkbox"/> Size:	No Charge	
<input type="checkbox"/> 12"Wx19"H ¹	<input type="checkbox"/> 17"Wx19"H	<input type="checkbox"/> 18"Wx23"H
<input type="checkbox"/> 12"Wx20"H ¹	<input type="checkbox"/> 17"Wx20"H	<input type="checkbox"/> 18"Wx24"H
<input type="checkbox"/> 14"Wx19"H	<input type="checkbox"/> 17"Wx21"H	<input type="checkbox"/> 18"Wx25"H
<input type="checkbox"/> 14"Wx20"H	<input type="checkbox"/> 17"Wx22"H	<input type="checkbox"/> 18"Wx26"H
<input type="checkbox"/> 16"Wx19"H	<input type="checkbox"/> 17"Wx23"H	<input type="checkbox"/> 18"Wx27"H
<input type="checkbox"/> 16"Wx20"H	<input type="checkbox"/> 17"Wx24"H	<input type="checkbox"/> 18"Wx28"H
<input type="checkbox"/> 16"Wx21"H	<input type="checkbox"/> 17"Wx25"H	<input type="checkbox"/> 19"Wx19"H
<input type="checkbox"/> 16"Wx22"H	<input type="checkbox"/> 17"Wx26"H	<input type="checkbox"/> 19"Wx20"H
<input type="checkbox"/> 16"Wx23"H	<input type="checkbox"/> 18"Wx19"H	<input type="checkbox"/> 19"Wx21"H
<input type="checkbox"/> 16"Wx24"H	<input type="checkbox"/> 18"Wx20"H	<input type="checkbox"/> 19"Wx22"H
<input type="checkbox"/> 16"Wx25"H	<input type="checkbox"/> 18"Wx21"H	<input type="checkbox"/> 19"Wx23"H
<input type="checkbox"/> 16"Wx26"H	<input type="checkbox"/> 18"Wx22"H	<input type="checkbox"/> 19"Wx24"H
<input type="checkbox"/> 19"Wx25"H	<input type="checkbox"/> 19"Wx26"H	<input type="checkbox"/> 20"Wx19"H
<input type="checkbox"/> 19"Wx26"H	<input type="checkbox"/> 20"Wx19"H	<input type="checkbox"/> 20"Wx20"H
<input type="checkbox"/> 20"Wx19"H	<input type="checkbox"/> 20"Wx20"H	<input type="checkbox"/> 20"Wx21"H
<input type="checkbox"/> 20"Wx20"H	<input type="checkbox"/> 20"Wx21"H	<input type="checkbox"/> 20"Wx22"H
<input type="checkbox"/> 20"Wx21"H	<input type="checkbox"/> 20"Wx22"H	<input type="checkbox"/> 20"Wx23"H
<input type="checkbox"/> 20"Wx22"H	<input type="checkbox"/> 20"Wx23"H	<input type="checkbox"/> 20"Wx24"H
<input type="checkbox"/> 20"Wx23"H	<input type="checkbox"/> 20"Wx24"H	<input type="checkbox"/> 20"Wx25"H
<input type="checkbox"/> 20"Wx24"H	<input type="checkbox"/> 20"Wx25"H	<input type="checkbox"/> 20"Wx26"H
<input type="checkbox"/> 20"Wx25"H	<input type="checkbox"/> 20"Wx26"H	<input type="checkbox"/> 20"Wx27"H
<input type="checkbox"/> 20"Wx26"H	<input type="checkbox"/> 20"Wx27"H	<input type="checkbox"/> 20"Wx28"H
<input type="checkbox"/> 21"Wx19"H ¹	<input type="checkbox"/> 21"Wx20"H ¹	<input type="checkbox"/> 21"Wx21"H ¹
<input type="checkbox"/> 21"Wx20"H ¹	<input type="checkbox"/> 21"Wx21"H ¹	<input type="checkbox"/> 21"Wx22"H ¹
<input type="checkbox"/> 21"Wx21"H ¹	<input type="checkbox"/> 21"Wx22"H ¹	<input type="checkbox"/> 21"Wx23"H ¹
<input type="checkbox"/> 21"Wx22"H ¹	<input type="checkbox"/> 21"Wx23"H ¹	<input type="checkbox"/> 21"Wx24"H ¹
<input type="checkbox"/> 21"Wx23"H ¹	<input type="checkbox"/> 21"Wx24"H ¹	<input type="checkbox"/> 21"Wx25"H ¹
<input type="checkbox"/> 21"Wx24"H ¹	<input type="checkbox"/> 21"Wx25"H ¹	<input type="checkbox"/> 21"Wx26"H ¹
<input type="checkbox"/> 21"Wx25"H ¹	<input type="checkbox"/> 21"Wx26"H ¹	<input type="checkbox"/> 21"Wx27"H ¹
<input type="checkbox"/> 21"Wx26"H ¹	<input type="checkbox"/> 21"Wx27"H ¹	<input type="checkbox"/> 21"Wx28"H ¹
<input type="checkbox"/> 22"Wx21"H	<input type="checkbox"/> 22"Wx22"H	<input type="checkbox"/> 22"Wx23"H
<input type="checkbox"/> 22"Wx22"H	<input type="checkbox"/> 22"Wx23"H	<input type="checkbox"/> 22"Wx24"H
<input type="checkbox"/> 22"Wx23"H	<input type="checkbox"/> 22"Wx24"H	<input type="checkbox"/> 22"Wx25"H
<input type="checkbox"/> 22"Wx24"H	<input type="checkbox"/> 22"Wx25"H	<input type="checkbox"/> 22"Wx26"H
<input type="checkbox"/> 22"Wx25"H	<input type="checkbox"/> 22"Wx26"H	<input type="checkbox"/> 22"Wx27"H
<input type="checkbox"/> 22"Wx26"H	<input type="checkbox"/> 22"Wx27"H	<input type="checkbox"/> 22"Wx28"H
<input type="checkbox"/> 23"Wx21"H ¹	<input type="checkbox"/> 23"Wx22"H ¹	<input type="checkbox"/> 23"Wx23"H ¹
<input type="checkbox"/> 23"Wx22"H ¹	<input type="checkbox"/> 23"Wx23"H ¹	<input type="checkbox"/> 23"Wx24"H ¹
<input type="checkbox"/> 23"Wx23"H ¹	<input type="checkbox"/> 23"Wx24"H ¹	<input type="checkbox"/> 23"Wx25"H ¹
<input type="checkbox"/> 23"Wx24"H ¹	<input type="checkbox"/> 23"Wx25"H ¹	<input type="checkbox"/> 23"Wx26"H ¹
<input type="checkbox"/> 23"Wx25"H ¹	<input type="checkbox"/> 23"Wx26"H ¹	<input type="checkbox"/> 23"Wx27"H ¹
<input type="checkbox"/> 23"Wx26"H ¹	<input type="checkbox"/> 23"Wx27"H ¹	<input type="checkbox"/> 23"Wx28"H ¹
<input type="checkbox"/> 24"Wx19"H	<input type="checkbox"/> 24"Wx20"H	<input type="checkbox"/> 24"Wx21"H
<input type="checkbox"/> 24"Wx20"H	<input type="checkbox"/> 24"Wx21"H	<input type="checkbox"/> 24"Wx22"H
<input type="checkbox"/> 24"Wx21"H	<input type="checkbox"/> 24"Wx22"H	<input type="checkbox"/> 24"Wx23"H
<input type="checkbox"/> 24"Wx22"H	<input type="checkbox"/> 24"Wx23"H	<input type="checkbox"/> 24"Wx24"H
<input type="checkbox"/> 24"Wx23"H	<input type="checkbox"/> 24"Wx24"H	<input type="checkbox"/> 24"Wx25"H
<input type="checkbox"/> 24"Wx24"H	<input type="checkbox"/> 24"Wx25"H	<input type="checkbox"/> 24"Wx26"H
<input type="checkbox"/> 24"Wx25"H	<input type="checkbox"/> 24"Wx26"H	
<input type="checkbox"/> 24"Wx26"H		
Back Width Package:		
<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W
<input type="checkbox"/> 25"W-28"W		

1. Not available with TRU-Comfort Back

TRU-Balance®3 Back Specifications

Back Color Shroud Option:

☐ ● Black in Black Standard

Back Color Options: \$100

<input type="checkbox"/> Red Corvette	<input type="checkbox"/> White Satin	<input type="checkbox"/> Deep Purple	<input type="checkbox"/> Grey Street	<input type="checkbox"/> Midnight Blue	<input type="checkbox"/> Yellow Blaze
<input type="checkbox"/> Electric Blue	<input type="checkbox"/> Pretty in Pink	<input type="checkbox"/> Lime Light	<input type="checkbox"/> Silver Bullet	<input type="checkbox"/> Orange Crush	

Additional TRU-Balance®3 Back Options



- | | |
|---|----------|
| <input type="checkbox"/> Large Lateral Positioning Wedge, pair | \$30 |
| Not available on Tilt Backs 16"-17"H w/ limited adjustability on 18"H. Not available on Recline Backs 19"-20"H w/ limited adjustability on 21"H. N/A with Back Cane Option | |
| <input type="checkbox"/> Aftermarket Back Interface | \$200 |
| Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Contact your Quantum specialist for details. Does not include brackets to mount aftermarket back. SOF is required on back widths less than 16". | |
| <input type="checkbox"/> Omit Back | (-\$100) |
| <input type="checkbox"/> Comfort Company Back <i>See Comfort Company order form</i> | Varies |
| Privacy Flap: Standard on recline seating systems. | |
| <input type="checkbox"/> 8.75"W | \$35 |
| <input type="checkbox"/> 12.75"W | \$35 |

Lap Belts

- | | |
|--|----------|
| <input type="checkbox"/> 38" Lap Belt | Standard |
| <input type="checkbox"/> 50" Lap Belt | Standard |
| <input type="checkbox"/> 60" Lap Belt | Standard |
| <input type="checkbox"/> 70" Lap Belt | Standard |
| <input type="checkbox"/> 80" Lap Belt | Standard |
| <input type="checkbox"/> 90" Lap Belt | Standard |
| <input type="checkbox"/> 60" Padded Push Button, center Pull | \$115 |
| <input type="checkbox"/> Large 13-19", Hip Padded, Push Button | \$92 |



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning




Description	Price	HCPSC Code
Laterals		
<input type="checkbox"/>  Stealth Thoracic Lateral <i>Swing-Away laterals may interfere with TB3 and 2-post, flip-up armrest. Not available with Back Cane Option</i>		
Bracket: <input type="checkbox"/> Swing-Away, Rail Mounts (pair)	\$450	E1028
Pads: <input type="checkbox"/> Curved Pad, 2" x 3" (pair)	\$200	E0956
<input type="checkbox"/> Curved Pad, 3" x 5" (pair)	\$200	E0956
<input type="checkbox"/> Curved Pad, 4" x 6" (pair)	\$200	E0956
Lateral Installation: <input type="checkbox"/> Install <input type="checkbox"/> Do not Install		
<input type="checkbox"/>  Stealth Pelvic/Thigh Guides <i>Pelvic/Thigh guides may interfere with the joystick receiver block and down posts.</i>		
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
Bracket: <input type="checkbox"/> Flip Down (each)	\$225	E1028
<input type="checkbox"/> Flip Down with Full Surface Contact Feature (each)	\$225	E1028
<input type="checkbox"/> Removable (each)	\$250	E1028
<input type="checkbox"/> Removable w/ 1" Offset (each)	\$325	E1028
Adapter: <input type="checkbox"/> 2"Height Adapter Plate (each)	\$25	K0108
<i>Height adapter is needed w/ TRU-Comfort seat cushions. Not available with the flip down bracket.</i>		
<input type="checkbox"/> 2 1/4" Adapter Plate (each) <i>Only available with the flip down full surface bracket</i>	\$50	K0108
Pads: <input type="checkbox"/> Short Pad, 4" x 4" (each)	\$50	E0956
<input type="checkbox"/> Medium Pad, 4" x 8" (each)	\$50	E0956
<input type="checkbox"/> Long Pad, 4" x 12" (each)	\$50	E0956
Lateral Installation: <input type="checkbox"/> Install <input type="checkbox"/> Do not Install		






Armrests		
<input type="checkbox"/> NEW!! TB3 Armrests Version 2 <i>May interfere with swing-away laterals. Choose height package below.</i>	No Charge	
<input type="checkbox"/> Down Post +/- up to 2" difference between seat and back width. Not Available on 14"W Seat/Back	\$175	
Left Armrest Height Package:	Right Armrest Height Package:	
<input type="checkbox"/> 8"-12" Height <input type="checkbox"/> 12"-16" Height	<input type="checkbox"/> 8"-12" Height <input type="checkbox"/> 12"-16" Height	
<input type="checkbox"/> 2-Post, Flip Back, Height Adjustable, Removable <i>Height range 9"-13"</i>	No Charge	
<i>Not available with a depth less than 16"; Required with flip forward back and laptrays.</i>		
<input type="checkbox"/> Pediatric, Removable, Height Adjustable <i>Height adjustment from 7.75" to 9.75" in 1/2" increments. Not available with full length armpads.</i>	No Charge	
<input type="checkbox"/> Quad Release Lever		
Left Armpad Length:	Right Armpad Length:	
<input type="checkbox"/> 14"L Full Length <input type="checkbox"/> 10"L Desk Length	<input type="checkbox"/> 14"L Full Length <input type="checkbox"/> 10"L Desk Length	
Left Armpad Style:	Right Armpad Style:	
<input type="checkbox"/> Straight <input type="checkbox"/> Waterfall	<input type="checkbox"/> Straight <input type="checkbox"/> Waterfall	
<input type="checkbox"/> Falcon, Gel, Full Length 3"W x 14"L \$92.50 ea.	<input type="checkbox"/> Falcon, Gel, Full Length 3"W x 14"L \$92.50 ea.	
<input type="checkbox"/> Mounting Plate, Full Length ¹ \$56.52 ea.	<input type="checkbox"/> Mounting Plate, Full Length ¹ \$56.52 ea.	
<input type="checkbox"/> Falcon, Gel, Desk Length 3"W x 10"L 70 ea.	<input type="checkbox"/> Falcon, Gel, Desk Length 3"W x 10"L \$70 ea.	
<input type="checkbox"/> Mounting Plate, Desk Length ¹ \$56.20 ea.	<input type="checkbox"/> Mounting Plate, Desk Length ¹ \$56.20 ea.	





1. Required only with 2-post, flip back, height adjustable, removable armrests or pediatric armrests. Interference with Back may occur.



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPSC Code
Additional Armrest Accessories		
<i>See positioning components order form for more options</i>		
<input type="checkbox"/>  Ottobock Arm Channels (each) <i>A 2" narrower back width is recommended to avoid interference w/ the back cushion.</i> <input type="checkbox"/> Medium (4.75" x 11.5") <i>(Desk length armrest required)</i> <input type="checkbox"/> Large (4.75" x 14") <i>(Full length armrest required)</i> <input type="checkbox"/> One Piece (4.75" x 18") <i>(Full length armrest required)</i> Side: <input type="checkbox"/> Left <input type="checkbox"/> Right	\$104	E2209
<input type="checkbox"/>  Ottobock Hand Pads (each) <i>Not available on same side as controller</i> Side: <input type="checkbox"/> Left <input type="checkbox"/> Right Size: <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$52	
<input type="checkbox"/>  Stealth Elbow Support (each) Side: <input type="checkbox"/> Bracket-Left (each) <input type="checkbox"/> Bracket-Right (each) Pad Size: <input type="checkbox"/> 3"W x 3"H (each) <input type="checkbox"/> 4"W x 4"H (each)	\$170 \$170 \$50 \$50	E1028 E1028 K0108 K0108









TRU-Balance®3 Complimentary Accessories		
<i>Choose two complimentary items in this section. Each additional item will be \$25 MSRP.</i>		
<input type="checkbox"/>  Cup Holder <i>Not available with 2-post, flip back or pediatric armrests or Back Cane Option</i>		
<input type="checkbox"/>  Personal Item Hook <i>Not available with Back Cane Option</i>		
<input type="checkbox"/>  Reflector Kit		
<input type="checkbox"/>  Cell Phone/Media Holder <i>Not available w/ 2-post, flip back or pediatric armrests or Back Cane Options</i>		
<input type="checkbox"/>  Backpack Holder, pair <i>Not available with Back Cane Option.</i>		

TRU-Balance®3 Premium Installed Accessories		
<input type="checkbox"/>  Trilogy Vent Tray <i>Not available with Back Cane option</i>	\$400	E1029
<input type="checkbox"/>  Compact Vent Tray	\$1,395	E1029
<input type="checkbox"/>  Push Handles (pair) <i>Not available with O2 holder, Back Cane Option, or vent tray. Required with attendant manual recline.</i>	\$175	
<input type="checkbox"/>  Accessory bar <i>Comes standard with the SCIM module, Trilogy vent tray and O2 holder. May be added to mount additional accessories to back of unit. Not available with Back Cane Option.</i>	\$100	

Additional Accessories		
<input type="checkbox"/> Quantum Backpack <i>May interfere with other rear accessories</i>	\$60	
<input type="checkbox"/> Weather Cover, Large, Black	\$158	







Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
TRU-Balanace®3 Premium Accessories		
<i>Premium accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur.</i>		
<input type="checkbox"/>  Clothing Guard, Pair	\$175	
<input type="checkbox"/>  Glove Box	\$100	
<input type="checkbox"/>  Hydration System	\$175	
<input type="checkbox"/>  Oxygen Holder	\$150	E2208
<input type="checkbox"/>  Oxygen Holder <i>for use w/ Back Cane Option</i>	\$150	E2208
<input type="checkbox"/>  Transfer Bar, Pair <i>(Quick Release)</i>	\$175	
<input type="checkbox"/>  Transfer Bar, Pair <i>(Bolt-On)</i>	\$175	
<input type="checkbox"/>  iPad Holder	\$375	
Synergy® Cushions		
<input type="checkbox"/> Matching Seat Size		
<input type="checkbox"/> Non-Matching Seat Size		
<i>Please refer to Synergy Cushion order forms for detailed pricing, sizes, codes and part #'s</i>		
Simplicity Cushion: General Use		
<input type="checkbox"/> 16" x 16" - 20" x 24"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
Solution® Cushion: Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608
Solution® 1 Cushion: Skin Protection		
<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604
Spectrum® Gel Cushion: Positioning		
<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606
Spectrum® Foam Cushion: Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608
Spectrum® Air Flat Cushion: Adjustable Skin Protection		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2623
Air Tube:		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
Spectrum® Air Contour Cushion: Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
Air Tube:		
<input type="checkbox"/> Left Side		
<input type="checkbox"/> Right Side		
TRU-Comfort™2 4-Way Stretch: Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™2 Air Mesh: Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™ Accessory:		
<input type="checkbox"/> Sport Seat Cushion Cover	\$100	
<i>Available as an additional cover option for both TRU-Comfort seat cushion styles. Only available in sizes listed.</i>		
<input type="checkbox"/> 16"W x 16"D		
<input type="checkbox"/> 18"W x 18"D		
<input type="checkbox"/> 20"W x 20"D		



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
Headrests		
Headrest Pads		
<input type="checkbox"/> Stealth 10" Comfort Plus	\$142	E0955
<input type="checkbox"/> Stealth 14" Comfort Plus	\$142	E0955
<input type="checkbox"/> Stealth Small Contoured (9.75"W x 3"H)	\$205	E0955
<input type="checkbox"/> Stealth Large Contoured (11"W x 6"H)	\$205	E0955
Headrest Mounts		
<input type="checkbox"/> Stealth Removable Headrest Bracket	\$210	E1028
<input type="checkbox"/> Stealth's Worlds Best Removable Headrest Bracket	\$210	E1028
<input type="checkbox"/> No Bracket		

Swing-Away Headrest Mounts		
Description	Price	HCPCS Code
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Right Side, for Comfort Series ¹ (Stealth: CPS9-R)	\$325.00	E1028
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Left Side, for Comfort Series ¹ (Stealth: CPS9-L)	\$325.00	E1028
<input type="checkbox"/>  Facial Lateral, 4.5" x 2", Left and Right Side, for Comfort Series ¹ (Stealth: CPS9-B)	\$425.00	E1028
<input type="checkbox"/>  Egg Switch Mount (Stealth: ESM) (each)	\$50.00	

1. Only available with the Stealth Products Comfort series headrests.

Center Mount Foot Platforms		
<i>All efforts will be made to accommodate legrest length selection, but minor changes may be needed to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of unit. Additional foot riggings can be found on page 11.</i>		
<input type="checkbox"/> Adjustable Center Mount Foot Platform Includes 10"W x 11"D Split Footplates	\$250	
<input type="checkbox"/> Center Mount Power Articulating Foot Platform Includes 10"W x 11"D Split Footplates	\$2,895	K0108
Extensions:		
<input type="checkbox"/> 10-16"L		
<input type="checkbox"/> 13-19"L		
Footplates:		
<input type="checkbox"/> 10"W x 11"D Split Footplates	Standard	
<input type="checkbox"/> 10"W x 8"D Split Footplates	No Charge	
Calf Supports:		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small ¹ (6"H x 6"W)	\$150	
<input type="checkbox"/> Large ¹ (8"H x 8"W, tapered to 6"W)	\$150	
<input type="checkbox"/> Calf Pad Covers (pair) (1/2" Soft Visco Foam)		
<input type="checkbox"/> Air Mesh	\$120	
<input type="checkbox"/> Reversed Dartex	\$120	
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W)	\$200	
<i>Only available on 20"W or greater.</i>		
<input type="checkbox"/> Calf Panel Covers (1/2" Soft Visco Foam)		
<input type="checkbox"/> Air Mesh	\$120	
<input type="checkbox"/> Reversed Dartex	\$120	

1. Comes standard with articulating foot platform.





Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
Center Mount Foot Platforms (cont.)		

All efforts will be made to accommodate legrest length selection, but minor changes may be needed to be made to meet shipping guidelines. Some adjustment may be needed upon receipt of unit.

Optional Modular Footplate Extension Components/Kits:

Assembly required. May require seating configuration adjustments to avoid front caster interference. Not available with 10"W x 8"D Split Footplates.

<input type="checkbox"/>		4" Width Extension Kit	\$50
<input type="checkbox"/>		2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/>		4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75
<input type="checkbox"/>		6-Piece Extension Stock Kit <i>(Includes all above pieces)</i>	\$95

Additional Optional Accessories:

<input type="checkbox"/>	Wheel Bumpers, pair	\$25
<input type="checkbox"/>	Footplate Connector	\$15
<input type="checkbox"/>	Footplate Cover, Reversed Dartex <i>(1/2" Soft Visco Foam)</i>	\$120








Not available with modular footplate extension components/kits.

Adjustable Center Mount Foot Platform Compatibility Matrix	TRU-Balance®3 Recline			TRU-Balance®3 Tilt			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
10"-16" Extension	!	!	✓	✓	✓	✓	✓	✓	✓
13"-19" Extension	!	!	!	!	!	!	!	!	!

- ☒ This rigging Configuration causes no interference issues at this seat-to-floor height.
- ☐ This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.

Additional Foot Rigging Accessories

See Positioning Components order form for more options

<input type="checkbox"/>		Therafin Flip-down Knee Abductor bracket (each)	\$172	E1028
		Pads:		
		<input type="checkbox"/> Small (3.5" x 3.25" x 4")	\$135	E0957
		<input type="checkbox"/> Small, Gel (3.5" x 3.25" x 4")	\$145	E0957
		<input type="checkbox"/> Medium (4" x 3.25" x 6")	\$135	E0957
		<input type="checkbox"/> Large (4.5" x 2.75" x 6")	\$135	E0957
<input type="checkbox"/>		POSAlinc Knee Adductor for Swing-Away Legrests (pair)	\$512	K0108
<input type="checkbox"/>		POSAlinc Rail Mounted Knee Adductors (pair)	\$574	K0108
<input type="checkbox"/>		Stealth Foot Boxes (pair)	\$475	K0108
		<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/>		Stealth Padded Shoe Holders (pair)	\$166	K0108
		<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/>		Ankle Huggers (pair)	\$156	K0108
		<input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/>		Bodypoint High Mount, Clamp-On Footplates (pair)		
		<input type="checkbox"/> Angle Adjustable: <i>(Swing-away legrest upper assembly and receiver is required)</i>		
		<input type="checkbox"/> Small, 4"x 7.75"	\$356	K0040
		<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	K0040
		<input type="checkbox"/> Multi-Axis:		
		<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	K0040



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPSC Code
Powered Elevating Legrests		
<i>All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.</i>		
<input type="checkbox"/> 70° Power Articulating Swing-Away Legrests 4" articulation range <i>Use code K0108 for an individual legrest.</i>	\$2,395	E1010
Rigging Selection:		
<input type="checkbox"/> Standard Matching Legrests (pair)	<input type="checkbox"/> Non-Matching Legrests (pair)	<input type="checkbox"/> Individual Legrests Left or Right Only
Operation Method:		
<input type="checkbox"/> Independent	<input type="checkbox"/> Combined	
Footplate Selection:		
<input type="checkbox"/> Angle Adjustable Footplates (pair) Includes heel loops.	\$150	K0040
Left Legrest Lower Extensions:		
<input type="checkbox"/> 11"-13"L	<input type="checkbox"/> 14.5"-18"L	<input type="checkbox"/> 17.5"-21"L
Right Legrest Lower Extensions:		
<input type="checkbox"/> 11"-13"L	<input type="checkbox"/> 14.5"-18"L	<input type="checkbox"/> 17.5"-21"L
<input type="checkbox"/> Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4"W x 6"D	<input type="checkbox"/> 5"W x 8"D	<input type="checkbox"/> 6"W x 11"D ¹
<input type="checkbox"/> 5"W x 6"D	<input type="checkbox"/> 6"W x 8"D ¹	<input type="checkbox"/> 8"W x 11"D ²
Accessories:		
<input type="checkbox"/> Small Calf Pad (pair) 5.35"W - 6"H	No Charge	
<input type="checkbox"/> Large Calf Pad (pair) 5.86"W x 8"H	No Charge	
<input type="checkbox"/> Gel Padded Knee Adductor Buttons (pair)	\$125	K0108

1. May cause interference on seat widths less than 18"

2. May cause interference on seat widths less than 20"

Power Elevating Legrest Compatibility Matrix	TRU-Balance®3 Recline			TRU-Balance®3 Tilt			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
Angle-Adjustable, 11" - 13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle-Adjustable, 14.5" - 18"	!	!	!	!	!	✓	!	!	✓
Angle-Adjustable, 17.5" - 21"	X	X	X	X	X	!	X	!	!

☒ This rigging Configuration causes no interference issues at this seat-to-floor height.☐ This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.☒ This rigging configuration does not work with this seat-to-floor height.

Non-Powered Legrests

<i>All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.</i>		
<input type="checkbox"/> 60-90° Swing-Away Legrests adjustable in 5° Increments	\$162	
Rigging Selection:		
<input type="checkbox"/> Standard Matching Legrests (pair)	<input type="checkbox"/> Non-Matching Legrests (pair)	<input type="checkbox"/> Individual Legrests Left or Right Only
Footplate Selection:		
<input type="checkbox"/> Angle-Adjustable Footplates (pair) Includes Heel Loops.	\$150	K0040
<input type="checkbox"/> Multi-Axis Footplates (pair) Includes Heel Loops.	\$495	K0040
<input type="checkbox"/> Standard, Extruded Footplates (pair)	No Charge	
Left Legrest:		
Lower Extensions:		
<input type="checkbox"/> 11"-13"L ¹	<input type="checkbox"/> 13"-16.5"L ¹	
<input type="checkbox"/> 11"-16.5"L ²	<input type="checkbox"/> 15.5"-19.5"L ²	
Right Legrest:		
Lower Extensions:		
<input type="checkbox"/> 11"-13"L ¹	<input type="checkbox"/> 13"-16.5"L ¹	
<input type="checkbox"/> 11"-16.5"L ²	<input type="checkbox"/> 15.5"-19.5"L ²	
<input type="checkbox"/> Angle Adjustable/Multi-Axis Footplate Size:		
<input type="checkbox"/> 4"W x 6"D	<input type="checkbox"/> 5"W x 8"D	<input type="checkbox"/> 6"W x 11"D ³
<input type="checkbox"/> 5"W x 6"D	<input type="checkbox"/> 6"W x 8"D ³	<input type="checkbox"/> 8"W x 11"D ⁴
<input type="checkbox"/> Standard, Extruded Footplate Size:		
<input type="checkbox"/> 7.25"W x 6"D	<input type="checkbox"/> 8.5"W x 6"D ³	
1. Not available with multi-axis footplates.	3. May cause interference on seat widths less than 18"	
2. Not available w/ Angle Adjustable Footplates	4. May cause interference on seat widths less than 20"	



Quantum Q6 Edge 2.0 with TRU-Balance®3 Power Positioning

Description	Price	HCPCS Code
Non-Powered Legrests (cont.)		

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

Accessories:

<input type="checkbox"/> Heel Loops (pair) <i>Standard on Angle Adjustable and Multi-Axis footplates</i>	\$38	E0951
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair)	\$80	K0108
<input type="checkbox"/> Therafin Multi-Axis Amputee Support (each) <input type="checkbox"/> Left <input type="checkbox"/> Right <i>Swing-Away legrest upper assembly and receiver is required</i>	\$230	E1020
<input type="checkbox"/> Amputee Support Gel Cover (each) <input type="checkbox"/> Left <input type="checkbox"/> Right	\$150	K0108

60-90° Swing-Away Legrest Compatibility Matrix	TRU-Balance®3 Recline			TRU-Balance®3 Tilt			TRU-Balance®3 Lift		
	Low	Medium	High	Low	Medium	High	Low	Medium	High
Standard, 13" - 16.5"	!	!	✓	!	✓	✓	✓	✓	✓
Standard, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!
Angle-Adjustable, 13" - 16.5"	!	!	✓	!	!	!	!	!	!
Angle-Adjustable, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!
Multi-Axis, 15.5" - 19.5"	X	!	!	!	!	!	!	!	!

✓ This rigging Configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may interfere with the ground (less than 1.75" clearance when fully extended) and will reduce obstacle climbing ability.

X This rigging configuration does not work with this seat-to-floor height.

☐ **90° Pediatric Swing-Away Legrests** No Charge

Rigging Selection:

☐ Standard Matching Legrests (pair) ☐ Non-Matching Legrests (pair) ☐ Individual Legrests *Left or Right Only*

Footplate Selection:

☐ Angle-Adjustable Footplates (pair) *Includes heel loops.*

Left Legrest:**Lower Extensions:**

☐ 7"-10"L

Angle Adjustable Footplate Size:

☐ 4"W x 6"D ☐ 5"W x 6"D

Right Legrest:**Lower Extensions:**

☐ 7"-10"L

☐ 4"W x 6"D

☐ 5"W x 6"D

\$150 K0040

Accessories:

<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Elevating Legrests	\$325	E0990
Rigging Selection:		
<input type="checkbox"/> Standard Matching Legrests (pair) <input type="checkbox"/> Non-Matching Legrests (pair) <input type="checkbox"/> Individual Legrests <i>Left or Right Only</i>		
Footplate Selection:		
<input type="checkbox"/> Angle-Adjustable Footplates (pair) <i>Includes Heel Loops.</i>		
<input type="checkbox"/> Standard, Extruded Footplates (pair)		
Left Legrest Lower Extensions:		
<input type="checkbox"/> 11"-14"L ¹ <input type="checkbox"/> 14"-17"L ¹		
<input type="checkbox"/> 12"-14.5"L ² <input type="checkbox"/> 14.5"-18"L ²		
Right Legrest Lower Extensions:		
<input type="checkbox"/> 11"-14"L ¹ <input type="checkbox"/> 14"-17"L ¹		
<input type="checkbox"/> 12"-14.5"L ² <input type="checkbox"/> 14.5"-18"L ²		
Left Legrest:		
<input type="checkbox"/> 4"W x 6"D <input type="checkbox"/> 5"W x 8"D <input type="checkbox"/> 6"W x 11"D ³		
<input type="checkbox"/> 5"W x 6"D <input type="checkbox"/> 6"W x 8"D ³ <input type="checkbox"/> 8"W x 11"D ⁴		
Right Legrest:		
<input type="checkbox"/> 4"W x 6"D <input type="checkbox"/> 5"W x 8"D <input type="checkbox"/> 6"W x 11"D ³		
<input type="checkbox"/> 5"W x 6"D <input type="checkbox"/> 6"W x 8"D ³ <input type="checkbox"/> 8"W x 11"D ⁴		
Standard, Extruded Footplate Size:		
Left Legrest:		
<input type="checkbox"/> 7.25"W x 6"D <input type="checkbox"/> 8.5"W x 6"D ³		
Left Legrest:		
<input type="checkbox"/> 7.25"W x 6"D <input type="checkbox"/> 8.5"W x 6"D ³		

\$150 K0040

No Charge

Accessories:

<input type="checkbox"/> Heel Loops (pair) <i>Standard on Angle Adjustable and Multi-Axis footplates</i>	\$38	E0951
<input type="checkbox"/> Wheel Bumpers (pair)	\$25	
<input type="checkbox"/> Gel Padded Adductor Buttons (pair)	\$125	K0108

1. Not available with angle adjustable footplates.

2. Not available with standard, extruded footplates.

3. May cause interference on seat widths less than 18".

4. May cause interference on seat widths less than 20".