

**Jazzy® 1170 XL Plus  
Group 3**  
Order Form / Price List  
Revised: June 2011



**Pride Mobility Products Corporation®**

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at [www.mypride.com](http://www.mypride.com)

<b>Account #:</b>	<b>Date:</b>	<b>Order No.</b>	<b>Fax No.</b>	<b>Phone No.</b>
<b>Bill To:</b>	<b>Ship To:</b>		<b>P.O. #:</b>	Discount:
<b>Address:</b>	<b>Address:</b>			Terms:
<b>Phone Contact:</b>	<b>Marked For:</b>			Freight:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				All pricing is MSRP. Taken By:

**BASE MODEL**

- Jazzy 1170 3VHD-C K0853 \$9,753  
*Comes standard with solid tires*
- Jazzy 1170 3VHD-SS K0852 \$9,753  
*Comes standard with solid tires*

**COLOR**

- Champagne NC
- Candy Apple Red \$52
- Viper Blue \$52
- Onyx Black \$52

**BATTERIES (Installed)**

- Group 24, 70 amp Gel E2363 \$352ea

**BATTERY CHARGER**

- 8 amp On-board Charger Std.

**NON-EXPANDABLE CONTROLS**

- VSI Integral 70 amp<sup>23</sup> NC

**EXPANDABLE CONTROLS**

**CONTROLLER OPTIONS**

- PG Controller \$695  
**All joysticks must be used in conjunction w/a PG controller**
- Harness for Expandable Controller \$350

**JOYSTICK OPTIONS**

- PG Electronics**
- 3-Key Remote Plus 70 amp<sup>13,14</sup> \$700  
(Select Box for Profile)
  - 1 Profile  Multiple Profile
- 7-Key Access. Remote Plus 70 amp<sup>11,13,17</sup> \$700
- Remote Plus 70 amp with Speed Dial<sup>13,17</sup> \$700
- Remote Plus 70 amp with Toggles<sup>13,17</sup> \$700
- Remote Plus 70 amp with Switch Jacks<sup>13,17</sup> \$700
- Microdrive, 70 amp<sup>13,14,17</sup> \$700
- Omni Plus<sup>18</sup> \$1721.63

Joystick includes electronics & junction boxes.

Pride has submitted for the establishment of a new HCPCS code for the Upgraded Proportional Joystick. Call for clarification.

**PROGRAMMER**

- PG Programmer \$423

**CONTROLLER MOUNT**

- Fixed, In-line Joystick Mount Std.
  - Left  Right
- Swing-away Joystick Mount E1028 \$210
  - 4" <sup>8</sup>  7"
  - Left  Right
- Multi-axis Joystick Mount K0108 \$100
  - Left  Right  
(Fixed, does not swing-away)

**Jazzy® 1170 XL Plus  
Group 3**  
Order Form / Price List  
Revised: June 2011

**JOYSTICK HANDLE OPTIONS**

**PG Controls**

- Std. Handle Std.
- Mushroom Handle E2323 \$71
- T-Bar E2323 \$107
- Chin Cup E2324<sup>24</sup> \$71
- Extended I-Handle E2323 \$71
- Soft Ball Handle E2323 \$50
- Medium Ball Handle E2323 \$50
- Goal Post Handle E2323 \$66
  - 3" W  4" W

**NEW SWING-AWAY LEGRESTS**

- 70° Swing-away footrests Style #8<sup>24</sup> \$162
- 80° Swing-away footrests Style #8<sup>24</sup> \$162
- 70° Tapered S/A footrests Style #8<sup>24</sup> \$162  
(N/C with power positioning systems)
- Quad Release Lever NC
- S/A Standard Extensions: NC
  - 11"-13" Ext. NC
  - 13"-16.5" Ext. Std.
  - 15.5"-19.5" Ext. NC
  - 18.5"-22" Ext. NC
- (Lengths over 17" may cause interference)
- Extruded Footplate Std.
- 7.25" x 6" <sup>10</sup> (Tapered legs require 16"W or greater seat)
- 8.5 x 6" <sup>9</sup> (Tapered legs require 18"W or greater seat)
- Angle Adjustable Footplates (pair) K0040 \$150**
  - 4"x6" <sup>19</sup> (Tapered legs require 14"W or greater seat)
  - 5"x6" <sup>20</sup> (Tapered legs require 16"W or greater seat)
  - 5"x8" <sup>20</sup> (Tapered legs require 16"W or greater seat)
  - 6"x8" <sup>21</sup> (Tapered legs require 18"W or greater seat)
  - 6"x11" <sup>21</sup> (Tapered legs require 18"W or greater seat)
  - 8"x11" <sup>22</sup> (Tapered legs require 22"W or greater seat)
- Wheel Bumpers \$25
- Multi-Axis Ext. (ball & socket):
  - 11"-16.5" Ext.
  - 15.5"-19.5" Ext.
  - 18.5"-22" Ext.
- (Lengths over 17" may cause interference)
- Multi-Axis Footplates (pair) K0040 \$495**
  - 4"x6" <sup>19</sup> (Tapered legs require 14"W or greater seat)
  - 5"x6" <sup>20</sup> (Tapered legs require 16"W or greater seat)
  - 5"x8" <sup>20</sup> (Tapered legs require 16"W or greater seat)
  - 6"x8" <sup>21</sup> (Tapered legs require 18"W or greater seat)
  - 6"x11" <sup>21</sup> (Tapered legs require 18"W or greater seat)
  - 8"x11" <sup>22</sup> (Tapered legs require 22"W or greater seat)
- Wheel Bumpers \$25

**FRONT RIGGINGS**

- HD Foot Platform Std.
- No Front Rigging NC

- Rigging Type (pair)** E0990 \$325
- Elevating Leg Rests
  - Select footplate size for above:
    - 7 1/4 x 6  8 1/2 x 6
  - Adjustable Calf Pads w/ ELRS (1 hole std.)
  - 1 hole  3 holes\*  6 holes\*
  - \*Potential for interference w/ shrouds
  - Extensions:**
    - 11" - 14" Ext.
    - 14" - 17" Ext.
    - 17" - 20" Ext.
  - (Lengths over 17" may cause interference)

- Footplates (pair)**
- Reinforced Footplate Std.
    - 7 1/4 x 6  8 1/2 x 6
  - 5"x8" Angle-Adj. Footplate<sup>1</sup> K0040 \$150
  - 6"x8" Angle-Adj. Footplate<sup>1</sup> K0040 \$150
  - 6"x11" Angle-Adj. Footplate<sup>1</sup> K0040 \$150
  - 8"x11" Angle-Adj. Footplate<sup>1</sup> K0040 \$150
  - Wheel Bumpers (Ang-adj. only) \$25
  - Extensions<sup>7</sup>**
    - 11"-14" Ext. NC
    - 14"-17" Ext. Std.
    - 17"-20" Ext<sup>3</sup> NC

- Rigging Options**
- Heel Loops (pair)<sup>6</sup> E0951 \$38
  - Gel Neoprene Footrest Sleeve (pair)<sup>12</sup> K0108 \$65
  - ELR / ALR Adductor Buttons (pair) K0108 \$100
  - POSAlinc Knee Adductor (pair) K0108 \$464 for S/A legrests
  - Residual limb support left E1020 \$300
  - Residual limb support right E1020 \$300

**Jazzy® 1170 XL Plus  
Group 3**  
Order Form / Price List  
Revised: June 2011

**SEATING<sup>5</sup>**

High-Back w/ Headrest **Std.**

- 26" W x 20" D
- 26" W x 22" D
- 26" W x 26" D

Gray Vinyl       Beige Vinyl

High-Back w/ Headrest  
& Solid Seat Pan<sup>9</sup> **Std.**

- 26" W x 20" D
- 26" W x 22" D
- 26" W x 26" D

Gray Vinyl       Beige Vinyl

**22" W AND ABOVE WILL NOT BE LTD RECLINE.**

Custom seat dimensions are available.  
Please contact the Custom Group for  
pricing and leadtime.

**ARMRESTS**

**For High-Back Seats**

- Flip-up, removable **Std.**
- Flip-up, height-adjustable **E0973** \$170  
removable

**Armpads**

- Desk Length, 10" **NC**
- Left **NC**       Right **NC**
- Full Length, 14" **NC**
- Left **Std.**       Right **Std.**

**CUSHIONS<sup>16</sup>**

- Spectrum Cushion **(Visco only)**
- 26 x 20 - 28 x 26 **E2608** \$580

**ACCESSORIES**

- Rear Basket<sup>4</sup> \$103
- Cane/Crutch Holder<sup>4</sup> **E2207** \$98
- Oxygen Holder **E2208** \$150
- Walker Holder<sup>4</sup> **K0108** \$92
- Flag Holder \$17
- Lap Belt **Std.**
- 60"    70"    80"    90"

**Footnotes:**

1. Standard with heel loops.
2. Standard w/ heel loops and angle-adjustable footplates, check size below.
3. May cause interference w/ front anti-tips w/ certain seat widths and lower leg lengths. Call Quantum Customer Service.
4. Includes mounting bracket and hardware.
5. Seat widths over 26" wide will be wider than the base of the unit.
6. Standard with ELRs.
7. Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
8. Can only mount one armrest accessory in combination w/ S/A joystick.
9. Includes black vinyl seat pan, back available in all listed colors.
10. Std. w/ one programmable drive profile.
11. Needed for lighting packages.
12. N/A on ELRs.
13. Avail. w/five individual drive profile.
14. Requires multi-axis joystick mount; N/C for mount.
15. N/A with waterfall armpads and dual post armrests.
16. Please refer to cushion order forms for pricing, codes, and sizes.
17. Needed for controlling power positioning systems.
18. Please consult specialty controls order form for a variety of secondary controls and mounting options.
19. Only avail. on 12" W or greater seats.
20. Only avail. on 14" W or greater seats.
21. Only avail. on 16" W or greater seats.
22. Only avail. on 20" W or greater seats.
23. Std. w/ one programmable drive profile.
24. HCPCS code E2323 when not used in conjunction with chin control.

**Quantum Custom Delivery Service**

Includes: liftgate service and delivery of your product(s) inside your facility  
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.