

**Jazzy® 1420 Series
Group 2**
Order Form / Price List
Revised: June 2011



Pride Mobility Products Corporation®
182 Susquehanna Ave., Exeter, PA 18643
Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

BASE MODELS

- Jazzy 1420 2HD-C K0825 \$8,700
- Jazzy 1420 2VHD-C K0827 \$8,700
- Jazzy 1420 2HD-SS K0824 \$8,700
- Jazzy 1420 2VHD-SS K0826 \$8,700

COLOR

- Champagne NC
- Candy Apple Red NC
- Viper Blue NC
- Onyx Black NC
- Sunburst Orange NC

BATTERIES (Installed)

- Group 24, 70 amp, Gel E2363 \$352ea

BATTERY CHARGER

- 8 amp Off-board Charger Std.

NON-EXPANDABLE CONTROLS

- PG VR2 70 amp¹ Std.
- PG 6-key VR2 70 amp^{1,4} \$200

EXPANDABLE CONTROLS

All joysticks must be used in conjunction w/ a PG controller

CONTROLLER OPTIONS

- PG Controller \$695
- All joysticks or enhanced display must be used in conjunction w/ a PG controller**
- Harness for Expandable Controller \$350

JOYSTICK OPTIONS

PG Electronics

- 3-Key Remote Plus 70 amp^{1,3} \$700
(Select Box for Profile)
- 1 Profile Multiple Profile
- Remote Plus 70 amp with Speed Dial^{3,5} \$700
- Remote Plus 70 amp with Toggles^{3,5} \$700
- Remote Plus 70 amp with Switch Jacks^{3,5} \$700
- Microdrive, 70 amp^{3,5,6} \$700
- Omni Plus¹⁷ \$1721.63
(Requires a secondary control to operate the power chair)

Joystick includes electronics & junction boxes.

Pride has submitted for the establishment of a new HCPCS code for the Upgraded Proportional Joystick. Call for clarification.

PROGRAMMER

- PG Programmer \$423

JOYSTICK MOUNT

- Fixed, In-Line Joystick Mount Std.
- Left Right
- Swing-away Joystick Mount¹⁸ E1028 \$210
- 4" ¹⁹ 7"
- Left Right
- Multi-axis Joystick Mount K0108 \$100
- Left Right
- (Fixed, does not swing-away)*
- Vertical Adj. Clamp⁸ K0108 \$95
- Provides vertical joystick adjustment w/o adjusting armrest height.

JOYSTICK HANDLE OPTIONS

PG Controls

- Std. Knob Std.
- Mushroom Handle E2323 \$71
- T-Bar E2323 \$107
- Chin Cup E2324¹⁵ \$71
- Extended I-Handle E2323 \$71
- Soft Ball Handle E2323 \$50
- Medium Ball Handle E2323 \$50
- Goal Post Handle E2323 \$66
- 3" W 4" W

Jazzy® 1420 Group 2

Order Form / Price List
Revised: June 2011

DRIVE WHEELS

<input type="checkbox"/> 16" Pneumatic (pair)	Std
<input type="checkbox"/> Flat-Free Inserts (pair)	\$213

CASTER WHEELS

<input type="checkbox"/> Pneumatic (pair)	Std.
<input type="checkbox"/> Flat-Free Inserts (pair)	\$105

SWING-AWAY LEGRESTS

<input type="checkbox"/> 70° Swing-away footrests Style #8	\$162
<input type="checkbox"/> 80° Swing-away footrests Style #8	\$162
<input type="checkbox"/> 70° Tapered S/A footrests Style #8	\$162

(N/C with power positioning systems)

<input type="checkbox"/> Quad Release Lever	NC
---	----

S/A Standard Extensions:

<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext.	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
<input type="checkbox"/> 18.5"-22" Ext.	NC

(Lengths over 17" may cause interference)

<input type="checkbox"/> Extruded Footplate	Std.
---	------

<input type="checkbox"/> 7.25" x 6" 10 (Tapered legs require 18"W or greater seat)	
--	--

<input type="checkbox"/> 8.5 x 6" 9 (Tapered legs require 20"W or greater seat)	
---	--

Angle Adjustable Footplates (pair) \$150

<input type="checkbox"/> 4"x6" 16 (Tapered legs require 16"W or greater seat)	
---	--

<input type="checkbox"/> 5"x6" 10 (Tapered legs require 18"W or greater seat)	
---	--

<input type="checkbox"/> 5"x8" 10 (Tapered legs require 18"W or greater seat)	
---	--

<input type="checkbox"/> 6"x8" 9 (Tapered legs require 20"W or greater seat)	
--	--

<input type="checkbox"/> 6"x11" 9 (Tapered legs require 20"W or greater seat)	
---	--

<input type="checkbox"/> 8"x11" 11 (Tapered legs require 24"W or greater seat)	
--	--

<input type="checkbox"/> Wheel Bumpers	\$25
--	------

Multi-Axis Ext. (ball & socket):

<input type="checkbox"/> 11"-16.5" Ext.	
---	--

<input type="checkbox"/> 15.5"-19.5" Ext.	
---	--

<input type="checkbox"/> 18.5"-22" Ext.	
---	--

(Lengths over 17" may cause interference)

Multi-Axis Footplates (pair) \$495

<input type="checkbox"/> 4"x6" 16 (Tapered legs require 16"W or greater seat)	
---	--

<input type="checkbox"/> 5"x6" 10 (Tapered legs require 18"W or greater seat)	
---	--

<input type="checkbox"/> 5"x8" 10 (Tapered legs require 18"W or greater seat)	
---	--

<input type="checkbox"/> 6"x8" 9 (Tapered legs require 20"W or greater seat)	
--	--

<input type="checkbox"/> 6"x11" 9 (Tapered legs require 20"W or greater seat)	
---	--

<input type="checkbox"/> 8"x11" 11 (Tapered legs require 24"W or greater seat)	
--	--

<input type="checkbox"/> Wheel Bumpers	\$25
--	------

FRONT RIGGINGS

Rigging Type (pair)

<input type="checkbox"/> Elevating Leg Rests ²⁰	E0990/K0195	\$325
--	-------------	-------

Select footplate size for above:

<input type="checkbox"/> 71/4 x 6	<input type="checkbox"/> 81/2 x 6
-----------------------------------	-----------------------------------

Adjustable Calf Pads w/ ELRS (1 hole std.)

<input type="checkbox"/> 1 hole	<input type="checkbox"/> 3 holes*	<input type="checkbox"/> 6 holes*
---------------------------------	-----------------------------------	-----------------------------------

*Potential for interference w/ shrouds

Extensions:

<input type="checkbox"/> 11" - 14" Ext.	
---	--

<input type="checkbox"/> 14" - 17" Ext.	
---	--

<input type="checkbox"/> 17" - 20" Ext.	
---	--

(Lengths over 17" may cause interference)

<input type="checkbox"/> HD Drop-Ins, 70 ^{o2,7,13}	\$420
---	-------

<input type="checkbox"/> HD Drop-Ins, 80 ^{o2,7,13}	\$420
---	-------

N/A w/ contoured or Specialty seating

Angle Extensions: NC

<input type="checkbox"/> 12" - 14.5" Ext.	
---	--

<input type="checkbox"/> 14.5" - 18" Ext.	
---	--

<input type="checkbox"/> 17.5" - 21" Ext.	
---	--

(Lengths over 17" may cause interference)

<input type="checkbox"/> HD Foot Platform (Substitute for S/A)	NC
--	----

<input type="checkbox"/> No Front Rigging	NC
---	----

Footplates (pair)

<input type="checkbox"/> Reinforced Footplate	Std.
---	------

<input type="checkbox"/> 7.25" x 6" ¹⁰	<input type="checkbox"/> 8.5 x 6" ⁹
---	--

<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ¹⁰	\$150
---	-------

<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ⁹	\$150
--	-------

<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ⁹	\$150
---	-------

<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ¹¹	\$150
--	-------

<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$25
--	------

Extensions:¹²

<input type="checkbox"/> 11"-14" Ext.	NC
---------------------------------------	----

<input type="checkbox"/> 14"-17" Ext.	Std.
---------------------------------------	------

<input type="checkbox"/> 17"-20" Ext.	NC
---------------------------------------	----

<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	\$326
---	-------

<input type="checkbox"/> Small	<input type="checkbox"/> Large
--------------------------------	--------------------------------

<input type="checkbox"/> High mount Clamp-on; Multi-axis	\$412
--	-------

<input type="checkbox"/> Small	<input type="checkbox"/> Large
--------------------------------	--------------------------------

<input type="checkbox"/> High mount footplatform ¹⁹	K0108	\$450
--	-------	-------

HD Foot Platform
Lower Extensions

<input type="checkbox"/> 10"-14" L

<input type="checkbox"/> 14"-17" L

<input type="checkbox"/> 17"-20" L

Rigging Options

<input type="checkbox"/> Heel Loops (pair) ¹³	E0951	\$38
--	-------	------

<input type="checkbox"/> Gel Neoprene Footrest	K0108	
--	-------	--

Sleeve (pair) ¹⁴	\$65
-----------------------------	------

<input type="checkbox"/> ELR / ALR Adductor Buttons (pair)	K0108	\$100
--	-------	-------

<input type="checkbox"/> POSAlinc Knee Adductor (pair) for S/A legrests	K0108	\$464
--	-------	-------

<input type="checkbox"/> Residual limb support left	E1020	\$300
---	-------	-------

<input type="checkbox"/> Residual limb support right	E1020	\$300
--	-------	-------

Footnotes for pages 1 & 2

- Std. w/ one programmable drive profile.
- May cause interference w/ front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- Needed for controlling power positioning systems.
- For use with power elevating seat.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- N/A w/ contoured seating.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and/or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- HCPCS code E2323 when not used in conjunction with chin control.
- Only avail. on 12" W or greater seats.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- N/A with waterfall armpads and dual post armrests.
- E0990 Code used for purchase in CBA. K0195 code used for capped rental.

Jazzy® 1420 Group 2

Order Form / Price List
Revised: June 2011

SEATING

High-Back w/ Headrest Std.

- 22" W x 20" D
- 24" W x 20" D
- 24" W x 22" D

Gray Vinyl Beige Vinyl

High-Back w/ Headrest
& Solid Seat Pan¹ Std.

- 22" W x 20" D
- 24" W x 20" D
- 24" W x 22" D

Gray Vinyl Beige Vinyl

22" W AND ABOVE WILL NOT BE LTD RECLINE.

Custom seat dimensions are available.
Please contact the Custom Group for pricing and leadtime.

ARMRESTS

For High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable removable⁴ \$170 E0973

Armpads

- Desk Length, 10" NC
- Left NC Right NC
- Full Length, 14" NC
- Left **Std.** Right **Std.**

POWER ELEVATING SEAT^{4,6,7}

(Max weight capacity 300lbs.)

- Power Elevating Seat \$975
- Through 6-Key VR2 NC
- Power Elevating Seat \$1375
- Through Remote Controller
- Lap Belt Std.

MANUAL RECLINER^{2,5}

Manual Recline High Back E1226 \$600

Manual Recline High Back¹
w/ Solid Seat Pan E1226 \$600
(180 degrees of Recline)
(Dimensions 18"W x 18"D)

- | | |
|-------------------------------------|--------------------------------------|
| <input type="checkbox"/> Gray Vinyl | <input type="checkbox"/> Beige Vinyl |
| Width: | Depth: |
| <input type="checkbox"/> 16" NC | <input type="checkbox"/> 16" NC |
| <input type="checkbox"/> 18" Std. | <input type="checkbox"/> 18" Std. |
| <input type="checkbox"/> 20" NC | <input type="checkbox"/> 20" NC |

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS³

- Simplicity Cushion

16 x 16 - 20 x 20	E2601	\$98
22 x 16 - 24 x 24	E2602	\$98
- Solution 1 Cushion

14 x 14 - 20 x 18	E2603	\$380
20 x 20 - 20 x 24	E2603	\$430
22 x 16 - 24 x 24	E2604	\$430
- Solution Cushion

10 x 10 - 20 x 18	E2607	\$380
20 x 20 - 20 x 24	E2607	\$430
22 x 16 - 24 x 24	E2608	\$430
- Spectrum Cushion
 - Gel

10 x 10 - 20 x 18	E2605	\$480
20 x 20 - 20 x 24	E2605	\$530
22 x 16 - 24 x 24	E2606	\$530
 - Visco

10 x 10 - 20 x 18	E2607	\$480
20 x 20 - 20 x 24	E2607	\$530
22 x 16 - 24 x 24	E2608	\$530
- Spectrum Air

14 x 12 - 20 x 24	E2622	\$480
22 x 16 - 24 x 24	E2623	\$530
- Spectrum Air Contour

14 x 12 - 20 x 24	E2624	\$499
22 x 16 - 24 x 24	E2625	\$549

Footnotes for pages 3:

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special measurements are not available.
3. See Synergy Cushion order forms for detailed pricing, sizes & codes.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.
6. Includes appropriate actuator module.
7. Max. weight capacity is 300 lbs.

Jazzy® 1420 Group 2

Order Form / Price List
Revised: June 2011

Synergy Seating Systems

Synergy Seating™ System **\$434**

WIDTH

Set Seat Width At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W

DEPTH

Set Seat Depth At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"W-24"W

Please consult your Quantum Specialist or contact the Custom Group before ordering

BACK CANES¹⁴

<input type="checkbox"/> Back Canes w/ push handles	Std.
<input type="checkbox"/> Back Canes w/out push handles	\$125
<input type="checkbox"/> Offset (Meso) \$195	<input type="checkbox"/> Inset (Endo) \$195
<input type="checkbox"/> 2" (1" per side)	<input type="checkbox"/> 4" (2" per side)
N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests.	
N/A w/ any reclining systems.	

Actual Back Cane Height is reflected in measurement

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

BACK OPTIONS

<input type="checkbox"/> Synergy Back up to 20" W	\$275
<input type="checkbox"/> Synergy Back 21-24" W	\$325
<input type="checkbox"/> Omit Back	

Back Height:¹

<input type="checkbox"/> 18"H	<input type="checkbox"/> 19"H	<input type="checkbox"/> 20"H	<input type="checkbox"/> 21"H
<input type="checkbox"/> 22"H	<input type="checkbox"/> 23"H	<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H
<input type="checkbox"/> Sling Back	\$98		

Back height will be the same as back cane height. Includes gap for cushion. Push Bar selection must match back width; center tube standard when selecting flip-up arms

BACK OPTIONS CONT.

TRU-Comfort Back^{9,13}
(Even sizes only 16" W - 22"W)

<input type="checkbox"/> Std.	<input type="checkbox"/> New Plus	
Back Height:	<input type="checkbox"/> 22"	<input type="checkbox"/> 25"
<input type="checkbox"/> Back up to 20" W	E2620	\$595
<input type="checkbox"/> 22" W	E2621	\$695
<input type="checkbox"/> TRU-Comfort Cushion		
16 x 16 - 20 x 18	E2601	\$380
20 x 20 - 20 x 22	E2601	\$380
22 x 16 - 22 x 24	E2602	\$480
<input type="checkbox"/> TRU-Comfort Plus Cushion		
16 x 16 - 20 x 18	E2603	\$380
20 x 20 - 20 x 22	E2603	\$380
22 x 16 - 22 x 24	E2604	\$480
<input type="checkbox"/> Omit Back		-\$100

See TRU-Comfort Back & Cushion order forms for specific sizes

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

ARMRESTS

<input type="checkbox"/> Flip-back, cane mounted ²	Std.
<input type="checkbox"/> 2-Post, flip-back, height ²	\$150ea
adjustable, removable armrests	E0973
(9"-13" of height adjustment in 1" increments)	
<input type="checkbox"/> Single Post height adj. ^{7,8}	E0973 \$75ea
quick height, removable	
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.	
<input type="checkbox"/> Quad Release Lever	NC
<input type="checkbox"/> Removable, height adjustable	\$75ea
heavy duty armrests ¹²	E0973
<input type="checkbox"/> 8"H-12"H	<input type="checkbox"/> 10"H-14"H
<input type="checkbox"/> Cantilever Armrests	SOF
Angle-Adjustable: - 30° through + 30°	
Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting aftermarket arm troughs or Universal Swivel w/ desk length arms.	

ARMPADS

<input type="checkbox"/> Full length, Straight, 2" X 14"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Straight, 2" X 10"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Full Length, Waterfall, 2" X 14" ⁴	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Waterfall, 2" X 10" ^{4,5}	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right

HEAD REST

<input type="checkbox"/> Stealth	E0955
<input type="checkbox"/> Large (11"W x 6" H)	<input type="checkbox"/> Small (9 3/4" W x 3" H) \$205
<input type="checkbox"/> Stealth Comfort Plus	E0955
<input type="checkbox"/> Small (10"W)	<input type="checkbox"/> Large (14" W) \$205
<input type="checkbox"/> Removable Head Rest Bracket	E1028 \$210
<input type="checkbox"/> Stealth, swing-away (Flip-Down)	E1028
Head Rest Bracket	\$210
Interfaces w/ any headrest that uses 1/2" square down-tube.	
<input type="checkbox"/> Medical Applications	E0955
<input type="checkbox"/> Large (11"W x 6" H)	<input type="checkbox"/> Small (9 3/4" W x 3" H) \$230
<input type="checkbox"/> Removable Head Rest Bracket	E1028 \$235
<input type="checkbox"/> Fixed, Offset Mounting Bracket	NC
<input type="checkbox"/> TRU-Comfort	E0955
<input type="checkbox"/> 8" W	<input type="checkbox"/> Leatherette \$205
<input type="checkbox"/> 10" W	<input type="checkbox"/> Plus \$205
<input type="checkbox"/> Cobra Removable Head Rest Bracket	E1028 \$210
<input type="checkbox"/> Removable Bracket for after market headrest	E1028 \$210
<input type="checkbox"/> No Headrest Required	

ACCESSORIES

<input type="checkbox"/> Rear Basket ^{10,11}		\$103
<input type="checkbox"/> Cane/Crutch Holder ^{10,11}	E2207	\$98
<input type="checkbox"/> Oxygen Holder ^{10,11}	E2208	\$103
<input type="checkbox"/> Oxygen Holder ^{10,11}	E2208	\$150
(Synergy Seat only)		
<input type="checkbox"/> Walker Holder ^{10,11}	K0108	\$92
<input type="checkbox"/> Flag Holder (SOFO on Power Positioning)		\$17
<input type="checkbox"/> Weather Cover, Large, Black		\$202
<input type="checkbox"/> Saddle bag (contoured seating only)		\$28
<input type="checkbox"/> Cup holder		\$17
<input type="checkbox"/> Left	<input type="checkbox"/> Right	
N/A on desk length arm with S/A mount or toggle box.		
<input type="checkbox"/> Lap Belt		Std.
<input type="checkbox"/> 50" ¹¹	<input type="checkbox"/> 60"	<input type="checkbox"/> 70"
<input type="checkbox"/> 80"	<input type="checkbox"/> 90"	

Footnotes for pages 4:

- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms.
- N/A w/ S/A joystick mount.
- Available w/ solid backs only.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Cane mounted, flip-back arms are N/A.
- Includes mounting bracket and hardware.
- N/A on manual recliner high-back.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Back height is measured from the top of a 3" cushion on the seat.
- Height is measured from seat pan to top of cane.

Lowest seat-to-floor height

Synergy Seat:	18.125"
High-back Seat:	22.5"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility. Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.