

**Jazzy® 600 Series  
Group 2**

Order Form / Price List  
Revised: June 2011



**Pride Mobility Products Corporation®**

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at [www.mypride.com](http://www.mypride.com)

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

**BASE MODEL**

- Jazzy 600 2S-C K0823 \$5,363
- Jazzy 600 2S-SS K0822 \$5,363

**JOYSTICK MOUNT**

- Fixed, In-Line Joystick Mount Std.
  - Left  Right
- Swing-away Joystick Mount<sup>16</sup> E1028 \$210
  - 4" <sup>18</sup>  7"
  - Left  Right
- Quick Disconnect, Joystick Mount<sup>20,22</sup> NC
  - Left  Right
- Multi-axis Joystick Mount K0108 \$100
  - Left  Right
  - (Fixed, does not swing-away)
- Vertical Adj. Clamp<sup>4,21</sup> K0108 \$95

**SWING-AWAY LEGRESTS**

- 70° Swing-away footrests Style #8 \$162
- 70° Tapered S/A footrests Style #8 \$162  
(N/C with power positioning systems)
- Quad Release Lever NC
- S/A Standard Extensions:** NC
  - 11"-13" Ext. NC
  - 13"-16.5" Ext. Std.
  - 15.5"-19.5" Ext. NC
  - (Lengths over 17" may cause interference)
- Extruded Footplate Std.
- 7.25" x 6" <sup>6</sup> (Tapered legs require 18"W or greater seat)
- 8.5 x 6" <sup>5</sup> (Tapered legs require 20"W or greater seat)

**COLORS**

- Metallic Colors:**
- Champagne NC
  - Candy Apple Red \$52
  - Viper Blue \$52
  - Onyx Black \$52
  - Pewter \$52
  - Silver \$52
  - Emerald Green \$52
  - Graphite Blue \$52
  - Sapphire Blue \$52
- Solid Colors:**
- Kelly Green \$52
  - Tangerine \$52

**JOYSTICK HANDLE OPTIONS**

- Std. Knob Std.
- Mushroom Handle E2323 \$71
- T-Bar E2323 \$107
- Chin Cup E2324<sup>9</sup> \$71
- Extended I-Handle E2323 \$71
- Soft Ball Handle E2323 \$50
- Medium Ball Handle E2323 \$50
- Goal Post Handle E2323 \$66
  - 3" W  4" W

**Angle Adjustable Footplates (pair) \$150**

- 4"x6" <sup>10</sup> (Tapered legs require 16"W or greater seat)
- 5"x6" <sup>6</sup> (Tapered legs require 18"W or greater seat)
- 5"x8" <sup>6</sup> (Tapered legs require 18"W or greater seat)
- 6"x8" <sup>5</sup> (Tapered legs require 20"W or greater seat)
- 6"x11" <sup>5</sup> (Tapered legs require 20"W or greater seat)
- 8"x11" <sup>7</sup> (Tapered legs require 24"W or greater seat)
- Wheel Bumpers \$25

**BATTERIES (2 required)**

- NF-22, 50 AGM (Installed) E2361 \$255
- NF-22, 50 Gel Cell (Installed) E2361 \$295

**BATTERY CHARGER**

- 5 amp off-board Charger NC

**Multi-Axis Ext. (ball & socket):**

- 11"-16.5" Ext.
- 15.5"-19.5" Ext.
- (Lengths over 17" may cause interference)
- Multi-Axis Footplates (pair) \$495**
- 4"x6" <sup>10</sup> (Tapered legs require 16"W or greater seat)
- 5"x6" <sup>6</sup> (Tapered legs require 18"W or greater seat)
- 5"x8" <sup>6</sup> (Tapered legs require 18"W or greater seat)
- 6"x8" <sup>5</sup> (Tapered legs require 20"W or greater seat)
- 6"x11" <sup>5</sup> (Tapered legs require 20"W or greater seat)
- 8"x11" <sup>7</sup> (Tapered legs require 24"W or greater seat)
- Wheel Bumpers \$25

**PROGRAMMER**

- PG Programmer \$423
- Q-Logic Programmer \$423

**NON-EXPANDABLE CONTROLS**

- PG 4-Key VR2 60 amp Std.
- NE 75 amp <sup>14</sup> NC
- Requires the selection of the 5 amp off-board charger
- NE + 6-Key 75 amp <sup>14,15</sup> \$200

Requires the selection of the 5 amp off-board charger

**DRIVE WHEELS**

- 14" Pneumatic (pair)<sup>17</sup> NC
- 14" Flat-Free (pair) \$213

**FRONT RIGGINGS**

- Foot Platform Std.
- Rigging Type (pair)**
- Elevating Leg Rests E0990 \$325
  - Select footplate size for above:
    - 71/4 x 6  81/2 x 6
  - Adjustable Calf Pads w/ ELRS (1 hole std.)
    - 1 hole  3 holes\*  6 holes\*
  - \*Potential for interference w/ shrouds
- Extensions:**
  - 11" - 14" Ext.
  - 14" - 17" Ext.
- (Lengths over 17" may cause interference)

## FRONT RIGGINGS CONT.

<input type="checkbox"/> HD Drop-Ins, 70° <sup>3,4,13</sup>	\$420
N/A w/ contoured or Specialty seating	
<b>Angle Extensions:</b>	NC
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	
(Lengths over 17" may cause interference)	

<input type="checkbox"/> Multi-Axis Ext. (ball & socket):	
<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	
(Lengths over 17" may cause interference)	

## Footplates (pair)

<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 71/4 x 6 <sup>6</sup>	
<input type="checkbox"/> 81/2 x 6 <sup>5</sup>	
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate <sup>6</sup>	\$150
<b>Extensions<sup>2</sup></b>	
<input type="checkbox"/> 11"-14" Ext. <sup>24</sup>	NC
<input type="checkbox"/> 14"-17" Ext.	Std.

<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	\$326
<input type="checkbox"/> Small	
<input type="checkbox"/> Large	
<input type="checkbox"/> High mount Clamp-on; Multi-axis	\$412
<input type="checkbox"/> Small	
<input type="checkbox"/> Large	

<input type="checkbox"/> High mount footplatform <sup>15</sup>	<b>K0108</b>	\$450
<input type="checkbox"/> Mid-wheel 6 foot platform		
Lower Extensions		
<input type="checkbox"/> 10"-14" L		
<input type="checkbox"/> 14"-17" L		

## Rigging Options

<input type="checkbox"/> Heel Loops (pair) <sup>13</sup>	<b>E0951</b>	\$38
<input type="checkbox"/> Gel Neoprene Footrest		
Sleeve (pair) <sup>2</sup>	<b>K0108</b>	\$65
<input type="checkbox"/> ELR / ALR Adductor Buttons (pair)	<b>K0108</b>	\$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair)		
for S/A legrests	<b>K0108</b>	\$464
<input type="checkbox"/> Residual limb support left	<b>E1020</b>	\$300
<input type="checkbox"/> Residual limb support right	<b>E1020</b>	\$300

## SEATING

<input type="checkbox"/> Med-Back	Std.
-----------------------------------	------

<input type="checkbox"/> Med-Back	
& Solid Seat Pan <sup>1</sup>	Std.
<input type="checkbox"/> 18" W x 18" D	
<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl

<input type="checkbox"/> Ltd. Recline, Gray Vinyl, High-Back	
w/ Headrest , <b>Depth Adjustable</b>	\$275
<input type="checkbox"/> 18" W x 18" - 20" D	
<input type="checkbox"/> 20" W x 18" - 20" D	

<input type="checkbox"/> Ltd. Recline High-Back w/ Headrest	
& Solid Seat Pan <sup>1</sup>	\$275
<input type="checkbox"/> 18" W x 18" D	
<input type="checkbox"/> 18" W x 20" D	
<input type="checkbox"/> Gray Vinyl	<input type="checkbox"/> Beige Vinyl

Custom seat dimensions are available.  
Please contact the Custom Group for pricing and leadtime.

## ARMRESTS

## For Medium and High-Back Seats

<input type="checkbox"/> Flip-up, removable		Std.
<input type="checkbox"/> Flip-up, height-adjustable	<b>E0973</b>	\$170
removable <sup>19</sup>		

## Armpads

<input type="checkbox"/> Desk Length, 10"		NC
<input type="checkbox"/> Left	NC	<input type="checkbox"/> Right
		NC
<input type="checkbox"/> Full Length, 14"		NC
<input type="checkbox"/> Left	Std.	<input type="checkbox"/> Right
		Std.

## Special Order Seat

W x D

Seat sizes can not exceed 22 W x 20D  
**\$450 MSRP**

MANUAL RECLINER <sup>5</sup>

<input type="checkbox"/> Manual Recline High Back	<b>E1226</b>	\$600
---	--------------	-------

<input type="checkbox"/> Manual Recline High Back <sup>1</sup>		
w/ Solid Seat Pan	<b>E1226</b>	\$600
(180 degrees of Recline)		
(Dimensions 18"W x 18"D)		

<input type="checkbox"/> Gray Vinyl		<input type="checkbox"/> Beige Vinyl
<b>Width:</b>		<b>Depth:</b>
<input type="checkbox"/> 16" NC		<input type="checkbox"/> 16" NC
<input type="checkbox"/> 18" Std.		<input type="checkbox"/> 18" Std.
<input type="checkbox"/> 20" NC		<input type="checkbox"/> 20" NC

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS<sup>17</sup>

<input type="checkbox"/> Simplicity Cushion		
16 x 16 - 20 x 20	<b>E2601</b>	\$98
22 x 16 - 24 x 24	<b>E2602</b>	\$98
<input type="checkbox"/> Solution 1 Cushion		
14 x 14 - 20 x 18	<b>E2603</b>	\$380
20 x 20 - 20 x 24	<b>E2603</b>	\$430
22 x 16 - 24 x 24	<b>E2604</b>	\$430
<input type="checkbox"/> Solution Cushion		
10 x 10 - 20 x 18	<b>E2607</b>	\$380
20 x 20 - 20 x 20	<b>E2607</b>	\$430
22 x 16 - 22 x 20	<b>E2608</b>	\$430
<input type="checkbox"/> Spectrum Cushion		
<input type="checkbox"/> Gel		
10 x 10 - 20 x 18	<b>E2605</b>	\$480
20 x 20 - 20 x 20	<b>E2605</b>	\$530
22 x 16 - 22 x 20	<b>E2606</b>	\$530
<input type="checkbox"/> Visco		
10 x 10 - 20 x 18	<b>E2607</b>	\$480
20 x 20 - 20 x 20	<b>E2607</b>	\$530
22 x 16 - 22 x 20	<b>E2608</b>	\$530
<input type="checkbox"/> Spectrum Air		
14 x 12 - 20 x 24	<b>E2622</b>	\$480
22 x 16 - 24 x 24	<b>E2623</b>	\$530
<input type="checkbox"/> Spectrum Air Contour		
14 x 12 - 20 x 24	<b>E2624</b>	\$499
22 x 16 - 24 x 24	<b>E2625</b>	\$549

## Footnotes for pages 1 &amp; 2:

- Includes black vinyl seat pan, back available in all listed colors..
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- May cause interference with front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Pediatric foot platform is N/A w/ 16"-18" & 18"-20" extensions.
- HPCPS code E2323 when not used in conjunction with chin control.
- Please note swing away footrests are standard when ordered with Power Positioning
- Only avail. on 12" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ five individual drive profiles.
- Needed for controlling power positioning
- Can only mount one armrest accessory in combination w/ S/A joystick.
- See Synergy Cushion order forms for detailed pricing, sizes & codes.
- N/A with waterfall armpads and dual post armrests.
- N/A on manual recliner.
- N/A w/ 2-post, desk length arms w/ a toggle.
- Provides vertical joystick adjustment w/o adjusting armrest height.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.

## Jazzy® 600 Series Group 2

Order Form / Price List  
Revised: June 2011

### SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) \$275

#### SYSTEM SPECIFICATIONS

##### Specialty Seating Kit Width

16" NC  
 18" Std.  
 20" NC

##### Specialty Seating Kit Depth

16" Std.  17" NC  
 18" NC  19" NC  
 20" NC

##### Specialty Seating Kit Back Cane Height

18" Std.  20" NC  
 22" NC

##### Back Selection

See Synergy Shape Back Order Form<sup>1,2,3</sup>

Synergy Back up to 20"W<sup>18</sup> \$275

Omit Back

See Positioning Components Order Form  
for additional backs

##### Back Height Solid Pan<sup>4</sup>

18" Std.  19" NC  
 20" NC  21" NC  
 22" NC  23" NC  
 24" NC  25" NC

##### Sling Back

16" wide \$98  
 18" wide \$98  
 20" wide \$98

##### Sling Back Height

(Hgt must match cane hgt.)

18" height NC  
 20" height NC  
 22" height NC

**Push Bar selection must match back width; center tube standard when selecting flip-up arms.**

#### ARMRESTS

2-Post, fixed, flip-back<sup>3,5</sup> Std.

Flip-back, cane mounted<sup>3,5</sup> \$95ea

Removable, adjustable,  
quick height<sup>6,19</sup> E0973 \$275ea  
11"-17" of height adj. 3" on top in 1/2" increments  
& 3" on bottom in 1" increments.

Removable, height adjustable<sup>19</sup> E0973 \$150ea  
heavy duty armrests

8"H-12"H  10"H-14"H

##### Armrest Pads

Full Length, straight, 2"x14" NC

Left  Right

Desk Length, straight, 2"x10" NC

Left  Right

Full Length, Waterfall, 2"x14"<sup>7</sup> NC

Left  Right

Desk Length, Waterfall, 2"x10"<sup>7,17</sup> NC

Left  Right

#### HEAD REST

Stealth E0955  
 Large (11"W x 6" H)  Small (9 3/4" W x 3" H) \$205

Stealth Comfort Plus E0955  
 Large (14"W)  Small (10" W) \$205

Removable Head Rest Bracket E1028 \$210

Stealth, swing-away (Flip-Down)  
Head Rest Bracket E1028 \$210

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications E0955  
 Large (11"W x 6" H)  Small (9 3/4" W x 3" H) \$230

Removable Head Rest Bracket E1028 \$235

### Synergy Seating Systems<sup>11</sup>

Synergy Seating™ Frame \$434

(On chair MSRP)

Std. Stem Mount NC

2" Extended Vertical Stem Mounts NC

**(Cannot achieve lowest Seat to Floor Height)**

Synergy Manual Recline Pkg.<sup>10,12,13,14,15,2</sup> \$700

Maximum weight 300 lbs. E1226

Synergy Manual Tilt Pkg.<sup>13,16</sup> \$550

(Max Dimensions 14"W x 14"D) K0108

Maximum weight 75 lbs.

**Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.**

#### WIDTH

##### Set Seat Width At:

10" NC  11" NC  
 12" NC  13" NC  
 14" NC  15" NC  
 16" NC  17" NC  
 18" Std.  19" NC  
 20" NC

##### Seat Width Package

10"W-14"W  14"W-18"W  
 16"W-20"W

##### Seat Width Package Manual Tilt

10"W-12"W  12"W-14"W  
 14"W-16"W  16"W-18"W

#### DEPTH

##### Set Seat Depth At:

10" NC  11" NC  
 12" NC  13" NC  
 14" NC  15" NC  
 16" Std.  17" NC  
 18" NC  19" NC  
 20" NC

##### Seat Depth Package

10"D-14"D  14"D-18"D  
 16"D-20"D

##### Seat Depth Package Manual Tilt

10"D-14"D  14"D-16"D  
 16"D-18"D

**Please consult your Quantum Specialist or contact the Custom Group before ordering**

#### BACK CANES<sup>2</sup>

Back Canes w/ push handles Std.  
 Back Canes w/out push handles \$125  
 Flip Forward<sup>21</sup> \$150  
 Offset (Meso) \$195  Inset (Endo) \$195  
 2" (1" per side)  4" (2" per side)  
**N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests. N/A w/ any reclining systems.**

**Actual Back Cane Height is reflected in measurement**

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

##### Back Cane Height

14" NC  15" NC  
 16" NC  17" NC  
 18" Std.  19" NC  
 20" NC  21" NC  
 22" NC  23" NC  
 24" NC  25" NC

#### Footnotes for pages 4:

- 14" & 22" width are N/A.
- Height is measured from seat pan to top of cane.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Available w/ solid backs only.
- N/A w/ manual or PES systems.
- N/A with TRU-Comfort
- Manual systems do not include a headrest.
- N/A w/ flip-back, cane mounted arms.
- Synergy Back not included.
- Includes removable, adjustable, quick-height armrests.
- N/A w/ Shape Back.
- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- N/A w/ S/A joystick mount.
- 20"D is N/A.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and waterfall armpads.

# Jazzy® 600 Series Group 2

Order Form / Price List

Revised: June 2011

## BACK OPTIONS

See Synergy Shape Back Order Form<sup>1</sup>

Synergy Back up to 20" W \$275

Omit Back

**Back Height:**<sup>2</sup>

18"H  19"H  20"H  21"H

22"H  23"H  24"H  25"H

Sling Back<sup>3</sup> \$98

*Back height will be the same as back cane height.*

*Includes gap for cushion.*

*Push Bar selection must match back width; center tube standard when selecting flip-up arms.*

TRU-Comfort Back<sup>14,1</sup>

(Even sizes only 16" W - 20" W)

Std.  New Plus

**Back Height:**  22"  25"

Back up to 20" W **E2620** \$595

TRU-Comfort Cushion

16 x 16 - 20 x 18 **E2601** \$380

20 x 20 - 20 x 22 **E2601** \$380

22 x 16 - 22 x 24 **E2602** \$480

TRU-Comfort Plus Cushion

16 x 16 - 20 x 18 **E2603** \$380

20 x 20 - 20 x 22 **E2603** \$380

22 x 16 - 22 x 24 **E2604** \$480

Omit Back -(\$100)

See TRU-Comfort Back & Cushion order forms for specific sizes

*TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs*

## ARMRESTS

Flip-back, cane mounted<sup>4,5,6</sup> Std.

2-Post, flip-back, height<sup>4,5</sup> **E0973** \$150ea adjustable, removable armrests

(9"-13" of height adjustment in 1" increments)

Single Post height adj.<sup>7,12</sup> **E0973** \$75ea quick height, removable

11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.

Quad Release Lever NC

Removable, height adjustable **E0973** \$75ea heavy duty armrests<sup>12</sup>

8"H-12"H  10"H-14"H

Pediatric removable<sup>12</sup> **E0973** \$75ea adjustable, quick height.

Match Shroud Color  Black

Quad Release Lever NC  
(7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)

Cantilever Armrests<sup>20</sup> SOF

Angle-Adjustable: - 30° through + 30°

## ARMPADS

Full length, Straight, 2" X 14" Std.

Left  Right

Desk Length, Straight, 2" X 10" Std.

Left  Right

Full Length, Waterfall, 2" X 14"<sup>8</sup> NC

Left  Right

Desk Length, Waterfall, 2" X 10"<sup>8,9</sup> NC

Left  Right

## HEAD REST

Stealth **E0955**

Large (11"W x 6" H)  Small (9 3/4" W x 3" H) \$205

Stealth Comfort Plus **E0955**

Large (14"W)  Small (10" W) \$205

Removable Head Rest Bracket **E1028** \$210

Stealth, swing-away (Flip-Down)

Head Rest Bracket **E1028** \$210

Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications **E0955**

Large (11"W x 6" H)  Small (9 3/4" W x 3" H) \$230

Removable Head Rest Bracket **E1028** \$235

Fixed, Offset Mounting Bracket NC

TRU-Comfort **E0955**

8" W  Leatherette \$205

10" W  Plus \$205

Removable Head Rest Bracket **E1028** \$210

Removable Bracket for **E1028** \$210 after market headrest

No Headrest Required

## ACCESSORIES

Rear Basket<sup>17,18</sup> \$103

Cane/Crutch Holder<sup>17,18</sup> **E2207** \$98

Oxygen Holder<sup>17,18</sup> **E2208** \$103

Oxygen Holder<sup>14,17,18,19</sup> **E2208** \$150

Stroller handle \$85

(Synergy Manual Tilt only)

Walker Holder<sup>17,18</sup> **K0108** \$92

Flag Holder \$17

Weather Cover, Medium, Black \$158

Saddle bag(contoured seating only) \$28

Cup holder \$17

Left  Right

N/A on desk length arm with S/A mount or toggle box.

Lap Belt Std.

Pediatric 38" N/A w/ limited recline seating

50"<sup>13</sup>  60"  70"  80"  90"

### Lowest seat-to-floor height

Specialty Seat: 17"

Synergy Seat: 17"

High-back Seat: 21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires).

Increased widths may increase seat-to-floor height. Contact customer service for details.

Listed measurements are in the lowest position. It does not include elevating systems.

### Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility  
Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.

- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ Synergy manual tilt or recline.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- S/A laterals interfere w/ flip-back arms.
- N/A w/ Synergy manual recline.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- N/A w/ S/A joystick mount.
- Available w/ solid backs only.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- Back height is measured from the top of a 3" cushion on the seat.
- Height is measured from seat pan to top of cane.
- Includes mounting bracket and hardware.
- N/A on manual recliner high-back
- Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
- Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting after market arm troughs or Universal Swivel w/ desk length arms.