

Jazzy® 600 Series Group 3

Order Form / Price List
Revised: June 2011



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at www.mypride.com

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:			P.O. #:
Address:	Address:			Discount:
Phone Contact:	Marked For:			Terms:
DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.				Freight:
				Taken By:

All pricing is MSRP.

BASE MODEL

- Jazzy 600 3S-C **K0849** \$5,363
- Jazzy 600 3S-SS **K0848** \$5,363
- Jazzy 600 3SP-SS **K0856** \$5,363
- Jazzy 600 3SP-C **K0857** \$5,363

COLORS

Metallic Colors:

- Champagne NC
- Candy Apple Red \$52
- Viper Blue \$52
- Onyx Black \$52
- Pewter \$52
- Silver \$52
- Emerald Green \$52
- Graphite Blue \$52
- Sapphire Blue \$52

Solid Colors:

- Kelly Green \$52
- Tangerine \$52

BATTERIES (2 required)

- NF-22, 50 AGM (Installed) **E2361** \$255
- NF-22, 50 Gel Cell (Installed) **E2361** \$295

BATTERY CHARGER

- 5 amp Off-board Charger Std.

PROGRAMMER

- PG Programmer \$423
- Q-logic Programmer \$423

NON-EXPANDABLE CONTROLS

- NE 75 amp⁵ Std.
- NE+ 6-Key 75 amp^{3,5} \$200
(Tilt through joystick or switch)

EXPANDABLE CONTROLS

All joysticks must be used in conjunction w/ a PG controller

CONTROLLER OPTIONS

- PG Controller **E2377** \$695
- Harness for Expandable Controller **E2313** \$350

All joysticks or enhanced display must be used in conjunction w/ a PG controller

JOYSTICK OPTIONS

Pride has submitted for the establishment of a new HCPCS code for the Upgraded Proportional Joystick. Call for clarification.

PG Electronics

- 3-Key Remote Plus 70 amp³ \$700
(Select Box for Profile)
 - 1 Profile
 - Multiple Profile
- 7-KeyAccess. Remote Plus 70 amp^{3,4} \$700
- Remote Plus 70 amp with Speed Dial ^{3,5} \$700
- Remote Plus 70 amp with Toggles^{3,5} \$700
- Remote Plus 70 amp with Switch Jacks ^{3,5} \$700
- Microdrive, 70 amp ^{3,5,6} \$700
- Omni Plus²³ \$1721.63

(Requires a secondary control to operate the power chair)

Joystick includes electronics & junction boxes.

JOYSTICK MOUNT

- Fixed, In-Line Joystick Mount Std.
 - Left
 - Right
- Swing-away Joystick Mount¹⁶ \$210
 - 4" ¹⁸
 - 7" **E1028**
 - Left
 - Right
- Quick Disconnect, Joystick Mount^{22,26} NC
 - Left
 - Right
- Multi-axis Joystick Mount \$100
 - Left
 - Right
 (Fixed, does not swing-away) **K0108**
- Vertical Adj. Clamp^{8,27} **K0108** \$95

JOYSTICK HANDLE OPTIONS

- Std. Knob Std.
- Mushroom Handle **E2323** \$71
- T-Bar **E2323** \$107
- Chin Cup **E2324⁴** \$71
- Extended I-Handle **E2323** \$71
- Soft Ball Handle **E2323** \$50
- Medium Ball Handle **E2323** \$50
- Goal Post Handle **E2323** \$66
 - 3" W
 - 4" W

DRIVE WHEELS

- 14" Pneumatic (pair)¹⁷ NC
- 14" Flat-Free (pair) \$213

Jazzy® 600 Series Group 3

Order Form / Price List
Revised: June 2011

NEW SWING-AWAY LEGRESTS

<input type="checkbox"/> 70° Swing-away footrests Style #8 ²⁹	\$162
<input type="checkbox"/> 70° Tapered S/A footrests Style #8 ²⁹	\$162
(N/C with power positioning systems)	
<input type="checkbox"/> Quad Release Lever	NC

<input type="checkbox"/> S/A Standard Extensions:	NC
<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext.	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
(Lengths over 17" may cause interference)	

<input type="checkbox"/> Extruded Footplate	Std.
<input type="checkbox"/> 7.25" x 6" ¹⁰ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8.5 x 6" ¹⁹ (Tapered legs require 20"W or greater seat)	

Angle Adjustable Footplates ³⁰ (pair) K0040 \$150

<input type="checkbox"/> 4"x6" ²⁸ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 24"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

<input type="checkbox"/> Multi-Axis Ext. (ball & socket)	
<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	
(Lengths over 17" may cause interference)	

Multi-Axis Footplates (pair): K0040 \$495

<input type="checkbox"/> 4"x6" ²⁸ (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x6" ¹⁰ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 5"x8" ¹⁰ (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x8" ⁹ (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 6"x11" ⁹ (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 8"x11" ¹¹ (Tapered legs require 24"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

<input type="checkbox"/> Foot Platform ²⁵	Std.
--	------

Rigging Type (pair) E0990 \$325

<input type="checkbox"/> Elevating Leg Rests	E0990	\$325
Select footplate size for above:		
<input type="checkbox"/> 71/4 x 6	<input type="checkbox"/> 81/2 x 6	
Adjustable Calf Pads w/ ELRS (1 hole std.)		
<input type="checkbox"/> 1 hole	<input type="checkbox"/> 3 holes*	<input type="checkbox"/> 6 holes*

*Potential for interference w/ shrouds

Extensions:

<input type="checkbox"/> 11" - 14" Ext.	
<input type="checkbox"/> 14" - 17" Ext.	
(Lengths over 17" may cause interference)	

FRONT RIGGINGS

<input type="checkbox"/> HD Drop-Ins, 70 ^{07,8,13}	\$ 420
N/A w/ contoured or Specialty seating	
Angle Extensions:	NC
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	
(Lengths over 17" may cause interference)	

<input type="checkbox"/> 90° Pediatric Swing-Away footrests Style #8 ²⁹	\$162
Extension Range 7" - 10"	
Footplates, Ang. Adj.: <input type="checkbox"/> 4"x6" ²	<input type="checkbox"/> 5"x6" ¹⁰

<input type="checkbox"/> Power-elevating, high-pivot ^{19,22}	E1010	\$1695
---	-------	--------

Legrests, (calf pads included)²¹

<input type="checkbox"/> Independent / Combined	<input type="checkbox"/> Combined ¹⁹
---	---

<input type="checkbox"/> Power-elevating,	
---	--

Articulating Legrests ^{19,20,22}	E1010	\$2395
Articulating Range: <input type="checkbox"/> 4" L		Std.
<input type="checkbox"/> 6" L (Drop-in only)		NC

Set at: <input type="checkbox"/> 70°	
--------------------------------------	--

<input type="checkbox"/> Independent / Combined ¹⁹	<input type="checkbox"/> Combined ¹⁹
---	---

<input type="checkbox"/> Swing-away Power-elevating,	
--	--

Articulating Legrests ^{14,19,20,24}	E1010	\$2395
Articulating Range: <input type="checkbox"/> 4" L		Std.
<input type="checkbox"/> 6" L (Drop-in only)		NC

Set at: <input type="checkbox"/> 70°	
--------------------------------------	--

<input type="checkbox"/> Quad Release Lever	NC
---	----

<input type="checkbox"/> Independent / Combined ¹⁹	<input type="checkbox"/> Combined ¹⁹
---	---

<input type="checkbox"/> Calf Pads Left & Right (6 x 5.35)	Std.
--	------

<input type="checkbox"/> Calf Pads Left & Right (8.00 x 5.86)	NC
---	----

<input type="checkbox"/> Power Articulating Foot Platform K0108	\$2,895
---	---------

<input type="checkbox"/> Synchronized	<input type="checkbox"/> Non-Synchronized
---------------------------------------	---

w/ split ang. adj. foot platform
-8" of Articulation

- Elevation range 90°-170°

Adjustment Packages	NC
---------------------	----

<input type="checkbox"/> 12"-14"	<input type="checkbox"/> 14"-16"
----------------------------------	----------------------------------

<input type="checkbox"/> 16"-18" (depth adj)	<input type="checkbox"/> 18"-20" (depth adj)
--	--

Select Foot platform for above: Adult Pediatric¹

(Lengths over 17" may cause interference)

18" - 20" Adj. Package requires a minimum STF of 20" to accommodate for ground clearance

<input type="checkbox"/> Connection Bracket	\$49
---	------

<input type="checkbox"/> Heel Cups	K0108	\$39
------------------------------------	-------	------

Footplates (pair)

<input type="checkbox"/> Reinforced Footplate	Std.
---	------

<input type="checkbox"/> 71/4 x 6 ¹⁰	<input type="checkbox"/> 81/2 x 6 ⁹
---	--

<input type="checkbox"/> 5"x8" Angle-Adj. Footplate ¹⁰	K0040	\$150
---	-------	-------

<input type="checkbox"/> 6"x8" Angle-Adj. Footplate ⁹	K0040	\$150
--	-------	-------

<input type="checkbox"/> 6"x11" Angle-Adj. Footplate ⁹	K0040	\$150
---	-------	-------

<input type="checkbox"/> 8"x11" Angle-Adj. Footplate ¹¹	K0040	\$150
--	-------	-------

<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$150
--	-------

Extensions¹²

<input type="checkbox"/> 11"-14" Ext. ²⁴ (\$125 w/ power legs)	NC
---	----

<input type="checkbox"/> 14"-17" Ext.	Std.
---------------------------------------	------

<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	K0040	\$326
---	-------	-------

<input type="checkbox"/> Small	<input type="checkbox"/> Large
--------------------------------	--------------------------------

<input type="checkbox"/> High mount Clamp-on; Multi-axis	K0040	\$412
--	-------	-------

<input type="checkbox"/> Small	<input type="checkbox"/> Large
--------------------------------	--------------------------------

<input type="checkbox"/> High mount footplatform ¹⁵	K0108	\$450
--	-------	-------

Mid-wheel 6 foot platform Lower Extensions

<input type="checkbox"/> 10"-14" L	
------------------------------------	--

<input type="checkbox"/> 14"-17" L	
------------------------------------	--

FRONT RIGGINGS CONT.

Rigging Options		
<input type="checkbox"/> Heel Loops (pair) ¹³	E0951	\$38
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	K0108	\$65
<input type="checkbox"/> ELR / ALR Adductor Buttons (pair)	K0108	\$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair) for S/A legrests	K0108	\$464
<input type="checkbox"/> Residual limb support left	E1020	\$300
<input type="checkbox"/> Residual limb support right	E1020	\$300

Footnotes for pages 1 & 2:

- Pediatric foot platform is N/A w/ 16"-18 & 18"-20" extensions.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- Needed for controlling power positioning systems.
- HPCPS code E2323 when not used in conjunction with chin control.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- May cause interference with front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and /or cushion.
- Std. w/ ang-adj. footplates.
- Please note swing away footrests are standard when ordered with Power Positioning.
- Std. calf pad is not depth adjustable.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- If ordered w/ reference guides or configuration packages discounts may apply.
- N/A with waterfall armpads and dual post armrests.
- Requires Remote Plus.
- Power legrests are only available on TRU-Balance and Synergy Seating Systems. Only the combined & independent motion can be selected on the Synergy Seat.
- Std. calf pad is not depth adjustable.
- N/A w/ waterfall armpads; only available on Quantum Seating Systems and only one option maybe mounted.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- 11" - 14" Ext. are N/A w/ power S/A legrests.
- Fixed bracket, does not elevate.
- N/A w/ 2-post, desk length arms w/ a toggle.
- Provides vertical joystick adjustment w/o adjusting armrest height.
- Only avail. on 12" W or greater seats.
- Please note swing away footrests are standard when ordered with Power Positioning.
- Standard w/ heel loops.

Jazzy® 600 Series Group 3

Order Form / Price List
Revised: June 2011

SEATING

Med-Back Std.

Med-Back
& Solid Seat Pan ¹ Std.
 18" W x 18" D
 Gray Vinyl Beige Vinyl

Ltd. Recline, Gray Vinyl, High-Back
w/ Headrest ² **Depth Adjustable**
 18" W x 18" - 20" D
 20" W x 18" - 20" D

Ltd. Recline High-Back w/ Headrest
& Solid Seat Pan ¹ \$275
 Gray Vinyl Beige Vinyl
 18" W x 18" D
 18" W x 20" D

Custom seat dimensions are available.
Please contact the Custom Group for pricing and leadtime.

Special Order Seat

W x D

Seat sizes can not exceed 22 W x 20D
\$450 MSRP

ARMRESTS

For Medium and High-Back Seats

Flip-up, removable Std.
 Flip-up, height-adjustable **E0973** \$170
removable⁴

Armpads

Desk Length, 10" NC
 Left NC Right NC
 Full Length, 14" NC
 Left Std. Right Std.

MANUAL RECLINER ⁵

Manual Recline High Back **E1226** \$600

Manual Recline High Back ¹ **E1226** \$600
w/ Solid Seat Pan
(180 degrees of Recline)
(Dimensions 18"W x 18"D)

Gray Vinyl Beige Vinyl

Width:		Depth:	
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 16"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 18"	Std.
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 20"	NC

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

CUSHIONS ³

Simplicity Cushion
16 x 16 - 20 x 20 **E2601** \$98
22 x 16 - 24 x 24 **E2602** \$98
 Solution 1 Cushion
14 x 14 - 20 x 18 **E2603** \$380
20 x 20 - 20 x 24 **E2603** \$430
22 x 16 - 24 x 24 **E2604** \$430

Solution Cushion
10 x 10 - 20 x 18 **E2607** \$380
20 x 20 - 20 x 20 **E2607** \$430
22 x 16 - 22 x 20 **E2608** \$430

Spectrum Cushion
 Gel
10 x 10 - 20 x 18 **E2605** \$480
20 x 20 - 20 x 20 **E2605** \$530
22 x 16 - 22 x 20 **E2606** \$530

Visco
10 x 10 - 20 x 18 **E2607** \$480
20 x 20 - 20 x 20 **E2607** \$530
22 x 16 - 22 x 20 **E2608** \$530

Spectrum Air
14 x 12 - 20 x 24 **K0734** \$480
22 x 16 - 24 x 24 **K0735** \$530
 Spectrum Air Contour
14 x 12 - 20 x 24 **K0736** \$499
22 x 16 - 24 x 24 **K0737** \$549

No Cushion Required

Footnotes for pages 3:

1. Includes black vinyl seat pan, back available in all listed colors.
2. Must choose 60" lapbelt.
3. See Synergy Cushion order forms for detailed pricing, sizes & codes.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.

Jazzy® 600 Series Group 3

Order Form / Price List

Revised: June 2011

SPECIALTY SEATING

 Specialty Seat Kit (limited adjustability) \$275

SYSTEM SPECIFICATIONS

Specialty Seating Kit Width

- 16" NC
 18" Std.
 20" NC

Specialty Seating Kit Depth

- 16" Std. 17" NC
 18" NC 19" NC
 20" NC

Specialty Seating Kit Back Cane Height

- 18" Std. 20" NC
 22" NC

Back Selection

 See Synergy Shape Back Order Form^{1,2,3}
 Solid Plate Back up to 20"W¹⁸ \$275

 Omit Back

 See Positioning Components Order Form for additional backs

Back Height Solid Pan⁴

- 18" Std. 19" NC
 20" NC 21" NC
 22" NC 23" NC
 24" NC 25" NC

 Sling Back

- 16" wide \$98
 18" wide \$98
 20" wide \$98

Sling Back Height

(Hgt must match cane hgt.)

- 18" height NC
 20" height NC
 22" height NC

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

ARMRESTS

 2-Post, fixed flip-back^{3,5} Std.

 Flip-back, cane mounted^{3,5} \$95ea

 Removable, adjustable, quick height^{6,19} E0973 \$275ea
 11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.

 Removable, height adjustable¹⁹ E0973 \$150ea
 heavy duty armrests

- 8"H-12"H 10"H-14"H

Armrest Pads

 Full Length, straight, 2"x14" NC

- Left Right

 Desk Length, straight, 2"x10" NC

- Left Right

 Full Length, Waterfall, 2"x14"⁷ NC

- Left Right

 Desk Length, Waterfall, 2"x10"^{7,17} NC

- Left Right

HEAD REST

- Stealth E0955
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
 Stealth Comfort Plus E0955
 Large (14"W) Small (10" W) \$205
 Removable Head Rest Bracket E1028 \$210

 Stealth, swing-away (Flip-Down)
 Head Rest Bracket E1028 \$210

Interfaces w/ any headrest that uses 1/2" square down-tube.

- Medical Applications E0955
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230
 Removable Head Rest Bracket E1028 \$235
 Fixed, Offset Mounting Bracket NC

Synergy Seating Systems¹¹

Synergy SeatingTM Frame \$434
(On chair MSRP)
 Std. Stem Mount NC

2" Extended Vertical Stem Mounts NC
(Cannot achieve lowest Seat to Floor Height)

Synergy Manual Recline Pkg.^{9,12,13,14,15,20} \$700
 Maximum weight 300 lbs. E1226

Synergy Manual Tilt Pkg.^{13,16} K0108 \$550
(Max Dimensions 14"W x 14"D)
 Maximum weight 75 lbs.

Synergy Manual Systems come standard with mercury drive lock-out and are in addition to the MSRP of the Synergy Seat.

WIDTH

Set Seat Width At:

- 10" NC 11" NC
 12" NC 13" NC
 14" NC 15" NC
 16" NC 17" NC
 18" Std. 19" NC
 20" NC

Seat Width Package

- 10"W-14"W 14"W-18"W
 16"W-20"W

Seat Width Package Manual Tilt

- 10"W-12"W 12"W-14"W
 14"W-16"W 16"W-18"W

DEPTH

Set Seat Depth At:

- 10" NC 11" NC
 12" NC 13" NC
 14" NC 15" NC
 16" Std. 17" NC
 18" NC 19" NC
 20" NC

Seat Depth Package

- 10"D-14"D 14"D-18"D
 16"D-20"D

Seat Depth Package Manual Tilt

- 10"D-14"D 14"D-16"D
 16"D-18"D

Please consult your Quantum Specialist or contact the Custom Group before ordering

BACK CANES²¹

- Back Canes w/ push handles Std.
 Back Canes w/out push handles \$125
 Flip Forward²² \$150
 Offset(Meso) \$195 Inset(Endo) \$195
 2" (1" per side) 4" (2" per side)
N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .
N/A w/ any reclining systems.

Actual Back Cane Height is reflected in measurement

Synergy Manual Recline comes standard with 22" non-height adjustable back canes.

Back Cane Height

- 14" NC 15" NC
 16" NC 17" NC
 18" Std. 19" NC
 20" NC 21" NC
 22" NC 23" NC
 24" NC 25" NC

Footnotes for pages 4:

- 14" & 22" width are N/A.
- Micro-drive, Shape Back & high-mount foot platform are N/A in the same configuration.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- Available w/ solid backs only.
- N/A w/ TRU-Comfort
- N/A w/ manual or PES systems.
- Manual systems do not include a headrest.
- N/A w/ flip-back, cane mounted arms.
- Synergy Back not included.
- Includes removable, adjustable, quick-height armrests.
- N/A w/ Shape Back.
- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- N/A w/ S/A joystick mount.
- 20"D is N/A.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- Solid plate back is N/A w/manual recline.
- Height is measured from seat pan to top of cane.
- N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and waterfall armpads..

Jazzy® 600 Series Group 3

Order Form / Price List
Revised: June 2011

BACK OPTIONS

See Synergy Shape Back Order Form¹
 Synergy Back up to 20" W \$275

Omit Back
Back Height:²

18"H 19"H 20"H 21"H
 22"H 23"H 24"H 25"H

Sling Back³ \$98

Back height will be the same as back cane height. Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

TRU-Comfort Back^{14,17}
(Even sizes only 16" W - 20"W)

Std. New Plus

Back Height: 22" 25"
 Back up to 20" W E2620 \$595

TRU-Comfort Cushion
16 x 16 - 20 x 18 E2601 \$380
20 x 20 - 20 x 22 E2601 \$380
22 x 16 - 22 x 24 E2602 \$480

TRU-Comfort Plus Cushion
16 x 16 - 20 x 18 E2603 \$380
20 x 20 - 20 x 22 E2603 \$380
22 x 16 - 22 x 24 E2604 \$480
 Omit Back -(\$100)

See TRU-Comfort Back & Cushion order forms for specific sizes

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

ARMRESTS

Flip-back, cane mounted^{4,5,6} Std.
 2-Post, flip-back, height^{4,5} E0973 \$150ea
adjustable, removable armrests

(9"-13" of height adjustment in 1" increments)
 Single Post height adj.^{7,11} E0973 \$75ea
quick height, removable
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.

Quad Release Lever NC
 Removable, height adjustable E0973 \$75ea
heavy duty armrests¹¹

8"H-12"H 10"H-14"H
 Pediatric removable¹¹ E0973 \$75ea
adjustable, quick height.

Match Shroud Color Black
 Quad Release Lever NC
(7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)

Cantilever Armrests SOF
Angle-Adjustable: - 30° through + 30°

Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting aftermarket arm troughs or Universal Swivel w/ desk length arms.

ARMPADS

Full length, Straight, 2" X 14" Std.
 Left Right
 Desk Length, Straight, 2" X 10" Std.
 Left Right
 Full Length, Waterfall, 2" X 14"¹⁸ NC
 Left Right
 Desk Length, Waterfall, 2" X 10"^{8,9} NC
 Left Right

HEAD REST

Stealth E0955
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
 Stealth Comfort Plus E0955
 Large (14"W) Small (10" W) \$205
 Removable Head Rest Bracket E1028 \$210
 Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$210
Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications E0955
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230
 Removable Head Rest Bracket E1028 \$235
 Fixed, Offset Mounting Bracket NC

TRU-Comfort E0955
 8" W Leatherette \$205
 10" W Plus \$205
 Cobra Removable Head Rest Bracket E1028 \$210
 Removable Bracket for E1028 \$210
after market headrest

TRU-BALANCE TILT

TRU-Balance Power Tilt E1002 \$4995
55 degree center of gravity (CG) power tilt seating system. 250 lbs. weight capacity

SYSTEM SPECIFICATIONS

Tower Selection
 Adjustable-angle Std.
Will raise minimum STF by 7/8".
 Fixed NC

Set Seat Width At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Width Package

14"W-18"W 16"W-20"W

Set Seat Depth At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	

Seat Depth Package

14"D-18"D 16"D-20"D

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

Actual Back Cane Height is reflected in measurement

Set Back Cane Height At:

<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Seat Overall Back Height¹³

Back heights of 22"-28" include a 6" lateral extension. Please select if required

with other back heights

<input type="checkbox"/> 6" Lateral Extensions	NC
<input type="checkbox"/> 14" NC	<input type="checkbox"/> 15" NC
<input type="checkbox"/> 16" NC.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Please consult your Quantum Specialist or contact the Custom Group before ordering

BACK CANES¹⁸

Back Canes w/ push handles Std.
 Back Canes w/out push handles \$125
 Flip Forward¹⁹ \$150
 Offset(Meso) \$195 Inset (Endo) \$195
 2" (1" per side) 4" (2" per side)
N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .
N/A w/ any reclining systems.

Back Type

See Synergy Shape Back Order Form
 Synergy Back up to 20" W Std.
 Attendant manual recline^{14,15,16} \$750
180 degrees of recline with 2 1/2" raised pivot point. Controlled through dual gas struts. 300 lb. weight capacity.

TRU-Comfort Back^{14,17}
(Even sizes only 16" W - 20"W)

Std. New Plus

Back Height: 22" 25"
 Back up to 20" W E2620 \$595
 TRU-Comfort Cushion
16 x 16 - 20 x 18 E2601 \$380
20 x 20 - 20 x 22 E2601 \$380
22 x 16 - 22 x 24 E2602 \$480
 TRU-Comfort Plus Cushion
16 x 16 - 20 x 18 E2603 \$380
20 x 20 - 20 x 22 E2603 \$380
22 x 16 - 22 x 24 E2604 \$480
 Omit Back -(\$100)

See TRU-Comfort Back & Cushion order forms for specific sizes

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

Footnotes for page 5:

- Manual Tilt w/ Shape back & oxygen tank holder requires the use of a pushbar.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ Synergy manual tilt or recline.
- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- S/A laterals interfere w/ flip-back arms.
- N/A w/ Synergy manual recline.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- N/A w/ S/A joystick mount.
- Available w/ solid backs only.
- Desk length armrests w/ an accessory limit the accessibility to the height adj.lever.
- Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
- N/A w/ flip-back, cane mounted arms.
- N/A w/ Shape Back.
- Includes price for removable, adjustable, quick-height arms & 2-post arms.
- Back height is measured from the top of a 3" cushion on the seat.
- Height is measured from seat pan to top of cane.
- N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and waterfall armpads. 02 Tank Holder is N/A when flip-down canes are selected w/ any TB Tilt.

Jazzy® 600 Series Group 3

Order Form / Price List
Revised: June 2011

ARMRESTS

- Flip-back, cane mounted^{1,2,3} Std.
Mounted on back canes with full width spreader bar.
- 2-Post, flip-back, hgt. adj, removable^{1,3} NC
(9"-13" of height adjustment in 1" increments)
- Single Post height adj., removable^{4,18} NC
 Quad Release Lever NC
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.
- Removable, single-post, height adjustable, heavy duty¹⁸ NC
 8"-12" Adjustment
 10"-14" Adjustment
- Pediatric, removable,¹⁸ NC
adjustable, quick height
(7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)
- Match Shroud Color Black NC
 Quad Release Lever NC
- Cantilever Armrests²¹ SOF
Angle-Adjustable: - 30° through + 30°

Armrest Pads

- Full Length, straight, 2"x14" NC
 Left Right
- Desk Length, straight, 2"x10" NC
 Left Right
- Full Length, Waterfall, 2"x14"¹⁵ NC
 Left Right
- Desk Length, Waterfall, 2"x10"^{5,6} NC
 Left Right

POWER POSITIONING ELECTRONICS

Select Electronics to be used

- NE (Through Toggle)
 NE+ 6-Key (Through Controller)

SINGLE POWER ELECTRONICS

Electronics for Tilt Only

- Through Toggle Box⁷ NC
 Left Right (mounts to side of armrest)
- Through 6-Key NE+(Tilt Only) E2310 \$1175
- Through Single Switch(1/8" Jack)
Tilt Module^{8,9,14} \$500
- Through Remote Plus Controller E2310 \$1175
- Through 2-Directional Control Box
w/ dual Switch Jacks(1/8" Jack)^{8,9,14,19} \$1175

* TB2 Tilt seat function plus power riggings

- Through Remote Plus Controller E2311 \$2375

HEAD REST

- Stealth E0955
 Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$205
- Stealth Comfort Plus E0955
 Large (14"W) Small (10" W) \$205
- Removable Head Rest Bracket E1028 \$210
- Stealth, swing-away (Flip-Down)
Head Rest Bracket E1028 \$210
- Interfaces w/ any headrest that uses 1/2" square down-tube.

Medical Applications E0955

- Large (11"W x 6" H) Small (9 3/4" W x 3" H) \$230
- Removable Head Rest Bracket E1028 \$235
- Fixed, Offset Mounting Bracket NC

TRU-Comfort E0955

- 8" W Leatherette \$205
 10" W Plus \$205
- Removable Head Rest Bracket E1028 \$210
- Removable Bracket for E1028 \$210
after market headrest

ACCESSORIES

- Rear Basket^{12,13,15} \$103
- Cane/Crutch Holder^{12,13,15} E2207 \$98
- Oxygen Holder^{12,13,15} E2208 \$103
- Oxygen Holder^{12,13,16,17,20} E2208 \$150
(TRU-Balance2 Tilt & Synergy Seat only)
- Stroller handle \$85
(Synergy Manual Tilt only)
- Walker Holder^{12,13,15} K0108 \$92
- Flag Holder (SOF on Power Positioning) \$17
- Weather Cover, Medium, Black \$158
- Saddle bag(contoured seating only) \$28
- Cup holder \$17
 Left Right
N/A on desk length arm with S/A mount or toggle box.
- Lap Belt Std.
 Pediatric 38" N/A w/ limited recline seating
- 50"¹³ 60" 70" 80" 90"

Footnotes for pages 6:

- Shape Back & w/ flip-back arms may interfere in the same configuration, call for details.
- N/A w/ Attendant manual recline.
- S/A laterals interfere w/ flip-back arms.
- Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
- Only avail. on removable, adjustable, quick-height arms & pediatric arms.
- N/A w/ S/A joystick mount.
- N/A w/ waterfall armpads.
- PG electronics only.
- Consult specialty control order form for a wide variety of switch types and sizes.
- Comes with swing-away / removable hardware; HCPCS code is E1028.
- These items are sold by Quantum Rehab in a package, reflecting a single MSRP.
- Includes mounting bracket and hardware.
- N/A on manual recliner high-back
- Does not include switch.
- Not available on tilt systems.
- N/A w/ flip-back arms.
- Manual Tilt w/ Shape Back and oxygen holder requires the use of a pushbar.
- Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
- 2 switches required, & Remote Plus.
- TB2 Tilt w/ TRU-Comfort & O² holder is SOF.
- Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting after market arm troughs or Universal Swivel w/ desk length arms.

Lowest seat-to-floor height

Specialty Seat:	17"
Synergy Seat:	17"
TB 2 Tilt:	18.5"
(w/ Ang-Adj towers)	
TB 2 Tilt:	17.625"
(w/ fixed towers towers)	
High-back Seat:	21"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility. Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.