

# Jazzy 614HD

Group 2 Heavy Duty  
Order Form  
450lb Weight Capacity



Pride Mobility Products Corporation  
182 Susquehanna Ave., Exeter, PA 18643  
Phone: (866)800-2002 • Fax: (866) 707-3422

Discount: \_\_\_\_\_  
Terms: \_\_\_\_\_  
Freight: \_\_\_\_\_  
Taken By: \_\_\_\_\_



Account Number: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Ship To Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.*

BASE MODEL			JOYSTICK HANDLES		
Description	Price	HCPCS Code	Description	Price	HCPCS Code
<input type="checkbox"/> Jazzy 614HD 2HD-C	\$6,295	K0825	<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Jazzy 614HD 2HD-SS	\$6,295	K0824	<input type="checkbox"/> DL Knob	\$20	E2323
<b>COLOR</b>			<input type="checkbox"/> Sponge Ball	\$35	E2323
<input type="checkbox"/> Jazzy Red			<input type="checkbox"/> Wrist Plate	\$45	E2323
<input type="checkbox"/> Jazzy Blue			<input type="checkbox"/> Round T-Bar (Chrome)	\$45	E2323
<b>BATTERY CHARGER</b>			<input type="checkbox"/> Flat T-Bar (Chrome)	\$45	E2323
Description	Price		<input type="checkbox"/> Straight T-Bar (Chrome)	\$45	E2323
<input type="checkbox"/> 5 Amp Off-board	Standard		<input type="checkbox"/> Finger Tube (Chrome)	\$45	E2323
<b>DRIVE WHEELS</b>			<b>CONTOURED SEATING</b>		
Description	Price		Description	Price	HCPCS Code
<input type="checkbox"/> 14" Flat-Free (pair)	Standard		<b>Depth Adjustable High Back:</b>		
<b>BASE MOUNTED FOOT PLATFORM</b>			<input type="checkbox"/> Limited Recline Depth Adjustable	No Charge	
Description	Price		High Back w/ Headrest <sup>1</sup>		
<input type="checkbox"/> Base Mounted Foot Platform	Standard		Color: <input type="checkbox"/> Black Vinyl		
<input type="checkbox"/> Omit Base Mounted Foot Platform	No Charge		Size: <input type="checkbox"/> 18"W x 18-20"D	<input type="checkbox"/> 20"W x 18-20"D	<input type="checkbox"/> 22"W x 20-22"D
<b>NON-EXPANDABLE CONTROLS</b>			<b>High Back with Solid Seat Pan:</b>		
Description	Price		<input type="checkbox"/> Limited Recline High-Back with Headrest & Solid Seat Pan	No Charge	
<input type="checkbox"/> Shark Controller, 50 Amp	No Charge		Color: <input type="checkbox"/> Black Vinyl		
<b>PROGRAMMING</b>			Size: <input type="checkbox"/> 16"W x 16"D	<input type="checkbox"/> 16"W x 18"D	<input type="checkbox"/> 16"W x 20"D
Description	Price		<input type="checkbox"/> 16"W x 22"D	<input type="checkbox"/> 18"W x 16"D	<input type="checkbox"/> 18"W x 18"D
<input type="checkbox"/> Dynamic Programmer w/ Flight Adapter	\$423		<input type="checkbox"/> 18"W x 20"D	<input type="checkbox"/> 18"W x 22"D	<input type="checkbox"/> 20"W x 16"D
<b>JOYSTICK MOUNT</b>			<input type="checkbox"/> 20"W x 18"D	<input type="checkbox"/> 20"W x 20"D	<input type="checkbox"/> 20"W x 22"D
Description	Price	HCPCS Code	<b>Armrests:</b>		
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>			<input type="checkbox"/> Flip-Up, Removable	No Charge	
<input type="checkbox"/> Fixed Inline <sup>1</sup>	Standard		<input type="checkbox"/> Flip-Up, Height-Adjustable, Removable	\$85 ea.	E0973
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right			<b>Armpads:</b>		
<input type="checkbox"/> Swing-away	\$210	E1028	Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length <sup>2</sup>		
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right			Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length <sup>2</sup>		
<b>Special Order Seat</b>			_____ "W x _____ "D		
			<b>Seat size can not exceed 24"W x 24"D</b>		
			**Contact your inside sales specialist for pricing**		
1. Not available with Q-logic EX Stand Alone Joystick.			1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part		
			2. Selection of desk length armpads on a depth-adjustable high back requires them to be added as an additional part.		

**SYNERGY® SEATING SYSTEMS**

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy Seating Frame	\$434	
<b>Seat Width:</b>		
<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"		
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19"		
<input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24"		
<b>Seat Width Package:</b>		
<input type="checkbox"/> 10"W-14"W <input type="checkbox"/> 14"W-18"W <input type="checkbox"/> 16"W-20"W <input type="checkbox"/> 20"W-22"W		
<b>Seat Depth:</b>		
<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"		
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19"		
<input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24"		
<b>Seat Depth Package:</b>		
<input type="checkbox"/> 10"D-14"D <input type="checkbox"/> 14"D-18"D <input type="checkbox"/> 15"D-20"D <input type="checkbox"/> 20"D-22"D		
<b>Back Cane Height:</b>		
<input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <sup>1</sup> <input type="checkbox"/> 17" <sup>1</sup>		
<input type="checkbox"/> 18" <input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21"		
<input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24" <input type="checkbox"/> 25"		

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

**SYNERGY® ARMRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> Flip-back, Cane Mounted <sup>1</sup>	No Charge	
<input type="checkbox"/> 2-post, Flip-back, Height-adjustable, Removable <sup>2,3</sup>	\$150 ea.	E0973
<input type="checkbox"/> Single-post, Height-adjustable, Removable <sup>3,4</sup>	\$75 ea.	E0973
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Removable, Single-post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
<input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/> Pediatric, Removable <sup>5</sup>	\$75 ea.	E0973
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Bariatric, Removable, Height Adjustable <sup>6</sup>	\$75 ea.	E0973
<input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/> Height Adjustable Cantilever Armrests <sup>7</sup>	\$195	
Angle-Adjustable: -30° through +30°		

1. Interferes with swing-away laterals. Not available on flip forward canes or manual recline.  
 2. Interferes with swing-away laterals. Height adjustable from 9"-13" in 1" increments.  
 3. Not available with a depth less than 16".  
 4. Height adjustable from 11"-17" in 1/2" increments  
 5. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armpads not available  
 6. Bariatric Armrests require Bariatric Armpads  
 7. Height adjustment from 11-17" depends on cane height. Not available with Flip Forward back canes or Shape Back.

**SYNERGY® BACK OPTIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy® Shape Back <i>Positioning Back</i>	Various	E2613/E2614
<i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>		
<input type="checkbox"/> Sling Back	\$98	
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	\$275/\$325 <sup>1</sup>	
<i>Overall Back</i> <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
<i>Height:</i> <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
<i>Choose Lateral Extensions<sup>2</sup>:</i>		
<input type="checkbox"/> 6" for 13"-16"H Back    \$25		
<input type="checkbox"/> 10" for 17"-18"H Back    \$25		
<input type="checkbox"/> 12" For 19"-22"H Back    \$25		
<input type="checkbox"/> 16" For 23"-25"H Back    \$25		
<input type="checkbox"/> TRU-Comfort™ Back <i>4-way Stretch Material Cover</i>		
<i>Width:</i> <input type="checkbox"/> 16"W <input type="checkbox"/> 18" W <input type="checkbox"/> 20"W    \$595    E2620		
<i>Overall Back</i> <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
<i>Height<sup>3</sup>:</i> <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
<i>Width:</i> <input type="checkbox"/> 17"W <input type="checkbox"/> 19"W    \$595    E2620		
<i>Overall Back</i> <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
<i>Height<sup>3</sup>:</i> <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
<i>Width:</i> <input type="checkbox"/> 22"W <input type="checkbox"/> 24"W    \$695    E2621		
<i>Overall Back</i> <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
<i>Height<sup>3</sup>:</i>		

**Synergy® Back Cane Options:**

<input type="checkbox"/> Standard Back Canes	Standard
<input type="checkbox"/> Omit Push Handles    \$125	
<input type="checkbox"/> Flip Forward Back Canes <sup>4</sup>	\$150
<input type="checkbox"/> Omit Push Handles    \$125	
<input type="checkbox"/> Offset (Meso) <sup>5</sup>	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	
<input type="checkbox"/> Inset (Endo)	\$195
<input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side)	

1. Price reflects back widths 21"W and wider.  
 2. Lateral extensions are required to mount lateral thoracic supports.  
 3. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".  
 4. Not available with O2 Holder, cane mounted or cantilever armrests. May interfere with Lateral Supports and Waterfall Armpads.  
 5. Dual Post Flip-back Armrests interfere with offset canes.

**SYNERGY® ARMPADS**

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Full Length (2"x14") <sup>1</sup>	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Desk Length (2"x10") <sup>1</sup>	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$92.50 ea.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	
\$56.69 ea.	
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$70.00 ea.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	
\$53.52 ea.	
<input type="checkbox"/> Flat, Bariatric <sup>2</sup>	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.  
 2. Bariatric Armrests require Bariatric Armpads

**LAP BELTS**

Description	Price
<input type="checkbox"/> 50" Lap Belt <sup>1</sup>	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="checkbox"/> 60" Padded Lap Belt <sup>1</sup>	\$115
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard

1. Only available with Synergy seating.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**HEADRESTS** **FOOT RIGGING ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
<b>TRU-Comfort™ Headrest<sup>1</sup>:</b>		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
<b>Stealth:</b>		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket <sup>2</sup> (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350)	\$210	
Used in Addition to Removable Headrest Bracket		
<b>Aftermarket Bracket:</b>		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

1. Measurement is taken from the inside of the headrest.

**NON-POWER LEGRESTS**

Description	Price
<input type="checkbox"/> <b>70° Swing-away Legrests<sup>1</sup></b>	<b>\$162</b>
<input type="checkbox"/> Standard Extruded Footplates	No Charge
<i>Standard Extruded Extensions:</i>	
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"	
<i>Standard Extruded Footplates:</i>	
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"	
<input type="checkbox"/> Angle-Adjustable Footplates <sup>2</sup>	\$150
<i>Angle Adjustable Extensions:</i>	
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"	
<i>Angle-Adjustable Footplates:</i>	
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"	
<input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>	
<input type="checkbox"/> Multi-Axis Footplates	\$495
<i>Multi-Axis Extensions:</i>	
<input type="checkbox"/> 15.5"-19.5"	
<i>Multi-Axis Footplates:</i>	
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"	
<input type="checkbox"/> 6" x 8" <sup>3</sup> <input type="checkbox"/> 6" x 11" <sup>3</sup> <input type="checkbox"/> 8" x 11" <sup>4</sup>	

1. Possible caster wheel interference when adjusted to less than 70 degrees.  
 2. Includes Heel Loops  
 3. Causes interference on seat widths less than 18".  
 4. Causes interference on seat widths less than 20".

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) <sup>1</sup>	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) <sup>2</sup>	\$65	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$512	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair)*	\$574	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Bodypoint Ankle Huggers (pair)*	\$134	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium		

1. Standard with angle adjustable footplates.  
 2. Available only with swing-away legrests.

**SEAT CUSHIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Simplicity Cushion</b> General Use		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$98	E2602
<b>Solutrum® 1 Cushion</b> Skin Protection:		
<input type="checkbox"/> 14" x 14" - 20" x 22"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$430	E2604
<b>Solutrum® Cushion</b> Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 22"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$430	E2608
<b>Spectrum® Gel Cushion</b> Positioning		
<input type="checkbox"/> 12" x 10" - 20" x 22"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$530	E2606
<b>Spectrum® Visco Cushion</b> Skin Protection & Positioning		
<input type="checkbox"/> 10" x 10" - 20" x 22"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$530	E2608
<b>Spectrum® Air</b> Adjustable Skin Protection		
<input type="checkbox"/> 14" x 12" - 20" x 22"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$530	E2623
<b>Spectrum® Air Contour</b> Adjustable Skin Protection & Positioning		
<input type="checkbox"/> 14" x 12" - 20" x 22"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 22" x 22"	\$549	E2625
<b>TRU-Comfort™ 2 4-Way Stretch</b> Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
<b>TRU-Comfort™ 2 Air Mesh</b> Skin Protection		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**LATERALS**

Description	Price	HCPCS Code
<input type="checkbox"/> <b>Stealth Thoracic Laterals (pair)*</b>		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away Rail Mounts	\$450	E1028
<i>For Use on Synergy Back and Tru-Comfort</i>		
<i>Pad:</i>		
<input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> <b>Stealth Pelvic/Thigh Laterals*</b>		
<i>Side:</i>		
<input type="checkbox"/> Left Side		<input type="checkbox"/> Right Side
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$250	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$225	E1028
<input type="checkbox"/> Flip Down, Full Surface Contact Bracket (ea.)	\$225	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$50	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$50	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$50	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) <sup>2</sup>	\$25	K0108
<input type="checkbox"/> 2 1/4" Depth/Height Adapter Plate (ea.) <sup>3</sup>	\$25	K0108

1. Swing-away laterals interfere with flip-up armrests.  
 2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.  
 3. Only available on the flip down bracket.

**BATTERY OPTIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> NF-22, 50-Amp Batteries		
<input type="checkbox"/> Maintenance Free AGM	\$295 ea.	E2361
<i>Installation:</i>		
<input type="checkbox"/> Batteries Installed		<input type="checkbox"/> Do Not Install

**SEAT TO FLOOR HEIGHT**

Seating System	Seat-To-Floor Height			
Contoured Seat	21.5"	22.5"	23.5"	24.5"
Synergy Seat (Standard Stem Mounts)	17.75"	18.75"	19.75"	

*Seat to floor height are +/- 1/8" with 18"W x 18"D with Flat-Free tires. Increased widths may increase lowest seat to floor heights.*

**NOTES:**

**MISCELLANEOUS POSITIONING COMPONENTS**

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket*	\$172	E1028
<i>Pad:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$135	E0957
<input type="checkbox"/> Stealth Foot Boxes (pair)*	\$475	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)*	\$166	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/> Otto Bock Arm Channels (ea.)*	\$104	E2209
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Otto Bock Flat Hand Pad(ea.)*	\$52	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Medium <input type="checkbox"/> Large		

*Please see positioning components order form for additional types, sizes and details.*

**ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> Rear Basket <sup>1</sup>	\$103	
<input type="checkbox"/> Cane/Crutch Holder <sup>1</sup>	\$98	E2207
<input type="checkbox"/> Weather Cover, Medium, Black	\$158	
<input type="checkbox"/> Cup Holder	\$17 <sup>2</sup>	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Saddle Bag	\$28	
<input type="checkbox"/> Flag Holder	\$17	
<input type="checkbox"/> Quantum Backpack <sup>3</sup>	\$60	
<input type="checkbox"/> Oxygen Holder, Contoured Seating <sup>1</sup>	\$103	E2208
<input type="checkbox"/> Oxygen Holder, Synergy Seat <sup>4</sup>	\$150	E2208
<input type="checkbox"/> Walker Holder <sup>4</sup>	\$92	
<input type="checkbox"/> Tablet Holder <sup>5</sup>	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder <sup>5</sup>	\$140	

1. Includes mounting hardware. Not available on Manual Recliner.  
 2. No charge on Synergy seating.  
 3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.  
 4. Not available with Synergy manual recline or manual tilt. Includes mounting hardware.  
 5. Only available with Synergy or Specialty seating.

*Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum*