

# Jazzy Select 14

Order Form / Price List

Revised: June 2011



Pride Mobility Products Corporation®

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

Account #:	Date:	Order No.
Bill To:	Ship To:	P.O. #:
Address:	Address:	Discount:
Phone Contact:	Marked For:	Terms:
<p><b>This product is available for purchase online at <a href="http://www.mypride.com">www.mypride.com</a></b> All pricing is MSRP.</p>		Freight:
<p>DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.</p>		Taken By:

## SPECIFIC CONFIGURATION

- Jazzy Select 14-C K0823 \$5,096  
Comes standard with Flat-free Tires
- Jazzy Select 14-SS K0822 \$5,096  
Comes standard with Flat-free Tires

## COLOR

- Jazzy Red NC
- Jazzy Blue NC

## BATTERIES INSTALLED

- NF-22, AGM E2361 \$255ea

## CHARGER

- 5 amp off-board Std.

## PROGRAMMER

- Dynamic Programmer \$423

## JOYSTICK MOUNT

- Fixed, In-Line Joystick Mount Std.
  - Left  Right
- Swing-away Joystick Mount E1028 \$210
  - Left  Right
  - 4"  7"

## FRONT RIGGINGS

- Foot Platform Std.

### Rigging Type (pair)

Select footplate size for S/A & ELRs:

- 71/4 x 6  81/2 x 6

- Swing-away footrests \$162

- 11"-13" Ext. NC

- 13"-16.5" Ext. Std.

- 15.5"-19.5" Ext. NC

- Elevating Leg Rests E0990/K0195<sup>14</sup> \$325

Adjustable Calf Pads w/ ELRS (1 hole std.)

- 1 hole  3 holes\*  6 holes\*

\*Potential for interference w/ shrouds

### Extensions<sup>4</sup>

- 11"-14" Ext. NC

- 14"-17" Ext. Std.

### Rigging Options

- Heel Loops (pair) E0951 \$38

## SEATING

- Med-Back Gray Vinyl Std.

- Med-Back, Gray Vinyl & Solid Seat Pan<sup>2</sup> NC  
 18" W x 18" D

- Ltd. Recline, Gray Vinyl, High-Back w/ Headrest, **Depth Adjustable**<sup>5</sup> \$275

- 18" W x 18" - 20" D
- 20" W x 18" - 20" D

## SOLID SEAT PAN SELECTION

- Ltd. Recline, Gray Vinyl, High-Back w/ Headrest & Solid Seat Pan<sup>2</sup> \$275
  - 18" W x 18" D
  - 18" W x 20" D
  - 20" W x 18" D
  - 20" W x 20" D

Custom seat dimensions are available. Please contact the Custom Group for pricing and leadtime.

## ARMRESTS

### For High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable E0973 \$170 removable

### Armpads

- Desk Length, 10"
  - Left NC  Right NC
- Full Length, 14"<sup>13</sup>
  - Left Std.  Right Std.

## Special Order Seat

\_\_\_\_\_ W x \_\_\_\_\_ D

Seat sizes can not exceed 20 W x 20 D  
\$450 MSRP

# Jazzy Select 14

Order Form / Price List

Revised: June 2011

## CUSHIONS<sup>3</sup>

<input type="checkbox"/> Simplicity Cushion			
16 x 16 - 20 x 20	<b>E2601</b>		\$98
<input type="checkbox"/> Solution 1 Cushion			
14 x 14 - 20 x 18	<b>E2603</b>		\$380
20 x 20	<b>E2603</b>		\$430
<input type="checkbox"/> Solution Cushion			
10 x 10 - 20 x 18	<b>E2607</b>		\$380
20 x 20	<b>E2607</b>		\$430
<input type="checkbox"/> Spectrum Cushion			
<input type="checkbox"/> Gel			
10 x 10 - 20 x 18	<b>E2605</b>		\$480
20 x 20	<b>E2605</b>		\$530
<input type="checkbox"/> Visco			
10 x 10 - 20 x 18	<b>E2607</b>		\$480
20 x 20 - 20 x 20	<b>E2607</b>		\$530
<input type="checkbox"/> Spectrum Air			
14 x 12 - 20 x 24	<b>E2622</b>		\$480
22 x 16 - 24 x 24	<b>E2623</b>		\$530
<input type="checkbox"/> Spectrum Air Contour			
14 x 12 - 20 x 24	<b>E2624</b>		\$499
22 x 16 - 24 x 24	<b>E2625</b>		\$549

## ACCESSORIES

<input type="checkbox"/> Rear Basket <sup>1</sup>		\$103
<input type="checkbox"/> Cane/Crutch Holder <sup>1</sup>	<b>E2207</b>	\$98
<input type="checkbox"/> Oxygen Holder <sup>1</sup>	<b>E2208</b>	\$103
<input type="checkbox"/> Walker Holder <sup>1</sup>	<b>K0108</b>	\$92
<input type="checkbox"/> Weather Cover, Medium, Black		\$158
<input type="checkbox"/> Cup Holder		\$17
<input type="checkbox"/> Right	<input type="checkbox"/> Left	
<input type="checkbox"/> Lap Belt		Std.
<input type="checkbox"/> 38"	<input type="checkbox"/> 50"	<input type="checkbox"/> 60"
<input type="checkbox"/> 70"	<input type="checkbox"/> 80"	<input type="checkbox"/> 90"

## SPECIALTY SEATING

Specialty Seat Kit (limited adjustability) **\$275**

### SYSTEM SPECIFICATIONS

#### Specialty Seating Kit Width

<input type="checkbox"/> 16" NC		<input type="checkbox"/> 20" NC
<input type="checkbox"/> 18" Std.		

#### Specialty Seating Kit Depth

<input type="checkbox"/> 16" Std.	<input type="checkbox"/> 17" NC
<input type="checkbox"/> 18" NC	

#### Specialty Seating Kit Back Cane Height

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 20" NC
<input type="checkbox"/> 22" NC	

## BACKS

See Synergy Shape Back Order Form<sup>6</sup>

Solid Plate Back up to 20"W **\$275**

Omit Back

See Positioning Components Order Form for additional backs

#### Back Height Solid Pan <sup>7</sup>

<input type="checkbox"/> 18" Std.	<input type="checkbox"/> 19" NC
<input type="checkbox"/> 20" NC	<input type="checkbox"/> 21" NC
<input type="checkbox"/> 22" NC	<input type="checkbox"/> 23" NC
<input type="checkbox"/> 24" NC	<input type="checkbox"/> 25" NC

Sling Back

<input type="checkbox"/> 16" wide	\$98
<input type="checkbox"/> 18" wide	\$98
<input type="checkbox"/> 20" wide	\$98

#### Sling Back Height

(Hgt must match cane hgt.)	
<input type="checkbox"/> 18" height	NC
<input type="checkbox"/> 20" height	NC
<input type="checkbox"/> 22" height	NC

## ARMRESTS

<input type="checkbox"/> 2-Post, fixed flip-back <sup>8,9</sup>	Std.
<input type="checkbox"/> Flip-back, cane mounted <sup>8,9</sup>	\$95ea
<input type="checkbox"/> Removable, adjustable, quick height <sup>10</sup>	<b>E0973</b> \$275ea
11"-17" of height adj 3" on top in 1/2" increments & 3" on bottom in 1" increments	
<input type="checkbox"/> Removable, height adjustable <sup>10</sup>	<b>E0973</b> \$150ea
heavy duty armrests	
<input type="checkbox"/> 8"H-12"H	<input type="checkbox"/> 10"H-14"H

### Armrest Pads

<input type="checkbox"/> Full Length, straight, 2"x14"	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, straight, 2"x10"	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Full Length, Waterfall, 2"x14" <sup>11</sup>	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Waterfall, 2"x10" <sup>11,12</sup>	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right

## HEAD REST

<input type="checkbox"/> Stealth	<b>E0955</b>
<input type="checkbox"/> Large (11"W x 6" H)	<input type="checkbox"/> Small (9 3/4" W x 3" H) \$205
<input type="checkbox"/> Stealth Comfort Plus	<b>E0955</b>
<input type="checkbox"/> Large (14"W)	<input type="checkbox"/> Small (10" W) \$205
<input type="checkbox"/> Removable Head Rest Bracket	<b>E1028</b> \$210
<input type="checkbox"/> Stealth, swing-away (Flip-Down) Head Rest Bracket	<b>E1028</b> \$210
Interfaces w/ any headrest that uses 1/2" square down-tube.	

<input type="checkbox"/> Medical Applications	<b>E0955</b>
<input type="checkbox"/> Large (11"W x 6" H)	
<input type="checkbox"/> Small (9 3/4" W x 3" H)	\$210
<input type="checkbox"/> Removable Head Rest Bracket	<b>E1028</b> \$235

### Footnotes for pages 1 & 2

- Includes mounting bracket and hardware.
- Additional leadtimes apply.
- See cushion order form for pricing and sizes.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured seating or rehab seating w/ cushions are selected. The overall length may vary approximately 3" depending on choice of seat and /or cushion.
- Must choose a 60" lapbelt.
- Shape Back and w/ flip-back arms may interfere in the same configuration, call for details.
- 20" D is N/A.
- Call for availability with Shape back.
- S/A laterals interfere w/ flip-back arms.
- Desk length armrests w/ an accessory limit the accessibility to the height adjustment lever.
- N/A w/ flip back, cane mounted arms.
- S/A joystick mount N/A.
- Full length arms required on Comfort & Depth Adjustable High-Back seating options.
- E0990 Code used for purchase in CBA. K0195 code used for capped rental.

# Jazzy Select 14

Order Form / Price List

Revised: June 2011

## Synergy Seating Systems

 Synergy Seating™ System \$434

### WIDTH

Set Seat Width At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Width Package

<input type="checkbox"/> 14"W-18"W	<input type="checkbox"/> 16"W-20"W
------------------------------------	------------------------------------

### DEPTH

Set Seat Depth At:

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC		

Seat Depth Package

<input type="checkbox"/> 14"D-18"D	<input type="checkbox"/> 16"D-20"D
------------------------------------	------------------------------------

Please consult your Quantum Specialist or contact the Custom Group before ordering

### BACK CANES<sup>1</sup>

<input type="checkbox"/> Back Canes w/ push handles	Std.
<input type="checkbox"/> Back Canes w/out push handles	\$125
<input type="checkbox"/> Flip Forward <sup>13</sup>	\$150
<input type="checkbox"/> Offset (Meso) \$195	<input type="checkbox"/> Inset (Endo) \$195
<input type="checkbox"/> 2" (1" per side)	<input type="checkbox"/> 4" (2" per side)
N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .	
N/A w/ any reclining systems.	

Actual Back Cane Height is reflected in measurement

Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

### BACK OPTIONS

- See Synergy Shape Back Order Form<sup>2,3</sup>
- Solid Plate Back up to 20" W \$275
- Omit Back
- Back Height:<sup>4</sup>
  - 18"H     19"H     20"H     21"H
  - 22"H     23"H     24"H     25"H
- Sling Back<sup>5</sup> \$98

Back height will be the same as back cane height.  
Includes gap for cushion.  
Push Bar selection must match back width; center tube standard when selecting flip-up arms.

- TRU-Comfort Back<sup>4,6</sup>
  - Std.     New Plus
  - Back Height:     22"     25"
  - Back up to 20" W **E2620** \$595
  - 22" W **E2621** \$695
- TRU-Comfort Cushion
  - 16 x 16 - 20 x 18 **E2601** \$380
  - 20 x 20 - 20 x 22 **E2601** \$380
  - 22 x 16 - 22 x 24 **E2602** \$480
- TRU-Comfort Plus Cushion
  - 16 x 16 - 20 x 18 **E2603** \$380
  - 20 x 20 - 20 x 22 **E2603** \$380
  - 22 x 16 - 22 x 24 **E2604** \$480
- Omit Back -(\$100)

See TRU-Comfort Back &amp; Cushion

order forms for specific sizes

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

### ARMRESTS

- Flip-back, cane mounted<sup>3,7</sup> Std.
  - 2-Post, flip-back, height<sup>3,7,8</sup> **E0973** \$150ea adjustable, removable armrests (9"-13" of height adjustment in 1" increments)
  - Single Post height adj.<sup>9,10</sup> **E0973** \$75ea quick height, removable
    - 11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.
    - Quad Release Lever NC (11"-17" of height adj in 3/8" increments)
  - Removable, height adjustable **E0973** \$75ea heavy duty armrests<sup>9</sup>
    - 8"H-12"H     10"H-14"H
  - Pediatric removable<sup>9</sup> **E0973** \$75ea adjustable, quick height.
    - Match Shroud Color     Black
  - Quad Release Lever NC (7 3/4" - 9 3/4" height adjustment in 1/2" increments from seat pan to top of armrest)
  - Cantilever Armrests<sup>14</sup> SOF
- Angle-Adjustable: - 30° through + 30°

### ARMPADS

- Full length, Straight, 2" X 14" Std.
  - Left     Right
- Desk Length, Straight, 2" X 10" Std.
  - Left     Right
- Full Length, Waterfall, 2" X 14"<sup>11</sup> NC
  - Left     Right
- Desk Length, Waterfall, 2" X 10"<sup>11,12</sup> NC
  - Left     Right

### HEAD REST<sup>5</sup>

- Stealth **E0955**
  - Large (11"W x 6" H)     Small (9 3/4" W x 3" H) \$205
- Stealth Comfort Plus **E0955**
  - Large (14"W)     Small (10" W) \$205
- Removable Head Rest Bracket **E1028** \$210
- Stealth, swing-away (Flip-Down) Head Rest Bracket **E1028** \$210  
Interfaces w/ any headrest that uses 1/2" square down-tube.
- Medical Applications **E0955**
  - Large (11"W x 6" H)     Small (9 3/4" W x 3" H) \$230
- Removable Head Rest Bracket **E1028** \$235
- Fixed, Offset Mounting Bracket NC
- TRU-Comfort **E0955**
  - 8" Wide     10" Wide \$205
- Cobra Removable Head Rest Bracket **E1028** \$210
- Removable Bracket for **E1028** \$210  
after market headrest

Footnotes for page 3:

1. Height is measured from seat pan to top of cane.
2. 14" & 22" widths are N/A.
3. Shape Back & w/ flip-back arms may interfere in the same configuration call for details.
4. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
5. Headrests are N/A w/ sling backs.
6. N/A w/ flip back, cane mounted arms.
7. S/A laterals interfere w/ flip-back arms.
8. May cause interference with certain dimension configurations. Call Customer Service.
9. Desk length armrests w/ an accessory limit the accessibility to the height adjustment lever.
10. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Available adj. range is 11" - 17".
11. Only avail. on removable, adjustable, quick-height arms & pediatric arms.
12. S/A joystick mount N/A.
13. N/A w/ TB power recline, Endo/Meso canes, Synergy Manual Systems, vent trays, bariatric armpads, occupied transit kit and TRU-Comfort or sling backs. May also interfere w/ lateral supports and waterfall armpads.
14. Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting after market arm troughs or Universal Swivel w/ desk length arms.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.