



**Outlander
Exterior Lift Series**

Order Form / Price List



Pride Mobility Products Corporation

182 Susquehanna Ave., Exeter, PA 18643

Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

Account Number: _____

Provider Name: _____

Contact: _____

Phone: _____

Fax: _____

Date: _____

PO#: _____

Marked For: _____

Ship To: _____

Address: _____

City: _____ State: _____ Zip: _____

Terms: _____

Freight: _____

Taken By: _____

Vehicle Make: _____ Model: _____ Year: _____

Please read all footnotes - May change vehicle application

CHOOSE LIFT, OPTIONS & ACCESSORIES

MODEL NAME	MSRP	WEIGHT(A)	WEIGHT(B)				HITCH TUBE	POWER OPTIONS		
↓	↓	↓	→	8lbs.	8.5lbs.	10lbs.	9lbs.	N/A	16.5lbs.	16.5lbs.
			CLASS II MSRP = \$0	CLASS III/IV MSRP = \$0	CLASS III/IV EXTENDER MSRP = \$213	UNIVERSAL RISER MSRP = \$89	VEHICLE HARD WIRE KIT ⁶ MSRP = \$0	COMBINATION HARD. WIRE & BATTERY PACK w/ CHARGER ⁴ MSRP = \$330	BATTERY PACK w/ CHARGER MSRP = \$320	
<input type="checkbox"/> OUTLANDER	\$2,050	92.5lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER DE	\$2,050	87.5lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER XL	\$2,050	102.5lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER LP 10"	\$2,050	73lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> OUTLANDER LP 14"	\$2,050	74lbs.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

SECUREMENT ACCESSORIES - (CHOOSE ONE)	MSRP	WEIGHT(C)	ADDITIONAL SECUREMENT ACCESSORIES	MSRP	WEIGHT(C)
<input type="checkbox"/> Scooter Lockdown Arm (LFTASMB1245)	N/C	8lbs.	<input type="checkbox"/> Scooter Lockdown Arm	\$250	8lbs.
<input type="checkbox"/> 3 Ratcheting Straps	N/C	14lbs.	<input type="checkbox"/> 3 Ratcheting Straps	\$250	14lbs.
<input type="checkbox"/> Powerchair Hold Down Arm (LFTASMB1273) ⁶	N/C	6lbs.	<input type="checkbox"/> Powerchair Hold Down Arm (LFTASMB1273) ⁶	\$250	6lbs.

Please read all footnotes - May change vehicle application

ACCESSORIES -	MSRP	WEIGHT(D)
<input type="checkbox"/> Low Profile Hitch Tube Class III (HTCASMB1022)	\$99.45	8.5lbs.
<input type="checkbox"/> Hitch Lowering Tube	\$90	10.5lbs.
<input type="checkbox"/> Swing Away Adaptor ²	\$600	33lbs.
<input type="checkbox"/> Black, Micro Scooter Weather Cover 61" L x 41.25" H	\$114	5lbs.
<input type="checkbox"/> Black, Medium Scooter Weather Cover 75.75" L x 46.5" H	\$158	6lbs.
<input type="checkbox"/> Black, Large Scooter Weather Cover 84" L x 49" H	\$202	7lbs.
<input type="checkbox"/> Black, Micro Power Chair Weather Cover 62" L x 46.25" H	\$114	5lbs.
<input type="checkbox"/> Black, Medium Power Chair Weather Cover 67" L x 51" H	\$158	6lbs.
<input type="checkbox"/> Black, Large Power Chair Weather Cover 76" L x 52" H	\$202	7lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube (Available as SOF) ¹	\$130	5lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 2" Rise (Available as SOF) ¹	\$130	5.5lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 4" Rise (Available as SOF) ¹	\$130	6lbs.
<input type="checkbox"/> Self-Leveling Pivot Tube with 6" Rise (Available as SOF) ¹	\$130	7lbs.

Footnotes:

1. 10 business day lead time applies.
2. For use when rear door or hatch needs to be opened, swings the exterior lift to one side to allow access w/o removing the lift.
3. Tongue Weight - "weight carrying rating only" - Do not use "weight distribution rating."
4. Must not exceed 380lbs. which is the max weight limit of the lift or the vehicle tongue weight rating (see footnote 3 above).
5. Only option available for the Outlander LP models.

CONFIRM VEHICLE APPLICATION

° Weights A _____ + B _____ + C _____ + D _____ = E _____

° E _____ + Mobility Device (including all accessory weights) _____ = F _____

° F _____ - Vehicle Tongue Weight⁴ _____ = G⁵ _____

(G must be less than Vehicle Tongue Weight⁴ for application to work properly)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.



OUTLANDER_LIFT_ORD_FORM-07-31-13