

Quantum® Q6000Z Series

Group 3 Standard
Order Form
300 lbs. weight capacity



ACCOUNT INFO

Account Number: _____
Provider Name: _____
Contact Name: _____
Phone #: _____
Fax #: _____
E-mail Address: _____

Date: _____
PO#: _____
Marked For: _____
3 Column Pricing on Quote/Order: Yes No
Quote #: _____

Ship To: _____
Address: _____
City: _____ State: _____ Zip: _____
Lift Gate Needed: Yes No

Taken By: _____

SET-TO-SPEC*

Please Supply Client Information:

Client Height: _____
Client Weight: _____

Mount Hip/Thigh Supports: Yes No
Mount Lateral Thoracic Supports: Yes No

Armrest Height from Pan: _____

Rigging Length (whole numbers only): _____

Seat to Floor Height (check one):

NOTES:

Seating System	Low	Middle	Tall
Contoured Seating	21"	22"	23"
Synergy® Seat (Standard Towers)	16.625"	17.625"	18.625"
Synergy® Seat (Extended Towers)	18.625"	19.625"	20.625"

Seat to floor height are +/- 1/8" with 18"W x 18"D with Flat-Free tires (add 3/8" for pneumatic tires). Increased widths may increase lowest seat to floor heights.

* While every effort will be taken to ensure dimensions will be preset to your specifications prior to shipping variations in options may effect final settings.

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP. Options with a circular checkbox (⊙) are NOT available through Quantum West.

BASE MODEL		
Description	Price	HCPCS Code
<input type="checkbox"/> Quantum Q6000Z 3S-C	\$6,395	K0849
<input type="checkbox"/> Quantum Q6000Z 3S-SS	\$6,395	K0848

COLORS	
Metallic Colors:	
<input type="checkbox"/> Metallic Red	<input type="checkbox"/> Metallic Black
<input type="checkbox"/> Metallic Blue	<input type="checkbox"/> Metallic Silver
<input type="checkbox"/> Sunburst Orange	

DRIVE WHEELS	
Description	Price
<input type="checkbox"/> 14" Flat-Free (pair)	Standard
<input type="checkbox"/> 14" Pneumatic (pair)	No Charge
<input type="checkbox"/> 4" Wide Flat-Free, Knobby, Extended-Traction Tires ¹	\$551
<input type="checkbox"/> 4" Wide Pneumatic, Knobby, Extended-Traction Tires ¹	\$551

1. Increases overall base width by 1.5" with flat-free, 1.8" with pneumatic.

MOTOR OPTIONS		
Description	Price	HCPCS Code
<input type="checkbox"/> Standard Motor Package	Standard	
<input type="checkbox"/> High Speed Motor Package	\$1,574	
<i>Programmable Speed Up To 8 MPH, 300 lbs. Weight Capacity. HS Motors decrease per-charge range.</i>		
<input type="checkbox"/> Accu-Trac Motor Technology	\$1,295	K0108
<i>Delivers Enhanced Traction and Tracking, Especially With Switched Drive Controls. Requires Q-Logic 2 EX.</i>		
<input type="checkbox"/> Community Use Package ¹	\$2,125	
<u>Package Includes:</u>		
- 300 lbs. Weight Capacity		
- 4-Pole High Speed Motor Package		
- 4" Wide Flat Free, Knobby, Extended Traction Tires		

1. Increases overall base width by 1.5".

BASE MOUNTED FOOT PLATFORM	
Description	Price
<input type="checkbox"/> Base Mounted Foot Platform Installed	Standard
<input type="checkbox"/> No Base Mounted Foot Platform	No Charge

BATTERY CHARGER	
Description	Price
<input type="checkbox"/> 8 Amp Off-board Charger	Standard

TRANSIT KIT	
Description	Price
<i>Refer to last page for Details Regarding Occupied and Unoccupied Transit Kits</i>	
<input type="checkbox"/> Occupied Transit Kit	\$350
<input type="checkbox"/> Unoccupied Transit Kit ¹	No Charge

1. Unoccupied Transit Kit with a Vent Tray may not allow lowest seat to floor height, vent tray unavailable w/ Occupied Transit Kit.

SUSPENSION OPTIONS	
Description	Price
<input type="checkbox"/> Standard Suspension	Standard
<input type="checkbox"/> Pediatric Strut Suspension	No Charge
<i>(For User Weights Less Than 120lbs)</i>	

LED LIGHT PACKAGE	
Description	Price
<input type="checkbox"/> Through Q-Logic 2 EX or NE+ 8-Key w/ Lights	\$995

NON-EXPANDABLE CONTROLS	
Description	Price
<input type="checkbox"/> 75 Amp NE	Standard
<input type="checkbox"/> 75 Amp NE+ 6-Key For Up To 2 Actuators Thru Joystick	\$200 <i>Standard on 3SP-C Package</i>
<input type="checkbox"/> 75 Amp NE+ 8-Key For Lights	\$400

EXPANDABLE CONTROLS	
Description	Price
<i>Q-Logic 2 EX Joystick or Enhanced Display must be ordered with a Q-Logic 2 EX Controller and Harness For Expandable Controls</i>	
<input type="checkbox"/> 100 Amp Q-Logic 2 EX Controller	\$695
<input type="checkbox"/> Expandable Controller Harness	\$350
<input type="checkbox"/> Q-Logic 2 EX Joystick	\$700
<input type="checkbox"/> Q-Logic 2 EX Joystick for Lights	\$700
<input type="checkbox"/> Stand Alone Joystick w/ External Charger Port	\$700

ARMOR COLOR OPTIONS:	
<input type="checkbox"/> Black Soft Touch (Standard color option)	
Colors: \$50	
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Sunburst Orange
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Pearl Pink

PROGRAMMING	
Description	Price
<input type="checkbox"/> Q-Logic Handheld Programmer	\$545
<input type="checkbox"/> Q-Logic PC Programming Cable with Disk ¹	\$750

1. Works with Q-Logic 2 EX only.

Options with a circular checkbox (○) are NOT available through Quantum West.

JOYSTICK HANDLES

Description	Price	HCPCS Code
<input type="checkbox"/> Standard Carrot Handle	Standard	
<input type="checkbox"/> Mushroom Handle	\$71	E2323
<input type="checkbox"/> T-Bar	\$107	E2323
<input type="checkbox"/> Chin Cup	\$71	E2324 ¹
<input type="checkbox"/> Extended I-Handle	\$71	E2323
<input type="checkbox"/> Soft Ball Handle	\$50	E2323
<input type="checkbox"/> Medium Ball Handle	\$50	E2323
<input type="checkbox"/> Bodypoint 3" Goal Post Handle	\$69	E2323
<input type="checkbox"/> Bodypoint 4" Goal Post Handle	\$69	E2323

1. Use HCPCS code E2323 when not used in conjunction with a chin control.

JOYSTICK MOUNTING BRACKETS

Description	Price	HCPCS Code
<i>Comes Standard with a Quick Release, Adjustable-Height Joystick Block</i>		
<input type="checkbox"/> Fixed Inline ¹	Standard	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Swing-away	\$210	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Retract4 Swing-Away Joystick Arm Kit	\$250	E1028
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Multi-axis	\$100	
Side: <input type="checkbox"/> Left <input type="checkbox"/> Right		

1. Not available with Q-Logic 2 EX Stand Alone Joystick

CONTOURED SEATING

Description	Price	HCPCS Code
Depth Adjustable High Back:		
<input type="checkbox"/> Limited Recline Depth Adjustable High Back w/ Headrest ¹	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 18"W x 18-20"D <input type="checkbox"/> 20"W x 18-20"D <input type="checkbox"/> 22"W x 20-22"D		

High Back with Solid Seat Pan:

<input type="checkbox"/> Limited Recline High-Back with Headrest & Solid Seat Pan	No Charge	
Color: <input type="checkbox"/> Black Vinyl		
Size: <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 16"W x 20"D <input type="checkbox"/> 16"W x 22"D <input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 18"W x 18"D <input type="checkbox"/> 18"W x 20"D <input type="checkbox"/> 18"W x 22"D <input type="checkbox"/> 20"W x 16"D <input type="checkbox"/> 20"W x 18"D <input type="checkbox"/> 20"W x 20"D <input type="checkbox"/> 20"W x 22"D		

Armrests:

<input type="checkbox"/> Flip-Up, Removable	No Charge	
<input type="checkbox"/> Flip-Up, Height-Adjustable, Removable	\$85 ea.	E0973
Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length ²		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length ²		

Special Order Seat

_____ "W x _____ "D

Seat size can not exceed 24"W x 24"D

Contact your inside sales specialist for pricing

1. Seat comes automatically with a 60" belt and full length armrests. Any other selection would need to be added as an additional part.
2. Selection of desk length armpads on a depth-adjustable high back requires them to be added as an additional part.

MANUAL RECLINER

Description	Price	HCPCS Code
<input type="checkbox"/> Manual Recliner High Back with Headrest	\$600	E1226
<input type="checkbox"/> Manual Recliner High Back with Headrest and Solid Seat Pan ¹	\$600	E1226
Upholstery:		
<input type="checkbox"/> Gray Vinyl		
Seat Width:		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
Seat Depth:		
<input type="checkbox"/> 16" <input type="checkbox"/> 18" <input type="checkbox"/> 20"		
Armpads:		
Left: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		
Right: <input type="checkbox"/> Full Length <input type="checkbox"/> Desk Length		

1. Includes black vinyl seatpan. Back available in listed colors.

SYNERGY® SEATING SYSTEMS

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy Seating Frame	\$434	
<input type="checkbox"/> Standard Stem Mounts		
<input type="checkbox"/> 2" Extended Vertical Stem Mounts ¹		
<input type="checkbox"/> Synergy Manual Recline Package ²	\$700	E1226
<input type="checkbox"/> Synergy Manual Tilt Package ³	\$550	K0108
<i>Max. Seat Size 14"W x 14"D, 75 lbs Weight Capacity</i>		

Seat Width:

<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19"
<input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24"

Seat Width Package:

<input type="checkbox"/> 10"W-14"W <input type="checkbox"/> 14"W-18"W <input type="checkbox"/> 16"W-20"W <input type="checkbox"/> 20"W-24"W

Seat Width Package (Power Elevate):

<input type="checkbox"/> 14"W-18"W <input type="checkbox"/> 18"W-22"W

Seat Depth:

<input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13" <input type="checkbox"/> 14"
<input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/> 19"
<input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23" <input type="checkbox"/> 24"

Seat Depth Package:

<input type="checkbox"/> 10"D-14"D <input type="checkbox"/> 14"D-18"D <input type="checkbox"/> 15"D-20"D <input type="checkbox"/> 20"D-24"D

Back Cane Height:

<input type="checkbox"/> 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" ⁴ <input type="checkbox"/> 17" ⁴ <input type="checkbox"/> 18"
<input type="checkbox"/> 19" <input type="checkbox"/> 20" <input type="checkbox"/> 21" <input type="checkbox"/> 22" <input type="checkbox"/> 23"
<input type="checkbox"/> 24" <input type="checkbox"/> 25"

1. Increases seat to floor height by 2". Lowest seat-to-floor height not available.
2. Not available with Flip Back Arms, Cantilever Arms, Shape Back or Power Elevate. Requires 22"H canes. Pivot point 2.5" above seat pan.
3. In i-Quantum this option is chosen under "Width Package Type"
4. Selection of a Synergy Back and cane mounted, flip up armrests disqualifies 16" and 17"H back canes.

Options with a circular checkbox (○) are NOT available through Quantum West.

SYNERGY® BACK OPTIONS

Description	Price	HCPCS Code
<input type="checkbox"/> Synergy® Shape Back Positioning Back <i>See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes</i>	Various	E2613/E2614
<input type="radio"/> Sling Back	\$98	
<input type="checkbox"/> Omit Back	No Charge	
<input type="checkbox"/> Synergy® Back with 1.5" Thick Back Cushion	\$275	
Overall Back <input type="checkbox"/> 18" H <input type="checkbox"/> 19" H <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H		
Height: <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Choose Lateral Extensions ¹ :		
<input type="checkbox"/> 6" for 13"-16" H Back	\$25	
<input type="checkbox"/> 10" for 17"-18" H Back	\$25	
<input type="checkbox"/> 12" For 19"-22" H Back	\$25	
<input type="checkbox"/> 16" For 23"-25" H Back	\$25	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width: <input type="checkbox"/> 16" W <input type="checkbox"/> 18" W <input type="checkbox"/> 20" W	\$595	E2620
Overall Back <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H		
Height ² : <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 17" W <input type="checkbox"/> 19" W	\$595	E2620
Overall Back <input type="checkbox"/> 22" H <input type="checkbox"/> 23" H <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H		
Height ² : <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Width: <input type="checkbox"/> 22" W <input type="checkbox"/> 24" W ³	\$695	E2621
Overall Back <input type="checkbox"/> 24" H <input type="checkbox"/> 25" H <input type="checkbox"/> 27" H <input type="checkbox"/> 28" H		
Height ² :		

Tru Comfort Color Shroud Option:

Black (Standard color option)

Metallic Colors:\$100

- Red
- Blue
- Sunburst Orange
- Silver

Synergy® Back Cane Options:

- Standard Back Canes Standard
- Omit Push Handles \$125
- Flip Forward Back Canes⁴ \$150
- Omit Push Handles \$125
- Offset (Meso)⁵ \$195
- 2" (1" Per Side) 4" (2" Per Side)
- Inset (Endo) \$195
- 2" (1" Per Side) 4" (2" Per Side)

1. Lateral extensions are required to mount lateral thoracic supports.
 2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".
 3. Available with offset (meso) brackets.
 4. Not available with Vent Tray, O2 Holder, cane mounted or cantilever armrests. Also not available w/ a Synergy manual recline or manual tilt. May interfere with Lateral Supports and Waterfall Armrests.
 5. Dual Post Flip-back Armrests interfere with offset canes.

SYNERGY® ARMRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> Flip-back, Cane Mounted ¹	No Charge	
<input type="checkbox"/> 2-post, Flip-back, Height-adjustable, Removable ^{2,3}	\$150 ea.	E0973
<input type="checkbox"/> Single-post, Height-adjustable, Removable ^{3,4}	\$75 ea.	E0973
<input type="checkbox"/> Quad Release Lever		
<input type="checkbox"/> Removable, Single-post, Height Adjustable, Heavy Duty	\$75 ea.	E0973
<input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/> Pediatric, Removable ⁵	\$75 ea.	E0973
<input type="checkbox"/> Quad Release Lever		
<input type="radio"/> Bariatric, Removable, Height Adjustable ⁶	\$75 ea.	E0973
<input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H		
<input type="checkbox"/> Height Adjustable Cantilever Armrests ⁷	\$195	
Angle-Adjustable: -30° through +30°		

1. Interferes with swing-away laterals. Not available on flip forward canes or manual recline.
 2. Interferes with swing-away laterals. Height adjustable from 9"-13" in 1" increments.
 3. Not available with a depth less than 16".
 4. Height adjustable from 11"-17" in 1/2" increments
 5. Height adjustable from 7.75"-9.75" in 1/2" increments. Full length armrests not available
 6. Bariatric Armrests require Bariatric Armrests
 7. Height adjustment from 11"-17" depends on cane height. Not available with Manual Recline, Flip Forward back canes & Shape Back.

SYNERGY® ARMPADS

Description	Price	HCPCS Code
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Full Length (2"x14") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Waterfall, Desk Length (2"x10") ¹	No Charge	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$92.50 ea.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$56.69 ea.	
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$70.00 ea.	
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$53.52 ea.	
<input type="radio"/> Flat, Bariatric ²	No Charge	
<input type="radio"/> Left <input type="radio"/> Right		

1. Not available w/ Cane Mounted. Flip Back Armrests combined with a Synergy Back.
 2. Bariatric Armrests require Bariatric Armpads

HEADRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> No Headrest Required		
TRU-Comfort™ Headrest¹:		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
Stealth:		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket ² (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket (FDM350)	\$210	
Used in Addition to Removable Headrest Bracket		
Aftermarket Bracket:		
<input type="checkbox"/> Stealth's World's Best Removable Headrest Hardware (TWB485)	\$210	E1028

1. Measurement is taken from the inside of the headrest.

NON-POWERED LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 60-90° Swing-away Legrests ¹	Standard	
(Adjustable in 5° increments)		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Standard Extruded Footplates:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates ²	\$150	K0040
Angle Adjustable Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Angle-Adjustable Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		
<input type="checkbox"/> Multi-Axis Footplates ²	\$495	K0040
Multi-Axis Extensions:		
<input type="checkbox"/> 15.5"-19.5"		
Multi-Axis Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴		

1. Possible caster wheel interference when adjusted to less than 70 degrees.
 2. Includes Heel Loops
 3. Causes interference on seat widths less than 18".
 4. Causes interference on seat widths less than 20".

Options with a circular checkbox (☉) are NOT available through Quantum West.


NON-POWERED LEGRESTS (CONT'D)

Description	Price	HCPCS Code
<input type="checkbox"/> High Mount, Clamp-On Footplates*		
<input type="checkbox"/> Angle Adjustable		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	
<input type="checkbox"/> Heavy Duty Drop-In Legrests¹	\$420	
<input type="checkbox"/> Angle Adjustable Footplates ²	\$150	K0040
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5"		<input type="checkbox"/> 14.5" -18"
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"
<input type="checkbox"/> 6" x 8" ³		<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 11" ³		<input type="checkbox"/> 8" x 11" ⁴
<input type="checkbox"/> Elevating Legrests²	\$325	E0990
<input type="checkbox"/> Standard Footplates		
<i>Standard Extruded Extension Length:</i>		
<input type="checkbox"/> 11"-14"		<input type="checkbox"/> 14"-17"
<i>Standard Extruded Footplate Size:</i>		
<input type="checkbox"/> 7.25" x 6"		<input type="checkbox"/> 8.5" x 6"
<input type="checkbox"/> Angle Adjustable Footplates ³	\$150	K0040
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 12"-14.5"		<input type="checkbox"/> 14.5" -18"
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"
<input type="checkbox"/> 6" x 8" ⁴		<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 11" ⁴		<input type="checkbox"/> 8" x 11" ⁵
<input type="checkbox"/> 90° Pediatric Swing-Away Legrests²	\$162	
<input type="checkbox"/> Angle Adjustable Footplates ³		
<i>Angle Adjustable Extension Length:</i>		
<input type="checkbox"/> 7"-10"		
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"

1. Only available with Synergy seating.
2. Extensions are measured from the seatpan to footplate.
3. Includes heel loops.
4. Causes interference on seat widths less than 18".
5. Causes interference on seat widths less than 20".

ADJUSTABLE CENTER MOUNT FOOT PLATFORM


Description	Price
<input type="checkbox"/> Adjustable Center Mount Foot Platform ¹	\$250
<i>Includes 10"W x 11"D Split Foot Platform</i>	
<i>Extensions:²</i>	
<input type="checkbox"/> 8"-12" Ext.	<input type="checkbox"/> 10"-14" Ext
<input type="checkbox"/> 10"-16" Ext.	<input type="checkbox"/> 13"-19" Ext ³
<i>Calf Supports⁴:</i>	
<input type="checkbox"/> Calf Pads (pair)	
<input type="checkbox"/> Small (6"H x 6"W)	\$150
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150
<input type="checkbox"/> Calf Panel (ea.)	
<input type="checkbox"/> Small (8"H x 18"W)	\$150
<input type="checkbox"/> Large (8"H x 24"W) ⁵	\$200
<i>Optional Modular Footplate Extension Components/Kits^{4,6}</i>	
<input type="checkbox"/> 4" Width Extension Kit	\$50
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50



1. Only available with Synergy seating.
2. Extension length is measured from top of the seat pan to footplate.
3. Not available on Recline only.
4. Choose in Rigging Accessory Section in iQuantum.
5. Only available on 20"W or greater.
6. Assembly required. May require seating configuration adjustments to avoid front caster interference.

ADJUSTABLE CENTER MOUNT FOOT PLATFORM (CONT'D)

Description	Price
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75
<input type="checkbox"/> 6-Piece Extension Stock Kit <i>(Includes all above pieces)</i>	\$95
<i>Additional Optional Accessories⁷:</i>	
<input type="checkbox"/> Footplate Connector	\$49
<input type="checkbox"/> Wheel Bumpers, pair	\$25



7. Assembly required. May require seating configuration adjustments to avoid front caster interference.

ELECTRONICS FOR POWER FOOT RIGGINGS

Description	Price	HCPCS Code
<input type="checkbox"/> Single Actuator Through Joystick <i>(AFP Only)</i>	\$1,175	E2310 ¹
<input type="checkbox"/> Dual Actuator Through NE+ 6-Key <i>(Power Legrests)</i>	\$2,375	E2311 ²
<input type="checkbox"/> Multiple Actuator Through Q-Logic 2 EX <i>(Power Legrests)</i>	\$2,375	E2311 ²

1. HCPCS used only with Power AFP.
2. HCPCS used only with Power Legrests


POWER ELEVATING LEGRESTS

Description	Price	HCPCS Code
<input type="checkbox"/> 70° Articulating Swing-Away Legrests, 4" Articulation Range¹	\$2,395	E1010
<input type="checkbox"/> Quad Release Lever		
<i>Calf Pads²:</i>		
<input type="checkbox"/> 6" x 5.35"		<input type="checkbox"/> 8" x 5.86"
<i>Angle Adjustable Extension Length³:</i>		
<input type="checkbox"/> 11"-13"		<input type="checkbox"/> 14.5"-18"
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"
<input type="checkbox"/> 6" x 8" ⁴		<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 11" ⁴		<input type="checkbox"/> 8" x 11" ⁵
<i>Operational Method:</i>		
<input type="checkbox"/> Independent		<input type="checkbox"/> Combined
<input type="checkbox"/> 70° Articulating Drop-In Legrests¹	\$2,395	E1010
<input type="checkbox"/> 4" Range		
<input type="checkbox"/> 6" Range		
<i>Calf Pads²:</i>		
<input type="checkbox"/> 6" x 5.35"		<input type="checkbox"/> 8" x 5.86"
<i>Angle Adjustable Extension Length³:</i>		
<input type="checkbox"/> 12"-14.5"		<input type="checkbox"/> 14.5"-18"
<i>Angle Adjustable Footplate Size:</i>		
<input type="checkbox"/> 4" x 6"		<input type="checkbox"/> 5" x 6"
<input type="checkbox"/> 6" x 8" ⁴		<input type="checkbox"/> 5" x 8"
<input type="checkbox"/> 6" x 11" ⁴		<input type="checkbox"/> 8" x 11" ⁵
<i>Operational Method:</i>		
<input type="checkbox"/> Independent		<input type="checkbox"/> Combined

1. Only available with Synergy seating.
2. Choose in Rigging Accessory Section in iQuantum.
3. Extensions are measured from top of seat pan to footplate.
4. Only available with seat width 18" or greater.
5. Only available with seat width 20" or greater.

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Power Articulating Foot Platform ¹	\$2895	K0108
<i>Includes 10"W x 11"D Split Foot Platform</i>		
<i>Extensions²:</i>		
<input type="checkbox"/> 10"-16" Ext.		<input type="checkbox"/> 13"-19" Ext.
<i>Calf Supports³:</i>		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) ⁴	\$200	
<i>Optional Modular Footplate Extension Components/Kits^{3,5}</i>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	





1. Only available with Synergy seating.
2. Extension length is measured from top of the seat pan to footplate.
3. Choose in Rigging Accessory Section in iQuantum.
4. Only available on 20"W or greater.
5. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

Options with a circular checkbox (⊙) are NOT available through Quantum West.

CENTER MOUNT POWER ARTICULATING FOOT PLATFORM (CONT'D)

Description	Price	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/> 6-Piece Extension Stock Kit <i>(Includes all above pieces)</i>	\$95	
<i>Additional Optional Accessories¹:</i>		
<input type="checkbox"/> Footplate Connector	\$49	
<input type="checkbox"/> Wheel Bumpers, pair	\$25	

1. Assembly required. May require seating configuration adjustments to avoid front caster interference.

FRONT RIGGING ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) ¹	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) ²	\$65	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$512	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair) *	\$574	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Bodypoint Ankle Huggers (pair)*	\$134	K0108
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Small <input type="checkbox"/> Medium		

1. Standard with Angle Adjustable or Multi-Axis footplates.
2. Only available with Swing-Away legrests.

SEAT CUSHIONS

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
Simplicity Cushion <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 24"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
Solution® 1 Cushion <i>Skin Protection:</i>		
<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604
Solution® Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608
Spectrum® Gel Cushion <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2606
Spectrum® Visco Cushion <i>Skin Protection & Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2608
Spectrum® Air <i>Adjustable Skin Protection</i>		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$480	E2622
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$530	E2623
Spectrum® Air Contour <i>Adjustable Skin Protection & Positioning</i>		
<input type="checkbox"/> 14" x 12" - 20" x 24"	\$499	E2624
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$549	E2625
TRU-Comfort™ 2 4-Way Stretch <i>Skin Protection</i>		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604
TRU-Comfort™ 2 Air Mesh <i>Skin Protection</i>		
<input type="checkbox"/> 16" x 14" - 20" x 20"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 22" x 24"	\$480	E2604

LATERALS

Description	Price	HCPCS Code
<input type="checkbox"/> Stealth Thoracic Laterals (pair)*		
<i>Bracket:</i>		
<input type="checkbox"/> Swing-Away Rail Mounts	\$450	E1028
<i>For Use on Synergy Back and Tru-Comfort</i>		
<i>Pad:</i>		
<input type="checkbox"/> Flat Pad, 3"x5", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Flat Pad, 4"x6", Leather-Like Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 2"x3", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 3"x5", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Curved Pad, 4"x6", 4-Way Stretch Cover	\$200	E0956
<input type="checkbox"/> Stealth Pelvic/Thigh Laterals*		
<i>Side:</i>		
<input type="checkbox"/> Left Side <input type="checkbox"/> Right Side		
<i>Bracket:</i>		
<input type="checkbox"/> Removable Bracket (ea.)	\$250	E1028
<input type="checkbox"/> Removable Bracket w/ 1" Offset (ea.)	\$325	E1028
<input type="checkbox"/> Flip Down Bracket (ea.)	\$225	E1028
<input type="checkbox"/> Flip Down, Full Surface Contact Bracket (ea.)	\$225	E1028
<i>Pad:</i>		
<input type="checkbox"/> Short Pad, 4"x4" (ea.)	\$50	E0956
<input type="checkbox"/> Medium Pad, 4"x8" (ea.)	\$50	E0956
<input type="checkbox"/> Long Pad, 4"x12" (ea.)	\$50	E0956
<i>Adapter:</i>		
<input type="checkbox"/> 2" Depth/Height Adapter Plate (ea.) ²	\$25	K0108
<input type="checkbox"/> 2 1/4" Depth/Height Adapter Plate (ea.) ³	\$25	K0108

1. Swing-away laterals interfere with flip-up armrests.
2. Height adapter plate is needed w/ Tru Comfort seat cushion. Not available on the flip down bracket.
3. Only available on the flip down bracket.

MISCELLANEOUS POSITIONING COMPONENTS

Description	Price	HCPCS Code
<input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket*	\$172	E1028
<i>Pad:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large	\$135	E0957
<input type="checkbox"/> Stealth Foot Boxes (pair)*	\$475	K0108
<input type="checkbox"/> Stealth Padded Shoe Holders (pair)*	\$166	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large		
<input type="checkbox"/> Otto Bock Arm Channels (ea.)*	\$104	E2209
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large		
<input type="checkbox"/> Otto Bock Flat Hand Pad(ea.)*	\$52	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<i>Size:</i>		
<input type="checkbox"/> Medium <input type="checkbox"/> Large		

Please see positioning components order form for additional types, sizes and details.

ACCESSORIES

Description	Price	HCPCS Code
<input type="checkbox"/> Rear Basket ¹	\$103	
<input type="checkbox"/> Cane/Crutch Holder ¹	\$98	E2207
<input type="checkbox"/> Weather Cover, Medium, Black	\$158	
<input type="checkbox"/> Cup Holder	\$17 ²	
<i>Side:</i>		
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Saddle Bag	\$28	
<input type="checkbox"/> Flag Holder	\$17	
<input type="checkbox"/> Quantum Backpack ³	\$60	
<input type="checkbox"/> Oxygen Holder, Contoured Seating ¹	\$103	E2208
<input type="checkbox"/> Oxygen Holder, Synergy Seat ⁴	\$150	E2208
<input type="checkbox"/> Walker Holder ⁴	\$92	
<input type="checkbox"/> Tablet Holder ⁵	\$315	
<input type="checkbox"/> Deluxe Cell Phone/Media Holder ⁵	\$140	

1. Includes mounting hardware. Not available on Manual Recliner.
2. No charge on Synergy seating.
3. Backpack may interfere w/ other rear accessories. Not available w/ contoured seating.
4. Not available with Synergy manual recline or manual tilt. Includes mounting hardware.
5. Only available with Synergy or Specialty seating.

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum

Options with a circular checkbox (☐) are NOT available through Quantum West.

LAPBELTS

Description	Price
<input type="checkbox"/> 38" Lap Belt ¹	Standard
<input type="checkbox"/> 50" Lap Belt ¹	Standard
<input type="checkbox"/> 60" Lap Belt	Standard
<input type="radio"/> 60" Non-Retractable (Sure-Lok) ²	Standard
<input type="checkbox"/> 60" Padded Lap Belt ³	\$115
<input type="checkbox"/> 70" Lap Belt	Standard
<input type="checkbox"/> 80" Lap Belt	Standard
<input type="checkbox"/> 90" Lap Belt	Standard
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width*	\$87
<input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width ³	\$92
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width*	\$112
<input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width*	\$116

1. Includes mounting hardware. Not available with Manual Recliner.
2. Only needed with the occupied transit kit.
3. Only available with Synergy seating.

BATTERIES

Description	Price	HCPCS Code
<input type="radio"/> Group 24, 70 Amp, Gel (Universal)	\$352 ea.	E2363
<input type="checkbox"/> Group 24, 70 Amp, Gel (MK)	\$352 ea.	E2363
<input type="checkbox"/> NF-22, 50-Amp Batteries		
<input type="checkbox"/> Maintenance Free AGM (MK)	\$295 ea.	E2361
<input type="checkbox"/> Maintenance Free Gel Cell (Universal)	\$295 ea.	E2361
<input type="checkbox"/> Maintenance Free Gel Cell (MK)	\$295 ea.	E2361
Installation:		
<input type="checkbox"/> Batteries Installed		
<input type="checkbox"/> Do Not Install		

TRANSIT KIT INFO

Q6000Z

THIS SECTION MUST BE READ

The Transit Securement Package, when used in conjunction with a wheelchair in the specific configuration of an 18"x 20" seat size with 25" high back canes, Full-length height-adjustable single post armrests, Synergy back 25"h, and style 8 swing away legrests and extruded footplates and heel-loops, conforms with the standard of ANSI/RESNA WC/Vol. 4, Section 19, for the purpose of enabling wheelchairs to be secured in certain types of motor vehicles. The Transit Securement Package, including manufacturer-installed front and rear securement brackets and pelvic belt anchoring brackets, has been crash tested in accordance with ANSI/RESNA WC/Vol. 4, Section 19, Frontal Impact Test, with a 168-lbs. (76.2 kg) surrogate occupant corresponding to a user weight range of 115 lbs. (52.16 kg) to 210 lbs. (95.25 kg). While the Transit Securement Package is compatible with additional seat sizes, optional components, and varied configurations, and provides designated securement points, Pride Mobility makes no claims that any other components or configurations have been tested beyond the above described ANSI/RESNA conforming configuration. Please refer to the Transit Securement Package supplement (part number INFMANU3658) for additional information.

RULES:

- All options are available when unoccupied.

Notes when selecting occupied transit kit:

- It is recommended that usable back height be 2" lower than shoulder height.
- Available on Synergy Seat up to 20"W x 20"D
- Power AFP, ALRs & ELRs are available on Synergy Seating
- All legrests must have heel loops
- Synergy manual seating is N/A
- TRU-Balance Power Positioning seating, options and accessories are N/A
- TRU-Comfort seat cushion is N/A
- Contoured seating is N/A
- Specialty seating is N/A
- Vent option is N/A
- Omit Back option is N/A

Items marked with a (*) should be chosen in the Positioning Components screen in iQuantum