

Quantum® R4000 Series

Group 3 Single Power & Multiple Power

Order Form

300 lbs weight capacity



Account Number: _____

Provider Name: _____

Contact: _____

Phone: _____

Fax: _____

Date: _____

PO#: _____

Marked For: _____

Ship To: _____

Address: _____

City: _____ State: _____ Zip: _____


Discount: _____

Terms: _____

Freight: _____

Taken By: _____

SET-TO-SPEC*



Please supply patient information for optimal setup prior to shipping.

Client Height: _____

Client Weight: _____

A. Seat Width: _____

B. Lateral Offset: _____

C. Hip/Thigh Lateral Offset: _____

D. Center of Gravity: *Optimized based on seat dimensions*

E. Seat Depth: _____

F. Back Height from Pan: _____

G. Lateral Height from Pan: _____

H. Armrest Height from Pan: _____

I. Rigging Length (whole numbers only): _____

J. Seat to Floor Height (circle one): _____

Clinician Notes: _____

| Seating System | Seat-To-Floor Heights | | | |
|---|-----------------------|---------|---------|---------|
| | 18.875" | 19.875" | 20.875" | 21.875" |
| TRU-Balance® Tilt or Tilt & Recline | 18.875" | 19.875" | 20.875" | 21.875" |
| Super-Low TRU-Balance® Tilt or Tilt & Recline | 17.875" | 18.875" | 19.875" | 19.875" |
| TRU-Balance® Recline | 16.75" | 17.625" | 18.5" | 18.5" |
| TRU-Balance® Lift | 19.625" | 20.5" | 21.375" | 21.375" |
| Super-Low TRU-Balance® Lift | 17.875" | 18.75" | 19.625" | 19.625" |

Seat to floor height are +/- 1/8" with 18"W x 18"D with Flat-Free tires and 8" casters (add 5/8" for 9" casters). Increased widths may increase lowest seat to floor heights.

* While every effort will be taken to ensure items will be preset to your specifications prior to shipping, variations in options may effect final settings.

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

HCPCS codes should not be considered as legal advice. All prices are MSRP.

BASE MODEL

| Description | Price | HCPCS Code |
|---|---------|------------|
| <input type="checkbox"/> Quantum R4000 3SP-SS | \$6,808 | K0856 |
| <input type="checkbox"/> Quantum R4000 3MP-SS | \$6,608 | K0861 |

COLOR

- Candy Apple Red
- Viper Blue
- Sunburst Orange
- Onyx Black
- Silver
- Champagne

DRIVE WHEELS

| Description | Price |
|--|-----------|
| <input type="checkbox"/> 14" Flat-Free <i>Rim Color:</i> <input type="checkbox"/> Black <input type="checkbox"/> Silver | Standard |
| <input type="checkbox"/> 14" Pneumatic <i>Rim Color:</i> <input type="checkbox"/> Black <input type="checkbox"/> Silver | No Charge |
| <input type="checkbox"/> 4" Wide Flat Free, Knobby, Extended-Traction Tires ¹ | \$551 |
| <input type="checkbox"/> 4" Wide Pneumatic, Knobby, Extended-Traction Tires ¹ | \$551 |

CASTER WHEELS

| | |
|--|----------|
| <input type="checkbox"/> 9" Pneumatic (pair) | Standard |
| <input type="checkbox"/> 9" Flat-Free (pair) | \$105 |
| <input type="checkbox"/> 8" Flat-Free (pair) | \$105 |

MOTOR OPTIONS

| Description | Price | HCPCS Code |
|---|----------|------------|
| <input type="checkbox"/> Standard Motor Package | Standard | |
| <input type="checkbox"/> High Speed H2 (6MPH) Motor Package <i>Programmable Speed Up To 6MPH, 300 lbs. Weight Capacity. HS Motors decrease per-charge range.</i> | \$1,574 | |
| <input type="checkbox"/> High Speed Hammer (8MPH) Motor Package ² <i>Programmable Speed Up To 8 MPH, 300 lbs. Weight Capacity. HS Motors decrease per-charge range.</i> | \$1,574 | |
| <input type="checkbox"/> Community Use Package ¹ <i>Package Includes:</i> - 300 lbs. Weight Capacity - 4-Pole High Speed (6MPH) Motor Package - 4" Wide Flat Free or Pneumatic Knobby, Extended Traction Tires | \$2,125 | |

LED LIGHT PACKAGE

| Description | Price |
|--|-------|
| <input type="checkbox"/> Through Q-Logic EX with Lights ³ | \$995 |
| <input type="checkbox"/> Through 9-Key VR2 with Lights | \$995 |

BATTERY CHARGER

| Description | Price |
|---|----------|
| <input type="checkbox"/> 8 Amp, Off-board Charger | Standard |

NON-EXPANDABLE CONTROLS

| Description | Price |
|--|--|
| <input type="checkbox"/> PG VR2 70 amp <input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile | No Charge |
| <input type="checkbox"/> PG 6-Key VR2 70 amp <input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile | \$200 <i>Standard on 3SP-SS Package</i> |
| <input type="checkbox"/> PG 9-Key VR2 w/ Lights 70 amp <input type="checkbox"/> 1 Profile <input type="checkbox"/> Multiple Profile | \$200 |

FOOTNOTES

- Increases overall base width by 1.5" with Flat Free, 1.8" with Pneumatic.
- Requires Q-Logic EX.
- LED light package with Tru Balance Lift requires SOF.
- All mounting hardware is included in kit. See Specialty Electronics order form for additional mounting options and secondary controls. Bluetooth mouse emulator requires a run plug harness for either static or power seating. Additional fees will apply.
- Use code E2323 when not used in conjunction with a chin control.
- Not available with Q-logic EX Stand Alone Joystick.

EXPANDABLE CONTROLS

Q-Logic EX Joystick or Enhanced Display must be ordered with a Q-Logic EX Controller and Harness For Expandable Controls

| Description | Price | HCPCS Code |
|---|------------|------------|
| <input type="checkbox"/> 100 Amp Q-Logic EX Controller | \$695 | E2377 |
| <input type="checkbox"/> Harness For Expandable Controller | \$350 | E2313 |
| <input type="checkbox"/> Q-Logic EX Joystick | \$700 | |
| <input type="checkbox"/> Q-Logic EX Joystick for Lights <i>Required For LED Light Package</i> | \$700 | |
| <input type="checkbox"/> Q-Logic EX Standalone Joystick | \$700 | E2373 |
| <input type="checkbox"/> Q-Logic EX Enhanced Display ⁴ <i>w/ Bluetooth Mouse & IR Capability (KITASMB1678)</i> Includes: - Q-Logic EX Enhanced Display - Flexible Gooseneck - Enhanced Display Mounting Bracket | \$1,884.93 | |

PROGRAMMING

| Description | Price |
|---|-------|
| <input type="checkbox"/> Q-Logic Handheld Programmer | \$545 |
| <input type="checkbox"/> Q-Logic PC Programming Cable with Disk | \$750 |
| <input type="checkbox"/> PG Handheld Programmer | \$423 |

JOYSTICK HANDLES

| Description | Price | HCPCS Code |
|--|----------|--------------------|
| <input type="checkbox"/> Standard Knob | Standard | |
| <input type="checkbox"/> Mushroom Handle | \$71 | E2323 |
| <input type="checkbox"/> T-Bar | \$107 | E2323 |
| <input type="checkbox"/> Chin Cup | \$71 | E2324 ⁵ |
| <input type="checkbox"/> Extended I-Handle | \$71 | E2323 |
| <input type="checkbox"/> Soft Ball Handle | \$50 | E2323 |
| <input type="checkbox"/> Medium Ball Handle | \$50 | E2323 |
| <input type="checkbox"/> Bodypoint 3" Goal Post Handle | \$66 | E2323 |
| <input type="checkbox"/> Bodypoint 4" Goal Post Handle | \$66 | E2323 |

JOYSTICK MOUNTING BRACKETS

| Description | Price | HCPCS Code |
|---|----------|------------|
| <input type="checkbox"/> Fixed Inline ⁶ <input type="checkbox"/> Left <input type="checkbox"/> Right | Standard | |
| <input type="checkbox"/> Swing-away <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> Inline <input type="checkbox"/> Offset <input type="checkbox"/> 4" Ext. <input type="checkbox"/> 7" Ext. | \$210 | E1028 |
| <input type="checkbox"/> Multi-axis Fixed <input type="checkbox"/> Left <input type="checkbox"/> Right | \$100 | K0108 |

TRU-BALANCE® POWER POSITIONING SYSTEMS

| Description | Price | HCPCS Code |
|---|-----------|------------|
| <input type="checkbox"/> TRU-Balance® Power Tilt (45°) <i>45 degree center of gravity (CG) power tilt seating system. 300lbs weight capacity.</i> | \$4995 | E1002 |
| <input type="checkbox"/> TRU-Balance® Power Recline (168°) <i>Full Drive lock-out at 20 degrees. 300 lbs. weight capacity</i> | \$5595 | E1004 |
| <input type="checkbox"/> TRU-Balance® Power Tilt (45°) & Recline (168°) <i>Full Drive lock-out at 20 degrees, 300 lbs. weight capacity 45 degree center of gravity (CG) power tilt seating system.</i> | \$9595 | E1007 |
| <input type="checkbox"/> TRU-Balance® Lift (10") ¹ <i>Must be selected in conjunction with TB Power Tilt or TB Power Tilt & Recline.</i> | \$2895 | |
| <input type="checkbox"/> TRU-Balance® Super Low Option ² <i>For Use With TRU-Balance Tilt, Tilt & Recline, or Lift To Achieve Lower Seat-To-Floor Height</i> | No Charge | |

SEATING SPECIFICATIONS

Seat Width:

14" 15" 16" 17" 18"

19" 20"

Seat Width Package:

14"W-18"W 16"W-20"W

Seat Depth:

14" 15" 16" 17" 18"

19" 20"

Seat Depth Package:

14"D-18"D 16"D-20"D

Back Cane Height:

14" 15" 16" 17" 18"

19" 20" 21" 22" 23"

24" 25"

Back Rest Height:

18" 19" 20" 21" 22"

23" 24" 25"

BACK OPTIONS

| Description | Price | HCPCS Code |
|---|----------|-------------|
| Tilt Back: | | |
| <input type="checkbox"/> Synergy Back with 1.5" thick back cushion <i>Choose Lateral Extension:</i> | Standard | |
| <input type="checkbox"/> 10" (For Back Heights 18") | \$25 | |
| <input type="checkbox"/> 12" (For Back Heights 19"-22") | \$25 | |
| <input type="checkbox"/> 16" (For Back Heights 23"-25") | \$25 | |
| <input type="checkbox"/> Synergy® Shape Back Positioning Back <i>See Synergy Shape Back order form for detailed pricing, sizes and codes</i> | Various | E2613/E2614 |
| <input type="checkbox"/> Omit Back | (-\$100) | |
| Tilt Back Options: | | |
| <input type="checkbox"/> Standard Back Canes | Standard | |
| <input type="checkbox"/> Omit Push Handles | \$125 | |
| <input type="checkbox"/> Flip Forward Canes ³ | \$150 | |
| <input type="checkbox"/> Omit Push Handles | \$125 | |
| <input type="checkbox"/> Offset (Meso) ⁴ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |
| <input type="checkbox"/> Inset (Endo) ⁴ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |
| <input type="checkbox"/> Attendant Manual Recline ⁵ <i>250lbs Weight Capacity, 180° of Recline</i> | \$750 | E1226 |
| Recline Back: | | |
| <input type="checkbox"/> Synergy Reclining Back <i>with 1.5" thick back cushion</i> | Standard | |
| <input type="checkbox"/> After Market Back Interface ⁶ | \$200 | |
| Recline Back Options: | | |
| <input type="checkbox"/> Offset (Meso) ⁷ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |
| <input type="checkbox"/> Inset (Endo) ⁷ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |

TRU-COMFORT™ SEATING

| Description | Price | HCPCS Code |
|---|-----------|------------|
| <i>Please refer to TRU-Comfort™ Cushion and TRU-Comfort Back order forms for detailed pricing, sizes and codes.</i> | | |
| TRU-Comfort™ Back Cushion: | | |
| <input type="checkbox"/> 16"W <input type="checkbox"/> 18" W <input type="checkbox"/> 20"W | \$595 | E2620 |
| Back Height⁸: | | |
| <input type="checkbox"/> 22" H <input type="checkbox"/> 25" H | No Charge | |
| Cover: | | |
| <input type="checkbox"/> Standard Allure leather-like cover | No Charge | |
| <input type="checkbox"/> Plus 4-way stretch material cover | No Charge | |
| TRU-Comfort™ Back Options: | | |
| <input type="checkbox"/> Offset (Meso) ⁹ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |
| <input type="checkbox"/> Inset (Endo) ⁹ | \$195 | |
| <input type="checkbox"/> 2" (1" Per Side) <input type="checkbox"/> 4" (2" Per Side) | | |
| <i>Please refer to TRU-Comfort™ Cushion and TRU-Comfort Back order forms for detailed pricing, sizes and codes.</i> | | |
| TRU-Comfort™ Seat Cushion General Use <i>Includes Allure leather-like cover</i> | | |
| <input type="checkbox"/> 16" x 16" - 20" x 20" | \$380 | E2601 |
| TRU-Comfort™ Plus Seat Cushion Skin Protection <i>Includes 4-way stretch material cover</i> | | |
| <input type="checkbox"/> 16" x 16" - 20" x 20" | \$380 | E2603 |

ARMRESTS

| Description | Price | HCPCS Code |
|---|-----------|------------|
| Tilt Armrests: | | |
| <input type="checkbox"/> Flip-back, Cane Mounted ¹⁰ | No Charge | |
| <input type="checkbox"/> 2-post, Flip-back, Height-adjustable, Removable ^{10,11} | No Charge | |
| <input type="checkbox"/> Single-post, Height-adjustable, Removable ¹² | No Charge | |
| <input type="checkbox"/> Quad Release Lever | | |
| <input type="checkbox"/> Removable, Single-post, Height Adjustable, Heavy Duty | No Charge | |
| <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H-14"H | | |
| <input type="checkbox"/> Pediatric, Removable ¹³ | No Charge | |
| Color: <input type="checkbox"/> Match Shroud <input type="checkbox"/> Black | | |
| <input type="checkbox"/> Quad Release Lever | | |
| <input type="checkbox"/> Height Adjustable Cantilever Armrests ¹⁴ | \$195 | |
| <i>Angle-Adjustable: - 30° through + 30°</i> | | |
| Recline Armrests: | | |
| <input type="checkbox"/> Reclining, Flip Back <i>TRU-Balance Power Recline Only</i> | | |
| <input type="checkbox"/> 8"H - 9"H <input type="checkbox"/> 10"H-11"H | | |
| <input type="checkbox"/> Height Adjustable Cantilever Armrests ¹⁴ | \$195 | |
| <i>Angle-Adjustable: - 30° through + 30°</i> | | |

FOOTNOTES

- 10" of scissor lift. Must select TB Tilt or TB Tilt & Recline in conjunction with lift. Includes 1/4 Speed Inhibit. Allowable dimensions 16" x 16" - 20" x 20". Requires Q-Logic EX. Not available with Manual Recline.
- Requires Group 34 batteries.
- Not available with Attendant Manual Recline, O2 Holder, Cane Mounted or Cantilever Arms. Tru Comfort Back requires SOF. May interfere with Lateral Supports and Waterfall Arm pads.
- Dual Post Flip-back Armrests interfere with Meso canes. Not available with Manual Recline.
- Not available with Flip Back Arms, Cantilever Arms or Shape Back. Requires Synergy Back. Pivot point 2.5" above seat pan. In iQuantum this option should be chosen under 'Depth Package Type' and also 'Back Cane Type'.
- Seat depth may be adjusted to accommodate certain after market backs. Contact your Quantum Specialist for details. Does not include brackets to mount After Market Back.
- Recline units require Cantilever Armrests. Back width cannot be less than 14"W or exceed 24"W.
- TRU-Comfort Back Rest Height is measured from the top of a 3" cushion.
- Dual Post Flip-back Armrests interfere with Meso brackets. Back width cannot be less than 16"W or exceed 22"W.
- Interferes with swing-away laterals.
- Height adjustable from 9"-13" from seat pan in 1" increments.
- Height adjustable from 11"-17" from seat pan in 1/2" increments.
- Height adjustable from 7.75"-9.75" from seat pan in 1/2" increments. Not available with full length arm pads.
- Height adjustment from 11-17" from seat pan depends on cane height. Not available with Manual Recline, Flip Forward back canes, Cane Clamp Pushbar, & Shape Back. Call for details when selecting aftermarket arm troughs.

ARMPADS

| Description | Price | HCPCS Code |
|--|-----------|------------|
| <input type="checkbox"/> Straight, Full Length (2"x14") | Standard | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Straight, Desk Length (2"x10") | No Charge | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Waterfall, Full Length (2"x14") | No Charge | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Waterfall, Desk Length (2"x10") | No Charge | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14") | \$150 pr. | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10") | \$115 pr. | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |

JOYSTICK MOUNTING BLOCKS

| Description | Price | HCPCS Code |
|---|-----------|------------|
| <input type="checkbox"/> Quick Disconnect ¹ | No Charge | |
| <input type="checkbox"/> Vertical Adjustment Clamp | \$95 | K0108 |
| <input type="checkbox"/> Figure 8 style (For use with waterfall arm pads) | No Charge | |
| <input type="checkbox"/> Figure 8 T-style (For use with recline & waterfall arm pads) | No Charge | |

POWER POSITIONING ELECTRONICS

| Description | Price | HCPCS Code |
|--|-----------|------------|
| Single Actuator Electronics | | |
| <input type="checkbox"/> Through Toggle Box (Tilt Only) | No Charge | |
| <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount | | |
| <input type="checkbox"/> Through Single Switch Module ² (Tilt Only) | \$500 | |
| <input type="checkbox"/> Through VR2 6-Key Controller (Tilt Only) | No Charge | |
| <input type="checkbox"/> Through Q-Logic EX Controller | \$1175 | E2310 |
| Multiple Actuator Electronics | | |
| <input type="checkbox"/> Through Q-Logic EX Controller | \$2375 | E2311 |
| <input type="checkbox"/> Q-Logic EX Through 4-way Toggle | \$2375 | |
| <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount | | |
| <input type="checkbox"/> Q-Logic EX Through 4-way Button Toggle | \$2375 | |
| <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount | | |
| <input type="checkbox"/> Q-Logic EX through Quad Push-button | \$2375 | |
| <input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount | | |

LAPBELTS

| Description | Price | HCPCS Code |
|---|----------|------------|
| <input type="checkbox"/> 38" Lap Belt | Standard | |
| <input type="checkbox"/> 50" Lap Belt | Standard | |
| <input type="checkbox"/> 60" Lap Belt | Standard | |
| <input type="checkbox"/> 70" Lap Belt | Standard | |
| <input type="checkbox"/> 80" Lap Belt | Standard | |
| <input type="checkbox"/> 90" Lap Belt | Standard | |
| <input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded Med, 9"-15" Hip Width | \$86 | |
| <input type="checkbox"/> Pelvic Positioning Belt, 2-Point Padded, Lg, 13"-19" Hip Width | \$91 | |
| <input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Med, 9"-15" Hip Width | \$112 | |
| <input type="checkbox"/> Pelvic Positioning Belt, 4-Point Padded, Lg, 13"-19" Hip Width | \$116 | |

HEADRESTS

| Description | Price | HCPCS Code |
|---|-------|------------|
| <input type="checkbox"/> No Headrest Required | | |
| TRU-Comfort™ Headrest³ | | |
| <input type="checkbox"/> 8" Wide Leatherette | \$205 | E0955 |
| <input type="checkbox"/> 8" TRU-Comfort Plus | \$205 | E0955 |
| <input type="checkbox"/> 10" Wide Leatherette | \$205 | E0955 |
| <input type="checkbox"/> 10" TRU-Comfort Plus | \$205 | E0955 |
| <input type="checkbox"/> Removable Headrest Bracket ⁴ | \$210 | E1028 |
| Stealth | | |
| <input type="checkbox"/> Small (9.75"W x 3"H) | \$205 | E0955 |
| <input type="checkbox"/> Large (11"W x 6"H) | \$205 | E0955 |
| <input type="checkbox"/> Pediatric Comfort Plus (6"W) | \$205 | E0955 |
| <input type="checkbox"/> Small Comfort Plus (10"W) | \$205 | E0955 |
| <input type="checkbox"/> Large Comfort Plus (14"W) | \$205 | E0955 |
| <input type="checkbox"/> Removable Headrest Bracket ⁴ | \$210 | E1028 |
| <input type="checkbox"/> Swing-away (Flip-down) Headrest Bracket ⁶ | \$210 | E1028 |

Used in Addition to Removable Headrest Bracket

HEADRESTS (cont.)

| Description | Price | HCPCS Code |
|--|-------|------------|
| Medical Applications | | |
| <input type="checkbox"/> Small (9.75"W x 3"H) | \$230 | E0955 |
| <input type="checkbox"/> Large (11"W x 6"H) | \$230 | E0955 |
| <input type="checkbox"/> Removable Headrest Bracket ⁴ | \$210 | E1028 |
| <input type="checkbox"/> Reclining Headrest (6.5"W x 7"H) | \$205 | E0955 |
| <input type="checkbox"/> Removable Bracket For Reclining Headrest ⁴ | \$210 | E1028 |
| Aftermarket Bracket | | |
| <input type="checkbox"/> Removable Bracket for Aftermarket Headrest | \$210 | E1028 |

LATERALS

| Description | Price | HCPCS Code |
|--|-------|------------|
| <input type="checkbox"/> Stealth Thoracic Laterals (pair) | | |
| <i>Bracket:</i> | | |
| <input type="checkbox"/> 125° Swing-Away Bracket ⁶ | \$450 | E1028 |
| (For Use with TRU-Comfort) | | |
| <input type="checkbox"/> Swing-Away, Rail Mount Bracket ⁶ | \$450 | E1028 |
| (For Use with Synergy Back) | | |
| <i>Pad:</i> | | |
| <input type="checkbox"/> Flat Pad, 3" x 5", Leather-Like Cover | \$200 | E0956 |
| <input type="checkbox"/> Flat Pad, 4" x 6", Leather-Like Cover | \$200 | E0956 |
| <input type="checkbox"/> Curved Pad, 2" x 3", 4-Way Stretch Cover | \$200 | E0956 |
| <input type="checkbox"/> Curved Pad, 3" x 5", 4-Way Stretch Cover | \$200 | E0956 |
| <input type="checkbox"/> Curved Pad, 4" x 6", 4-Way Stretch Cover | \$200 | E0956 |
| <input type="checkbox"/> Stealth Pelvic/Thigh Laterals | | |
| <i>Side:</i> | | |
| <input type="checkbox"/> Left Side <input type="checkbox"/> Right Side | | |
| <i>Bracket:</i> | | |
| <input type="checkbox"/> Removable Bracket (ea.) | \$250 | E1028 |
| <input type="checkbox"/> Flip Down Bracket (ea.) | \$225 | E1028 |
| <i>Pad:</i> | | |
| <input type="checkbox"/> Short Pad, 4" x 4" (ea.) | \$90 | E0956 |
| <input type="checkbox"/> Medium Pad, 4" x 8" (ea.) | \$90 | E0956 |
| <input type="checkbox"/> Long Pad, 4" x 12" (ea.) | \$90 | E0956 |

NON-POWER LEGRESTS

| Description | Price | HCPCS Code |
|--|-----------|------------|
| <input type="checkbox"/> Swing-away Legrests⁷ | | |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° | Standard | |
| <input type="checkbox"/> Straight <input type="checkbox"/> Tapered (70° only) | | |
| <input type="checkbox"/> Standard Extruded Footplates | No Charge | |
| <i>Standard Extruded Extensions:</i> | | |
| <input type="checkbox"/> 11"-13" <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5" | | |
| <i>Standard Extruded Footplates:</i> | | |
| <input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6" | | |
| <input type="checkbox"/> Angle-Adjustable Footplates ⁸ | \$150 | K0040 |
| <i>Angle Adjustable Extensions:</i> | | |
| <input type="checkbox"/> 11"-13" <input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5" | | |
| <i>Angle-Adjustable Footplates:</i> | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ⁹ <input type="checkbox"/> 6" x 11" ⁹ <input type="checkbox"/> 8" x 11" ¹⁰ | | |
| <input type="checkbox"/> Multi-Axis Footplates ⁸ | \$495 | K0040 |
| <i>Multi-Axis Extensions:</i> | | |
| <input type="checkbox"/> 11"-16.5" <input type="checkbox"/> 15.5"-19.5" | | |
| <i>Multi-Axis Footplates:</i> | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ⁹ <input type="checkbox"/> 6" x 11" ⁹ <input type="checkbox"/> 8" x 11" ¹⁰ | | |
| <input type="checkbox"/> High Mount, Clamp-On Footrest | | |
| <input type="checkbox"/> Angle Adjustable | \$326 | K0040 |
| <input type="checkbox"/> Small, 4" x 7.75" <input type="checkbox"/> Large, 4.25" x 9.25" | | |
| <input type="checkbox"/> Multi-Axis | \$412 | K0040 |
| <input type="checkbox"/> Small, 4" x 7.75" <input type="checkbox"/> Large, 4.25" x 9.25" | | |

FOOTNOTES

- Not available with waterfall arm pads.
- Single switch not included. See specialty controls order form for a list of available switches.
- TRU-Comfort™ Headrest is measured at the front of the pad.
- When headrest pad is chosen i-Quantum automatically selects Removable Headrest Bracket.
- Not available with power recline.
- Swing-away laterals interfere with flip-back armrests.
- Extensions measured from top of seatpan to footplate.
- Includes heel loops.
- Requires seat width 18" or greater.
- Requires seat width 20" or greater.

NON-POWER LEGRESTS (cont.)

| Description | Price | HCPCS Code |
|---|------------------|--------------|
| <input type="checkbox"/> Heavy Duty Drop-In Legrests¹ | \$420 | |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° | | |
| <input type="checkbox"/> Angle Adjustable Footplates² | \$150 | K0040 |
| Angle Adjustable Extension Length: | | |
| <input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5"-18" | | |
| Angle Adjustable Footplate Size: | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴ | | |
| <input type="checkbox"/> Elevating Legrests¹ | \$325 | E0990 |
| <input type="checkbox"/> Standard Footplates | | |
| Standard Extruded Extension Length: | | |
| <input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17" | | |
| Standard Extruded Footplate Size: | | |
| <input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6" | | |
| <input type="checkbox"/> Angle Adjustable Footplates² | \$150 | K0040 |
| Angle Adjustable Extension Length: | | |
| <input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5"-18" | | |
| Angle Adjustable Footplate Size: | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴ | | |
| <input type="checkbox"/> Calf Pad Adjustment | | |
| <input type="checkbox"/> 1 Hole <input type="checkbox"/> 3 Hole <input type="checkbox"/> 6 Hole | | |
| <input type="checkbox"/> 90° Pediatric Swing-Away Legrests¹ | No Charge | |
| <input type="checkbox"/> Angle Adjustable Footplates² | | |
| Angle Adjustable Extension Length: | | |
| <input type="checkbox"/> 7"-10" | | |
| Angle Adjustable Footplate Size: | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" | | |

POWER ELEVATING LEGRESTS

| Description | Price | HCPCS Code |
|---|----------------|--------------|
| <input type="checkbox"/> Articulating Swing-Away Legrests, 4" Articulation Range | \$2,395 | E1010 |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° | | |
| <input type="checkbox"/> Quad Release Lever | | |
| Calf Pads ⁵ : | | |
| <input type="checkbox"/> 6" x 5.35" <input type="checkbox"/> 8" x 5.86" | | |
| Angle Adjustable Extension Length ¹ : | | |
| <input type="checkbox"/> 11"-13" <input type="checkbox"/> 14.5"-18" | | |
| Angle Adjustable Footplate Size: | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴ | | |
| Operational Method: | | |
| <input type="checkbox"/> Independent <input type="checkbox"/> Combined | | |
| <input type="checkbox"/> Articulating Drop-In Legrests | \$2,395 | E1010 |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° | | |
| <input type="checkbox"/> 4" Range <input type="checkbox"/> 6" Range | | |
| Calf Pads ⁵ : | | |
| <input type="checkbox"/> 6" x 5.35" <input type="checkbox"/> 8" x 5.86" | | |
| Angle Adjustable Extension Length ¹ : | | |
| <input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5"-18" | | |
| Angle Adjustable Footplate Size: | | |
| <input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8" | | |
| <input type="checkbox"/> 6" x 8" ³ <input type="checkbox"/> 6" x 11" ³ <input type="checkbox"/> 8" x 11" ⁴ | | |
| Operational Method: | | |
| <input type="checkbox"/> Independent <input type="checkbox"/> Combined | | |

FOOT RIGGING ACCESSORIES

| Description | Price | HCPCS Code |
|---|--------------|--------------|
| <input type="checkbox"/> Wheel Bumpers⁶ | \$25 | |
| <input type="checkbox"/> Heel Loops (pair)⁷ | \$38 | E0951 |
| <input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair)⁸ | \$65 | K0108 |
| <input type="checkbox"/> ELR/ALR Adductor Buttons (pair) | \$100 | K0108 |
| <input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair) | \$464 | K0108 |
| <input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair) | \$464 | K0108 |
| <input type="checkbox"/> Residual Limb Support | \$300 | E1020 |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| <input type="checkbox"/> Ankle Huggers (pair) | \$126 | K0108 |
| <input type="checkbox"/> Small <input type="checkbox"/> Medium | | |

SYNERGY® CUSHIONS

| Description | Price | HCPCS Code |
|--|--------------|--------------|
| <input type="checkbox"/> No Cushion Required | | |
| <i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i> | | |
| <input type="checkbox"/> Simplicity Cushion <i>General Use</i> | | |
| <input type="checkbox"/> 16" x 16" - 20" x 20" | \$98 | E2601 |
| <input type="checkbox"/> Solution® 1 Cushion <i>Skin Protection</i> | | |
| <input type="checkbox"/> 14" x 14" - 20" x 18" | \$380 | E2603 |
| <input type="checkbox"/> 20" x 20" - 20" x 20" | \$430 | E2603 |
| <input type="checkbox"/> Solution® Cushion <i>Skin Protection & Positioning</i> | | |
| <input type="checkbox"/> 10" x 10" - 20" x 18" | \$380 | E2607 |
| <input type="checkbox"/> 20" x 20" - 20" x 20" | \$430 | E2607 |
| <input type="checkbox"/> Spectrum® Gel Cushion <i>Positioning</i> | | |
| <input type="checkbox"/> 12" x 10" - 20" x 19" | \$480 | E2605 |
| <input type="checkbox"/> 20" x 20" - 20" x 20" | \$530 | E2605 |
| <input type="checkbox"/> Spectrum® Visco Cushion <i>Skin Protection & Positioning</i> | | |
| <input type="checkbox"/> 10" x 10" - 20" x 18" | \$480 | E2607 |
| <input type="checkbox"/> 20" x 20" - 20" x 20" | \$530 | E2607 |
| <input type="checkbox"/> Spectrum Air <i>Adjustable Skin Protection</i> | | |
| <input type="checkbox"/> 14" x 12" - 20" x 20" | \$480 | E2622 |
| <input type="checkbox"/> Spectrum Air Contour <i>Adjustable Skin Protection & Positioning</i> | | |
| <input type="checkbox"/> 14" x 12" - 20" x 20" | \$499 | E2624 |

MISCELLANEOUS POSITIONING COMPONENTS

| Description | Price | HCPCS Code |
|---|--------------|--------------|
| <input type="checkbox"/> Therafin Flip-down Knee Abductor Bracket | \$190 | E1028 |
| Pad: | | |
| <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large | \$125 | E0957 |
| <input type="checkbox"/> Stealth Shoe Boxes (pair) | \$475 | K0108 |
| <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large | | |
| <input type="checkbox"/> Stealth Padded Shoe Holders (pair) | \$166 | K0108 |
| <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Large | | |
| <input type="checkbox"/> Otto Bock Arm Channels (ea.) | \$102 | E2209 |
| Side: | | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| Size: | | |
| <input type="checkbox"/> Small <input type="checkbox"/> Medium <input type="checkbox"/> Large | | |
| <input type="checkbox"/> Otto Bock Flat Hand Pad (ea.) | \$51 | |
| Side: | | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |
| Size: | | |
| <input type="checkbox"/> Medium <input type="checkbox"/> Large | | |

Please see positioning components order form for additional types, sizes and details.

ACCESSORIES

| Description | Price | HCPCS Code |
|---|--------------|--------------|
| <input type="checkbox"/> Quantum Backpack⁹ | \$60 | |
| <input type="checkbox"/> Oxygen Holder¹⁰ | \$150 | E2208 |
| <input type="checkbox"/> Recline Accessory Mounting Bracket¹¹ | \$50 | |
| <input type="checkbox"/> Weather Cover, Medium, Black | \$158 | |
| <input type="checkbox"/> Cup Holder | \$17 | |
| Side: | | |
| <input type="checkbox"/> Left <input type="checkbox"/> Right | | |

BATTERY OPTIONS

| Description | Price | HCPCS Code |
|--|------------------|--------------|
| <input type="checkbox"/> Group 24, 70 Amp, Gel (Universal) | \$352 ea. | E2363 |
| <input type="checkbox"/> Group 24, 70 Amp, Gel (MK) | \$352 ea. | E2363 |
| <input type="checkbox"/> Group 34, 60-Amp Gel¹² | \$330 ea. | K0108 |
| Installation: | | |
| <input type="checkbox"/> Batteries Installed <input type="checkbox"/> Do Not Install | | |

FOOTNOTES

- Extensions measured from top of seatpan to footplate.
- Includes heel loops.
- Requires seat width 18" or greater.
- Requires seat width 20" or greater.
- Calf pads are selected in the Rigging Accessory section of i-Quantum.
- Standard on Power Elevating Legrests.
- Standard with angle adjustable footplates.
- Available only with swing-away legrests.
- Backpack may interfere w/ other rear accessories.
- Oxygen holder on Tilt with Tru Comfort back or on Power Recline requires SOF.
- 1" adaptor tube allows cane mounted accessories to mount to power recline back or siderail.
- Required with Super Low Option.