

# Tru-Balance® 2 Power Positioning Order Form



Discount: \_\_\_\_\_  
Terms: \_\_\_\_\_  
Freight: \_\_\_\_\_  
Taken By: \_\_\_\_\_

Account Number: \_\_\_\_\_  
Provider Name: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
E-Mail: \_\_\_\_\_  
Date: \_\_\_\_\_  
PO#: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Ship To Address: \_\_\_\_\_  
City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice. All prices are MSRP.*

## BASES: SEATING SPECIFICATIONS:

- 250lb Weight Capacity. Max Dimensions 20"W x 20"D**  
 J6
- 300lb Weight Capacity. Max Dimensions 20"W x 20"D**  
 Q610
- 300lb Weight Capacity. Max Dimensions 22"W x 22"D**  
 Rival
- 300lb Weight Capacity. Max Dimensions 22"W x 22"D**  
 Q6 Edge
- 300lb Weight Capacity. Max Dimensions 24"W x 24"D**  
 Q6000Z

**\*\*For more detailed information, please refer to the base order forms\*\***

**For Discontinued Bases, Contact Inside Sales Specialist**

- Seat Width:**  
 14"  15"  16"  17"  18"  
 19"  20"  21"  22"  23"  24"
- Seat Width Package:**  
 14"W-18"W  16"W-20"W  20"W-24"W
- Seat Depth:**  
 14"  15"  16"  17"  18"  
 19"  20"  21"  22"  23"  24"
- Seat Depth Package:**  
 14"D-18"D  15"D-20"D  20"D-24"D
- Back Cane Height:**  
 14"  15"  16"<sup>1</sup>  17"<sup>1</sup>  
 18"  19"  20"  21"  
 22"  23"  24"  25"

1. Selection of a Synergy® Back & Flip Up Armrests disqualified 16" or 17"H back canes.

## TRU-BALANCE® 2 SEATING SYSTEMS BACK OPTIONS

Description	Price	HCPCS Code
List base serial # for off chair seating		
S/N: _____		
<input type="checkbox"/> TRU-Balance® 2 Power Tilt (55°) <i>55 degree center of gravity (CG) power tilt seating system.</i>	\$4,995	E1002
<input type="checkbox"/> TRU-Balance® 2 Power Recline (168°) <i>168 degrees of power recline with full drive lock-out set at 20 degrees</i>	\$5,595	E1004
<input type="checkbox"/> TRU-Balance® 2 Power Tilt (55°) & Recline (168°) <i>55 degree center of gravity (CG) power tilt seating system with 168 degrees of power recline Full drive lock-out at 20 degrees</i>	\$9,595	E1007
<input type="checkbox"/> TRU-Balance® Lift (10") <sup>1</sup> <i>Must be selected in conjunction with TB Power Tilt or TB Power Tilt &amp; Recline</i>	\$2,895	
<input type="checkbox"/> TRU-Balance® 6" Power Elevate <sup>2</sup> <i>Includes 1/4 Speed Inhibit.</i>	\$1,895	
<input type="checkbox"/> TRU-Balance® 7" Power Elevate <sup>3</sup> <i>Includes 1/4 Speed Inhibit.</i>	\$1,895	
<input type="checkbox"/> TRU-Balance® 400lb Power Tilt (55°) <sup>4</sup> <i>55 degree center of gravity (CG) power tilt seating system. 400lb wt cap.</i>	\$4,995	E1002

- 10" of scissor Lift, must select TB Tilt or TB Tilt & Recline in conjunction with lift, 1/4 Speed Inhibit, 300 lbs. max. wt., Allowable dimensions 16" x 16" - 20" x 20" on most bases, Requires Q-Logic 2 EX, Not available with Manual Recline or Power Recline Only.
- Only available on Q6 Edge. Maximum seat size 20"W x 20"D. Available with all power positioning except 10" Lift. Increases STF Height. Includes 1/4 speed inhibit.
- Only available on Q6000Z and Q6 Edge. Maximum seat size 20"W x 20"D. Available with all power positioning except 10" Lift. Increases STF Height. Includes 1/4 speed inhibit.
- Q6000Z HD only.

Description	Price	HCPCS Code
-------------	-------	------------

- Tilt Back:**
- Synergy® Shape Back *Positioning Back* Various E2613/E2614  
*See Synergy Shape Back Order Form For Detailed Pricing, Sizes and Codes*
- Omit Back (-\$100)
- Synergy® Back with 1.5" Thick Back Cushion No Charge
- Overall Back*  18" H  19" H  20" H  21" H  
*Height:*  22" H  23" H  24" H  25" H
- Choose Lateral Extensions<sup>1</sup>:*
- 6" for 13"-16"H Back \$25  
 10" for 17"-18"H Back \$25  
 12" For 19"-22"H Back \$25  
 16" For 23"-25"H Back \$25
- TRU-Comfort™ Back *4-way Stretch Material Cover*
- Width:*  16"W  18" W  20"W \$595 E2620
- Overall Back*  20" H  21" H  22" H  23" H  
*Height<sup>2</sup>:*  24" H  25" H  27" H  28" H
- Width:*  17"W  19"W \$595 E2620
- Overall Back*  22" H  23" H  24" H  25" H  
*Height<sup>2</sup>:*  27" H  28" H
- Width:*  22"W  24"W \$695 E2621
- Overall Back*  24" H  25" H  27" H  28" H  
*Height<sup>2</sup>:*

1. Lateral extensions are required to mount lateral thoracic supports.
2. Overall Back height measurement is measured from the seat pan to the top of the back, +/- 1/2".

**BACK OPTIONS (CONT'D)**

Description	Price	HCPCS Code
<b>Tru Comfort Color Shroud Option:</b>		
<input type="checkbox"/> Black (Standard color option)		
<b>Metallic Colors:\$100</b>		
<input type="checkbox"/> Candy Apple Red		<input type="checkbox"/> Emerald Green
<input type="checkbox"/> Viper Blue		<input type="checkbox"/> Sapphire Blue
<input type="checkbox"/> Sunburst Orange		<input type="checkbox"/> Tru Blue <sup>1</sup>
<input type="checkbox"/> Pearl Pink		<input type="checkbox"/> Orion Yellow <sup>1</sup>
<input type="checkbox"/> Silver		
<b>Tilt Back Options:</b>		
<input type="checkbox"/> Standard Back Canes	Standard	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Flip Forward Back Canes <sup>2</sup>	\$150	
<input type="checkbox"/> Omit Push Handles	\$125	
<input type="checkbox"/> Offset (Meso) <sup>3</sup>	\$195	
<input type="checkbox"/> 2" (1" Per Side)		<input type="checkbox"/> 4" (2" Per Side)
<input type="checkbox"/> Inset (Endo) <sup>3</sup>	\$195	
<input type="checkbox"/> 2" (1" Per Side)		<input type="checkbox"/> 4" (2" Per Side)
<input type="checkbox"/> Attendant Manual Recline <sup>4</sup>	\$750	E1226

**Power Recline Back:**

<input type="checkbox"/> Synergy® Reclining Back with 1.5" Thick Back Cushion	Standard	
<input type="checkbox"/> TRU-Comfort™ Back 4-way Stretch Material Cover		
Width: <input type="checkbox"/> 16"W <input type="checkbox"/> 17" W <sup>5</sup> <input type="checkbox"/> 18"W	\$595/\$795	E2620
<input type="checkbox"/> 19"W <sup>5</sup> <input type="checkbox"/> 20"W		
<input type="checkbox"/> 22"W	\$695/\$795	E2621
Overall Back Height <sup>6</sup> : <input type="checkbox"/> 20" H <input type="checkbox"/> 21" H <input type="checkbox"/> 22"H <input type="checkbox"/> 23"H		
<input type="checkbox"/> 24"H <input type="checkbox"/> 25"H <input type="checkbox"/> 27"H <input type="checkbox"/> 28"H		
<input type="checkbox"/> After-Market Back Interface <sup>7</sup>	\$200	

**Recline Back Options:**

<input type="checkbox"/> Offset (Meso) <sup>3</sup>	\$195
<input type="checkbox"/> 2" (1" Per Side)	<input type="checkbox"/> 4" (2" Per Side)
<input type="checkbox"/> Inset (Endo) <sup>3</sup>	\$195
<input type="checkbox"/> 2" (1" Per Side)	<input type="checkbox"/> 4" (2" Per Side)

**Recline Color Shroud Option:**

<input type="checkbox"/> Black (Standard color option)	
<b>Metallic Colors:N/C</b>	
<input type="checkbox"/> Candy Apple Red	<input type="checkbox"/> Emerald Green
<input type="checkbox"/> Viper Blue	<input type="checkbox"/> Sapphire Blue
<input type="checkbox"/> Sunburst Orange	<input type="checkbox"/> Tru Blue <sup>1</sup>
<input type="checkbox"/> Pearl Pink	<input type="checkbox"/> Orion Yellow <sup>1</sup>
<input type="checkbox"/> Silver	

- SOF required.
- Not available with Attendant Manual Recline, Vent Tray, or O2 Holder or 400LB tilt. May interfere with Lateral Supports and Waterfall Armpads.
- Dual Post Flip-back Armrests interfere with offset canes. Not available with Manual Recline.
- Not available with Flip Back Arms, Cantilever Arms, Shape Back, Tru Comfort back, power elevating seat or 10" Lift. Requires 22"H canes. Pivot point 2.5" above seat pan.
- SOF is required with no extended lead time
- Overall Back height measurement is measured from the seat pan to the top of the back. All back heights except 25" and 28" are SOF with no extended lead time.
- Overall Seat depth may need to be adjusted to accommodate certain backs. Contact your Quantum Specialist for details. Does not include brackets to mount after market back.

**JOYSTICK RECEIVER MOUNT**

Description	Price
<input type="checkbox"/> Universal Joystick Receiver Mount	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Omit Universal Joystick Receiver Mount	No Charge

**SYNERGY® ARMRESTS**

Description	Price
<b>Tilt Armrests:</b>	
<input type="checkbox"/> Flip Back, Cane Mounted Armrests <sup>1</sup>	No Charge
<input type="checkbox"/> 2-Post, Flip-back, Height Adjustable, Removable <sup>2</sup>	No Charge
<input type="checkbox"/> Single Post Height Adjustable Removable <sup>3</sup>	No Charge
<input type="checkbox"/> Quad Release Lever	
<input type="checkbox"/> Removable, Single Post, Height Adjustable, Heavy Duty	No Charge
Height: <input type="checkbox"/> 8"H - 12"H <input type="checkbox"/> 10"H - 14"H	
<input type="checkbox"/> Pediatric, Removable <sup>4</sup>	No Charge
<input type="checkbox"/> Quad Release Lever	
<input type="checkbox"/> Height Adjustable Cantilever Armrests <sup>5</sup>	\$195
	Angle Adjustable: -30° through +30°
<b>Power Recline Armrests:</b>	
<input type="checkbox"/> Reclining, Flip Back	Standard
Height: <input type="checkbox"/> 8"H - 9"H <input type="checkbox"/> 10"H - 11"H	
<input type="checkbox"/> Height Adjustable Cantilever Armrests <sup>1</sup>	\$195
	Angle Adjustable: -30° through +30°

**SYNERGY® ARMPADS**

Description	Price
<input type="checkbox"/> Straight, Full Length (2"x14")	Standard
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Straight, Desk Length (2"x10")	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Full Length (2"x14") <sup>1</sup>	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Waterfall, Desk Length (2"x10") <sup>1</sup>	No Charge
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Falcon Gel, Full Length, (3.5"x14")*	\$80.00 ea.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Full Length Gel Armpad Mounting Plate (Required)*	\$56.69 ea.
<input type="checkbox"/> Falcon Gel, Desk Length, (3.5"x10")*	\$62.50 ea.
<input type="checkbox"/> Left <input type="checkbox"/> Right	
<input type="checkbox"/> Desk Length Gel Armpad Mounting Plate (Required)*	\$53.52 ea.

1. Not available w/ Cane Mounted, Flip Back Armrests combined with a Synergy Back.

**POWER POSITIONING ELECTRONICS**

Description	Price	HCPCS Code
Choose Power Positioning Electronics From The Section That Matches The Electronics On The Base.		
<u>Q-Logic 2 Electronics</u>		
<b>Single Actuator System:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Q-Logic 2 EX Controller		
<input type="checkbox"/> Single Actuator Module (AM1)	\$1,175	E2310 <sup>1</sup>
<input type="checkbox"/> Multiple Actuator Module (AAM) <sup>2</sup>	\$2,375	E2311 <sup>1</sup>
<b>Multiple Actuator System:</b>		
<input type="checkbox"/> Through Q-Logic 2 EX Controller	\$2,375	E2311 <sup>1</sup>
<input type="checkbox"/> Through 4-Way Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through 4-Way Button Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Quad Pushbutton	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		

1. Applies to covered actuators only.  
2. For adding up to 4 additional actuator functions in the field.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**POWER POSITIONING ELECTRONICS (CONT'D)      NON-POWER LEGRESTS**

Description	Price	HCPCS Code
<u>Remote Plus Electronics</u>		
<b>Single Actuator System:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Remote+ Controller	\$1,175	E2310 <sup>1</sup>
<b>Multiple Actuator System:</b>		
<input type="checkbox"/> Through Remote+ Controller	\$2,375	E2311 <sup>1</sup>
<input type="checkbox"/> Through 4-Way Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through 4-Way Button Toggle	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through Quad Pushbutton	\$2,375	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<u>NE+ Electronics</u>		
<b>Single Actuator System:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through NE+ 6-Key Joystick	\$1,175	E2310 <sup>1</sup>
<b>Two Actuator System:</b>		
<input type="checkbox"/> Through NE+ 6-Key Controller	\$2,375	E2311 <sup>1</sup>
<u>NE Electronics</u>		
<b>Tilt Only:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<u>VR2 Electronics</u>		
<b>Tilt Only:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through VR2 6-Key Joystick	No Charge	
<u>VSI Electronics</u>		
<b>Tilt Only:</b>		
<input type="checkbox"/> Through Toggle Box (Tilt Only)	No Charge	
<input type="checkbox"/> Left Mount <input type="checkbox"/> Right Mount		
<input type="checkbox"/> Through VSI 6-Key Joystick		

1. Applies to covered actuators only.

**HEADRESTS**

Description	Price	HCPCS Code
-------------	-------	------------

<input type="checkbox"/> No Headrest Required		
<b>TRU-Comfort™ Headrest<sup>1</sup>:</b>		
<input type="checkbox"/> 8" W Leatherette	\$205	E0955
<input type="checkbox"/> 8"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> 10" W Leatherette	\$205	E0955
<input type="checkbox"/> 10"W TRU-Comfort Plus	\$205	E0955
<input type="checkbox"/> Removable Headrest Bracket	\$210	E1028
<b>Stealth:</b>		
<input type="checkbox"/> Small (9.75"W x 3"H) (STL-OTSM250)	\$205	E0955
<input type="checkbox"/> Large (11"W x 6"H) (STL-OTLG150)	\$205	E0955
<input type="checkbox"/> Pediatric Comfort Plus (6"W) (CPN350)	\$142	E0955
<input type="checkbox"/> Small Comfort Plus (10"W) (CP250)	\$142	E0955
<input type="checkbox"/> Large Comfort Plus (14"W) (CP150)	\$142	E0955
<input type="checkbox"/> Removable Headrest Bracket <sup>2</sup> (HMO475P)	\$210	E1028
<input type="checkbox"/> Swing-away (Flip-down) Headrest	\$210	
Bracket <sup>2</sup> (FDM350)		
Used in Addition to Removable Headrest Bracket		
<b>Medical Applications:</b>		
<input type="checkbox"/> Removable Bracket For Reclining Headrest	\$210	E1028
<b>Aftermarket Bracket:</b>		
<input type="checkbox"/> Stealth's World's Best Removable Headrest	\$210	E1028
Hardware (TWB485)		

1.. Measurement is taken from the inside of the headrest.  
2. Not available on power recline.

Description	Price	HCPCS Code
<input type="checkbox"/> <b>60-90° Swing-away Legrests<sup>1</sup></b>	<b>Standard</b>	
(Adjustable in 5° increments)		
<input type="checkbox"/> Standard Extruded Footplates	No Charge	
Standard Extruded Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Standard Extruded Footplates:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle-Adjustable Footplates <sup>2</sup>	\$150	K0040 <sup>3</sup>
Angle Adjustable Extensions:		
<input type="checkbox"/> 13"-16.5" <input type="checkbox"/> 15.5"-19.5"		
Angle-Adjustable Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>4</sup> <input type="checkbox"/> 6" x 11" <sup>4</sup> <input type="checkbox"/> 8" x 11" <sup>5</sup>		
<input type="checkbox"/> Multi-Axis Footplates	\$495	K0040 <sup>3</sup>
Multi-Axis Extensions:		
<input type="checkbox"/> 15.5"-19.5"		
Multi-Axis Footplates:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>4</sup> <input type="checkbox"/> 6" x 11" <sup>4</sup> <input type="checkbox"/> 8" x 11" <sup>5</sup>		
<input type="checkbox"/> <b>90° Pediatric Swing-Away Legrests<sup>6</sup></b>	<b>\$162</b>	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>		
Angle Adjustable Extension Length:		
<input type="checkbox"/> 7"-10"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6"		
<input type="checkbox"/> <b>Heavy Duty Drop-In Legrests<sup>6</sup></b>	<b>\$420</b>	
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$150	K0040 <sup>3</sup>
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>4</sup> <input type="checkbox"/> 6" x 11" <sup>4</sup> <input type="checkbox"/> 8" x 11" <sup>5</sup>		
<input type="checkbox"/> <b>Elevating Legrests<sup>6</sup></b>	\$325	E0990
<input type="checkbox"/> Standard Footplates		
Standard Extruded Extension Length:		
<input type="checkbox"/> 11"-14" <input type="checkbox"/> 14"-17"		
Standard Extruded Footplate Size:		
<input type="checkbox"/> 7.25" x 6" <input type="checkbox"/> 8.5" x 6"		
<input type="checkbox"/> Angle Adjustable Footplates <sup>2</sup>	\$150	K0040 <sup>3</sup>
Angle Adjustable Extension Length:		
<input type="checkbox"/> 12"-14.5" <input type="checkbox"/> 14.5" -18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>4</sup> <input type="checkbox"/> 6" x 11" <sup>4</sup> <input type="checkbox"/> 8" x 11" <sup>5</sup>		
<input type="checkbox"/> <b>High Mount, Clamp-On Footplates*</b>		
<input type="checkbox"/> Angle Adjustable		K0040 <sup>3</sup>
<input type="checkbox"/> Small, 4" x 7.75"	\$356	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$356	
<input type="checkbox"/> Multi-Axis		K0040 <sup>3</sup>
<input type="checkbox"/> Small, 4" x 7.75"	\$456	
<input type="checkbox"/> Large, 4.25" x 9.25"	\$460	

1. Possible caster wheel interference when adjusted to less than 70 degrees.  
2. Includes Heel Loops  
3. K0040 is not separately reimbursable on Group 2 units.  
4. Causes interference on seat widths less than 18".  
5. Causes interference on seat widths less than 20".  
6. Extensions are measured from the seatpan to footplate.

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

**ADJUSTABLE CENTER MOUNT FOOT PLATFORM**

Description	Price
<input type="checkbox"/> Adjustable Center Mount Foot Platform <sup>1</sup> <i>Includes 10"W x 11"D Split Foot Platform</i>	\$250
Extensions: <sup>2</sup> <input type="checkbox"/> 10"-16" Ext. <input type="checkbox"/> 13"-19" Ext.	
Calf Supports: <sup>3</sup>	
<input type="checkbox"/> Calf Pads (pair)	
<input type="checkbox"/> Small (6"H x 6"W)	\$150
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)	\$150
<input type="checkbox"/> Calf Panel (ea.)	
<input type="checkbox"/> Small (8"H x 18"W)	\$150
<input type="checkbox"/> Large (8"H x 24"W) <sup>4</sup>	\$200
Optional Modular Footplate Extension Components/Kits <sup>3,5</sup>	
<input type="checkbox"/> 4" Width Extension Kit	\$50
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75
<input type="checkbox"/> 6-Piece Extension Stock Kit <i>(Includes all above pieces)</i>	\$95
Additional Optional Accessories <sup>3,5</sup> :	
<input type="checkbox"/> Footplate Connector	\$49
<input type="checkbox"/> Wheel Bumpers, pair	\$25

1. Not available on 614 series.
2. Extension length is measured from top of the seat pan to footplate.
3. Choose in Rigging Accessory Section in iQuantum.
4. Only available on 20"W or greater.
5. Assembly required. May require seating configuration adjustments to avoid front caster interference.

**POWER ELEVATING LEGRESTS**

Description	Price	HCPCS Code
<input type="checkbox"/> 70° Articulating Swing-Away Legrests, 4" Articulation Range <sup>1</sup>	\$2,395	E1010
<input type="checkbox"/> Quad Release Lever		
Calf Pads <sup>2</sup> :		
<input type="checkbox"/> 6" x 5.35" <input type="checkbox"/> 8" x 5.86"		
Angle Adjustable Extension Length <sup>3</sup> :		
<input type="checkbox"/> 11"-13" <input type="checkbox"/> 14.5"-18"		
Angle Adjustable Footplate Size:		
<input type="checkbox"/> 4" x 6" <input type="checkbox"/> 5" x 6" <input type="checkbox"/> 5" x 8"		
<input type="checkbox"/> 6" x 8" <sup>4</sup> <input type="checkbox"/> 6" x 11" <sup>4</sup> <input type="checkbox"/> 8" x 11" <sup>5</sup>		
Operational Method:		
<input type="checkbox"/> Independent <input type="checkbox"/> Combined		

1. Not available on J6 or Q610 Series
2. Choose in Rigging Accessory Section in iQuantum.
3. Extensions are measured from top of seat pan to footplate.
4. Only available with seat width 18" or greater.
5. Only available with seat width 20" or greater.

**CENTER MOUNT ARTICULATING FOOT PLATFORM**

Description	Price	HCPCS Code
<input type="checkbox"/> Center Mount Articulating Foot Platform	\$2895	K0108
<i>Includes 10"W x 11"D Split Foot Platform</i>		
Extensions <sup>1</sup> : <input type="checkbox"/> 10"-16" Ext. <input type="checkbox"/> 13"-19" Ext.		
Calf Supports <sup>2</sup> :		
<input type="checkbox"/> Calf Pads (pair)		
<input type="checkbox"/> Small (6"H x 6"W)		
<input type="checkbox"/> Large (8"H x 8"W, Tapered to 6"W)		
<input type="checkbox"/> Calf Panel (ea.)		
<input type="checkbox"/> Small (8"H x 18"W)	\$150	
<input type="checkbox"/> Large (8"H x 24"W) <sup>3</sup>	\$200	
Optional Modular Footplate Extension Components/Kits <sup>2,4</sup>		
<input type="checkbox"/> 4" Width Extension Kit	\$50	
<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	\$50	
<input type="checkbox"/> 4" Width Extension with 2" Bell Shaped, Length Extension Kit	\$75	
<input type="checkbox"/> 6-Piece Extension Stock Kit <i>(Includes all above pieces)</i>	\$95	
Additional Optional Accessories <sup>4</sup> :		
<input type="checkbox"/> Footplate Connector	\$49	
<input type="checkbox"/> Wheel Bumpers, pair	\$25	

1. Extension length is measured from top of the seat pan to footplate.
2. Choose in Rigging Accessory Section in iQuantum.
3. Only available on 20"W or greater.
4. Assembly required. May require seating configuration adjustments to avoid front caster interference.

**FOOT RIGGING ACCESSORIES**

Description	Price	HCPCS Code
<input type="checkbox"/> Wheel Bumpers	\$25	
<input type="checkbox"/> Heel Loops (pair) <sup>1</sup>	\$38	E0951
<input type="checkbox"/> Gel Neoprene Footrest Sleeve (pair) <sup>2</sup>	\$65	K0108
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	\$125	K0108
<input type="checkbox"/> POSAlinc Knee Adductor for Swing-Away Legrests (pair)*	\$512	K0108
<input type="checkbox"/> POSAlinc Rail Mounted Knee Adductor (pair)*	\$574	K0108
<input type="checkbox"/> Residual Limb Support	\$300	E1020
<input type="checkbox"/> Left <input type="checkbox"/> Right		
<input type="checkbox"/> Bodypoint Ankle Huggers (pair)*	\$134	K0108
<input type="checkbox"/> Small <input type="checkbox"/> Medium		

1. Standard with angle adjustable footplates.
2. Available only with swing-away legrests.

**SYNERGY® SEAT CUSHIONS**

Description	Price	HCPCS Code
<input type="checkbox"/> No Cushion Required		
<i>Please refer to Synergy order forms for detailed pricing, sizes and codes.</i>		
<b>Simplicity Cushion</b> <i>General Use</i>		
<input type="checkbox"/> 16" x 16" - 20" x 20"	\$98	E2601
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$98	E2602
<b>Solution@ 1 Cushion</b> <i>Skin Protection:</i>		
<input type="checkbox"/> 14" x 14" - 20" x 24"	\$380	E2603
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2604
<b>Solution@ Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$380	E2607
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2608
<b>Spectrum@ Gel Cushion</b> <i>Positioning</i>		
<input type="checkbox"/> 12" x 10" - 20" x 24"	\$480	E2605
<input type="checkbox"/> 22" x 16" - 24" x 24"	\$430	E2606
<b>Spectrum@ Visco Cushion</b> <i>Skin Protection &amp; Positioning</i>		
<input type="checkbox"/> 10" x 10" - 20" x 24"	\$480	E2607
<input type="checkbox"/> 22" x 16" - 32" x 28"	\$530	E2608

Items marked with a (\*) should be chosen in the Positioning Components screen in iQuantum

