



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum® 4Front® 2 Order Form with Contoured Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

QUANTUM
4FRONT²



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 3S-C.....\$8,870
Part: 4FRONT 2 3S-C. HCPCS: K0849
Captain seating standard.

Quantum 4Front 2 3SP-C.....\$8,559
Part: 4FRONT 2 3SP-C. HCPCS: K0857
Contoured seat (captain), single power option for power adjustable seat height.

Quantum 4Front 2 3S-SS.....\$8,870
Part: 4FRONT 2 3S-SS. HCPCS: K0848
Solid seat pan seating standard.

Quantum 4Front 2 3SP-SS.....\$8,559
Part: 4FRONT 2 3SP-SS. HCPCS: K0856
Solid seat pan seating standard, single power option for power adjustable seat height.

Account Number: _____ Date: _____ Quote Name/Marked for: _____

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1a. Standard Shrouds

 Red Corvette Part: DGN2003028-PRD1	 Deep Purple Part: DGN2003028-PPR5	 Midnight Blue Part: DGN2003028-PBU10	 Green Machine Part: DGN2003028-PGN11
 Back in Black Part: DGN2003028-PBK1	 Raspberry Beret Part: DGN2003028-PPK4	 Electric Blue Part: DGN2003028-PBU9	 Military Green Part: DGN2003028-PGN8
 Night Train Part: DGN2003028-PBK3	 Pretty in Pink Part: DGN2003028-PPK3	 Blue Jean Part: DGN2003028-PBU11	
 Grey Street Part: DGN2003028-PSL3	 Evolution Orange Part: DGN2003028-POR7	 Lemon Crush Part: DGN2003028-PYL4	

2.1b. Hydrographic Shrouds

 Military Digital Camo \$263 Part: DGN2003028-WCM7	 Faded Glory \$263 Part: DGN2003028-WFL2
---	---

3. POWER ADJUSTABLE SEAT HEIGHT

Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

10" Power Lift for Captain Seat	\$3,725
Part: PTOINDV3979. HCPCS: E2298/E2300^	
For use with K0856 and K0857 base code.	

4. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera	\$368
Part: MEC1910365	

1. Optional Camera Mounting Bracket

Backup Camera Captain's Headrest Bracket	\$21	Backup Camera Track Mount Bracket	\$27
Part: MEC1901231		Part: MEC2104227	

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket	\$105
Part: MEC2106123	



5. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color.

5.1. Drive Wheels

3"x14" Flat Free Black Tires with Accent	Standard	4"x14" Knobby, Flat Free Black Tires	\$580
Part: PTOINDV3525		Part: XRF000368	
3"x14" Pneumatic Black Tires with Accent	No Charge	4"x14" Knobby, Pneumatic Black Tires	\$580
Part: WHL1704059 x2		Part: XRF000367	

5.2. Drive Wheel Accent Color

Not available with knobby tires.


Back in Black	Standard	Pretty in Pink	No Charge	Green Machine	No Charge
Part: DGN164245-PBK1 x2		Part: DGN164245-PPK3 x2		Part: DGN164245-PGN11 x2	
Red Corvette	No Charge	Evolution Orange	No Charge	Military Green	No Charge
Part: DGN164245-PRD1 x2		Part: DGN164245-POR7 x2		Part: DGN164245-PGN8 x2	
Night Train	No Charge	Midnight Blue	No Charge	Military Digital Camo	\$58 ea.
Part: DGN164245-PBK3 x2		Part: DGN164245-PBU10 x2		Part: DGN164245-WCM7 x2	
Grey Street	No Charge	Electric Blue	No Charge	Faded Glory	\$58 ea.
Part: DGN164245-PSL3 x2		Part: DGN164245-PBU9 x2		Part: DGN164245-WFL2 x2	
Deep Purple	No Charge	Blue Jean	No Charge		
Part: DGN164245-PPR5 x2		Part: DGN164245-PBU11 x2			
Raspberry Beret	No Charge	Lemon Crush	No Charge		
Part: DGN164245-PPK4 x2		Part: DGN164245-PYL4 x2			

6. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Operated Fender Lights	No Charge	Omit Fender Lights	No Charge
Part: DGN2004686		Part: DGN1700032	

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.



7. BATTERIES

Optional. The powerchair requires 2 flat top batteries and will be installed if ordered.

Group 24, AGM Cell Battery (MK)..... \$541 ea.
Part: BAT1712207. HCPCS: E2363

Group 34, Gel Cell Battery (MK) \$509 ea.
Part: BATGEL1006. HCPCS: E2359
Not available with Power Lift.

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight
Part: PTOINDV3837

201 lb-300 lb User Weight
Part: PTOINDV3838

9. JOYSTICK

REQUIRED. Select a type of drive control from the following sections.

9.1a. Q-Logic 3 Non-Expandable Capacity Joystick

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.

QLNES 10 Key Joystick \$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

9.1b. Q-Logic Expandable Capacity Controls

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity Joysticks

Q-Logic 3 Expandable Capacity Joystick\$840
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick\$840
for Full Light Package
Part: CTL167725
Required with full light packages.

Q-Logic 3e Joystick\$735
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls. Includes magnetic charger.



2. Q-Logic 3 Expandable Capacity Required Components

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$730
Part: XRF000002 HCPCS: E2377

Q-Logic 3 Expandable Capacity Controller Harness..... \$458
Part: XRF000264. HCPCS: E2313

9.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635

9.3. Q-Logic Joystick Armor

REQUIRED with Q-Logic 3 Expandable Joystick selections above. Not available with Q-Logic 3e and QLNES Joystick. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds.

Red Corvette..... \$63
Part: DGN124336-PRD1

Back in Black.....\$63
Part: DGN124336-PBK1

Grey Street..... \$63
Part: DGN124336-PSL3

Deep Purple.....\$63
Part: DGN124336-PPR5

Raspberry Beret..... \$63
Part: DGN124336-PPK4

Pretty in Pink.....\$63
Part: DGN124336-PPK3

Evolution Orange.....\$63
Part: DGN124336-POR7

Midnight Blue.....\$63
Part: DGN124336-PBU10

Electric Blue \$63
Part: DGN124336-PBU9

Blue Jean \$63
Part: DGN124336-PBU11

Lemon Crush.....\$63
Part: DGN124336-PYL4

Green Machine\$63
Part: DGN124336-PGN11

Military Green \$63
Part: DGN124336-PGN8

10. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in the image on the previous page. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C)..... \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F)..... \$80

Part: ACCASMB1856. HCPCS: E2323


Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323


Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



11. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket Standard
Q-Logic Part: FRMASMB10617Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1028**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201043. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201044. HCPCS: E1032**11.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.**11.2. Joystick Mount Location**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**12. CONTOURED SEATING****REQUIRED.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.**12.1. Occupied Transit, WC-19**WC-19 Occupied Transit Option for Contoured Seat \$368
Part: KIT122609To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**12.2a. Q-Captain Seat**16"W x 16"D No Charge
Part: SET200514518"W x 20"D No Charge
Part: SET200515120"W x 22"D No Charge
Part: SET200590916"W x 18"D No Charge
Part: SET200514618"W x 22"D No Charge
Part: SET200590822"W x 18"D No Charge
Part: SET200515516"W x 20"D No Charge
Part: SET200514720"W x 16"D No Charge
Part: SET200515222"W x 20"D No Charge
Part: SET200515618"W x 16"D No Charge
Part: SET200514920"W x 18"D No Charge
Part: SET200515322"W x 22"D No Charge
Part: SET200515718"W x 18"D No Charge
Part: SET200515020"W x 20"D No Charge
Part: SET2005154**12.2b. Q-Captain Seat with Solid Seat Pan**16"W x 16"D No Charge
Part: SET200539418"W x 20"D No Charge
Part: SET200539920"W x 22"D No Charge
Part: SET200591116"W x 18"D No Charge
Part: SET200539518"W x 22"D No Charge
Part: SET200591022"W x 18"D No Charge
Part: SET200540316"W x 20"D No Charge
Part: SET200539620"W x 16"D No Charge
Part: SET200540022"W x 20"D No Charge
Part: SET200540418"W x 16"D No Charge
Part: SET200539720"W x 18"D No Charge
Part: SET200540122"W x 22"D No Charge
Part: SET200540518"W x 18"D No Charge
Part: SET200539820"W x 20"D No Charge
Part: SET2005402**12.3. Seat to Floor Height**

Power adjustable seat height option is fixed at lowest seat to floor height.


Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

	Cushion	SSP
High	22.5"	18.5"
Med	21.5"	17.5"
Q-Captain Seat to Floor Heights		

Contoured seating section continued on following pages.

12. CONTOURED SEATING CONTINUED**12.4. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39

Part: POP1909006-S

14"W and smaller

Small Wedge Positioning Kit\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit.....\$105

Part: POP1909006-L

14"W and smaller

Large Wedge Positioning Kit.....\$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

FRG Top Foam

Reference part: XRF000381

No Top Foam

Reference part: XRF000382

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

Stealth Polyurethane Liner

Reference part: XRF000384

No Liner

Reference part: XRF000385

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Mesh Cover

Reference part: XRF000387

Stealth Tek Cover

Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

Contoured seating section continued on following pages.

9. POWER POSITIONING SEATING CONTINUED**6. Optional Base Insert**

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

QTY

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

QTY

Additional Stealth Tek Cover

Reference Part: XRF000394

QTY

Additional Mesh Cover

Reference Part: XRF000393

QTY

12.5. Armrests**Optional.** Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Q Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat arm pads.

Full Length, Left Armrest \$90 ea.

Part: MEC2201304. HCPCS: K0020

Full Length, Right Armrest \$90 ea.

Part: MEC2201303. HCPCS: K0020

1b. Drop-In, Height Adjustable, Armrests (T-Arm)

Includes flat arm pads.

☐ ☐ Full Length, Drop-In Armrests \$79 ea.
L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973☐ ☐ Desk Length, Drop-In Armrests \$79 ea.
L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973**1c. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973☐ ☐ Desk Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973**1d. Bariatric, Height Adjustable, Removable Armrests**

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H \$90 ea.
L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973☐ ☐ Full Length, Bariatric Armrests, 10"-14"H \$90 ea.
L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973**1e. Heavy Duty, 8"-12" Height Adjustable, Removable Armrests**

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.
L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H \$90 ea.
L R Left Part: XRF000255. Right part: XRF000256. HCPCS: E0973**1e. Heavy Duty, 10"-14" Height Adjustable, Removable Armrests**

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
L R Left Part: XRF000257. Right part: XRF000258. HCPCS: E0973☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14" H \$90 ea.
L R Left part: XRF000259. Right part: XRF000260. HCPCS: E0973**1f. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight arm pads.

☐ ☐ Pediatric, Removable Armrests \$79 ea.
L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973**12.6. Q-Captain Armrest Spacer****Optional.**

Armrest Spacer \$35

Part: XRF000370

1" spacer for Q-Captain, Flip-Up, Height Adjustable Armrests only.

12.7. Armrest Height**REQUIRED with Bariatric, Heavy Duty, and Pediatric armrests.****2. Left Armrest Height**

7.75" H

PTOINDV3722

8.25" H

PTOINDV3724

10" H

PTOINDV3537

12" H

PTOINDV3539

14" H

PTOINDV3541

16" H

PTOINDV3543

8" H

PTOINDV3535

9" H

PTOINDV3536

11" H

PTOINDV3538

13" H

PTOINDV3540

15" H

PTOINDV3542

3. Right Armrest Height

7.75" H

PTOINDV3721

8.25" H

PTOINDV3723

10" H

PTOINDV3452

12" H

PTOINDV3454

14" H

PTOINDV3456

16" H

PTOINDV3458

8" H

PTOINDV3450

9" H

PTOINDV3451

11" H

PTOINDV3453

13" H

PTOINDV3455

15" H

PTOINDV3457

Contoured seating section continued on following pages.

12. CONTOURED SEATING CONTINUED**12.8. Armpads****REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.****1a. Straight Armpads for Single Post Armrests**Straight, Full Length Left and Right No Charge
Part: ARM134275Straight, Full Length Left, Desk Length Right No Charge
Part: ARM134281Straight, Desk Length Left, Full Length Right No Charge
Part: ARM134280Straight, Desk Length Left and Right No Charge
Part: ARM134254**1b. Waterfall Armpads for Single Post Armrests**Waterfall, Full Length Left and Right No Charge
Part: ARM134285Waterfall, Full Length Left, Desk Length Right No Charge
Part: ARM134289Waterfall, Desk Length Left, Full Length Right No Charge
Part: ARM134288Waterfall, Desk Length Left and Right No Charge
Part: ARM134282**1c. Falcon Gel Armpads for Single Post Armrests**Falcon Gel, Full Length Left and Right \$311.22
Part: ARM136438Falcon Gel, Full Length Left, Desk Length Right \$287.93
Part: ARM136442Falcon Gel, Desk Length Left, Full Length Right \$287.93
Part: ARM136439Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**1d. Armpads for Bariatric or Heavy Duty Armrests**Bariatric Armpads No Charge
Part: UPSASMB1027
Only for use with Bariatric Armrests.Straight Armpads No Charge
Part: UPSASMB1021
Only for use with Heavy Duty ArmrestsWaterfall Armpads No Charge
Part: UPSASMB1024
Only for use with Heavy Duty Armrests**12.9. Quad Release Lever****Optional. Available with Single Post and Pediatric armrests.**Quad Release Lever, Pair No Charge
Part: FRMASMB10517**13. LAP BELTS****REQUIRED except with depth adjustable high back. Select a lap belt below.**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**14. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**14.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform TypeBasic Center Mount No Charge
Part: XRF000323**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates No Charge
Part: MEC1906675Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906680Long Extension, 5"x9" Footplates No Charge
Part: MEC1906677Medium Lower Extension, 6.25"W Tapered 5"W x 12"D... No Charge
Part: MEC2104853Long Lower Extension, 6.25"W Tapered 5"W x 12"D..... No Charge
Part: MEC2104854**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Length Extension.**4. Basic Center Mount Additional Accessories**Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following page.

14. LEGRESTS CONTINUED**14.1b. Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform Assembly

Quick Adjustable Center Mount Foot Platform\$263
 Part: XRF000311. HCPCS: K0040

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Pediatric Extension, 5"x8" Footplates
 Part: RIG2404015

Medium Extension, 7"x8" Footplates
 Part: RIG2404027

Long Extension, 7"x8" Footplates
 Part: RIG2404028

Medium Extension, 5"x6" Footplates
 Part: RIG150206

Medium Extension, 7"x11" Footplates
 Part: RIG2404029

Long Extension, 7"x11" Footplates
 Part: RIG2404030

3a. Center Mount Calf Pads

Small Calf Pads, Quick Adjustable Center Mount Only, \$158
 6"Wx6"H (pair)
 Part: XRF000156

Large Calf Pads, Quick Adjustable Center Mount Only, \$158
 8" taper to 6"Wx8"H (pair)
 Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel for Quick Adjustable, 18"Wx8"H \$184
 Part: FRMASMB7110012

Large Calf Panel for Quick Adjustable, 24"Wx8"H \$237
 Part: RIG123825
 Only available on widths of 20" or greater.

4. Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
 Part: RIG152519
 Only available with the 5x6 footplates.

Outside Edge \$69
 Part: RIG2404922
 Available with 7" width footplates.

Width Extension for 8" \$69
 Part: RIG2404919

Front and Rear Cap \$16
 Part: RIG2404923

Width Extension for 11" \$69
 Part: RIG2404920

Available with 5"x8" footplates.

Outside Edge \$69
 Part: RIG2404921

Available with 5"x8" footplates.

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
 Part: ST-CSCS0775-Q

Footplate Riser \$75
 Part: XRF000348

Large Foam Calf Pad Cover (pair) \$83
 Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Small Foam Panel Cover (pair) \$176
 Part: ST-CSCP0920-Q

Footplate Riser Calf Panel \$50
 Part: FRMASMB7110011

Large Foam Panel Cover (pair) \$282
 Part: ST-CSCP0927-Q

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$16
 Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-AFPG5008E-Q. HCPCS: K0108

7"x11" Tapered Gel Footplate Covers \$270
 Part: ST-AFPG7008T-Q. HCPCS: K0108

7"x11" Tapered Gel Footplate Cover for Width Extension \$293
 Part: ST-AFPG7011T-Q. HCPCS: K0108

14.1c. 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Swing-Away Legrest Upper

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest \$79 ea.
 Pair part: XRF000082. Left: XRF000083.
 Right: XRF000084
 Not available with 11"-13" extensions.

☐ P ☐ L ☐ R 60°-90° Swing-Away Legrest, Short \$79 ea.
 Pair part: XRF000085. Left: XRF000086.
 Right: XRF000087.
 For use only with 11"-13" lower extensions.

2. Lower Extension

☐ P ☐ L ☐ R 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ P ☐ L ☐ R 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ P ☐ L ☐ R 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

Legrests section continued on following page.

14. LEGRESTS CONTINUED**3b. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ P ☐ L ☐ R	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ P ☐ L ☐ R	6"Wx11"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ P ☐ L ☐ R	6"Wx8"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ P ☐ L ☐ R	8"Wx11"D ² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea. Part: INDPART2822. HCPCS: E1020	Therafin 32387 Amputee Support Gel Cover \$155 ea. Part: INDPART2823. HCPCS: K0108
--	---

14.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ P ☐ L ☐ R	Elevating Legrest Assembly \$171 ea. Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990 Select lower extension and footplates
--	--

2. Lower Extension

☐ P ☐ L ☐ R	11"-14" Length Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced only	☐ P ☐ L ☐ R	14"-17" Length Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced only
☐ P ☐ L ☐ R	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103 For angle adjustable only	☐ P ☐ L ☐ R	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only

3a. Reinforced Cast Footplate(s)

☐ P ☐ L ☐ R	7.25"Wx6"D Reinforced Footplate..... Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	☐ P ☐ L ☐ R	8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
--	--	--	--

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ P ☐ L ☐ R	5"Wx8"D with Heel Loops..... \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	☐ P ☐ L ☐ R	6"Wx11"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
☐ P ☐ L ☐ R	6"Wx8"D ¹ with Heel Loops..... \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951	☐ P ☐ L ☐ R	8"Wx11"D ² with Heel Loops..... \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25 Part: INDPART2242. HCPCS: E0953
--

14.2. Legrest Accessories

Heel Loops (pair)..... \$53 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi axis footplates when noted.	Wheel Bumper..... \$37 Part: XRF000100
--	---

14.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

Account Number: _____ Date: _____ Quote Name/Marked for: _____

15. POWER LIFT ELECTRONICS

REQUIRED with 12" Power Lift for Captain Seat with iLevel.

Single Actuator Control through Q-Logic Joystick..... 2,060
Part: ELE2009405. HCPCS: E2310

Single Actuator Through Controller Interface (AAM)\$3,040
Part: XRF000068. HCPCS: E2310

16. SEAT MOUNTED LED LIGHTS

Optional. Requires QLNES 10 Key Joystick, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

17. ACCESSORIES

Optional.

Cup Holder, Q-Captain No Charge
Part: ACC2006355

Front Mounted Cell Phone/ Media Holder, X-Grip, Q-Captain..... \$74
Part: ACC2103691

Rear View Mirror, Q-Captain \$69
Part: ACC2103692

Oxygen Holder.....\$118
Part: ACC2006353. HCPCS: E2208
Not available with Power Lift option.

Oxygen Holder..... \$180
Part: ACC2400456. HCPCS: E2208
Select with Power Lift option.

Cane/Crutch Holder\$110
Part: ACC2006352. HCPCS: E2207
Not available with Power Lift option.

Cane/Crutch Holder \$250
Part: ACC2400455. HCPCS: E2207
Select with Power Lift option.

Rear Basket\$118
Part: ACC2006351
Not available with Power Lift option.

Rear Basket\$124
Part: ACC2400454
Select with Power Lift option.

Saddle Bag.....\$40
Part: DWR1020D010

Large, Black, Weather Cover \$202
Part: ACCCOVR1021

Quantum Backpack \$168
Part: ACC1905853

Grip Lap Board (20005)..... \$38.95
Part: ACC152992

Grip Drink Holder (20085) \$15.50
Part: ACC152993

Reflector Kit \$53
Part: ACC124670

Clothing Guard (Pair)..... \$237
Part: ACC125007
Available with Flip-Up armrests only. Not available with Elevating Legrests.

Transfer Handles\$105 ea.
Part: MEC2205680. HCPCS: E0705
Available with Flip-Up armrests only. Not available with Elevating Legrests.

Privacy Bag.....\$84
Part: ACC2005723

Personal Item Pouch, Rail Mounted.....\$84
Part: ACC2101017

Personal Item Pouch, Arm Mounted\$84
Part: ACC2103694

Glove Box, Rail Mounted \$158
Part: ACC125011

Glove Box, Arm Mounted..... \$158
Part: ACC2103693

Remote Stop Switch \$1,295
Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

18. PROGRAMMER

Optional.

Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976

Q-Logic PC Programming Disk with Cable \$445
Part: ELEASMB5216

19. BRANCH PATCHES

Optional.

Air Force Patch \$21
Part: ACC140240

Army Patch \$21
Part: ACC140241

Coast Guard Patch \$21
Part: ACC140239

Marines Patch..... \$21
Part: ACC140243

Navy Patch \$21
Part: ACC140244

Purple Heart Patch \$21
Part: ACC140238



Military Patches

20. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**20.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

_____ QTY UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware <i>Part: ST-UTTP5-Q. HCPCS: E1028/E1034^</i>	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) <i>Part: POS140428. HCPCS: E1028/E1034^</i>
_____ QTY UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware <i>Part: ST-UTTP7-Q. HCPCS: E1028/E1034^</i>	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) <i>Part: POS140429. HCPCS: E1028/E1034^</i>

2. UniLink Pelvic/Thigh Pads

_____ QTY UniLink, 4" x 4" Pad \$106 ea. <i>Part: ST-UTP544056100-Q. HCPCS: E0953/E0956</i>	_____ QTY UniLink, 4" x 12" Pad \$106 ea. <i>Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956</i>
_____ QTY UniLink, 4" x 8" Pad..... \$106 ea. <i>Part: ST-UTP548056100-Q. HCPCS: E0953/E0956</i>	

3. UniLink Pelvic/Thigh Additional Links

_____ QTY UniLink, 1.5" Link \$81 ea. <i>Part: ST-UIL001150A000-Q</i>	_____ QTY UniLink, 2" Link..... \$81 ea. <i>Part: ST-UIL001200A000-Q</i>
---	--

4. UniLink Pelvic/Thigh Guide Installation

Install on Left	Install on Right	Do Not Install
<i>Part: PTOINDV3585</i>	<i>Part: PTOINDV3586</i>	<i>Part: PTOINDV3719</i>

20.2. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube\$199 ea. Style Arms (ARMS-ES-L) <i>Part: POS124064. HCPCS: E1028</i>	Right, Swing-Away, Elbow Support Bracket, Tube.....\$199 ea. Style Arms (ARMS-ES-R) <i>Part: POS124066. HCPCS: E1028</i>
--	--

2. Elbow Support Pad

_____ QTY 3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea. <i>Part: POS124065. HCPCS: K0108</i>	_____ QTY 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea. <i>Part: POS124067. HCPCS: K0108</i>
---	--

20.3. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539 (SPLE107) (pair) <i>Part: INDPART2786. HCPCS: E0954</i>	Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539 (SPLE109) (pair) <i>Part: INDPART2788. HCPCS: E0954</i>
Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539 (SPLE108) (pair) <i>Part: INDPART2787. HCPCS: E0954</i>	

20.4. Stealth Positioning Shoe Holder

_____ QTY Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea. <i>Part: ST-IFSH4211-Q. HCPCS: E0951+E0952</i>	_____ QTY Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea. <i>Part: ST-IFSH4213-Q. HCPCS: E0951+E0952</i>
_____ QTY Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea. <i>Part: ST-IFSH4212-Q. HCPCS: E0951+E0952</i>	

21. LAP TRAY

21.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	M	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.