



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum® 4Front® 2 Order Form with TRU Balance® 4 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

QUANTUM
4FRONT²



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 3SP-SS\$8,559

Part: 4FRONT 2 3SP-SS. HCPCS: K0856

Single power option.

Quantum 4Front 2 3MP-SS\$8,373

Part: 4FRONT 2 3MP-SS. HCPCS: K0861


Multiple actuator power option.

2. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**2.1a. Standard Shrouds**

 Red Corvette Part: DGN2003028-PRD1	 Deep Purple Part: DGN2003028-PPR5	 Midnight Blue Part: DGN2003028-PBU10	 Green Machine Part: DGN2003028-PGN11
 Back in Black Part: DGN2003028-PBK1	 Raspberry Beret Part: DGN2003028-PPK4	 Electric Blue Part: DGN2003028-PBU9	 Military Green Part: DGN2003028-PGN8
 Night Train Part: DGN2003028-PBK3	 Pretty in Pink Part: DGN2003028-PPK3	 Blue Jean Part: DGN2003028-PBU11	
 Grey Street Part: DGN2003028-PSL3	 Evolution Orange Part: DGN2003028-POR7	 Lemon Crush Part: DGN2003028-PYL4	

2.1b. Hydrographic Shrouds**Military Digital Camo** \$263
Part: DGN2003028-WCM7**Faded Glory** \$263
Part: DGN2003028-WFL2**3. BASE MOUNTED LIGHTS****Optional.** Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.**Switch Operated Fender Lights** No Charge
Part: DGN2004686**Omit Fender Lights** No Charge
Part: DGN1700032

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.



Fender Lights

4. QUANTUM BACKUP CAMERA**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.**Quantum Backup Camera** \$368
Part: MEC1910365**1. Optional Mounting Bracket****Backup Camera Headrest Bracket** \$21
Part: MEC1910373**Backup Camera Track Mount Bracket** \$27
Part: MEC2104227

Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket**Backup Camera Q-Logic Joystick Mount Bracket** \$105
Part: MEC2106123

Backup Camera

5. DRIVE WHEELS**REQUIRED.** Select a type of drive wheels and accent color.**5.1. Drive Wheels****3"x14" Flat Free Black Tires with Accent** Standard
Part: PTOINDV3525**4"x14" Knobby, Flat Free Black Tires** \$580
Part: XRF000368**3"x14" Pneumatic Black Tires with Accent** No Charge
Part: WHL1704059 x2**4"x14" Knobby, Pneumatic Black Tires** \$580
Part: XRF000367**5.2. Drive Wheel Accent Color**

Not available with knobby tires.

Back in Black Standard
Part: DGN164245-PBK1 x2**Raspberry Beret** No Charge
Part: DGN164245-PPK4 x2**Green Machine** No Charge
Part: DGN164245-PGN11 x2**Red Corvette** No Charge
Part: DGN164245-PRD1 x2**Pretty in Pink** No Charge
Part: DGN164245-PPK3 x2**Military Green** No Charge
Part: DGN164245-PGN8 x2**Cherry Bomb** No Charge
Part: DGN164245-PRD4 x2**Evolution Orange** No Charge
Part: DGN164245-POR7 x2**Sugar Plum** No Charge
Part: DGN164245-PPR7 x2**Night Train** No Charge
Part: DGN164245-PBK3 x2**Midnight Blue** No Charge
Part: DGN164245-PBU10 x2**Military Digital Camo** \$58 ea.
Part: DGN164245-WCM7 x2**Grey Street** No Charge
Part: DGN164245-PSL3 x2**Electric Blue** No Charge
Part: DGN164245-PBU9 x2**Faded Glory** \$58 ea.
Part: DGN164245-WFL2 x2**White Satin** No Charge
Part: DGN164245-PWH2 x2**Blue Jean** No Charge
Part: DGN164245-PBU11 x2**Deep Purple** No Charge
Part: DGN164245-PPR5 x2**Lemon Crush** No Charge
Part: DGN164245-PYL4 x2**6. BATTERIES****Optional.** The powerchair requires 2 flat top batteries and will be installed if ordered.**Group 24, AGM Battery (MK)** \$541 ea.
Part: BAT1712207 x2. HCPCS: E2363

7. PATIENT'S WEIGHT**REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 200 lb User Weight
Part: PTOINDV3837201 lb-300 lb User Weight
Part: PTOINDV3838**8. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**8.1. Q-Logic 3 Expandable Capacity Controls**

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity JoysticksQ-Logic 3 Expandable Capacity Joystick\$840
Part: CTL167724Q-Logic 3 Expandable Capacity Joystick\$840
for Full Light Package
Part: CTL167725
Required with full light packages.Q-Logic 3e Joystick\$735
Part: CTL1910600


Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic
Part: XRF000010**2. Q-Logic 3 Expandable Capacity Required Components**

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Capacity Controller Harness..... \$458
Part: XRF000264. HCPCS: E2313**8.2. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.Black, Soft TouchStandard
Part: DGN124307Red Corvette..... \$63
Part: DGN124336-PRD1Back in Black..... \$63
Part: DGN124336-PBK1Grey Street..... \$63
Part: DGN124336-PSL3Deep Purple..... \$63
Part: DGN124336-PPR5Raspberry Beret..... \$63
Part: DGN124336-PPK4Pretty in Pink \$63
Part: DGN124336-PPK3Evolution Orange..... \$63
Part: DGN124336-POR7Midnight Blue \$63
Part: DGN124336-PBU10Electric Blue \$63
Part: DGN124336-PBU9Blue Jean \$63
Part: DGN124336-PBU11Lemon Crush \$63
Part: DGN124336-PYL4Green Machine \$63
Part: DGN124336-PGN11Military Green \$63
Part: DGN124336-PGN8Military Digital Camo..... \$63
Part: DGN124336-WCM7Faded Glory \$63
Part: DGN124336-WFL2**9. JOYSTICK HANDLES****Optional.** The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308**10. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**10.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket Standard
Part: FRMASMB10617**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**10.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032


Joystick mounting brackets section continued on following page.

10. JOYSTICK MOUNTING BRACKETS CONTINUED**10.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
 Part: XRF000003. HCPCS: E1028
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.



Retract4

10.2. Joystick Mount Location

☐ **Left Side Mount**
 Part: PTOINDV1225

☐ **Right Side Mount**
 Part: PTOINDV1226



Multi-Axis Joystick Mount

11. TRU BALANCE 4 SEATING

Contains memory seating and positional feedback. Maximum seat size is 22"W x 22"D.

11.1. TRU Balance 4 Configuration**REQUIRED.**

TRU-Balance 4 Power Tilt Only \$8,595
 Part: PTOINDV4014. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 4 Power Tilt and Recline \$12,790
 Part: PTOINDV3954. HCPCS: E1007
 50° of posterior tilt and 165° of recline.

11.2. 10" Power Lift

Requires combined selection of TRU-Balance 4 Power Tilt and Recline and Articulating Foot Platform. Requires selection in 12.3 Anterior Tilt. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

Power Lift \$3,725
 Part: PTOINDV3973. HCPCS: E2298/E2300^
 10" scissor lift with positional feedback.

11.3. Anterior Tilt

Only available with Power Lift. Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

1. Anterior Tilt Angle

Only available with Power Lift.

10° Anterior Tilt Angle No Charge
 Part: MEC2009356

30° Anterior Tilt Angle \$2,595
 Part: MEC2009400. HCPCS: K0108
 Requires chest support and knee block.

20° Anterior Tilt Angle \$1,995
 Part: MEC2009399. HCPCS: K0108
 Requires chest support and knee block.

0° Anterior Tilt No Charge
 Part: PTOINDV3980

11.4. Knee Block Assembly**REQUIRED with 20° or 30° anterior tilt. Optional with 10° anterior tilt.**

Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$314
 Part: ST-KBATB4ATC-1-Q. HCPCS: K0108

Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) .. \$354
 Part: ST-KBATB4ATC-U-Q. HCPCS: K0108

11.5. Occupied Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option \$289
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

11.6. Seat Width**REQUIRED with seat selection. 12"-15" widths not available with Power Lift. Maximum seat width is 22". Part number is SETWDTH plus number following # below.**

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.7. Seat Width Package**REQUIRED.**

12"-16" Width Package
 Part: XRF000006
 Available with Power Tilt Only.

16"-20" Width Package
 Part: XRF000008

20"-24" Width Package
 Part: XRF000158
 Maximum width is 22".

11.8. Seat Depth**REQUIRED with seat selection. 12"-13" depth not available with Power Lift. Maximum seat depth is 22". Part number is SETDPHT plus number following # below.**

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU Balance 4 Seating section continued on following pages.

11. TRU BALANCE 4 SEATING CONTINUED**11.9. Seat Depth Package**

REQUIRED with seat selection.

1a. Power Tilt Depth Package12"-16" Depth Package
Part: XRF00001116"-20" Depth Package
Part: XRF00001320"-24" Depth Package
Part: XRF000159
Maximum depth is 22".**1b. Power Lift Depth Package**16"-22" Depth Package
Part: XRF000309


Low	Med	High
Without Lift		
16.75"	17.25"	17.75"
With Lift		
17.75"	18.25"	18.75"
TRU Balance 4 Seat to Floor Heights		

11.10. Seat to Floor Height

REQUIRED.

Low Seat to Floor
Part: PTOINDV3976Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560**11.9. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit\$72
Part: POP1909006-M
14"W and smaller

Medium Wedge Positioning Kit\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit\$105
Part: POP1909006-L
14"W and smaller

Large Wedge Positioning Kit\$105
Part: POP1712340-L
15"W and larger

TRU Balance 4 Seating section continued on following pages.

11. TRU BALANCE 4 SEATING CONTINUED**3. Cushion Top Foam**Visco Top Foam
Reference part: XRF000380No Top Foam
Reference part: XRF000382FRG Top Foam
Reference part: XRF000381**4. Cushion Liner**Stealth Urethane Liner
Reference part: XRF000383No Liner
Reference part: XRF000385Stealth Polyurethane Liner
Reference part: XRF000384**5. Cushion Cover**CoolCore Cover
Reference part: XRF000386Stealth Tek Cover
Reference part: XRF000388Mesh Cover
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					Cover
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69**7. Additional Cushion Liner**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____**8. Additional Cushion Cover**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____Additional Mesh Cover
Reference Part: XRF000393
QTY _____**11.12. Back Type Selection****REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.Redesigned, Static Type Back
Part: XRF000014

Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
Select cane type, back cane height, and back cushion configuration below.Redesigned, Flip Forward Back \$184
Part: XRF000016

Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

11.13. Back Cane Selection (Back Cane Option Only)**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**Angle Adjustable Back Cane No Charge
Part: MEC1402662" Endomorph (1" per side) Back Cane \$74
Part: XRF000269Flip Forward Back Cane \$53
Part: XRF0002734" Endomorph (2" per side) Back Cane \$74
Part: XRF000271**2. Back Cane Height**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**11.14. Back Cushion Configuration****REQUIRED.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**Sport Back No Charge
Part: XRF000017
Complete following width package and size.After-market Back Interface \$535
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.Omit Back Pan/Cushion No Charge
Part: PTOINDV1287
Select a following width package.

TRU Balance 4 seating section continued on following page.

11. TRU BALANCE 4 SEATING CONTINUED**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), 12" to 20" Width..... \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width..... \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$568
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$568
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$568
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$643
 Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$643
 Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
 Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
 Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$643
 Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$643
 Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
 Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$643
 Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$643
 Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$643
 Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$643
 Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$703
 Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
 Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
 Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$703
 Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$703
 Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$703
 Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$703
 Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$703
 Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$703
 Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$703
 Part: QB2AL20EX. HCPCS: E2620

Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL20EXX. HCPCS: E2620

11.15. TRU Balance 4 Back Shroud

REQUIRED. The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette
 Part: XRF000019

Grey Street
 Part: XRF000023

Pretty in Pink
 Part: XRF000027

Electric Blue
 Part: XRF000030

Green Machine
 Part: XRF000033

Back in Black
 Part: XRF000021

Deep Purple
 Part: XRF000025

Evolution Orange
 Part: XRF000028

Blue Jean
 Part: XRF000031

Military Green
 Part: XRF000034

Night Train
 Part: XRF000022

Raspberry Beret
 Part: XRF000026

Midnight Blue
 Part: XRF000029

Lemon Crush
 Part: XRF000032

11.16. Lateral Positioning Wedge

Optional. Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge \$16 ea.
 Part: SET123465 x2

TRU Balance 4 seating section continued on following page.

11. TRU BALANCE 4 SEATING CONTINUED**11.17. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003716"W-20"W Package
Part: XRF00003818"W-22"W Package
Part: XRF000039**11.18. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.12" Back Width
Part: PTOINDV356215" Back Width
Part: PTOINDV385617" Back Width
Part: PTOINDV356519" Back Width
Part: PTOINDV356721" Back Width
Part: PTOINDV356914" Back Width
Part: PTOINDV356316" Back Width
Part: PTOINDV356418" Back Width
Part: PTOINDV356620" Back Width
Part: PTOINDV356822" Back Width
Part: PTOINDV3570**11.19. Overall Back Height****Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.**19" Height
Part: SETHGHT103222" Height
Part: SETHGHT103525" Height
Part: SETHGHT103828" Height
Part: SETHGHT104131" Height
Part: SETHGHT110220" Height
Part: SETHGHT103323" Height
Part: SETHGHT103626" Height
Part: SETHGHT103929" Height
Part: SETHGHT104232" Height
Part: SETHGHT110321" Height
Part: SETHGHT103424" Height
Part: SETHGHT103727" Height
Part: SETHGHT104030" Height
Part: SETHGHT1043

		Tilt Only Width											Recline Width										
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		
Overall Height	14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"	
	15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"	
	16"	S	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	-	16"	
	17"	S	T	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	-	17"	
	18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	-	-	-	-	-	18"	
	19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	19"	
	20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	20"	
	21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	21"	
	22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	22"	
	23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	23"	
	24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	24"	
	25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	25"	
	26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	26"	
	27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	27"	
	28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	28"	
	29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	29"	
	30"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	30"	
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	31"		
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	32"		

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

11.20. Back Height Extension**Optional.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.Back Height Extension \$48
Part: XRF000041**12. ARMRESTS****REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**12.1. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and armrest height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

L

R

Flip-Back, 8"-16"H Armrests Standard
 Left part: ARM1806662. Right part: ARM1806663

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

L

R

Flip-Back Armrest Downpost \$92 ea.
 Left part: XRF000265. Right part: XRF000266

12.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

L

R

Full Length, Two-Post Armrest for Power Seating No Charge
 Left part: XRF000047. Right part: XRF000048

L

R

Desk Length, Two-Post Armrest for Power Seating.... No Charge
 Left part: XRF000049. Right part: XRF000050

Armrest section continued on following page.

12. ARMRESTS CONTINUED**12.1c. Pediatric, Removable, Height Adjustable Armrests**Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

L	R	Pediatric, Removable Armrests \$79 ea.
		Left Part: ARM142876. Right part: ARM142877

12.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1a. Flip-Back and Two-Post Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

12.3. Armrest Size

Select a left (L) and right (R) armrest size. This also affects armpad selection.

L	R	10" Length Left part: XRF000057. Right part: XRF000123	L	R	14" Length Left part: XRF000058. Right part: XRF000124	L	R	17" Length Left part: XRF000059. Right part: XRF000125
----------	----------	--	----------	----------	--	----------	----------	--

13. ARMPADS**REQUIRED.****13.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge	Interface for Aftermarket Pad No Charge
Part: XRF000060	Part: XRF000308
Not available with Flip-back, 8"-16" H armrests.	Only available with Flip-back, 8"-16" armrests.

13.1b. Foam Armpads

L	R	Straight, Armpad Standard	L	R	Waterfall, Armpad No Charge
		Left part: XRF000061. Right part: XRF000126			Left part: XRF000062. Right part: XRF000127
					Not available with Flip-back, 8"-16" armrests.

13.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195	Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806807	Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171	Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806813	Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22	Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136438	Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93	Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136442	Part: ARM136437

14. JOYSTICK BRACKETS**REQUIRED.**

Standard Fixed Receiver Standard	Stealth Height Adjustable Receiver..... \$236
Part: XRF000063	Part: XRF000064. HCPCS: K0108

15. HEADRESTS**Optional.** Select a headrest pad and bracket.**15.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q, HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q, HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q, HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q, HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q, HCPCS: E1028**15.2. Q Cam Lock****Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

16. LAP BELTS**REQUIRED.** Select a lap belt below.**16.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**16.1b. Optional Lap Belts**Stealth Positioning 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning 60" Single Pull, Push Button, Extra Large.....\$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable).....\$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.**17. LEGRESTS**

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

17.1a. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable FootplatesShort Extension, 5"x9" Footplates\$180
Part: XRF000312, HCPCS: K0040Medium Extension, 5"x9" Footplates\$180
Part: XRF000313, HCPCS: K0040Long Extension, 5"x9" Footplates\$180
Part: XRF000314, HCPCS: K0040**3. Basic Center Mount Calf Pads**Small Calf Pads, Basic Center Mount.....\$158
Part: RIG2104369Large Calf Pads, Basic Center Mount\$158
Part: RIG2202997

Only available with Long Extension.

4. Basic Center Mount Additional AccessoryBasic Center Mount Footplate Connector\$16
Part: RIG2104240**17.1b. Center Mount Foot Platform****1. Foot Platform Assembly**Dual Actuator, Power Articulating Foot Platform.....\$3,595
Part: RIG2009246, HCPCS: E1012+K0040
Positional feedback with independent control of articulation.Quick Adjustable Center Mount Foot Platform.....\$263
Part: XRF000311, HCPCS: K0040
Not available with power lift.

Legrest Compatibility Articulating Foot Platform TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Medium Length Extension	11.5"-14.5"	11.5"-15.5"	11.5"-16.5"
Long Length Extension	N/A	14.5"-15.5"	14.5"-16.5"

Legrest Compatibility Quick Adjustable CM TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Pediatric Length Extension	8.5"-12.5"	8.5"-12.5"	8.5"-12.5"
Medium Length Extension	11"-15"	11"-16"	11"-17"
Long Length Extension	14"-15"	14"-16"	14"-17"

Legrests section continued on following page.

17. LEGRESTS CONTINUED**2. Lower Extension and Angle Adjustable Footplates**

Pediatric Extension, 5"x8" Footplates
Part: RIG2404015
Not available with power articulating platform.
Medium Extension, 5"x6" Footplates
Part: RIG150206

Medium Extension, 7"x8" Footplates
Part: RIG2404027
Medium Extension, 7"x11" Footplates
Part: RIG2404029

Long Extension, 7"x8" Footplates
Part: RIG2404028
Long Extension, 7"x11" Footplates
Part: RIG2404030

3. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8" \$69
Part: RIG2404919

Outside Edge \$69
Part: RIG2404922

Width Extension for 11" \$69
Part: RIG2404920

Available with 7" width footplates.

Outside Edge \$69
Part: RIG2404921

Front and Rear Cap \$16
Part: RIG2404923

Available with 5"x8" footplates.

Available with 5"x8" footplates.

5. Center Mount Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133717

Footplate Connector \$16
Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108

7"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG7008T-Q. HCPCS: K0108

7"x11" Tapered Gel Footplate Cover for Width Extension \$293
Part: ST-AFPG7011T-Q. HCPCS: K0108

17.2. Legrest Accessories

Wheel Bumper \$37
Part: XRF000100

17.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)
10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)

18. POWER POSITIONING ELECTRONICS**REQUIRED.****18.1a. Q-Logic 3 Controls through Joystick**

Single Actuator Through Controller and Interface \$2,163
Part: XRF000335. HCPCS: E2310

Single Actuator Through Controller \$3,040
Interface (AAM)

Multiple Actuator Through Controller and Interface \$3,040
Part: XRF000336. HCPCS: E2311

Part: XRF000068. HCPCS: E2310

18.1b. Q-Logic 3 Controls with iAccess Switch Control

Single Actuator Control Through Flat Button \$2,547
iAccess and Controller and Interface
Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000074. HCPCS: E2310

Single Actuator Control Through Toggle Button \$2,547
iAccess and Controller and Interface
Part: XRF000338. HCPCS: E2310

Multiple Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface
Part: XRF000339. HCPCS: E2311

Single Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000070. HCPCS: E2310

Multiple Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface
Part: XRF000340. HCPCS: E2311

Power positioning electronics section continued on following page.

18. POWER POSITIONING ELECTRONICS CONTINUED**18.2. iAccess Mount****REQUIRED with iAccess options.****1. iAccess Mount Selection**

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$26 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26 Part: MEC2001392
Under Arm Mount, Tube Style Armrest No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26 Part: MEC2001391

2. Under Arm Mount Side

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

19. ATTENDANT CONTROL**Optional.** Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.


Attendant Control \$525 Part: XRF000078. HCPCS: E2331 Select this option if no charge option criteria is not met.	Attendant Control with Tilt, Recline, Alternative Controls and Q-Logic Display or SCIM Module No Charge Part: XRF000079
---	--

20. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045 Part: XRF000130 Not available with Base Mounted, LED, Full Light Package.

21. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

Accessory Bar \$132 Part: XRF000128 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.	Oxygen Holder \$210 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Backpack Holder (Pair) (A) \$69 Part: ACC147774	Personal Item Hook \$69 Part: ACC125008
Cane/Crutch Holder \$105 Part: ACC152515. HCPCS: E2207	Privacy Bag \$42 Part: ACC2005723
Clothing Guard (Pair) \$237 Part: ACC125007	Push Handles (pair) \$195 Part: ACC2005027 Not available with oxygen holder or vent tray.
Compact Vent Tray \$1,465 Part: ACC133889. HCPCS: E1029	Q-Logic EX/NE Programmer \$573 Part: ELEASMB6976
Cup Holder No Charge Part: ACC125003	Quantum Backpack \$168 Part: ACC1905853
Flag Kit \$69 Part: ACC1700443	Rear View Mirror \$69 Part: ACC155833
Flip-Back Armrest, Side Track Accessory Bracket ... \$11 Part: MEC1809178	Reflector Kit \$53 Part: ACC124670
Front Mounted, Cell Phone/ Media Holder, \$69 X-Grip (E) Part: ACC155835	Stealth Large Tablet Holder \$625.28 Part: ACC141818
Glove Box (C) \$158 Part: ACC125011	Transfer Handle (D) \$105 ea. Part: MEC1809930 x2. HCPCS: E0705
Grip Drink Holder (20085) \$15.50 Part: ACC152993	Personal Item Pouch \$84 Part: ACC2101017
Grip Lap Board (20005) \$38.95 Part: ACC152992	Depth Adjustable Vent Tray \$595 Part: ACC2302500. HCPCS: E1029 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Hydration System (B) \$237 Part: ACC125010	Remote Stop Switch \$1,295 Part: XRF000379
Large, Black, Weather Cover \$202 Part: ACCCOVR1021	

**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea. Part: ACC2407046	Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea. Part: ACC2405496
--	---

22. BRANCH PATCHES

Optional.

Air Force Patch\$21
Part: ACC140240

Army Patch\$21
Part: ACC140241

Coast Guard Patch\$21
Part: ACC140239

Marines Patch\$21
Part: ACC140243

Navy Patch\$21
Part: ACC140244

Purple Heart Patch\$21
Part: ACC140238

**23. POSITIONING COMPONENTS**Optional. More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable & Removable Hardware \$265 ea.
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable & Removable Hardware \$276 ea.
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q

UniLink, 2" Link \$81 ea.
Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
Part: PTOINDV3585

Install on Right
Part: PTOINDV3586

Do Not Install
Part: PTOINDV3719

**23.2. Stealth Thoracic Laterals**

Swing-away laterals may interfere with flip-back and 2-post, flip up armrests.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664
Part: ST-ULS4042000-Q. HCPCS: E1034
Only available with Redesigned back assemblies.

Swing-away, UniLink, 125" w/Track, pair (ULPK404108B00) ... \$664
Part: ST-ULPK404108B00-Q. HCPCS: E1034
Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034
Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332
Part: ST-ULPK404103L00-Q. HCPCS: E1034
Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332
Part: ST-ULPK404103R00-Q. HCPCS: E1034
Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250
Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
Part: PTOINDV3499

Do Not Install
Part: PTOINDV3500

23.4. Chest Bar Assembly

Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.

Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$368
Part: ST-CBA-TB4-S. HCPCS: K0108
Available with 14" to 20" back widths.

Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$368
Part: ST-CBA-TB4-M. HCPCS: K0108
Available with 17" to 22" back widths.

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.5. Stealth Positioning Chest Harness**

Upper Classic, Medium (IFPS1103) \$154 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

**23.6. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor\$97 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor\$97 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

**23.7. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

Facial Lateral for Comfort Series..... \$318 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Lateral for Comfort Series..... \$318 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028

Facial Laterals for Comfort Series \$376 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

**2. Facial Lateral Pads and Accessories**


Egg Switch Mount (ESM) (A)\$57 ea.
 Part: INDPART2815

Facial Pad, 4"x2" (900-Q).....\$61 ea.
 Part: ST-900-Q. HCPCS: K0108

Facial Pad, 6"x2" (920-Q).....\$61 ea.
 Part: ST-920-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C)..... \$276 ea.
 Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest

**23.8. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support..... \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q


Right, Swing-Away, Elbow Support..... \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket.....\$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$57 ea.
 Part: POS124067. HCPCS: K0108

**23.9. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

23.10. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Quantum Lap Tray****Mounting bracket is required with lap tray selection.** These trays are compatible with Flip-Back armrests.**1a. Black Lap Tray**

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

24.2. Simple Solutions Lap Tray**Mounting bracket is required with lap tray selection.** See chart below for compatibility.**1. Lap Tray**

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

:Compatible:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**25.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,980

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN137728

Grey Street \$63

Part: DGN137729-PSL3

Midnight Blue \$63

Part: DGN137729-PBU10

Red Corvette..... \$63

Part: DGN137729-PRD1

Deep Purple \$63

Part: DGN137729-PPR5

Electric Blue \$63

Part: DGN137729-PBU9

Military Green \$63

Part: DGN137729-PGN8

Back in Black..... \$63

Part: DGN137729-PBK1

Pretty in Pink \$63

Part: DGN137729-PPK3

Lemon Crush..... \$63

Part: DGN137729-PYL4

Military Digital Camo..... \$63

Part: DGN137729-WCM7

Night Train..... \$63

Part: DGN137729-PBK3

Evolution Orange..... \$63

Part: DGN137729-POR7

Green Machine \$63

Part: DGN137729-PGN11

Faded Glory \$63

Part: DGN137729-WFL2

25.2. Q-Logic SCIM Module**Optional.****SCIM Module Kit for Q-Logic 3**.....\$525

Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 Expandable Capacity Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ**.....\$840

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

25.4. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch Head Array (IDHBT200-1)** \$4,774

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

25.5. Sip and Puff Drive**1. Module****Sip and Puff Module for Q-Logic 3**.....\$1,990

Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series,** \$318**Left Mount (SP625-CL)**

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$318**Head Array (SP625-U)**

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, \$318**Right Mount (SP622-CR)**

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Sip & Puff Module

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)**.....\$35

Part: ST-SP625-S12-Q



Stealth Sip & Puff

25.6. Egg Switches**Red Egg Switch (ESRED)**..... \$87

Part: SWTMCRO1077

Blue Egg Switch (ESBLU)..... \$87

Part: SWTMCRO1079

Black Egg Switch (ESBLK)..... \$87

Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$87

Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$87

Part: SWTMCRO1080



Egg Switch

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

25.7. Switch Adapter

2 to 1 Switch Adapter \$15

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520

Part: CTL168057

