



Quantum® 4Front® 2 Order Form with TRU Balance® Flex

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____
 A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 2 3S-SS..... \$8,870

Part: 4FRONT 2 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1a. Standard Shrouds

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN2003028-PRD1 |  Deep Purple
Part: DGN2003028-PPR5 |  Midnight Blue
Part: DGN2003028-PBU10 |  Green Machine
Part: DGN2003028-PGN11 |
|  Back in Black
Part: DGN2003028-PBK1 |  Raspberry Beret
Part: DGN2003028-PPK4 |  Electric Blue
Part: DGN2003028-PBU9 |  Military Green
Part: DGN2003028-PGN8 |
|  Night Train
Part: DGN2003028-PBK3 |  Pretty in Pink
Part: DGN2003028-PPK3 |  Blue Jean
Part: DGN2003028-PBU11 | |
|  Grey Street
Part: DGN2003028-PSL3 |  Evolution Orange
Part: DGN2003028-POR7 |  Lemon Crush
Part: DGN2003028-PYL4 | |

2.1b. Hydrographic Shrouds

- | | |
|---|---|
|  Military Digital Camo \$263
Part: DGN2003028-WCM7 |  Faded Glory \$263
Part: DGN2003028-WFL2 |
|---|---|

3. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

- | | |
|--|---|
| Switch Operated Accent Lights No Charge
Part: DGN2004686
Accent lights operate by pressing included single switch button. Accent lighting is not available with Single Actuator Control through Toggle Box. | Omit Accent Lights No Charge
Part: DGN1700032 |
|--|---|



4. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

- Quantum Backup Camera** \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

- | | |
|--|---|
| Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only. | Backup Camera Track Mount Bracket \$27
Part: MEC2104227 |
|--|---|

2. Optional Display Mounting Bracket

- Backup Camera Q-Logic Joystick Mount Bracket** \$105
Part: MEC2106123



5. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color.

5.1. Drive Wheels

- | | |
|--|--|
| 3"x14" Flat Free Black Tires with Accent Standard
Part: PTOINDV3525 | 4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368 |
| 3"x14" Pneumatic Black Tires with Accent No Charge
Part: WHL1704059 x2 | 4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367 |

5.2. Drive Wheel Accent Color

Not available with knobby tires.

- | | | |
|---|--|--|
| Back in Black Standard
Part: DGN164245-PBK1 x2 | Pretty in Pink No Charge
Part: DGN164245-PPK3 x2 | Green Machine No Charge
Part: DGN164245-PGN11 x2 |
| Red Corvette No Charge
Part: DGN164245-PRD1 x2 | Evolution Orange No Charge
Part: DGN164245-POR7 x2 | Military Green No Charge
Part: DGN164245-PGN8 x2 |
| Night Train No Charge
Part: DGN164245-PBK3 x2 | Midnight Blue No Charge
Part: DGN164245-PBU10 x2 | Military Digital Camo \$58 ea.
Part: DGN164245-WCM7 x2 |
| Grey Street No Charge
Part: DGN164245-PSL3 x2 | Electric Blue No Charge
Part: DGN164245-PBU9 x2 | Faded Glory \$58 ea.
Part: DGN164245-WFL2 x2 |
| Deep Purple No Charge
Part: DGN164245-PPR5 x2 | Blue Jean No Charge
Part: DGN164245-PBU11 x2 | |
| Raspberry Beret No Charge
Part: DGN164245-PPK4 x2 | Lemon Crush No Charge
Part: DGN164245-PYL4 x2 | |

6. BATTERIES

Optional. The powerchair requires 2 flat top batteries and will be installed if ordered.

- | | |
|--|--|
| Group 24, AGM Cell Battery (MK) \$541 ea.
Part: BAT1712207. HCPCS: E2363 | Group 34, Gel Cell Battery (MK) \$509 ea.
Part: BATGEL1006. HCPCS: E2359 |
|--|--|

7. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight
Part: PTOINDV3837

201 lb-300 lb User Weight
Part: PTOINDV3838

8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

8.1.a. Omit Joystick

No Joystick, Setup for Q-Logic Electronics
Part: XRF000010
Requires selection of Enhanced Display

8.1.b. Q-Logic 3 Non-Expandable Capacity Controls

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.

QLNES 10 Key Joystick \$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

8.1.c. Q-Logic 3 Expandable Capacity Controls

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity Joysticks

Q-Logic 3 Expandable Capacity Joystick\$840
Part: CTL167724

Q-Logic 3e Joystick\$735
Part: CTL1910600

Q-Logic 3 Expandable Capacity Joystick\$840
for Full Light Package
Part: CTL167725
Required with full light packages.

Can operate full light package. Not available with alternative or specialty controls.



2. Q-Logic 3 Expandable Capacity Required Components

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$730
Part: XRF000002 HCPCS: E2377

Q-Logic 3 Expandable Capacity Controller Harness..... \$458
Part: XRF000264. HCPCS: E2313

8.2. Q-Logic Magnetic Charger

Available with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635

8.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds.
Red Corvette.....\$63
Part: DGN124336-PRD1
Back in Black.....\$63
Part: DGN124336-PBK1
Grey Street.....\$63
Part: DGN124336-PSL3
Deep Purple.....\$63
Part: DGN124336-PPR5

Raspberry Beret..... \$63
Part: DGN124336-PPK4
Pretty in Pink.....\$63
Part: DGN124336-PPK3
Evolution Orange.....\$63
Part: DGN124336-POR7
Midnight Blue\$63
Part: DGN124336-PBU10
Electric Blue\$63
Part: DGN124336-PBU9

Blue Jean \$63
Part: DGN124336-PBU11
Lemon Crush.....\$63
Part: DGN124336-PYL4
Green Machine\$63
Part: DGN124336-PGN11
Military Green\$63
Part: DGN124336-PGN8

9. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323
Chin Cup (C)..... \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)
T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323
Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323
Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323
Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323
Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323
Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



10. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

10.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket Standard
 Part: FRMASMB10617



10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
 Part: ACC1908434. HCPCS: E1032



10.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245
 Part: ACC2201041. HCPCS: E1032

Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201042. HCPCS: E1032

10.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount \$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.



10.2. Joystick Mount Location

Left Side Mount
 Part: PTOINDV1225

Right Side Mount
 Part: PTOINDV1226

11. TRU BALANCE FLEX SEATING SYSTEM

Select one option from each section. Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.

11.1. TRU Balance Flex Frame Selection

TRU Balance Flex Seating System No Charge
 Part: PTOINDV3853

11.2. Occupied Transit Kit

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

WC-19 Occupied Transit Option \$289
 Part: KIT152047

11.3. Width

Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.4. Depth

Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

High	18.5" at front with 4° posterior angle
Med	17.5" at front with 2° posterior angle
TRU Balance Flex Seat to Floor Heights	


11.5. Seat to Floor Height

Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

11.6. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



TRU Balance Flex Seating System section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

1. Cushion Type

TRU-Comfort 2 SPP
 Reference part: XRF000116. *See cushion order form.*
 Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
 Reference part: XRF000118. *See cushion order form.*
 Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
 Reference part: XRF000113. *See cushion order form.*
 Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
 Reference part: XRF000112. *See cushion order form.*
 Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
 Reference part: XRF000117. *See cushion order form.*
 Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
 Reference part: XRF000122. *See cushion order form.*
 Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
 Reference part: XRF000114. *See cushion order form.*
 Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
 Reference part: XRF000115. *See cushion order form.*
 Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
 Reference part: XRF000111. *See cushion order form.*
 Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39 Part: POP1909006-S 14"W and smaller	Medium Wedge Positioning Kit\$72 Part: POP1909006-M 14"W and smaller	Large Wedge Positioning Kit \$105 Part: POP1909006-L 14"W and smaller
Small Wedge Positioning Kit\$39 Part: POP1712340-S 15"W and larger	Medium Wedge Positioning Kit\$72 Part: POP1712340-M 15"W and larger	Large Wedge Positioning Kit \$105 Part: POP1712340-L 15"W and larger

3. Cushion Top Foam

Visco Top Foam Reference part: XRF000380	No Top Foam Reference part: XRF000382
FRG Top Foam Reference part: XRF000381	

4. Cushion Liner

Stealth Urethane Liner Reference part: XRF000383	No Liner Reference part: XRF000385
Stealth Polyurethane Liner Reference part: XRF000384	

5. Cushion Cover

CoolCore Cover Reference part: XRF000386	Stealth Tek Cover Reference part: XRF000388
Mesh Cover Reference part: XRF000387	

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	✓	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389
 MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

<u> </u> Additional Stealth Urethane Liner Reference Part: XRF000390	<u> </u> Additional Stealth Polyurethane Liner Reference Part: XRF000391
--	--

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

<u> </u> Additional CoolCore Cover Reference Part: XRF000392	<u> </u> Additional Stealth Tek Cover Reference Part: XRF000394
<u> </u> Additional Mesh Cover Reference Part: XRF000393	

TRU Balance Flex Seating System section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

11.7. TRU Balance Flex Back Cane Selection

- Angle Adjustable Back Cane..... No Charge
Part: MEC140266
- Flip Forward Back Cane \$184
Part: XRF000274
Not available with oxygen holder, or cane mounted and cantilever armrests.
- 2" Mesomorph (1" per side) Back Cane..... \$205
Part: MEC1804736
- 2" Endomorph (1" per side) Back Cane \$205
Part: XRF000270

- 4" Endomorph (2" per side) Back Cane..... \$205
Part: XRF000272
- Manual Recline Back \$788
Part: FRMASMB12910. HCPCS: E1226
Select a back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

11.8. TRU Balance Flex Back Cane Height Selection

16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

11.9. TRU Balance Flex Back Cushion Options

Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Synergy Back Selection

- Synergy Back..... \$289
Part: XRF000363

1b. Sport Back Selection

- Sport Back Cushion, 21" Wide or Less..... \$394
Part: XRF000361
- Sport Back Cushion, 22" Wide or More..... \$604
Part: XRF000362

1c. TRU-Comfort Back Selection

- TRU-Comfort Back, Back Cane Option, Up to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620
Back size based on seat width, back cane type, and overall height option selections.
- TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621
Back size based on seat width, back cane type, and overall height option selections.

		Back Width										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	13"	-	-	S	-	-	-	-	-	-	-	-
	14"	P	T	-	S	P	-	-	-	-	-	-
	15"	S	P	T	-	S	P	-	-	-	-	-
	16"	S	P	T	-	S	P	T	P	T	P	T
	17"	P	T	S	S	P	T	S	P	T	S	P
	18"	-	S	S	S	S	P	T	S	P	T	S
	19"	-	-	-	S	S	P	T	S	P	T	S
	20"	-	-	-	S	S	P	T	S	P	T	S
	21"	-	-	S	S	S	P	T	S	P	T	S
	22"	-	-	S	S	S	P	T	S	P	T	S
	23"	-	-	S	S	S	P	T	S	P	T	S
24"	-	-	S	S	S	P	T	S	P	T	S	
25"	-	-	S	S	S	P	T	S	P	T	S	

S Synergy Back
 P Sport Back
 T TRU-Comfort Back

TRU Balance Flex Seat Back Availability Chart

TRU Balance Flex Seating System section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

1d. ADI Back, AL Series, 4-Point, Active Contour

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|---|---|
| Active Contour, 16" Wide x 19" or 20" Overall Height \$622
Part: ALBR16SA4. HCPCS: E2613 | Active Contour, 18" Wide x 23" or 24" Overall Height \$622
Part: ALBR18XA4. HCPCS: E2613 |
| Active Contour, 16" Wide x 21" or 22" Overall Height \$622
Part: ALBR16TA4. HCPCS: E2613 | Active Contour, 20" Wide x 19" or 20" Overall Height \$622
Part: ALBR20SA4. HCPCS: E2613 |
| Active Contour, 16" Wide x 23" or 24" Overall Height \$622
Part: ALBR16XA4. HCPCS: E2613 | Active Contour, 20" Wide x 21" or 22" Overall Height \$622
Part: ALBR20TA4. HCPCS: E2613 |
| Active Contour, 18" Wide x 19" or 20" Overall Height \$622
Part: ALBR18SA4. HCPCS: E2613 | Active Contour, 20" Wide x 23" or 24" Overall Height \$622
Part: ALBR20XA4. HCPCS: E2613 |
| Active Contour, 18" Wide x 21" or 22" Overall Height \$622
Part: ALBR18TA4. HCPCS: E2613 | |

1e. ADI Back, AL Series, 4-Point, Deep Contour

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|--|---|
| Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR16SD4. HCPCS: E2615 | Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR18TD4. HCPCS: E2615 |
| Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735
Part: ALBR16TD4. HCPCS: E2615 | Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
Part: ALBR20SD4. HCPCS: E2615 |
| Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735
Part: ALBR16XD4. HCPCS: E2615 | Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
Part: ALBR20XD4. HCPCS: E2615 |
| Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735
Part: ALBR18SD4. HCPCS: E2615 | |

1f. ADI Back, AL Series, 4-Point, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

- | | |
|--|--|
| Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
Part: ALBR18S64. HCPCS: E2615 | Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
Part: ALBR18X64. HCPCS: E2615 |
| Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
Part: ALBR18T64. HCPCS: E2615 | |

1h. Omit Back Selection

- Omit Back..... No Charge
Part: PTOINDV1287

11.10. Back Width Package

REQUIRED with Sport and TRU-Comfort back selections.

- | | |
|---|---|
| 16"-20" Back Width Package
Part: XRF000038 | 18"-22" Back Width Package
Part: XRF000039 |
|---|---|

11.11. Overall Back Height

- | | | | |
|--|--|--|--|
| 16" Overall Back Height
Part: SETHGHT1029 | 19" Overall Back Height
Part: SETHGHT1032 | 22" Overall Back Height
Part: SETHGHT1035 | 25" Overall Back Height
Part: SETHGHT1038 |
| 17" Overall Back Height
Part: SETHGHT1030 | 20" Overall Back Height
Part: SETHGHT1033 | 23" Overall Back Height
Part: SETHGHT1036 | |
| 18" Overall Back Height
Part: SETHGHT1031 | 21" Overall Back Height
Part: SETHGHT1034 | 24" Overall Back Height
Part: SETHGHT1037 | |

11.12. TRU Balance Flex Armrests

Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Includes flat armpads.

- | | |
|--|---|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Drop-In Armrests..... No Charge
L R Left part: ARM1803253. Right part: ARM156742 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743 |
|--|---|

1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

- | | |
|--|---|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Flip Back Armrests..... \$63
L R Left part: XRF000371. Right part: XRF000372 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Flip Back Armrests \$63
L R Left part: XRF000373. Right part: XRF000374 |
|--|---|

1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

- | | |
|--|---|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Two-Post Armrests..... \$158 ea.
L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Two-Post Armrests \$158 ea.
L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973 |
|--|---|

1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

- | | |
|---|---|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973 |
|---|---|

TRU Balance Flex Seating System section continued on following pages.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED

11.13. TRU Balance Flex Armpads

Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads

Straight, Full Length Left and Right No Charge
Part: ARM134275

Straight, Desk Length Left, Full Length Right No Charge
Part: ARM134280

Straight, Full Length Left, Desk Length Right No Charge
Part: ARM134281

Straight, Desk Length Left and Right No Charge
Part: ARM134254

1b. Waterfall Armpads

Waterfall, Full Length Left and Right No Charge
Part: ARM134285

Waterfall, Desk Length Left, Full Length Right No Charge
Part: ARM134288

Waterfall, Full Length Left, Desk Length Right No Charge
Part: ARM134289

Waterfall, Desk Length Left and Right No Charge
Part: ARM134282

1c. Falcon Gel Armpads

Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right\$287.93
Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right\$287.93
Part: ARM136442

Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437

12. HEADRESTS

Optional. Select a headrest pad and bracket.

12.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q, HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q, HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179
Part: ST-CP10CCFPM-Q, HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q, HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink\$311
Part: ST-UHB2D33200-Q, HCPCS: E1033

Multi-Axis Removable Headrest Bracket\$174
Part: ST-HMO475-17M-Q, HCPCS: E1033

12.2. Q Cam Lock

Optional.

Q Cam Lock \$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

13. LAP BELTS

REQUIRED. Select a lap belt below.

13.1a. Standard Lap Belts

38" Lap Belt
Part: ACCBELT1006

50" Lap Belt
Part: ACCBELT1000

60" Lap Belt
Part: ACCBELT1002

70" Lap Belt
Part: ACCBELT1003

80" Lap Belt
Part: ACCBELT1004

90" Lap Belt
Part: ACCBELT1001

13.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350)\$181
Part: ST-IFRB5350

Stealth Positioning 60" Single Pull, Push Button, Extra Large \$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

Stealth Positioning 60" Single Pull, Push Button, Medium \$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

60" Sure-Lok Lap Belt (Non-retractable) \$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.

Stealth Positioning 60" Single Pull, Push Button, Large \$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

14. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

14.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
 Part: XRF000323

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates No Charge Part: MEC1906675	Medium Lower Extension, 6.25"W Tapered 5"W x 12"D... No Charge Part: MEC2104853
Medium Extension, 5"x9" Footplates No Charge Part: MEC1906680	Long Lower Extension, 6.25"W Tapered 5"W x 12"D..... No Charge Part: MEC2104854
Long Extension, 5"x9" Footplates No Charge Part: MEC1906677	

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158 Part: RIG2104369	Large Calf Pads, Basic Center Mount \$158 Part: RIG2202997 Only available with Long Length Extension.
--	---

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

14.1b. Center Mount Foot Platform

Select legrest length in following section.

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

1. Foot Platform Assembly

Quick Adjustable Center Mount Foot Platform \$263
 Part: XRF000311. HCPCS: K0040

2. Extension and Angle Adjustable Footplates

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 7"x8" Footplates Part: RIG2404027	Long Extension, 7"x8" Footplates Part: RIG2404028
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Medium Extension, 7"x11" Footplates Part: RIG2404029	Long Extension, 7"x11" Footplates Part: RIG2404030

3a. Center Mount Calf Support

Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, Quick Adjustable Center Mount Only, \$158 6"Wx6"H (pair) Part: XRF000156	Large Calf Pads, Quick Adjustable Center Mount Only, \$158 8" taper to 6"Wx8"H (pair) Part: XRF000157
---	---

3b. Center Mount Calf Panel

Small Calf Panel for Quick Adjustable, 18"Wx8"H \$184 Part: FRMASMB7110012	Large Calf Panel for Quick Adjustable, 24"Wx8"H \$237 Part: RIG123825 Only available on widths of 20" or greater.
---	---

4. Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$69 Part: RIG2404922 Available with 7" width footplates.
Width Extension for 8" \$69 Part: RIG2404919 Available with 5"x8" footplates.	Front and Rear Cap \$16 Part: RIG2404923 Available with 5"x8" footplates.
Outside Edge \$69 Part: RIG2404921 Available with 5"x8" footplates.	

Legrests section continued on following pages.

14. LEGRESTS CONTINUED

5. Additional Accessories

- Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q
- Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q
- Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q
- Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q
- Footplate Connector \$16
Part: RIGASMB7120023

- Footplate Riser \$75
Part: XRF000348
Not available with Long Extension. Select legrest length between 4" to 7" at 1" increments.
- Footplate Riser Calf Panel \$50
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

- 5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108
- 7"x8" Tapered Gel Footplate Covers \$270
Part: ST-AFPG7008T-Q. HCPCS: K0108

- 5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108
- 7"x11" Tapered Gel Footplate Cover \$293
Part: ST-AFPG7011T-Q. HCPCS: K0108

14.1c. 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. 60°-90° Swing-Away Legrest Upper

- P L R 60°-90° Swing-Away Legrest \$81 ea.
Pair part: XRF000082. Left: XRF000083.
Right: XRF000084
Not available with 11"-13" extensions.

- P L R 60°-90° Swing-Away Legrest, Short \$81 ea.
Pair part: XRF000085. Left: XRF000086.
Right: XRF000087.
For use only with 11"-13" lower extensions.

2. Swing-Away Legrests Lower Extension

- P L R 11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.
- P L R 13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.
- P L R 15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

Legrest Compatibility Matrix 60°-90° Swing-Away Legrest TRU Balance Flex	Seat to Floor Option		
	Low	Middle	High
Extruded/Angle A., 11"-13"	Full	Full	Full
Extruded/Angle A., 13"-16.5"	13"	13"-14"	13"-15"
Extruded/Angle A., 15.5"-19.5"	N/A	N/A	N/A
Multi-Axis, 11"-16.5"	Full	Full	Full
Multi-Axis, 15.5"-19.5"	N/A	N/A	N/A

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

- P L R 5"Wx8"D with Heel Loops \$90 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951
- P L R 6"Wx8"D¹ with Heel Loops \$90 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

- P L R 6"Wx11"D¹ with Heel Loops \$90 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951
- P L R 8"Wx11"D² with Heel Loops \$90 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- Therafin 31515 Multi-Axis Amputee Support \$230 ea.
Part: INDPART2822. HCPCS: E1020

- Therafin 32387 Amputee Support Gel Cover \$155 ea.
Part: INDPART2823. HCPCS: K0108

14.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

- P L R Elevating Legrest Assembly \$171 ea.
Pair part: FRMASMB6483. Left: RIG138012
Right: RIG138013. HCPCS: E0990
Select lower extension and footplates

Legrest Compatibility Matrix Elevating Legrest TRU Balance Flex	Seat to Floor Option		
	Low	Middle	High
Reinforced, 11"-14"	11"-13"	Full	Full
Reinforced, 14"-17"	N/A	14"	14"-15"
Angle Adjustable, 12"-16.5"	12"-13"	12.5"-14"	12.5"-15"
Angle Adjustable, 15.5"-19.5"	N/A	N/A	N/A

2a. Lower Extension

- P L R 11"-14" Length
Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
For reinforced only
- P L R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
Pair part: XRF000101. Left: XRF000102. Right: XRF000103
For angle adjustable only

- P L R 14"-17" Length
Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
For reinforced only
- P L R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
Pair part: XRF000104. Left: XRF000105. Right: XRF000106
For angle adjustable only

Legrests section continued on following pages.

14. LEGRESTS CONTINUED

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342
Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R 7.25"Wx6"D Reinforced Footplate..... Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R 8.5"Wx6"D Reinforced Footplate Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D with Heel Loops..... \$90 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

P L R 6"Wx11"D¹ with Heel Loops..... \$90 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

P L R 6"Wx8"D¹ with Heel Loops \$90 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

P L R 8"Wx11"D² with Heel Loops..... \$90 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
Part: INDPART2242. HCPCS: E0953

14.2. Legrest Accessories

Heel Loops (pair)..... \$53
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$37
Part: XRF000100

14.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

15. SEAT MOUNTED LED LIGHTS

Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130
Not available with Base Mounted, LED, Full Light Package.

16. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Flex Seat Cup Holder..... No Charge
Part: ACCASMB2634
Not available on TRU Balance Flex seating without two-post armrests.
Oxygen Holder for Flex Seat \$210
Part: ACC125006. HCPCS: E2208
Not available with Manual Recline.
Walker Holder \$109
Flex part: FRMASMB7031
Cane/Crutch Holder \$105
Flex part: ACC152515. HCPCS: E2207
Flag Holder \$55
Flex part: ACCASMB2415.
Saddle Bag.....\$40
Part: DWR1020D010

Stealth Cell Phone/Media Device Holder \$215
Part: ACC7120006
Stealth Large Tablet Holder \$625.28
Part: ACC141818
Quantum Backpack \$168
Part: ACC1905853
Grip Lap Board (20005).....\$38.95
Part: ACC152992
Grip Drink Holder (20085) \$15.50
Part: ACC152993
Accessory Bar\$132
Part: XRF000128

Accessories section continued on following page.

16. ACCESSORIES CONTINUED

Large, Black, Weather Cover \$202 Part: ACCCOVR1021	Transfer Handle \$105 ea. Part: MEC2205680. HCPCS: E0705
Privacy Bag \$42 Part: ACC2005723	Remote Stop Switch \$1,295 Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea. QTY Part: ACC2407046	Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea. QTY Part: ACC2405496
--	---

17. BRANCH PATCHES

Optional.

Air Force Patch\$21 Part: ACC140240	Coast Guard Patch\$21 Part: ACC140239	Navy Patch\$21 Part: ACC140244
Army Patch\$21 Part: ACC140241	Marines Patch\$21 Part: ACC140243	Purple Heart Patch\$21 Part: ACC140238



18. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

18.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$562 Part: INDPART2896. HCPCS: E1034	Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1034 Only available with back cane and omit back configurations.
---	---

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$250 Part: INDPART2663. HCPCS: E0956	Flat Pads, 3"x5", pair (ULP846050320) \$250 Coming soon! Part: ST-ULP846050320-Q. HCPCS: E0956
Curved Pads, 4"x6", pair (TWBL-4X6) \$250 Part: INDPART2664. HCPCS: E0956	Flat Pads, 4"x6", pair (ULP835050320) \$250 Coming soon! Part: ST-ULP835050320-Q. HCPCS: E0956

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

4. Lateral Extensions

REQUIRED with Synergy Back and thoracic lateral selections.

Lateral Extension for Synergy Back (pair) \$28 ea. Part: XRF000150

18.2. Stealth UniLink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable & Removable Hardware \$265 ea. QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable & Removable Hardware \$276 ea. QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Left \$358.87 ea. Part: ST-UTT8A116O2-L-Q. HCPCS: E1028/E1034^
Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^	UniLink Multi-Angle Thigh Guide, Right \$358.87 ea. Part: ST-UTT8A116O2-R-Q. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad \$106 ea. QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad \$106 ea. QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. QTY Part: ST-UIL001150A000-Q	UniLink, 2" Link \$81 ea. QTY Part: ST-UIL001200A000-Q
---	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

Positioning components section continued on following page.

18. POSITIONING COMPONENTS CONTINUED

18.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960
 Upper Classic, Large (IFPS1104) \$154 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960
 Upper Contour, Large (IFPS2104) \$154 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

18.4. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor\$97 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor\$97 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960



Chest Strap

18.5. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
 Part: ST-SAR2A-R-Q. HCPCS: E1028
 UniLink Facial Lateral, Swing Away, Left \$318
 Part: ST-SAR2A-L-Q. HCPCS: E1028
 UniLink Facial Lateral, Dual Swing Away\$492
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
 UniLink Facial Lateral, Dual Swing Away\$434
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028




Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and Accessories

_____ Egg Switch Mount (ESM) (A)\$57 ea.
 QTY Part: INDPART2815
 _____ Figure 8 Bracket with 13" Rod..... \$289 ea.
 QTY and Egg Switch Mount (B)
 Part: ST-F8675-13-Q
 _____ Facial Pad, 4".5x2", Cool Core (900M-Q)\$61 ea.
 QTY Part: ST-900M-Q. HCPCS: K0108

_____ Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea.
 QTY Part: ST-920M-Q. HCPCS: K0108
 _____ Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
 QTY Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

18.6. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket\$199 ea.
 Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket.....\$199 ea.
 Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

_____ 3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
 QTY Part: POS124065. HCPCS: K0108

_____ 4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
 QTY Part: POS124067. HCPCS: K0108

18.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954
 Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

18.8. Stealth Padded Shoe Holder

_____ Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
 QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952
 _____ Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
 QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

_____ Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

19. LAP TRAY

19.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q, HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
 Part: ST-SSLXL010001171-Q, HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
 Part: ST-SSLXM010001171-Q, HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	L	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	L	M	L	L
Drop-In, T-Arm	S	S	S	S	S	S	M	M	L	M	L	L
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	L	M	L	L
Cantilever	S	S	S	S	S	M	M	M	L	M	L	L
Q-Captain, Flip-Up						S	M		L			
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	M	L	M	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
 :Compatible :Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.