



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Jazzy® 1450 Order Form with Synergy® Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

1. BASE MODEL

REQUIRED. Select a group 2 power base. 600 lb. weight capacity with standard 14" flat free tires, 4.5 mph motors, and 8 amp off board charger. Dimensions of base without legrests: 29.25"W x 38.25"L.

Jazzy 1450 2VHD-SS..... \$6,729
Part: J1450 2VHD-SS. HCPCS: K0826


2. BASE COLOR


REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

Candy Apple Red
Part: SHDASMB6703




Viper Blue
Part: SHDASMB6704



3. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.Group 24, Gel Cell, Interceptor 2 Battery.....\$399 ea.
Part: BAT2100315. HCPCS: E2363**4. ELECTRONICS****REQUIRED.**VR2 6-Key Joystick Standard
Part: CTLDC1498**5. JOYSTICK HANDLES****Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323Chin Cup (C) \$90
Part: ACCASMB1861.
HCPCS: E2323. (E2324 when chin control)T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308**6. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**6.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for VR2..... Standard
Part: FRMASMB9532**6.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, VR2 \$319
Part: ACC1908439. HCPCS: E1032**6.1c. Basic Swing Away Joystick Mount**Basic Swing Away, VR2, Left Side \$245
Part: ACC2201051. HCPCS: E1032Basic Swing Away, VR2, Right Side \$245
Part: ACC2201052. HCPCS: E1032**6.2. Joystick Mount Location**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**7. SYNERGY SEATING SYSTEM**

Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 32"W x 28"D.

7.1. Synergy Frame Selection**REQUIRED.**Synergy Seating System \$840
Part: XRF000261**7.2. Seat Width****REQUIRED.** Maximum seat width is 32". Part number is SETWDTH plus number following # below.

20" W	21" W	22" W	23" W	24" W	25" W	26" W	27" W	28" W	29" W	30" W	31" W	32" W
(#1010)	(#1011)	(#1012)	(#1013)	(#1014)	(#1056)	(#1054)	(#1057)	(#1055)	(#1058)	(#1059)	(#1060)	(#1061)

7.3. Seat Width Package**REQUIRED.**20"-24" Width Package
Part: FRMASMB712001426"-28" Width Package
Part: FRMASMB712001930"-32" Width Package
Part: FRMASMB712002124"-26" Width Package
Part: FRMASMB712001528"-30" Width Package
Part: FRMASMB7120020**7.4. Seat Depth****REQUIRED.** Maximum seat depth is 28". Part number is SETDPTH plus number following # below.

20" D	21" D	22" D	23" D	24" D	25" D	26" D	27" D	28" D
(#1010)	(#1011)	(#1012)	(#1013)	(#1014)	(#1060)	(#1061)	(#1062)	(#1063)

7.5. Seat Depth Package**REQUIRED.**15"-20" Depth Package
Part: FRMASMB773220"-24" Depth Package
Part: FRMASMB773324"-28" Depth Package
Part: FRMASMB6605

Synergy seating continued on following page.

7. SYNERGY SEATING SYSTEM CONTINUED**7.6. Seat to Floor****REQUIRED.**

Lowest Seat to Floor
Part: PTOINDV3558


Middle Seat to Floor
Part: PTOINDV3559

Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
Quantum 1450		
16.5"	17.5"	18.5"
Synergy Seating Seat to Floor Heights		

7.7. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**

Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP

Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP

Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP

Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP

Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P

Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP

Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP

Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G

Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit.....\$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit.....\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit.....\$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

Synergy seating continued on following page.

7. SYNERGY SEATING SYSTEM CONTINUED**6. Optional Base Insert**

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

QTY _____

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

QTY _____

Additional Stealth Tek Cover

Reference Part: XRF000394

QTY _____

Additional Mesh Cover

Reference Part: XRF000393

QTY _____

7.8. Back Cane Selection**REQUIRED.**

Synergy Back Cane No Charge

Part: FRMASMB15845

Select a back cane height. Back width will be the same as seat width.

2" Endomorph (1" per side) Back Cane \$205

Part: FRMASMB12115

Select a cane height. Back width will be 2" smaller than seat width.

2" Mesomorph (1" per side) Back Cane \$205

Part: FRMASMB12067

Select a cane height. Back width will be 2" larger than seat width.

4" Endomorph (2" per side) Back Cane \$205

Part: FRMASMB12139

Select a cane height. Back width will be 4" smaller than seat width.

7.9. Back Cane Height**REQUIRED EXCEPT if a TRU-Comfort back is desired.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes

Part: FRMASMB15832

17" High Canes

Part: FRMASMB15835

20" High Canes

Part: FRMASMB15838

23" High Canes

Part: FRMASMB15841

15" High Canes

Part: FRMASMB15833

18" High Canes

Part: FRMASMB15836

21" High Canes

Part: FRMASMB15839

24" High Canes

Part: FRMASMB15842

16" High Canes

Part: FRMASMB15834

19" High Canes

Part: FRMASMB15837

22" High Canes

Part: FRMASMB15840

25" High Canes

Part: FRMASMB15843

7.10. Back Cushion**REQUIRED.** Select one type of back.**1a. Synergy Back Selection**

Synergy Back Cushion \$289

Part: XRF000363

1b. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744

Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort 4 Back, Back Cane Option, 22" Width \$855

Part: TC4V22W32WBC. HCPCS: E2621

1c. Omit Back Selection

Omit Back No Charge

Part: PTOINDV1287

		Back Width																
		16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"	27"	28"	29"	30"	31"	32"
Overall Height	17"	S	S	S	S	S	-	STC	-	S	-	-	-	-	-	-	-	-
	18"	S	S	S	S	S	-	STC	-	S	-	-	-	-	-	-	-	-
	19"	S	-	S	S	S	S	STC	S	S	-	S	S	S	S	-	S	-
	20"	S	-	STC	S	S	S	STC	S	S	-	S	S	S	S	-	S	-
	21"	S	S	STC	S	S	S	STC	S	S	S	-	S	S	-	S	-	-
	22"	S	S	STC	STC	STC	S	STC	S	STC	S	S	S	S	-	S	-	-
	23"	STC	STC	S	STC	STC	S	STC	S	S	S	S	S	S	-	-	S	S
	24"	STC	STC	STC	S	STC	S	S	S	S	S	S	S	S	-	-	S	S
	25"	S	S	STC	S	STC	S	S	S	S	S	S	S	S	-	-	S	S

S Synergy STC Synergy/TRU-Comfort 4 Back

Synergy and TRU-Comfort Back Availability Chart

Synergy seating continued on following page.

7. SYNERGY SEATING SYSTEM CONTINUED**7.11. Overall Back Height****REQUIRED.** Refer to the Synergy and TRU-Comfort back availability chart for stock order back size ranges.

17" Overall Back Height Part: XRF000139	20" Overall Back Height Part: XRF000142	23" Overall Back Height Part: XRF000145	27" Overall Back Height Part: XRF000201
18" Overall Back Height Part: XRF000140	21" Overall Back Height Part: XRF000143	24" Overall Back Height Part: XRF000146	28" Overall Back Height Part: XRF000202
19" Overall Back Height Part: XRF000141	22" Overall Back Height Part: XRF000144	25" Overall Back Height Part: XRF000147	

7.12. Lateral Extensions**REQUIRED with Synergy Back and thoracic lateral selections.**Lateral Extension for Synergy Back (pair)..... \$58
Part: XRF000150**7.13. Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Bariatric, Single Post, Heavy Duty, Height Adjustable, Removable Armrests**

May interfere with swing-away laterals.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H No Charge	☐ ☐ Full Length, Bariatric Armrests, 10"-14"H No Charge
☐ ☐ Left part: FRMASMB6728. Right part: FRMASMB6729	☐ ☐ Left part: FRMASMB6730. Right part: FRMASMB6731.

1b. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.	☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H \$90 ea.
☐ ☐ Left part: XRF000249. Right part: XRF000246. HCPCS: E0973	☐ ☐ Left part: XRF000247. Right part: XRF000248. HCPCS: E0973

1c. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H \$90 ea.	☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
☐ ☐ Left part: XRF000251. Right part: XRF000252. HCPCS: E0973	☐ ☐ Left part: XRF000251. Right part: XRF000252. HCPCS: E0973

1d. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests \$103 ea.	☐ ☐ Desk Length, Cantilever Armrests \$103 ea.
☐ ☐ Left part: FRMASMB16188. Right part: FRMASMB16189	☐ ☐ Left part: FRMASMB16190. Right part: FRMASMB16191

7.14. Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Bariatric Armpads No Charge	Waterfall Armpads No Charge
Part: UPSASMB1025	Part: XRF000149
Only available with bariatric armrests.	May cause interference with seat back. Not available with bariatric armrests.
Straight Armpads No Charge	Omit Synergy Armpads No Charge
Part: XRF000148	Part: PTOINDV1373
Not available with bariatric armrests.	

7.15. Joystick Receiver Bracket**REQUIRED.**

Standard Fixed Receiver Standard	Stealth Height Adjustable Receiver \$236
Part: XRF000063	Part: XRF000064. HCPCS: K0108

8. HEADRESTS**Optional.** Select a headrest pad and bracket.**8.1a. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$179	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955	Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$179	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955	Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink \$311	Multi-Axis Removable Headrest Bracket \$174
Part: ST-UHB2D33200-Q. HCPCS: E1033	Part: ST-HMO475-17M-Q. HCPCS: E1033

8.2. Flip-down Headrest BracketStealth Flip-down Headrest Bracket (FDM350) \$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

8.3. Q Cam Lock**Optional.**Q Cam Lock \$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

9. LAP BELTS**REQUIRED.** Select a lap belt below.**9.1a. Standard Lap Belts****60" Lap Belt**

Part: ACCBELT1002

70" Lap Belt

Part: ACCBELT1003

80" Lap Belt

Part: ACCBELT1004

90" Lap Belt

Part: ACCBELT1001

9.1b. Optional Lap Belts**Stealth Positioning 60" Single Pull, Push Button, Medium.....\$50**
(7.5") Padded, 1.5" Belt (IFPB5213)

Part: ST-IFPB5213-Q

Stealth Positioning 60" Single Pull, Push Button, Extra Large....\$50
Padded (11.5"), 2" Belt (IFPB5315)

Part: ST-IFPB5315-Q

Stealth Positioning 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)

Part: ST-IFPB5314-Q

10. BASE MOUNTED FOOT PLATFORM**REQUIRED.** Select a foot platform option. A base mounted foot platform is not available with center mount foot platforms. May interfere with seat mounted legrests.**Heavy Duty, 17"Wx13"D, Foot Platform..... Standard**

Part: FRMASMB17955

Omit Base Mounted Foot Platform No Charge

Part: PTOINDV3674

Heavy Duty, 17"Wx16"D, Foot Platform..... No Charge

Part: FRMASMB17956

11. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**11.1a. Manually Adjustable Center Mount Foot Platform****1. Foot Platform Assembly****Center Mount Foot Platform.....\$263**

Part: XRF000155

2. Footplates**17"Wx13"D One Piece Foot Platform**

Part: RIG123289

17"Wx16"D One Piece Foot Platform

Part: RIGASMB7120011

3. Center Mount Calf Support

Includes 1/2" of soft visco foam.

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$158

Part: XRF000157

Large Calf Panel for Manually Adjustable, 24"Wx8"H \$237

Part: RIG123825

4. Additional Accessories**Reversed Dartex Calf Pad Cover, Large (pair) \$158**

Part: RIG133719

Large Reversed Dartex Panel Cover (pair) \$210

Part: RIG134744

Large Foam Calf Pad Cover (pair) \$83

Part: ST-CSCS0985-Q

Large Foam Panel Cover (pair)..... \$282

Part: ST-CSCP0927-Q

11.1b. Heavy Duty, Drop-In, Swing-Away Legrests

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection☐ P ☐ L ☐ R**70° HD Drop-in for 8.5"-12" Lower Extension No Charge**

Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707

☐ P ☐ L ☐ R**70° HD Drop-in for 13"-18" or 16"-21" Lower No Charge**

Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709

2. Lower Extension☐ P ☐ L ☐ R**8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**

Pair part: XRF000101. Left: XRF000102. Right: XRF000103.

☐ P ☐ L ☐ R**16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length**

Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393

☐ P ☐ L ☐ R**13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**

Pair part: XRF000104. Left: XRF000105. Right: XRF000106

3. Angle Adjustable Footplate(s) with Heel Loops☐ P ☐ L ☐ R**5"Wx8"D Footplates \$90 ea.**Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: E0951☐ P ☐ L ☐ R**6"Wx11"D \$90 ea.**Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.
HCPCS: E0951☐ P ☐ L ☐ R**6"Wx8"D \$90 ea.**Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.
HCPCS: E0951☐ P ☐ L ☐ R**8"Wx11"D \$90 ea.**Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.
HCPCS: E0951**11.2. Legrest Accessories****Wheel Bumper..... \$37**

Part: XRF000100

Legrests continued on following page.

11. LEGRESTS CONTINUED**11.3. Set Legrest Extension Length**

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)	18" L (#3492)
11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	17.5" L (#3491)	

12. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be mounted by the provider.

Synergy Seat Cup Holder No Charge Part: ACCASMB2634	Stealth Large Tablet Holder \$420 Part: ACC7120005
Oxygen Holder for Synergy Seat \$158 Not available with Synergy Manual Recline. Part: FRMASMB7030. HCPCS: E2208	Quantum Backpack \$168 Part: ACC1905853
Walker Holder \$109 Synergy part: FRMASMB7031	Grip Lap Board (20005) \$38.95 Part: ACC152992
Cane/Crutch Holder \$116 Synergy part: FRMASMB4822. HCPCS: E2207	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Rear Basket \$124 Synergy part: FRMASMB4821	Accessory Bar \$132 Part: XRF000128
Flag Holder \$35 Synergy part: ACCASMB2415	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
Saddle Bag \$40 Part: DWR1020D010	Privacy Bag \$42 Part: ACC1907570
Stealth Cell Phone/Media Device Holder \$215 Part: ACC7120006	Personal Item Pouch \$84 Part: ACC2101017

13. BRANCH PATCHES

Optional.

Air Force Patch \$21 Part: ACC140240	Coast Guard Patch \$21 Part: ACC140239	Navy Patch \$21 Part: ACC140244
Army Patch \$21 Part: ACC140241	Marines Patch \$21 Part: ACC140243	Purple Heart Patch \$21 Part: ACC140238

**14. POSITIONING COMPONENTS**

Optional. More options are available on the [Positioning Components Order Form](#).

14.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBTLTC) \$562 Part: INDPART2896. HCPCS: E1034	Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1034 Only available with back cane and omit back configurations.
--	---

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$250 Part: INDPART2663. HCPCS: E0956	Curved Pads, 4"x6", pair (TWBL-4X6) \$250 Part: INDPART2664. HCPCS: E0956
--	--

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

14.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink Flip Down Hip/Thigh Support, Left \$332 ea. Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^	UniLink Flip Down Hip/Thigh Support, Right \$332 ea. Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^
Flip Down with LEFT Full Surface Contact \$281 Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact \$281 Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956	UniLink, 4" x 12" Pad \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956
Qty _____	Qty _____
UniLink, 4" x 8" Pad \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956	
Qty _____	

Positioning components section continued on following page.

14. POSITIONING COMPONENTS CONTINUED**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link \$81 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link..... \$81 ea. Part: ST-UIL001200A000-Q
---	--

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

14.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$154 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$154 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$154 ea. Part: ST-IFPS2104-Q. HCPCS: E0960

**14.5. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor (IFCS4113) \$97 ea. Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large w/ Belt Anchor (IFCS4114) \$97 ea. Part: ST-IFCS4114-Q. HCPCS: E0960
---	--

**14.6. Stealth Swing-Away Headrest Mounts**


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Stealth, Swing-Away, Headrest Mount Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318 Part: ST-SAR2A-R-Q. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away \$492 2 Upper & 1 Lower Left Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
UniLink Facial Lateral, Swing Away, Left \$318 Part: ST-SAR2A-L-Q. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away \$434 2 Upper Part: ST-SAR4A-DBU-Q. HCPCS: E1028
UniLink Facial Lateral, Dual Swing Away \$492 2 Upper & 1 Lower Right Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028	

**2. Stealth, Swing-Away, Headrest Mount, Facial Lateral Pads and Accessories**

Egg Switch Mount (ESM) (A) \$57 ea. Part: INDPART2815	Swing-Away Joystick Mount,..... \$276 ea. 18" Rod (C) Part: ST-SUS18-Q Allows mounting of facial pads and attaches to the headrest
Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea. Part: ST-900M-Q. HCPCS: K0108	
Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea. Part: ST-920M-Q. HCPCS: K0108	

**14.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket \$199 ea. Tube Style Arms (ARMS-ES-L) Part: POS124064. HCPCS: E1028	Right, Swing-Away, Elbow Support Bracket..... \$199 ea. Tube Style Arms (ARMS-ES-R) Part: POS124066. HCPCS: E1028
---	---

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea. Part: POS124065. HCPCS: K0108	4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea. Part: POS124067. HCPCS: K0108
---	--

14.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539 (SPLE107) (pair) Part: INDPART2786. HCPCS: E0954	Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539 (SPLE109) (pair) Part: INDPART2788. HCPCS: E0954
Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539 (SPLE108) (pair) Part: INDPART2787. HCPCS: E0954	

14.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea. Part: ST-IFSH4211-Q. HCPCS: E0951+E0952	Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea. Part: ST-IFSH4213-Q. HCPCS: E0951+E0952
Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea. Part: ST-IFSH4212-Q. HCPCS: E0951+E0952	