

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum Outback Order Form with Contoured Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

FDA Class II
Medical Device*

(TRU Balance Power Positioning shown above)

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. 350 lb. weight capacity with Q-Captain seating, 15" pneumatic drive wheels, base mounted lights, and 8 amp off board charger. Dimensions of base without legrests: 28.82"Wx42"L.

Quantum Outback \$19,595

Part: Q OUTBACK. HCPCS: K0899 - Power mobility device, not coded by PDAC.

*The Quantum® Outback is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

2. BASE COLOR

REQUIRED. Select one shroud color.

Cadet Base Shroud Color No Charge Part: DGN2301319-TGY3-PCL2	Desert Camouflage Base Shroud Color No Charge Part: DGN2301319-WCM11-PCL2
Tan Base Shroud Color No Charge Part: DGN2301319-TTN1-PCL2	Tropic Lightning Camouflage Base Shroud Color No Charge Part: DGN2301319-WCM12-PCL2
Black Camouflage Base Shroud Color No Charge Part: DGN2301319-WCM9-PCL2	Jungle Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM13-PCL2
White Camouflage Base Shroud Color No Charge Part: DGN2301319-WCM10-PCL2	

3. MOTORS

REQUIRED.

5 MPH (High Torque) Motors No Charge Part: MOT2205971	6 MPH (High Speed) Motors No Charge Part: MOT2300107
--	---

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

10" Power Lift for Captain Seat \$3,725 Part: XRF000353. HCPCS: E2298/E2300^ 10" of scissor lift. 22"W x 22"D maximum seat size.	
--	--

5. OMIT MOUNTED LIGHTS

Optional.

Omit Base Mounted Lights No Charge Part: DGN2400803	
--	--

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368 Part: MEC1910365	
---	--

1. Optional Camera Mounting Bracket

Backup Camera Captain's Headrest Bracket \$21 Part: MEC1901231	Backup Camera Track Mount Bracket \$27 Part: MEC2104227
---	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105 Part: MEC2106123	
--	--



7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

Group 27, Interceptor 2, Gel Battery Assembly \$1,038 Part: BAT2401900	
---	--

8. JOYSTICK

REQUIRED. Select a type of drive control from the following sections.

8.1a. Omit Joystick

No Joystick, Setup for Q-Logic Electronics
Part: XRF000010
Requires selection of Enhanced Display

8.1b. Q-Logic 3e Joystick

Q-Logic 3e Joystick No Charge Part: CTL1910600 Can operate full light package. Not available with alternative or specialty controls.	
--	--

8.1c. Expandable Controls

Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module \$730 Part: XRF000002. HCPCS: E2377	Q-Logic 3 Expandable Controller Harness \$458 Part: XRF000264. HCPCS: E2313
---	--

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick \$840 Part: CTL167724	Q-Logic 3 EX Joystick for Full Light Package \$840 Part: CTL167725 Required with full light packages.
--	---



9. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in the image on the previous page. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
Medium Ball Handle (E) \$69 Part: ACCASMB1859. HCPCS: E2323	Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
Soft Ball Handle (F) \$80 Part: ACCASMB1856. HCPCS: E2323	

**10. JOYSTICK MOUNTING BRACKETS**

Optional. Select a mount to be used with a joystick selected in previous section.


10.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket Standard Part: FRMASMB10617	
--	--


10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319 Part: ACC1908434. HCPCS: E1032	
---	--

**10.1c. Basic Swing Away Joystick Mount**

Basic Swing Away, Q-Logic, Left Side \$245 Part: ACC2201041. HCPCS: E1032	Basic Swing Away, Q-Logic, Right Side \$245 Part: ACC2201042. HCPCS: E1032
--	---

**10.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368 Part: XRF000003. HCPCS: E1032 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Multi-Axis Mount \$126 Part: XRF000004 Allows height adjustment of joystick.
--	--

**10.2. Joystick Mount Location**

Left Side Mount Part: PTOINDV1225	Right Side Mount Part: PTOINDV1226
--------------------------------------	---------------------------------------

11. CONTOURED SEATING

REQUIRED.

11.1. Occupied Transit, WC-19

WC-19 Occupied Transit Option for Contoured Seat \$368
--

Coming soon. Part: KIT122609

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

11.2a. Quantum Captain Seat

16"W x 16"D No Charge Part: SET2005145	18"W x 20"D No Charge Part: SET2005151	20"W x 22"D No Charge Part: SET2005909
16"W x 18"D No Charge Part: SET2005146	18"W x 22"D No Charge Part: SET2005908	22"W x 18"D No Charge Part: SET2005155
16"W x 20"D No Charge Part: SET2005147	20"W x 16"D No Charge Part: SET2005152	22"W x 20"D No Charge Part: SET2005156
18"W x 16"D No Charge Part: SET2005149	20"W x 18"D No Charge Part: SET2005153	22"W x 22"D No Charge Part: SET2005157
18"W x 18"D No Charge Part: SET2005150	20"W x 20"D No Charge Part: SET2005154	

Contoured seating section continued on following pages.

11. CONTOURED SEATING CONTINUED**11.2b. Quantum Captain Seat with Solid Seat Pan**16"W x 16"D No Charge
Part: SET200539416"W x 18"D No Charge
Part: SET200539516"W x 20"D No Charge
Part: SET200539618"W x 16"D No Charge
Part: SET200539718"W x 18"D No Charge
Part: SET200539818"W x 20"D No Charge
Part: SET200539918"W x 22"D No Charge
Part: SET200591020"W x 16"D No Charge
Part: SET200540020"W x 18"D No Charge
Part: SET200540120"W x 20"D No Charge
Part: SET200540220"W x 22"D No Charge
Part: SET200591122"W x 18"D No Charge
Part: SET200540322"W x 20"D No Charge
Part: SET200540422"W x 22"D No Charge
Part: SET2005405**11.3. Seat to Floor Height**


Power adjustable seat height option is fixed at lowest seat to floor height.

Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
19.5"	20.5"	21.5"
Q Captain Seat to Floor Heights		

11.4. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit\$72
Part: POP1909006-M
14"W and smaller

Medium Wedge Positioning Kit\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit\$105
Part: POP1909006-L
14"W and smaller

Large Wedge Positioning Kit\$105
Part: POP1712340-L
15"W and larger

Contoured seating section continued on following pages.

11. CONTOURED SEATING CONTINUED**3. Cushion Top Foam**

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner				Cover	
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

11.5. Armrests

REQUIRED. Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Q-Captain, Flip-Up, Height Adjustable Armrests

Includes full length, flat arm pads.

Full Length, Left Armrest \$90 ea.
Part: MEC2201304. HCPCS: K0020

Full Length, Right Armrest..... \$90 ea.
Part: MEC2201303. HCPCS: K0020

1b. Drop-In, Height Adjustable, Armrests (T-Arm)

Includes flat arm pads.

☐ ☐ Full Length, Drop-In Armrests..... \$79 ea.
L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973

☐ ☐ Desk Length, Drop-In Armrests \$79 ea.
L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973

1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests \$79 ea.
L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973

1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H \$90 ea.
L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973

☐ ☐ Full Length, Bariatric Armrests, 10"-14"H \$90 ea.
L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973

1e. Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.
L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea.
L R Left Part: XRF000255. Right part: XRF000256. HCPCS: E0973

1f. Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.
L R Left part: XRF000257. Right part: XRF000258. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14" H \$90 ea.
L R Left part: XRF000259. Right part: XRF000260. HCPCS: E0973

1g. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 8.25". Includes desk length, straight arm pads.

☐ ☐ Pediatric, Removable Armrests \$79 ea.
L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

Contoured seating section continued on following pages.

11. CONTOURED SEATING CONTINUED**11.6. Armrest Height****REQUIRED with Bariatric, Heavy Duty, and Pediatric armrests.****2. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

3. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

11.7. Armpads**REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.****1a. Straight Armpads for Single Post Armrests**

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right No Charge Part: ARM134254

1b. Waterfall Armpads for Single Post Armrests

Waterfall, Full Length Left and Right No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads for Single Post Armrests

Falcon Gel, Full Length Left and Right \$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right \$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right \$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437

1d. Armpads for Bariatric or Heavy Duty Armrests

Bariatric Armpads No Charge Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
Straight Armpads No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests	

11.8. Quad Release Lever**Optional. Available with Single Post and Pediatric armrests.**

Quad Release Lever, Pair No Charge Part: FRMASMB10517	
--	--

12. LAP BELTS**REQUIRED. Select a lap belt below.****12.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

12.1b. Optional Lap Belts

Stealth Positioning 60" Single Pull, Push Button, Medium \$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning 60" Single Pull, Push Button, Extra Large \$50 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning 60" Single Pull, Push Button, Large \$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	60" Sure-Lok Lap Belt (Non-retractable) \$120 Required with ISO 7176-19 Compliant Transit Option. Coming soon. Part: ACCBELT1016

13. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

13.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount.....No Charge
Part: XRF000323

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" FootplatesNo Charge Long Extension, 5"x9" FootplatesNo Charge
Part: MEC1906675 Part: MEC1906677

Medium Extension, 5"x9" FootplatesNo Charge
Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158 Large Calf Pads, Basic Center Mount \$158
Part: RIG2104369 Part: RIG2202997

Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

13.1b. Heavy Duty, Drop-In Legrests

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

☐ P ☐ L ☐ R	70° HD Drop-in for 8.5"-12" Lower Extension No Charge Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707	☐ P ☐ L ☐ R	70° HD Drop-in 13"-18" or 16"-21" Lower No Charge Extension Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
--	--	--	---

2. Lower Extension

☐ P ☐ L ☐ R	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103.	☐ P ☐ L ☐ R	16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393
--	--	--	---

☐ P ☐ L ☐ R	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106
--	--

3. Angle Adjustable Footplate(s) with Heel Loops

☐ P ☐ L ☐ R	6"Wx8"D ¹ \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795. HCPCS: K0040+E0951	☐ P ☐ L ☐ R	8"Wx11"D ² \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799. HCPCS: K0040+E0951
--	---	--	--

☐ P ☐ L ☐ R	6"Wx11"D ¹ \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797. HCPCS: K0040+E0951
--	--

May cause interference on seat widths less than 18".

May cause interference on seat widths less than 20".

13.2. Legrest Accessories

Heel Loops (pair)..... \$53 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi axis footplates when noted.	Wheel Bumper..... \$37 Part: XRF000100
--	---

13.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

Legrests section continued on following page.

14. POWER LIFT ELECTRONICS**REQUIRED with 12" Power Lift for Captain Seat with iLevel.**

Single Actuator Through Controller and Interface \$2,163
 Part: XRF000335. HCPCS: E2310

15. ACCESSORIES**Optional.**

Cup Holder, Q-Captain No Charge Part: ACC2006355	Reflector Kit \$53 Part: ACC124670
Front Mounted Cell Phone/ Media Holder, X-Grip, Q-Captain \$74 Part: ACC2103691	Clothing Guard (Pair) \$237 Part: ACC125007 Available with Flip-Up armrests only. Not available with Elevating Legrests.
Rear View Mirror, Q-Captain \$69 Part: ACC2103692	Transfer Handles \$105 ea. Part: MEC1809930. HCPCS: E0705 Available with Flip-Up armrests only. Not available with Elevating Legrests.
Oxygen Holder \$118 Part: ACC2006353. HCPCS: E2208	Privacy Bag \$84 Part: ACC2005723
Cane/Crutch Holder \$110 Part: ACC2006352. HCPCS: E2207	Personal Item Pouch, Rail Mounted \$84 Part: ACC2101017
Rear Basket \$118 Part: ACC2006351	Personal Item Pouch, Arm Mounted \$84 Part: ACC2103694
Saddle Bag \$40 Part: DWR1020D010	Glove Box, Rail Mounted \$158 Part: ACC125011
Large, Black, Weather Cover \$202 Part: ACCCOVR1021	Glove Box, Arm Mounted \$158 Part: ACC2103693
Quantum Backpack \$168 Part: ACC1905853	Q-Logic EX/NE Programmer \$573 Part: ELEASMB6976
Grip Lap Board (20005) \$38.95 Part: ACC152992	Q-Logic PC Programming Disk with Cable \$445 Part: ELEASMB5216
Grip Drink Holder (20085) \$15.50 Part: ACC152993	

Aftermarket Accessory Mounting Bracket**The Quantum Outback includes 2 aftermarket accessory mounting bracket, slotted at no charge. Up to 5 additional brackets may be ordered below.**

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
 Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
 Part: ACC2405496

16. BRANCH PATCHES**Optional.**

Air Force Patch \$21 Part: ACC140240	Marines Patch \$21 Part: ACC140243
Army Patch \$21 Part: ACC140241	Navy Patch \$21 Part: ACC140244
Coast Guard Patch \$21 Part: ACC140239	Purple Heart Patch \$21 Part: ACC140238

**17. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**17.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea. Removable Mounting Hardware Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028/E1034^
UniLink, 7" Height Adjustable, \$276 ea. Removable Mounting Hardware Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	Flip Down with RIGHT Full Surface Contact \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link \$81 ea. Part: ST-UIL001200A000-Q
---	---

Positioning components section continued on following pages.

17. POSITIONING COMPONENTS CONTINUED

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
Part: PTOINDV3585

Install on Right
Part: PTOINDV3586

Do Not Install
Part: PTOINDV3719

17.2. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket, Tube\$199 ea.
Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube\$199 ea.
Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108
QTY

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
Part: POS124067. HCPCS: K0108
QTY

17.3. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

17.4. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211) \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952
QTY

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952
QTY

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952
QTY

18. LAP TRAY

18.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						S	M	ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	ML	ML	L	L	

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

 :Compatible

 :Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.