



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum R-TRAK Order Form with TRU Balance 3 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 7 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx34.8"L

R-TRAK 3MP-SS\$8,373

Part: R-TRAK 3MP-SS. HCPCS: K0861

Multiple actuator power option.

R-TRAK 3SP-SS\$8,559

Part: R-TRAK 3SP-SS. HCPCS: K0856

Single actuator power option.

R-TRAK 3S-SS.....\$8,870

Part: R-TRAK 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Account Number: _____ Date: _____ Quote Name/Marked for: _____

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

 Red Corvette Part: DGN2106050-PRD1	 Deep Purple Part: DGN2106050-PPR5	 Midnight Blue Part: DGN2106050-PBU10	 Green Machine Part: DGN2106050-PGN11
 Back in Black Part: DGN2106050-PBK1	 Raspberry Beret Part: DGN2106050-PPK4	 Electric Blue Part: DGN2106050-PBU9	 Military Green Part: DGN2106050-PGN8
 Night Train Part: DGN2106050-PBK3	 Pretty in Pink Part: DGN2106050-PPK3	 Blue Jean Part: DGN2106050-PBU11	
 Grey Street Part: DGN2106050-PSL3	 Evolution Orange Part: DGN2106050-POR7	 Lemon Crush Part: DGN2106050-PYL4	

3. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Operated Front Caster and Rear Fender Lights.... No Charge
Part: DGN2106044

Omit Front Caster and Rear Fender Lights..... No Charge
Part: DGN2106043

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

4. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21 Part: MEC1910373 Works with .5" x .5" vertical supports only.	Backup Camera Track Mount Bracket \$27 Part: MEC2104227
--	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



5. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color.

5.1. Drive Wheels

3"x14" Flat Free Black Tires with Accent No Charge Part: WHL1704058 x2	4"x14" Knobby, Flat Free Black Tires \$580 Part: XRF000368
3"x14" Pneumatic Black Tires with Accent Standard Part: PTOINDV1243	4"x14" Knobby, Pneumatic Black Tires \$580 Part: XRF000367

5.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in Black.....Standard Part: DGN164245-PBK1 x2	Raspberry Beret.....No Charge Part: DGN164245-PPK4 x2	Blue JeanNo Charge Part: DGN164245-PBU11 x2
Red Corvette.....No Charge Part: DGN164245-PRD1 x2	Pretty in PinkNo Charge Part: DGN164245-PPK3 x2	Lemon CrushNo Charge Part: DGN164245-PYL4 x2
Night Train.....No Charge Part: DGN164245-PBK3 x2	Evolution Orange.....No Charge Part: DGN164245-POR7 x2	Green MachineNo Charge Part: DGN164245-PGN11 x2
Grey StreetNo Charge Part: DGN164245-PSL3 x2	Midnight BlueNo Charge Part: DGN164245-PBU10 x2	Military GreenNo Charge Part: DGN164245-PGN8 x2
Deep PurpleNo Charge Part: DGN164245-PPR5 x2	Electric BlueNo Charge Part: DGN164245-PBU9 x2	

6. BATTERIES

Optional. The powerchair requires 2 flat top batteries and will be installed if ordered.

Group 24, Interceptor 2 Battery \$399 ea.
Part: BAT2100315. HCPCS: E2363

7. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

7.1a. Omit Joystick

No Joystick, Setup for Q-Logic Electronics

Part: XRF000010

Requires selection of Enhanced Display

7.1b. Non-Expandable Joystick

QLNES 6 Key Joystick Standard

Part: CTL2205379

Controls up to 2 actuators.

QLNES 10 Key Joystick \$420

Part: CTL2205380

Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

7.1c. Q-Logic 3 Expandable Capacity Controls

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity Joysticks

Q-Logic 3 EX Joystick\$840

Part: CTL167724

Q-Logic 3 EX Joystick for Full Light Package\$840

Part: CTL167725

Required with full light packages.

Q-Logic 3e Joystick\$735

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

**2. Q-Logic 3 Expandable Capacity Required Components**

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$730

Part: XRF000002 HCPCS: E2377

Q-Logic 3 Expandable Capacity Controller Harness..... \$458

Part: XRF000264. HCPCS: E2313

7.2. Q-Logic Magnetic Charger

Available with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60

Part: ELE2103635

7.3. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN124307

Also matches Night Train shrouds.

Red Corvette.....\$63

Part: DGN124336-PRD1

Back in Black.....\$63

Part: DGN124336-PBK1

Grey Street.....\$63

Part: DGN124336-PSL3

Deep Purple.....\$63

Part: DGN124336-PPR5

Raspberry Beret..... \$63

Part: DGN124336-PPK4

Pretty in Pink.....\$63

Part: DGN124336-PPK3

Evolution Orange.....\$63

Part: DGN124336-POR7

Midnight Blue.....\$63

Part: DGN124336-PBU10

Electric Blue.....\$63

Part: DGN124336-PBU9

Blue Jean..... \$63

Part: DGN124336-PBU11

Lemon Crush.....\$63

Part: DGN124336-PYL4

Green Machine.....\$63

Part: DGN124336-PGN11

Military Green.....\$63

Part: DGN124336-PGN8

8. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90

Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$90

Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$126

Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E) \$69

Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F) \$80

Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90

Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$104

Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104

Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge

Part: DGN124308




9. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**9.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket Standard
Part: FRMASMB10617**9.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**9.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032
Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**9.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.**9.2. Joystick Mount Location**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**10. TRU BALANCE 3® SEATING - QUANTUM INTELLIGENT SYSTEMS™****10.1. TRU Balance 3 Configuration****REQUIRED.** Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.TRU Balance 3 Static Seat \$2,410
Part: PTOINDV3442
85° to 105° of limited recline (-5/ +15) via gas strut within back cover.TRU Balance 3 Power Recline Only \$8,234
Part: PTOINDV3444. HCPCS: E1004
85° to 175° of power recline.TRU Balance 3 Power Tilt Only \$8,202
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.TRU Balance 3 Power Tilt & Recline \$12,615
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.**10.2. 10" Power Lift****Optional.** Select one option below. ^Please verify the appropriate billing code to be used for this item with all third-party payors.TRU Balance 3 Power 10" Lift \$3,725
Part: PTOINDV3775. HCPCS: E2298/E2300
10" of scissor lift. 22"W x 22"D maximum seat size.**10.3. Occupied Transit****Optional.** Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.Occupied Transit \$289
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**10.4. Seat Width****REQUIRED with seat selection.** 15" width is only available with non-matching back width. Minimum seat width is 12". Maximum seat width is 22".

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.5. Seat Width Package**REQUIRED with seat selection.**12"-16" Width Package
Part: XRF000006
Not available with Power Lift.14"-16" Width Package
Part: XRF000007
Only available with Power Lift.16"-20" Width Package
Part: XRF00000820"-24" Width Package
Part: XRF000158
Maximum width is 22".**10.6. Seat Depth****REQUIRED with seat selection.** Minimum seat depth is 12". Maximum seat depth is 22".

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.7. Seat to Floor Height**REQUIRED.**Lowest Seat to Floor
Part: PTOINDV3558

*Not available with Tilt Only or Tilt and Recline configurations.

Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560


*Not available with Lift configurations.

Low	Med	High
Static Seat/ Recline Only		
16"	17"	18"
Any Power Lift		
17.25"	18.25"	N/A
Tilt/Tilt & Recline		
N/A	18.5"	19.5"
TRU Balance 3 Seat to Floor Heights		

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**10.8. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$39

Part: POP1909006-S

14"W and smaller

Small Wedge Positioning Kit.....\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit.....\$105

Part: POP1909006-L

14"W and smaller

Large Wedge Positioning Kit.....\$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

FRG Top Foam

Reference part: XRF000381

No Top Foam

Reference part: XRF000382

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

Stealth Polyurethane Liner

Reference part: XRF000384

No Liner

Reference part: XRF000385

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Mesh Cover

Reference part: XRF000387

Stealth Tek Cover

Reference part: XRF000388

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**7. Additional Cushion Liner**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 _____ QTY Reference Part: XRF000390

Additional Stealth Polyurethane Liner
 _____ QTY Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 _____ QTY Reference Part: XRF000392

Additional Stealth Tek Cover
 _____ QTY Reference Part: XRF000394

Additional Mesh Cover
 _____ QTY Reference Part: XRF000393

10.9. Back Type Selection**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
 Part: XRF000014
 12" and 24" back width will use the original TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
 Part: XRF000015
 Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184
 Part: XRF000016
 12" and 24" back width will use the original TRU Balance 3 Back assembly. Not available with recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

10.10. Back Cane Selection**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable Back Cane..... No Charge
 Part: MEC140266

Flip Forward Back Cane \$53
 Part: XRF000273

2" Endomorph (1" per side) Back Cane..... \$74
 Part: XRF000269

4" Endomorph (2" per side) Back Cane..... \$74
 Part: XRF000271

2. Back Cane Height

16"H Cane
 Part: FRMASMB15834

17"H Cane
 Part: FRMASMB15835

20"H Cane
 Part: FRMASMB15838

21"H Cane
 Part: FRMASMB15839

24"H Cane
 Part: FRMASMB15842

25"H Cane
 Part: FRMASMB15843

10.11. Back Cushion Configuration**REQUIRED.** Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**

Sport Back..... No Charge
 Part: XRF000017
 Complete following width package and size.

Omit Back Pan/Cushion..... No Charge
 Part: PTOINDV1287
 Select a following width package.

After-market Back Interface \$535
 Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" Width \$744
 Part: SET123622. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width..... \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
 Part: TC4V10W21WBC. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$855
 Part: TC4V22W32WBC. HCPCS: E2621
 Complete following width package and size.

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**1c. ADI Back, AL Series, 4-Point, Active Contour**

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622 Part: ALBR16SA4. HCPCS: E2613	Active Contour, 18" Wide x 23" or 24" Overall Height \$622 Part: ALBR18XA4. HCPCS: E2613
Active Contour, 16" Wide x 21" or 22" Overall Height \$622 Part: ALBR16TA4. HCPCS: E2613	Active Contour, 20" Wide x 19" or 20" Overall Height \$622 Part: ALBR20SA4. HCPCS: E2613
Active Contour, 16" Wide x 23" or 24" Overall Height \$622 Part: ALBR16XA4. HCPCS: E2613	Active Contour, 20" Wide x 21" or 22" Overall Height \$622 Part: ALBR20TA4. HCPCS: E2613
Active Contour, 18" Wide x 19" or 20" Overall Height \$622 Part: ALBR18SA4. HCPCS: E2613	Active Contour, 20" Wide x 23" or 24" Overall Height \$622 Part: ALBR20XA4. HCPCS: E2613
Active Contour, 18" Wide x 21" or 22" Overall Height \$622 Part: ALBR18TA4. HCPCS: E2613	

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height \$735 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 18" Wide x 21" or 22" Overall Height \$735 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height \$735 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height \$735 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height \$735 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height \$735 Part: ALBR20XD4. HCPCS: E2615
Deep Contour, 18" Wide x 19" or 20" Overall Height \$735 Part: ALBR18SD4. HCPCS: E2615	

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height \$822 Part: ALBR18S64. HCPCS: E2615	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height \$822 Part: ALBR18X64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height \$822 Part: ALBR18T64. HCPCS: E2615	

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL18AX. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18AXX. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL16AX. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16AXX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20AXX. HCPCS: E2613

1g. Q-Back 2.0, AL Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED**1h. Q-Back 2.0, AL Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
Part: QB2AL16ES. HCPCS: E2620Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
Part: QB2AL16ET. HCPCS: E2620Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620**1i. Q-Back 2.0, AFT Series, Active Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613**1j. Q-Back 2.0, AFT Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT16DT. HCPCS: E2620AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620**1k. Q-Back 2.0, AFT Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AFT16EXX. HCPCS: E2620**10.12. Lateral Positioning Wedge****Optional.**Lateral Positioning Wedge, pair \$16 ea.
Part: SET123465 x2

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

10.13. TRU Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.Red Corvette
Part: XRF000019Grey Street
Part: XRF000023Pretty in Pink
Part: XRF000027Electric Blue
Part: XRF000030Green Machine
Part: XRF000033Back in Black
Part: XRF000021Deep Purple
Part: XRF000025Evolution Orange
Part: XRF000028Blue Jean
Part: XRF000031Military Green
Part: XRF000034Night Train
Part: XRF000022Raspberry Beret
Part: XRF000026Midnight Blue
Part: XRF000029Lemon Crush
Part: XRF000032**10.14. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003718"W-22"W Package
Part: XRF00003910"W-14"W Package
Part: XRF00003616"W-20"W Package
Part: XRF00003820"W-24"W Package
Part: XRF000040**10.15. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.12" W
(#3562)14" W
(#3563)15" W
(#3856)16" W
(#3564)17" W
(#3565)18" W
(#3566)19" W
(#3567)20" W
(#3568)21" W
(#3569)22" W
(#3570)24" W
(#3571)

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

		Tilt or Static Width												Recline Width											
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"	
	15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"	
	16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	16"	
	17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	17"	
	18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	18"	
	19"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	T	19"
	20"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	T	20"
	21"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	21"
	22"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	22"
	23"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	23"
	24"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	24"
	25"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	25"
	26"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	26"
27"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	27"	
28"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	28"	
29"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	29"	
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	30"	
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	31"	
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	32"	

S

Sport Back

Q

Q-Back

T

TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

10.16. Overall Back Height

Static Seat and Power Tilt available height range: 14"-29"H. *Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.*

Power Recline available height range: 19"-23"H. *Odd heights include 6" gap and even heights include 7" gap. See the back availability chart for stock order back size ranges.*

14" Height Part: SETHGHT1027	18" Height Part: SETHGHT1031	22" Height Part: SETHGHT1035	26" Height Part: SETHGHT1039	30" Height Part: SETHGHT1043
15" Height Part: SETHGHT1028	19" Height Part: SETHGHT1032	23" Height Part: SETHGHT1036	27" Height Part: SETHGHT1040	31" Height Part: SETHGHT1102
16" Height Part: SETHGHT1029	20" Height Part: SETHGHT1033	24" Height Part: SETHGHT1037	28" Height Part: SETHGHT1041	32" Height Part: SETHGHT1103
17" Height Part: SETHGHT1030	21" Height Part: SETHGHT1034	25" Height Part: SETHGHT1038	29" Height Part: SETHGHT1042	

10.17. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$48
Part: XRF000351

10.18. Back Height Extension

Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension \$48
Part: XRF000041

11. TRU BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

11.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and armrest height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ ☐ Flip-Back, 8"-16"H Armrests Standard	☐ ☐ Flip-Back, 8"-16"H Armrests for Static Seat \$289
☐ ☐ Left part: ARM1806662. Right part: ARM1806663	☐ ☐ Left part: XRF000320. Right part: XRF000321
☐ ☐ Not available with TRU Balance 3 Static Seat.	

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ ☐ Flip-Back Armrest Downpost \$92 ea.
☐ ☐ Left part: XRF000265. Right part: XRF000266

11.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Only available with TRU Balance 3 Static Seat. Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests No Charge	☐ ☐ Desk Length, Drop-In Armrests No Charge
☐ ☐ Left part: ARM1803253. Right part: ARM156742	☐ ☐ Left part: ARM1803254. Right part: ARM156743

TRU Balance 3 armrest section continued on following page.

11. TRU BALANCE 3 ARMRESTS CONTINUED**11.1c. Two-Post, Height Adjustable, Removable Armrests**

Select the armrest and the left and right armrest heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ ☐ Full Length, Two-Post Armrest for Power Seating No Charge
 Left part: XRF000047. Right part: XRF000048

☐ ☐ Desk Length, Two-Post Armrest for Power Seating.... No Charge
 Left part: XRF000049. Right part: XRF000050

1b. Two-Post Armrest for Static Seating

☐ ☐ Full Length, Two-Post Armrest for Static Seating.....\$158 ea.
 Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrest for Static Seating\$158 ea.
 Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

11.1d. Heavy Duty, Removable Armrests

HD removable armrests require armrest height and armrest size selections. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$79 ea.
 Heavy Duty (Sold as pair)
 Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$79 ea.
 Heavy Duty (Sold as pair)
 Part: XRF000056. HCPCS w/ Static Seat: E0973

11.1e. Pediatric, Removable, Height Adjustable ArmrestsHeight adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armrests.****1. Pediatric Armrest Selection**

☐ ☐ Pediatric, Removable Armrests \$79 ea.
 Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

2. Pediatric, Left Armrest Height

7.75" H
PTOINDV3722

8.25" H
PTOINDV3724

3. Pediatric, Right Armrest Height

7.75" H
PTOINDV3721

8.25" H
PTOINDV3723

11.2. Armrest Height

Do not select with Pediatric Armrests. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H
PTOINDV3535

10" H
PTOINDV3537

12" H
PTOINDV3539

14" H
PTOINDV3541

16" H
PTOINDV3543

9" H
PTOINDV3536

11" H
PTOINDV3538

13" H
PTOINDV3540

15" H
PTOINDV3542

2. Right Armrest Height

8" H
PTOINDV3450

10" H
PTOINDV3452

12" H
PTOINDV3454

14" H
PTOINDV3456

16" H
PTOINDV3458

9" H
PTOINDV3451

11" H
PTOINDV3453

13" H
PTOINDV3455

15" H
PTOINDV3457

11.3. Armrest Size

Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ ☐ 10" Length
 Left part: XRF000057. Right part:
 XRF000123

☐ ☐ 14" Length
 Left part: XRF000058. Right part:
 XRF000124

☐ ☐ 17" Length
 Left part: XRF000059. Right part:
 XRF000125
 Only available with Flip-Back, 8"-16" H
 Armrests and straight arm pads.

12. TRU BALANCE 3 ARMPADS**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.**12.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge
 Part: XRF000060
 Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge
 Part: XRF000308
 Only available with Flip-back, 8"-16" armrests.

12.1b. Foam Armpads

☐ ☐ Straight, Armpad Standard
 Left part: XRF000061. Right part: XRF000126

☐ ☐ Waterfall, Armpad No Charge
 Left part: XRF000062. Right part: XRF000127
 Not available with Flip-back, 8"-16" armrests.

TRU Balance 3 armpads section continued on following page.

Account Number: _____ Date: _____ Quote Name/Marked for: _____

12. TRU BALANCE 3 ARMPADS CONTINUED

12.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195
Part: ARM1806807

Falcon Gel, Full Length Left, Desk Length Right\$171
Part: ARM1806813

Falcon Gel, Desk Length Left, Full Length Right\$171
Part: ARM1806810

Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438

Falcon Gel, Full Length Left, Desk Length Right\$287.93
Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right\$287.93
Part: ARM136439

Falcon Gel, Desk Length Left and Right\$265.02
Part: ARM136437

13. TRU BALANCE 3 JOYSTICK BRACKETS

REQUIRED.

Standard Fixed Receiver Standard
Part: XRF000063

Stealth Height Adjustable Receiver \$236
Part: XRF000064. HCPCS: K0108

14. TRU BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

14.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM)\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033

Multi-Axis Removable Headrest Bracket\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033

14.2. Q Cam Lock

Optional.

Q Cam Lock\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

15. LAP BELTS

REQUIRED. Select a lap belt below.

15.1a. Standard Lap Belts

38" Lap Belt
Part: ACCBELT1006

50" Lap Belt
Part: ACCBELT1000

60" Lap Belt
Part: ACCBELT1002

70" Lap Belt
Part: ACCBELT1003

80" Lap Belt
Part: ACCBELT1004

90" Lap Belt
Part: ACCBELT1001

15.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350)\$181
Part: ST-IFRB5350

Stealth Positioning 60" Single Pull, Push Button, Medium\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

Stealth Positioning 60" Single Pull, Push Button, Extra Large\$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable) \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

16. POWER LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)
Part: XRF000108

Combined Legrests (Requires 2 actuator control)
Part: XRF000109

Articulating Foot Platform (Requires 1 Actuator Control)
Part: XRF000110

17. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

17.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates No Charge
Part: MEC1906675

Long Extension, 5"x9" Footplates No Charge
Part: MEC1906677

Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997

Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

17.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311. HCPCS: K0040

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577

Medium Extension, 5"x8" Footplates
Part: RIG2404016

Long Extension, 5"x8" Footplates
Part: RIG2404017

Pediatric Extension, 5"x8" Footplates
Part: RIG2404015

Medium Extension, 5"x11" Footplates
Part: RIG2404018

Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium Extension, 5"x6" Footplates
Part: RIG150206

Long Extension, 5"x6" Footplates
Part: RIG151505

3a. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

Small Calf Pads, Quick Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156

Large Calf Pads, Quick Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

Calf pads are standard with a power articulating foot platform.

Small Calf Panel, 18"Wx8"H..... \$184
Part: FRMASMB7110012

Large Calf Panel, 24"Wx8"H \$237
Part: RIG123825

Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.

Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.

Width Extension for 8" \$69
Part: RIG2404919

Front and Rear Cap \$16
Part: RIG2404923

Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" footplates.

Legrests section continued on following pages.

17. LEGRESTS CONTINUED**5. Additional Accessories**

Small Foam Calf Pad Cover (pair).....	\$67
Part: ST-CSCS0775-Q	
Large Foam Calf Pad Cover (pair)	\$83
Part: ST-CSCS0985-Q	
Small Foam Panel Cover (pair)	\$176
Part: ST-CSCP0920-Q	
Large Foam Panel Cover (pair).....	\$282
Part: ST-CSCP0927-Q	
Footplate Connector	\$16
Part: RIGASMB7120023	

Footplate Riser	\$75
Part: XRF000348	
Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.	
Footplate Riser Calf Panel	\$50
Part: FRMASMB7110011	
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers.....	\$222
Part: ST-FPG5008T-Q. HCPCS: K0108	
5"x11" Tapered Gel Footplate Covers.....	\$270
Part: ST-AFPG5011T-Q. HCPCS: K0108	

5"x8" Tapered Gel Footplate Covers for Width Extension	\$247
Part: ST-AFPG5008E-Q. HCPCS: K0108	
5"x11" Tapered Gel Footplate Cover for Width Extension	\$270
Part: ST-AFPG5011E-Q. HCPCS: K0108	

17.1c. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair	\$2,516
Part: XRF000341. HCPCS: E1010+K0040+E0951	

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Static Seat, Recline Only			
8.5"-14.5" Lower Length	8.5"-12.5"	8.5"-13.5"	8.5"-14.5"
12.5"-21" Lower Length	N/A	12.5"-13.5"	12.5"-14.5"

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
12.5"-21" Lower Length	12.5"-15.5"	12.5"-15.5"	12.5"-17"

2. Lower Extension

8.5" to 14.5" Lower Length	12.5" to 21" Lower Length
Part: XRF000344	Part: XRF000347

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair	5"Wx8"D, Pair	6"Wx11"D ¹ , Pair
Part: FRMASMB7448	Part: FRMASMB7450	Part: FRMASMB7452
5"Wx6"D, Pair	6"Wx8"D ¹ , Pair	8"Wx11"D ² , Pair
Part: FRMASMB7449	Part: FRMASMB7451	Part: FRMASMB7453

17.1d. 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

☐ P ☐ L ☐ R	60°-90° Swing-Away Legrest, Power Seating..... No Charge
	Pair part: RIG120062. Left: FRMASMB17558.
	Right: FRMASMB17559
	Available with Power Positioning only. Not available with 11"-13" extensions.
☐ P ☐ L ☐ R	60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge
	Pair part: RIG136481. Left: RIG132607.
	Right: RIG132609.
	Available with Power Positioning only. For use only with 11"-13" lower extensions.
☐ P ☐ L ☐ R	60°-90° Swing-Away Legrest, Static Seat \$86 ea.
	Pair part: XRF000082. Left: XRF000083.
	Right: XRF000084
	Available with Static Seat only. Not available with 11"-13" extensions.
☐ P ☐ L ☐ R	60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
	Pair part: XRF000085. Left: XRF000086.
	Right: XRF000087.
	Available with Static Seat only. For use only with 11"-13" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

17. LEGRESTS CONTINUED**2. Swing-Away Legrests Lower Extension**

11"-13" Length
 P L R Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

13"-16.5" Length
 P L R Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

15.5"-19.5" Length
 P L R Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

5"Wx8"D with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

6"Wx11"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

6"Wx8"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

8"Wx11"D² with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. 60°-90° Swing-Away Legrest Accessories

Gel Neoprene Footrest Sleeve (pair) \$84
 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

17.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

Elevating Legrest Assembly \$171 ea.
 P L R Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

11"-14" Length
 P L R Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 P L R Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

14"-17" Length
 P L R Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 P L R Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

7.25"Wx6"D Reinforced Footplate..... Standard
 P L R Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

8.5"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

5"Wx8"D with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

6"Wx11"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

6"Wx8"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

8"Wx11"D² with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

 This rigging configuration causes no interference issues at this seat-to-floor height.
 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

17. LEGRESTS CONTINUED**17.2. Legrest Accessories**

Heel Loops (pair).....\$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper.....\$37
 Part: XRF000100

17.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options.

18.1a. Actuator Control through Joystick**1. Actuator Control Selection**

Single Actuator Through Controller.....\$2,163
 Interface
 Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
 Interface
 Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and \$3,040
 Interface (AAM)
 Part: XRF000068. HCPCS: E2310

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations.

4-Function Push Button
 Part: CTL123218

4-Way Toggle
 Part: ELEASMB5925

4-Way Button Toggle
 Part: ELEASMB5926

See section 21.2 for switch mounting.

3. User Controls Side Selection

Required with 4-Way Toggle and 4-Way Button Toggle.

Mount User Controls on Left
 Part: PTOINDV3996

Mount User Controls on Right
 Part: PTOINDV3997

**18.1b. Actuator Control with iAccess Switch**

Single Actuator Control Through Flat Button\$2,547
 iAccess and Controller and Interface
 Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button\$2,547
 iAccess and Controller and Interface
 Part: XRF000338. HCPCS: E2310

Single Actuator Control Through Flat Button.....\$3,460
 iAccess and Controller and Interface (AAM)
 Part: XRF000070. HCPCS: E2310

Single Actuator Control Through Toggle Button\$3,460
 iAccess and Controller and Interface (AAM)
 Part: XRF000074. HCPCS: E2310

Multiple Actuator Control Through Flat Button\$3,460
 iAccess and Controller and Interface
 Part: XRF000339. HCPCS: E2311

Multiple Actuator Control Through Toggle Button.....\$3,460
 iAccess and Controller and Interface
 Part: XRF000340. HCPCS: E2311

18.2. 4-Function Push Button/iAccess Mount

REQUIRED with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge
 Part: MEC1907500

Under Arm Mount, Flip-Back Armrest..... No Charge
 Part: MEC1906531

Under Arm Mount, Tube Style Armrest No Charge
 Part: MEC1906530

Attendant Mount.....\$26
 Part: MEC1906564

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard.....\$26
 Part: MEC2001392

Under Arm Mount, Tube Style Armrest, plus iAccess Guard\$26
 Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount
 Part: PTOINDV3678

Right Side Under Arm Mount
 Part: PTOINDV3679

19. TRU BALANCE 3 ATTENDANT CONTROL

Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525
 Part: XRF000078. HCPCS: E2331
 Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
 Controls and Q-Logic Display or SCIM Module
 Part: XRF000079

20. SEAT MOUNTED LED LIGHTS

Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045
 Part: XRF000130
 Not available with Base Mounted, LED, Full Light Package.

21. PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU Balance 3 Cane/Crutch Holder \$105
 Part: ACC152515. HCPCS: E2207

Quantum Backpack \$168
 Part: ACC1905853

Grip Lap Board (20005) \$38.95
 Part: ACC152992

Grip Drink Holder (20085) \$15.50
 Part: ACC152993

TRU Balance 3 Hydration System (B) \$237
 Part: ACC125010

TRU Balance 3 Transfer Handle (D) \$105 ea.
 Part: MEC1809930. HCPCS: E0705
 Not available with Heavy Duty armrests.

TRU Balance 3 Cup Holder No Charge
 Part: ACC125003
 For use with TRU Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest) No Charge
 Part: ACCASMB2634
 For use with tube-style armrests.

TRU Balance 3 Personal Item Hook \$69
 Part: ACC125008
 Not available with back cane option.

TRU Balance 3 Reflector Kit \$53
 Part: ACC124670

TRU Balance 3 Backpack Holder (Pair) (A) \$69
 Part: ACC147774
 Not available with back cane option.

TRU Balance 3 Front Mounted \$69
 Cell Phone/ Media Holder, X-Grip (E)
 Part: ACC155835
 Not available with Heavy Duty armrests.

TRU Balance 3 Rear View Mirror \$69
 Part: ACC155833
 Not available with Heavy Duty armrests.

TRU Balance 3 Oxygen Holder \$210
 Part: ACC125006. HCPCS: E2208
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

TRU Balance 3 Clothing Guard (Pair) \$237
 Part: ACC125007

Not available with Heavy Duty armrests.

TRU Balance 3 Glove Box (C) \$158
 Part: ACC125011

Stealth Large Tablet Holder \$625.28
 Part: ACC141818

Not available with Heavy Duty armrests.

Compact Vent Tray \$1,465
 Part: ACC133889. HCPCS: E1029

Push Handles (pair) \$195
 Part: ACC2005027
 Not available with oxygen holder, back cane option, or vent tray.

Accessory Bar \$132
 Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the oxygen holder. Not available with back cane option.

Large, Black, Weather Cover \$202
 Part: ACCCOVR1021

Q-Logic EX/NE Programmer \$573
 Part: ELEASMB6976

Flip-Back Armrest, Side Track Accessory Bracket ... \$11
 Part: MEC1809178


Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.

Privacy Bag \$42
 Part: ACC2005723

Personal Item Pouch \$84
 Part: ACC2101017

Depth Adjustable Vent Tray \$595
 Part: ACC2302500. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOIN-DV3769.

Remote Stop Switch \$1,295
 Part: XRF000379



Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
 Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
 Part: ACC2405496

22. BRANCH PATCHES

Optional.

Air Force Patch \$21
 Part: ACC140240

Coast Guard Patch \$21
 Part: ACC140239

Navy Patch \$21
 Part: ACC140244

Army Patch \$21
 Part: ACC140241

Marines Patch \$21
 Part: ACC140243

Purple Heart Patch \$21
 Part: ACC140238



23. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664

Part: ST-ULS4042000-Q. HCPCS: E1034

Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664

Part: ST-ULPK404108B00-Q. HCPCS: E1034

Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664

Part: ST-ULPK404104B00-Q. HCPCS: E1034

Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332

Part: ST-ULPK404103L00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332

Part: ST-ULPK404103R00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250

Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250

Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals

Part: PTOINDV3499

Do Not Install

Part: PTOINDV3500

23.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.

Removable Mounting Hardware

Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.

Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.

Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.

Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.

Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.

Feature (TWBADD-L & LPHW407)

Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.

Feature (TWBADD-R & LPHW407)

Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad..... \$106 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$106 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.

Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.

Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left

Part: PTOINDV3585

Install on Right

Part: PTOINDV3586

Do Not Install

Part: PTOINDV3719

23.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.

Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.

Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.

Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.

Part: ST-IFPS2104-Q. HCPCS: E0960

23.4. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor \$97 ea.

(IFCS4113)

Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$97 ea.

(IFCS4114)

Part: ST-IFCS4114-Q. HCPCS: E0960



UniLink Pad & Hardware



Contoured and Classic Chest Harness



Chest Strap

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.5. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right..... \$318
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left \$318
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Left


Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$434
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and AccessoriesEgg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea.
Part: ST-920M-Q. HCPCS: K0108Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

23.6. Stealth Elbow Support**1. Elbow Support Bracket**Left, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Left, Fixed Elbow Support Bracket \$199 ea.
TRU Balance 3 Arms (ARMS-ESFTB3-L)
Part: POS133465Left, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-QRight, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Right, Fixed Elbow Support Bracket \$199 ea.
TRU Balance 3 Arms (ARMS-ESFTB3-R)
Part: POS133466Right, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108

Elbow Supports

23.7. Stealth Foot BoxesStealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954**23.8. Stealth Positioning Shoe Holder**Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Quantum Lap Tray****Mounting bracket is required with lap tray selection.** These trays are compatible with Flip-Back armrests.**1a. Black Lap Tray**

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

24.2. Simple Solutions Lap Tray**Mounting bracket is required with lap tray selection.** See chart below for compatibility.**1. Lap Tray**

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**25.1. Enhanced Display**Q-Logic 3 Enhanced Display Kit.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Grey Street \$63 Part: DGN137729-PSL3	Evolution Orange..... \$63 Part: DGN137729-POR7	Lemon Crush..... \$63 Part: DGN137729-PYL4
Red Corvette..... \$63 Part: DGN137729-PRD1	Deep Purple \$63 Part: DGN137729-PPR5	Midnight Blue \$63 Part: DGN137729-PBU10	Green Machine \$63 Part: DGN137729-PGN11
Back in Black..... \$63 Part: DGN137729-PBK1	Raspberry Beret..... \$63 Part: DGN137729-PPK4	Electric Blue \$63 Part: DGN137729-PBU9	Military Green \$63 Part: DGN137729-PGN8
Night Train..... \$63 Part: DGN137729-PBK3	Pretty in Pink \$63 Part: DGN137729-PPK3	Blue Jean \$63 Part: DGN137729-PBU11	

25.2. Q-Logic SCIM Module**Optional.**SCIM Module Kit for Q-Logic 3.....\$552
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 Expandable Capacity Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJQ-Logic 3 EX SAJ\$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

25.4. Head Arrayi-Drive 4.0 Bluetooth Tri-array Series,..... \$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330Junior, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330Adult, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330

Tri-Array Series

25.5. Sip and Puff Drive**1. Module**Sip and Puff Module for Q-Logic 3.....\$1,990
Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.**Stealth Sip and Puff for Comfort Series, \$318
Left (SP625M-CL)
Part: ST-SP625M-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.Stealth Sip and Puff for Stealth i-Drive, \$318
Head Array (SP625-U)
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.Stealth Sip and Puff for Comfort Series, \$318
Right (SP625M-CR)
Part: ST-SP625M-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Sip & Puff Module



Stealth Sip & Puff

3. Sip & Puff AccessoriesSip & Puff Straws, 12 pack (SP625-S12).....\$35
Part: ST-SP625-S12-Q

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

25.6. Egg Switches

<u>QTY</u>	Red Egg Switch (ESRED)..... \$87 Part: SWTMCRO1077	<u>QTY</u>	Green Egg Switch (ESGRN) \$87 Part: SWTMCRO1078
<u>QTY</u>	Blue Egg Switch (ESBLU) \$87 Part: SWTMCRO1079	<u>QTY</u>	Yellow Egg Switch (ESYEL) \$87 Part: SWTMCRO1080
<u>QTY</u>	Black Egg Switch (ESBLK) \$87 Part: SWTMCRO1081		



25.7. Switch Adapter

<u>QTY</u>	2 to 1 Switch Adapter..... \$15 Part: DWR1361H037 Connects 2 single switches to a stereo jack to act as a dual switch.
------------	--

25.8. Environmental Controls

<u>QTY</u>	Q-Logic 3 ECU Module Kit \$2,520 Part: CTL168057
------------	---