



QUANTUM®
#1 FOR REHAB POWER

Quantum Rehab®, 401 York Ave., Duryea, PA 18642
Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Quantum Outback Order Form with Q Static Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



FDA Class II
Medical Device*

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, option numbers, and availability are subject to change without notice. Prices and option numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT DETAILS

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

BASE MODEL

REQUIRED. 300 lb. weight capacity with TRU-Balance 3 positioning, 15" pneumatic drive wheels, and 10 amp off board charger. Dimensions of base without legrests: 28.82"Wx42"L.

OPTION NUMBER	DESCRIPTION	NOTES	HCPCS	MSRP
Q OUTBACK	Quantum Outback	*The Quantum® Outback is an FDA Class II Medical Device designed to aid individuals with mobility impairments.	K0899	\$19,595

Account Number:

Date:

Quote Name/Marked for:

DRIVE MOTORS**Required.**

10037	5 mph Motors	High torque motors.	No Charge
10034	6 mph Drive Motors	High speed motors.	No Charge

BASE COLOR**REQUIRED. Select one.**

10026	Cadet Base Color		No Charge
10027	Tan Base Color		No Charge
10028	White Camouflage Base Color		No Charge
10029	Desert Camouflage Base Color		No Charge
10032	Black Camouflage Base Color		No Charge

BASE MOUNTED LIGHTS**Optional.**

10040	Omit Base Mounted Lights	Accent lights operate by pressing included single switch button.	No Charge
-------	--------------------------	--	-----------

QUANTUM BACKUP CAMERA**Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.**

10041	Quantum Backup Camera		\$375
-------	-----------------------	--	-------

OPTIONAL CAMERA MOUNTING BRACKET

10042	Backup Camera Headrest Bracket	Works with .5" x .5" vertical supports only.	\$21
10043	Backup Camera Track Mount Bracket		\$28

OPTIONAL DISPLAY MOUNTING BRACKET

10044	Backup Camera Display, Joystick Mount Bracket	Only available with Q Logic EX Joystick.	\$107
-------	---	--	-------

**DRIVE WHEELS****REQUIRED.**

10045	Knobby, Flat Free Wheel	4"x14" Tires	No Charge
10046	Knobby, Pneumatic Wheel	4"x14" Tires	No Charge

BATTERIES**Optional. The battery assembly includes 2 batteries and will be installed if ordered.**

10075	Group 27, Interceptor 2, Gel Battery Assembly		K0108	\$1,038
-------	---	--	-------	---------

ELECTRONICS**REQUIRED. Select one drive control. The No Joystick option is used to omit the joystick and utilize a specialty control, like those on the Alternative Controls and Electronics Order Form.**

10081	Q-Logic 3e Joystick	Not available with alternative or specialty controls. Includes Q-Logic Magnetic Charger below. Operates Full Lights.	No Charge
10079	Q-Logic 3 Expandable Capacity Joystick	(EX) Includes XLR USB charger.	\$857
10080	Q-Logic 3 Expandable Capacity Joystick for Full Light Package	(EX) Includes XLR USB charger. Operates Full Lights.	\$857
XRF000010	No Joystick, Setup for Q-Logic	Requires specialty controls selection.	No Charge

The Q-Logic 3 Controller Module and Controller Harness options below are required with (EX) electronics selections above. If a Q-Logic specialty control is selected elsewhere, Customer Service will quote the order with these items.

XRF000002	Q-Logic 3 Controller Module	E2377	\$745
XRF000264	Q-Logic 3 Expandable Controller Harness	E2313	\$467



Account Number:

Date:



Quote Name/Marked for:

JOYSTICK HANDLES**Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.**

ACCASMB1860	Mushroom Handle	E2323	\$155
ACCASMB1861	Chin Cup	E2323/ E2324	\$155
ACCASMB1858	T-Bar	E2323	\$195
ACCASMB1859	Medium Ball Handle	E2323	\$155
ACCASMB1856	Soft Ball Handle	E2323	\$174
ACCASMB1857	Extended I-Handle	E2323	\$174
PLSKNCP1057	Bodypoint 3" Goal Post Handle	PC101B	E2323 \$179
PLSKNCP1058	Bodypoint 4" Goal Post Handle	PC102B	E2323 \$179

**JOYSTICK MOUNTING BRACKETS****Optional. Select a mount to be used with a joystick selected in previous section.**

10109	Fixed, In-line, Bracket	No Charge
10110	Swing-away, Adjustable, Joystick Bracket	Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. E1032 \$325
10111	Basic Swing-away Joystick Bracket	E1032 \$245
XRF000003	Retract4 Mount	Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. E1032 \$375
XRF000004	Multi-Axis Mount	Allows height adjustment of joystick. \$129

**JOYSTICK MOUNT LOCATION**

PTOINDV1225	Mount Controller on Left	No Charge
PTOINDV1226	Mount Controller on Right	No Charge

Q STATIC SEAT**Q STATIC SEAT CONFIGURATION****REQUIRED. Maximum seat size is 22"W x 22"D. 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.**

10112	Q Static Basic	Excludes Power Lift, Center Spine, Flip Forward Back, Aftermarket Back, Flip Back Armrests, AFP, Power ELR, and Custom options.	No Charge
10114	Q Static Positioning		\$2,458

10" POWER LIFT**Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.**

10116	Q Static 10" Lift	Available with Q Static Positioning. 10" of scissor lift with iLevel System. Requires seat function selection above.	E2298/ E2300^ \$3,725
-------	-------------------	--	--------------------------

OCCUPIED TRANSIT

Optional. Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU Balance 3 armrests and 60" Sure-lok lap belt. Only available with swingaway legrests and center mount foot platforms. Omitting a backrest selection requires a provider to facilitate an aftermarket backrest that best fits the consumer's needs. The third-party backrest if applicable must be occupied transit approved per ISO/RESNA standards. When modifying a power chair, including the addition of third-party componentry, to accommodate the medical needs of the power chair occupant, a risk assessment in conformance with ISO 14971, as outlined in ISO 7176, should be performed.

10120	Occupied Transit	ISO occupied transit.	\$295
-------	------------------	-----------------------	-------

Account Number:

Date:

Quote Name/Marked for:

SEAT WIDTH**Required.**

SETWDTH1004	14" Seat Width	No Charge	SETWDTH1009	19" Seat Width	No Charge
SETWDTH1005	15" Seat Width	No Charge	SETWDTH1010	20" Seat Width	No Charge
SETWDTH1006	16" Seat Width	No Charge	SETWDTH1011	21" Seat Width	No Charge
SETWDTH1007	17" Seat Width	No Charge	SETWDTH1012	22" Seat Width	No Charge
SETWDTH1008	18" Seat Width	No Charge			

SEAT WIDTH PACKAGE**REQUIRED with seat selection. Maximum width is 22".**

XRF000007	14"-16" Width Package	No Charge	10385	20"-22" Width Package	No Charge
XRF000008	16"-20" Width Package	No Charge			

SEAT DEPTH**REQUIRED with seat selection. Maximum seat depth is 22".**

SETDPTH1004	14" Seat Depth	No Charge	SETDPTH1009	19" Seat Depth	No Charge
SETDPTH1005	15" Seat Depth	No Charge	SETDPTH1010	20" Seat Depth	No Charge
SETDPTH1006	16" Seat Depth	No Charge	SETDPTH1011	21" Seat Depth	No Charge
SETDPTH1007	17" Seat Depth	No Charge	SETDPTH1012	22" Seat Depth	No Charge
SETDPTH1008	18" Seat Depth	No Charge			


SEAT TO FLOOR HEIGHT**REQUIRED.**

PTOINDV3558	Lowest Seat to Floor	No Charge
PTOINDV3559	Middle Seat to Floor	No Charge
PTOINDV3560	Highest Seat to Floor	No Charge

Low	Med	High
Power Lift		
19.5"	20.5"	21.5"
Q Static		
20"	21"	22"
Seat to Floor		

STEALTH SEAT CUSHION

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete Stealth Cushion Order Forms. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**CUSHION TYPE**

Ref: XRF000116	TRU-Comfort 2 SPP Seat Cushion	Part series: TC2-###XXYY. HCPCS: E2607/E2608. MSRP: \$441-\$651.
Ref: XRF000118	Essence SPP Seat Cushion	Part series: ESS-###XXYY. HCPCS: E2607/E2608. MSRP: \$441-\$630.
Ref: XRF000113	Solution SPP Seat Cushion	Part series: SOL-###XXYY. HCPCS: E2607/E2608. MSRP: \$414-\$533.
Ref: XRF000112	Spectrum Foam SPP Seat Cushion	Part series: SPCF-###XXYY. HCPCS: E2607/E2608. MSRP: \$507-\$709.
Ref: XRF000117	Spectrum Gel SPP Seat Cushion	Part series: SPCG-###XXYY. HCPCS: E2607/E2608. MSRP: \$507-\$602.
Ref: XRF000122	Premiere P Seat Cushion	Part series: PREM-###XXYY. HCPCS: E2605/E2606. MSRP: \$358-\$603.
Ref: XRF000114	Glacial SP Seat Cushion	Part series: GLA-###XXYY. HCPCS: E2603/E2604. MSRP: \$364-\$564.
Ref: XRF000115	Zen SP Seat Cushion	Part series: ZEN-###XXYY. HCPCS: E2603/E2604. MSRP: \$348-\$547.
Ref: XRF000111	Simplicity G Seat Cushion	Part series: SIM-###XXYY. HCPCS: E2601/E2602. MSRP: \$122-\$182.

Account Number:

Date:

Quote Name/Marked for:

ESSENCE CUSHION WEDGE KIT**Only available with Essence cushion selection above.**

POP1909006-S	Small Wedge Positioning Kit	14"W and smaller
POP1712340-S	Small Wedge Positioning Kit	15"W and larger
POP1909006-M	Medium Wedge Positioning Kit	14"W and smaller
POP1712340-M	Medium Wedge Positioning Kit	15"W and larger
POP1909006-L	Large Wedge Positioning Kit	14"W and smaller
POP1712340-L	Large Wedge Positioning Kit	15"W and larger

CUSHION TOP FOAM

Ref: XRF000380	Visco Top Foam
Ref: XRF000381	FRG Top Foam
Ref: XRF000382	No Top Foam

CUSHION LINER

Ref: XRF000383	Stealth Urethane Liner
Ref: XRF000384	Stealth Polyurethane Liner
Ref: XRF000385	No Liner

CUSHION COVER

Ref: XRF000386	CoolCore Cover
Ref: XRF000387	Mesh Cover
Ref: XRF000388	Stealth Tek Cover

OPTIONAL BASE INSERT**Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.**

Ref: XRF000389	ABS Flat Rigidizer	MSRP: \$48-\$69
----------------	--------------------	-----------------

ADDITIONAL CUSHION LINER**Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.**

Ref: XRF000390	Additional Stealth Urethane Liner	Quantity: ____
Ref: XRF000391	Additional Stealth Polyurethane Liner	Quantity: ____

ADDITIONAL CUSHION COVER**Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.**

Ref: XRF000392	Additional CoolCore Cover	Quantity: ____
Ref: XRF000393	Additional Mesh Cover	Quantity: ____
Ref: XRF000394	Additional Stealth Tek Cover	Quantity: ____

BACK TYPE**REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.**

XRF000014	Center Seat Spine	Not available with Q Static Basic. Select a back cushion configuration in the subsequent section.	No Charge
XRF000015	Back Cane Option	Select cane type, back cane height, and back cushion configuration below.	No Charge
XRF000016	Flip Forward Back	Not available with Q Static Basic. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.	\$188

BACK CANE SELECTION**REQUIRED when Back Cane Option is selected above. Select type and height of back canes.****BACK CANE TYPE**

10123	Angle Adjustable Back Cane		No Charge
XRF000274	Flip Forward Back Cane	Available with Q Static Basic	\$188
XRF000273	Flip Forward Back Cane	Available with Q Static Positioning	\$54
XRF000270	2" Endomorph (1" per side) Back Cane	Available with Q Static Basic	\$209
XRF000272	4" Endomorph (2" per side) Back Cane	Available with Q Static Basic	\$209

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	-	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

Account Number:

Date:

Quote Name/Marked for:

XRF000269	2" Endomorph (1" per side) Back Cane	Available with Q Static Positioning	\$75
XRF000271	4" Endomorph (2" per side) Back Cane	Available with Q Static Positioning	\$75

BACK CANE HEIGHT

10127	16" Back Cane Height	No Charge
10128	17" Back Cane Height	No Charge
10131	20" Back Cane Height	No Charge
10132	21" Back Cane Height	No Charge
10135	24" Back Cane Height	No Charge
10136	25" Back Cane Height	No Charge

BACK CUSHION CONFIGURATION

REQUIRED. Select one type of back configuration. See the back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

SPORT BACK/OMIT/AFTER-MARKET INTERFACE

XRF000361	Sport Back Cushion, 21" Wide or less	Available with Q Static Basic. Complete following width package and size.	\$402
XRF000362	Sport Back Cushion, 22" Wide or more	Available with Q Static Basic. Complete following width package and size.	\$616
XRF000017	Sport Back	Available with Q Static Positioning. Complete following width package and size.	No Charge
PTOINDV1287	Omit Back Pan/Cushion	Select a following width package.	No Charge
ST-SSAPT4-Q	After-market Back Interface	Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.	\$535

TRU-COMFORT BACK

TC4V10W21W	TRU-Comfort Back, Less than 22" W	Center Spine or Flip Forward Back compatibility. Complete following width package and size.	E2620	\$759
TC4V10W21WBC	TRU-Comfort Back, Back Cane Option, Less than 22" W	Back Cane compatibility. Complete following width package and size.	E2620	\$759
TC4V22W32W	TRU-Comfort Back, 22" W or Greater	Center Spine or Flip Forward Back compatibility. Complete following width package and size.	E2621	\$872
TC4V22W32WBC	TRU-Comfort Back, Back Cane Option, 22" W or Greater	Back Cane compatibility. Complete following width package and size.	E2621	\$872

		Back Width											
		10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
Overall Height	14"	P	T	P	T	P	T	-	-	-	-	-	-
	15"	P	T	P	T	P	T	-	-	-	-	-	-
	16"	P	T	P	T	P	T	P	T	P	T	P	T
	17"	P	T	P	T	P	T	P	T	P	T	P	T
	18"	P	P	P	T	P	T	P	T	P	T	P	T
	19"	P	P	P	T	P	T	P	T	P	T	P	T
	20"	-	-	P	T	P	T	P	T	P	T	P	T
	21"	-	-	P	T	P	T	P	T	P	T	P	T
	22"	-	-	P	T	P	T	P	T	P	T	P	T
	23"	-	-	P	T	P	T	P	T	P	T	P	T
	24"	-	-	P	T	P	T	P	T	P	T	P	T
	25"	-	-	P	T	P	T	P	T	P	T	P	T
	26"	-	-	P	T	P	T	P	T	P	T	P	T
	27"	-	-	P	T	P	T	P	T	P	T	P	T
	28"	-	-	P	T	P	T	P	T	P	T	P	T
	29"	-	-	P	T	P	T	P	T	P	T	P	T
		P Sport Back T TRU-Comfort Back											
		Seat Back Availability Chart											

Account Number:

Date:

Quote Name/Marked for:

Q-BACK 2.0, AL SERIES, ACTIVE CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AL14AS	Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL14AT	Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL16AS	Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL16AT	Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL16AX	Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL16AXX	Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	E2613	\$703
QB2AL18AS	Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL18AT	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL18AX	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL18AXX	Active Contour, 18"Wx25-26"H Tilt/27"-28" H Recline	E2613	\$703
QB2AL20AS	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL20AT	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL20AX	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL20AXX	Active Contour, 20"Wx25-26"H Tilt/27"-28" H Recline	E2613	\$703

Q-BACK 2.0, AL SERIES, DEEP CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AL14DS	Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL14DT	Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AL16DS	Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL16DT	Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AL16DX	Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$643
QB2AL16DXX	Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703
QB2AL18DS	Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL18DT	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AL18DX	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$643
QB2AL18DXX	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703
QB2AL20DS	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL20DT	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AL20DX	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$643
QB2AL20DXX	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703

Q-BACK 2.0, AL SERIES, EXTRA DEEP CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AL14ES	Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$703
QB2AL16ES	Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$703
QB2AL16ET	Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$703
QB2AL16EX	Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$703
QB2AL18ES	Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$703
QB2AL18ET	Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$703
QB2AL18EX	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$703
QB2AL18EXX	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703
QB2AL20ES	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$703
QB2AL20ET	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$703
QB2AL20EX	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	E2620	\$703
QB2AL20EXX	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703

Account Number:

Date:

Quote Name/Marked for:

Q-BACK 2.0, AFT SERIES, ACTIVE CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AFT16AS	AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AFT18AS	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568

Q-BACK 2.0, AFT SERIES, DEEP CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AFT16DS	AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AFT16DT	AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AFT18DS	AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AFT18DT	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AFT20DT	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643

Q-BACK 2.0, AFT SERIES, EXTRA DEEP CONTOUR**Select a back type and the specific overall back height in subsequent section.**

QB2AFT16EXX	AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	E2620	\$703
-------------	--	-------	-------

BACK WIDTH PACKAGE**REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.**

PTOINDV3468	Match Seat Width Selections	No Charge
XRF000036	10"-14" Back Width Package	No Charge
XRF000037	14"-16" Back Width Package	No Charge
XRF000038	16"-20" Back Width Package	No Charge
XRF000039	18"-22" Back Width Package	No Charge
XRF000040	20"-24" Back Width Package	No Charge

BACK WIDTH**REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.**

PTOINDV3561	Set Back Width at 10"W	No Charge	PTOINDV3566	Set Back Width at 18"W	No Charge
PTOINDV3562	Set Back Width at 12"W	No Charge	PTOINDV3567	Set Back Width at 19"W	No Charge
PTOINDV3563	Set Back Width at 14"W	No Charge	PTOINDV3568	Set Back Width at 20"W	No Charge
PTOINDV3856	Set Back Width at 15"W	No Charge	PTOINDV3569	Set Back Width at 21"W	No Charge
PTOINDV3564	Set Back Width at 16"W	No Charge	PTOINDV3570	Set Back Width at 22"W	No Charge
PTOINDV3565	Set Back Width at 17"W	No Charge			

OVERALL BACK HEIGHT**REQUIRED with Sport or TRU Comfort back selection. QStatic odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. See the back availability chart for stock order back size ranges.**

SETHGHT1027	Set Back Height to 14"H	No Charge	SETHGHT1035	Set Back Height to 22"H	No Charge
SETHGHT1028	Set Back Height to 15"H	No Charge	SETHGHT1036	Set Back Height to 23"H	No Charge
SETHGHT1029	Set Back Height to 16"H	No Charge	SETHGHT1037	Set Back Height to 24"H	No Charge
SETHGHT1030	Set Back Height to 17"H	No Charge	SETHGHT1038	Set Back Height to 25"H	No Charge
SETHGHT1031	Set Back Height to 18"H	No Charge	SETHGHT1039	Set Back Height to 26"H	No Charge
SETHGHT1032	Set Back Height to 19"H	No Charge	SETHGHT1040	Set Back Height to 27"H	No Charge
SETHGHT1033	Set Back Height to 20"H	No Charge	SETHGHT1041	Set Back Height to 28"H	No Charge
SETHGHT1034	Set Back Height to 21"H	No Charge	SETHGHT1042	Set Back Height to 29"H	No Charge

PRIVACY FLAP**Optional.**

XRF000351	Privacy Flap	\$49
-----------	--------------	------

OPTIONAL BACK HEIGHT EXTENSION**Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.**

XRF000041	Back Height Extension	\$49
-----------	-----------------------	------

Account Number:

Date:

Quote Name/Marked for:

ARMRESTS**REQUIRED. Select one type of armrest.**

10153	Drop-In, H Adj, Removable, T-Arm Armrests (n)	Includes straight armpads. Select armrest height and size.	No Charge
10152	Q Captain, Flip-Up, Height Adjustable, Full Length	Includes full length, straight armpads. Select armrest height.	K0020 \$180
10157	Pediatric, Removable Armrests	Includes desk length, straight armpads. Select height adjustment from 7.75" to 8.25".	E0973 \$162
10165	Cantilever, Height Adjustable Armrests	Not available with flip forward back canes. Select armrest height (11"-17"), size, and armpad.	\$105
10156	Two-Post, H Adj, Removable Armrest (u)	Select armrest height, size, and armpad.	E0973 \$322
10158	8"-12"H, Single Post, H Adj, Remov, HD Armrests	Select armrest height, size, and armpad.	E0973 \$161
10159	10"-14"H, Single Post, H Adj, Remov, HD Armrests	Select armrest height, size, and armpad.	E0973 \$161
10150	Flip-Back, 8"-16"H Armrests (u)	Not available with Back Cane Option. Select armrest height, size, and armpad (Not available with waterfall armpads)	\$590

ARMREST SUPPORT**Optional with Flip-Back, 8"-16"H Armrests.**

10323	Flip-Back Armrest Downposts (pair)	Not available less than 16" width seat or back.	\$188
-------	------------------------------------	---	-------

ARMREST HEIGHT**REQUIRED.**

LEFT ARMREST HEIGHT			RIGHT ARMREST HEIGHT		
PTOINDV3722	7.75" Left Armrest Height	No Charge	PTOINDV3721	7.75" Right Armrest Height	No Charge
PTOINDV3535	8" Left Armrest Height	No Charge	PTOINDV3450	8" Right Armrest Height	No Charge
PTOINDV3724	8.25" Left Armrest Height	No Charge	PTOINDV3723	8.25" Right Armrest Height	No Charge
PTOINDV3536	9" Left Armrest Height	No Charge	PTOINDV3451	9" Right Armrest Height	No Charge
PTOINDV3537	10" Left Armrest Height	No Charge	PTOINDV3452	10" Right Armrest Height	No Charge
PTOINDV3538	11" Left Armrest Height	No Charge	PTOINDV3453	11" Right Armrest Height	No Charge
PTOINDV3539	12" Left Armrest Height	No Charge	PTOINDV3454	12" Right Armrest Height	No Charge
PTOINDV3540	13" Left Armrest Height	No Charge	PTOINDV3455	13" Right Armrest Height	No Charge
PTOINDV3541	14" Left Armrest Height	No Charge	PTOINDV3456	14" Right Armrest Height	No Charge
PTOINDV3542	15" Left Armrest Height	No Charge	PTOINDV3457	15" Right Armrest Height	No Charge
PTOINDV3543	16" Left Armrest Height	No Charge	PTOINDV3458	16" Right Armrest Height	No Charge

ARMREST SIZE**REQUIRED.**

XRF000057	10" Armrest Length, Left	No Charge
XRF000123	10" Armrest Length, Right	No Charge
XRF000058	14" Armrest Length, Left	No Charge
XRF000124	14" Armrest Length, Right	No Charge
XRF000059	17" Armrest Length, Left	Only available with Flip-Back, 8"-16" H Armrests and straight arm pads. No Charge
XRF000125	17" Armrest Length, Right	Only available with Flip-Back, 8"-16" H Armrests and straight arm pads. No Charge

ARMPADS**REQUIRED.**

XRF000060	No Armpad Selection	No Charge
XRF000308	Interface for Aftermarket Pad	No Charge
XRF000061	Straight, Armpad, Left	No Charge
XRF000126	Straight, Armpad, Right	No Charge
XRF000062	Waterfall, Armpad, Left	Not available with Flip-back armrests. No Charge
XRF000127	Waterfall, Armpad, Right	Not available with Flip-back armrests. No Charge

Account Number:

Date:

Quote Name/Marked for:

FALCON GEL ARMPADS FOR FLIP-BACK ARMRESTS**Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.**

ARM1806807	Falcon Gel, Full Length Left and Right	\$199
ARM1806813	Falcon Gel, Full Length Left, Desk Length Right	\$174
ARM1806810	Falcon Gel, Desk Length Left, Full Length Right	\$174
ARM1806806	Falcon Gel, Desk Length Left and Right	\$150

FALCON GEL ARMPADS FOR TUBE STYLE ARMRESTS**Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.**

ARM136438	Falcon Gel, Full Length Left and Right	\$317
ARM136442	Falcon Gel, Full Length Left, Desk Length Right	\$294
ARM136439	Falcon Gel, Desk Length Left, Full Length Right	\$294
ARM136437	Falcon Gel, Desk Length Left and Right	\$270

JOYSTICK BRACKET RECEIVER**REQUIRED.**

XRF000063	Standard Fixed Receiver	No Charge
XRF000064	Stealth Height Adjustable Receiver	K0108 \$236

HEADREST**Optional.****HEADREST PADS**

CP08CCSPM	Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)	E0955	\$179
CP10CCSPM	Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)	E0955	\$179
CP14CCSPM	Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)	E0955	\$179
CP06CCSPM	Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)	E0955	\$179

Standard	\$50
TRU-Comfort/Sport	
14"-21"H	22"-29"H
Q-Back 2	
19"-24"H	25"-26"H
Q Cam Lock Avail.	
(overall height)	

HEADREST BRACKETS

ST-UHB2D33200-Q	Comfort Plus Headrest Mounting Bracket, UniLink	E1033	\$311
ST-HM0475-17M-Q	Multi-Axis Removable Headrest Bracket	E1033	\$174

Q CAM LOCK

XRF000376	Q Cam Lock	\$51
-----------	------------	------

LAP/ PELVIC BELTS**REQUIRED. Select one lap/pelvic belt below.**

10167	38" Lap Belt		No Charge
10168	50" Lap Belt		No Charge
10169	60" Lap Belt		No Charge
10170	70" Lap Belt		No Charge
10171	80" Lap Belt		No Charge
10172	90" Lap Belt		No Charge
ST-IFRB5350	Retractable Lap Belt	Stealth Positioning IFRB5350.	\$181
ST-IFPB5213-Q	60" Single Pull, 1.5", Push Button, M Pads (7.5")	Stealth Positioning IFPB5213.	\$50
ST-IFPB5314-Q	60" Single Pull, 2", Push Button, L Pads (9.5")	Stealth Positioning IFPB5214.	\$50
ST-IFPB5315-Q	60" Single Pull, 2", Push Button, XL Pads (11.5")	Stealth Positioning IFPB5215.	\$50
ACCBELT1016	60" Sure-Lok Lap Belt (Non-retractable)	Required with Compliant Transit Option.	\$120

LEGREST CONTROLS**Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.**

XRF000108	Independent Legrests	Adds 2 actuator control requirement.	No Charge
XRF000109	Combined Legrests	Adds 2 actuator control requirement.	No Charge
XRF000110	Articulating Foot Platform	Adds 1 actuator control requirement.	No Charge

Account Number:

Date:

Quote Name/Marked for:

FRONT RIGGINGS

Optional. Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

BASIC CENTER MOUNT**Basic Center Mount Foot Platform Assembly**

XRF000323	Basic Center Mount	Select extension and footplate. Calf pads available depending on configuration.	No Charge
-----------	--------------------	---	-----------

Basic Center Mount Extension and Angle Adjustable Footplates

10189	Short Extension, 5"x9" Footplates	Not available with calf pads.	No Charge
10190	Medium Extension, 5"x9" Footplates	Requires small calf pads with power tilt configurations.	No Charge
10191	Long Extension, 5"x9" Footplates	Requires small or large calf pads with power tilt configurations.	No Charge

Basic Center Mount Calf Pads

PTOINDV3768	No Calf Pads	No Charge	Short	Med.	Long
10194	Small Calf Pads, Basic Center Mount	\$161	8-10	10-14	14-19
10195	Large Calf Pads, Basic Center Mount	\$161	Basic Center Mount Extension Range (in.)		

Basic Center Mount Additional Accessory

10196	Basic Center Mount Footplate Connector	Optional.	\$16
-------	--	-----------	------

CENTER MOUNT FOOT PLATFORM**Center Mount Foot Platform Type**

PTOINDV3754	Power Articulating Foot Platform	Select extension with footplates and calf support. Other items are optional.	E1012+ K0040	\$3,145	Medium	Long
					11"-17"	14"-20"
					Power AFP Extension Range	

Extension and Angle Adjustable Footplates

10370	Medium Extension, 5"x8" Footplates	No Charge
10371	Medium Extension, 5"x11" Footplates	No Charge
10372	Long Extension, 5"x8" Footplates	No Charge
10373	Long Extension, 5"x11" Footplates	No Charge

Center Mount Calf Pads

10206	Small Calf Pads, AFP Only, 6"Wx6"H (pair)	Small or large calf pad standard with AFP.	No Charge
10357	Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair)	Small or large calf pad standard with AFP.	No Charge

Center Mount Calf Panel

10207	Small Calf Panel, 18"Wx8"H		\$188
10208	Large Calf Panel, 24"Wx8"H	Only available on widths of 20" or greater.	\$242

Optional Center Mount Modular Footplate Extension Components

Assembly required. May require seating configuration adjustments to avoid front caster interference.

10209	Width Extension	\$70
10210	Outside Edge	\$70
10211	Front and Rear Cap	\$16

Optional Center Mount Additional Accessories

ST-CSCS0775-Q	Small Foam Calf Pad Cover (pair)		\$67
ST-CSCS0985-Q	Large Foam Calf Pad Cover (pair)		\$83
ST-CSCP0920-Q	Small Foam Panel Cover (pair)		\$176
ST-CSCP0927-Q	Large Foam Panel Cover (pair)		\$282
10212	Footplate Connector		\$16
XRF000348	Footplate Riser	Not available with Long Extension. Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments.	\$77
10213	Footplate Riser Calf Panel	Only available and required with Footplate Riser selection and extension lengths of 6" to 9".	\$51

Account Number:

Date:

Quote Name/Marked for:

Optional Stealth Gel Footplate Covers

ST-FPG5008T-Q	5"x8" Tapered Gel Footplate Covers	\$222
ST-AFPG5011T-Q	5"x11" Tapered Gel Footplate Covers	\$270
ST-AFPG5008E-Q	5"x8" Tapered Gel Footplate Covers for Width Extension	\$247
ST-AFPG5011E-Q	5"x11" Tapered Gel Footplate Cover for Width Extension	\$270

POWER ELEVATING LEGRESTS**Power Elevating Legrest Upper Assembly****Requires power legrest control selection found in earlier section.**

XRF000341	Power Elevating Legrests, Pair	Choose extension and pair footplates.	E1010+ K0040+ E0951	\$2,516
XRF000342	Power Elevating Legrests, Left	Choose extension and left footplate.	K0108	\$1,258
XRF000343	Power Elevating Legrests, Right	Choose extension and right footplate.	K0108	\$1,258

Legrest Length Matrix	Seat to Floor Option			Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-12.5"	8.5"-13.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
12.5"-21" Lower Length	N/A	12.5"-13.5"	12.5"-14.5"	12.5"-21" Lower Length	12.5"-15.5"	12.5"-15.5"	12.5"-17"

Lower Extension

XRF000344	8.5" to 14.5" Lower Length	No Charge
XRF000347	12.5" to 21" Lower Length	No Charge

Angle Adjustable Footplate(s)**The 6"W footplates may cause interference with seat widths less than 16". The 8"W footplates may cause interference with seat widths less than 19".**

10216	4"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10358	4"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10364	4"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge
10217	5"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10359	5"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10365	5"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge
10218	5"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10360	5"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10366	5"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge
10219	6"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10361	6"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10367	6"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge
10220	6"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10362	6"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10368	6"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge
10221	8"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10363	8"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Lt	No Charge
10369	8"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Rt	No Charge

Account Number:

Date:

Quote Name/Marked for:

HEAVY DUTY, DROP-IN, SWING-AWAY LEGRESTS**70° Drop-In Legrest**

10269	70° Heavy Duty, Drop-In, Swing-Away Legrests, Pair	No Charge
10270	70° Heavy Duty, Drop-In, Swing-Away Legrest, Left	No Charge
10271	70° Heavy Duty, Drop-In, Swing-Away Legrest, Right	No Charge

Lower Extension**Rigging configurations may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability when noted below.**

10272	8.5"-12" Ext Length, HD Drop-Ins, Pair	No Charge
10273	8.5"-12" Ext Length, HD Drop-Ins, Left	No Charge
10274	8.5"-12" Ext Length, HD Drop-Ins, Right	No Charge
10275	13"-18" Ext Length, HD Drop-Ins, Pair	No Charge
10276	13"-18" Ext Length, HD Drop-Ins, Left	No Charge
10277	13"-18" Ext Length, HD Drop-Ins, Right	No Charge
10382	16"-21" Ext Length, HD Drop-Ins, Pair	No Charge
10383	16"-21" Ext Length, HD Drop-Ins, Left	No Charge
10384	16"-21" Ext Length, HD Drop-Ins, Right	No Charge

Angle Adjustable Footplate(s)

10278	6"Wx8"D AAdj Ftpts, Heel Loops, HD Drop-In, Pair	K0040+ E0951	\$184
10279	6"Wx8"D AAdj Ftpts, Heel Loops, HD Drop-In, Left	K0040+ E0951	\$92
10280	6"Wx8"D AAdj Ftpts, Heel Loops, HD Drop-In, Right	K0040+ E0951	\$92
10281	6"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Pair	K0040+ E0951	\$184
10282	6"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Left	K0040+ E0951	\$92
10283	6"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Right	K0040+ E0951	\$92
10284	8"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Pair	K0040+ E0951	\$184
10285	8"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Left	K0040+ E0951	\$92
10286	8"Wx11"D AAdj Ftpts, Heel Loops, HD Drop-In, Right	K0040+ E0951	\$92

SET LEGREST EXTENSION LENGTH**REQUIRED with center mounts and matching left and right rigging configurations. Select one length.****CENTERMOUNT EXTENSION LENGTH**

PTOINDV3989	2" L Legrest Length	No Charge	PTOINDV3689	12.5" L Legrest Length	No Charge
PTOINDV3990	3" L Legrest Length	No Charge	PTOINDV3690	13" L Legrest Length	No Charge
PTOINDV3991	4" L Legrest Length	No Charge	PTOINDV3691	13.5" L Legrest Length	No Charge
PTOINDV3992	5" L Legrest Length	No Charge	PTOINDV3692	14" L Legrest Length	No Charge
PTOINDV3993	6" L Legrest Length	No Charge	PTOINDV3693	14.5" L Legrest Length	No Charge
PTOINDV3994	7" L Legrest Length	No Charge	PTOINDV3694	15" L Legrest Length	No Charge
PTOINDV3680	8" L Legrest Length	No Charge	PTOINDV3695	15.5" L Legrest Length	No Charge
PTOINDV3681	8.5" L Legrest Length	No Charge	PTOINDV3696	16" L Legrest Length	No Charge
PTOINDV3682	9" L Legrest Length	No Charge	PTOINDV3697	16.5" L Legrest Length	No Charge
PTOINDV3683	9.5" L Legrest Length	No Charge	PTOINDV3698	17" L Legrest Length	No Charge
PTOINDV3684	10" L Legrest Length	No Charge	PTOINDV3699	17.5" L Legrest Length	No Charge
PTOINDV3685	10.5" L Legrest Length	No Charge	PTOINDV3700	18" L Legrest Length	No Charge
PTOINDV3686	11" L Legrest Length	No Charge	PTOINDV3701	18.5" L Legrest Length	No Charge
PTOINDV3687	11.5" L Legrest Length	No Charge	PTOINDV3702	19" L Legrest Length	No Charge
PTOINDV3688	12" L Legrest Length	No Charge			

Account Number:

Date:

Quote Name/Marked for:

FRONT RIGGING EXTENSION LENGTH

PTOINDV3470	7" L Legrest Length	No Charge	PTOINDV3486	15" L Legrest Length	No Charge
PTOINDV3471	7.5" L Legrest Length	No Charge	PTOINDV3487	15.5" L Legrest Length	No Charge
PTOINDV3472	8" L Legrest Length	No Charge	PTOINDV3488	16" L Legrest Length	No Charge
PTOINDV3473	8.5" L Legrest Length	No Charge	PTOINDV3489	16.5" L Legrest Length	No Charge
PTOINDV3474	9" L Legrest Length	No Charge	PTOINDV3490	17" L Legrest Length	No Charge
PTOINDV3475	9.5" L Legrest Length	No Charge	PTOINDV3491	17.5" L Legrest Length	No Charge
PTOINDV3476	10" L Legrest Length	No Charge	PTOINDV3492	18" L Legrest Length	No Charge
PTOINDV3477	10.5" L Legrest Length	No Charge	PTOINDV3493	18.5" L Legrest Length	No Charge
PTOINDV3478	11" L Legrest Length	No Charge	PTOINDV3494	19" L Legrest Length	No Charge
PTOINDV3479	11.5" L Legrest Length	No Charge	PTOINDV3495	19.5" L Legrest Length	No Charge
PTOINDV3480	12" L Legrest Length	No Charge	PTOINDV3496	20" L Legrest Length	No Charge
PTOINDV3481	12.5" L Legrest Length	No Charge	PTOINDV3497	20.5" L Legrest Length	No Charge
PTOINDV3482	13" L Legrest Length	No Charge	PTOINDV3498	21" L Legrest Length	No Charge
PTOINDV3483	13.5" L Legrest Length	No Charge	PTOINDV3600	21.5" L Legrest Length	No Charge
PTOINDV3484	14" L Legrest Length	No Charge	PTOINDV3669	22" L Legrest Length	No Charge
PTOINDV3485	14.5" L Legrest Length	No Charge	PTOINDV3670	22.5" L Legrest Length	No Charge

SEATING ELECTRONICS**REQUIRED with power seating or legrest options. Select one type of power actuator control.****ACTUATOR CONTROL THROUGH JOYSTICK**

XRF000335	Single Actuator Through Controller and Interface	E2310	\$2,206
XRF000336	Multiple Actuator Through Controller and Interface	E2311	\$3,101

ATTENDANT CONTROL**Optional. Mounting bracket is included.**


XRF000078	Attendant Control	E2331	\$536
-----------	-------------------	-------	-------

ACCESSORIES**Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.**

XRF000128	Accessory Bar	Horizontal rail attached to rear of seat. Included with certain accessories when noted.	\$135.00	Backpack Holder	
10375	Accessory Hook	Coming soon.	\$29.00	Clothing Guard	
10326	Backpack Holder		\$70.00	Glove Box	
10327	Cane/Crutch Holder		E2207 \$107.00	Phone Holder	
10331	Cell Phone/ Media Holder, X-Grip	Not available with Heavy Duty armrests.	\$70.00	Hydration System	
10332	Clothing Guard		\$242.00	Personal Item Hook	
10334	Cup Holder		No Charge		
10337	Flag Kit		\$70.00		
10338	Flip-Back Armrest, Side Track Accessory Bracket		\$11.00		
10340	Glove Box		\$161.00		
10342	Hydration System		\$242.00		
10343	Large, Black, Weather Cover		\$202.00		
10344	Oxygen Holder	Includes accessory bar, a credit will be applied PTOINDV3769.	E2208 \$214.00		
10348	Personal Item Hook		\$70.00		

Account Number: Date: Quote Name/Marked for:

10349	Personal Item Pouch		\$86.00
10350	Privacy Bag		\$43.00
10351	Push Handles	Not available with oxygen holder or vent tray.	\$199.00
10352	Quantum Backpack		\$171.00
10353	Rear View Mirror		\$70.00
10354	Reflector Kit		\$54.00
XRF000379	Remote Stop Switch		\$1,321.00
ACC7120006	Stealth Cell Phone/Media Device Holder		\$255.00
ACC141818	Stealth Large Tablet Holder		\$625.28
XRF000396	Transfer Handles (pair)	E0705	\$214.00
10324	Aftermarket Accessory Mounting Bracket, Non-slotted	Quantity: __ Select a quantity up to 5.	\$54.00
10325	Aftermarket Accessory Mounting Bracket, Slotted	Quantity: __ Select a quantity up to 5.	\$54.00



BRANCH PATCHES

Optional.

ACC140240	Air Force Patch	\$21
ACC140241	Army Patch	\$21
ACC140239	Coast Guard Patch	\$21
ACC140243	Marines Patch	\$21
ACC140244	Navy Patch	\$21
ACC140238	Purple Heart Patch	\$21



POSITIONING COMPONENTS

Optional. More options are available on the Positioning Components Order Form. Please provide quantity where applicable.

STEALTH THORACIC LATERALS

Swing-away laterals may interfere with flip-back and 2-post, flip up armrests.

Stealth Thoracic Lateral Bracket

ST-ULS4042000-Q	Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)	Only available with Redesigned back assemblies.	E1034	\$664.00
ST-ULPK404108B00-Q	Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00)	Only available with Sport or TRU Comfort backs.	E1034	\$664.00
ST-ULPK404104B00-Q	Swing-away, UniLink, Back Cane, pair (ULPK404104B00)	Only available with back cane and omit back configurations.	E1034	\$664.00
ST-ULPK404103L00-Q	Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00)	Only available with Active Contour Q-Backs.	E1034	\$332.00
ST-ULPK404103R00-Q	Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00)	Only available with Active Contour Q-Backs.	E1034	\$332.00

Stealth Thoracic Lateral Pads

INDPART2663	Curved Pads, 3"x5", pair (TWBL-3X5)	E0956	\$250.00
INDPART2664	Curved Pads, 4"x6", pair (TWBL-4X6)	E0956	\$250.00
ST-ULP846050320-Q	Flat Pads, 3"x5", pair (ULP846050320)	E0956	\$250.00
ST-ULP835050320-Q	Flat Pads, 4"x6", pair (ULP835050320)	E0956	\$250.00

Stealth Thoracic Lateral Installation

PTOINDV3499	Install Laterals	No Charge
PTOINDV3500	Do Not Install	No Charge

Account Number:

Date:

Quote Name/Marked for:

STEALTH UNILINK PELVIC/THIGH GUIDE**UniLink Pelvic/Thigh Hardware****^ Please verify the appropriate billing code depending upon the placement of the support pad.**

ST-UTTP5-Q	UniLink, 5" Height Adjustable & Removable Hardware	Quantity: ____	E1028/ E1034^	\$265.00
ST-UTTP7-Q	UniLink, 7" Height Adjustable & Removable Hardware	Quantity: ____	E1028/ E1034^	\$276.00
ST-UTT8A11602-L-Q	UniLink Multi-Angle Thigh Guide, Left		E1028/ E1034^	\$358.87
ST-UTT8A11602-R-Q	UniLink Multi-Angle Thigh Guide, Right		E1028/ E1034^	\$358.87
ST-UTS25100-L-Q	UniLink Flip Down Hip/Thigh Support, Left		E1028/ E1034^	\$332.00
ST-UTS25100-R-Q	UniLink Flip Down Hip/Thigh Support, Right		E1028/ E1034^	\$332.00
POS140428	Flip Down with LEFT Full Surface Contact Feature (TWBADD-L & LPHW407)		E1028/ E1034^	\$281.00
POS140429	Flip Down with RIGHT Full Surface Contact Feature (TWBADD-R & LPHW407)		E1028/ E1034^	\$281.00

UniLink Pelvic/Thigh Pads**^ Please verify the appropriate billing code to be used with all third-party payors.**

ST-UTP544056100-Q	UniLink, 4" x 4" Pad	Quantity: ____	E0953/ E0956^	\$106.00
ST-UTP548056100-Q	UniLink, 4" x 8" Pad	Quantity: ____	E0953/ E0956^	\$106.00
ST-UTP54D056100-Q	UniLink, 4" x 12" Pad	Quantity: ____	E0953/ E0956^	\$106.00

UniLink Pelvic/Thigh Additional Links

ST-UIL001150A000-Q	UniLink, 1.5" Link	Quantity: ____		\$81.00
ST-UIL001200A000-Q	UniLink, 2" Link	Quantity: ____		\$81.00

UniLink Pelvic/Thigh Guide Installation

PTOINDV3585	Install Left Lateral Pelvic/ Thigh Bracket		No Charge
PTOINDV3586	Install Right Lateral Pelvic/ Thigh Bracket		No Charge
PTOINDV3719	Do Not Install Lateral Pelvic/ Thigh Bracket		No Charge

STEALTH POSITIONING CHEST HARNESS

ST-IFPS1103-Q	Upper Classic, Medium (IFPS1103)	E0960	\$154
ST-IFPS1104-Q	Upper Classic, Large (IFPS1104)	E0960	\$154
ST-IFPS2103-Q	Upper Contour, Medium (IFPS2103)	E0960	\$154
ST-IFPS2104-Q	Upper Contour, Large (IFPS2104)	E0960	\$154

**STEALTH POSITIONING CHEST STRAPS**

ST-IFCS4113-Q	Chest Strap, Medium w/ Belt Anchor (IFCS4113)	E0960	\$97
ST-IFCS4114-Q	Chest Strap, Large w/ Belt Anchor (IFCS4114)	E0960	\$97

**STEALTH SWING-AWAY HEADREST MOUNTS****These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.****Facial Lateral Hardware**

ST-SAR2A-L-Q	UniLink Facial Lateral, Swing Away, Left	E1028	\$318
ST-SAR2A-R-Q	UniLink Facial Lateral, Swing Away, Right	E1028	\$318
ST-SAR4A-DBU-LL-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper & 1 Lower Left	E1028	\$492
ST-SAR4A-DBU-RL-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper & 1 Lower Right	E1028	\$492
ST-SAR4A-DBU-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper	E1028	\$434



Account Number:

Date:

Quote Name/Marked for:

Facial Lateral Pads and Accessories

INDPART2815	Egg Switch Mount (ESM)	Quantity: ____		\$57
ST-900M-Q	Facial Pad, 4.5"x2", Cool Core (900M-Q)	Quantity: ____	K0108	\$61
ST-920M-Q	Facial Pad, 6"x2", Cool Core (920M-Q)	Quantity: ____	K0108	\$61
ST-SUS18-Q	Swing-Away Joystick Mount, 18" Rod	Quantity: ____		\$276
Allows mounting of facial pads and attaches to the headrest.				

STEALTH ELBOW SUPPORT**Elbow Support Bracket**

POS124064	Left, Swing-Away, Elbow Support Bracket Tube Style Arms (ARMS-ES-L)	E1028	\$199
POS124066	Right, Swing-Away, Elbow Support Bracket Tube Style Arms (ARMS-ES-R)	E1028	\$199
ST-UETB4021FL-Q	Left, UniLink, Fixed, Elbow Support Bracket for Flip-Back Armrests		\$211
ST-UETB4021FR-Q	Right, UniLink, Fixed, Elbow Support Bracket for Flip-Back Armrests		\$211

Elbow Support Pad

POS124065	3"Wx3"H Elbow Pad (ARMS-ES-3X3)	Quantity: ____	K0108	\$57
POS124067	4"Wx4"H Elbow Pad (ARMS-ES-4X4)	Quantity: ____	K0108	\$57

STEALTH ARM STOP**Arm Stop Hardware and Pad**

ST-ASF6X2C-S-Q	Stealth Arm Stop	Quantity: ____		\$351
----------------	------------------	----------------	--	-------

Arm Stop Location

PTOINDV3936	Mount Left Arm Stop On Inside	No Charge
PTOINDV3937	Mount Right Arm Stop On Inside	No Charge
PTOINDV3938	Mount Left Arm Stop On Outside	No Charge
PTOINDV3939	Mount Right Arm Stop On Outside	No Charge

STEALTH FOOT BOXES

INDPART2786	Stealth Small Foot Boxes, 5"Wx8"Dx7"H	E0954	\$539
INDPART2787	Stealth Medium Foot Boxes, 6"Wx10"Dx9"H	E0954	\$539
INDPART2788	Stealth Large Foot Boxes, 7"Wx12"Dx11"H	E0954	\$539

STEALTH POSITIONING SHOE HOLDER

ST-IFSH4211-Q	Small Shoe Holder, 3/4" Webbing (IFSH4211)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87
ST-IFSH4212-Q	Medium Shoe Holder, 3/4" Webbing (IFSH4212)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87
ST-IFSH4213-Q	Large Shoe Holder, 3/4" Webbing (IFSH4213)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87



Account Number:

Date:

Quote Name/Marked for:

QUANTUM LAP TRAY**Optional. Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.****BLACK QUANTUM LAP TRAY**

MEC1905328-11B	Black Lap Tray, 11" W	For use with 12"-13" back width.	E0950	\$204
MEC1905328-13B	Black Lap Tray, 13" W	For use with 14"-15" back width.	E0950	\$204
MEC1905328-15B	Black Lap Tray, 15" W	For use with 16"-17" back width.	E0950	\$204
MEC1905328-17B	Black Lap Tray, 17" W	For use with 18"-19" back width.	E0950	\$204
MEC1905328-19B	Black Lap Tray, 19" W	For use with 20"-22" back width.	E0950	\$204

CLEAR QUANTUM LAP TRAY

MEC1905328-11C	Clear Lap Tray, 11" W	For use with 12"-13" back width.	E0950	\$204
MEC1905328-13C	Clear Lap Tray, 13" W	For use with 14"-15" back width.	E0950	\$204
MEC1905328-15C	Clear Lap Tray, 15" W	For use with 16"-17" back width.	E0950	\$204
MEC1905328-17C	Clear Lap Tray, 17" W	For use with 18"-19" back width.	E0950	\$204
MEC1905328-19C	Clear Lap Tray, 19" W	For use with 20"-22" back width.	E0950	\$204

QUANTUM LAP TRAY MOUNTING HARDWARE

MEC2104886	Mounting Bracket, Chest Bar/Lap Tray			\$77
------------	--------------------------------------	--	--	------

SIMPLE SOLUTIONS LAP TRAY**Optional. Mounting bracket is required with lap tray selection. See chart below for compatibility.****SIMPLE SOLUTIONS LAP TRAY**

ST-SSLXS010001171-Q	Simple Solutions Small Tray with Cut Out	E0950	\$309
ST-SSLXM010001171-Q	Simple Solutions Medium Tray with Cut Out	E0950	\$309
ST-SSLXL010001171-Q	Simple Solutions Large Tray with Cut Out	E0950	\$309

SIMPLE SOLUTIONS LAP TRAY MOUNTING HARDWARE

ST-SSLXHD10000060-Q	Swing Away Mounting Hardware		\$302
---------------------	------------------------------	--	-------

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■ :Compatible ■ :Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

