



Quantum Rehab®, 401 York Ave., Duryea, PA 18642
Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2 Order Form with Q Static Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____


Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, option numbers, and availability are subject to change without notice. Prices and option numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT DETAILS

Weight: _____ Height: _____



A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L.

OPTION NUMBER	DESCRIPTION	NOTES	HCPCS	MSRP
Q6 EDGE 2.0 3S-SS	Quantum Q6 Edge 2.0 3S-SS	Solid seat pan standard with no power option.	K0848	\$8,704
Q6 EDGE 2.0 3SP-SS	Quantum Q6 Edge 2.0 3SP-SS	Single power option.	K0856	\$8,399

Account Number:

Date:

Quote Name/Marked for:

DRIVE MOTORS**REQUIRED. Select one.**

10034	6 mph Drive Motors		No Charge
10035	6 mph Accu-Trac Motor Technology	Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.	K0108 \$1,387.00

BASE COLOR**REQUIRED. Select one.**

10001	 Red Corvette Base Color	No Charge	10006	 Evolution Orange Base Color	No Charge
10003	 Back in Black Base Color	No Charge	10009	 Midnight Blue Base Color	No Charge
10004	 Night Train Base Color	No Charge	10002	 Electric Blue Base Color	No Charge
10008	 Grey Street Base Color	No Charge	10023	 Lemon Crush Base Color	No Charge
10007	 Deep Purple Base Color	No Charge	10010	 Green Machine Base Color	No Charge
10005	 Pretty in Pink Base Color	No Charge	10011	 Military Green Base Color	No Charge

BASE MOUNTED LIGHTS**Optional. Select one. Seat Mounted, LED, Full Light Package option can be found later on this form.**

10038	Accent Lights	Accent lights operate by pressing included single switch button.	No Charge
XRF000177	Switch Operated Accent Lights without Lift System	Accent lights operate by pressing included single switch button.	\$268.00
XRF000129	Base Mounted, LED, Full Light Package	Requires QLNES 10 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.	\$1,066.00

**QUANTUM BACKUP CAMERA****Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.**

10041	Quantum Backup Camera		\$375
OPTIONAL CAMERA MOUNTING BRACKET			
10042	Backup Camera Headrest Bracket	Works with .5" x .5" vertical supports only.	\$21
10043	Backup Camera Track Mount Bracket		\$28
OPTIONAL DISPLAY MOUNTING BRACKET			
10044	Backup Camera Display, Joystick Mount Bracket	Only available with Q Logic EX Joystick.	\$107

**DRIVE WHEELS****REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**

XRF000356	3" x 14" Pneumatic Wheel	Black tires.	No Charge
XRF000355	3" x 14" Flat Free Wheel	Black tires.	No Charge
XRF000367	Knobby, Pneumatic Wheel	4"x14" Black Tires.	\$592
XRF000368	Knobby, Flat Free Wheel	4"x14" Black Tires.	\$592

Account Number:

Date:

Quote Name/Marked for:

DRIVE WHEEL ACCENT COLOR**Optional. Not available with knobby tires.**

XRF000179	Back in Black Wheel Accent Color	\$64.00	XRF000186	Evolution Orange Wheel Accent Color	\$64.00
XRF000180	Red Corvette Wheel Accent Color	\$64.00	XRF000187	Midnight Blue Wheel Accent Color	\$64.00
XRF000181	Night Train Wheel Accent Color	\$64.00	XRF000188	Electric Blue Wheel Accent Color	\$64.00
XRF000182	Grey Street Wheel Accent Color	\$64.00	XRF000189	Lemon Crush Wheel Accent Color	\$64.00
XRF000184	Deep Purple Wheel Accent Color	\$64.00	XRF000190	Green Machine Wheel Accent Color	\$64.00
XRF000185	Pretty in Pink Wheel Accent Color	\$64.00	XRF000191	Military Green Wheel Accent Color	\$64.00

BATTERIES**Optional. The powerchair requires 2 batteries and will be installed if ordered.**

10067	NF-22, Gel Cell, Interceptor 2 Battery	MSRP is each battery.	E2361	\$362
10068	NF-22, AGM, MK Battery	MSRP is each battery.	E2361	\$441
10069	NF-22, Gel Cell, MK Battery	MSRP is each battery.	E2361	\$441

ELECTRONICS**REQUIRED. Select one drive control. The No Joystick option is used to omit the joystick and utilize a specialty control, like those on the Alternative Controls and Electronics Order Form.****NON-EXPANDABLE CONTROLS****Non-expandable do not allow supplementary drive controllers.**

10077	QLNES 6 Key Joystick	Controls up to 2 actuators.	No Charge
10078	QLNES 10 Key Joystick	Controls up to 2 actuators. Operates Full Lights.	\$429

EXPANDABLE CONTROLS**Q-Logic 3 controls require the controller module and expandable harness. If a Q-Logic joystick or specialty control is selected elsewhere, Customer Service will quote the order with these items.**

XRF000377	Q-Logic 3e Joystick	Includes 10082 Q-Logic Magnetic Charger below. Operates Full Lights.	\$735
10079	Q-Logic 3 Expandable Capacity Joystick	(EX) Includes XLR USB charger.	\$857
10080	Q-Logic 3 Expandable Capacity Joystick for Full Light Package	(EX) Includes XLR USB charger. Operates Full Lights.	\$857
XRF000010	No Joystick, Setup for Q-Logic	Requires specialty controls selection.	No Charge

The Q-Logic 3 Controller Module and Controller Harness options below are required with Expandable electronics selections above.

XRF000002	Q-Logic 3 Controller Module	E2377	\$745
XRF000264	Q-Logic 3 Expandable Controller Harness	E2313	\$467

Q-LOGIC MAGNETIC CHARGER FOR QLNES

10082	Magnetic to XLR Charger Adapter with Wall Holder	Optional with QLNES.	\$61
-------	--	----------------------	------



QLNES Joystick



Q-Logic Joystick

Account Number:

Date:

Quote Name/Marked for:

Q-LOGIC 3 EX JOYSTICK ARMOR**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Black Soft Touch matches Night Train color.**



10084	Black, Soft Touch, Joystick Armor	No Charge	10091	Evolution Orange Joystick Armor	\$64.00
10085	Red Corvette Joystick Armor	\$64.00	10092	Midnight Blue Joystick Armor	\$64.00
10086	Back in Black Joystick Armor	\$64.00	10093	Electric Blue Joystick Armor	\$64.00
10087	Grey Street Joystick Armor	\$64.00	10095	Lemon Crush Joystick Armor	\$64.00
10088	Deep Purple Joystick Armor	\$64.00	10096	Green Machine Joystick Armor	\$64.00
10090	Pretty in Pink Joystick Armor	\$64.00	10097	Military Green Joystick Armor	\$64.00

JOYSTICK HANDLES**Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.**

ACCASMB1860	Mushroom Handle	E2323	\$155
ACCASMB1861	Chin Cup	E2323/ E2324	\$155
ACCASMB1858	T-Bar	E2323	\$195
ACCASMB1859	Medium Ball Handle	E2323	\$155
ACCASMB1856	Soft Ball Handle	E2323	\$174
ACCASMB1857	Extended I-Handle	E2323	\$174
PLSKNCP1057	Bodypoint 3" Goal Post Handle	PC101B	E2323 \$179
PLSKNCP1058	Bodypoint 4" Goal Post Handle	PC102B	E2323 \$179

**JOYSTICK MOUNTING BRACKETS****Optional. Select a mount to be used with a joystick selected in previous section.**

10109	Fixed, In-line, Bracket	No Charge
10110	Swing-away, Adjustable, Joystick Bracket	Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. E1032 \$325
10111	Basic Swing-away Joystick Bracket	E1032 \$245
XRF000003	Retract4 Mount	Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. E1032 \$375
XRF000004	Multi-Axis Mount	Allows height adjustment of joystick. \$129

**JOYSTICK MOUNT LOCATION**

PTOINDV1225	Mount Controller on Left	No Charge
PTOINDV1226	Mount Controller on Right	No Charge

Q STATIC SEAT**Q STATIC SEAT CONFIGURATION****REQUIRED. Maximum seat size is 22"W x 22"D. 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.**

10112	Q Static Basic	Excludes Power Lift, Center Spine, Flip Forward Back, Aftermarket Back, Flip Back Armrests, AFP, Power ELR, and Custom options.	No Charge
10114	Q Static Positioning		\$2,458

12" POWER LIFT**Optional. ^Please verify the appropriate billing code to be used for this item with all third-party payors.**

10115	Q Static 12" Lift	Available with Q Static Positioning. 12" of scissor lift with iLevel System. Requires seat function selection above.	E2298/ E2300^	\$3,725
-------	-------------------	--	------------------	---------

Account Number:

Date:

Quote Name/Marked for:

OCCUPIED TRANSIT

Optional. The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Omitting a backrest selection requires a provider to facilitate an aftermarket backrest that best fits the consumer's needs. The third-party backrest If applicable must be occupied transit approved per ISO/RESNA standards. When modifying a power chair, including the addition of third-party componentry, to accommodate the medical needs of the power chair occupant, a risk assessment in conformance with ISO 14971, as outlined in ISO 7176, should be performed.

10120	Occupied Transit	WC-19 occupied transit. To determine if your chair is WC-19 compliant (occupied), please refer to this website: https://wc-transportation-safety.umtri.umich.edu/	\$295
-------	------------------	---	-------

SEAT WIDTH**Required. Maximum seat width is 22".**

SETWDTH1000	10" Seat Width	No Charge	SETWDTH1007	17" Seat Width	No Charge
SETWDTH1001	11" Seat Width	No Charge	SETWDTH1008	18" Seat Width	No Charge
SETWDTH1002	12" Seat Width	No Charge	SETWDTH1009	19" Seat Width	No Charge
SETWDTH1003	13" Seat Width	No Charge	SETWDTH1010	20" Seat Width	No Charge
SETWDTH1004	14" Seat Width	No Charge	SETWDTH1011	21" Seat Width	No Charge
SETWDTH1005	15" Seat Width	No Charge	SETWDTH1012	22" Seat Width	No Charge
SETWDTH1006	16" Seat Width	No Charge			

SEAT WIDTH PACKAGE**REQUIRED with seat selection. Maximum width is 22".**

10376	10"-12" Width Package	Only available with back canes and 10" to 14" seat depth.	No Charge
10377	12"-14" Width Package	Only available with back canes and 12" to 20" seat depth.	No Charge
XRF000006	12"-16" Width Package		No Charge
XRF000007	14"-16" Width Package		No Charge
XRF000008	16"-20" Width Package		No Charge
10385	20"-22" Width Package		No Charge

SEAT DEPTH**REQUIRED with seat selection. Maximum seat depth is 22".**

SETDPTH1000	10" Seat Depth	No Charge	SETDPTH1007	17" Seat Depth	No Charge
SETDPTH1001	11" Seat Depth	No Charge	SETDPTH1008	18" Seat Depth	No Charge
SETDPTH1002	12" Seat Depth	No Charge	SETDPTH1009	19" Seat Depth	No Charge
SETDPTH1003	13" Seat Depth	No Charge	SETDPTH1010	20" Seat Depth	No Charge
SETDPTH1004	14" Seat Depth	No Charge	SETDPTH1011	21" Seat Depth	No Charge
SETDPTH1005	15" Seat Depth	No Charge	SETDPTH1012	22" Seat Depth	No Charge
SETDPTH1006	16" Seat Depth	No Charge			

SEAT TO FLOOR HEIGHT**REQUIRED.**

PTOINDV3558	Lowest Seat to Floor	No Charge
PTOINDV3559	Middle Seat to Floor	No Charge
PTOINDV3560	Highest Seat to Floor	No Charge

Low	Med	High
16.5"	17.5"	18.5"
Seat to Floor		

Account Number:

Date:

Quote Name/Marked for:

STEALTH SEAT CUSHION

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete Stealth Cushion Order Forms. Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

Cushion Model	####-###XXYY	Width	Depth
Cushion Top Foam V=Visco Top Foam: 1" thick, XXSoft viscoelastic foam topper. F=FRG Top Foam: Under 22"W, blue XXSoft FRG foam topper. 22"W and over, green MedSoft FRG foam topper. N=No Top Foam	Cushion Liner U=Urethane Liner: 2Way Stretch. P=Polyurethane Liner: 4Way Stretch. N=No Liner	Cushion Cover C=CoolCore: Chemical free, 4Way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. M=Mesh: Superior air circulation, non-slip. T=Stealth Tek: Waterproof with stretch capability.	

CUSHION TYPE

Ref: XRF000116	TRU-Comfort 2 SPP Seat Cushion	Part series: TC2-###XXYY. HCPCS: E2607/E2608. MSRP: \$441-\$651.
Ref: XRF000118	Essence SPP Seat Cushion	Part series: ESS-###XXYY. HCPCS: E2607/E2608. MSRP: \$441-\$630.
Ref: XRF000113	Solution SPP Seat Cushion	Part series: SOL-###XXYY. HCPCS: E2607/E2608. MSRP: \$414-\$533.
Ref: XRF000112	Spectrum Foam SPP Seat Cushion	Part series: SPCF-###XXYY. HCPCS: E2607/E2608. MSRP: \$507-\$709.
Ref: XRF000117	Spectrum Gel SPP Seat Cushion	Part series: SPCG-###XXYY. HCPCS: E2607/E2608. MSRP: \$507-\$602.
Ref: XRF000122	Premiere P Seat Cushion	Part series: PREM-###XXYY. HCPCS: E2605/E2606. MSRP: \$358-\$603.
Ref: XRF000114	Glacial SP Seat Cushion	Part series: GLA-###XXYY. HCPCS: E2603/E2604. MSRP: \$364-\$564.
Ref: XRF000115	Zen SP Seat Cushion	Part series: ZEN-###XXYY. HCPCS: E2603/E2604. MSRP: \$348-\$547.
Ref: XRF000111	Simplicity G Seat Cushion	Part series: SIM-###XXYY. HCPCS: E2601/E2602. MSRP: \$122-\$182.

ESSENCE CUSHION WEDGE KIT

Only available with Essence cushion selection above.

POP1909006-S	Small Wedge Positioning Kit	14"W and smaller	\$39
POP1712340-S	Small Wedge Positioning Kit	15"W and larger	\$39
POP1909006-M	Medium Wedge Positioning Kit	14"W and smaller	\$72
POP1712340-M	Medium Wedge Positioning Kit	15"W and larger	\$72
POP1909006-L	Large Wedge Positioning Kit	14"W and smaller	\$105
POP1712340-L	Large Wedge Positioning Kit	15"W and larger	\$105

CUSHION TOP FOAM

Ref: XRF000380	Visco Top Foam
Ref: XRF000381	FRG Top Foam
Ref: XRF000382	No Top Foam

CUSHION LINER

Ref: XRF000383	Stealth Urethane Liner
Ref: XRF000384	Stealth Polyurethane Liner
Ref: XRF000385	No Liner

CUSHION COVER

Ref: XRF000386	CoolCore Cover
Ref: XRF000387	Mesh Cover
Ref: XRF000388	Stealth Tek Cover

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

OPTIONAL BASE INSERT

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

Ref: XRF000389	ABS Flat Rigidizer	MSRP: \$48-\$69
----------------	--------------------	-----------------

Account Number:

Date:

Quote Name/Marked for:

ADDITIONAL CUSHION LINER**Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.**

Ref: XRF000390 Additional Stealth Urethane Liner Quantity: _____

Ref: XRF000391 Additional Stealth Polyurethane Liner Quantity: _____

ADDITIONAL CUSHION COVER**Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.**

Ref: XRF000392 Additional CoolCore Cover Quantity: _____

Ref: XRF000393 Additional Mesh Cover Quantity: _____

Ref: XRF000394 Additional Stealth Tek Cover Quantity: _____

BACK TYPE**REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.**

XRF000014	Center Seat Spine	Not available with Q Static Basic. Select a back cushion configuration in the subsequent section.	No Charge
XRF000015	Back Cane Option	Select cane type, back cane height, and back cushion configuration below.	No Charge
XRF000016	Flip Forward Back	Not available with Q Static Basic. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.	\$188

BACK CANE SELECTION**REQUIRED when Back Cane Option is selected above. Select type and height of back canes.****BACK CANE TYPE**

10123	Angle Adjustable Back Cane		No Charge
XRF000274	Flip Forward Back Cane	Available with Q Static Basic	\$188
XRF000273	Flip Forward Back Cane	Available with Q Static Positioning	\$54
XRF000270	2" Endomorph (1" per side) Back Cane	Available with Q Static Basic	\$209
XRF000272	4" Endomorph (2" per side) Back Cane	Available with Q Static Basic	\$209
XRF000269	2" Endomorph (1" per side) Back Cane	Available with Q Static Positioning	\$75
XRF000271	4" Endomorph (2" per side) Back Cane	Available with Q Static Positioning	\$75

BACK CANE HEIGHT

10127	16" Back Cane Height		No Charge
10128	17" Back Cane Height		No Charge
10131	20" Back Cane Height		No Charge
10132	21" Back Cane Height		No Charge
10135	24" Back Cane Height		No Charge
10136	25" Back Cane Height		No Charge

BACK CUSHION CONFIGURATION**REQUIRED. Select one type of back configuration. See the back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.****SPORT BACK/OMIT/AFTER-MARKET INTERFACE**

XRF000361	Sport Back Cushion, 21" Wide or less	Available with Q Static Basic. Complete following width package and size.	\$402
XRF000362	Sport Back Cushion, 22" Wide or more	Available with Q Static Basic. Complete following width package and size.	\$616
XRF000017	Sport Back	Available with Q Static Positioning. Complete following width package and size.	No Charge
PTOINDV1287	Omit Back Pan/Cushion	Select a following width package.	No Charge
ST-SSAPT4-Q	After-market Back Interface	Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.	\$535

Account Number:

Date:

Quote Name/Marked for:

TRU-COMFORT BACK

TC4V10W21W	TRU-Comfort Back, Less than 22" W	Center Spine or Flip Forward Back compatibility. Complete following width package and size.	E2620	\$759
TC4V10W21WBC	TRU-Comfort Back, Back Cane Option, Less than 22" W	Back Cane compatibility. Complete following width package and size.	E2620	\$759
TC4V22W32W	TRU-Comfort Back, 22" W or Greater	Center Spine or Flip Forward Back compatibility. Complete following width package and size.	E2621	\$872
TC4V22W32WBC	TRU-Comfort Back, Back Cane Option, 22" W or Greater	Back Cane compatibility. Complete following width package and size.	E2621	\$872

Overall Height	Back Width											
	10"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
14"	P	T	P	T	P	T	-	-	-	-	-	-
15"	P	T	P	T	P	T	-	-	-	-	-	-
16"	P	T	P	T	P	T	P	T	P	T	P	T
17"	P	T	P	T	P	T	P	T	P	T	P	T
18"	P	P	P	T	P	T	P	T	P	T	P	T
19"	P	P	P	T	P	T	P	T	P	T	P	T
20"	-	-	P	T	P	T	P	T	P	T	P	T
21"	-	-	P	T	P	T	P	T	P	T	P	T
22"	-	-	P	T	P	T	P	T	P	T	P	T
23"	-	-	P	T	P	T	P	T	P	T	P	T
24"	-	-	P	T	P	T	P	T	P	T	P	T
25"	-	-	P	T	P	T	P	T	P	T	P	T
26"	-	-	P	T	P	T	P	T	P	T	P	T
27"	-	-	P	T	P	T	P	T	P	T	P	T
28"	-	-	P	T	P	T	P	T	P	T	P	T
29"	-	-	P	T	P	T	P	T	P	T	P	T

P Sport Back
 T TRU-Comfort Back

Seat Back Availability Chart

Q-BACK 2.0, AL SERIES, ACTIVE CONTOUR

Select a back type and the specific overall back height in subsequent section.

QB2AL14AS	Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL14AT	Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL16AS	Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL16AT	Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL16AX	Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL16AXX	Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline	E2613	\$703
QB2AL18AS	Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL18AT	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL18AX	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL18AXX	Active Contour, 18"Wx25-26"H Tilt/27"-28" H Recline	E2613	\$703
QB2AL20AS	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline	E2613	\$568
QB2AL20AT	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline	E2613	\$568
QB2AL20AX	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline	E2613	\$568
QB2AL20AXX	Active Contour, 20"Wx25-26"H Tilt/27"-28" H Recline	E2613	\$703

Q-BACK 2.0, AL SERIES, DEEP CONTOUR

Select a back type and the specific overall back height in subsequent section.

QB2AL14DS	Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL14DT	Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643
QB2AL16DS	Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline	E2620	\$643
QB2AL16DT	Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline	E2620	\$643

Account Number:		Date:	Quote Name/Marked for:		
QB2AL16DX	Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$643	
QB2AL16DXX	Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
QB2AL18DS	Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$643	
QB2AL18DT	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$643	
QB2AL18DX	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$643	
QB2AL18DXX	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
QB2AL20DS	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$643	
QB2AL20DT	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$643	
QB2AL20DX	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$643	
QB2AL20DXX	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
Q-BACK 2.0, AL SERIES, EXTRA DEEP CONTOUR					
Select a back type and the specific overall back height in subsequent section.					
QB2AL14ES	Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$703	
QB2AL16ES	Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$703	
QB2AL16ET	Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$703	
QB2AL16EX	Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$703	
QB2AL18ES	Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$703	
QB2AL18ET	Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$703	
QB2AL18EX	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$703	
QB2AL18EXX	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
QB2AL20ES	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$703	
QB2AL20ET	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$703	
QB2AL20EX	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline		E2620	\$703	
QB2AL20EXX	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
Q-BACK 2.0, AFT SERIES, ACTIVE CONTOUR					
Select a back type and the specific overall back height in subsequent section.					
QB2AFT16AS	AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline		E2613	\$568	
QB2AFT18AS	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline		E2613	\$568	
Q-BACK 2.0, AFT SERIES, DEEP CONTOUR					
Select a back type and the specific overall back height in subsequent section.					
QB2AFT16DS	AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$643	
QB2AFT16DT	AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$643	
QB2AFT18DS	AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline		E2620	\$643	
QB2AFT18DT	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$643	
QB2AFT20DT	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline		E2620	\$643	
Q-BACK 2.0, AFT SERIES, EXTRA DEEP CONTOUR					
Select a back type and the specific overall back height in subsequent section.					
QB2AFT16EXX	AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline		E2620	\$703	
BACK SHROUD					
The back shroud is not available when Omit Back Pan/Cushion is selected.					
XRF000019	Red Corvette Back Shroud	No Charge	XRF000030	Electric Blue Back Shroud	No Charge
XRF000021	Back in Black Back Shroud	No Charge	XRF000032	Lemon Crush Back Shroud	No Charge
XRF000022	Night Train Back Shroud	No Charge	XRF000033	Green Machine Back Shroud	No Charge
XRF000023	Grey Street Back Shroud	No Charge	XRF000034	Military Green Back Shroud	No Charge
XRF000025	Deep Purple Back Shroud	No Charge			
XRF000027	Pretty in Pink Back Shroud	No Charge			
XRF000028	Evolution Orange Back Shroud	No Charge			
XRF000029	Midnight Blue Back Shroud	No Charge			

Account Number:

Date:

Quote Name/Marked for:

BACK WIDTH PACKAGE**REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.**

PTOINDV3468	Match Seat Width Selections	No Charge
XRF000036	10"-14" Back Width Package	No Charge
XRF000037	14"-16" Back Width Package	No Charge
XRF000038	16"-20" Back Width Package	No Charge
XRF000039	18"-22" Back Width Package	No Charge
XRF000040	20"-24" Back Width Package	No Charge

BACK WIDTH**REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.**

PTOINDV3561	Set Back Width at 10"W	No Charge	PTOINDV3566	Set Back Width at 18"W	No Charge
PTOINDV3562	Set Back Width at 12"W	No Charge	PTOINDV3567	Set Back Width at 19"W	No Charge
PTOINDV3563	Set Back Width at 14"W	No Charge	PTOINDV3568	Set Back Width at 20"W	No Charge
PTOINDV3856	Set Back Width at 15"W	No Charge	PTOINDV3569	Set Back Width at 21"W	No Charge
PTOINDV3564	Set Back Width at 16"W	No Charge	PTOINDV3570	Set Back Width at 22"W	No Charge
PTOINDV3565	Set Back Width at 17"W	No Charge			

OVERALL BACK HEIGHT**REQUIRED with Sport or TRU Comfort back selection. QStatic odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. See the back availability chart for stock order back size ranges.**

SETHGHT1027	Set Back Height to 14"H	No Charge	SETHGHT1035	Set Back Height to 22"H	No Charge
SETHGHT1028	Set Back Height to 15"H	No Charge	SETHGHT1036	Set Back Height to 23"H	No Charge
SETHGHT1029	Set Back Height to 16"H	No Charge	SETHGHT1037	Set Back Height to 24"H	No Charge
SETHGHT1030	Set Back Height to 17"H	No Charge	SETHGHT1038	Set Back Height to 25"H	No Charge
SETHGHT1031	Set Back Height to 18"H	No Charge	SETHGHT1039	Set Back Height to 26"H	No Charge
SETHGHT1032	Set Back Height to 19"H	No Charge	SETHGHT1040	Set Back Height to 27"H	No Charge
SETHGHT1033	Set Back Height to 20"H	No Charge	SETHGHT1041	Set Back Height to 28"H	No Charge
SETHGHT1034	Set Back Height to 21"H	No Charge	SETHGHT1042	Set Back Height to 29"H	No Charge

PRIVACY FLAP**Optional.**

XRF000351	Privacy Flap	\$49
-----------	--------------	------

OPTIONAL BACK HEIGHT EXTENSION**Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.**

XRF000041	Back Height Extension	\$49
-----------	-----------------------	------

ARMRESTS**REQUIRED. Select one type of armrest.**

10153	Drop-In, H Adj, Removable, T-Arm Armrests (n)	Includes straight armpads. Select armrest height and size.	No Charge
10152	Q Captain, Flip-Up, Height Adjustable, Full Length	Includes full length, straight armpads. Select armrest height.	K0020 \$180
10157	Pediatric, Removable Armrests	Includes desk length, straight armpads. Select height adjustment from 7.75" to 8.25".	E0973 \$162
10165	Cantilever, Height Adjustable Armrests	Not available with flip forward back canes. Select armrest height (11"-17"), size, and armpad.	\$105
10156	Two-Post, H Adj, Removable Armrest (u)	Select armrest height, size, and armpad.	E0973 \$322
10158	8"-12"H, Single Post, H Adj, Remov, HD Armrests	Select armrest height, size, and armpad.	E0973 \$161
10159	10"-14"H, Single Post, H Adj, Remov, HD Armrests	Select armrest height, size, and armpad.	E0973 \$161
10150	Flip-Back, 8"-16"H Armrests (u)	Not available with Back Cane Option. Select armrest height, size, and armpad (Not available with waterfall armpads)	\$590

Account Number:

Date:

Quote Name/Marked for:

ARMREST SUPPORT**Optional with Flip-Back, 8"-16"H Armrests.**

10323	Flip-Back Armrest Downposts (pair)	Not available less than 16" width seat or back.	\$188
-------	------------------------------------	---	-------

ARMREST HEIGHT**REQUIRED.****LEFT ARMREST HEIGHT**

PTOINDV3722	7.75" Left Armrest Height	No Charge
PTOINDV3535	8" Left Armrest Height	No Charge
PTOINDV3724	8.25" Left Armrest Height	No Charge
PTOINDV3536	9" Left Armrest Height	No Charge
PTOINDV3537	10" Left Armrest Height	No Charge
PTOINDV3538	11" Left Armrest Height	No Charge
PTOINDV3539	12" Left Armrest Height	No Charge
PTOINDV3540	13" Left Armrest Height	No Charge
PTOINDV3541	14" Left Armrest Height	No Charge
PTOINDV3542	15" Left Armrest Height	No Charge
PTOINDV3543	16" Left Armrest Height	No Charge

RIGHT ARMREST HEIGHT

PTOINDV3721	7.75" Right Armrest Height	No Charge
PTOINDV3450	8" Right Armrest Height	No Charge
PTOINDV3723	8.25" Right Armrest Height	No Charge
PTOINDV3451	9" Right Armrest Height	No Charge
PTOINDV3452	10" Right Armrest Height	No Charge
PTOINDV3453	11" Right Armrest Height	No Charge
PTOINDV3454	12" Right Armrest Height	No Charge
PTOINDV3455	13" Right Armrest Height	No Charge
PTOINDV3456	14" Right Armrest Height	No Charge
PTOINDV3457	15" Right Armrest Height	No Charge
PTOINDV3458	16" Right Armrest Height	No Charge

ARMREST SIZE**REQUIRED.**

XRF000057	10" Armrest Length, Left	No Charge
XRF000123	10" Armrest Length, Right	No Charge
XRF000058	14" Armrest Length, Left	No Charge
XRF000124	14" Armrest Length, Right	No Charge
XRF000059	17" Armrest Length, Left	Only available with Flip-Back, 8"-16" H Armrests and straight arm pads. No Charge
XRF000125	17" Armrest Length, Right	Only available with Flip-Back, 8"-16" H Armrests and straight arm pads. No Charge

ARMPADS**REQUIRED.**

XRF000060	No Armpad Selection	No Charge
XRF000308	Interface for Aftermarket Pad	No Charge
XRF000061	Straight, Armpad, Left	No Charge
XRF000126	Straight, Armpad, Right	No Charge
XRF000062	Waterfall, Armpad, Left	Not available with Flip-back, Q Captain, or Pediatric armrests. No Charge
XRF000127	Waterfall, Armpad, Right	Not available with Flip-back, Q Captain, or Pediatric armrests. No Charge

FALCON GEL ARMPADS FOR FLIP-BACK ARMRESTS**Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.**

ARM1806807	Falcon Gel, Full Length Left and Right	\$199.00
ARM1806813	Falcon Gel, Full Length Left, Desk Length Right	\$174.00
ARM1806810	Falcon Gel, Desk Length Left, Full Length Right	\$174.00
ARM1806806	Falcon Gel, Desk Length Left and Right	\$150.00

FALCON GEL ARMPADS FOR TUBE STYLE ARMRESTS**Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.**

ARM136438	Falcon Gel, Full Length Left and Right	\$317.00
ARM136442	Falcon Gel, Full Length Left, Desk Length Right	\$294.00
ARM136439	Falcon Gel, Desk Length Left, Full Length Right	\$294.00
ARM136437	Falcon Gel, Desk Length Left and Right	\$270.00

Account Number:

Date:

Quote Name/Marked for:

JOYSTICK BRACKET RECEIVER**REQUIRED.**

XRF000063	Standard Fixed Receiver	No Charge
XRF000064	Stealth Height Adjustable Receiver	K0108 \$236

HEADREST**Optional.****HEADREST PADS**

CP08CCSPM	Stealth 8" Comfort Plus, Cool Core (CP08CCSPM)	E0955	\$179
CP10CCSPM	Stealth 10" Comfort Plus, Cool Core (CP10CCSPM)	E0955	\$179
CP14CCSPM	Stealth 14" Comfort Plus, Cool Core (CP14CCSPM)	E0955	\$179
CP06CCSPM	Stealth 6", Comfort Plus, Cool Core (CP06CCSPM)	E0955	\$179

HEADREST BRACKETS

ST-UHB2D33200-Q	Comfort Plus Headrest Mounting Bracket, UniLink	E1033	\$311
ST-HMO475-17M-Q	Multi-Axis Removable Headrest Bracket	E1033	\$174

Q CAM LOCK

XRF000376	Q Cam Lock	\$51
-----------	------------	------

Standard	\$50
Tilt/Static/Back Cane	
TRU-Comfort/Sport	
14"-21"H 22"-29"H	
Q-Back 2	
19"-24"H 25"-26"H	
Power Recline Back	
TRU-Comfort/Sport	
19"-21"H 22"-32"H	
Q-Back 2	
22"-25"H 26"-29"H	
Q Cam Lock Avail. (overall height)	

LAP/ PELVIC BELTS**REQUIRED. Select one lap/pelvic belt below.**

10167	38" Lap Belt	No Charge
10168	50" Lap Belt	No Charge
10169	60" Lap Belt	No Charge
10170	70" Lap Belt	No Charge
10171	80" Lap Belt	No Charge
10172	90" Lap Belt	No Charge
ST-IFRB5350	Retractable Lap Belt	Stealth Positioning IFRB5350. \$181
ST-IFPB5213-Q	60" Single Pull, 1.5", Push Button, M Pads (7.5")	Stealth Positioning IFPB5213. \$50
ST-IFPB5314-Q	60" Single Pull, 2", Push Button, L Pads (9.5")	Stealth Positioning IFPB5214. \$50
ST-IFPB5315-Q	60" Single Pull, 2", Push Button, XL Pads (11.5")	Stealth Positioning IFPB5215. \$50
ACCBELT1016	60" Sure-Lok Lap Belt (Non-retractable)	Required with Compliant Transit Option. \$120

LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

XRF000108	Independent Legrests	Adds 2 actuator control requirement.	No Charge
XRF000109	Combined Legrests	Adds 2 actuator control requirement.	No Charge
XRF000110	Articulating Foot Platform	Adds 1 actuator control requirement.	No Charge

BASE MOUNTED FOOT PLATFORM

Optional. May interfere with seat mounted legrests. Not available with center mount foot platforms.

10185	Base Mounted Foot Platform	No Charge
10186	Base Mounted, Bell Shaped Foot Platform	\$252
PTOINDV3674	Omit Base Mounted Foot Platform	No Charge

FRONT RIGGINGS

Optional. Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

BASIC CENTER MOUNT**Basic Center Mount Foot Platform Assembly**

XRF000323	Basic Center Mount	Select extension and footplate. Calf pads available depending on configuration.	No Charge
-----------	--------------------	---	-----------

Account Number:

Date:

Quote Name/Marked for:

Basic Center Mount Extension and Angle Adjustable Footplates

XRF000312	Short Extension, 5"x9" Footplates	Not available with calf pads.	K0040	\$184
XRF000313	Medium Extension, 5"x9" Footplates	Requires small calf pads with power tilt configurations.	K0040	\$184
XRF000314	Long Extension, 5"x9" Footplates	Requires small or large calf pads with power tilt configurations.	K0040	\$184

Basic Center Mount Calf Pads

PTOINDV3768	No Calf Pads	No Charge	Short	Med.	Long
10194	Small Calf Pads, Basic Center Mount	\$161	8-10	10-14	14-19
10195	Large Calf Pads, Basic Center Mount	\$161	Basic Center Mount Extension Range (in.)		

Basic Center Mount Additional Accessory

10196	Basic Center Mount Footplate Connector	Optional.	\$16
-------	--	-----------	------

CENTER MOUNT FOOT PLATFORM**Center Mount Foot Platform Type**

PTOINDV3754	Power Articulating Foot Platform	Select extension with footplates and calf support. Other items are optional.	E1012+K0040	\$3,145	Power AFP Extension Range		
XRF000311	Quick Adjustable Center Mount Foot Platform	Select extension with footplates. Other items are optional.	K0040	\$268	Ped.	Med.	Long
					8.5-12.5	11-16.5	14-20
Extension and Angle Adjustable Footplates					Quick Adj Center Mnt Extension Range (in)		
10201	Pediatric Extension. 5"x8" Footplates	Not available with AFP		No Charge			

Extension and Angle Adjustable Footplates

10201	Pediatric Extension, 5"x8" Footplates	Not available with AFP.	No Charge
10370	Medium Extension, 5"x8" Footplates		No Charge
10371	Medium Extension, 5"x11" Footplates		No Charge
10372	Long Extension, 5"x8" Footplates		No Charge
10373	Long Extension, 5"x11" Footplates		No Charge

Center Mount Calf Pads

10206	Small Calf Pads, AFP Only, 6"Wx6"H (pair)	Small or large calf pad standard with AFP.	No Charge
10357	Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair)	Small or large calf pad standard with AFP.	No Charge
XRF000156	Small Calf Pads, Adjustable Center Mount, 6"Wx6"H (pair)	Available with adjustable center mount.	\$161
XRF000157	Large Calf Pads, Adjustable Center Mount, 8" taper to 6"Wx8"H (pair)	Available with adjustable center mount.	\$161

Center Mount Calf Panel

10207	Small Calf Panel, 18"Wx8"H	\$188
10208	Large Calf Panel, 24"Wx8"H	Only available on widths of 20" or greater. \$242

Optional Center Mount Modular Footplate Extension Components**Assembly required. May require seating configuration adjustments to avoid front caster interference.**

10209	Width Extension	\$70
10210	Outside Edge	\$70
10211	Front and Rear Cap	\$16

Optional Center Mount Additional Accessories

ST-CSCS0775-Q	Small Foam Calf Pad Cover (pair)	\$67
ST-CSCS0985-Q	Large Foam Calf Pad Cover (pair)	\$83
ST-CSCP0920-Q	Small Foam Panel Cover (pair)	\$176
ST-CSCP0927-Q	Large Foam Panel Cover (pair)	\$282
10212	Footplate Connector	\$16
XRF000348	Footplate Riser	Not available with Long Extension. Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments. \$77
10213	Footplate Riser Calf Panel	Only available and required with Footplate Riser selection and extension lengths of 6" to 9". \$51

Account Number:

Date:

Quote Name/Marked for:

Optional Stealth Gel Footplate Covers

ST-FPG5008T-Q	5"x8" Tapered Gel Footplate Covers	\$222
ST-AFPG5011T-Q	5"x11" Tapered Gel Footplate Covers	\$270
ST-AFPG5008E-Q	5"x8" Tapered Gel Footplate Covers for Width Extension	\$247
ST-AFPG5011E-Q	5"x11" Tapered Gel Footplate Cover for Width Extension	\$270

SWING-AWAY LEGRESTS**60°-90° Swing-Away Legrest Upper**

XRF000082	60°-90° Swing-Away Legrest, Static Seat, Pair	Select pair ext. plus pair footplate or amputee support(s).	\$176
XRF000085	60°-90° Swing-Away Legrest, Short, Static Seat, P	Select pair ext. plus pair footplate or amputee support(s).	\$176

Lower Extension

Rigging configurations may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability when noted below.

XRF000088	11"-13" Length, Pair	See note above with recline plus low STF.	No Charge
XRF000094	13"-16.5" Length, Pair	See note above with recline plus low STF.	No Charge
XRF000097	15.5"-19.5" Length, Pair	See note above.	No Charge

Angle Adjustable Footplate(s)

10222	5"Wx8"D, Ang Adj Ftpts, Heel Loops, S/A, Pair	K0040+ E0951	\$184
10225	6"Wx8"D, Ang Adj Ftpts, Heel Loops, S/A, Pair	K0040+ E0951	\$184
10228	6"Wx11"D, Ang Adj Ftpts, Heel Loops, S/A, Pair	K0040+ E0951	\$184
10231	8"Wx11"D, Ang Adj Ftpts, Heel Loops, S/A, Pair	K0040+ E0951	\$184

Amputee Supports with Swing-Away Upper

INDPART2822	Therafin 31515 Multi-Axis Amputee Support	Quantity: _____	E1020	\$259
INDPART2823	Therafin 32387 Amputee Support Gel Cover	Quantity: _____	K0108	\$183

POWER ELEVATING LEGRESTS**Power Elevating Legrest Upper Assembly**

Requires power legrest control selection found in earlier section.

XRF000341	Power Elevating Legrests, Pair	Choose extension and pair footplates.	E1010+ K0040+ E0951	\$2,516
-----------	--------------------------------	---------------------------------------	---------------------	---------

Lower Extension

XRF000344	8.5" to 14.5" Lower Length	No Charge
XRF000347	12.5" to 21" Lower Length	No Charge

Angle Adjustable Footplate(s)

The 6"W footplates may cause interference with seat widths less than 16". The 8"W footplates may cause interference with seat widths less than 19".

10216	4"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10217	5"Wx6"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10218	5"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10219	6"Wx8"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10220	6"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge
10221	8"Wx11"D, AA Ftplts, Heel Loops, Wheel Bumpers, Pr	No Charge

Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High
Power ELR: Q Static Seat			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14	12.5"-15"	12.5"-16"

Account Number:

Date:

Quote Name/Marked for:

ELEVATING LEGRESTS**An upper assembly, one extension or limb support, and one type of footplate is required to configure elevating legrests.****Elevating Legrest Upper Assembly**

10234	Elevating Legrest Assembly, Pair	E0990	\$342
-------	----------------------------------	-------	-------

Elevating Legrest Lower Extension

XRF000344	8.5" to 14.5" Lower Length	No Charge
-----------	----------------------------	-----------

XRF000347	12.5" to 21" Lower Length	No Charge
-----------	---------------------------	-----------

Elevating Legrest Residual Limb Support

FRMASMB5301	Left Residual Limb Support with Calf Pad	E1020	\$342
-------------	--	-------	-------

FRMASMB5300	Right Residual Limb Support with Calf Pad	E1020	\$342
-------------	---	-------	-------

Elevating Legrest Angle Adjustable Footplate(s)

10255	5"Wx8"D, Ang Adj Ftpts, Heel Loops, ELR, Pair	K0040+ E0951	\$184
-------	---	--------------	-------

10258	6"Wx8"D, Ang Adj Ftpts, Heel Loops, ELR, Pair	K0040+ E0951	\$184
-------	---	--------------	-------

10261	6"Wx11"D, Ang Adj Ftpts, Heel Loops, ELR, Pair	K0040+ E0951	\$184
-------	--	--------------	-------

10264	8"Wx11"D, Ang Adj Ftpts, Heel Loops, ELR, Pair	K0040+ E0951	\$184
-------	--	--------------	-------

Elevating Legrest Accessories

INDPART2242	Falcon Gel Padded Knee Adductor Buttons (KBE) (Pr)	E0953	\$134
-------------	--	-------	-------

LEGREST ACCESSORIES

10267	Heel Loops (Pair)	Standard on angle adjustable.	E0951	\$54
-------	-------------------	-------------------------------	-------	------

XRF000100	Wheel Bumper(s)	Pair of bumpers.	\$38
-----------	-----------------	------------------	------

10268	Wheel Bumper	One bumper.	\$19
-------	--------------	-------------	------

SET LEGREST EXTENSION LENGTH**REQUIRED with center mounts and matching left and right rigging configurations. Select one length.****CENTERMOUNT EXTENSION LENGTH**

PTOINDV3989	2" L Legrest Length	No Charge	PTOINDV3689	12.5" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	------------------------	-----------

PTOINDV3990	3" L Legrest Length	No Charge	PTOINDV3690	13" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	----------------------	-----------

PTOINDV3991	4" L Legrest Length	No Charge	PTOINDV3691	13.5" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	------------------------	-----------

PTOINDV3992	5" L Legrest Length	No Charge	PTOINDV3692	14" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	----------------------	-----------

PTOINDV3993	6" L Legrest Length	No Charge	PTOINDV3693	14.5" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	------------------------	-----------

PTOINDV3994	7" L Legrest Length	No Charge	PTOINDV3694	15" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	----------------------	-----------

PTOINDV3680	8" L Legrest Length	No Charge	PTOINDV3695	15.5" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	------------------------	-----------

PTOINDV3681	8.5" L Legrest Length	No Charge	PTOINDV3696	16" L Legrest Length	No Charge
-------------	-----------------------	-----------	-------------	----------------------	-----------

PTOINDV3682	9" L Legrest Length	No Charge	PTOINDV3697	16.5" L Legrest Length	No Charge
-------------	---------------------	-----------	-------------	------------------------	-----------

PTOINDV3683	9.5" L Legrest Length	No Charge	PTOINDV3698	17" L Legrest Length	No Charge
-------------	-----------------------	-----------	-------------	----------------------	-----------

PTOINDV3684	10" L Legrest Length	No Charge	PTOINDV3699	17.5" L Legrest Length	No Charge
-------------	----------------------	-----------	-------------	------------------------	-----------

PTOINDV3685	10.5" L Legrest Length	No Charge	PTOINDV3700	18" L Legrest Length	No Charge
-------------	------------------------	-----------	-------------	----------------------	-----------

PTOINDV3686	11" L Legrest Length	No Charge	PTOINDV3701	18.5" L Legrest Length	No Charge
-------------	----------------------	-----------	-------------	------------------------	-----------

PTOINDV3687	11.5" L Legrest Length	No Charge	PTOINDV3702	19" L Legrest Length	No Charge
-------------	------------------------	-----------	-------------	----------------------	-----------

PTOINDV3688	12" L Legrest Length	No Charge			
-------------	----------------------	-----------	--	--	--

Account Number:

Date:

Quote Name/Marked for:

FRONT RIGGING EXTENSION LENGTH

PTOINDV3470	7" L Legrest Length	No Charge	PTOINDV3486	15" L Legrest Length	No Charge
PTOINDV3471	7.5" L Legrest Length	No Charge	PTOINDV3487	15.5" L Legrest Length	No Charge
PTOINDV3472	8" L Legrest Length	No Charge	PTOINDV3488	16" L Legrest Length	No Charge
PTOINDV3473	8.5" L Legrest Length	No Charge	PTOINDV3489	16.5" L Legrest Length	No Charge
PTOINDV3474	9" L Legrest Length	No Charge	PTOINDV3490	17" L Legrest Length	No Charge
PTOINDV3475	9.5" L Legrest Length	No Charge	PTOINDV3491	17.5" L Legrest Length	No Charge
PTOINDV3476	10" L Legrest Length	No Charge	PTOINDV3492	18" L Legrest Length	No Charge
PTOINDV3477	10.5" L Legrest Length	No Charge	PTOINDV3493	18.5" L Legrest Length	No Charge
PTOINDV3478	11" L Legrest Length	No Charge	PTOINDV3494	19" L Legrest Length	No Charge
PTOINDV3479	11.5" L Legrest Length	No Charge	PTOINDV3495	19.5" L Legrest Length	No Charge
PTOINDV3480	12" L Legrest Length	No Charge	PTOINDV3496	20" L Legrest Length	No Charge
PTOINDV3481	12.5" L Legrest Length	No Charge	PTOINDV3497	20.5" L Legrest Length	No Charge
PTOINDV3482	13" L Legrest Length	No Charge	PTOINDV3498	21" L Legrest Length	No Charge
PTOINDV3483	13.5" L Legrest Length	No Charge	PTOINDV3600	21.5" L Legrest Length	No Charge
PTOINDV3484	14" L Legrest Length	No Charge	PTOINDV3669	22" L Legrest Length	No Charge
PTOINDV3485	14.5" L Legrest Length	No Charge	PTOINDV3670	22.5" L Legrest Length	No Charge

SEATING ELECTRONICS**REQUIRED with power seating or legrest options. Select one type of power actuator control.****ACTUATOR CONTROL THROUGH JOYSTICK**

XRF000335	Single Actuator Through Controller and Interface	E2310	\$2,206
XRF000336	Multiple Actuator Through Controller and Interface	E2311	\$3,101


ATTENDANT CONTROL**Optional. Mounting bracket is included.**

XRF000078	Attendant Control	E2331	\$536
-----------	-------------------	-------	-------

SEAT MOUNTED LED LIGHTS**Optional. Requires QLNE5 10 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.**

XRF000130	Seat Mounted, LED, Full Light Package	Not available with Base Mounted, LED, Full Light Package.	\$1,045
-----------	---------------------------------------	---	---------

ACCESSORIES**Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.**

XRF000128	Accessory Bar	Horizontal rail attached to rear of seat. Included with certain accessories when noted.	\$135	Backpack Holder	
10375	Accessory Hook	Coming soon.	\$29	Clothing Guard	
10326	Backpack Holder		\$70	Glove Box	
10327	Cane/Crutch Holder		E2207 \$107	Phone Holder	
10331	Cell Phone/ Media Holder, X-Grip	Not available with Heavy Duty armrests.	\$70	Hydration System	
10332	Clothing Guard		\$242	Personal Item Hook	
10334	Cup Holder		No Charge		
10337	Flag Kit		\$70		
10338	Flip-Back Armrest, Side Track Accessory Bracket		\$11		
10340	Glove Box		\$161		
10342	Hydration System		\$242		
10343	Large, Black, Weather Cover		\$202		
10344	Oxygen Holder	Includes accessory bar, a credit will be applied PTOINDV3769.	E2208 \$214		
10348	Personal Item Hook		\$70		

Account Number: Date: Quote Name/Marked for:

10349	Personal Item Pouch		\$86
10350	Privacy Bag		\$43
10351	Push Handles	Not available with oxygen holder or vent tray.	\$199
10352	Quantum Backpack		\$171
10353	Rear View Mirror		\$70
10354	Reflector Kit		\$54
XRF000379	Remote Stop Switch		\$1,321
ACC7120006	Stealth Cell Phone/Media Device Holder		\$255
ACC141818	Stealth Large Tablet Holder		\$625
XRF000396	Transfer Handles (pair)	E0705	\$214
10324	Aftermarket Accessory Mounting Bracket, Non-slotted	Quantity: __ Select a quantity up to 5.	\$54
10325	Aftermarket Accessory Mounting Bracket, Slotted	Quantity: __ Select a quantity up to 5.	\$54



BRANCH PATCHES

Optional.

ACC140240	Air Force Patch	\$21
ACC140241	Army Patch	\$21
ACC140239	Coast Guard Patch	\$21
ACC140243	Marines Patch	\$21
ACC140244	Navy Patch	\$21
ACC140238	Purple Heart Patch	\$21



Branch Patches

POSITIONING COMPONENTS

Optional. More options are available on the Positioning Components Order Form. Please provide quantity where applicable.

STEALTH THORACIC LATERALS

Swing-away laterals may interfere with flip-back and 2-post, flip up armrests.

Stealth Thoracic Lateral Bracket

ST-ULS4042000-Q	Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)	Only available with Redesigned back assemblies.	E1034	\$664
ST-ULPK404108B00-Q	Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00)	Only available with Sport or TRU Comfort backs.	E1034	\$664
ST-ULPK404104B00-Q	Swing-away, UniLink, Back Cane, pair (ULPK404104B00)	Only available with back cane and omit back configurations.	E1034	\$664
ST-ULPK404103L00-Q	Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00)	Only available with Active Contour Q-Backs.	E1034	\$332
ST-ULPK404103R00-Q	Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00)	Only available with Active Contour Q-Backs.	E1034	\$332

Stealth Thoracic Lateral Pads

INDPART2663	Curved Pads, 3"x5", pair (TWBL-3X5)	E0956	\$250
INDPART2664	Curved Pads, 4"x6", pair (TWBL-4X6)	E0956	\$250
ST-ULP846050320-Q	Flat Pads, 3"x5", pair (ULP846050320)	E0956	\$250
ST-ULP835050320-Q	Flat Pads, 4"x6", pair (ULP835050320)	E0956	\$250

Stealth Thoracic Lateral Installation

PTOINDV3499	Install Laterals	No Charge
PTOINDV3500	Do Not Install	No Charge

Account Number:

Date:

Quote Name/Marked for:

STEALTH UNILINK PELVIC/THIGH GUIDE**UniLink Pelvic/Thigh Hardware****^ Please verify the appropriate billing code depending upon the placement of the support pad.**

ST-UTTP5-Q	UniLink, 5" Height Adjustable & Removable Hardware	Quantity: ____	E1028/ E1034^	\$265
ST-UTTP7-Q	UniLink, 7" Height Adjustable & Removable Hardware	Quantity: ____	E1028/ E1034^	\$276
ST-UTT8A11602-L-Q	UniLink Multi-Angle Thigh Guide, Left		E1028/ E1034^	\$359
ST-UTT8A11602-R-Q	UniLink Multi-Angle Thigh Guide, Right		E1028/ E1034^	\$359
ST-UTS25100-L-Q	UniLink Flip Down Hip/Thigh Support, Left		E1028/ E1034^	\$332
ST-UTS25100-R-Q	UniLink Flip Down Hip/Thigh Support, Right		E1028/ E1034^	\$332
POS140428	Flip Down with LEFT Full Surface Contact Feature (TWBADD-L & LPHW407)		E1028/ E1034^	\$281
POS140429	Flip Down with RIGHT Full Surface Contact Feature (TWBADD-R & LPHW407)		E1028/ E1034^	\$281

UniLink Pelvic/Thigh Pads**^ Please verify the appropriate billing code to be used with all third-party payors.**

ST-UTP544056100-Q	UniLink, 4" x 4" Pad	Quantity: ____	E0953/ E0956^	\$106
ST-UTP548056100-Q	UniLink, 4" x 8" Pad	Quantity: ____	E0953/ E0956^	\$106
ST-UTP54D056100-Q	UniLink, 4" x 12" Pad	Quantity: ____	E0953/ E0956^	\$106

UniLink Pelvic/Thigh Additional Links

ST-UIL001150A000-Q	UniLink, 1.5" Link	Quantity: ____		\$81
ST-UIL001200A000-Q	UniLink, 2" Link	Quantity: ____		\$81

UniLink Pelvic/Thigh Guide Installation

PTOINDV3585	Install Left Lateral Pelvic/ Thigh Bracket		No Charge
PTOINDV3586	Install Right Lateral Pelvic/ Thigh Bracket		No Charge
PTOINDV3719	Do Not Install Lateral Pelvic/ Thigh Bracket		No Charge

STEALTH POSITIONING CHEST HARNESS

ST-IFPS1103-Q	Upper Classic, Medium (IFPS1103)	E0960	\$154
ST-IFPS1104-Q	Upper Classic, Large (IFPS1104)	E0960	\$154
ST-IFPS2103-Q	Upper Contour, Medium (IFPS2103)	E0960	\$154
ST-IFPS2104-Q	Upper Contour, Large (IFPS2104)	E0960	\$154

**STEALTH POSITIONING CHEST STRAPS**

ST-IFCS4113-Q	Chest Strap, Medium w/ Belt Anchor (IFCS4113)	E0960	\$97
ST-IFCS4114-Q	Chest Strap, Large w/ Belt Anchor (IFCS4114)	E0960	\$97

**STEALTH SWING-AWAY HEADREST MOUNTS****These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.****Facial Lateral Hardware**

ST-SAR2A-L-Q	UniLink Facial Lateral, Swing Away, Left	E1028	\$318
ST-SAR2A-R-Q	UniLink Facial Lateral, Swing Away, Right	E1028	\$318
ST-SAR4A-DBU-LL-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper & 1 Lower Left	E1028	\$492
ST-SAR4A-DBU-RL-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper & 1 Lower Right	E1028	\$492
ST-SAR4A-DBU-Q	Unilink Facial Laterals, Dual Swing Away, 2 Upper	E1028	\$434



Account Number:

Date:

Quote Name/Marked for:

Facial Lateral Pads and Accessories

INDPART2815	Egg Switch Mount (ESM)	Quantity: ____	\$57
ST-900M-Q	Facial Pad, 4.5"x2", Cool Core (900M-Q)	Quantity: ____ K0108	\$61
ST-920M-Q	Facial Pad, 6"x2", Cool Core (920M-Q)	Quantity: ____ K0108	\$61
ST-SUS18-Q	Swing-Away Joystick Mount, 18" Rod	Quantity: ____	\$276

Allows mounting of facial pads and attaches to the headrest.

**Headrest Accessories****STEALTH ELBOW SUPPORT****Elbow Support Bracket**

POS124064	Left, Swing-Away, Elbow Support Bracket Tube Style Arms (ARMS-ES-L)	E1028	\$199
POS124066	Right, Swing-Away, Elbow Support Bracket Tube Style Arms (ARMS-ES-R)	E1028	\$199
ST-UETB4021FL-Q	Left, UniLink, Fixed, Elbow Support Bracket for Flip-Back Armrests		\$211
ST-UETB4021FR-Q	Right, UniLink, Fixed, Elbow Support Bracket for Flip-Back Armrests		\$211

**Elbow Support****Elbow Support Pad**

POS124065	3"Wx3"H Elbow Pad (ARMS-ES-3X3)	Quantity: ____ K0108	\$57
POS124067	4"Wx4"H Elbow Pad (ARMS-ES-4X4)	Quantity: ____ K0108	\$57

STEALTH ARM STOP**Arm Stop Hardware and Pad**

ST-ASF6X2C-S-Q	Stealth Arm Stop	Quantity: ____	\$351
----------------	------------------	----------------	-------

Arm Stop Location

PTOINDV3936	Mount Left Arm Stop On Inside	No Charge
PTOINDV3937	Mount Right Arm Stop On Inside	No Charge
PTOINDV3938	Mount Left Arm Stop On Outside	No Charge
PTOINDV3939	Mount Right Arm Stop On Outside	No Charge

STEALTH FOOT BOXES

INDPART2786	Stealth Small Foot Boxes, 5"Wx8"Dx7"H	E0954	\$539
INDPART2787	Stealth Medium Foot Boxes, 6"Wx10"Dx9"H	E0954	\$539
INDPART2788	Stealth Large Foot Boxes, 7"Wx12"Dx11"H	E0954	\$539

STEALTH POSITIONING SHOE HOLDER

ST-IFSH4211-Q	Small Shoe Holder, 3/4" Webbing (IFSH4211)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87
ST-IFSH4212-Q	Medium Shoe Holder, 3/4" Webbing (IFSH4212)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87
ST-IFSH4213-Q	Large Shoe Holder, 3/4" Webbing (IFSH4213)	Quantity: ____ Part and MSRP is 1 holder, 2 required for pair footplates.	E0951+ E0952	\$87

QUANTUM LAP TRAY**Optional. Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.****BLACK QUANTUM LAP TRAY**

MEC1905328-11B	Black Lap Tray, 11" W	For use with 12"-13" back width.	E0950	\$204
MEC1905328-13B	Black Lap Tray, 13" W	For use with 14"-15" back width.	E0950	\$204
MEC1905328-15B	Black Lap Tray, 15" W	For use with 16"-17" back width.	E0950	\$204
MEC1905328-17B	Black Lap Tray, 17" W	For use with 18"-19" back width.	E0950	\$204
MEC1905328-19B	Black Lap Tray, 19" W	For use with 20"-22" back width.	E0950	\$204

CLEAR QUANTUM LAP TRAY

MEC1905328-11C	Clear Lap Tray, 11" W	For use with 12"-13" back width.	E0950	\$204
MEC1905328-13C	Clear Lap Tray, 13" W	For use with 14"-15" back width.	E0950	\$204

Account Number:	Date:	Quote Name/Marked for:		
MEC1905328-15C	Clear Lap Tray, 15" W	For use with 16"-17" back width.	E0950	\$204
MEC1905328-17C	Clear Lap Tray, 17" W	For use with 18"-19" back width.	E0950	\$204
MEC1905328-19C	Clear Lap Tray, 19" W	For use with 20"-22" back width.	E0950	\$204
QUANTUM LAP TRAY MOUNTING HARDWARE				
MEC2104886	Mounting Bracket, Chest Bar/Lap Tray			\$77

SIMPLE SOLUTIONS LAP TRAY

Optional. Mounting bracket is required with lap tray selection. See chart below for compatibility.

SIMPLE SOLUTIONS LAP TRAY

ST-SSLXS010001171-Q	Simple Solutions Small Tray with Cut Out	E0950	\$309
ST-SSLXM010001171-Q	Simple Solutions Medium Tray with Cut Out	E0950	\$309
ST-SSLXL010001171-Q	Simple Solutions Large Tray with Cut Out	E0950	\$309

SIMPLE SOLUTIONS LAP TRAY MOUNTING HARDWARE

ST-SSLXHD10000060-Q	Swing Away Mounting Hardware		\$302
---------------------	------------------------------	--	-------

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible ■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.