

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Quantum 4Front™ Order Form with TRU-Balance 3 Positioning**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____


Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

**INTRODUCTION**

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24.2"Wx39.9"L

Quantum 4Front 3MP-SS..... \$7,895

Part: 4FRONT 3MP-SS. HCPCS: K0861

Multiple actuator power option.

Quantum 4Front 3SP-SS..... \$7,895

Part: 4FRONT 3SP-SS. HCPCS: K0856

Single power option.

Quantum 4Front 3S-SS..... \$7,895

Part: 4FRONT 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

2.1a. Standard Shrouds

 Red Corvette Part: DGN1700109-PRD1	 White Satin Part: DGN1700109-PWT2	 Midnight Blue Part: DGN1700109-PBU10	 Military Green Part: DGN1700109-PGN8
 Cherry Bomb Part: DGN1700109-PRD4	 Deep Purple Part: DGN1700109-PPR5	 Electric Blue Part: DGN1700109-PBU9	 Sugar Plum Part: DGN1700109-PPR7
 Back in Black Part: DGN1700109-PBK1	 Raspberry Beret Part: DGN1700109-PPK4	 Blue Jean Part: DGN1700109-PBU11	
 Night Train Part: DGN1700109-PBK3	 Pretty in Pink Part: DGN1700109-PPK3	 Lemon Crush Part: DGN1700109-PYL4	
 Grey Street Part: DGN1700109-PSL3	 Evolution Orange Part: DGN1700109-POR7	 Green Machine Part: DGN1700109-PGN11	

2.1b. Hydrographic Shrouds

☐  Military Digital Camo \$250 Part: DGN1700109-WCM7	☐  Faded Glory \$250 Part: DGN1700109-WFL2	☐  Carbon Fiber \$250 Part: DGN1700109-WCB1
--	--	---

3. NO FRONT ANTI-TIP OPTION

Optional. With this option, maximum speed will be 5 mph and 2.7 mph when elevated, maximum weight will be 200 lbs., and seat towers will be fixed in lowest position.

No Front Anti-Tips	No Charge
Part: PTOINDV3851	

4. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Operated Fender Lights..... No Charge Part: DGN1700033 Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.	Base Mounted, LED, Full Light Package \$995 Part: ELE1700035 Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.
Omit Fender Lights	No Charge
Part: DGN1700032	

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera	\$350
Part: MEC1910365	
Requires Q-Logic 3 electronics.	

1. Optional Mounting Bracket

Backup Camera Headrest Bracket.....	\$20
Part: MEC1910373	
Works with .5" x .5" vertical supports only.	



6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color.

6.1. Drive Wheels

3"x14" Flat Free Black Tires with Accent..... Standard Part: PTOINDV3525	4"x14" Flat Free Knobby Black Tires..... \$275.50 ea. Part: WHL2006822 x2
3"x14" Pneumatic Black Tires with Accent..... No Charge Part: WHL1704059 x2	4"x14" Pneumatic Knobby Black Tires..... \$275.50 ea. Part: WHL2006647 x2

6.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in BlackStandard Part: DGN164245-PBK1 x2	Raspberry BeretNo Charge Part: DGN164245-PPK4 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red CorvetteNo Charge Part: DGN164245-PRD1 x2	Pretty in PinkNo Charge Part: DGN164245-PPK3 x2	Military GreenNo Charge Part: DGN164245-PGN8 x2
Cherry BombNo Charge Part: DGN164245-PRD4 x2	Evolution OrangeNo Charge Part: DGN164245-POR7 x2	Sugar PlumNo Charge Part: DGN164245-PPR7 x2
Night TrainNo Charge Part: DGN164245-PBK3 x2	Midnight BlueNo Charge Part: DGN164245-PBU10 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Grey StreetNo Charge Part: DGN164245-PSL3 x2	Electric BlueNo Charge Part: DGN164245-PBU9 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
White SatinNo Charge Part: DGN164245-PWH2 x2	Blue JeanNo Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
Deep PurpleNo Charge Part: DGN164245-PPR5 x2	Lemon CrushNo Charge Part: DGN164245-PYL4 x2	

Account Number: _____ Date: _____ Marked for: _____

7. BATTERIES

Optional. The powerchair requires 2 flat top batteries and will be installed if ordered.

Group 34, AGM Cell Battery (Interceptor)..... \$350 ea.
Part: BAT1708797. HCPCS: E2359

Group 24, AGM Cell Battery (MK)..... \$495 ea.
Part: BAT1712207. HCPCS: E2363

Group 34, Gel Cell Battery (MK) \$465 ea.
Part: BATGEL1006. HCPCS: E2359

Not available with TRU-Balance 3 Power Lift.

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight
Part: PTOINDV3837

201 lb-300 lb User Weight
Part: PTOINDV3838

9. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

9.1a. Omit Joystick

No Joystick, Setup for Q-Logic Electronics
Part: XRF000010

9.1b. Q-Logic 3 Non-Expandable Capacity Controls

Q-Logic 3 non-expandable capacity allows control up to one actuator.

Q-Logic 3 Non-Expandable Capacity Joystick Standard
Part: CTL167724-NE

Q-Logic 3 Non-Expandable Capacity Joystick No Charge
for Full Light Package
Part: CTL167725-NE

Q-Logic 3e Non-Expandable Capacity Joystick No Charge
Part: CTL1910600-NE

Can operate full light package. Not available with alternative or specialty controls.

9.1c. Q-Logic 3 Expandable Capacity Controls

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Expandable Capacity Joysticks

Q-Logic 3 Expandable Capacity Joystick \$800
Part: CTL167724

Q-Logic 3 Expandable Capacity Joystick \$800
for Full Light Package
Part: CTL167725
Required with full light packages.

Q-Logic 3e Joystick \$700
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

2. Q-Logic 3 Expandable Capacity Required Components

Q-Logic 3 expandable capacity controls require the controller module and expandable harness below. Customer Service will include these items in a quote when expandable controls or specialty electronics are selected on this form.

Q-Logic 3 Expandable Capacity Controller Module..... \$695
Part: XRF000002 HCPCS: E2377

Q-Logic 3 Expandable Capacity Controller Harness..... \$425
Part: XRF000264. HCPCS: E2313

9.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307

Non-matching Joystick Armor \$50
Part: PTOINDV3422

Select color below

Match Joystick Armor to Base Color \$50
Part: PTOINDV3421

Part below will be selected with system configurator.

Red Corvette
Part: DGN124336-PRD1

Grey Street
Part: DGN124336-PSL3

Pretty in Pink
Part: DGN124336-PPK3

Blue Jean
Part: DGN124336-PBU11

Sugar Plum
Part: DGN124336-PPR7

Cherry Bomb
Part: DGN124336-PRD4

White Satin
Part: DGN124336-PWH2

Evolution Orange
Part: DGN124336-POR7

Lemon Crush
Part: DGN124336-PYL4

Military Digital Camo
Part: DGN124336-WCM7

Back in Black
Part: DGN124336-PBK1

Deep Purple
Part: DGN124336-PPR5

Midnight Blue
Part: DGN124336-PBU10

Green Machine
Part: DGN124336-PGN11

Faded Glory
Part: DGN124336-WFL2

⚡ Night Train*
* See Black, Soft Touch

Raspberry Beret
Part: DGN124336-PPK4

Electric Blue
Part: DGN124336-PBU9

Military Green
Part: DGN124336-PGN8

Carbon Fiber
Part: DGN124336-WCB1

10. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$85
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C) \$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$120
Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E) \$65
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F) \$65
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$85
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle \$99
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$99
Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308

**11. JOYSTICK MOUNTING BRACKETS**

Optional. Select a mount to be used with a joystick selected in previous section.

11.1a. Fixed, In-line Joystick Bracket


Fixed, In-line, Bracket Standard
Part: FRMASMB10617

11.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$280
Part: ACC1908434. HCPCS: E1028

Swing-away, NE/NE+ \$280
Part: ACC1908437. HCPCS: E1028


**11.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$350
Part: XRF000003. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Stealth Swing-Away/Flip Down Bracket \$294.25
Part: XRF000005. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

Multi-Axis Mount \$100
Part: XRF000004
Allows height adjustment of joystick.

**11.2. Joystick Mount Location**

☐ Left Side Mount
Part: PTOINDV1225

☐ Right Side Mount
Part: PTOINDV1226

12. TRU-BALANCE 3® SEATING - QUANTUM INTELLIGENT SYSTEMS™**12.1. TRU-Balance 3 Configuration**

REQUIRED. Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.

TRU-Balance 3 Static Seat \$2,295
Part: PTOINDV3442
20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Recline Only \$7,595
Part: PTOINDV3444. HCPCS: E1004
175° of power recline.

TRU-Balance 3 Power Tilt Only \$6,995
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Tilt & Recline \$11,895
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

12.2. 10" Power Lift

Optional. Select one option below.

TRU-Balance 3 Power 10" Lift \$2,600
Part: PTOINDV3775. HCPCS: E2300
10" of scissor lift. 22"W x 22"D maximum seat size.

12.3. Occupied Transit

Optional. Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU-Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.

Occupied Transit \$275
Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

12.4. Seat Dimensions - Width

REQUIRED with seat selection. 15" width is only available with non-matching back width. Minimum seat width is 12". Maximum seat width is 22".

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**12.5. Seat Dimensions - Width Package****REQUIRED with seat selection.****12"-16" Width Package**

Part: XRF000006

Not available with Power Lift.

14"-16" Width Package

Part: XRF000007

Only available with Power Lift.

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

Maximum width is 22".

12.6. Seat Dimensions - Depth**REQUIRED with seat selection.** Minimum seat depth is 12". Maximum seat depth is 22".**12" D**
(#1002)**13" D**
(#1003)**14" D**
(#1004)**15" D**
(#1005)**16" D**
(#1006)**17" D**
(#1007)**18" D**
(#1008)**19" D**
(#1009)**20" D**
(#1010)**21" D**
(#1011)**22" D**
(#1012)**12.7. Seat Dimensions - Depth Package****REQUIRED with seat selection.****12"-16" Depth Package**

Part: XRF000011

Not available with Power Lift.

14"-16" Depth Package

Part: XRF000012

Only available with Power Lift.

16"-20" Depth Package

Part: XRF000013

20"-24" Depth Package

Part: XRF000159

Maximum depth is 22".

12.8. Seat to Floor Height**REQUIRED.****Lowest Seat to Floor**

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559

Highest Seat to Floor

Part: PTOINDV3560

Low	Med	High
Static Seat/ Recline Only		
15"	16"	17"
Tilt/Tilt & Recline/Lift		
17.75"	18.75"	19.75"
TRU-Balance 3 Seat to Floor Heights		

12.9. Matching Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.


Stealth Simplicity G Cushion

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.

Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size

Reference part: XRF000111. MSRP: \$110/\$165


**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size

Reference part: XRF000115. MSRP: \$315/\$425

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size

Reference part: XRF000112. MSRP: \$460/\$530


**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size

Reference part: XRF000113. MSRP: \$375/\$430


**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608

TRU-Comfort 2 SPP, Matching Seat Size

Reference part: XRF000116. MSRP: \$400/\$500


**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

Glacial SP, Matching Seat Size

Reference part: XRF000114. MSRP: \$330/\$440


**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size

Reference part: XRF000122. MSRP: \$325/\$475

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-ESSSPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size

Reference part: XRF000118. MSRP: \$400/\$500

Essence SPP, Small Wedge, Matching Seat Size

Reference part: XRF000119. MSRP: \$435/\$535

Essence SPP, Medium Wedge, Matching Seat Size

Reference part: XRF000120. MSRP: \$465/\$565

Essence SPP, Large Wedge, Matching Seat Size

Reference part: XRF000121. MSRP: \$495/\$595


**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCGSPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size

Reference part: XRF000117. MSRP: \$460/\$530

**12.10. Back Type Selection****REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.**Redesigned, Static Type Back**

Part: XRF000014

12" and 24" back width will use the TRU-Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: XRF000015

Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$150
Part: XRF000016

12" and 24" back width will use the TRU-Balance 3 Back assembly. Not available with recline configurations. Only available with 2-post, flip back or pediatric armrests. Select a back cushion type in the subsequent section.

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**12.11. Back Cane Selection****REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable Back Cane No Charge Part: MEC140266	2" Endomorph (1" per side) Back Cane \$70 Part: XRF000269
Flip Forward Back Cane \$25 Part: XRF000273	4" Endomorph (2" per side) Back Cane \$70 Part: XRF000271

2. Back Cane Height

16"H Cane Part: FRMASMB15834	20"H Cane Part: FRMASMB15838	24"H Cane Part: FRMASMB15842
17"H Cane Part: FRMASMB15835	21"H Cane Part: FRMASMB15839	25"H Cane Part: FRMASMB15843

12.12. Back Cushion Configuration**REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions.** Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**

Sport Back No Charge Part: XRF000017 Complete following width package and size.	After-market Back Interface \$495 Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs. A Special Order is required on back widths less than 12".
Omit Back Pan/Cushion No Charge Part: PTOINDV1287 Select a following width package.	

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" Width \$675 Part: SET123622. HCPCS: E2620 Complete following width package and size.	TRU-Comfort Back (Static/Flip Forward), 22" Width \$775 Part: TC4V22W32W. HCPCS: E2621 Complete following width package and size.
TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$675 Part: TC4V10W21W. HCPCS: E2620 Complete following width package and size.	TRU-Comfort Back, Back Cane Option, 22" Width \$775 Part: TC4V22W32WBC. HCPCS: E2621 Complete following width package and size.
TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$675 Part: TC4V10W21WBC. HCPCS: E2620 Complete following width package and size.	

1c. ADI Back, AL Series, 4-Point, Active Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height \$610 Part: ALBR14LA4. HCPCS: E2613	Active Contour, 18" Wide x 21" or 22" Overall Height \$655 Part: ALBR18TA4. HCPCS: E2613
Active Contour, 14" Wide x 16" or 17" Overall Height \$625 Part: ALBR14MA4. HCPCS: E2613	Active Contour, 18" Wide x 23" or 24" Overall Height \$670 Part: ALBR18XA4. HCPCS: E2613
Active Contour, 16" Wide x 19" or 20" Overall Height \$640 Part: ALBR16SA4. HCPCS: E2613	Active Contour, 20" Wide x 19" or 20" Overall Height \$640 Part: ALBR20SA4. HCPCS: E2613
Active Contour, 16" Wide x 21" or 22" Overall Height \$655 Part: ALBR16TA4. HCPCS: E2613	Active Contour, 20" Wide x 21" or 22" Overall Height \$655 Part: ALBR20TA4. HCPCS: E2613
Active Contour, 16" Wide x 23" or 24" Overall Height \$670 Part: ALBR16XA4. HCPCS: E2613	Active Contour, 20" Wide x 23" or 24" Overall Height \$670 Part: ALBR20XA4. HCPCS: E2613
Active Contour, 18" Wide x 19" or 20" Overall Height \$640 Part: ALBR18SA4. HCPCS: E2613	

1d. ADI Back, AL Series, 4-Point, Deep Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height \$685 Part: ALBR14LD4. HCPCS: E2615	Deep Contour, 18" Wide x 19" or 20" Overall Height \$715 Part: ALBR18SD4. HCPCS: E2615
Deep Contour, 14" Wide x 16" or 17" Overall Height \$700 Part: ALBR14MD4. HCPCS: E2615	Deep Contour, 18" Wide x 21" or 22" Overall Height \$730 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 19" or 20" Overall Height \$715 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height \$715 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height \$730 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 21" or 22" Overall Height \$730 Part: ALBR20TD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height \$745 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height \$745 Part: ALBR20XD4. HCPCS: E2615

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**1e. ADI Back, AL Series, 4-Point, Extra Deep Contour**

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height \$805 Part: ALBR16T64. HCPCS: E2615	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$820 Part: ALBR18X64. HCPCS: E2615
Extra Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$820 Part: ALBR16X64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 19" or 20" Overall Height..... \$790 Part: ALBR20S64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$790 Part: ALBR18S64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 21" or 22" Overall Height..... \$805 Part: ALBR20T64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$805 Part: ALBR18T64. HCPCS: E2615	Extra Deep Contour, 20" Wide x 23" or 24" Overall Height..... \$820 Part: ALBR20X64. HCPCS: E2615

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$525 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$525 Part: QB2AL18AS. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$525 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$525 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$525 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$525 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$525 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$525 Part: QB2AL20AT. HCPCS: E2613

1g. Q-Back 2.0, AL Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$595 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$595 Part: QB2AL18DS. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$595 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$595 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$595 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$595 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$595 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$595 Part: QB2AL20DT. HCPCS: E2620

12.13. Lateral Positioning Wedge

Optional.

Lateral Positioning Wedge, pair \$30 Part: SET123465 Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

12.14. TRU-Balance 3 Back Shroud

REQUIRED. The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Blue Jean Part: XRF000031	Sugar Plum Part: XRF000035
Cherry Bomb Part: XRF000020	White Satin Part: XRF000024	Evolution Orange Part: XRF000028	Lemon Crush Part: XRF000032	
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Midnight Blue Part: XRF000029	Green Machine Part: XRF000033	
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Electric Blue Part: XRF000030	Military Green Part: XRF000034	

12.15. TRU-Balance 3 Power Tilt Side Shrouds

REQUIRED with TRU-Balance 3 Power Tilt selections. Reference colors shown on page 2.

Red Corvette Part: DGN1709093-PRD1	Grey Street Part: DGN1709093-PSL3	Pretty in Pink Part: DGN1709093-PPK3	Blue Jean Part: DGN1709093-PBU11	Sugar Plum Part: DGN1709093-PPR7
Cherry Bomb Part: DGN1709093-PRD4	White Satin Part: DGN1709093-PWT2	Evolution Orange Part: DGN1709093-POR7	Lemon Crush Part: DGN1709093-PYL4	
Back in Black Part: DGN1709093-PBK1	Deep Purple Part: DGN1709093-PPR5	Midnight Blue Part: DGN1709093-PBU10	Green Machine Part: DGN1709093-PGN11	
Night Train Part: DGN1709093-PBK3	Raspberry Beret Part: DGN1709093-PPK4	Electric Blue Part: DGN1709093-PBU9	Military Green Part: DGN1709093-PGN8	

TRU-Balance 3 seating section continued on following page.

12. TRU-BALANCE 3 SEATING CONTINUED**12.16. TRU-Balance 3 Power Lift Side Shrouds****REQUIRED** with TRU-Balance 3, 10" Power Lift selection. The power lift side shrouds are located at the base of the TRU-Balance 3 system and is available in matching shroud colors. Reference colors shown on page 2.

Red Corvette Part: DGN1800070	Grey Street Part: DGN1800071	Pretty in Pink Part: DGN1800068	Blue Jean Part: DGN1905060	Sugar Plum Part: DGN1905063
Cherry Bomb Part: DGN1905045	White Satin Part: DGN1800073	Evolution Orange Part: DGN1905051	Lemon Crush Part: DGN1905054	
Back in Black Part: DGN1800061	Deep Purple Part: DGN1800069	Midnight Blue Part: DGN1800063	Green Machine Part: DGN1905057	
Night Train Part: DGN1800062	Raspberry Beret Part: DGN1905048	Electric Blue Part: DGN1800064	Military Green Part: DGN1800066	

12.17. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections Part: PTOINDV3468	14"W-16"W Package Part: XRF000037	18"W-22"W Package Part: XRF000039
10"W-14"W Package Part: XRF000036	16"W-20"W Package Part: XRF000038	20"W-24"W Package Part: XRF000040

12.18. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	14" W (#3563)	15" W (#3856)	16" W (#3564)	17" W (#3565)	18" W (#3566)	19" W (#3567)	20" W (#3568)	21" W (#3569)	22" W (#3570)	24" W (#3571)
-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------	-------------------------

12.19. Overall Back Height**Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion. Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.**

14" Height Part: SETHGHT1027	18" Height Part: SETHGHT1031	22" Height Part: SETHGHT1035	26" Height Part: SETHGHT1039	30" Height Part: SETHGHT1043
15" Height Part: SETHGHT1028	19" Height Part: SETHGHT1032	23" Height Part: SETHGHT1036	27" Height Part: SETHGHT1040	31" Height Part: SETHGHT1102
16" Height Part: SETHGHT1029	20" Height Part: SETHGHT1033	24" Height Part: SETHGHT1037	28" Height Part: SETHGHT1041	32" Height Part: SETHGHT1103
17" Height Part: SETHGHT1030	21" Height Part: SETHGHT1034	25" Height Part: SETHGHT1038	29" Height Part: SETHGHT1042	

S Sport Back**T** TRU-Comfort Back**ST** Sport or TRU-Comfort Back**SQT** Sport, Q-Back, or TRU-Comfort Back

	Tilt or Static Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
	15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
	16"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	-	-	-	-	-	-	-	-	-		
	17"	S	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	-	-	-	-	-	-	-	-	-		
	18"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	-	-	-	-	-	-	-	-	-		
	19"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST	S	ST	ST	ST	ST	ST	ST	ST	ST	S		
	20"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST	S	ST	ST	ST	ST	ST	ST	ST	ST	S		
	21"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST	-	ST	ST	ST	ST	ST	ST	ST	ST	ST		
	22"	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST	ST	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST		
	23"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	ST	ST	SQT	ST	SQT	ST	SQT	ST	ST		
	24"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	S	-	ST	ST	SQT	ST	SQT	ST	SQT	ST		
	25"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	S	-	ST	ST	SQT	ST	SQT	ST	SQT	ST		
	26"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	ST	ST	ST	ST	ST	ST	ST	ST	ST		
	27"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	ST	ST	ST	ST	ST	ST	ST	ST	ST		
	28"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	ST	ST	ST	ST	ST	ST	ST	ST	S		
	29"	-	ST	ST	ST	ST	ST	ST	ST	ST	ST	ST	-	ST	ST	ST	ST	ST	ST	ST	ST	ST		
30"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	ST	ST	ST	ST			
31"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	ST	ST	ST	ST			
32"	-	-	-	-	-	-	-	-	-	-	-	-	ST	ST	ST	ST	ST	ST	ST	ST	ST			

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.20. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap, 9.75" \$35 Part: DWR1044E456B	Privacy Flap, 12.75" \$35 Part: DWR1044E456A
---	--

12.21. Back Height Extension**Optional. Not available with back cane option, 12" back width, or 24" back width.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension \$44.95 Part: XRF000041

13. TRU-BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

13.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

L R Flip-Back, 8"-16"H Armrests Standard
Left part: ARM1806662. Right part: ARM1806663

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

L R Flip-Back Armrest Downpost \$87.50 ea.
Left part: XRF000265. Right part: XRF000266

13.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

L R Full Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000047. Right part: XRF000048

L R Desk Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000049. Right part: XRF000050

1b. Two-Post Armrest for Static Seating

L R Full Length, Two-Post Armrest for Static Seating \$150 ea.
Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

L R Desk Length, Two-Post Armrest for Static Seating \$150 ea.
Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

13.1c. Heavy Duty, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold as pair)
Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$75 ea.
Heavy Duty (Sold as pair)
Part: XRF000056. HCPCS w/ Static Seat: E0973

13.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H
PTOINDV3535

10" H
PTOINDV3537

12" H
PTOINDV3539

14" H
PTOINDV3541

16" H
PTOINDV3543

9" H
PTOINDV3536

11" H
PTOINDV3538

13" H
PTOINDV3540

15" H
PTOINDV3542

2. Right Armrest Height

8" H
PTOINDV3450

10" H
PTOINDV3452

12" H
PTOINDV3454

14" H
PTOINDV3456

16" H
PTOINDV3458

9" H
PTOINDV3451

11" H
PTOINDV3453

13" H
PTOINDV3455

15" H
PTOINDV3457

13.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

L R 10" Length
Left part: XRF000057. Right part:
XRF000123

L R 14" Length
Left part: XRF000058. Right part:
XRF000124

L R 17" Length
Left part: XRF000059. Right part:
XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads.

14. TRU-BALANCE 3 ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

14.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection No Charge
Part: XRF000060
Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge
Part: XRF000308
Only available with Flip-back, 8"-16" armrests.

14.1b. Foam Armpads

L R Straight, Armpad Standard
Left part: XRF000061. Right part: XRF000126

L R Waterfall, Armpad No Charge
Left part: XRF000062. Right part: XRF000127
Not available with Flip-back, 8"-16" armrests.

TRU-Balance 3 armpads section continued on following page.

14. TRU-BALANCE 3 ARMPADS CONTINUED**14.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back ArmrestsFalcon Gel, Full Length Left and Right \$185
Part: ARM1806807Falcon Gel, Full Length Left, Desk Length Right \$162.50
Part: ARM1806813Falcon Gel, Desk Length Left, Full Length Right \$162.50
Part: ARM1806810Falcon Gel, Desk Length Left and Right \$140
Part: ARM1806806**b. Falcon Gel Armpads for Tube Style Armrests**Falcon Gel, Full Length Left and Right \$296.40
Part: ARM136438Falcon Gel, Full Length Left, Desk Length Right \$274.22
Part: ARM136442Falcon Gel, Desk Length Left, Full Length Right \$274.22
Part: ARM136439Falcon Gel, Desk Length Left and Right \$252.40
Part: ARM136437**15. TRU-BALANCE 3 ARMREST BRACKETS****15.1. Joystick Bracket****REQUIRED.**Standard Fixed Receiver Standard
Part: XRF000063Stealth Height Adjustable Receiver \$214
Part: XRF000064. HCPCS: K0108**16. TRU-BALANCE 3 HEADRESTS****Optional.** Select a headrest pad and bracket.**16.1a. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$167
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$167
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$167
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$167
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink \$294
Part: ST-UHB2D33200-Q. HCPCS: E1028**16.1b. Headrests, Comfort Plus****1. Headrest Pads**Stealth 8" Comfort Plus (CP450) \$151.94
Part: POS148484. HCPCS: E0955Stealth 10" Comfort Plus (CP250) \$151.94
Part: SETHEAD1057. HCPCS: E0955Stealth 14" Comfort Plus (CP150) \$151.94
Part: INDPART2327. HCPCS: E0955Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$151.94
Part: SETHEAD1069. HCPCS: E0955**2. Headrest Brackets**Multi-Axis Removable Headrest Bracket (HMO475P) \$160.50
Part: HDWASMB1084. HCPCS: E1028**17. LAP BELTS****REQUIRED.** Select a lap belt below.**17.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**17.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350) \$175
Part: ST-IFRB5350Stealth Positioning 60" Single Pull, Push Button, Medium \$86
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning 60" Single Pull, Push Button, Large \$86
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning 60" Single Pull, Push Button, Extra Large \$86
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable) \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

18. TRU-BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)
Part: XRF000108

Articulating Foot Platform (Requires 1 Actuator Control)
Part: XRF000110

Combined Legrests (Requires 2 actuator control)
Part: XRF000109

19. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

19.1a. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Power Articulating Foot Platform \$2,995
Part: PTOINDV3754. HCPCS: E1012+K0040

Manually Adjustable Center Mount Foot Platform... \$250
Part: XRF000155. HCPCS: K0040

**2. Lower Extension and Angle Adjustable Footplates**

8"-12" Extension, 5"x6" Footplates

Part: RIG152577

Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates

Part: RIG150206

10"-16" Extension, 5"x8" Footplates with Width Extension

Part: MEC1707469

10"-16" Extension, 5"x12" Footplates with Width Extension
Part: MEC1707465

13"-19" Extension, 5"x6" Footplates

Part: RIG151505

13"-19" Extension, 5"x8" Footplates with Width Extension

Part: MEC1707474

13"-19" Extension, 5"x12" Footplates with Width Extension
Part: MEC1707468

Center Mount Foot Platform Compatibility Matrix for TRU-Balance® 3	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
10"-16" Length Extension	!	✓	✓	✓	✓	✓	✓	✓	✓
13"-19" Length Extension	!	!	!	!	!	!	!	!	✓



No interference issues at this seat-to-floor height.



May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

Small Calf Pads, Manually Adjustable Center Mount Only, \$150
6"Wx6"H (pair)
Part: XRF000156

Large Calf Pads, Manually Adjustable Center Mount Only, \$150
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel for AFP, 18"Wx8"H \$150
Part: FRMASMB7110007

Small Calf Panel for Manually Adjustable, 18"Wx8"H \$150
Part: FRMASMB7110012

Large Calf Panel for AFP, 24"Wx8"H \$200
Part: RIG123824
Only available on widths of 20" or greater.

Large Calf Panel for Manually Adjustable, 24"Wx8"H \$200
Part: RIG123825
Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

4" Outside Edge (pair) \$50
Part: RIG152510
Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension \$50
Part: RIG152519

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

6" Outside Edge (pair) \$50
Part: RIG152516
Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG133717

Large Reversed Dartex Panel Cover (pair) \$150
Part: RIG134744

Large Reversed Dartex Calf Pad Cover (pair) \$150
Part: RIG133719

Footplate Connector \$15
Part: RIGASMB7120023

Small Reversed Dartex Panel Cover (pair) \$150
Part: RIG134673

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$205
Part: ST-FPG5008T-Q

5"x11" Tapered Gel Footplate Cover for Width Extension \$250
Part: ST-FPG7011T-Q

5"x8" Tapered Gel Footplate Covers for Width Extension \$225
Part: ST-FPG6508T-Q

Legrests section continued on following pages.

19. LEGRESTS CONTINUED**19.1b. 70° Power Articulating, Swing Away Legrests**

4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

P L R 70° Articulating, Swing-Away.....\$1,197.50 ea.
Legrest with 14.5"-18" Extension
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.
Right Part: FRMASMB9075. Pair HCPCS: E1010+K0040+E0951

P L R 70° Articulating, Swing-Away.....\$1,197.50 ea.
Legrest with 11"-13" Extension
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.
Right Part: FRMASMB11367. Pair HCPCS: E1010+K0040+E0951

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Angle Adjustable, 11"-13"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!	!	!	!	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 4"Wx6"D
Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.

P L R 5"Wx8"D
Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.

P L R 6"Wx11"D¹
Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.

P L R 5"Wx6"D
Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.

P L R 6"Wx8"D¹
Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.

P L R 8"Wx11"D²
Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest Accessories

Calf Pad (pair) 6"Wx6"H..... No Charge
Part: MEC153735

Falcon Gel Padded Knee Adductor Buttons (KBE)(pair)\$125
Part: INDPART2242. HCPCS: E0953

19.1c. 60°-90° Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

P L R 60°-90° Swing-Away Legrest, Power Seating..... No Charge
Pair part: RIG120062. Left: FRMASMB17558.
Right: FRMASMB17559
Available with Power Positioning only. Not available with 11"-13" extensions.

P L R 60°-90° Swing-Away Legrest, Static Seat \$81 ea.
Pair part: XRF000082. Left: XRF000083.
Right: XRF000084
Available with Static Seat only. Not available with 11"-13" extensions.

P L R 60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge
Pair part: RIG136481. Left: RIG132607.
Right: RIG132609.
Available with Power Positioning only. For use only with 11"-13" lower extensions.

P L R 60°-90° Swing-Away Legrest, Short, Static Seat \$81 ea.
Pair part: XRF000085. Left: XRF000086.
Right: XRF000087.
Available with Static Seat only. For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Upper

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

70° Swing-Away Legrest (pair).....\$162
Part: FRMASMB7841
Not available with 11"-13" extensions.

2. Swing-Away Legrests Lower Extension

P L R 11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.
N/A with Multi-axis.

P L R 13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.
N/A with Multi-Axis.

P L R 11"-16.5" Length
Pair part: XRF000091.
Left: XRF000092.
Right: XRF000093.
Multi-Axis only.

P L R 15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	!
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	!
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	!
11"-13" extension causes no interference at all seat-to-floor heights.									

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Extruded Footplate(s)

P L R 7.25"Wx6"D Extruded Footplate..... Standard
Part pair: FRMASMB7847. Left: FRMASMB7861
Right: FRMASMB7862

P L R 8.5"Wx6"D Extruded Footplate Standard
Part pair: FRMASMB7848. Left: FRMASMB7863
Right: FRMASMB7864
May cause interference on seat widths less than 18."

Legrests section continued on following pages.

19. LEGRESTS CONTINUED**3b. Angle Adjustable Footplate(s) and Heel Loops**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐	☐	☐	4"Wx6"D \$85 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040+E0951	☐	☐	☐	6"Wx8"D¹ \$85 ea. Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563 HCPCS: K0040+E0951
☐	☐	☐	5"Wx6"D \$85 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040+E0951	☐	☐	☐	6"Wx11"D¹ \$85 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040+E0951
☐	☐	☐	5"Wx8"D \$85 ea. Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 HCPCS: K0040+E0951	☐	☐	☐	8"Wx11"D² \$85 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040+E0951

3c. Multi-Axis Footplate(s) and Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐	☐	☐	4"Wx6"D \$247.50 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 HCPCS: K0040+E0951	☐	☐	☐	6"Wx8"D¹ \$247.50 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 HCPCS: K0040+E0951
☐	☐	☐	5"Wx6"D \$247.50 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 HCPCS: K0040+E0951	☐	☐	☐	6"Wx11"D¹ \$247.50 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 HCPCS: K0040+E0951
☐	☐	☐	5"Wx8"D \$247.50 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 HCPCS: K0040+E0951	☐	☐	☐	8"Wx11"D² \$247.50 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 HCPCS: K0040+E0951

4. 60°-90° Swing-Away Legrest Accessories

Gel Neoprene Footrest Sleeve (pair) \$80
Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
QTY Part: INDPART2823. HCPCS: K0108

19.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ **Elevating Legrest Assembly** \$162.50 ea.
Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013.
HCPCS: E0990
Select lower extension and footplate

2. Lower Extension

☐ ☐ ☐ **11"-14" Length**
Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
For reinforced only

☐ ☐ ☐ **8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length**
Pair part: XRF000101. Left: XRF000102. Right: XRF000103
For angle adjustable only

☐ ☐ ☐ **14"-17" Length**
Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
For reinforced only

☐ ☐ ☐ **13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length**
Pair part: XRF000104. Left: XRF000105. Right: XRF000106
For angle adjustable only

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ **7.25"Wx6"D Reinforced Footplate** Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ **8.5"Wx6"D Reinforced Footplate** Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18".

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". ³4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ **4"Wx6"D** \$85 ea.
Pair parts: FRMASMB8812 x2.
HCPCS: K0040

☐ ☐ ☐ **5"Wx6"D** \$85 ea.
Pair parts: FRMASMB8811 x2.
HCPCS: K0040

☐ ☐ ☐ **5"Wx8"D with Heel Loops** \$85 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx8"D¹ with Heel Loops** \$85 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

☐ ☐ ☐ **6"Wx11"D¹ with Heel Loops** \$85 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

☐ ☐ ☐ **8"Wx11"D² with Heel Loops** \$85 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

19. LEGRESTS CONTINUED**4. Elevating Legrest Accessories**

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$125
 Part: INDPART2242. HCPCS: E0953

19.2. Legrest Accessories

Heel Loops (pair) \$50 Wheel Bumper \$12.50 ea.
 Part: FRMASMB7873. HCPCS: E0951 Part: XRF000100
 Standard on angle adjustable and multi axis footplates when noted.

19.3. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

20. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options.

20.1a. Q-Logic 3 Controls through Joystick

Single Actuator Through Q-Logic 3 \$1,895 Multiple Actuator Through Q-Logic 3 \$2,895
 Controller and Interface (AM1) Controller and Interface (AAM)
 Part: XRF000067. HCPCS: E2310 Part: XRF000069. HCPCS: E2311

Single Actuator Through Q-Logic 3 \$2,895
 Controller and Interface (AAM)
 Part: XRF000068. HCPCS: E2310

Optional Switch Control

Not available with AM1 electronics. Switch will be mounted to same side of joystick.

4-Function Push Button 4-Way Toggle 4-Way Button Toggle
 Part: CTL123218 Part: ELEASMB5925 Part: ELEASMB5926
 See section 21.2 for switch mounting.

**20.1b. Q-Logic 3 Controls with iAccess Switch Control**

iAccess control is not available with combined or independent power legrests.

Single Actuator Control, Through Flat Button \$2,295 iAccess and Q-Logic 3 Controller and Interface (AM1) Part: XRF000070. HCPCS: E2310	Single Actuator Control, Through Toggle Button \$2,295 iAccess and Q-Logic 3 Controller and Interface (AM1) Part: XRF000073. HCPCS: E2310
Single Actuator Control, Through Flat Button \$3,295 iAccess and Q-Logic 3 Controller and Interface (AAM) Part: XRF000071. HCPCS: E2310	Single Actuator Control, Through Toggle Button \$3,295 iAccess and Q-Logic 3 Controller and Interface (AAM) Part: XRF000074 HCPCS: E2310
Multiple Actuator Control Through Flat Button \$3,295 iAccess and Q-Logic 3 Controller and Interface (AAM) Part: XRF000072. HCPCS: E2311	Multiple Actuator Control Through Toggle Button \$3,295 iAccess and Q-Logic 3 Controller and Interface (AAM) Part: XRF000075. HCPCS: E2311

20.2. 4-Function Push Button/iAccess Mount

REQUIRED with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU-Balance 3 HD armrests.

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$25 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$25 Part: MEC2001400
Under Arm Mount, Tube Style Armrest No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$25 Part: MEC2001401

1. Under Arm Mount Side Selection

Left Side Under Arm Mount Right Side Under Arm Mount
 Part: PTOINDV3678 Part: PTOINDV3679

21. TRU-BALANCE 3 ATTENDANT CONTROL

Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.


Attendant Control \$500 Attendant Control with Tilt, Recline, Alternative No Charge
 Part: XRF000078. HCPCS: E2331 Controls and Q-Logic Display or SCIM Module
 Select this option if no charge option criteria is not met. Part: XRF000079

22. SEAT MOUNTED LED LIGHTS**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$995

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

**23. PROVIDER INSTALLED ACCESSORIES****Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.TRU-Balance 3 Cane/Crutch Holder \$100
Part: ACC152515. HCPCS: E2207Quantum Backpack \$150
Part: ACC1905853Grip Lap Board (20005)..... \$38.95
Part: ACC152992Grip Drink Holder (20085) \$15.50
Part: ACC152993TRU-Balance 3 Hydration System (B) \$200
Part: ACC125010TRU-Balance 3 Transfer Handle (D) \$100 ea.
Part: MEC1809930. HCPCS: E0705
Not available with Heavy Duty armrests.TRU-Balance 3 Cup Holder..... No Charge
Part: ACC125003
For use with TRU-Balance 3 (track style) armrests.Cup Holder (Tube-style Armrest)..... No Charge
Part: ACCASMB2634
For use with tube-style armrests.TRU-Balance 3 Personal Item Hook \$50
Part: ACC125008
Not available with back cane option.TRU-Balance 3 Reflector Kit..... \$50
Part: ACC124670TRU-Balance 3 Backpack Holder (Pair) (A) \$50
Part: ACC147774
Not available with back cane option.TRU-Balance 3 Front Mounted \$50
Cell Phone/ Media Holder, X-Grip (E)
Part: ACC155835
Not available with Heavy Duty armrests.TRU-Balance 3 Rear View Mirror \$50
Part: ACC155833
Not available with Heavy Duty armrests.TRU-Balance 3 Oxygen Holder \$200
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOIN-DV3769.TRU-Balance 3 Clothing Guard (Pair) \$200
Part: ACC125007
Not available with Heavy Duty armrests.TRU-Balance 3 Glove Box (C) \$125
Part: ACC125011Stealth Large Tablet Holder \$375
Part: ACC141818
Not available with Heavy Duty armrests.Trilogy Vent Tray \$400
Part: ACC136791. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769. Contact custom for models other than 100 and 200.Compact Vent Tray \$1,395
Part: ACC133889. HCPCS: E1029Push Handles (pair)..... \$175
Part: ACC2005027
Not available with oxygen holder, back cane option, or vent tray.Accessory Bar \$100
Part: ACC1909634
Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.Large, Black, Weather Cover \$202
Part: ACCCOVR1021Q-Logic EX/NE Programmer \$545
Part: ELEASMB6976Flip-Back Armrest, Side Track Accessory Bracket ... \$10
Part: MEC1809178
Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.Privacy Bag..... \$39.95
Part: ACC1907570VOCSN Vent Tray for Power Tilt, Seat Mount \$895
Part: ACC2000162. HCPCS: E1029VOCSN Vent Tray for Power Lift, Seat Mount \$895
Part: ACC2000159. HCPCS: E1029VOCSN Vent Tray for Accessory Bar Mount..... \$895
Part: ACC2000166. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769.**24. BRANCH PATCHES****Optional.**Air Force Patch \$20
Part: ACC140240Army Patch \$20
Part: ACC140241Coast Guard Patch \$20
Part: ACC140239Marines Patch..... \$20
Part: ACC140243Navy Patch \$20
Part: ACC140244Purple Heart Patch \$20
Part: ACC140238

25. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**25.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$614
 Part: ST-ULS4042000-Q. HCPCS: E1028
 Only available with Redesigned back assemblies.

Swing-away, Rail Mounts, pair (TWBTLTC) \$472.50
 Part: INDPART2677. HCPCS: E1028
 Not available with Redesign back assemblies.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$210
 Part: INDPART2663. HCPCS: E0956

Flat Pads, 3"x5", pair (TWBLP3X5) \$210
 Part: INDPART2676. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$210
 Part: INDPART2664. HCPCS: E0956

Flat Pads, 4"x6", pair (TWBLP4X6) \$210
 Part: INDPART2675. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

25.2. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea.
 Part: INDPART2804. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$236.25 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Removable Bracket (LP430)..... \$262.50 ea.
 Part: INDPART2805. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

Removable Bracket with 1" Offset (LP430 & LHW-121)... \$400 ea.
 Part: INDPART2854. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$26.25
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) \$53.50
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.
 Part: INDPART2801. HCPCS: E0953/E0956

Long Pad, 4"x12" (TWBADD-4X12) \$75 ea.
 Part: INDPART2803. HCPCS: E0953/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.
 Part: INDPART2802. HCPCS: E0953/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

25.3. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware... \$245 ea.
 Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable & Removable Hardware... \$255 ea.
 Part: ST-UTTP7-Q. HCPCS: E1028

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$95 ea.
 Part: ST-UTP544056100-Q. HCPCS: E0956

UniLink, 4" x 12" Pad \$125 ea.
 Part: ST-UTP54D056100-Q. HCPCS: E0956

UniLink, 4" x 8" Pad..... \$105 ea.
 Part: ST-UTP548056100-Q. HCPCS: E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$75 ea.
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$75 ea.
 Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719



Positioning components section continued on following page.

25. POSITIONING COMPONENTS CONTINUED**25.4. Stealth Stealth Positioning Chest Harness**

Upper Classic, Medium (IFPS1103) \$142 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$142 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960

25.5. Stealth Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor \$90 ea.
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960

25.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

Facial Lateral for Comfort Series.....\$294.25 ea.
Right Side (CPS9-R-HW)
Part: INDPART2647. HCPCS: E1028

Facial Lateral for Comfort Series.....\$294.25 ea.
Left Side (CPS9-L-HW)
Part: INDPART2648. HCPCS: E1028

Facial Laterals for Comfort Series \$347.75 ea.
Left and Right Side (CPS9-B-HW)
Part: INDPART2649. HCPCS: E1028

2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A) \$52.50 ea.
Part: INDPART2815

Figure 8 Bracket with 13" Rod..... \$267.50 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.
Part: ST-900-Q. HCPCS: K0108

Part: ST-900-Q. HCPCS: K0108

Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.
Part: ST-920-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$250 ea.
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

25.7. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support..... \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$195
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$183.75 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket.....\$195
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea.
Part: POS124067. HCPCS: K0108

25.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$498.75
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$498.75
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

25.9. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder 3.75"Wx6.25"D \$266
(SPLE111) (pair)
Part: INDPART2783. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, 4"Wx8.25"D \$266
(SPLE112) (pair)
Part: INDPART2784. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, 5"Wx10.75"D \$266
(SPLE114) (pair)
Part: INDPART2785. HCPCS: E0951+E0952




Contoured and Classic Chest Harness



Chest Strap



Left and Right Facial Lateral with Pads



Headrest Accessories



Elbow Supports

26. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**26.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,884.93

Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard **Non-matching Shroud** \$50 **Match Shroud to Base Color** \$50
Part: DGN137728 Select color below Part: See color below**Red Corvette**
Part: DGN137729-PRD1**Grey Street**
Part: DGN137729-PSL3**Pretty in Pink**
Part: DGN137729-PPK3**Electric Blue**
Part: DGN137729-PBU9**Green Machine**
Part: DGN137729-PGN11**Back in Black**
Part: DGN137729-PBK1**White Satin**
Part: DGN137729-PWH2**Evolution Orange**
Part: DGN137729-POR7**Lemon Crush**
Part: DGN137729-PYL4**Military Green**
Part: DGN137729-PGN8**Night Train**
Part: DGN137729-PBK3**Deep Purple**
Part: DGN137729-PPR5**Midnight Blue**
Part: DGN137729-PBU10

Enhanced Display

26.2. Q-Logic SCIM Module**SCIM Module Kit for Q-Logic 3**.....\$525

Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 Expandable Capacity Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

26.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ**.....\$800

Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

26.4. Head Array**Adult, i-Drive Tri-Array Series, 3-Switch** \$4,414
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
(IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.




Adult Pro-Series

26.5. Sip and Puff Drive**1. Module****Sip and Puff Module for Q-Logic 3**.....\$1,895

Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.



Sip & Puff Module

2. Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series,**\$294.25
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive,\$294.25
Head Array (SP625-U)

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series,\$294.25
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)**.....\$32.10
Part: ST-SP625-S12-Q**Sip & Puff Bus Cable**.....\$350
Part: DWR1111H046**26.6. Sip and Puff Drive****Red Egg Switch (ESRED)**..... \$80.85
Part: SWTMCRO1077**Blue Egg Switch (ESBLU)**..... \$80.85
Part: SWTMCRO1079**Black Egg Switch (ESBLK)**..... \$80.85
Part: SWTMCRO1081**Green Egg Switch (ESGRN)**..... \$80.85
Part: SWTMCRO1078**Yellow Egg Switch (ESYEL)**..... \$80.85
Part: SWTMCRO1080

Egg Switch

Alternative controls and electronics section continued on following page.

26. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

26.7. Switch Adapter

2 to 1 Switch Adapter \$13.39

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

26.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,400

Part: CTL168057

27. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

Secondary Multiplier (SM) \$221.21

Part: XRF000107

Arm Mounts		1 Arm			
		0 Back	1 Back	2 Back	3+ Back
Arm Mounts	Joystick (Hand Control)	NA	SM	SM	SOF
	Enhanced Display	NA	SM	SM	SOF
	iAccess	NA	SM	SM	SOF
	USB Charger	NA	SM	SM	SOF
Back Mounts	SAJ	NA	SM	SM	SOF
	Attendant Control	NA	SM	SM	SOF
	SCIM	NA	SM	SM	SOF
	Sip n Puff	NA	SM	SM	SOF
Back Mounts	ECU	NA	SM	SM	SOF
	ECU	NA	SM	SM	SOF
	ECU	NA	SM	SM	SOF
	ECU	NA	SM	SM	SOF

Power Positioning Chart

Back Mounts		1 Arm			
		0 Back	1 Back	2 Back	3+ Back
Back Mounts	0 Back	NA	SM	SOF	SOF
	1 Back	NA	SM	SOF	SOF
	2 Back	SM	SM	SOF	SOF
	3+ Back	SOF	SOF	SOF	SOF

Static Seat Chart

Q-Logic 3 Multiplier Charts

