



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with Pediatric TRU Balance® 3 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 2.0 3S-SS\$8,533

Part: Q6 EDGE 2.0 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 2.0 3SP-SS\$8,234

Part: Q6 EDGE 2.0 3SP-SS. HCPCS: K0856

Single power option.

☐ Quantum Q6 Edge 2.0 3MP-SS\$8,055

Part: Q6 EDGE 2.0 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. DRIVE MOTORS

REQUIRED.

- ☐ 6 mph Drive Motors Standard
Part: MOT147433
- ☐ Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN153201-PRD1 |  Night Train
Part: DGN153201-PBK3 |  Deep Purple
Part: DGN153201-PPR5 |  Lemon Crush
Part: DGN153201-PYL4 |
|  Electric Blue
Part: DGN153201-PBU9 |  Pretty in Pink
Part: DGN153201-PPK3 |  Grey Street
Part: DGN153201-PSL3 |  Green Machine
Part: DGN153201-PGN11 |
|  Back in Black
Part: DGN153201-PBK1 |  Evolution Orange
Part: DGN153201-POR7 |  Midnight Blue
Part: DGN153201-PBU10 |  Military Green
Part: DGN153201-PGN8 |

4. BASE MOUNTED LIGHTS

Optional. Select only one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

4.1a. Fender Lights


Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights (without iLevel 12" Lift System) \$263
Part: XRF000177

Switch Fender Lights, Setup for iLevel 12" Lift System ... No Charge
Part: DGN154117

4.1b. Base Mounted Light

Base Mounted, LED, Full Light Package \$1045
Part: XRF000129,
Requires Q-Logic 3 Joystick for Full Light Package.



5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$363
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



6. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

- | | |
|---|---|
| 3"x14" Flat Free Black Tires (No Accent Color) Standard
Part: PTOINDV3525 (pair) | 3"x14" Pneumatic Black Tires (No Accent Color) No Charge
Part: WHL150034 x2 |
| 3"x14" Flat Free Black Tires for Wheel Accents No Charge
Part: WHL1611297 x2 | 3"x14" Pneumatic Black Tires for Wheel Accents No Charge
Part: WHL1611296 x2 |
| 4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368 | 4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367 |

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

- | | | |
|--|---|---|
| Back in Black Accent Color \$63
Part: XRF000179 | Deep Purple Accent Color \$63
Part: XRF000184 | Electric Blue Accent Color \$63
Part: XRF000188 |
| Red Corvette Accent Color \$63
Part: XRF000180 | Pretty in Pink Accent Color \$63
Part: XRF000185 | Lemon Crush Accent Color \$63
Part: XRF000189 |
| Night Train Accent Color \$63
Part: XRF000181 | Evolution Orange Accent Color \$63
Part: XRF000186 | Green Machine Accent Color \$63
Part: XRF000190 |
| Grey Street Accent Color \$63
Part: XRF000182 | Midnight Blue Accent Color \$63
Part: XRF000187 | Military Green Accent Color \$63
Part: XRF000191 |

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor Battery..... \$355 ea.
Part: BAT1704613. HCPCS: E2361NF-22, Gel Cell, MK Battery \$432 ea.
Part: BATGEL1008. HCPCS: E2361NF-22, AGM, MK Battery..... \$432 ea.
Part: BAT121346. HCPCS: E2361**8. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**8.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key JoystickStandard

Part: CTL2205379

Controls up to 2 actuators.

QLNES 10 Key Joystick\$420

Part: CTL2205380

Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**8.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 Controller Module\$730

Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$458

Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick\$840

Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights\$840

Part: CTL167725

Q-Logic 3e Joystick\$735

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic

Part: XRF000010

**8.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60

Part: ELE2103635

8.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.

Black, Soft Touch Standard

Part: DGN124307

Also matches Night Train shrouds.

Deep Purple.....\$63

Part: DGN124336-PPR5

Electric Blue\$63

Part: DGN124336-PBU9

Red Corvette.....\$63

Part: DGN124336-PRD1

Pretty in Pink\$63

Part: DGN124336-PPK3

Lemon Crush.....\$63

Part: DGN124336-PYL4

Back in Black.....\$63

Part: DGN124336-PBK1

Evolution Orange.....\$63

Part: DGN124336-POR7

Green Machine\$63

Part: DGN124336-PGN11

Grey Street.....\$63

Part: DGN124336-PSL3

Midnight Blue\$63

Part: DGN124336-PBU10

Military Green\$63

Part: DGN124336-PGN8

9. JOYSTICK HANDLES**Optional.** The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.☐ Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323☐ Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323☐ Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)☐ Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323☐ T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323☐ Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323☐ Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308☐ Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323

10. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**10.1a. Fixed, In-line Joystick Bracket**
☐ Fixed, In-line, Bracket for Q-Logic 3 Standard
 Part: FRMASMB10617
10.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.



 Swing-away, Q-Logic 3 \$319
 Part: ACC1908434. HCPCS: E1032
10.1c. Basic Swing Away Joystick Mount
 Basic Swing Away, Q-Logic, Left Side \$245
 Part: ACC2201041. HCPCS: E1032

 Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201042. HCPCS: E1032
10.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

 Retractable Mount \$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

 Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.
10.2. Joystick Mount Location
 Left Side Mount
 Part: PTOINDV1225

 Right Side Mount
 Part: PTOINDV1226
**11. TRU BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™****REQUIRED.** Maximum seat size is 22"W x 22"D.**11.1. TRU Balance 3 Configuration**
 TRU Balance 3 Static Seat \$2,410
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

 TRU Balance 3 Power Recline Only \$8,243
 Part: PTOINDV3444. HCPCS: E1004
 175° of power recline.

 TRU Balance 3 Power Tilt Only \$8,202
 Part: PTOINDV3443. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

 TRU Balance 3 Power Tilt & Recline \$12,615
 Part: PTOINDV3446. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.
11.2. 12" Power Lift**Optional.** Select one option below. ^Please verify the appropriate billing code to be used for this item with all third-party payors.
 TRU Balance 3, 12" Power Lift with iLevel® \$3,725
 Part: PTOINDV3773. HCPCS: E2298/E2300^
 12" of scissor lift with iLevel® System. Not available with high speed motors or manual recline.
**11.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Only available with 16"-20" width and depth packages. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>
 WC-19 Occupied Transit Option \$289
 Part: KIT152047
11.4. Width**REQUIRED with seat selection.** Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)
------------------	------------------	------------------	------------------	------------------

11.5. Width Package**REQUIRED with seat selection.**
 12"-16" Width Package
 Part: XRF000006.

 14"-16" Width Package
 Part: XRF000007

 16"-20" Width Package
 Part: XRF000008
11.6. Depth**REQUIRED with seat selection.** Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**11.7. Depth Package****REQUIRED with seat selection.****12"-16" Depth Package**

Part: XRF000011

14"-16" Depth Package

Part: XRF000012

16"-20" Depth Package

Part: XRF000013

11.8. Seat to Floor Height**REQUIRED.** Select a seat to floor height compatible with the seat functions previously selected.**16.5" Seat to Floor**

Part: PTOINDV3459

Not available with Power Tilt or
Power Lift Configurations.**17.5" Seat to Floor**

Part: PTOINDV3460

18.5" Seat to Floor


Part: PTOINDV3463

19.5" Seat to Floor

Part: PTOINDV3505

Not available with Static Seat or
Power Recline only.**11.9. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39

Part: POP1909006-S

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Large Wedge Positioning Kit..... \$105

Part: POP1909006-L

14"W and smaller

Small Wedge Positioning Kit\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit..... \$105

Part: POP1712340-L

15"W and larger

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**3. Cushion Top Foam**

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390

Additional Stealth Polyurethane Liner
Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392

Additional Stealth Tek Cover
Reference Part: XRF000394

Additional Mesh Cover
Reference Part: XRF000393

11.10. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Available with Static Seat and Power Tilt (1 actuator). Only available with 2-post, flip back, or pediatric armrests. Select a back cushion type in the subsequent section.

11.11. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane..... No Charge
Part: MEC140266

Flip Forward Back Cane \$53
Part: XRF000273

2" Endomorph (1" per side) Back Cane..... \$74
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... \$74
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

11.12. Back Cushion Configuration

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back..... No Charge
Part: XRF000017
Complete following width package and size.

Omit Back Pan/Cushion..... No Charge
Part: PTOINDV1287
Select a following width package.

After-market Back Interface \$535
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armpads and arm through. Does not include brackets to mount after-market backs.

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), 12" Width \$744
 Part: SET123622. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$744
 Part: TC4V10W21WBC. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$855
 Part: TC4V22W32WBC. HCPCS: E2621
 Complete following width package and size.

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
 Part: ALBR16SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622
 Part: ALBR16TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622
 Part: ALBR16XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622
 Part: ALBR18SA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$622
 Part: ALBR18TA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$622
 Part: ALBR18XA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622
 Part: ALBR20SA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622
 Part: ALBR20TA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622
 Part: ALBR20XA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height \$735
 Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height \$735
 Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height \$735
 Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height \$735
 Part: ALBR18SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height \$735
 Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
 Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height \$822
 Part: ALBR18S64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height \$822
 Part: ALBR18T64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height \$822
 Part: ALBR18X64. HCPCS: E2615

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
 Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AL Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1h. Q-Back 2.0, AL Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$703 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL18ET. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620

1i. Q-Back 2.0, AFT Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

1j. Q-Back 2.0, AFT Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT18DS. HCPCS: E2620	

1k. Q-Back 2.0, AFT Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AFT16EXX. HCPCS: E2620	
---	--

11.13. TRU Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Evolution Orange Part: XRF000028	Lemon Crush Part: XRF000032
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Midnight Blue Part: XRF000029	Green Machine Part: XRF000033
Night Train Part: XRF000022	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Military Green Part: XRF000034

11.14. Lateral Positioning Wedge**Optional.**

Lateral Positioning Wedge, pair \$32 Part: SET123465	
---	--

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**11.15. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections Part: PTOINDV3468	14"W-16"W Package Part: XRF000037
10"W-14"W Package Part: XRF000036	16"W-20"W Package Part: XRF000038
	18"W-22"W Package Part: XRF000039

11.16. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

12" W (#3562)	17" W (#3565)
14" W (#3563)	18" W (#3566)
16" W (#3564)	

	Tilt or Static Width						Recline Width					
	12"	14"	15"	16"	17"	18"	12"	14"	15"	16"	17"	18"
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-
	15"	T	-	-	-	-	-	-	-	-	-	-
	16"	S	S	T	S	T	S	T	S	T	-	-
	17"	S	S	T	S	T	S	T	S	T	-	-
	18"	-	S	T	S	T	S	T	S	T	-	-
	19"	-	S	QT	S	T	S	QT	S	T	S	QT
	20"	-	S	QT	S	T	S	QT	S	T	S	QT
	21"	-	S	QT	S	T	S	QT	S	T	S	QT
	22"	-	S	QT	S	T	S	QT	S	T	S	QT
	23"	-	S	T	S	T	S	QT	S	T	S	QT
	24"	-	S	T	S	T	S	QT	S	T	S	QT
	25"	-	S	T	S	T	S	QT	S	T	S	QT
	26"	-	S	T	S	T	S	QT	S	T	S	QT
	27"	-	S	T	S	T	S	QT	S	T	S	QT
	28"	-	S	T	S	T	S	QT	S	T	S	QT
	29"	-	S	T	S	T	S	QT	S	T	S	QT

S Sport Back
 Q Q-Back
 T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

11.17. Overall Back Height

Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.

Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	18" Height Part: SETHGHT1031	22" Height Part: SETHGHT1035	26" Height Part: SETHGHT1039
15" Height Part: SETHGHT1028	19" Height Part: SETHGHT1032	23" Height Part: SETHGHT1036	27" Height Part: SETHGHT1040
16" Height Part: SETHGHT1029	20" Height Part: SETHGHT1033	24" Height Part: SETHGHT1037	28" Height Part: SETHGHT1041
17" Height Part: SETHGHT1030	21" Height Part: SETHGHT1034	25" Height Part: SETHGHT1038	29" Height Part: SETHGHT1042

11.18. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap	\$48
Part: XRF000351	

11.19. Back Height Extension**Optional.** Not available with back cane option or 12" back width. Adds 2" of seat back pan height.

Back Height Extension	\$48
Part: XRF000041	

12. TRU BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**12.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

1. Flip-Back Armrest Selection

Flip-Back, 8"-16"H Armrests Standard Left part: ARM1806662. Right part: ARM1806663 Not available with TRU Balance 3 Static Seat.	Flip-Back, 8"-16"H Armrests for Static Seat \$289 Left part: XRF000320. Right part: XRF000321
--	--

2. Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

Flip-Back Armrest Downpost \$92 ea. Left part: XRF000265. Right part: XRF000266
--

3. Flip-Back, Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

TRU Balance 3 armrests section continued on following page.

12. TRU BALANCE 3 ARMRESTS CONTINUED**4. Flip-Back, Right Armrest Height**

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

12.1b. Pediatric, Removable, Height Adjustable ArmrestsHeight adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.****1. Pediatric Armrest Selection**

L R Pediatric, Removable Armrests \$79 ea.
 Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

2. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

3. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

12.2. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

L R 10" Length
 Left part: XRF000057. Right part: XRF000123

L R 14" Length
 Left part: XRF000058. Right part: XRF000124

13. TRU BALANCE 3 ARMPADS**REQUIRED.** Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.**13.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge
 Part: XRF000060
 Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge
 Part: XRF000308
 Only available with Flip-back, 8"-16" armrests.

13.1b. Foam Armpads

L R Straight, Armpad Standard
 Left part: XRF000061. Right part: XRF000126

L R Waterfall, Armpad No Charge
 Left part: XRF000062. Right part: XRF000127
 Not available with Flip-back, 8"-16" armrests.

13.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195
 Part: ARM1806807

Falcon Gel, Full Length Left, Desk Length Right.....\$171
 Part: ARM1806813

Falcon Gel, Desk Length Left, Full Length Right.....\$171
 Part: ARM1806810

Falcon Gel, Desk Length Left and Right\$147
 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22
 Part: ARM136438

Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
 Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
 Part: ARM136439

Falcon Gel, Desk Length Left and Right \$265.02
 Part: ARM136437

14. TRU BALANCE 3 ARMREST BRACKETS**14.1. Joystick Bracket****REQUIRED.**

Standard Fixed ReceiverStandard
 Part: XRF000063

Stealth Height Adjustable Receiver.....\$236
 Part: XRF000064. HCPCS: K0108



15. TRU BALANCE 3 HEADRESTS**Optional.** Select a headrest pad and bracket.**15.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Bracket**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**15.2. Q Cam Lock****Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

16. LAP BELTS**REQUIRED.** Select a lap belt below.**16.1a. Standard Lap Belts**38" Lap Belt 50" Lap Belt 60" Lap Belt
Part: ACCBELT1006 Part: ACCBELT1000 Part: ACCBELT1002**16.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350).....\$181
Part: ST-IFRB5350Stealth Positioning 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-QStealth Positioning 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-QStealth Positioning 60" Single Pull, Push Button, Extra Large....\$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q60" Sure-Lok Lap Belt (Non-retractable).....\$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016**17. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**17.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit base mounted foot platform No Charge
Part: PTOINDV3674Bell Shaped, Base Mounted Foot Platform \$247
Part: RIGASMB7120041**17.1b. Basic Center Mount Foot Platform**

Select legrest length in following section.

1. Center Mount Foot Platform TypeBasic Center Mount..... No Charge
Part: XRF000323**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**Short Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Length Extension.**4. Basic Center Mount Additional Accessories**Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

Legrests section continued on following pages.

17. LEGRESTS CONTINUED**17.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform TypePower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311. HCPCS: K0040**2. Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577Medium Extension, 5"x11" Footplates
Part: RIG2404018Pediatric Extension, 5"x8" Footplates
Part: RIG2404015Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Quick Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Quick Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**Small Calf Panel 18"Wx8"H \$184
Part: FRMASMB7110012**4. Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" footplates.

5. Additional AccessoriesSmall Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-QFootplate Riser \$75
Part: XRF000348Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-QFootplate Riser Calf Panel \$50
Part: FRMASMB7110011Footplate Connector \$16
Part: RIGASMB7120023

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108**17.1d. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ P	☐ L	☐ R	Elevating Legrest Assembly \$171 ea. Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990 Select lower extension and footplates
----------------------------	----------------------------	----------------------------	--

2a. Lower Extension

☐ P	☐ L	☐ R	11"-14" Length Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030 For reinforced only
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	14"-17" Length Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035 For reinforced only
----------------------------	----------------------------	----------------------------	---

☐ P	☐ L	☐ R	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103 For angle adjustable only
----------------------------	----------------------------	----------------------------	--

☐ P	☐ L	☐ R	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106 For angle adjustable only
----------------------------	----------------------------	----------------------------	---

Legrests section continued on following pages.

17. LEGRESTS CONTINUED

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12.5"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342
Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R 7.25"Wx6"D Reinforced Footplate..... Standard
Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R 8.5"Wx6"D Reinforced Footplate Standard
Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D with Heel Loops..... \$90 ea.
Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
HCPCS: K0040+E0951

P L R 6"Wx11"D¹ with Heel Loops..... \$90 ea.
Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
HCPCS: K0040+E0951

P L R 6"Wx8"D¹ with Heel Loops..... \$90 ea.
Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
HCPCS: K0040+E0951

P L R 8"Wx11"D² with Heel Loops..... \$90 ea.
Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
Part: INDPART2242. HCPCS: E0953

17.2. Legrest Accessories

Heel Loops (pair)..... \$53
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$19 ea.
Part: XRF000100

17.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options. Select one option from section 19.1 or 19.3.**18.1a. Tilt Only through Toggle****Only available with QLNES Joystick.**

Single Actuator Control, Power Tilt Only, Through Toggle No Charge

Part: ELE1905907

Select toggle side. Will require extended lead time when selected with Fender Lights.

1. Toggle Mount Side SelectionLeft Side Toggle
Part: PTOINDV1839Right Side Toggle
Part: PTOINDV1838**18.1b. Controls through Joystick**Single Actuator Through Controller and \$2,163
Interface

Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
Interface

Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and \$3,040
Interface (AAM)

Part: XRF000068. HCPCS: E2310

1. Optional Switch Control**Only available with Q-Logic 3 electronics.** Not available with single actuator through controller and interface.4-Function Push Button
Part: CTL123218

See section 20.2 for switch mounting.

4-Way Toggle
Part: ELEASMB59254-Way Button Toggle
Part: ELEASMB5926**2. User Controls Side Selection****Required with 4-Way Toggle and 4-Way Button Toggle.**Mount User Controls on Left
Part: PTOINDV3996Mount User Controls on Right
Part: PTOINDV3997**18.1c. Q-Logic 3 Controls with iAccess Switch Control****Only available with Q-Logic 3 electronics.** iAccess control is not available with combined or independent power legrests.Single Actuator Control Through Flat Button \$2,547
iAccess and Controller and Interface

Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000074. HCPCS: E2310

Single Actuator Control Through Toggle Button \$2,547
iAccess and Controller and Interface

Part: XRF000338. HCPCS: E2310

Multiple Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface

Part: XRF000339. HCPCS: E2311

Single Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000070. HCPCS: E2310

Multiple Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface

Part: XRF000340. HCPCS: E2311

18.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-function push button and iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.Joystick Mount No Charge
Part: MEC1907500Attendant Mount \$26
Part: MEC1906564Under Arm Mount, Flip-Back Armrest No Charge
Part: MEC1906531Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26
Part: MEC2001392Under Arm Mount, Tube Style Armrest No Charge
Part: MEC1906530Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26
Part: MEC2001391**1. Under Arm Mount Side Selection**Left Side Under Arm Mount
Part: PTOINDV3678Right Side Under Arm Mount
Part: PTOINDV3679**19. TRU BALANCE 3 ATTENDANT CONTROL****Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.Attendant Control \$525
Part: XRF000078. HCPCS: E2331


Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module

Part: XRF000079


20. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.



21. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.TRU Balance 3 Cane/Crutch Holder \$105
Part: ACC152515. HCPCS: E2207.Quantum Backpack \$168
Part: ACC1905853Grip Lap Board (20005) \$38.95
Part: ACC152992Grip Drink Holder (20085) \$15.50
Part: ACC152993TRU Balance 3 Hydration System (B) \$237
Part: ACC125010TRU Balance 3 Transfer Handle (D) \$105 ea.
Part: MEC2205680. HCPCS: E0705
Not available with Heavy Duty armrests.TRU Balance 3 Cup Holder No Charge
Part: ACC125003
For use with TRU Balance 3 (track style) armrests.Cup Holder (Tube-style Armrest) No Charge
Part: ACCASMB2634
For use with tube-style armrests.TRU Balance 3 Personal Item Hook \$69
Part: ACC125008
Not available with back cane option.TRU Balance 3 Reflector Kit \$53
Part: ACC124670TRU Balance 3 Backpack Holder (Pair) (A) \$69
Part: ACC147774
Not available with back cane option.TRU Balance 3 Front Mounted \$69
Cell Phone/ Media Holder, X-Grip (E)
Part: ACC155835
Not available with Heavy Duty armrests.TRU Balance 3 Rear View Mirror \$69
Part: ACC155833
Not available with Heavy Duty armrests.TRU Balance 3 Oxygen Holder \$210
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOIN-DV3769.TRU Balance 3 Clothing Guard (Pair) \$237
Part: ACC125007

Not available with Heavy Duty armrests.

TRU Balance 3 Glove Box (C) \$158
Part: ACC125011Stealth Large Tablet Holder \$625.28
Part: ACC141818
Not available with Heavy Duty armrests.Compact Vent Tray \$1,465
Part: ACC133889. HCPCS: E1029Push Handles (pair) \$195
Part: ACC2005027
Not available with oxygen holder, back cane option, or vent tray. Required with attendant manual recline.Accessory Bar \$132
Part: XRF000128
Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.Large, Black, Weather Cover \$202
Part: ACCCOVR1021Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Privacy Bag \$42
Part: ACC20005723Flip-Back Armrest, Side Track Accessory Bracket ... \$11
Part: MEC1809178
Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nutsPersonal Item Pouch \$84
Part: ACC2101017Depth Adjustable Vent Tray \$595
Part: ACC2302500. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769.Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496**22. BRANCH PATCHES****Optional.**Air Force Patch \$21
Part: ACC140240Army Patch \$21
Part: ACC140241Coast Guard Patch \$21
Part: ACC140239Marines Patch \$21
Part: ACC140243Navy Patch \$21
Part: ACC140244Purple Heart Patch \$21
Part: ACC140238

23. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#).**23.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664

Part: ST-ULS4042000-Q. HCPCS: E1034

Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664

Part: ST-ULPK404108B00-Q. HCPCS: E1034

Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664

Part: ST-ULPK404104B00-Q. HCPCS: E1034

Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332

Part: ST-ULPK404103L00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332

Part: ST-ULPK404103R00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250

Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250

Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals

Part: PTOINDV3499

Do Not Install

Part: PTOINDV3500

23.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.

Removable Mounting Hardware

Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.

Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.

Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left \$332 ea.

Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.

Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)

Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact \$281 ea.

Feature (TWBADD-R & LPHW407)

Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad \$106 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$106 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.

Part: ST-UIL001150A000-Q

UniLink, 2" Link \$81 ea.

Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left

Part: PTOINDV3585

Install on Right

Part: PTOINDV3586

Do Not Install

Part: PTOINDV3719

23.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.

Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.

Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$154 ea.

Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.

Part: ST-IFPS2104-Q. HCPCS: E0960



Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.5. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor\$97 ea.
(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor\$97 ea.
(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960

**23.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318
Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$492
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$492
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$434
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028


**2. Facial Lateral Pads and Accessories**

Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815

Facial Pad, 4"x2", Cool Core (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea.
Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

**23.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket.....\$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$57 ea.
Part: POS124067. HCPCS: K0108

23.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

23.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

24.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SM	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**25.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Night Train\$63 Part: DGN137729-PBK3	Pretty in Pink\$63 Part: DGN137729-PPK3	Lemon Crush\$63 Part: DGN137729-PYL4
Red Corvette\$63 Part: DGN137729-PRD1	Grey Street\$63 Part: DGN137729-PSL3	Midnight Blue\$63 Part: DGN137729-PBU10	Green Machine\$63 Part: DGN137729-PGN11
Back in Black\$63 Part: DGN137729-PBK1	Deep Purple\$63 Part: DGN137729-PPR5	Electric Blue\$63 Part: DGN137729-PBU9	Military Green\$63 Part: DGN137729-PGN8

**25.2. SCIM Module****Optional.**☐ **SCIM Module for Q-Logic 3**.....\$525
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

**25.3. Q-Logic 3 SAJ****Q-Logic 3 EX SAJ**.....\$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.

**25.4. Head Array****i-Drive 4.0 Bluetooth Tri-array Series**,.....\$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330**Adult, Curved, i-Drive Tri-Array with Mounting Hardware**.....\$4,474
Part: ST-IDHBT400-U-Q. HCPCS: E2330**i-Drive 4.0 Bluetooth Pro Series**,.....\$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330**Junior, Curved, i-Drive Tri-Array with Mounting Hardware**.....\$4,474
Part: ST-IDHBT450-U-Q. HCPCS: E2330**25.5. Sip and Puff Drive****1. Sip and Puff Module****Sip and Puff Module for Q-Logic 3**.....\$1,990
Part: CTL168056. HCPCS: E2325
Includes module and mounting hardware.**2. Sip and Puff Straw and Mount****The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series**,.....\$318
Left (SP625M-CL)
Part: ST-SP625M-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.**Stealth Sip and Puff for Stealth i-Drive**,.....\$318
Head Array (SP625-U)
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.**Stealth Sip and Puff for Comfort Series**,.....\$318
Right (SP625M-CR)
Part: ST-SP625M-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.**3. Sip & Puff Accessories****Sip & Puff Straws, 12 pack (SP625-S12)**.....\$35
Part: ST-SP625-S12-Q**25.6. Egg Switches****Red Egg Switch (ESRED)**.....\$87
Part: SWTMCRO1077**Green Egg Switch (ESGRN)**.....\$87
Part: SWTMCRO1078**Blue Egg Switch (ESBLU)**.....\$87
Part: SWTMCRO1079**Yellow Egg Switch (ESYEL)**.....\$87
Part: SWTMCRO1080**Black Egg Switch (ESBLK)**.....\$87
Part: SWTMCRO1081

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

25.7. Switch Adapter

2 to 1 Switch Adapter \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
 Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

