



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with Synergy Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 2.0 3S-SS \$8,533
Part: Q6 EDGE 2.0 3S-SS. HCPCS: K0848
Solid seat pan standard with no power option.

Quantum Q6 Edge 2.0 3SP-SS \$8,234
Part: Q6 EDGE 2.0 3SP-SS. HCPCS: K0856
Single power option for Power Elevating Seat.

2. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard
Part: MOT147433

Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

Account Number: _____ Date: _____ Quote Name: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

 Red Corvette Part: DGN153201-PRD1	 Night Train Part: DGN153201-PBK3	 Deep Purple Part: DGN153201-PPR5	 Lemon Crush Part: DGN153201-PYL4
 Electric Blue Part: DGN153201-PBU9	 Pretty in Pink Part: DGN153201-PPK3	 Grey Street Part: DGN153201-PSL3	 Green Machine Part: DGN153201-PGN11
 Back in Black Part: DGN153201-PBK1	 Evolution Orange Part: DGN153201-POR7	 Midnight Blue Part: DGN153201-PBU10	 Military Green Part: DGN153201-PGN8

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height is 15.875" with 6" stroke and 16.875" with 7" stroke. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. *Please verify the appropriate billing code to be used for this item with all third-party payors.

Power Elevating Seat with 7" Stroke \$1,990 Part: XRF000131. HCPCS: E2298/E2300^	Power Elevating Seat with 6" Stroke \$1,990 Part: XRF000132. HCPCS: E2298/E2300^
---	---

5. BASE MOUNTED LIGHTS

Optional. Select only one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

5.1a. Fender Lights

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights (without iLevel 12" Lift System) \$263
Part: XRF000177



5.1b. Base Mounted Light

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129
Requires Q-Logic 3 Joystick for Full Light Package.

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21 Part: MEC1910373 Works with .5" x .5" vertical supports only.	Backup Camera Track Mount Bracket \$27 Part: MEC2104227
--	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



7. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color if applicable. Knobby tires increase overall base width by 1.5".

7.1. Drive Wheels

3"x14" Flat Free Black Tires (No Accent Color)..... Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires (No Accent Color)..... No Charge Part: WHL150034 x2
3"x14" Flat Free Black Tires for Wheel Accents..... No Charge Part: WHL1611297 x2	3"x14" Pneumatic Black Tires for Wheel Accents..... No Charge Part: WHL1611296 x2
4"x14" Knobby, Flat Free Black Tires..... \$580 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires..... \$580 Part: XRF000367

7.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black Accent Color..... \$63 Part: XRF000179	Deep Purple Accent Color..... \$63 Part: XRF000184	Electric Blue Accent Color \$63 Part: XRF000188
Red Corvette Accent Color \$63 Part: XRF000180	Pretty in Pink Accent Color \$63 Part: XRF000185	Lemon Crush Accent Color..... \$63 Part: XRF000189
Night Train Accent Color \$63 Part: XRF000181	Evolution Orange Accent Color \$63 Part: XRF000186	Green Machine Accent Color \$63 Part: XRF000190
Grey Street Accent Color..... \$63 Part: XRF000182	Midnight Blue Accent Color \$63 Part: XRF000187	Military Green Accent Color \$63 Part: XRF000191

8. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea. Part: BAT2004376. HCPCS: E2361	NF-22, Gel Cell, MK Battery \$432 ea. Part: BATGEL1008. HCPCS: E2361
NF-22, AGM, MK Battery..... \$432 ea. Part: BAT121346. HCPCS: E2361	

9. SEAT TOWERS**Optional.** Not available with Power Adjustable Seat Height.

Extended Seat Towers No Charge

Part: MEC120301

Increases seat to floor height.

10. ELECTRONICS**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**10.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard

Part: CTL2205379

Controls up to 2 actuators.

QLNES 10 Key Joystick \$420

Part: CTL2205380

Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

**10.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module \$730

Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness \$458

Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick \$840

Part: CTL167724

Q-Logic 3e Joystick \$735

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

Q-Logic 3 EX Joystick for Full Lights \$840

Part: CTL167725

No Joystick, Setup for Q-Logic

Part: XRF000010

Requires selection of Enhanced Display

**10.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60

Part: ELE2103635

10.3. Q-Logic 3 Joystick Armor**REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch Standard

Part: DGN124307

Also matches Night Train shrouds.

Deep Purple \$63

Part: DGN124336-PPR5

Electric Blue \$63

Part: DGN124336-PBU9

Red Corvette \$63

Part: DGN124336-PRD1

Pretty in Pink \$63

Part: DGN124336-PPK3

Lemon Crush \$63

Part: DGN124336-PYL4

Back in Black \$63

Part: DGN124336-PBK1

Evolution Orange \$63

Part: DGN124336-POR7

Green Machine \$63

Part: DGN124336-PGN11

Grey Street \$63

Part: DGN124336-PSL3

Midnight Blue \$63

Part: DGN124336-PBU10

Military Green \$63

Part: DGN124336-PGN8

11. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90

Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G) \$90

Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C) \$90

Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Bodypoint 3" Goal Post Handle \$104

Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D) \$126

Part: ACCASMB1858. HCPCS: E2323

Bodypoint 4" Goal Post Handle \$104

Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E) \$69

Part: ACCASMB1859. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge

Part: DGN124308


Soft Ball Handle (F) \$80

Part: ACCASMB1856. HCPCS: E2323



12. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**12.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**12.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.


Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**12.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**12.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032Multi-Axis Mount \$126
Part: XRF000004

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Allows height adjustment of joystick.

**12.2. Joystick Mount Location**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**13. SYNERGY SEATING SYSTEM****Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.**13.1. Synergy Frame Selection****REQUIRED.**Synergy Seating System \$456
Part: PTOINDV2995**13.2. Synergy Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>☐ WC-19 Occupied Transit Option \$289
Part: KITASMB1830**13.3. Seat Width****REQUIRED.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.4. Seat Width Package**REQUIRED.** Maximum seat width is 22".10"-14" Width Package
Part: FRMASMB1801114"-18" Width Package
Part: FRMASMB1801216"-20" Width Package
Part: FRMASMB1801320"-24" Width Package
Part: FRMASMB18014**13.5. Seat Depth****REQUIRED.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.6. Seat Depth Package**REQUIRED.** Maximum seat depth is 22".10"-14" Depth Package
Part: FRMASMB698114"-18" Depth Package
Part: FRMASMB1782315"-20" Depth Package
Part: FRMASMB773220"-24" Depth Package
Part: FRMASMB7733


Low	Med	High
Synergy Seat		
15.875"	16.875"	17.875"
Extended Towers		
Add 3 inches to above		
Seat to Floor Height		

13.7. Synergy Seat to Floor**REQUIRED.**Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**13.8. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$39

Part: POP1909006-S

14"W and smaller

Small Wedge Positioning Kit.....\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit.....\$105

Part: POP1909006-L

14"W and smaller

Large Wedge Positioning Kit.....\$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

FRG Top Foam

Reference part: XRF000381

No Top Foam

Reference part: XRF000382

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

Stealth Polyurethane Liner

Reference part: XRF000384

No Liner

Reference part: XRF000385

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Mesh Cover

Reference part: XRF000387

Stealth Tek Cover

Reference part: XRF000388

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**7. Additional Cushion Liner**

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

 QTY Reference Part: XRF000390

Additional Stealth Polyurethane Liner

 QTY Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

 QTY Reference Part: XRF000392

Additional Stealth Tek Cover

 QTY Reference Part: XRF000394

Additional Mesh Cover

 QTY Reference Part: XRF000393

13.9. Synergy Seating Back Cane Selection**REQUIRED.**

Synergy Back Cane No Charge
 Part: FRMASMB15845

TRU-Comfort Back Canes No Charge
 Part: XRF000133

Flip Forward Back Cane \$184
 Part: FRMASMB15846
 Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$205
 Part: FRMASMB12067

2" Endomorph (1" per side) Back Cane \$205
 Part: FRMASMB12115

4" Endomorph (2" per side) Back Cane \$205
 Part: FRMASMB12139

Manual Recline Back \$788
 Part: FRMASMB12910. HCPCS: E1226
 Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.

13.10. Synergy Seating Back Cane Height**REQUIRED with Synergy seating EXCEPT if a TRU-Comfort back is desired.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes
 Part: FRMASMB15832

17" High Canes
 Part: FRMASMB15835

20" High Canes
 Part: FRMASMB15838

23" High Canes
 Part: FRMASMB15841

15" High Canes
 Part: FRMASMB15833

18" High Canes
 Part: FRMASMB15836

21" High Canes
 Part: FRMASMB15839

24" High Canes
 Part: FRMASMB15842

16" High Canes
 Part: FRMASMB15834

19" High Canes
 Part: FRMASMB15837

22" High Canes
 Part: FRMASMB15840

25" High Canes
 Part: FRMASMB15843

13.11. Synergy Seating Back Cushion**REQUIRED.** Select one type of back.**1a. Synergy Back Selection**

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

Synergy Back \$289
 Part: XRF000363

1b. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744
 Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort 4 Back, Back Cane Option, 22" Width \$855
 Part: TC4V22W32WBC. HCPCS: E2621

1c. ADI Back, AL Series, 4-Point, Active Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
 Part: ALBR16SA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$622
 Part: ALBR18XA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622
 Part: ALBR16TA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622
 Part: ALBR20SA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622
 Part: ALBR16XA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622
 Part: ALBR20TA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622
 Part: ALBR18SA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622
 Part: ALBR20XA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$622
 Part: ALBR18TA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height \$735
 Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height \$735
 Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height \$735
 Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
 Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height \$735
 Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
 Part: ALBR20XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height \$735
 Part: ALBR18SD4. HCPCS: E2615

Synergy seating section continued on following pages.

13. SYNERGY SEATING SYSTEM CONTINUED**1e. ADI Back, AL Series, 4-Point, Extra Deep Contour**

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
Part: ALBR18S64. HCPCS: E2620

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
Part: ALBR18X64. HCPCS: E2620

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
Part: ALBR18T64. HCPCS: E2620

1g. Omit Back Selection

Omit Back..... No Charge
Part: PTOINDV1287

13.12. Synergy Seating Overall Back Height**REQUIRED.**

12" Overall Back Height
Part: XRF000134

19" Overall Back Height
Part: XRF000141

13" Overall Back Height
Part: XRF000135

20" Overall Back Height
Part: XRF000142

14" Overall Back Height
Part: XRF000136

21" Overall Back Height
Part: XRF000143

15" Overall Back Height
Part: XRF000137

22" Overall Back Height
Part: XRF000144

16" Overall Back Height
Part: XRF000138

23" Overall Back Height
Part: XRF000145

17" Overall Back Height
Part: XRF000139

24" Overall Back Height
Part: XRF000146

18" Overall Back Height
Part: XRF000140

25" Overall Back Height
Part: XRF000147

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	-	-	-	-	-	-	-	-
	16"	S	S	-	-	-	-	-	-	-	TC
	17"	-	S	S	S	S	S	S	S	-	S TC
	18"	-	S	S	S	S	S	S	S	-	S TC
	19"	-	S	S	S	-	S	S	S	S	S TC
	20"	-	S	S	S	-	S TC	S	S	S	S TC
	21"	-	S	S	S	S	S TC	S	S	S	S TC
	22"	-	S	S	S TC	S TC	S TC	S	S TC	S	S TC
	23"	-	S	S	S TC	S TC	S	S	S TC	S	S TC
	24"	-	S	S	S	S	S TC	S	S TC	S	S
	25"	-	S	S	S	S	S TC	S	S TC	S	S

☐ Synergy ☒ TRU-Comfort 4

Synergy and TRU-Comfort Back Availability Chart

13.13. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... \$63
L R Left Part: FRMASMB7038. Right part: FRMASMB7039

☐ ☐ Desk Length, Flip Back Armrests \$63
L R Left part: FRMASMB7040. Right part: FRMASMB7041

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests..... \$158 ea.
L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests \$158 ea.
L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.
L R Left Part: XRF000245. Right part: XRF000246. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea.
L R Left part: XRF000247. Right part: XRF000248. HCPCS: E0973

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.
L R Left Part: XRF000249. Right part: XRF000250. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
L R Left part: XRF000251. Right part: XRF000252. HCPCS: E0973

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

☐ ☐ Pediatric, Removable Armrests \$79 ea.
L R Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests \$103 ea.
L R Left Part: FRMASMB16188. Right part: FRMASMB16189

☐ ☐ Desk Length, Cantilever Armrests..... \$103 ea.
L R Left part: FRMASMB16190. Right part: FRMASMB16191

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$79 ea.
L R Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests \$79 ea.
L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

1h. Drop-In, Height Adjustable Armrests (T-Arm)

Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests..... No Charge
L R Left part: ARM1803253. Right part: ARM156742

☐ ☐ Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743.

Synergy seating section continued on following pages.

Account Number: _____ Date: _____ Quote Name: _____

13. SYNERGY SEATING SYSTEM CONTINUED

13.14. Synergy Armpads

REQUIRED with armrest selection above. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard Waterfall Armpads..... No Charge
 Part: XRF000148 Part: XRF000149

No Armpad Selection No Charge
 Part: PTOINDV1373

13.15. Joystick Bracket Receiver

REQUIRED.

☐ Standard Fixed Receiver Standard ☐ Stealth Height Adjustable Receiver..... \$236
 Part: FRMASMB16062 Part: XRF000064. HCPCS: K0108

14. HEADRESTS

Optional. Select a headrest pad and bracket.

14.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179 Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
 Part: ST-CP08CCFPM-Q. HCPCS: E0955 Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179 Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
 Part: ST-CP10CCFPM-Q. HCPCS: E0955 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311 Multi-Axis Removable Headrest Bracket.....\$174
 Part: ST-UHB2D33200-Q. HCPCS: E1033 Part: ST-HMO475-17M-Q. HCPCS: E1033

14.2. Flip-down Headrest Bracket

Stealth Flip-down Headrest Bracket (FDM350)..... \$170
 Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

14.3. Q Cam Lock

Optional.

Q Cam Lock..... \$50
 Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

15. LAP BELTS

REQUIRED. Select one lap belt below.

15.1a. Standard Lap Belts

☐ 38" Lap Belt ☐ 50" Lap Belt ☐ 60" Lap Belt ☐ 70" Lap Belt ☐ 80" Lap Belt ☐ 90" Lap Belt
 Part: ACCBELT1006 Part: ACCBELT1000 Part: ACCBELT1002 Part: ACCBELT1003 Part: ACCBELT1004 Part: ACCBELT1001

15.1b. Optional Lap Belts

☐ Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50 ☐ Stealth Positioning 60" Single Pull, Push Button, Extra Large.... \$50
 (7.5") Padded, 1.5" Belt (IFPB5213) Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5213-Q Part: ST-IFPB5315-Q

☐ Stealth Positioning 60" Single Pull, Push Button, Large \$50 ☐ 60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Padded (9.5"), 2" Belt (IFPB5314) Required with ISO 7176-19 Compliant Transit Option.
 Part: ST-IFPB5314-Q Part: ACCBELT1016

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1a. Base Mounted Foot Platform

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Omit base mounted foot platform No Charge
 Part: FRMASMB16077 Part: PTOINDV3674

Bell Shaped, Base Mounted Foot Platform \$247
 Part: RIGASMB7120041

16.1b. Basic Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Basic Center Mount..... No Charge
 Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

Legrests section continued on following page.

16. LEGRESTS CONTINUED**2. Basic Center Mount Extension and Angle Adjustable Footplates**

Short Extension, 5"x9" Footplates \$180 Part: XRF000312. HCPCS: K0040	Long Extension, 5"x9" Footplates \$180 Part: XRF000314. HCPCS: K0040
Medium Extension, 5"x9" Footplates \$180 Part: XRF000313. HCPCS: K0040	

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$158 Part: RIG2104369	Large Calf Pads, Basic Center Mount \$158 Part: RIG2202997 Only available with Long Length Extension.
---	---

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16 Part: RIG2104240

16.1c. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Foot Platform Assembly

Power Articulating Foot Platform a\$3,145 Part: PTOINDV3754. HCPCS: E1012+K0040
--

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

2. AFP Lower Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Medium Extension, 5"x6" Footplates Part: RIG150206	Medium Extension, 5"x11" Footplates Part: RIG2404018	Long Extension, 5"x8" Footplates Part: RIG2404017
Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x6" Footplates Part: RIG151505	Long Extension, 5"x11" Footplates Part: RIG2404019

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Part: RIGASMB7120018	Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard Part: RIGASMB7120019
--	--

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H \$184 Part: FRMASMB7110012	Large Calf Panel, 24"Wx8"H \$237 Part: RIG123825 Only available on widths of 20" or greater.
--	--

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$69 Part: RIG2404921 Available with 5"x8" footplates.
Width Extension for 8" \$69 Part: RIG2404919	Front and Rear Cap \$16 Part: RIG2404923 Available with 5"x8" footplates.
Width Extension for 11" \$69 Part: RIG2404920	

5. Additional Accessories

Small Foam Calf Pad Cover (pair) \$67 Part: ST-CSCS0775-Q	Large Foam Panel Cover (pair) \$282 Part: ST-CSCP0927-Q
Large Foam Calf Pad Cover (pair) \$83 Part: ST-CSCS0985-Q	Footplate Connector \$16 Part: RIGASMB7120023
Small Foam Panel Cover (pair) \$176 Part: ST-CSCP0920-Q	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension \$247 Part: ST-AFPG5008E-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers \$270 Part: ST-AFPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension \$270 Part: ST-AFPG5011E-Q. HCPCS: K0108

16.1d. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$86 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084 Not available with 11"-13" extensions.	☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087 For use only with 11"-13" lower extensions.
--	--

Legrests section continued on following pages.

16. LEGRESTS CONTINUED

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Extruded, 13" - 16.5"	!	✓	✓
Extruded, 15.5" - 19.5"	!	!	!
Angle Adjustable, 13" - 16.5"	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	!	!	!
Multi-Axis, 15.5" - 19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
☒ This rigging configuration does not work with this seat-to-floor height.

2. Lower Extension

P L R

11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

P L R

13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

P L R

15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

P L R

5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R

6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R

6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R

8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

16.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R

Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

P L R

11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

P L R

14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

P L R

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

P L R

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R

7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R

8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**3b. Angle Adjustable Footplate(s)**

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

16.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$37
 Part: XRF000100

16.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

17. SYNERGY POWER OPTION ELECTRONICS

REQUIRED for Synergy Seat with Power Legrests or Elevating Seat. Select the desired control type for the seating system.

Single Actuator Through Controller Interface \$2,163
 Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and Interface \$3,040
 Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and Interface (AAM) \$3,040
 Part: XRF000068. HCPCS: E2310

18. SEAT MOUNTED LED LIGHTS

Optional. Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045
 Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Synergy Seat Cup Holder No Charge Part: ACCASMB2634	Grip Lap Board (20005)..... \$38.95 Part: ACC152992
Oxygen Holder for Synergy Seat..... \$158 Not available with Synergy Manual Recline. Part: FRMASMB7030. HCPCS: E2208	Grip Drink Holder (20085) \$15.50 Part: ACC152993
Walker Holder \$109 Synergy part: FRMASMB7031	Accessory Bar \$132 Part: XRF000128
Cane/Crutch Holder \$105 Synergy part: FRMASMB4822. HCPCS: E2207	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
Rear Basket \$124 Synergy part: FRMASMB4821	Privacy Bag \$42 Part: ACC1907570
Flag Holder \$55 Synergy/Flex part: ACCASMB2415.	Q-Logic EX/NE Programmer \$573 Part: ELEASMB6976
Saddle Bag \$40 Part: DWR1020D010	Q-Logic PC Programming Disk with Cable \$445 Part: ELEASMB5216
Stealth Cell Phone/Media Device Holder \$215 Part: ACC7120006	Personal Item Pouch \$84 Part: ACC2101017
Stealth Large Tablet Holder \$625.28 Part: ST-URPKTAB0000-Q	Transfer Handle \$105 ea. Part: MEC2205680 x2. HCPCS: E0705 Available with Flip-Up armrests only.
Quantum Backpack \$168 Part: ACC1905853	Remote Stop Switch \$1,295 Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
 QTY Part: ACC2407046

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
 QTY Part: ACC2405496

20. BRANCH PATCHES**Optional.**

Air Force Patch \$21 Part: ACC140240	Coast Guard Patch \$21 Part: ACC140239	Navy Patch \$21 Part: ACC140244
Army Patch \$21 Part: ACC140241	Marines Patch \$21 Part: ACC140243	Purple Heart Patch \$21 Part: ACC140238

**21. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#).**21.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$562 Part: INDPART2896. HCPCS: E1034	Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1034 Only available with back cane and omit back configurations.
Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664 Part: ST-ULPK404108B00-Q. HCPCS: E1034 Only available with TRU Comfort backs.	

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250 Part: INDPART2663. HCPCS: E0956	Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250 Part: INDPART2664. HCPCS: E0956
---	---

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

4. Lateral Extensions**REQUIRED with Synergy Back and thoracic lateral selections.**

Lateral Extension for Synergy Back (pair)..... \$56 Part: XRF000150
--

Positioning components section continued on following page.

21. POSITIONING COMPONENTS CONTINUED**21.2. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**


^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^Flip Down with LEFT Full Surface Contact \$281
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^Flip Down with RIGHT Full Surface Contact..... \$281
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**

^ Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^**3. UniLink Pelvic/Thigh Additional Links**UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-QUniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q**4. UniLink Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**21.4. Stealth Positioning Chest Harness**Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960**21.5. Stealth Positioning Chest Straps**Chest Strap, Medium w/ Belt Anchor (IFCS4113).. \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor (IFCS4114)..... \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960**21.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028**2. Facial Lateral Pads and Accessories**Egg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea.
Part: ST-920M-Q. HCPCS: K0108Swing-Away Joystick Mount,..... \$276 ea.
18" Rod (C)Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest**21.7. Stealth Elbow Support****1. Elbow Support Bracket**Left, Swing-Away, Elbow Support Bracket \$199 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Right, Swing-Away, Elbow Support Bracket..... \$199 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108

Positioning components section continued on following page.

21. POSITIONING COMPONENTS CONTINUED**21.8. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

**21.9. Stealth Positioning Shoe Holder**

Small Shoe Holder, \$87 ea.
QTY 3/4" Webbing (IFSH4211)
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, \$87 ea.
QTY 3/4" Webbing (IFSH4212)
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, \$87 ea.
QTY 3/4" Webbing (IFSH4213)
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

**22. LAP TRAY**

Optional.

22.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	M	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	M	ML	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	M	ML	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	M	ML	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						M		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	M	ML	ML	ML	L	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

23. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below or to refer to the specialty controls order form.

23.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN137728

Red Corvette..... \$63
Part: DGN137729-PRD1

Back in Black..... \$63
Part: DGN137729-PBK1

Night Train..... \$63
Part: DGN137729-PBK3

Grey Street..... \$63
Part: DGN137729-PSL3

Deep Purple..... \$63
Part: DGN137729-PPR5

Pretty in Pink \$63
Part: DGN137729-PPK3

Midnight Blue \$63
Part: DGN137729-PBU10

Electric Blue \$63
Part: DGN137729-PBU9

Lemon Crush \$63
Part: DGN137729-PYL4

Green Machine \$63
Part: DGN137729-PGN11

Military Green \$63
Part: DGN137729-PGN8



Alternative controls and electronics section continued on following page.

23. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

23.2. Q-Logic 3 SAJ

- ☐ Q-Logic 3 EX SAJ \$840
 Part: CTL143980
 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

23.3. Head Array

- | | |
|--|--|
| i-Drive 4.0 Bluetooth Tri-array Series, \$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330 | Adult, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330 |
| i-Drive 4.0 Bluetooth Pro Series, \$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330 | Junior, Curved, i-Drive Tri-Array with \$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330 |



Adult Pro-Series

23.4. Egg Switches

- | | |
|--|---|
| Red Egg Switch (ESRED) \$87 ea.
Part: SWTMCRO1077 | Green Egg Switch (ESGRN) \$87 ea.
Part: SWTMCRO1078 |
| Blue Egg Switch (ESBLU) \$87 ea.
Part: SWTMCRO1079 | Yellow Egg Switch (ESYEL) \$87 ea.
Part: SWTMCRO1080 |
| Black Egg Switch (ESBLK) \$87 ea.
Part: SWTMCRO1081 | |



Egg Switch

23.5. Switch Adapter

- 2 to 1 Switch Adapter \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

23.6. Environmental Controls

- Q-Logic 3 ECU Module Kit \$2,520
 Part: CTL168057