



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with TRU Balance® 3 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 2.0 3S-SS\$8,533

Part: Q6 EDGE 2.0 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 2.0 3SP-SS\$8,234

Part: Q6 EDGE 2.0 3SP-SS. HCPCS: K0856

Single power option.

☐ Quantum Q6 Edge 2.0 3MP-SS\$8,055

Part: Q6 EDGE 2.0 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. DRIVE MOTORS

REQUIRED.

- ☐ 6 mph Drive Motors Standard
Part: MOT147433
- ☐ Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- | | | | |
|---|--|--|--|
|  Red Corvette
Part: DGN153201-PRD1 |  Night Train
Part: DGN153201-PBK3 |  Deep Purple
Part: DGN153201-PPR5 |  Lemon Crush
Part: DGN153201-PYL4 |
|  Electric Blue
Part: DGN153201-PBU9 |  Pretty in Pink
Part: DGN153201-PPK3 |  Grey Street
Part: DGN153201-PSL3 |  Green Machine
Part: DGN153201-PGN11 |
|  Back in Black
Part: DGN153201-PBK1 |  Evolution Orange
Part: DGN153201-POR7 |  Midnight Blue
Part: DGN153201-PBU10 |  Military Green
Part: DGN153201-PGN8 |

4. BASE MOUNTED LIGHTS

Optional. Select only one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

4.1a. Fender Lights

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Switch Fender Lights (without iLevel 12" Lift System) \$263
Part: XRF000177


Switch Fender Lights, Setup for iLevel 12" Lift System ... No Charge
Part: DGN154117

4.1b. Base Mounted Light

Base Mounted, LED, Full Light Package \$1045
Part: XRF000129,
Requires Q-Logic 3 Joystick for Full Light Package.



Base Mounted Lights



Fender Lights

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$363
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket..... \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

6. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires (No Accent Color)..... Standard
Part: PTOINDV3525 (pair)

3"x14" Flat Free Black Tires for Wheel Accents..... No Charge
Part: WHL1611297 x2

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

3"x14" Pneumatic Black Tires (No Accent Color)..... No Charge
Part: WHL150034 x2

3"x14" Pneumatic Black Tires for Wheel Accents..... No Charge
Part: WHL1611296 x2

4"x14" Knobby, Pneumatic Black Tires..... \$580
Part: XRF000367

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black Accent Color..... \$63
Part: XRF000179

Red Corvette Accent Color \$63
Part: XRF000180

Night Train Accent Color \$63
Part: XRF000181

Grey Street Accent Color..... \$63
Part: XRF000182

Deep Purple Accent Color..... \$63
Part: XRF000184

Pretty in Pink Accent Color \$63
Part: XRF000185

Evolution Orange Accent Color \$63
Part: XRF000186

Midnight Blue Accent Color \$63
Part: XRF000187

Electric Blue Accent Color \$63
Part: XRF000188

Lemon Crush Accent Color..... \$63
Part: XRF000189

Green Machine Accent Color \$63
Part: XRF000190

Military Green Accent Color \$63
Part: XRF000191

7. BATTERIES**Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376. HCPCS: E2361NF-22, Gel Cell, MK Battery \$432 ea.
Part: BATGEL1008. HCPCS: E2361NF-22, AGM, MK Battery \$432 ea.
Part: BAT121346. HCPCS: E2361**8. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**8.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick \$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

QLNES Joystick

8.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 Controller Module \$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness \$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**Q-Logic 3e Joystick \$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.Q-Logic 3 EX Joystick for Full Lights \$840
Part: CTL167725Q-Logic 3 EX Joystick \$840
Part: CTL167724No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display

Q-Logic 3 Joysticks

8.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635**8.3. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds.Deep Purple \$63
Part: DGN124336-PPR5Electric Blue \$63
Part: DGN124336-PBU9Red Corvette \$63
Part: DGN124336-PRD1Pretty in Pink \$63
Part: DGN124336-PPK3Lemon Crush \$63
Part: DGN124336-PYL4Back in Black \$63
Part: DGN124336-PBK1Evolution Orange \$63
Part: DGN124336-POR7Green Machine \$63
Part: DGN124336-PGN11Grey Street \$63
Part: DGN124336-PSL3Midnight Blue \$63
Part: DGN124336-PBU10Military Green \$63
Part: DGN124336-PGN8**9. JOYSTICK HANDLES****Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.☐ Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323☐ Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323☐ Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)☐ Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323☐ T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323☐ Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323☐ Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308☐ Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323

Joystick Handles

10. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**10.1a. Fixed, In-line Joystick Bracket**
☐ Fixed, In-line, Bracket for Q-Logic 3 Standard
 Part: FRMASMB10617
10.1.b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.



 Swing-away, Q-Logic 3 \$294
 Part: ACC1908434. HCPCS: E1032
10.1c. Basic Swing Away Joystick Mount
 Basic Swing Away, Q-Logic, Left Side \$245
 Part: ACC2201041. HCPCS: E1032

 Basic Swing Away, Q-Logic, Right Side \$245
 Part: ACC2201042. HCPCS: E1032
10.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

 Retract4 Mount \$368
 Part: XRF000003. HCPCS: E1032
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

 Multi-Axis Mount \$126
 Part: XRF000004
 Allows height adjustment of joystick.
10.2. Joystick Mount Location
 Left Side Mount
 Part: PTOINDV1225

 Right Side Mount
 Part: PTOINDV1226
**11. TRU BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™****REQUIRED.** Maximum seat size is 22"W x 22"D.**11.1. TRU Balance 3 Configuration**
 TRU Balance 3 Static Seat \$2,410
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

 TRU Balance 3 Power Recline Only \$8,243
 Part: PTOINDV3444. HCPCS: E1004
 175° of power recline.

 TRU Balance 3 Power Tilt Only \$8,202
 Part: PTOINDV3443. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

 TRU Balance 3 Power Tilt & Recline \$12,615
 Part: PTOINDV3446. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.
11.2. 12" Power Lift**Optional.** Select one option below. ^Please verify the appropriate billing code to be used for this item with all third-party payors.
 TRU Balance 3, 12" Power Lift with iLevel® \$3,725
 Part: PTOINDV3773. HCPCS: E2298/E2300^
 12" of scissor lift with iLevel® System.
**11.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Only available with 16"-20" width and depth packages. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>
 WC-19 Occupied Transit Option \$289
 Part: KIT152047
11.4. Width**REQUIRED with seat selection.** Maximum seat width is 22". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.5. Width Package**REQUIRED with seat selection.**
 12"-16" Width Package
 Part: XRF000006

 14"-16" Width Package
 Part: XRF000007

 16"-20" Width Package
 Part: XRF000008

 20"-24" Width Package
 Part: XRF000158
 Maximum width is 22".
11.6. Depth**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------


TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**11.7. Depth Package****REQUIRED with seat selection.****12"-16" Depth Package**
Part: XRF000011**14"-16" Depth Package**
Part: XRF000012**16"-20" Depth Package**
Part: XRF000013**20"-24" Depth Package**
Part: XRF000159
Maximum depth is 22".

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
17.5"	18.5"	19.5"
TRU Balance 3 Seat to Floor Heights		

11.8. Seat to Floor Height**REQUIRED.****Lowest Seat to Floor**
Part: PTOINDV3558**Middle Seat to Floor**
Part: PTOINDV3559**Highest Seat to Floor**
Part: PTOINDV3560**11.9. Stealth Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608**Essence SPP**Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608**Solution SPP**Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608**Spectrum Foam SPP**Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608**Spectrum Gel SPP**Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608**Premiere P**Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606**Glacial SP**Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604**Zen SP**Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604**Simplicity G**Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$39
Part: POP1909006-S
14"W and smaller**Small Wedge Positioning Kit.....\$39**
Part: POP1712340-S
15"W and larger**Medium Wedge Positioning Kit.....\$72**
Part: POP1909006-M
14"W and smaller**Medium Wedge Positioning Kit.....\$72**
Part: POP1712340-M
15"W and larger**Large Wedge Positioning Kit.....\$105**
Part: POP1909006-L
14"W and smaller**Large Wedge Positioning Kit.....\$105**
Part: POP1712340-L
15"W and larger

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**3. Cushion Top Foam**

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390

Additional Stealth Polyurethane Liner
Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392

Additional Stealth Tek Cover
Reference Part: XRF000394

Additional Mesh Cover
Reference Part: XRF000393

11.10. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back \$184
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Not available with power recline. Only available with 2-post or pediatric armrests. Select a back cushion type in the subsequent section.

11.11. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane No Charge
Part: MEC140266

2" Endomorph (1" per side) Back Cane \$74
Part: XRF000269

Flip Forward Back Cane \$53
Part: XRF000273

4" Endomorph (2" per side) Back Cane \$74
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

11.12. Back Cushion Configuration

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
Part: XRF000017
Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287
Select a following width package.

After-market Back Interface \$535
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall armrests and arm through. Does not include brackets to mount after-market backs.

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**1b. TRU-Comfort Back**

TRU-Comfort Back (Static/Flip Forward), 12" Width \$744
 Part: SET123622. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width..... \$744
 Part: TC4V10W21W. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
 Part: TC4V10W21WBC. HCPCS: E2620
 Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
 Part: TC4V22W32W. HCPCS: E2621
 Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$855
 Part: TC4V22W32WBC. HCPCS: E2621
 Complete following width package and size.

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
 Part: ALBR16SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622
 Part: ALBR16TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622
 Part: ALBR16XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622
 Part: ALBR18SA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$622
 Part: ALBR18TA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$622
 Part: ALBR18XA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622
 Part: ALBR20SA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622
 Part: ALBR20TA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622
 Part: ALBR20XA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735
 Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735
 Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735
 Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735
 Part: ALBR18SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735
 Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
 Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
 Part: ALBR18S64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
 Part: ALBR18T64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
 Part: ALBR18X64. HCPCS: E2615

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
 Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
 Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
 Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$568
 Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$568
 Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$568
 Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$568
 Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703
 Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AL Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1h. Q-Back 2.0, AL Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline..... \$703 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL18ET. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620

1i. Q-Back 2.0, AFT Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

1j. Q-Back 2.0, AFT Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT18DS. HCPCS: E2620	

1k. Q-Back 2.0, AFT Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AFT16EXX. HCPCS: E2620	
---	--

11.13. TRU Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Evolution Orange Part: XRF000028	Lemon Crush Part: XRF000032
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Midnight Blue Part: XRF000029	Green Machine Part: XRF000033
Night Train Part: XRF000022	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Military Green Part: XRF000034

11.14. Lateral Positioning Wedge**Optional.**

Lateral Positioning Wedge, pair \$32 Part: SET123465	
---	--

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

TRU Balance 3 Seating section continued on following pages.

11. TRU BALANCE 3 SEATING CONTINUED**11.15. Back Width Package****REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003718"W-22"W Package
Part: XRF00003910"W-14"W Package
Part: XRF00003616"W-20"W Package
Part: XRF00003820"W-24"W Package
Part: XRF000040

	Tilt or Static Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"		
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"	
	15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"	
	16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	16"	
	17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	17"	
	18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	18"	
	19"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	19"
	20"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	20"
	21"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	21"
	22"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	22"
	23"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	23"
	24"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	24"
	25"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	25"
	26"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	26"
	27"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	27"
	28"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	28"
	29"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	29"
	30"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	30"
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	31"	
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	32"	

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

11.16. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.12" Back Width
Part: PTOINDV356216" Back Width
Part: PTOINDV356419" Back Width
Part: PTOINDV356722" Back Width
Part: PTOINDV357014" Back Width
Part: PTOINDV356317" Back Width
Part: PTOINDV356520" Back Width
Part: PTOINDV356824" Back Width
Part: PTOINDV357115" Back Width
Part: PTOINDV385618" Back Width
Part: PTOINDV356621" Back Width
Part: PTOINDV3569**11.17. Overall Back Height**

Static Seat and Power Tilt available height range: 14"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion.

Power Recline available height range: 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height
Part: SETHGHT102718" Height
Part: SETHGHT103122" Height
Part: SETHGHT103526" Height
Part: SETHGHT103930" Height
Part: SETHGHT104315" Height
Part: SETHGHT102819" Height
Part: SETHGHT103223" Height
Part: SETHGHT103627" Height
Part: SETHGHT104031" Height
Part: SETHGHT110216" Height
Part: SETHGHT102920" Height
Part: SETHGHT103324" Height
Part: SETHGHT103728" Height
Part: SETHGHT104132" Height
Part: SETHGHT110317" Height
Part: SETHGHT103021" Height
Part: SETHGHT103425" Height
Part: SETHGHT103829" Height
Part: SETHGHT1042**11.18. Power Tilt Privacy Flap****Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.Privacy Flap, 9.75" \$48
Part: XRF000351**11.19. Back Height Extension****Optional. Not available with back cane option, 12" back width, or 24" back width.** Adds 2" of seat back pan height. 25"H selection includes a back height extension.Back Height Extension \$48
Part: XRF000041

12. TRU BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**12.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

- | | |
|--|--|
| ☐ ☐ Flip-Back, 8"-16"H Armrests Standard
L R Left part: ARM1806662. Right part: ARM1806663
Not available with TRU Balance 3 Static Seat. | ☐ ☐ Flip-Back, 8"-16"H Armrests for Static Seat \$289
L R Left part: XRF000320. Right part: XRF000321 |
|--|--|

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

- | |
|--|
| ☐ ☐ Flip-Back Armrest Downpost \$92 ea.
L R Left part: XRF000265. Right part: XRF000266 |
|--|

12.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests**Only available with TRU Balance 3 Static Seat.** Includes flat armpads.

- | | |
|---|---|
| ☐ ☐ Full Length, Drop-In Armrests No Charge
L R Left part: ARM1803253. Right part: ARM156742 | ☐ ☐ Desk Length, Drop-In Armrests No Charge
L R Left part: ARM1803254. Right part: ARM156743 |
|---|---|

12.1c. Two-Post, Flip-Back, Height Adjustable, Removable Armrests

Select a left and right armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes. Not available with a depth less than 16". Required with flip forward back and lap trays.

1a. Two-Post Armrest for Power Seating

- | | |
|--|--|
| ☐ ☐ Full Length, Two-Post Armrest for Power Seating No Charge
L R Left part: XRF000047. Right part: XRF000048 | ☐ ☐ Desk Length, Two-Post Armrest for Power Seating No Charge
L R Left part: XRF000049. Right part: XRF000050 |
|--|--|

1b. Two-Post Armrest for Static Seating

- | | |
|---|---|
| ☐ ☐ Full Length, Two-Post Armrest for Static Seating \$158 ea.
L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 | ☐ ☐ Desk Length, Two-Post Armrest for Static Seating \$158 ea.
L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973 |
|---|---|

12.1d. Heavy Duty, Single Post, Height Adjustable, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

- | | |
|--|---|
| 8"-12"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000055. HCPCS w/ Static Seat: E0973 | 10"-14"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000056. HCPCS w/ Static Seat: E0973 |
|--|---|

12.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1. Left Armrest Height**

- | | | | | |
|---------------------|----------------------|----------------------|----------------------|----------------------|
| 8" H
PTOINDV3535 | 10" H
PTOINDV3537 | 12" H
PTOINDV3539 | 14" H
PTOINDV3541 | 16" H
PTOINDV3543 |
| 9" H
PTOINDV3536 | 11" H
PTOINDV3538 | 13" H
PTOINDV3540 | 15" H
PTOINDV3542 | |

2. Right Armrest Height

- | | | | | |
|---------------------|----------------------|----------------------|----------------------|----------------------|
| 8" H
PTOINDV3450 | 10" H
PTOINDV3452 | 12" H
PTOINDV3454 | 14" H
PTOINDV3456 | 16" H
PTOINDV3458 |
| 9" H
PTOINDV3451 | 11" H
PTOINDV3453 | 13" H
PTOINDV3455 | 15" H
PTOINDV3457 | |

12.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

- | | | |
|---|---|--|
| ☐ ☐ 10" Length
L R Left part: XRF000057. Right part: XRF000123 | ☐ ☐ 14" Length
L R Left part: XRF000058. Right part: XRF000124 | ☐ ☐ 17" Length
L R Left part: XRF000059. Right part: XRF000125
Only available with Flip-Back, 8"-16" H Armrests and straight arm pads. |
|---|---|--|

13. TRU BALANCE 3 ARMPADS**REQUIRED.** Please note that waterfall and gel armpads may cause interference with the seat back.**13.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

- | | |
|--|---|
| No Armpad Selection No Charge
Part: XRF000060
Not available with Flip-back, 8"-16" H armrests. | Interface for Aftermarket Pad No Charge
Part: XRF000308
Only available with Flip-back, 8"-16" armrests. |
|--|---|

13.1b. Foam Armpads

- | | |
|--|--|
| ☐ ☐ Straight, Armpad Standard
L R Left part: XRF000061. Right part: XRF000126 | ☐ ☐ Waterfall, Armpad No Charge
L R Left part: XRF000062. Right part: XRF000127
Not available with Flip-back, 8"-16" armrests. |
|--|--|

TRU Balance 3 armpads section continued on following page.

12. TRU BALANCE 3 ARMRESTS CONTINUED**13.1c. Gel Armpads**

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back ArmrestsFalcon Gel, Full Length Left and Right\$195
Part: ARM1806807Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806810Falcon Gel, Full Length Left, Desk Length Right.....\$171
Part: ARM1806813Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806806**b. Falcon Gel Armpads for Tube Style Armrests**Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437**14. TRU BALANCE 3 ARMREST BRACKETS****14.1. Joystick Bracket****REQUIRED.**Standard Fixed ReceiverStandard
Part: XRF000063Stealth Height Adjustable Receiver.....\$236
Part: XRF000064. HCPCS: K0108**15. TRU BALANCE 3 HEADRESTS****Optional.** Select a headrest pad and bracket.**15.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**15.2. Q Cam Lock****Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

16. LAP BELTS**REQUIRED.** Select a lap belt below.**16.1a. Standard Lap Belts**38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**16.1b. Optional Lap Belts**Stealth Retractable Lap Belt (IFRB5350)..... \$182
Part: ST-IFRB5350Stealth Positioning 60" Single Pull, Push Button, Extra Large....\$50
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-QStealth Positioning 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016Stealth Positioning 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q**17. TRU BALANCE 3 LEGREST CONTROLS****Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.☐ Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108☐ Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110☐ Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109

18. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

18.1a. Base Mounted Foot Platform

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform	Standard	Omit base mounted foot platform	No Charge
Part: FRMASMB16077		Part: PTOINDV3674	
Bell Shaped, Base Mounted Foot Platform	\$247		
Part: RIGASMB7120041			

18.1b. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount	No Charge
Part: XRF000323	

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates	\$180	Long Extension, 5"x9" Footplates	\$180
Part: XRF000312. HCPCS: K0040		Part: XRF000314. HCPCS: K0040	
Medium Extension, 5"x9" Footplates	\$180		
Part: XRF000313. HCPCS: K0040			

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount	\$158	Large Calf Pads, Basic Center Mount	\$158
Part: RIG2104369		Part: RIG2202997	
		Only available with Long Length Extension.	

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector	\$16
Part: RIG2104240	

18.1c. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform	\$3,145	Quick Adjustable Center Mount Foot Platform	\$263
Part: PTOINDV3754. HCPCS: E1012+K0040		Part: XRF000311. HCPCS: K0040	

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x11" Footplates Part: RIG2404018
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Long Extension, 5"x6" Footplates Part: RIG151505
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x8" Footplates Part: RIG2404017
Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x11" Footplates Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pad

Small Calf Pads, AFP Only, 6"Wx6"H (pair)	Standard	Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair)	Standard
Part: RIGASMB7120018		Part: RIGASMB7120019	
Small Calf Pads, Quick Adjustable Center Mount Only, 6"Wx6"H (pair)	\$158	Large Calf Pads, Quick Adjustable Center Mount Only, 8" taper to 6"Wx8"H (pair)	\$158
Part: XRF000156		Part: XRF000157	

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H	\$184	Large Calf Panel, 24"Wx8"H	\$237
Part: FRMASMB7110012		Part: RIG123825	
		Only available on widths of 20" or greater.	

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension	\$69	Outside Edge	\$69
Part: RIG152519		Part: RIG2404921	
Only available with the 5x6 footplates.		Available with 5"x8" footplates.	
Width Extension for 8"	\$69	Front and Rear Cap	\$16
Part: RIG2404919		Part: RIG2404923	
Width Extension for 11"	\$69	Available with 5"x8" footplates.	
Part: RIG2404920			

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**5. Additional Accessories**

Small Foam Calf Pad Cover (pair) \$67
 Part: ST-CSCS0775-Q

Large Foam Calf Pad Cover (pair) \$83
 Part: ST-CSCS0985-Q

Small Foam Panel Cover (pair) \$176
 Part: ST-CSCP0920-Q

Large Foam Panel Cover (pair) \$282.21
 Part: ST-CSCP0927-Q

Footplate Connector \$16
 Part: RIGASMB7120023

Footplate Riser \$75
 Part: XRF000348
 Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel \$50
 Part: FRMASMB7110011
 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270
 Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
 Part: ST-AFPG5011E-Q. HCPCS: K0108

18.1d. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,516
 Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix Center Mount	Seat to Floor Option			Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14"	12.5"-15"	12.5"-16"	14.5"-21" Lower Length	12.5"-15"	12.5"-16"	12.5"-17"

2. Lower Extension

8.5" to 14.5" Lower Length
 Part: XRF000344

12.5" to 21.5" Lower Length
 Part: XRF000347

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair
 Part: FRMASMB7448

5"Wx8"D, Pair
 Part: FRMASMB7450

6"Wx11"D¹, Pair
 Part: FRMASMB7452

5"Wx6"D, Pair
 Part: FRMASMB7449

6"Wx8"D¹, Pair
 Part: FRMASMB7451

8"Wx11"D², Pair
 Part: FRMASMB7453

18.1e. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

P L R 60°-90° Swing-Away Legrest, Power Seating No Charge
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Available with Power Positioning only. Not available with 11"-13" extensions.

P L R 60°-90° Swing-Away Legrest, Short, Power Seating No Charge
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 Available with Power Positioning only. For use only with 11"-13" lower extensions.

P L R 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
 Pair part: XRF000082. Left: XRF000083.
 Right: XRF000084
 Available with Static Seat only. Not available with 11"-13" extensions.

P L R 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
 Pair part: XRF000085. Left: XRF000086.
 Right: XRF000087.
 Available with Static Seat only. For use only with 11"-13" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
11"-13" extension causes no interference at all seat-to-floor heights.									

✓ This rigging configuration causes no interference issues at this seat-to-floor height.

! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**2. Swing-Away Legrests Lower Extension**

11"-13" Length
 P L R Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

13"-16.5" Length
 P L R Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

15.5"-19.5" Length
 P L R Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

5"Wx8"D with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

6"Wx11"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

6"Wx8"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

8"Wx11"D² with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

18.1f. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

Elevating Legrest Assembly \$171 ea.
 P L R Pair part: FRMASMB483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

11"-14" Length
 P L R Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

14"-17" Length
 P L R Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 P L R Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 P L R Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12.5"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

7.25"Wx6"D Reinforced Footplate..... Standard
 P L R Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

8.5"Wx6"D Reinforced Footplate Standard
 P L R Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

5"Wx8"D with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

6"Wx11"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

6"Wx8"D¹ with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

8"Wx11"D² with Heel Loops..... \$90 ea.
 P L R Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

18.2. Legrest Accessories

Heel Loops (pair)..... \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$19 ea.
 Part: XRF000100

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**18.3. Set Legrest Extension Length****REQUIRED** with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.**1a. Centermount Extension Length**

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

19. TRU BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options.**19.1a. Tilt Only through Toggle**

Only available with QLNE5 Joystick.

Single Actuator Control, Power Tilt Only, Through Toggle No Charge
Part: ELE1905907

Select toggle side. Will require extended lead time when selected with Fender Lights.

1. Toggle Mount Side SelectionLeft Side Toggle
Part: PTOINDV1839Right Side Toggle
Part: PTOINDV1838**19.1b. Actuator Control through Joystick**Single Actuator Through Controller and \$2,163
Interface
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and \$3,040
Interface
Part: XRF000336. HCPCS: E2311Single Actuator Through Controller and \$3,040
Interface (AAM)

Part: XRF000068. HCPCS: E2310

1. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button

Part: CTL123218

4-Way Toggle

Part: ELEASMB5925

4-Way Button Toggle

Part: ELEASMB5926

See section 20.2 for switch mounting.

2. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

19.1c. Actuator Control with iAccess Switch**1. Actuator Control Selection**Single Actuator Control Through Flat Button \$2,547
iAccess and Controller and Interface

Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button \$2,547
iAccess and Controller and Interface

Part: XRF000338. HCPCS: E2310

Single Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000070. HCPCS: E2310

Single Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000074. HCPCS: E2310

Multiple Actuator Control Through Flat Button \$3,460
iAccess and Controller and Interface

Part: XRF000339. HCPCS: E2311

Multiple Actuator Control Through Toggle Button \$3,460
iAccess and Controller and Interface

Part: XRF000340. HCPCS: E2311

19.2. 4-Function Push Button/iAccess Mount**REQUIRED** with 4-function push button and iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$26

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26

Part: MEC2001392

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26

Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679



20. TRU BALANCE 3 ATTENDANT CONTROL

Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525

Part: XRF000078. HCPCS: E2331

Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module

Part: XRF000079

21. SEAT MOUNTED LED LIGHTS

Optional. Requires QLINES 10 Key Joystick, Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

**22. ACCESSORIES**

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU Balance 3 Cane/Crutch Holder \$105

Part: ACC152515. HCPCS: E2207

Quantum Backpack \$168

Part: ACC1905853

Grip Lap Board (20005) \$38.95

Part: ACC152992

Grip Drink Holder (20085) \$15.50

Part: ACC152993

TRU Balance 3 Hydration System (B) \$237

Part: ACC125010

TRU Balance 3 Transfer Handle (D) \$105 ea.

Part: MEC2205680. HCPCS: E0705

Not available with Heavy Duty armrests.

TRU Balance 3 Cup Holder No Charge

Part: ACC125003

For use with TRU Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest) No Charge

Part: ACCASMB2634

For use with tube-style armrests.

TRU Balance 3 Personal Item Hook \$69

Part: ACC125008

Not available with back cane option.

TRU Balance 3 Reflector Kit \$53

Part: ACC124670

TRU Balance 3 Backpack Holder (Pair) (A) \$69

Part: ACC147774

Not available with back cane option.

TRU Balance 3 Front Mounted \$69

Cell Phone/ Media Holder, X-Grip (E)

Part: ACC155835

Not available with Heavy Duty armrests.

TRU Balance 3 Rear View Mirror \$69

Part: ACC155833

Not available with Heavy Duty armrests.

TRU Balance 3 Oxygen Holder \$210

Part: ACC125006. HCPCS: E2208

Includes accessory bar, a credit will be applied PTOIN-DV3769.

TRU Balance 3 Clothing Guard (Pair) \$237

Part: ACC125007

Not available with Heavy Duty armrests.

TRU Balance 3 Glove Box (C) \$158

Part: ACC125011

Stealth Large Tablet Holder \$625.28

Part: ACC141818

Not available with Heavy Duty armrests.

Compact Vent Tray \$1,465

Part: ACC133889. HCPCS: E1029

Push Handles (pair) \$195

Part: ACC2005027

Not available with oxygen holder, back cane option, or vent tray.

Accessory Bar \$132

Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.

Large, Black, Weather Cover \$202

Part: ACCCOVR1021

Q-Logic EX/NE Programmer \$573

Part: ELEASMB6976

Privacy Bag \$420

Part: ACC20005723

Flip-Back Armrest, Side Track Accessory Bracket ... \$11

Part: MEC1809178

Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.

Personal Item Pouch \$84

Part: ACC2101017


Depth Adjustable Vent Tray \$595

Part: ACC2302500. HCPCS: E1029

Includes accessory bar, a credit will be applied PTOIN-DV3769.

Remote Stop Switch \$1,295

Part: XRF000379

**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.

Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.

Part: ACC2405496

QTY

23. BRANCH PATCHES

Optional.

Air Force Patch \$21

Part: ACC140240

Coast Guard Patch \$21

Part: ACC140239

Navy Patch \$21

Part: ACC140244

Army Patch \$21

Part: ACC140241

Marines Patch \$21

Part: ACC140243

Purple Heart Patch \$21

Part: ACC140238



24. POSITIONING COMPONENTSOptional. More options are available on the [Positioning Components Order Form](#).**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664

Part: ST-ULS4042000-Q. HCPCS: E1034

Only available with Redesigned back assemblies.

Lateral Hardware, Q-Back Right \$281

Part: ST-TWBLT-QL-R. HCPCS: E1034

Lateral Hardware, Q-Back Left \$281

Part: ST-TWBLT-QL-L. HCPCS: E1034

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664

Part: ST-ULPK404108B00-Q. HCPCS: E1034

Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664

Part: ST-ULPK404104B00-Q. HCPCS: E1034

Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332

Part: ST-ULPK404103L00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332

Part: ST-ULPK404103R00-Q. HCPCS: E1034

Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$250

Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$250

Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals

Part: PTOINDV3499

Do Not Install

Part: PTOINDV3500

24.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.

Removable Mounting Hardware
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.

Removable Mounting Hardware
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.

Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.

Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.

Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.

Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.

Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.

Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad..... \$106 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$106 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.

Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.

Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left

Part: PTOINDV3585

Install on Right

Part: PTOINDV3586

Do Not Install

Part: PTOINDV3719

24.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea.

Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154 ea.

Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154 ea.

Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154 ea.

Part: ST-IFPS2104-Q. HCPCS: E0960

24.4. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor \$97 ea.

(IFCS4113)
Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$97 ea.

(IFCS4114)
Part: ST-IFCS4114-Q. HCPCS: E0960

Positioning components section continued on following pages.

24. POSITIONING COMPONENTS CONTINUED**24.5. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

UniLink Facial Lateral, Swing Away, Right..... \$318
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left \$318
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$492
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$492
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$434
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028



Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A)\$57 ea.
 Part: INDPART2815

Facial Pad, 4".5x2", Cool Core (900M-Q)\$61 ea.
 Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea.
 Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
 Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

24.6. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028


Left, Fixed Elbow Support Bracket \$199 ea.
 TRU Balance 3 Arms (ARMS-ESFTB3-L)
 Part: POS133465

Left, UniLink, Fixed, Elbow Support Bracket..... \$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support..... \$199 ea.
 Bracket Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, Fixed Elbow Support Bracket \$199 ea.
 TRU Balance 3 Arms (ARMS-ESFTB3-R)
 Part: POS133466

Right, UniLink, Fixed, Elbow Support Bracket.....\$211
 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q



Elbow Supports

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
 Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$57 ea.
 Part: POS124067. HCPCS: K0108

24.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

24.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

25. LAP TRAY**25.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

25.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

26. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**26.1. Enhanced Display**Q-Logic 3 Enhanced Display Kit.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Night Train.....\$63 Part: DGN137729-PBK3	Pretty in Pink\$63 Part: DGN137729-PPK3	Lemon Crush.....\$63 Part: DGN137729-PYL4
Red Corvette.....\$63 Part: DGN137729-PRD1	Grey Street.....\$63 Part: DGN137729-PSL3	Midnight Blue\$63 Part: DGN137729-PBU10	Green Machine\$63 Part: DGN137729-PGN11
Back in Black.....\$63 Part: DGN137729-PBK1	Deep Purple.....\$63 Part: DGN137729-PPR5	Electric Blue\$63 Part: DGN137729-PBU9	Military Green\$63 Part: DGN137729-PGN8

**26.2. SCIM Module****Optional.**☐ SCIM Module for Q-Logic 3\$525
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.

**26.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ\$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.

**26.4. Head Array**i-Drive 4.0 Bluetooth Tri-array Series,.....\$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330Adult, Curved, i-Drive Tri-Array with\$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330i-Drive 4.0 Bluetooth Pro Series,.....\$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330Junior, Curved, i-Drive Tri-Array with\$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330**26.5. Sip and Puff Drive****1. Sip and Puff Module**Sip and Puff Module for Q-Logic 3\$1,990
Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.**Stealth Sip and Puff for Comfort Series,\$318
Left (SP625M-CL)
Part: ST-SP625M-CL-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.Stealth Sip and Puff for Stealth i-Drive,\$318
Head Array (SP625-U)
Part: ST-SP625-U-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.Stealth Sip and Puff for Comfort Series,\$318
Right (SP625M-CR)
Part: ST-SP625M-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.**3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12).....\$35
Part: ST-SP625-S12-Q**26.6. Egg Switches**Red Egg Switch (ESRED).....\$87
Part: SWTMCRO1077Green Egg Switch (ESGRN)\$87
Part: SWTMCRO1078Blue Egg Switch (ESBLU)\$87
Part: SWTMCRO1079Yellow Egg Switch (ESYEL)\$87
Part: SWTMCRO1080Black Egg Switch (ESBLK)\$87
Part: SWTMCRO1081

Alternative controls and electronics section continued on following page.

26. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

26.7. Switch Adapter

2 to 1 Switch Adapter \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

26.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
 Part: CTL168057

26.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

