



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 Order Form with TRU Balance® Flex Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



(Synergy seating shown above)

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 2.0 3S-SS \$8,533

Part: Q6 EDGE 2.0 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

2. DRIVE MOTORS**REQUIRED.**

6 mph Drive Motors Standard
Part: MOT147433

Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

 Red Corvette
Part: DGN153201-PRD1

 Night Train
Part: DGN153201-PBK3

 Deep Purple
Part: DGN153201-PPR5

 Lemon Crush
Part: DGN153201-PYL4

 Electric Blue
Part: DGN153201-PBU9

 Pretty in Pink
Part: DGN153201-PPK3

 Grey Street
Part: DGN153201-PSL3

 Green Machine
Part: DGN153201-PGN11

 Back in Black
Part: DGN153201-PBK1

 Evolution Orange
Part: DGN153201-POR7

 Midnight Blue
Part: DGN153201-PBU10

 Military Green
Part: DGN153201-PGN8

4. BASE MOUNTED LIGHTS

Optional. Select only one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

4.1a. Fender Lights

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.


Switch Fender Lights (without iLevel 12" Lift System) \$263
Part: XRF000177

4.1b. Base Mounted Light

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129
Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

5. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic 3 electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21
Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

6. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color if applicable. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires (No Accent Color) Standard
Part: PTOINDV3525 (pair)

3"x14" Flat Free Black Tires for Wheel Accents No Charge
Part: WHL1611297 x2

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

3"x14" Pneumatic Black Tires (No Accent Color) No Charge
Part: WHL150034 x2

3"x14" Pneumatic Black Tires for Wheel Accents No Charge
Part: WHL1611296 x2

4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black Accent Color \$63
Part: XRF000179

Red Corvette Accent Color \$63
Part: XRF000180

Night Train Accent Color \$63
Part: XRF000181

Grey Street Accent Color \$63
Part: XRF000182

Deep Purple Accent Color \$63
Part: XRF000184

Pretty in Pink Accent Color \$63
Part: XRF000185

Evolution Orange Accent Color \$63
Part: XRF000186

Midnight Blue Accent Color \$63
Part: XRF000187

Electric Blue Accent Color \$63
Part: XRF000188

Lemon Crush Accent Color \$63
Part: XRF000189

Green Machine Accent Color \$63
Part: XRF000190

Military Green Accent Color \$63
Part: XRF000191

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376. HCPCS: E2361

NF-22, AGM, MK Battery \$432 ea.
Part: BAT121346. HCPCS: E2361

NF-22, Gel Cell, MK Battery \$432 ea.
Part: BATGEL1008. HCPCS: E2361

8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

8.1a. Non-Expandable Controls

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key JoystickStandard

Part: CTL2205379

Controls up to 2 actuators.

QLNES 10 Key Joystick\$420

Part: CTL2205380

Controls up to 2 actuators through joystick and seat mounted, LED, full light package.



QLNES Joystick

8.1b. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730

Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$447

Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick\$840

Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights\$840

Part: CTL167725

Q-Logic 3e Joystick\$735

Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic

Part: XRF000010

Requires selection of Enhanced Display



Q-Logic 3 Joysticks

8.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder\$60

Part: ELE2103635

8.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft TouchStandard

Part: DGN124307

Also matches Night Train shrouds.

Deep Purple\$63

Part: DGN124336-PPR5

Electric Blue\$63

Part: DGN124336-PBU9

Red Corvette\$63

Part: DGN124336-PRD1

Pretty in Pink\$63

Part: DGN124336-PPK3

Lemon Crush\$63

Part: DGN124336-PYL4

Back in Black\$63

Part: DGN124336-PBK1

Evolution Orange\$63

Part: DGN124336-POR7

Green Machine\$63

Part: DGN124336-PGN11

Grey Street\$63

Part: DGN124336-PSL3

Midnight Blue\$63

Part: DGN124336-PBU10

Military Green\$63

Part: DGN124336-PGN8

9. JOYSTICK HANDLES

Optional. The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)\$90

Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G)\$90

Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C)\$90

Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Bodypoint 3" Goal Post Handle\$104

Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D)\$126

Part: ACCASMB1858. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$104

Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E)\$69

Part: ACCASMB1859. HCPCS: E2323

Q-Logic 3 Carrot Handle (A)No Charge


Part: DGN124308

Soft Ball Handle (F)\$80


Part: ACCASMB1856. HCPCS: E2323



Joystick Handles

10. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**10.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**10.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**10.1d. Specialty Joystick Mounts**


Select a bracket and side to mount it.

Retract4 Mount \$368
Part: XRF000003. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.**10.2. Joystick Mount Location**☐ Left Side Mount
Part: PTOINDV1225☐ Right Side Mount
Part: PTOINDV1226**11. TRU BALANCE FLEX SEATING SYSTEM****Required.** Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D.**11.1. TRU Balance Flex Frame Selection**TRU Balance Flex Seating System No Charge
Part: PTOINDV3853**11.2. TRU Balance Flex Seat Dimension - Width****REQUIRED.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

11.3. TRU Balance Flex Seat Dimension - Depth**REQUIRED.** Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

11.4. TRU Balance Flex Seat to Floor**REQUIRED.**16.5" at front and rear
Part: PTOINDV345917.5" at front with 2° posterior angle
Part: PTOINDV346018.25" at front with 4° posterior angle
Part: PTOINDV3805**11.5. Stealth Seat Cushion****Optional.** Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

TRU Balance Flex seating section continued on following page.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39

Part: POP1909006-S

14"W and smaller

Medium Wedge Positioning Kit.....\$72

Part: POP1909006-M

14"W and smaller

Large Wedge Positioning Kit..... \$105

Part: POP1909006-L

14"W and smaller

Small Wedge Positioning Kit\$39

Part: POP1712340-S

15"W and larger

Medium Wedge Positioning Kit.....\$72

Part: POP1712340-M

15"W and larger

Large Wedge Positioning Kit..... \$105

Part: POP1712340-L

15"W and larger

3. Cushion Top Foam**Visco Top Foam**

Reference part: XRF000380

No Top Foam

Reference part: XRF000382

FRG Top Foam

Reference part: XRF000381

4. Cushion Liner**Stealth Urethane Liner**

Reference part: XRF000383

No Liner

Reference part: XRF000385

Stealth Polyurethane Liner

Reference part: XRF000384

5. Cushion Cover**CoolCore Cover**

Reference part: XRF000386

Stealth Tek Cover

Reference part: XRF000388

Mesh Cover

Reference part: XRF000387

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer

Reference part: XRF000389

MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner

Reference Part: XRF000390

QTY

Additional Stealth Polyurethane Liner

Reference Part: XRF000391

QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover

Reference Part: XRF000392

QTY

Additional Stealth Tek Cover

Reference Part: XRF000394

QTY

Additional Mesh Cover

Reference Part: XRF000393

QTY

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner					
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	-	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

TRU Balance Flex seating section continued on following page.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**11.6. TRU Balance Flex Back Cane Selection****REQUIRED.**

Angle Adjustable Back Cane..... No Charge
Part: MEC140266

Flip Forward Back Cane \$184
Part: XRF000274
Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane..... \$205
Part: MEC1804736

2" Endomorph (1" per side) Back Cane..... \$205
Part: XRF000270

4" Endomorph (2" per side) Back Cane..... \$205
Part: XRF000272

Manual Recline Back \$788
Part: FRMASMB12910. HCPCS: E1226
Select Synergy Back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

11.7. TRU Balance Flex Back Cane Height Selection**REQUIRED.**

14" High Canes Part: FRMASMB15832	17" High Canes Part: FRMASMB15835	20" High Canes Part: FRMASMB15838	23" High Canes Part: FRMASMB15841
15" High Canes Part: FRMASMB15833	18" High Canes Part: FRMASMB15836	21" High Canes Part: FRMASMB15839	24" High Canes Part: FRMASMB15842
16" High Canes Part: FRMASMB15834	19" High Canes Part: FRMASMB15837	22" High Canes Part: FRMASMB15840	25" High Canes Part: FRMASMB15843

11.8. TRU Balance Flex Back Cushion Options

REQUIRED. Select one type of back. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1a. Synergy Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

Synergy Back..... \$289
Part: XRF000363

1b. Sport Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

Sport Back Cushion, 21" Wide or Less..... \$394
Part: XRF000361

Sport Back Cushion, 22" Wide or More..... \$604
Part: XRF000362

1c. TRU-Comfort Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width..... \$744
Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621

1d. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622 Part: ALBR16SA4. HCPCS: E2613	Active Contour, 18" Wide x 23" or 24" Overall Height \$622 Part: ALBR18XA4. HCPCS: E2613
Active Contour, 16" Wide x 21" or 22" Overall Height \$622 Part: ALBR16TA4. HCPCS: E2613	Active Contour, 20" Wide x 19" or 20" Overall Height \$622 Part: ALBR20SA4. HCPCS: E2613
Active Contour, 16" Wide x 23" or 24" Overall Height \$622 Part: ALBR16XA4. HCPCS: E2613	Active Contour, 20" Wide x 21" or 22" Overall Height \$622 Part: ALBR20TA4. HCPCS: E2613
Active Contour, 18" Wide x 19" or 20" Overall Height \$622 Part: ALBR18SA4. HCPCS: E2613	Active Contour, 20" Wide x 23" or 24" Overall Height \$622 Part: ALBR20XA4. HCPCS: E2613
Active Contour, 18" Wide x 21" or 22" Overall Height \$622 Part: ALBR18TA4. HCPCS: E2613	

1e. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height \$735 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height \$735 Part: ALBR20XD4. HCPCS: E2615
Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735 Part: ALBR18SD4. HCPCS: E2615	

1f. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822 Part: ALBR18S64. HCPCS: E2615	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822 Part: ALBR18X64. HCPCS: E2615
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822 Part: ALBR18T64. HCPCS: E2615	

TRU Balance Flex seating section continued on following page.

11. TRU BALANCE FLEX SEATING SYSTEM CONTINUED**1h. Omit Back Selection**

Omit Back..... No Charge
Part: PTOINDV1287

11.9. Back Width Package

REQUIRED with Sport and TRU-Comfort Back selection.

16"W-20"W Package 18"W-22"W Package
Part: XRF000038 Part: XRF000039

11.10. Overall Back Height

16" Overall Back Height Part: SETHGHT1029	23" Overall Back Height Part: SETHGHT1036
17" Overall Back Height Part: SETHGHT1030	24" Overall Back Height Part: SETHGHT1037
18" Overall Back Height Part: SETHGHT1031	25" Overall Back Height Part: SETHGHT1038
19" Overall Back Height Part: SETHGHT1032	26" Height Part: SETHGHT1039
20" Overall Back Height Part: SETHGHT1033	27" Height Part: SETHGHT1040
21" Overall Back Height Part: SETHGHT1034	28" Height Part: SETHGHT1041
22" Overall Back Height Part: SETHGHT1035	29" Height Part: SETHGHT1042

		Back Width									
		14"	15"	16"	17"	18"	19"	20"	21"	22"	
Overall Height	16"	P	T	P	T	P	T	P	T	P	T
	17"	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT
	18"	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT
	19"	P	T	SPT	SPT	P	T	SPT	SPT	SPT	SPT
	20"	P	T	SPT	SPT	P	T	SPT	SPT	SPT	SPT
	21"	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT
	22"	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT
	23"	P	T	SPT	SPT	SPT	SPT	SPT	SPT	SPT	SPT
	24"	P	T	P	T	SPT	S	P	SPT	SPT	S
	25"	P	T	P	T	SPT	S	P	SPT	SPT	S
	26"	P	T	P	T	P	T	P	T	T	P
	27"	P	T	P	T	P	T	P	T	T	P
	28"	P	T	P	T	P	T	P	T	T	P
	29"	P	T	P	T	P	T	P	T	T	P

☐ Synergy Back ☐ Sport Back ☐ TRU-Comfort Back

TRU Balance Flex Seat Back Availability Chart

11.11. TRU Balance Flex Armrests

REQUIRED. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. TRU Balance Flex Seat, Drop-In Armrests (T-Arm)

☐ Full Length, Drop-In Armrests..... No Charge L R Left part: ARM1803253. Right part: ARM156742	☐ Desk Length, Drop-In Armrests No Charge L R Left part: ARM1803254. Right part: ARM156743
---	--

1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ Full Length, Flip Back Armrests..... \$63 L R Left part: XRF000371. Right part: XRF000372	☐ Desk Length, Flip Back Armrests \$63 L R Left part: XRF000373. Right part: XRF000374
---	--

1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ Full Length, Two-Post Armrests.....\$158 ea. L R Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	☐ Desk Length, Two-Post Armrests\$158 ea. L R Left part: XRF000053. Right part: XRF000054. HCPCS: E0973
--	---

1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ Full Length, Single Post Armrests \$79 ea. L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973	☐ Desk Length, Single Post Armrests \$79 ea. L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973
--	--

11.12. TRU Balance Flex Armpads

REQUIRED with armrest selection above, except with drop-in armrests. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

1a. Straight Armpads

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134275
Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right..... No Charge Part: ARM134254

1b. Waterfall Armpads

Waterfall, Full Length Left and Right..... No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	Waterfall, Desk Length Left and Right No Charge Part: ARM134282

1c. Falcon Gel Armpads

Falcon Gel, Full Length Left and Right.....\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right.....\$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437

12. HEADRESTS**Optional.** Select a headrest pad and bracket.**12.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955**2. Headrest Brackets**Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033**12.2. Flip-down Headrest Bracket**Stealth Flip-down Headrest Bracket (FDM350)..... \$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

12.3. Q Cam Lock**Optional.**Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

13. LAP BELTS**REQUIRED.** Select one lap belt below.38" Lap Belt
Part: ACCBELT100650" Lap Belt
Part: ACCBELT100060" Lap Belt
Part: ACCBELT100270" Lap Belt
Part: ACCBELT100380" Lap Belt
Part: ACCBELT100490" Lap Belt
Part: ACCBELT1001**14. LEGRESTS****Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**14.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077Omit base mounted foot platform No Charge
Part: PTOINDV3674Bell Shaped, Base Mounted Foot Platform \$247
Part: RIGASMB7120041**14.1b. Basic Center Mount Foot Platform**

Select legrest length in following section. See legrest length matrix for compatibility.

1. Basic Center Mount Foot Platform AssemblyBasic Center Mount..... No Charge
Part: XRF000323**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**Short Extension, 5"x9" Footplates\$180
Part: XRF000312. HCPCS: K0040Long Extension, 5"x9" Footplates\$180
Part: XRF000314. HCPCS: K0040Medium Extension, 5"x9" Footplates\$180
Part: XRF000313. HCPCS: K0040**3. Basic Center Mount Center Mount Calf Support**

Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997

Only available with Long Length Extension.

4. Basic Center Mount Additional AccessoriesBasic Center Mount Footplate Connector \$16
Part: RIG2104240

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

Legrest section continued on following pages.

14. LEGRESTS CONTINUED**14.1c. Center Mount Foot Platform**

Select legrest length in following section.

1. Foot Platform Assembly

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

Quick Adjustable Center Mount Foot Platform \$263
 Part: XRF000311. HCPCS: K0040

2. Lower Extension and Angle Adjustable Footplates

Pediatric Extension, 5"x6" Footplates
 Part: RIG152577

Medium Extension, 5"x8" Footplates
 Part: RIG2404016

Long Extension, 5"x8" Footplates
 Part: RIG2404017

Pediatric Extension, 5"x8" Footplates
 Part: RIG2404015

Medium Extension, 5"x11" Footplates
 Part: RIG2404018

Long Extension, 5"x11" Footplates
 Part: RIG2404019

Medium Extension, 5"x6" Footplates
 Part: RIG150206

Long Extension, 5"x6" Footplates
 Part: RIG151505

3. Center Mount Calf Support

Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H \$158
 Part: XRF000156

Small Calf Panel, 18"Wx8"H \$184
 Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$158
 Part: XRF000157

Large Calf Panel, 24"Wx8"H \$237
 Part: RIG123825
 Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
 Part: RIG152519
 Only available with the 5x6 footplates.

Outside Edge \$69
 Part: RIG2404921
 Available with 5"x8" footplates.

Width Extension for 8" \$69
 Part: RIG2404919

Front and Rear Cap \$16
 Part: RIG2404923
 Available with 5"x8" footplates.

Width Extension for 11" \$69
 Part: RIG2404920

5. Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
 Part: ST-CSCS0775-Q

Footplate Riser \$75
 Part: XRF000348

Large Foam Calf Pad Cover (pair) \$83
 Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Small Foam Panel Cover (pair) \$176
 Part: ST-CSCP0920-Q

Footplate Riser Calf Panel \$50
 Part: FRMASMB7110011

Large Foam Panel Cover (pair) \$282
 Part: ST-CSCP0927-Q

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$16
 Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
 Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
 Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270
 Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
 Part: ST-AFPG5011E-Q. HCPCS: K0108

14.1d. Swing-Away Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
 Pair part: XRF000082. Left: XRF000083. Right: XRF000084
 Not available with 11"-13" extensions.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
 Pair part: XRF000085. Left: XRF000086. Right: XRF000087
 For use only with 11"-13" lower extensions.

2. Lower Extension

☐ ☐ ☐ 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ ☐ ☐ 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ ☐ ☐ 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

Legrest section continued on following pages.

14. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

14.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

P L R 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

P L R 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

P L R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

P L R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R 7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18".

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

14.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$19 ea.
 Part: XRF000100

14.3 Set Legrest Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

Legrest section continued on following pages.

14. LEGRESTS CONTINUED**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

15. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

16. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Flex Seat Cup Holder No Charge

Part: ACCASMB2634

Not available on TRU Balance Flex seating without two-post armrests.

Oxygen Holder for Flex Seat \$210

Part: ACC125006. HCPCS: E2208

Not available with Manual Recline.

Cane/Crutch Holder \$105

Flex part: ACC152515. HCPCS: E2207

Flag Holder \$35

Flex part: ACCASMB2415

Saddle Bag \$40

Part: DWR1020D010

Stealth Cell Phone/Media Device Holder \$215

Part: ACC7120006

Stealth Large Tablet Holder \$625.28

Part: ACC7120005

Quantum Backpack \$168

Part: ACC1905853

Grip Lap Board (20005)..... \$38.95

Part: ACC152992

Grip Drink Holder (20085) \$15.50

Part: ACC152993

Accessory Bar \$132

Part: XRF000128

Large, Black, Weather Cover \$202

Part: ACCCOVR1021

Privacy Bag \$42

Part: ACC20005723

Transfer Handle \$105 ea.

Part: MEC2205680. HCPCS: E0705

Q-Logic EX/NE Programmer \$573

Part: ELEASMB6976

Q-Logic PC Programming Disk with Cable \$445

Part: ELEASMB5216

Personal Item Pouch \$84

Part: ACC2101017

Remote Stop Switch \$1,295

Part: XRF000379

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.

Part: ACC2407046

QTY

Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.

Part: ACC2405496

QTY

17. BRANCH PATCHES**Optional.**Air Force Patch \$21
Part: ACC140240Coast Guard Patch \$21
Part: ACC140239Navy Patch \$21
Part: ACC140244Army Patch \$21
Part: ACC140241Marines Patch \$21
Part: ACC140243Purple Heart Patch \$21
Part: ACC140238

Military Patches

18. POSITIONING COMPONENTSOptional. More options are available on the [Positioning Components Order Form](#).**18.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBSTL) \$562
Part: INDPART2896. HCPCS: E1034Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034
Only available with back cane and omit back configurations.**2. Pads**Curved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6) \$250
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**4. Lateral Extensions****REQUIRED with Synergy Back and thoracic lateral selections.**Lateral Extension for Synergy Back (pair) \$29 ea.
Part: XRF000150**18.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP5-Q. HCPCS: E1028/E1034^Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^UniLink, 7" Height Adjustable, \$276 ea.
Removable Mounting Hardware
QTY Part: ST-UTTP7-Q. HCPCS: E1028/E1034^Flip Down with RIGHT Full Surface Contact \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Right \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 12" Pad \$106 ea.
QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 8" Pad \$106 ea.
QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^**3. UniLink Pelvic/Thigh Additional Links**UniLink, 1.5" Link \$81 ea.
QTY Part: ST-UIL001150A000-QUniLink, 2" Link \$81 ea.
QTY Part: ST-UIL001200A000-Q**4. UniLink Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**18.4. Stealth Positioning Chest Harness**Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103) \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960**18.5. Stealth Positioning Chest Straps**Chest Strap, Medium w/ Belt Anchor (IFCS4113) ... \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor (IFCS4114) \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960**18.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

Positioning components section continued on following pages.

18. POSITIONING COMPONENTS CONTINUED**2. Facial Lateral Pads and Accessories**


QTY Egg Switch Mount (ESM) (A)\$57 ea.
Part: INDPART2815

QTY Facial Pad, 4.5"x2", Cool Core (900M-Q)\$61 ea.
Part: ST-900M-Q. HCPCS: K0108

QTY Facial Pad, 6"x2", Cool Core (920M-Q)\$61 ea.
Part: ST-920M-Q. HCPCS: K0108

QTY Swing-Away Joystick Mount,\$276 ea.
18" Rod (C)

Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest

**Headrest Accessories****18.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket\$199 ea.
Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket\$199 ea.
Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

QTY 3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea.
Part: POS124065. HCPCS: K0108

QTY 4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$57 ea.
Part: POS124067. HCPCS: K0108

18.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H\$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H\$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H\$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

**Stealth Foot Boxes****18.9. Stealth Positioning Shoe Holder**

QTY Small Shoe Holder,\$87 ea.
3/4" Webbing (IFSH4211)
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

QTY Large Shoe Holder,\$87 ea.
3/4" Webbing (IFSH4213)
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

QTY Medium Shoe Holder,\$87 ea.
3/4" Webbing (IFSH4212)
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

**Stealth Shoe Holders**

19. LAP TRAY

Optional.

19.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q, HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q, HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q, HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware.....\$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	ML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.