



Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 2.0 X Order Form with Contoured and Synergy Seating

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a power base. This power base has a 300 lb. weight capacity with standard 6 mph programmable motors, 8 amp off board charger, fender lights and USB mobile device charger. Dimensions of base without legrests: 24"W x 35.5"L

Quantum Q6 Edge 2.0 X 3S-SS\$10,195

Part: Q6EDGE 2.0X 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge 2.0 X 3S-C\$10,195

Part: Q6EDGE 2.0X 3S-C. HCPCS: K0849

Contoured Seat (Captain) Standard.

Quantum Q6 Edge 2.0 X 3SP-SS\$10,195

Part: Q6EDGE 2.0X 3SP-SS. HCPCS: K0856

Single power option for Power Elevating Seat for Synergy Seat.

2. BASE COLOR**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.**2.1. Match Colors**

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

2.2. Standard Shrouds

 Red Corvette Part: DGN153201-PRD1	 White Satin Part: DGN153201-PWT2	 Grey Street Part: DGN153201-PSL3	 Military Green Part: DGN153201-PGN8
 Electric Blue Part: DGN153201-PBU9	 Pretty in Pink Part: DGN153201-PPK3	 Midnight Blue Part: DGN153201-PBU10	To Be Determined Part: N/A
 Back in Black Part: DGN153201-PBK1	 Evolution Orange Part: DGN1801101-POR7	 Lemon Crush Part: DGN1801101-PYL4	
 Night Train Part: DGN153201-PBK3	 Deep Purple Part: DGN153201-PPR5	 Green Machine Part: DGN1801101-PGN11	

3. DRIVE MOTORS**REQUIRED.**

6 mph Drive Motors Standard Part: MOT147433	Accu-Trac Motor Technology \$1,295 Part: MOT147881. HCPCS: K0108 Delivers enhanced traction and tracking, especially with switched drive controls.
--	--

4. BASE MOUNTED FOOT PLATFORM**REQUIRED.** May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Bell Shaped, Base Mounted Foot Platform \$235 Part: RIGASMB7120041	

5. BATTERIES**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor Battery.....\$310 ea. Part: BAT1704613. HCPCS: E2361	NF-22, AGM, MK Battery..... \$395 ea. Part: BAT121346. HCPCS: E2361
No batteries. Provider will assume responsibility of batteries. Part: Configuration selection	NF-22, Gel Cell, MK Battery \$395 ea. Part: BATGEL1008. HCPCS: E2361

6. DRIVE WHEELS**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**6.1. Drive Wheels**

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair). Wheel w/ accent part: WHL1611297 x qty 2	3"x14" Pneumatic Black Tires No Charge Part: WHL150034 x2. Wheel w/ accent part: WHL1611296 x qty 2
4"x14" Knobby, Flat Free Black Tires..... No Charge Part: WHL147601 x qty 2	4"x14" Knobby, Pneumatic Black Tires..... No Charge Part: WHL147599 x qty 2

6.2. Drive Wheel Accent Color**Optional.** Sold in Pairs (parts below will be quantity of 2). Not available with knobby tires.

Match to Base Color Part: Configuration selection, see specific part below		Non-matching Color Part: Configuration selection, see specific part below	
Back in Black..... \$24.98 ea. Part: DGN164245-PBK1 x2	Evolution Orange..... \$24.98 ea. Part: DGN164245-POR7 x2	Lemon Crush..... \$24.98 ea. Part: DGN164245-PYL4 x2	
Red Corvette..... \$24.98 ea. Part: DGN164245-PRD1 x2	Night Train..... \$24.98 ea. Part: DGN164245-PBK3 x2	Green Machine \$24.98 ea. Part: DGN164245-PGN11 x2	
Electric Blue \$24.98 ea. Part: DGN164245-PBU9 x2	Deep Purple \$24.98 ea. Part: DGN164245-PPR5 x2	Military Green \$24.98 ea. Part: DGN164245-PGN8 x2	
White Satin..... \$24.98 ea. Part: DGN164245-PWH2 x2	Grey Street..... \$24.98 ea. Part: DGN164245-PSL3 x2	To Be Determined Follow-up required.	
Pretty in Pink..... \$24.98 ea. Part: DGN164245-PPK3 x2	Midnight Blue \$24.98 ea. Part: DGN164245-PBU10 x2		

7. SEAT TOWERS/POWER ADJUSTABLE SEAT HEIGHT**REQUIRED.** Select one option below.**7.1a. Power Base Seat Towers**

Standard Seat Towers Standard Part: Configuration selection, no part number	Extended Seat Towers No Charge Part: MEC120301 Increases seat to floor height.
--	--

7.1b. Power Adjustable Seat Height

The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

Power Elevating Seat with 6" Stroke \$1,895 Contoured Seat Part: FRMASMB16309. HCPCS: E2300 Synergy Seat Part: FRMASMB16311. HCPCS: E2300	Power Elevating Seat with 7" Stroke \$1,895 Contoured Seat Part: FRMASMB16310. HCPCS: E2300 Synergy Seat Part: FRMASMB16312. HCPCS: E2300
---	---

8. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

8.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller ModuleStandard
Part: See XRef

Q-Logic 3 Expandable Controller Harness.....Standard
Part: HAR144005

2. Q-Logic 3 Joystick

Q-Logic 3 EX JoystickStandard
Part: CTL167724

No Joystick, Setup for Q-Logic
Part: Configuration selection



Q-Logic 3 EX Joystick

8.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic joystick selections above. Reference colors shown in section 2.

Black, Soft TouchStandard
Part: DGN124307

Non-matching Armor\$50
Part: PTOINDV3422

Match Armor to Base Color\$50
Part: PTOINDV3421

Select color below

Red Corvette
Part: DGN124336-PRD1

Night Train
See Black, Soft Touch

Evolution Orange
Part: DGN124336-POR7

Midnight Blue
Part: DGN124336-PBU10

Military Green
Part: DGN124336-PGN8

Electric Blue
Part: DGN124336-PBU9

White Satin
Part: DGN124336-PWH2

Deep Purple
Part: DGN124336-PPR5

Lemon Crush
Part: DGN124336-PYL4

Back in Black
Part: DGN124336-PBK1

Pretty in Pink
Part: DGN124336-PPK3

Grey Street
Part: DGN124336-PSL3

Green Machine
Part: DGN124336-PGN11

9. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in the image on the previous page. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B)\$85
Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G)\$85
Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C)\$85
Part: ACCASMB1861. HCPCS: E2323. (E2324 when chin control)

Bodypoint 3" Goal Post Handle\$99
Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D)\$120
Part: ACCASMB1858. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$99
Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E)\$65
Part: ACCASMB1859. HCPCS: E2323

Q-Logic 3 Carrot Handle (A)No Charge
Part: DGN124308

Soft Ball Handle (F)\$65
Part: ACCASMB1856. HCPCS: E2323



Joystick Handles

10. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.


Omit Joystick Mounting Bracket

10.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, BracketStandard
Q-Logic Part: FRMASMB10617

☐ Left Side Mount
Part: PTOINDV1225

☐ Right Side Mount
Part: PTOINDV1226



Swing Away Joystick Offset Mount (Right)

10.1b. Swing-Away (Retractable) Joystick Bracket


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

☐ Swing-away, Inline, Left Side Mount\$245
Part: See XRef. HCPCS: E1028

☐ Swing-away, 1" Offset, Left Side Mount\$245
Part: See XRef. HCPCS: E1028

☐ Swing-away, Inline, Right Side Mount\$245
Part: See XRef. HCPCS: E1028

☐ Swing-away, 1" Offset, Right Side Mount\$245
Part: See XRef. HCPCS: E1028



Swing Away Joystick Inline Mount (Right)

10.1c. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount\$350
Part: MEC138860. HCPCS: E1028


Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-axis Mount\$100
Q-Logic Part: ACC131216. HCPCS: E1028

SAJ Part: ACC142740. HCPCS: E1028
Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket\$294.25
Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028

SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted. Not available with contoured seating.



Retract4



Multi-Axis Joystick Mount

Specialty Mount Side

☐ Left Side Mount
Part: PTOINDV1225

☐ Right Side Mount
Part: PTOINDV1226

Account Number: _____ Date: _____ Marked for: _____

11. CONTOURED SEATING

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

11.1a. Captain Style Seating

1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline

Comfort seat, full length armrests must be selected with a depth adjustable high back. Includes 60" lap belt.

16"Wx16"-18"D No Charge	18"Wx18"-20"D \$150	20"Wx18"-20"D \$150
Part: SET131843	Part: SET131844	Part: SET131845

1b. Depth Adjustable High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt No Charge	22"Wx20"-22"D with 70" Lap Belt No Charge
Part: SETASMB6895	Part: SETASMB6897
20"Wx18"-20"D with 60" Lap Belt No Charge	
Part: SETASMB6896	

1c. High Back with Limited Recline

16"Wx16"D No Charge	16"Wx22"D No Charge	20"Wx16"D No Charge
Part: SETASMB2770	Part: SETASMB4887	Part: SETASMB2780
16"Wx18"D No Charge	18"Wx16"D No Charge	20"Wx22"D No Charge
Part: SETASMB2771	Part: SETASMB2775	Part: SETASMB4889
16"Wx20"D No Charge	18"Wx22"D No Charge	
Part: SETASMB2772	Part: SETASMB4888	

1d. Medium Back

16"Wx16" No Charge	16"Wx18" No Charge	18"Wx16" No Charge	18"Wx18" No Charge
Part: SETASMB2488	Part: SETASMB2560	Part: SETASMB2159	Part: SETASMB2551

11.1b. Contoured Seat with Solid Seat Pan

1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline and Solid Seat Pan

Comfort Seat, full length armrests must be selected with a depth adjustable high back.

16"Wx16"-18"D with 60" Lap Belt \$150	18"Wx18"-20"D with 60" Lap Belt \$150
Part: SET145814	Part: SET145815

1b. High Back with Limited Recline and Solid Seat Pan

16"Wx16"D No Charge	18"Wx16"D No Charge	20"Wx16"D No Charge
Part: SETASMB3178	Part: SETASMB3181	Part: SETASMB3183
16"Wx18"D No Charge	18"Wx18"D No Charge	20"Wx18"D No Charge
Part: SETASMB3179	Part: SETASMB2848	Part: SETASMB3184
16"Wx20"D No Charge	18"Wx20"D No Charge	20"Wx20"D No Charge
Part: SETASMB3180	Part: SETASMB3182	Part: SETASMB3185
16"Wx22"D No Charge	18"Wx22"D No Charge	20"Wx22"D No Charge
Part: SETASMB5006	Part: SETASMB5007	Part: SETASMB5008

1c. Medium Back with Solid Seat Pan

16"Wx16" No Charge	16"Wx18" No Charge	18"Wx16" No Charge	18"Wx18" No Charge
Part: SETASMB3112	Part: SETASMB3113	Part: SETASMB2492	Part: SETASMB3114

11.2. Occupied Transit, WC-19

WC-19 Occupied Transit Option for Contoured Seat \$350
Part: KIT122609

Only available with Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

11.3. Seat to Floor Height

Required with any contoured seating selection.

Lowest Seat to Floor Part: PTOINDV3558	Middle Seat to Floor Part: PTOINDV3559 Not available with Comfort high back seat.	Highest Seat to Floor Part: PTOINDV3560
---	---	--

11.4. Set Depth Selection for Comfort Seat

Required with Comfort Seat selection above.

14" Depth Part: SETDPTH1004	16" Depth Part: SETDPTH1006	18" Depth Part: SETDPTH1008	20" Depth Part: SETDPTH1010
--------------------------------	--------------------------------	--------------------------------	--------------------------------

11.5. Back Height Selection for Comfort Seat

Required with Comfort Seat selection above.

Set Back Height to Lowest Position Part: SETHGHT1100	Set Back Height to Highest Position Part: SETHGHT1101
---	--

11.6. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

Contoured seating section continued on following page.

11. CONTOURED SEATING CONTINUED**11.7. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests****REQUIRED with High Back Seat.** Select a left and right armrest or to omit them.

Full Length, Left Armrest.....\$85 Part: FRMASMB8614. HCPCS: E0973	Comfort Seat, Right, Full Length, Armrest.....\$85 Part: FRMASMB16315. HCPCS: E0973
Desk Length, Left Armrest.....\$85 Part: FRMASMB8612. HCPCS: E0973	Comfort Seat, Left, Full Length, Armrest.....\$85 Part: FRMASMB16316. HCPCS: E0973
Full Length, Right Armrest.....\$85 Part: FRMASMB8615. HCPCS: E0973	No Armrest.....Standard Part: N/A
Desk Length, Right Armrest.....\$85 Part: FRMASMB8613. HCPCS: E0973	

12. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.**12.1. Synergy Frame Selection**

Synergy Seating Frame.....\$434 Part: PTOINDV3746 Standard with K0848

12.2. Synergy Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

☐ WC-19 Occupied Transit Option.....\$275 Part: KITASMB1830

12.3. Seat Dimension - Width**REQUIRED with Synergy seating.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.4. Seat Dimension - Width Package**REQUIRED with Synergy seating.** Maximum seat width is 22".

10"-14" Width Part: FRMASMB18011	14"-18" Width Part: FRMASMB18012	16"-20" Width Part: FRMASMB18013	20"-24" Width Part: FRMASMB18014 (Maximum Width is 22")
-------------------------------------	-------------------------------------	-------------------------------------	---

12.5. Seat Dimension - Depth**REQUIRED with Synergy seating.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.6. Seat Dimension - Depth Package**REQUIRED with Synergy seating.** Maximum seat depth is 22".

10"-14" Depth Package Part: FRMASMB6981	14"-18" Depth Package Part: FRMASMB6982	15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package Part: FRMASMB7733 (Maximum Depth is 22")
--	--	--	--

12.7. Synergy Seat to Floor**REQUIRED with Synergy seating.**

15.875" Part: PTOINDV3652	16.875" Part: PTOINDV3653	17.875" Part: PTOINDV3654
------------------------------	------------------------------	------------------------------

12.8. Synergy Back Options**REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Cane Selection**

Synergy Back Cane.....No Charge Part: FRMASMB15844 Select a back cane height.	2" Endomorph (1" per side) Back Cane.....\$195 Part: FRMASMB12115 Select a cane height.
Flip Forward Back Cane.....\$150 Part: FRMASMB15846 Select a back cane height. Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.	4" Endomorph (2" per side) Back Cane.....\$195 Part: FRMASMB12139 Select a cane height.
2" Mesomorph (1" per side) Back Cane.....\$195 Part: FRMASMB12067 Select a cane height.	Manual Recline Back.....\$750 Part: FRMASMB12910. HCPCS: E1226 Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, Shape Back, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.
4" Mesomorph (2" per side) Back Cane.....\$195 Part: FRMASMB12091 Select a cane height.	

Synergy seating section continued on following page.

12. SYNERGY SEATING SYSTEM CONTINUED**2. Synergy Back Cane Height**

Omit this selection if a TRU-Comfort back is desired. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes
Part: FRMASMB15832

15" High Canes
Part: FRMASMB15833

16" High Canes
Part: FRMASMB15834

17" High Canes
Part: FRMASMB15835

18" High Canes
Part: FRMASMB15836

19" High Canes
Part: FRMASMB15837

20" High Canes
Part: FRMASMB15838

21" High Canes
Part: FRMASMB15839

22" High Canes
Part: FRMASMB15840

23" High Canes
Part: FRMASMB15841

24" High Canes
Part: FRMASMB15842

25" High Canes
Part: FRMASMB15843

3. Back Selection

Synergy Back..... \$275
Part: See XRef

TRU-Comfort, 20" Width or less..... \$675
Part: See XRef. HCPCS: E2620

TRU-Comfort, 22" Width or greater..... \$775
Part: See XRef. HCPCS: E2621

Sling Back..... \$98
Part: See XRef

Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back..... No Charge
Part: PTOINDV1287

4. Synergy Overall Back Height

Parts: See XRef, configuration selection with back selection.

12" High 15" High 18" High 21" High 24" High

13" High 16" High 19" High 22" High 25" High

14" High 17" High 20" High 23" High

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	-	-	-	-	-	-	-	-
	16"	S	S	-	-	-	-	-	-	-	-
	17"	-	S	S	S	S	S	S	S	S	S
	18"	-	S	S	S	S	S	S	S	S	S
	19"	-	S	S	S	-	S	S	S	S	S
	20"	-	S	S	S TC	-	S TC	S	S TC	S	S
	21"	-	S	S	S TC	S	S TC	S	S TC	S	S
	22"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	23"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	24"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
	25"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC

S Synergy Back S TC Synergy or TRU-Comfort Back

Synergy and TRU-Comfort Back Availability Chart

12.9. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... Standard
Left Part: FRMASMB7038. Right part: FRMASMB7039

☐ ☐ Desk Length, Flip Back Armrests Standard
Left part: FRMASMB7040. Right part: FRMASMB7041

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$150 ea.
Left Part: FRMASMB7062. Right part: FRMASMB7063.
HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests\$150 ea.
Left part: FRMASMB7060. Right part: FRMASMB7061.
HCPCS: E0973

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$75 ea.
Left Part: FRMASMB9982. Right part: FRMASMB9983.
HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$75 ea.
Left part: FRMASMB9984. Right part: FRMASMB9985
HCPCS: E0973

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$75 ea.
Left Part: FRMASMB9986. Right part: FRMASMB9987
HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H..... \$75 ea.
Left part: FRMASMB9988. Right part: FRMASMB9989
HCPCS: E0973

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

☐ ☐ Pediatric, Removable Armrests \$75 ea.
Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests \$97.50 ea.
Left Part: FRMASMB16188. Right part: FRMASMB16189

☐ ☐ Desk Length, Cantilever Armrests..... \$97.50 ea.
Left part: FRMASMB16190. Right part: FRMASMB16191

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$75 ea.
Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS:E0973

☐ ☐ Desk Length, Single Post Armrests \$75 ea.
Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

Synergy seating section continued on following page.

12. SYNERGY SEATING SYSTEM CONTINUED**12.10. Synergy Armpads****REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard
 Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
 Desk/Full part: UPSASMB1023

Waterfall Armpads..... No Charge
 Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
 Desk/Full part: UPSASMB1026

No Armpad Selection No Charge
 Part: PTOINDV1373

12.11. Joystick Bracket Receiver**REQUIRED.**

Standard Fixed Receiver Standard
 Stealth Height Adjustable Receiver \$214
 Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

**13. HEADRESTS****Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****13.1. Headrest Pads**

Stealth 8" Comfort Plus (CP450)..... \$151.94
 Part: POS148484. HCPCS: E0955


Stealth 10" Comfort Plus (CP250) (A) \$151.94
 Part: SETHEAD1057. HCPCS: E0955

Stealth 14" Comfort Plus (CP150) \$151.94
 Part: INDPART2327. HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) \$151.94
 (CPN350)
 Part: SETHEAD1069. HCPCS: E0955

Stealth Small Contoured (9.75"Wx3"H) \$64.20
 (STL-OTSM250)
 Part: SETHEAD1052. HCPCS: E0955

Stealth Large Contoured (11"Wx6"H) (B) \$80.25
 (STL-OTLG150)
 Part: SETHEAD1053. HCPCS: E0955

**13.2. Headrest Brackets**

Stealth's Worlds Best Removable \$220.50
 Headrest Bracket (TWB485-1) (C)
 Part: FRMASMB15999. HCPCS: E1028
 Pediatric, Removable, Offset (HMO475P-17) \$160.50
 Part: HDWASMB1166. HCPCS: E1028

Stealth Multi-Axis Removable (D)..... \$160.50
 Headrest Bracket (HMO475P)
 Part: HDWASMB1084. HCPCS: E1028

13.3. Flip-down Headrest Bracket

Stealth Flip-down Headrest Bracket (FDM350)..... \$224.70
 Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent, headrest post modification necessary to allow the headrest bracket to flip.

14. LAP BELTS**REQUIRED except with select contoured seats with included lap belt.** Select one lap belt below.**14.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

14.1b. Optional Lap Belts

Not available with contoured seating.

i-Fit 60" Single Pull, Push Button, Medium \$86
 (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q

i-Fit 60" Single Pull, Push Button, Large \$86
 Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q

i-Fit 60" Single Pull, Push Button, Extra Large \$86
 Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Required with ISO 7176-19 Compliant Transit Option.
 Part: ACCBELT1016

15. SYNERGY POWER OPTION ELECTRONICS**REQUIRED for Synergy Seat with Power Legrests or Elevating Seat.** Select the desired control type for the seating system.

One Actuator Control through through \$1,895
 Q-Logic 3 Electronics
 Part: See XRef. HCPCS for Articulating foot platform: E2310

Multiple Actuator Control through Q-Logic 3 \$2,895
 Electronics
 Part: See XRef. HCPCS for power legrests: E2311

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured seating.**

1. Foot Platform Assembly

Power Articulating Foot Platform \$2,995
Part: PTOINDV3754. HCPCS: E1012+K0040

Manually Adjustable Center Mount No Charge
Foot Platform
Part: See XRef. HCPCS: K0040

**2. Lower Extension and Angle Adjustable Footplates**

8"-12" Extension, 5"x6" Footplates

Part: RIG152577

Not available with power articulating platform.

8"-12" Extension, 5"x8" Footplates

Part: RIG152578

Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates

Part: RIG150206

10"-16" Extension, 5"x8" Footplates

Part: RIG150203

10"-16" Extension, 5"x11" Footplates

Part: RIG150212

13"-19" Extension, 5"x6" Footplates

Part: RIG151505

13"-19" Extension, 5"x8" Footplates

Part: RIG151506

13"-19" Extension, 5"x11" Footplates

Part: RIG151508

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	! May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
13"-19" Length Extension	!	!	!	

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H \$150*
Part: RIGASMB7120018

Small Calf Panel 18"Wx8"H \$150
Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H \$150*
Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H \$200
Part: RIG123824

Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$50
Part: RIG152517

4" Outside Edge (pair) \$50
Part: RIG152510

Only available with 6"D footplates and 8"D footplate extension.

Width Extension for 11"D Footplates \$50
Part: RIG152518

2" Bell-Shaped Length Extension \$50
Part: RIG152519

Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

6" Outside Edge (pair) \$50
Part: RIG152516

Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150
Small part: RIG133717. Large part: RIG133719

Wheel Bumpers (pair) \$12.50 ea.
Part: WHLASMB1578 x2

Reversed Dartex Panel Cover (pair) \$150
Small part: RIG134673. Large part: RIG134744

Footplate Connector \$15
Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$205
Part: ST-FPG5008T-Q. HCPCS: K0108

6.5"x8" Tapered Gel Footplate Covers for Width Extension \$225
Part: ST-FPG6508T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$225
Part: ST-FPG5011T-Q. HCPCS: K0108

7"x11" Tapered Gel Footplate Cover for Width Extension \$250
Part: ST-FPG7011T-Q. HCPCS: K0108

16.2. 70° Power Articulating, Swing Away Legrests

Not available with contoured seating. 4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

P L R

70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 14.5"-18" Extension
Pair Part: FRMASMB9073. Left Part: FRMASMB9074.
Right Part: FRMASMB9075. HCPCS: E1010+K0040+E0951

P L R

70° Articulating, Swing-Away \$1,197.50 ea.
Legrest with 11"-13" Extension
Pair Part: FRMASMB11444. Left Part: FRMASMB11366.
Right Part: FRMASMB11367. HCPCS: E1010+K0040+E0951

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
Angle Adjustable, 11" - 13"	✓	✓	✓	! This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 14.5" - 18"	!	!	!	

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**16.2. 70° Power Articulating, Swing Away Legrests continued****2. Angle Adjustable Footplate with Heel Loop and Wheel Bumper**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐	4"Wx6"D Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.	☐ ☐ ☐	5"Wx8"D Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.	☐ ☐ ☐	6"Wx11"D¹ Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.
☐ ☐ ☐	5"Wx6"D Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.	☐ ☐ ☐	6"Wx8"D¹ Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.	☐ ☐ ☐	8"Wx11"D² Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest Accessories

Calf Pad (pair) 6"Wx6"H..... No Charge
Part: MEC153735

Falcon Gel Padded Knee Adductor Buttons (KBE)(pair)\$125
Part: INDPART2242. HCPCS: E0953^
^Please verify the appropriate billing code to be used for this item with all third-party payors.

16.3. Pediatric, Swing-Away Legrests

Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Pediatric Legrests Upper

☐ ☐ ☐	90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension)\$81 ea. Pair part: FRMASMB8808. Left: FRMASMB8809 Right: FRMASMB8810
--	--

2. Angle Adjustable Footplate(s)

☐ ☐ ☐	4"Wx6"D Angle Adjustable\$85 ea. Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560 HCPCS: K0040	☐ ☐ ☐	5"Wx6"D Angle Adjustable\$85 ea. Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561 HCPCS: K0040
--	--	--	--

3. Accessories

Heel Loops (pair).....\$50
Part: FRMASMB7873. HCPCS: E0951

☐ ☐ ☐ **Wheel Bumper(s)\$12.50 ea.**
Pair part: WHL169508. Left: WHL169510. Right: WHL169511.

16.4. Swing-Away Legrests

Adjustable in 5° increments. Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

☐ ☐ ☐	60°-90° Swing-Away Legrest..... No Charge Pair part: RIG120062. Left: FRMASMB17558. Right: FRMASMB17559 Not available with 11"-13" extensions.	☐ ☐ ☐	60°-90° Swing-Away Legrest, Short No Charge Pair part: RIG136481. Left: RIG132607. Right: RIG132609. For use only with 11"-13" lower extensions.
--	--	--	--

1b. 70° Swing-Away Legrest Upper

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

70° Swing-Away Legrest (pair)\$162
Pair part: FRMASMB7841
Not available with 11"-13" extensions.

70° Tapered, Swing-Away Legrest (pair)\$162
Pair part: FRMASMB7843
Not available with 11"-13" extensions.

70° Tapered, Swing-Away Legrest, Short (pair)\$162
Pair part: FRMASMB8052
For use only with 11"-13" lower extensions.

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat				
	Low	Medium	High		
Extruded, 13" - 16.5"	!	✓	✓	✓	This rigging configuration causes no interference issues at this seat-to-floor height.
Extruded, 15.5" - 19.5"	!	!	!	!	This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 13" - 16.5"	!	✓	✓	✓	This rigging configuration does not work with this seat-to-floor height.
Angle Adjustable, 15.5" - 19.5"	!	!	!	!	
Multi-Axis, 15.5" - 19.5"	!	!	!	!	

2. Lower Extension

☐ ☐ ☐	11"-13" Length Part: See XRef N/A with Multi-axis.	☐ ☐ ☐	11"-16.5" Length Part: See XRef N/A with Angle Adjustable.	☐ ☐ ☐	13"-16.5" Length Part: See XRef N/A with Multi-Axis.	☐ ☐ ☐	15.5"-19.5" Length Part: See XRef
--	---	--	---	--	---	--	---

3a. Extruded Footplate(s)

☐ ☐ ☐	7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	☐ ☐ ☐	8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
--	---	--	---

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**16.4. 60°-90° Swing-Away Legrests continued****3b. Angle Adjustable Footplate(s) with Heel Loops**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	4"Wx6"D \$85 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040+E0951	P	L	R	6"Wx8"D ¹ \$85 ea. Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563 HCPCS: K0040+E0951
P	L	R	5"Wx6"D \$85 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040+E0951	P	L	R	6"Wx11"D ¹ \$85 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040+E0951
P	L	R	5"Wx8"D \$85 ea. Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 HCPCS: K0040+E0951	P	L	R	8"Wx11"D ² \$85 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040+E0951

3c. Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	4"Wx6"D \$247.50 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 HCPCS: K0040+E0951	P	L	R	6"Wx8"D ¹ \$247.50 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 HCPCS: K0040+E0951
P	L	R	5"Wx6"D \$247.50 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 HCPCS: K0040+E0951	P	L	R	6"Wx11"D ¹ \$247.50 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 HCPCS: K0040+E0951
P	L	R	5"Wx8"D \$247.50 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 HCPCS: K0040+E0951	P	L	R	8"Wx11"D ² \$247.50 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 HCPCS: K0040+E0951

4. 60°-90° Swing-Away Legrest Accessories

Heel Loops (pair) \$50
Part: FRMASMB7873. HCPCS: E0951
Standard on angle adjustable and multi-axis footplates.

Wheel Bumper \$12.50 ea.
Part: See XRef

Gel Neoprene Footrest Sleeve (pair) \$80
Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
Part: INDPART2823. HCPCS: K0108

16.5. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P	L	R	Elevating Legrest Assembly \$162.50 ea. Pair part: FRMASMB6483. Left: RIG138012 Right: RIG138013. HCPCS: E0990 Select lower extension and footplates
---	---	---	---

2. Lower Extension

P	L	R	11"-14" Length Pair: FRMASMB4755 Left: RIG138026 Right: RIG138030 For reinforced only	P	L	R	14"-17" Length Pair: FRMASMB4756 Left: RIG138033 Right: RIG138035 For reinforced only
P	L	R	12"-14.5" Length Part: See XRef For angle adjustable only	P	L	R	14.5"-18" Length Part: See XRef For angle adjustable only

3a. Reinforced Cast Footplate(s)

P	L	R	7.25"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056	P	L	R	8.5"Wx6"D Reinforced Footplate Standard Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062 May cause interference on seat widths less than 18."
---	---	---	---	---	---	---	--

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20". ³4"x6" and 5"x6" footplates are only available in pairs.

P	L	R	4"Wx6"D \$85 ea. Pair parts: FRMASMB8812 x2. HCPCS: K0040	P	L	R	6"Wx8"D ¹ with Heel Loops \$85 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951
P	L	R	5"Wx6"D \$85 ea. Pair parts: FRMASMB8811 x2. HCPCS: K0040	P	L	R	6"Wx11"D ¹ with Heel Loops \$85 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
P	L	R	5"Wx8"D with Heel Loops \$85 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951	P	L	R	8"Wx11"D ² with Heel Loops \$85 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12"-14.5"	✓	✓	✓
Angle Adjustable, 14.5"-18"	!	!	!

 This rigging configuration causes no interference issues at this seat-to-floor height.

 This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

Legrests section continued on following page.

16. LEGRESTS CONTINUED**4. Elevating Legrest Accessories**

Heel Loops (pair).....\$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable footplates.

Wheel Bumper.....\$12.50 ea.
 Part: See XRef
 Available on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair)\$125
 Part: INDPART2242. HCPCS: E0953^
 ^Please verify the appropriate billing code to be used for this item with all third-party payors.

17. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	


18. SEAT CUSHION

Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

18.1. Seat Cushion Model Selection**Stealth Simplicity G Cushion**


General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
 Part series: ST-SIMGXXYY

Simplicity G, 10"W to 21" W, E2601.....\$110
 Simplicity G, 22"W, E2602.....\$165

**Stealth Spectrum Foam SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
 Part series: ST-SPCFSPXXYY

Spectrum Foam SPP, 16"W to 20" W, E2607.....\$460
 Spectrum Foam SPP, 22"W, E2608.....\$530

**Stealth Solution SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
 Part series: ST-SOLSPPXXYY

Solution SPP, 10"W to 20" W, E2607.....\$375
 Solution SPP, 22"W, E2608.....\$430

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.
 Part series: ST-GLSPXXYY

Glacial SP, 10"W to 21"W, E2603.....\$330
 Glacial SP, 22"W, E2604.....\$440

**Stealth Zen SP Cushion**


Skin protection with Coolcore, 4-way stretch, breathable cover.
 Part series: ST-ZENSPXXYY

Zen SP, 10"W to 21" W, E2603.....\$315
 Zen SP, 22"W, E2604.....\$425

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
 Part series: ST-TC2SPPXXYY

TRU-Comfort 2 SPP, 16"W to 20" W, E2607.....\$400
 TRU-Comfort 2 SPP, 22"W, E2608.....\$500

**Stealth Spectrum Gel SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
 Part series: ST-SPCGSPPXXYY

Spectrum Gel SPP, 14"W to 20" W, E2607.....\$460
 Spectrum Gel SPP, 22"W, E2608.....\$530

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
 Part series: ST-ESSSPPXXYY

Essence SPP, No Wedge, 16"-20" Wide, E2607.....\$400
 Essence SPP, No Wedge, 22" Wide, E2608.....\$500
 Essence SPP, Small Wedge, 16"-20" Wide, E2607....\$435
 Essence SPP, Small Wedge, 22" Wide, E2608.....\$535
 Essence SPP, Medium Wedge, 16"-20" Wide, E2607..\$465
 Essence SPP, Medium Wedge, 22" Wide, E2608.....\$565
 Essence SPP, Large Wedge, 16"-20" Wide, E2607....\$495
 Essence SPP, Large Wedge, 22" Wide, E2608.....\$595

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-PREMPXXYY

Premiere P, 16"W to 20" W, E2605.....\$325
 Premiere P, 22" W, E2606.....\$475

**18.2. Seat Cushion Dimensions**

REQUIRED with cushion model selection. Select to match or pick non-matching dimensions. See respective cushion order form for available sizes.

Match Seat Width and Depth

Cushion Width

10" W	14" W	16" W	18" W	20" W	22" W
12" W	15" W	17" W	19" W	21" W	

Cushion Depth

10" D	14" D	16" D	18" D	20" D	22" D
12" D	15" D	17" D	19" D	21" D	

18.3. Seat Cushion Substitute/Add Fluid Proof Liners

Fluid Proof Urethane Liner.....\$50-\$70
 Part: URL1709070-XXYY
 Not available with Simplicity, Essence, or TRU-Comfort 2.

Fluid Proof Polyurethane Liner.....\$55-\$75
 Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY
 Not available with Simplicity.

Seat cushion section continued on following page.

18. SEAT CUSHION CONTINUED**18.4. Seat Cushion Substitute Cover**

An additional cover may be ordered or the standard cover may be substituted with a cover by checking the option in subsection 2. The ""prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover..... No Charge
Part: *prefix-TBLKXXYY
Not available with Simplicity, TRU-Comfort 2, Essence, or Premiere.

Mesh Cover No Charge
Part: *prefix-MBLKXXYY
Not available with Simplicity or TRU-Comfort 2.

19. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder..... \$20
Part: ACCASMB2287. Part with 12" lift: ACCASMB2288

Synergy Seat Cup Holder..... No Charge
Part: ACCASMB2634

Oxygen Holder for Contoured Seat \$118
Part: ACCASMB2292. HCPCS: E2208

Oxygen Holder for Synergy Seat..... \$150
Part: FRMASMB7030. HCPCS: E2208
Not available with Synergy Manual Recline.

Walker Holder \$103
Contoured seat part: ACCASMB1912. Synergy part: FRMASMB7031

Cane/Crutch Holder \$110
Contoured seat part: ACCASMB2294. HCPCS: E2207
Synergy part: FRMASMB4822. HCPCS: E2207

Cell Phone Holder for 12" Lift Captain Seat..... \$70
Part: ACC1611300

Rear Basket \$118
Contoured seat part: ACCASMB1914
Synergy part: FRMASMB4821

Flag Holder \$25
Contoured seat part: ACCASMB1930
Synergy part: ACCASMB2415.

Saddle Bag \$35
Part: DWR1020D010

Stealth Cell Phone/Media Device Holder \$195
Part: ACC7120006

Stealth Large Tablet Holder \$375
Part: ACC7120005

Quantum Backpack \$150
Part: ACC1905853

Grip Lap Board (20005)..... \$38.95
Part: ACC152992

Grip Drink Holder (20085) \$15.50
Part: ACC152993

Accessory Bar \$100
12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975.
20"W-22"W part: MEC146018

Large, Black, Weather Cover \$202
Part: ACCCOVR1021

20. BRANCH PATCHES

Optional.

Air Force Patch \$20
Part: ACC140240

Army Patch \$20
Part: ACC140241

Coast Guard Patch \$20
Part: ACC140239

Marines Patch \$20
Part: ACC140243

Navy Patch \$20
Part: ACC140244

Purple Heart Patch \$20
Part: ACC140238



Military Patches

21. ELECTRONICS ACCESSORIES

Optional.

☐ Q-Logic EX/NE Programmer \$545
Part: ELEASMB6976

☐ Q-Logic PC Programming Disk with Cable \$750
Part: ELEASMB5216

22. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Only available with the Synergy Seating system. Select a type of specialty controls below or to refer to the specialty controls order form.

22.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$1,720
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Non-matching Shroud \$50 Match Shroud to Base Color \$50
Part: DGN137728 Select color below Part: See color below

Red Corvette
Part: DGN137729-PRD1

Night Train
Part: DGN137729-PBK3

Evolution Orange
Part: DGN137729-POR7

Midnight Blue
Part: DGN137729-PBU10

Green Machine
Part: DGN137729-PGN11

Electric Blue
Part: DGN137729-PBU9

White Satin
Part: DGN137729-PWH2

Deep Purple
Part: DGN137729-PPR5

Lemon Crush
Part: DGN137729-PYL4

Military Green
Part: DGN137729-PGN8

Back in Black
Part: DGN137729-PBK1

Pretty in Pink
Part: DGN137729-PPK3

Grey Street
Part: DGN137729-PSL3



Enhanced Display

Alternative controls and electronics section continued on following page.

22. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**22.2. SCIM Module****Optional.**

- ☐ SCIM Module for Q-Logic 3 \$525
Part: CTL168055

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

22.3. Q-Logic 3 SAJ

- ☐ Q-Logic 3 EX SAJ \$700
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

22.4. Head Array

Adult, i-Drive Tri-Array Series, 3-Switch \$4,414
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461
Head Array (IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

22.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 \$1,895

Part: CTL168056. HCPCS: E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and Mount

Stealth Sip and Puff for Comfort Series, \$294.25
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series, \$294.25
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$294.25
Head Array (SP625-U)

Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Sip & Puff Module



Stealth Sip & Puff

22.6. Egg Switches

Red Egg Switch (ESRED) \$80.85 ea.
Part: SWTMCRO1077

Blue Egg Switch (ESBLU) \$80.85 ea.
Part: SWTMCRO1079

Black Egg Switch (ESBLK) \$80.85 ea.
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$80.85 ea.
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$80.85 ea.
Part: SWTMCRO1080



Egg Switch

22.7. Switch Adapter

2 to 1 Switch Adapter \$13.39
Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

22.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,400
Part: CTL168057

23. POSITIONING COMPONENTS

Optional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).

23.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$472.50
Part: INDPART2896. HCPCS: E1028

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$210
Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$210
Part: INDPART2664. HCPCS: E0956

3. Installation

Install Laterals
Part: PTOINDV3499

Do Not Install
Part: PTOINDV3500

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, \$245 ea.
 Removable Mounting Hardware
 Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable, \$255 ea.
 Removable Mounting Hardware
 Part: ST-UTTP7-Q

2. UniLink Pelvic/Thigh Pads

^APlease verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$95 ea.
 Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad \$125 ea.
 Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad \$105 ea.
 Part: ST-UTP548056100-Q

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link \$75 ea.
 Part: ST-UIL00150A000-Q

UniLink, 2" Link \$75 ea.
 Part: ST-UIL001200A00-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

23.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea.
 Part: INDPART2804. HCPCS: E1028

Flip Down with LEFT Full Surface Contact \$236.25 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Removable Bracket (LP430) \$262.50 ea.
 Part: INDPART2805. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact \$236.25 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

Removable Bracket with 1" Offset \$400 ea.
 (LP430 & LHW-121)
 Part: INDPART2854. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122) \$26.25
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129) \$53.50
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

^APlease verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4) \$75 ea.
 Part: INDPART2801. HCPCS: E0953~/E0956

Long Pad, 4"x12" (TWBADD-4X12) \$75 ea.
 Part: INDPART2803. HCPCS: E0953~/E0956

Medium Pad, 4"x8" (TWBADD-4X8) \$75 ea.
 Part: INDPART2802. HCPCS: E0953~/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

23.4. Stealth i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103) \$142 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$142 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$142 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

**23.5. Stealth i-Fit Chest Straps**

Chest Strap, Medium w/ Belt Anchor \$90 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$90 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

**23.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

Facial Lateral for Comfort Series \$294.25 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Laterals for Comfort Series \$347.75 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

Facial Lateral for Comfort Series \$294.25 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028



Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED

2. Facial Lateral Pads and Accessories

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25
Part: POS150532
Required with facial laterals to mount to comfort headrest brackets.

Egg Switch Mount (ESM) (A) \$52.50 ea.
Part: INDPART2815


Figure 8 Bracket with 13" Rod \$267.50 ea.
and Egg Switch Mount (B)
Part: ST-F8675-13-Q

Facial Pad, 4"x2" (900-Q) \$53.50 ea.
Part: ST-900-Q. HCPCS: K0108

Facial Pad, 6"x2" (920-Q) \$53.50 ea.
Part: ST-920-Q. HCPCS: K0108

Swing-Away Joystick Mount, \$250 ea.
18" Rod (C)
Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest



23.7. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support Bracket, Tube \$183.75 ea.
Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube \$183.75 ea.
Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$52.50 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$52.50 ea.
Part: POS124067. HCPCS: K0108

23.8. Stealth Foot Boxes

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$498.75
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$498.75
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$498.75
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954^



23.9. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder, \$266
3.75"Wx6.25"D (SPLE111) (pair)
Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, \$266
5"Wx10.75"D (SPLE114) (pair)
Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, \$266
4"Wx8.25"D (SPLE112) (pair)
Part: INDPART2784. HCPCS: E0951+E0952



24. TRU-BALANCE 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA) No Charge
Part: Configuration selection

Secondary Multiplier (SM) \$221.21
Part: HAR144004

Special Order Configuration (SOF)
Part: Configuration selection

Arm Mounts	Back Mounts			
		1 Arm	2 Arm	3+ Arm
Joystick (Hand Control) Enhanced Display iAccess USB Charger SAJ	Attendant Control SCIM Sip n Puff ECU	0 Back	NA	SM
		1 Back	NA	SM
		2 Back	SM	SM
		3+ Back	SOF	SOF

Q-Logic 3 Multiplier with Synergy Seat Chart

25. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)

Quantum® Rehab

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

©2019 Quantum Rehab - A Pride Mobility Products Corporation company. All rights reserved.

