

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Q6 Edge® 2.0 X Order Form with TRU-Balance® 3 Positioning**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Quote Name/Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

Q6 edge®
2.0x**INTRODUCTION**

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. This power base has a 300 lb. weight capacity, standard 6 mph programmable motors, 8 amp off board charger, fender lights, 2 multipliers, and USB mobile device charger. Power positioning includes an attendant control and iAccess interface. Dimensions of base without legrests: 24"W x 35.5"L.

Quantum Q6 Edge 2.0 X 3S-SS\$10,195

Part: Q6EDGE 2.0X 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge 2.0 X 3MP-SS\$10,195

Part: Q6EDGE 2.0X 3MP-SS. HCPCS: K0861

Multiple actuator power option.

Quantum Q6 Edge 2.0 X 3SP-SS\$10,195

Part: Q6EDGE 2.0X 3SP-SS. HCPCS: K0856

Single power option.

Account Number: _____ Date: _____ Quote Name: _____

2. DRIVE MOTORS

REQUIRED.

- ☐ 6 mph Drive Motors Standard
Part: MOT147433
- ☐ Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

- | | | | |
|--|---|---|---|
|  Red Corvette
Part: DGN153201-PRD1 |  Night Train
Part: DGN153201-PBK3 |  Deep Purple
Part: DGN153201-PPR5 |  Lemon Crush
Part: DGN153201-PYL4 |
|  Electric Blue
Part: DGN153201-PBU9 |  Pretty in Pink
Part: DGN153201-PPK3 |  Grey Street
Part: DGN153201-PSL3 |  Green Machine
Part: DGN153201-PGN11 |
|  Back in Black
Part: DGN153201-PBK1 |  Evolution Orange
Part: DGN153201-POR7 |  Midnight Blue
Part: DGN153201-PBU10 |  Military Green
Part: DGN153201-PGN8 |

4. QUANTUM BACKUP CAMERA

Optional The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21
Part: MEC1910373
Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123

Backup Camera NE/NE+ Joystick Mount Bracket \$105
Part: MEC2106132



5. WHEEL OPTIONS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

5.1. Drive Wheels

- | | |
|---|---|
| 3"x14" Flat Free Black Tires (No Accent Color) Standard
Part: PTOINDV3525 (pair) | 3"x14" Pneumatic Black Tires (No Accent Color) No Charge
Part: WHL150034 x2 |
| 3"x14" Flat Free Black Tires for Wheel Accents No Charge
Part: WHL1611297 x2 | 3"x14" Pneumatic Black Tires for Wheel Accents No Charge
Part: WHL1611296 x2 |
| 4"x14" Knobby, Flat Free Black Tires (No Accent Color) \$290 ea.
Part: WHL2006647 x2 | 4"x14" Knobby, Pneumatic Black Tires (No Accent Color) \$290 ea.
Part: WHL2006822 x2 |

5.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

- | | | |
|--|---|---|
| Back in Black Accent Color \$63
Part: XRF000179 | Deep Purple Accent Color \$63
Part: XRF000184 | Electric Blue Accent Color \$63
Part: XRF000188 |
| Red Corvette Accent Color \$63
Part: XRF000180 | Pretty in Pink Accent Color \$63
Part: XRF000185 | Lemon Crush Accent Color \$63
Part: XRF000189 |
| Night Train Accent Color \$63
Part: XRF000181 | Evolution Orange Accent Color \$63
Part: XRF000186 | Green Machine Accent Color \$63
Part: XRF000190 |
| Grey Street Accent Color \$63
Part: XRF000182 | Midnight Blue Accent Color \$63
Part: XRF000187 | Military Green Accent Color \$63
Part: XRF000191 |

6. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

- NF-22, Gel Cell, Interceptor Battery \$341 ea.
Part: BAT2004376. HCPCS: E2361
- NF-22, AGM, MK Battery \$415 ea.
Part: BAT121346. HCPCS: E2361
- NF-22, Gel Cell, MK Battery \$415 ea.
Part: BATGEL1008. HCPCS: E2361

7. ELECTRONICS

The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

7.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Joystick**REQUIRED.**

Q-Logic 3 EX JoystickStandard
Part: XRF000241

No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display.

Q-Logic 3e Joystick No Charge
Part: CTL1910600

Can operate full light package. Not available with alternative or specialty controls.



Q-Logic 3 Joysticks

2. Q-Logic 3 Required Components**REQUIRED with 3SP-SS or 3MP-SS.**

Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness.....\$447
Part: XRF000264. HCPCS: E2313

7.2. Q-Logic 3 Joystick Armor

REQUIRED with Q-Logic 3 Joystick selections above. Reference colors shown on page 2.

Black, Soft TouchStandard
Part: DGN124307
Also matches Night Train shrouds.

Red Corvette..... \$63
Part: DGN124336-PRD1

Back in Black.....\$63
Part: DGN124336-PBK1

Grey Street..... \$63
Part: DGN124336-PSL3

Deep Purple.....\$63
Part: DGN124336-PPR5

Pretty in Pink.....\$63
Part: DGN124336-PPK3

Evolution Orange.....\$63
Part: DGN124336-POR7

Midnight Blue\$63
Part: DGN124336-PBU10

Electric Blue\$63
Part: DGN124336-PBU9

Lemon Crush.....\$63
Part: DGN124336-PYL4

Green Machine\$63
Part: DGN124336-PGN11

Military Green \$63
Part: DGN124336-PGN8

8. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Chin Cup (C)..... \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323

Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323

Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Bodypoint 3" Goal Post Handle\$104
Part: PLSKNCP1057. HCPCS: E2323

Bodypoint 4" Goal Post Handle\$104
Part: PLSKNCP1058. HCPCS: E2323

Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308



Joystick Handles

9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

9.1a. Fixed, In-line Joystick Bracket

☐ Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617

9.1b. Swing-Away (Retractable) Joystick Bracket


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3.....\$294
Part: ACC1908434. HCPCS: E1028

9.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side\$245
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side\$245
Part: ACC2201042. HCPCS: E1028



Swing Away Joystick

Joystick mounting brackets section continued on following page.

9. JOYSTICK MOUNTING BRACKETS CONTINUED**9.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368

Part: XRF000003. HCPCS: E1028

Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126


Part: XRF000004

Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket \$250

Part: XRF000005. HCPCS: E1028

Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.

**9.2. Joystick Mount Location**☐ **Left Side Mount**
Part: PTOINDV1225☐ **Right Side Mount**
Part: PTOINDV1226**10. TRU-BALANCE® 3 SEATING - QUANTUM INTELLIGENT SYSTEMS™****REQUIRED.** Maximum seat size is 22"W x 22"D.**10.1. TRU-Balance 3 Configuration****TRU-Balance 3 Static Seat** \$2,410

Part: PTOINDV3747

20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU-Balance 3 Power Tilt Only \$7,944

Part: PTOINDV3748. HCPCS: E1002

50° Center of Gravity (CG) power tilt seating system.

TRU-Balance 3 Power Recline Only \$7,975

Part: PTOINDV3749. HCPCS: E1004

175° of power recline.

TRU-Balance 3 Power Tilt & Recline \$12,490

Part: PTOINDV3750. HCPCS: E1007

50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

10.2. 12" Power Lift**Optional.** Select one option below.**TRU-Balance 3, 12" Power Lift with iLevel®** \$2,730

Part: PTOINDV3751. HCPCS: E2300

12" of scissor lift with iLevel® System. Not available with high speed motors.

**10.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option. Only available with 16"-20" width and depth packages. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**WC-19 Occupied Transit Option** \$289

Part: KIT152047

10.4. Width**REQUIRED with seat selection.** Maximum seat width is 22". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W	13" W	14" W	15" W	16" W	17" W	18" W	19" W	20" W	21" W	22" W
(#1002)	(#1003)	(#1004)	(#1005)	(#1006)	(#1007)	(#1008)	(#1009)	(#1010)	(#1011)	(#1012)

10.5. Width Package**REQUIRED with seat selection.****12"-16" Width Package**

Part: XRF000006

14"-16" Width Package

Part: XRF000007

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

Maximum width is 22".

10.6. Depth**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D	13" D	14" D	15" D	16" D	17" D	18" D	19" D	20" D	21" D	22" D
(#1002)	(#1003)	(#1004)	(#1005)	(#1006)	(#1007)	(#1008)	(#1009)	(#1010)	(#1011)	(#1012)

10.7. Depth Package**REQUIRED with seat selection.****12"-16" Depth Package**

Part: XRF000011

14"-16" Depth Package

Part: XRF000012

16"-20" Depth Package

Part: XRF000013

20"-24" Depth Package

Part: XRF000159

Maximum depth is 22".

10.8. Seat to Floor Height**REQUIRED.** Select a seat to floor height compatible with the seat functions previously selected.**Lowest Seat to Floor**

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559

Highest Seat to Floor

Part: PTOINDV3560

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
17.5"	18.5"	19.5"
TRU-Balance 3 Seat to Floor Heights		

TRU-Balance 3 Seating section continued on following pages.


10. TRU-BALANCE 3 SEATING CONTINUED**10.9. Matching Seat Cushion**

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

Stealth TRU-Comfort 2 SPP Cushion


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY. HCPCS: E2607/E2608

TRU-Comfort 2 SPP, Matching Seat Size
Reference part: XRF000116. MSRP: \$420/\$525

**Stealth Solution SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SOLSPPXXYY. HCPCS: E2607/E2608

Solution SPP, Matching Seat Size
Reference part: XRF000113. MSRP: \$394/\$452

**Stealth Spectrum Foam SPP Cushion**


Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SPCFSPXXYY. HCPCS: E2607/E2608

Spectrum Foam SPP, Matching Seat Size
Reference part: XRF000112. MSRP: \$483/\$557

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. HCPCS: E2605/E2606

Premiere P, Matching Seat Size
Reference part: XRF000122. MSRP: \$341/\$499

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPXXYY. HCPCS: E2603/E2604

Zen SP, Matching Seat Size
Reference part: XRF000115. MSRP: \$331/\$446

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Includes urethane liner.
Part series: ST-SPCGSPPXXYY. HCPCS: E2607/E2608

Spectrum Gel SPP, Matching Seat Size
Reference part: XRF000117. MSRP: \$483/\$557

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPPXXYY. HCPCS: E2607/E2608

Essence SPP, No Wedge, Matching Seat Size
Reference part: XRF000118. MSRP: \$420/\$525

Essence SPP, Small Wedge, Matching Seat Size
Reference part: XRF000119. MSRP: \$457/\$562


Essence SPP, Medium Wedge, Matching Seat Size
Reference part: XRF000120. MSRP: \$488/\$593

Essence SPP, Large Wedge, Matching Seat Size
Reference part: XRF000121. MSRP: \$520/\$625

**Stealth Glacial SP Cushion**


Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. HCPCS: E2603/E2604

Glacial SP, Matching Seat Size
Reference part: XRF000114. MSRP: \$347/\$462

**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. HCPCS: E2601/E2602

Simplicity G, Matching Seat Size
Reference part: XRF000111. MSRP: \$116/\$173

**10.10. Back Type Selection**

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back

Part: XRF000014

12" and 24" back width will use the TRU-Balance 3 Back assembly.
Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: XRF000015

Available only with TRU-Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184

Part: XRF000016

12" and 24" back width will use the TRU-Balance 3 Back assembly. Not available with power recline. Available with Static Seat and Power Tilt (1 actuator). Only available with 2-post or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

10.11. Back Cane Selection (Back Cane Option Only)

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable No Charge
Part: MEC140266

Flip Forward \$53
Part: XRF000273

2" Endomorph (1" per side) Back Cane \$74
Part: XRF000269

4" Endomorph (2" per side) Back Cane..... \$74
Part: XRF000271

2. Back Cane Height

16"H Cane
Part: FRMASMB15834

17"H Cane
Part: FRMASMB15835

20"H Cane
Part: FRMASMB15838

21"H Cane
Part: FRMASMB15839

24"H Cane
Part: FRMASMB15842

25"H Cane
Part: FRMASMB15843

TRU-Balance 3 Seating section continued on following pages.

10. TRU-BALANCE 3 SEATING CONTINUED**10.12. Back Cushion Configuration - Type**

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Quantum Backs

Sport Back No Charge
Part: XRF000017

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 12" Width \$709
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$709
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$709
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$814
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$814
Part: TC4V22W32WBC. HCPCS: E2621

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface \$210
Part: XRF000018. Select a set width in following section, range of 16"-24". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.

1b. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$540.75
Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$540.75
Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$540.75
Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$540.75
Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$540.75
Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$540.75
Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$540.75
Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AL20AXX. HCPCS: E2613

1c. Q-Back 2.0, AL Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$612.85
Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$612.85
Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$612.85
Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$612.85
Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$612.85
Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$612.85
Part: QB2AL20DXX. HCPCS: E2620

TRU-Balance 3 Seating section continued on following page.

10. TRU-BALANCE 3 SEATING CONTINUED**1d. Q-Back 2.0, AL Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL14ES. HCPCS: E2620Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL14ET. HCPCS: E2620Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline.. \$669.50
Part: QB2AL16ES. HCPCS: E2620Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL16ET. HCPCS: E2620Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline.. \$669.50
Part: QB2AL16EX. HCPCS: E2620Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline.. \$669.50
Part: QB2AL16EXX. HCPCS: E2620Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL18ES. HCPCS: E2620Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline.. \$669.50
Part: QB2AL18ET. HCPCS: E2620Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline . \$669.50
Part: QB2AL18EX. HCPCS: E2620Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline . \$669.50
Part: QB2AL18EXX. HCPCS: E2620Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline . \$669.50
Part: QB2AL20ES. HCPCS: E2620Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline . \$669.50
Part: QB2AL20ET. HCPCS: E2620Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline . \$669.50
Part: QB2AL20EX. HCPCS: E2620Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline . \$669.50
Part: QB2AL20EXX. HCPCS: E2620**1i. Q-Back 2.0, AFT Series, Active Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AFT16AS. HCPCS: E2613AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$540.75
Part: QB2AFT18AS. HCPCS: E2613**1j. Q-Back 2.0, AFT Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline... \$612.85
Part: QB2AFT16DS. HCPCS: E2620AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline... \$612.85
Part: QB2AFT18DT. HCPCS: E2620AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline... \$612.85
Part: QB2AFT16DT. HCPCS: E2620AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline .. \$612.85
Part: QB2AFT20DT. HCPCS: E2620AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline .. \$612.85
Part: QB2AFT18DS. HCPCS: E2620**1k. Q-Back 2.0, AFT Series, Extra Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$669.50
Part: QB2AFT16EXX. HCPCS: E2620**10.13. TRU-Balance 3 Back Shroud****REQUIRED.** The back shroud is a feature of the TRU-Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.Red Corvette
Part: XRF000019Grey Street
Part: XRF000023Pretty in Pink
Part: XRF000027Electric Blue
Part: XRF000030Green Machine
Part: XRF000033Back in Black
Part: XRF000021Deep Purple
Part: XRF000025Evolution Orange
Part: XRF000028Blue Jean
Part: XRF000031Military Green
Part: XRF000034Night Train
Part: XRF000022Raspberry Beret
Part: XRF000026Midnight Blue
Part: XRF000029Lemon Crush
Part: XRF000032**10.14. Lateral Positioning Wedge****Optional.**

Lateral Positioning Wedge, pair \$32

Part: SET123465

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

10.15. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.Match Seat Width Selections
Part: PTOINDV346814"W-16"W Package
Part: XRF00003718"W-22"W Package
Part: XRF00003910"W-14"W Package
Part: XRF00003616"W-20"W Package
Part: XRF00003820"W-24"W Package
Part: XRF000040**10.16. Back Width****REQUIRED** except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.12" Back Width
Part: PTOINDV356216" Back Width
Part: PTOINDV356419" Back Width
Part: PTOINDV356722" Back Width
Part: PTOINDV357014" Back Width
Part: PTOINDV356317" Back Width
Part: PTOINDV356520" Back Width
Part: PTOINDV356824" Back Width
Part: PTOINDV357115" Back Width
Part: PTOINDV385618" Back Width
Part: PTOINDV356621" Back Width
Part: PTOINDV3569

TRU-Balance 3 Seating section continued on following page.

10. TRU-BALANCE 3 SEATING CONTINUED**10.17. Overall Back Height**

Static Seat and Power Tilt available height range: 14"-25"H. Odd heights include 4" gap and even heights include 3" gap between seat pan and bottom of back cushion.
Power Recline available height range: 19"-28"H. Odd heights include 6" gap and even heights include 7" gap. See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	18" Height Part: SETHGHT1031	22" Height Part: SETHGHT1035	26" Height Part: SETHGHT1039	30" Height Part: SETHGHT1043
15" Height Part: SETHGHT1028	19" Height Part: SETHGHT1032	23" Height Part: SETHGHT1036	27" Height Part: SETHGHT1040	31" Height Part: SETHGHT1102
16" Height Part: SETHGHT1029	20" Height Part: SETHGHT1033	24" Height Part: SETHGHT1037	28" Height Part: SETHGHT1041	32" Height Part: SETHGHT1103
17" Height Part: SETHGHT1030	21" Height Part: SETHGHT1034	25" Height Part: SETHGHT1038	29" Height Part: SETHGHT1042	

	Tilt or Static Width												Recline Width												
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	24"			
Overall Height	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"		
	15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"		
	16"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	16"	
	17"	S	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	17"	
	18"	-	S	T	S	T	S	T	S	T	S	T	S	-	-	-	-	-	-	-	-	-	-	18"	
	19"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	T	19"
	20"	-	S	T	S	T	S	T	S	T	S	T	S	S	S	T	S	T	S	T	S	T	S	T	20"
	21"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	21"
	22"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	22"
	23"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	23"
	24"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	24"
	25"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	25"
	26"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	26"
	27"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	27"
	28"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	28"
	29"	-	S	T	S	T	S	T	S	T	S	T	S	-	S	T	S	T	S	T	S	T	S	T	29"
	30"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	S	30"
	31"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	S	31"
	32"	-	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	S	T	S	32"

S

 Sport Back

Q

 Q-Back

T

 TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

10.18. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$48
 Part: XRF000351

10.19. Back Height Extension

Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension \$48
 Part: XRF000041

11. TRU-BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

11.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

L R Flip-Back, 8"-16"H Armrests Standard Left part: ARM1806662. Right part: ARM1806663 Not available with TRU-Balance 3 Static Seat.	L R Flip-Back, 8"-16"H Armrests for Static Seat \$289 Left part: XRF000320. Right part: XRF000321
---	---

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

L R Flip-Back Armrest Downpost \$92 ea.
 Left part: XRF000265. Right part: XRF000266

TRU-Balance 3 Armrests section continued on following page.

11. TRU-BALANCE 3 ARMRESTS CONTINUED**11.1b. Two-Post, Flip-Back, Height Adjustable, Removable Armrests**

Select a left and right armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes. Not available with a depth less than 16". Required with flip forward back and lap trays.

1a. Two-Post Armrest for Power Seating

☐ ☐ Full Length, Two-Post Armrest for Power Seating No Charge
Left part: XRF000047. Right part: XRF000048

☐ ☐ Desk Length, Two-Post Armrest for Power Seating.... No Charge
Left part: XRF000049. Right part: XRF000050

1b. Two-Post Armrest for Static Seating

☐ ☐ Full Length, Two-Post Armrest for Static Seating.....\$158 ea.
Left part: XRF000051. Right part: XRF000052. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrest for Static Seating\$158 ea.
Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

11.1c. Heavy Duty, Single Post, Height Adjustable, Removable Armrests

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$79 ea.
Heavy Duty (Sold in pairs)
Part: XRF000056. HCPCS w/ Static Seat: E0973

11.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H
PTOINDV3535

10" H
PTOINDV3537

12" H
PTOINDV3539

14" H
PTOINDV3541

16" H
PTOINDV3543

9" H
PTOINDV3536

11" H
PTOINDV3538

13" H
PTOINDV3540

15" H
PTOINDV3542

2. Right Armrest Height

8" H
PTOINDV3450

10" H
PTOINDV3452

12" H
PTOINDV3454

14" H
PTOINDV3456

16" H
PTOINDV3458

9" H
PTOINDV3451

11" H
PTOINDV3453

13" H
PTOINDV3455

15" H
PTOINDV3457

11.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ ☐ 10" Length
Left part: XRF000057. Right part:
XRF000123

☐ ☐ 14" Length
Left part: XRF000058. Right part:
XRF000124

☐ ☐ 17" Length
Left part: XRF000059. Right part:
XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads.

12. TRU-BALANCE 3 ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

12.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge
Part: XRF000060
Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge
Part: XRF000308
Only available with Flip-back, 8"-16" armrests.

12.1b. Foam Armpads

☐ ☐ Straight, Armpad Standard
Left part: XRF000061. Right part: XRF000126

☐ ☐ Waterfall, Armpad No Charge
Left part: XRF000062. Right part: XRF000127
Not available with Flip-back, 8"-16" armrests.

12.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right.....\$195
Part: ARM1806807

Falcon Gel, Desk Length Left, Full Length Right.....\$171
Part: ARM1806810

Falcon Gel, Full Length Left, Desk Length Right.....\$171
Part: ARM1806813

Falcon Gel, Desk Length Left and Right\$147
Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right.....\$311.22
Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442

Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437

Account Number: _____ Date: _____ Quote Name: _____

13. TRU-BALANCE 3 ARMREST BRACKETS

13.1. Joystick Bracket

REQUIRED.

Standard Fixed Receiver Standard
Part: XRF000063

Stealth Height Adjustable Receiver \$225
Part: XRF000064. HCPCS: K0108



14. TRU-BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

14.1a. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM) \$170.25
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM) \$170.25
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM) \$170.25
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM) \$170.25
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink \$296.50
Part: ST-UHB2D33200-Q. HCPCS: E1028

14.1b. Headrests, Comfort Plus

1. Headrest Pads

Stealth 8" Comfort Plus (CP450) \$170.25
Part: POS148484. HCPCS: E0955

Stealth 10" Comfort Plus (CP250) \$170.25
Part: SETHEAD1057. HCPCS: E0955

Stealth 14" Comfort Plus (CP150) \$170.25
Part: INDPART2327. HCPCS: E0955

Pediatric Comfort Plus (6"Wx4.5"H) (CPN350) \$170.25
Part: SETHEAD1069. HCPCS: E0955

2. Headrest Brackets

Multi-Axis Removable Headrest Bracket (HMO475P) \$165.32
Part: HDWASMB1084. HCPCS: E1028

15. LAP BELTS

REQUIRED. Select a lap belt below.

15.1a. Standard Lap Belts

38" Lap Belt
Part: ACCBELT1006

50" Lap Belt
Part: ACCBELT1000

60" Lap Belt
Part: ACCBELT1002

70" Lap Belt
Part: ACCBELT1003

80" Lap Belt
Part: ACCBELT1004

90" Lap Belt
Part: ACCBELT1001

15.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350) \$172.53
Part: ST-IFRB5350

Stealth Positioning 60" Single Pull, Push Button, Medium \$47.38
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning 60" Single Pull, Push Button, Large \$47.38
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

Stealth Positioning 60" Single Pull, Push Button, Extra Large \$47.38
Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable) \$120
Required with ISO 7176-19 Compliant Transit Option.
Part: ACCBELT1016

16. TRU-BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

☐ Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108

☐ Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110

☐ Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109

17. BASE MOUNTED FOOT PLATFORM

REQUIRED. May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077

Omit Base Mounted Foot Platform No Charge
Part: PTOINDV3674

Bell Shaped, Base Mounted Foot Platform \$247
Part: RIGASMB7120041

18. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

18.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount No Charge
Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040

Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997

Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

18.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040

Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311. HCPCS: K0040

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility. Do not select with Basic Center Mount.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577

Medium Extension, 5"x11" Footplates
Part: RIG150212

Pediatric Extension, 5"x8" Footplates
Part: RIG152578

Long Extension, 5"x6" Footplates
Part: RIG151505

Medium Extension, 5"x6" Footplates
Part: RIG150206

Long Extension, 5"x8" Footplates
Part: RIG151506

Medium Extension, 5"x8" Footplates
Part: RIG150203

Long Extension, 5"x11" Footplates
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

Small Calf Pads, Quick Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156

Large Calf Pads, Quick Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel for AFP, 18"Wx8"H \$184
Part: FRMASMB7110007

Small Calf Panel for Quick Adjustable, 18"Wx8"H \$184
Part: FRMASMB7110012

Large Calf Panel for AFP, 24"Wx8"H \$237
Part: RIG123824
Only available on widths of 20" or greater.

Large Calf Panel for Quick Adjustable, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$69
Part: RIG152517

4" Outside Edge (pair) \$69
Part: RIG152510

Width Extension for 11"D Footplates \$69
Part: RIG152518

Only available with 6"D footplates and 8"D footplate extension.

2" Bell-Shaped Length Extension \$69
Part: RIG152519

6" Outside Edge (pair) \$69
Part: RIG152516

Available with 8"D and 11"D footplates and 11"D footplate extension.

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**5. Additional Accessories**

Small Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133717

Large Reversed Dartex Calf Pad Cover (pair) \$210
Part: RIG133719

Small Reversed Dartex Panel Cover (pair) \$210
Part: RIG134673

Large Reversed Dartex Panel Cover (pair) \$210
Part: RIG134744

Small Foam Calf Pad Cover (pair) \$63.74
Part: ST-CSCS0775-Q

Large Foam Calf Pad Cover (pair) \$78.72
Part: ST-CSCS0985-Q

Small Foam Panel Cover (pair) \$167.57
Part: ST-CSCP0920-Q

Large Foam Panel Cover (pair) \$268.77
Part: ST-CSCP0927-Q

Footplate Connector \$16
Part: RIGASMB7120023

Footplate Riser \$75
Part: XRF000348

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel \$50
Part: FRMASMB7110011

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$211.15
Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$231.75
Part: ST-FPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension ... \$231.75
Part: ST-FPG6508T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension ... \$257.50
Part: ST-FPG7011T-Q. HCPCS: K0108

18.1c. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,516
Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix Center Mount	Seat to Floor Option			Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High		Low	Middle	High
Power ELR: Static Seat, Recline Only				Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"	8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14"	12.5"-15"	12.5"-16"	14.5"-21" Lower Length	12.5"-15"	12.5"-16"	12.5"-17"

2. Lower Extension

8.5" to 14.5" Lower Length
Part: XRF000344

14.5" to 18.5" Lower Length
Part: XRF000345

13.5" to 19.5" Lower Length
Part: XRF000346

14.5" to 21.5" Lower Length
Part: XRF000347

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

¹The 6"W footplates may cause interference with seat widths less than 16". ²The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair
Part: FRMASMB7448

5"Wx8"D, Pair
Part: FRMASMB7450

6"Wx11"D¹, Pair
Part: FRMASMB7452

5"Wx6"D, Pair
Part: FRMASMB7449

6"Wx8"D¹, Pair
Part: FRMASMB7451

8"Wx11"D², Pair
Part: FRMASMB7453

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**18.1d. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1a. 60°-90° Swing-Away Legrest Upper

P L R 60°-90° Swing-Away Legrest, Power Seating..... No Charge
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Not available with 11"-13" extensions.

P L R 60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 For use only with 11"-13" lower extensions.

1b. 70° Swing-Away Legrest Upper

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

70° Swing-Away Legrest (pair) \$172

Part: FRMASMB7841

Not available with 11"-13" extensions.

2. Swing-Away Legrests Lower Extension

P L R 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.
 N/A with Multi-axis.

P L R 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.
 N/A with Multi-Axis.

P L R 11"-16.5" Length
 Pair part: XRF000091.
 Left: XRF000092.
 Right: XRF000093.
 Multi-Axis only.

P L R 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3a. Swing-Away Legrests, Extruded Footplate(s)

P L R 7.25"Wx6"D Extruded Footplate..... Standard
 Part pair: FRMASMB7847. Left: FRMASMB7861
 Right: FRMASMB7862

P L R 8.5"Wx6"D Extruded Footplate Standard
 Part pair: FRMASMB7848. Left: FRMASMB7863
 Right: FRMASMB7864
 May cause interference on seat widths less than 18."

3b. Swing-Away Legrests, Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 4"Wx6"D \$90 ea.
 Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ \$90 ea.
 Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563
 HCPCS: K0040+E0951

P L R 5"Wx6"D..... \$90 ea.
 Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ \$90 ea.
 Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564
 HCPCS: K0040+E0951

P L R 5"Wx8"D \$90 ea.
 Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562
 HCPCS: K0040+E0951

P L R 8"Wx11"D² \$90 ea.
 Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565
 HCPCS: K0040+E0951

3c. Swing-Away Legrests, Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 4"Wx6"D \$260 ea.
 Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ \$260 ea.
 Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555
 HCPCS: K0040+E0951

P L R 5"Wx6"D..... \$260 ea.
 Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ \$260 ea.
 Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558
 HCPCS: K0040+E0951

P L R 5"Wx8"D \$260 ea.
 Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554
 HCPCS: K0040+E0951

P L R 8"Wx11"D² \$260 ea.
 Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559
 HCPCS: K0040+E0951

4. Swing-Away Legrest Accessories

Gel Neoprene Footrest Sleeve (pair) \$84
 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

QTY

QTY

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**18.1e. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D \$90 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx8"D¹ with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx6"D \$90 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx11"D¹ with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx8"D with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

18.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$19 ea.
 Part: XRF000100

18.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

19. TRU-BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options.**19.1a. Actuator Control through Joystick**

Single Actuator Through Controller and\$2,163

Interface
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and \$3,040
Interface

Part: XRF000336. HCPCS: E2311

1. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations. 4-Way switches will be mounted to same side of joystick.

4-Function Push Button

Part: CTL123218

See section 20.2 for switch mounting.

4-Way Toggle

Part: ELEASMB5925

4-Way Button Toggle

Part: ELEASMB5926

2. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

**19.1b. Actuator Control with iAccess Switch****1. Actuator Control Selection**

Single Actuator Through Flat Button iAccess,\$2,547

Controller, and Interface

Part: XRF000337. HCPCS: E2310

Single Actuator Through Toggle Button iAccess,\$2,547

Controller, and Interface

Part: XRF000338. HCPCS: E2310

Multiple Actuator Through Flat Button iAccess,\$3,460
Controller, and Interface

Part: XRF000339. HCPCS: E2311

Multiple Actuator Through Toggle Button iAccess,\$3,460
Controller, and Interface

Part: XRF000340. HCPCS: E2311

20. TRU-BALANCE 3 ATTENDANT CONTROL**Optional.** The attendant control is standard with power positioning configurations.

Attendant Control (Static Seat) \$525

Part: XRF000078. HCPCS: E2331

Attendant Control (Power Seating) Standard

Part: XRF000079

Select this option if no charge option criteria is not met.

21. SEAT MOUNTED LED LIGHTS**Optional.** Requires NE+ 8 Key, Q-Logic 3 Joystick for Full Light Package, or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package\$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

**22. ACCESSORIES****Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU-Balance 3 Cane/Crutch Holder \$105

Part: ACC152515. HCPCS: E2207

Quantum Backpack \$168

Part: ACC1905853

Grip Lap Board (20005)..... \$38.95

Part: ACC152992

Grip Drink Holder (20085) \$15.50

Part: ACC152993

TRU-Balance 3 Hydration System (B) \$237

Part: ACC125010

TRU-Balance 3 Transfer Handle (D) \$105 ea.

Part: MEC2205680. HCPCS: E0705

Not available with Heavy Duty armrests.

TRU-Balance 3 Cup Holder..... No Charge

Part: ACC125003

For use with TRU-Balance 3 (track style) armrests.

Cup Holder (Tube-style Armrest)..... No Charge

Part: ACCASMB2634

For use with tube-style armrests.

TRU-Balance 3 Personal Item Hook \$69

Part: ACC125008

Not available with back cane option.

TRU-Balance 3 Reflector Kit..... \$53

Part: ACC124670

TRU-Balance 3 Backpack Holder (Pair) (A) \$69

Part: ACC147774

Not available with back cane option.

TRU-Balance 3 Front Mounted \$69

Cell Phone/ Media Holder, X-Grip (E)

Part: ACC155835

Not available with Heavy Duty armrests.

TRU-Balance 3 Rear View Mirror \$69

Part: ACC155833

Not available with Heavy Duty armrests.

TRU-Balance 3 Oxygen Holder \$210

Part: ACC125006. HCPCS: E2208

Includes accessory bar, a credit will be applied PTOIN-DV3769.

TRU-Balance 3 Clothing Guard (Pair) \$237

Part: ACC125007

Not available with Heavy Duty armrests.


TRU-Balance 3 Glove Box (C) \$158

Part: ACC125011

Stealth Large Tablet Holder\$400

Part: ACC141818

Not available with Heavy Duty armrests.



Accessories section continued on following page.

22. ACCESSORIES CONTINUED

Trilogy Vent Tray \$420
 Part: ACC136791. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOINDV3769. Contact custom for models other than 100 and 200.

Compact Vent Tray \$1,465
 Part: ACC133889. HCPCS: E1029

Push Handles (pair) \$195
 Part: ACC2005027
 Not available with oxygen holder, back cane option, or vent tray.

Accessory Bar \$132
 Part: XRF000128
 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.

Large, Black, Weather Cover \$202
 Part: ACCCOVR1021

Flip-Back Armrest, Side Track Accessory Bracket \$11
 Part: MEC1809178
 Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.

Q-Logic EX/NE Programmer \$573
 Part: ELEASMB6976

VOCSN Vent Tray for Accessory Bar Mount \$940
 Part: ACC2000166. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOINDV3769.

Personal Item Pouch \$84
 Part: ACC2101017

Depth Adjustable Vent Tray \$595
 Coming Soon. Part: ACC2302500. HCPCS: E1029
 Includes accessory bar, a credit will be applied PTOINDV3769.

23. BRANCH PATCHES

Optional.

Air Force Patch \$21
 Part: ACC140240

Army Patch \$21
 Part: ACC140241

Coast Guard Patch \$21
 Part: ACC140239

Marines Patch \$21
 Part: ACC140243

Navy Patch \$21
 Part: ACC140244

Purple Heart Patch \$21
 Part: ACC140238

**24. POSITIONING COMPONENTS**Optional. More options are available on the [Positioning Components Order Form](#).**24.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q) \$632
 Part: ST-ULS4042000-Q. HCPCS: E1028
 Only available with Redesigned back assemblies.

Swing-away, Rail Mounts, pair (TWBTLTC) \$535
 Part: INDPART2677. HCPCS: E1028
 Not available with Redesign back assemblies.

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$238
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$238
 Part: INDPART2664. HCPCS: E0956

Flat Pads, 3"x5", pair (TWBLP3X5) \$238
 Part: INDPART2676. HCPCS: E0956

Flat Pads, 4"x6", pair (TWBLP4X6) \$238
 Part: INDPART2675. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

24.2. Stealth UniLink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable & Removable Hardware . \$252.21 ea.
 QTY Part: ST-UTTP5-Q. HCPCS: E1028

UniLink, 7" Height Adjustable & Removable Hardware . \$262.50 ea.
 QTY Part: ST-UTTP7-Q. HCPCS: E1028

2. UniLink Pelvic/Thigh Pads

UniLink, 4" x 4" Pad \$101.25 ea.
 QTY Part: ST-UTP544056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 8" Pad \$101.25 ea.
 QTY Part: ST-UTP548056100-Q. HCPCS: E0953/E0956

UniLink, 4" x 12" Pad \$101.25 ea.
 QTY Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$77.20 ea.
 QTY Part: ST-UIL001150A000-Q

UniLink, 2" Link \$77.20 ea.
 QTY Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719



Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**24.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW)\$267.50 ea.
 QTY Part: INDPART2804. HCPCS: E1028

Removable Bracket (LP430)..... \$270.38 ea.
 QTY Part: INDPART2805. HCPCS: E1028

Removable Bracket with 1" Offset (LP430 & LHW-121)\$363.61 ea.
 QTY Part: INDPART2854. HCPCS: E1028

Flip Down with LEFT Full Surface Contact\$267.75 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028

Flip Down with RIGHT Full Surface Contact.....\$267.75 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$27.04
 Part: INDPART2821. HCPCS: K0108
 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.

2 1/4" Adapter Plate (LHW-129)\$54.08
 Part: POS140430. HCPCS: K0108
 Only available with the flip down full surface bracket.

3. Pelvic/Thigh Guide Pad

Short Pad, 4"x4" (TWBADD-4X4)..... \$77.25 ea.
 QTY Part: INDPART2801. HCPCS: E0953/E0956

Medium Pad, 4"x8" (TWBADD-4X8)..... \$77.25 ea.
 QTY Part: INDPART2802. HCPCS: E0953/E0956

Long Pad, 4"x12" (TWBADD-4X12)..... \$77.25 ea.
 QTY Part: INDPART2803. HCPCS: E0953/E0956

4. Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

24.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$146.26 ea.
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$146.26 ea.
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$146.26 ea.
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$146.26 ea.
 Part: ST-IFPS2104-Q. HCPCS: E0960

24.5. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor \$92.70 ea.
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor \$92.70 ea.
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960

24.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware

Facial Lateral for Comfort Series.....\$303.08 ea.
 Right Side (CPS9-R-HW)
 Part: INDPART2647. HCPCS: E1028

Facial Lateral for Comfort Series.....\$303.08 ea.
 Left Side (CPS9-L-HW)
 Part: INDPART2648. HCPCS: E1028

Facial Laterals for Comfort Series \$358.18 ea.
 Left and Right Side (CPS9-B-HW)
 Part: INDPART2649. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Right.....\$303.08
 Part: ST-SAR2A-R-Q. HCPCS: E1028

UniLink Facial Lateral, Swing Away, Left\$303.08
 Part: ST-SAR2A-L-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$468.39
 2 Upper & 1 Lower Right
 Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$468.39
 2 Upper & 1 Lower Left
 Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away\$413.29
 2 Upper
 Part: ST-SAR4A-DBU-Q. HCPCS: E1028

2. Facial Lateral Pads and Accessories


Egg Switch Mount (ESM) (A) \$54.08 ea.
 QTY Part: INDPART2815

Figure 8 Bracket with 13" Rod..... \$275.53 ea.
 QTY and Egg Switch Mount (B)
 Part: ST-F8675-13-Q

Facial Pad, 4".5x2", Cool Core (900M-Q) . \$57.68 ea.
 QTY Part: ST-900M-Q. HCPCS: K0108

Facial Pad, 6"x2", Cool Core (920M-Q)....\$57.68 ea.
 QTY Part: ST-920M-Q. HCPCS: K0108

Swing-Away Joystick Mount, 18" Rod (C) \$262.50 ea.
 QTY Part: ST-SUS18-Q
 Allows mounting of facial pads and attaches to the headrest

**24.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube\$189.26 ea.
 Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Left, UniLink, Fixed, Elbow Support Bracket for Flip-Back, . \$200.85
 8"-16"H Armrests (UETB4021FL-Q)
 Part: ST-UETB4021FL-Q

Right, Swing-Away, Elbow Support Bracket, Tube\$189.26 ea.
 Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

Right, UniLink, Fixed, Elbow Support Bracket for Flip-Back, \$200.85
 8"-16"H Armrests (UETB4021FR-Q)
 Part: ST-UETB4021FR-Q

Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$54.08 ea.
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$54.08 ea.
 QTY Part: POS124067. HCPCS: K0108

24.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H\$513.71
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H\$513.71
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H\$513.71
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

24.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)\$82.40 ea.
 QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213)\$82.40 ea.
 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212)\$82.40 ea.
 QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

25. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

25.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit\$1,980
 Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard
 Part: CTL167726

Pretty in Pink\$63
 Part: DGN137729-PPK3

Lemon Crush\$63
 Part: DGN137729-PYL4

Red Corvette\$63
 Part: DGN137729-PRD1

Evolution Orange\$63
 Part: DGN137729-POR7

Green Machine\$63
 Part: DGN137729-PGN11

Electric Blue\$63
 Part: DGN137729-PBU9

Deep Purple\$63
 Part: DGN137729-PPR5

Military Green\$63
 Part: DGN137729-PGN8

Back in Black\$63
 Part: DGN137729-PBK1

Grey Street\$63
 Part: DGN137729-PSL3

Night Train\$63
 Part: DGN137729-PBK3

Midnight Blue\$63
 Part: DGN137729-PBU10

25.2. SCIM Module

Optional.

☐ SCIM Module for Q-Logic 3\$525
 Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

25.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ\$840
 Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

25.4. Head Array

i-Drive 4.0 Bluetooth Tri-array Series,\$4,546.42
 3 Switch with Egg Switch
 Part: ST-IDHBT200-U-Q. HCPCS: E2330

Junior, Curved, i-Drive Tri-Array with\$4,546.42
 Mounting Hardware
 Part: ST-IDHBT450-U-Q. HCPCS: E2330

i-Drive 4.0 Bluetooth Pro Series,\$4,594.83
 3 Switch with Egg Switch
 Part: ST-IDHBT300-U-Q. HCPCS: E2330

Adult, i-Drive Pro Series Head Array\$4,594.83
 (IDHBT300-1)
 Part: ST-IDHBT300-1-Q. HCPCS: E2330

Adult, Curved, i-Drive Tri-Array with\$4,546.42
 Mounting Hardware
 Part: ST-IDHBT400-U-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

Alternative controls and electronics section continued on following page.

25. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**25.5. Sip and Puff Drive****1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3.....\$1,990
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.

**2. Sip and Puff Straw and Mount**

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series,\$303.08
 Left (SP625M-CL)
 Part: ST-SP625M-CL-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive,\$303.08
 Head Array (SP625-U)
 Part: ST-SP625-U-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Stealth Sip and Puff for Comfort Series,\$303.08
 Right (SP625M-CR)
 Part: ST-SP625M-CR-Q. HCPCS: E2326
 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12)..... \$33.06
 Part: ST-SP625-S12-Q

Sip & Puff Bus Cable\$447
 Part: DWR1111H046

25.6. Egg Switches

Red Egg Switch (ESRED)..... \$83.28
 Part: SWTMCRO1077

Green Egg Switch (ESGRN) \$83.28
 Part: SWTMCRO1078

Blue Egg Switch (ESBLU) \$83.28
 Part: SWTMCRO1079

Yellow Egg Switch (ESYEL) \$83.28
 Part: SWTMCRO1080

Black Egg Switch (ESBLK) \$83.28
 Part: SWTMCRO1081

**25.7. Switch Adapter**

2 to 1 Switch Adapter..... \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

25.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
 Part: CTL168057

25.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

