



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Order Form with Contoured Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

- | | |
|---|---|
| ☐ Quantum Q6 Edge 3 3S-SS\$8,870
Part: Q6 EDGE 3 3S-SS. HCPCS: K0848
Solid seat pan standard with no power option. | ☐ Quantum Q6 Edge 3 3S-C\$8,870
Part: Q6 EDGE 3 3S-C. HCPCS: K0849
Contoured seat (captain) standard. |
| ☐ Quantum Q6 Edge 3 3SP-SS\$8,559
Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856
Solid seat pan with single power option for power adjustable seat height. | ☐ Quantum Q6 Edge 3 3SP-C\$8,559
Part: Q6 EDGE 3 3SP-C. HCPCS: K0857
Contoured seat (captain), single power option for power adjustable seat height. |

2. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7 mph, 4 Pole Drive Motors \$1,653
Part: MOT1902455

Requires Q-Logic 3 electronics.

7 mph, 4 Pole, Accu-Trac Motor Technology \$1,360
Part: MOT2008501. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette
Part: DGN1801101-PRD1

 Back in Black
Part: DGN1801101-PBK1

 Night Train
Part: DGN1801101-PBK3

 Grey Street
Part: DGN1801101-PSL3

 Deep Purple
Part: DGN1801101-PPR5

 Raspberry Beret
Part: DGN1801101-PPK4

 Pretty in Pink
Part: DGN1801101-PPK3

 Evolution Orange
Part: DGN1801101-POR7

 Midnight Blue
Part: DGN1801101-PBU10

 Electric Blue
Part: DGN1801101-PBU9

 Blue Jean
Part: DGN1801101-PBU11

 Lemon Crush
Part: DGN1801101-PYL4

 Green Machine
Part: DGN1801101-PGN11

 Military Green
Part: DGN1801101-PGN8

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$263
Part: DGN1801101-WCM7

☐  Faded Glory \$263
Part: DGN1801101-WFL2

☐  Carbon Fiber \$263
Part: DGN1801101-WCB1

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. *Please verify the appropriate billing code to be used for this item with all third-party payors.

12" Power Lift for Captain Seat with iLevel \$3,725
Part: PTOINDV3857. HCPCS: E2298/E2300^
For use with K0856 and K0857 base code.

5. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights No Charge
Part: DGN154117

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129

Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21
Part: MEC1901231

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

7. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

7.1. Drive Wheels

3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair)

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2

4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367

Drive wheels section continued on following page.

7. DRIVE WHEELS CONTINUED**7.2. Drive Wheel Accent Color**

Not available with knobby tires.

Back in Black.....Standard
Part: DGN164245-PBK1 x2Red Corvette.....No Charge
Part: DGN164245-PRD1 x2Night Train.....No Charge
Part: DGN164245-PBK3 x2Grey Street.....No Charge
Part: DGN164245-PSL3 x2Deep Purple.....No Charge
Part: DGN164245-PPR5 x2Raspberry Beret.....No Charge
Part: DGN164245-PPK4 x2Pretty in Pink.....No Charge
Part: DGN164245-PPK3 x2Evolution Orange.....No Charge
Part: DGN164245-POR7 x2Midnight Blue.....No Charge
Part: DGN164245-PBU10 x2Electric Blue.....No Charge
Part: DGN164245-PBU9 x2Blue Jean.....No Charge
Part: DGN164245-PBU11 x2Lemon Crush.....No Charge
Part: DGN164245-PYL4 x2Green Machine.....No Charge
Part: DGN164245-PGN11 x2Military Green.....No Charge
Part: DGN164245-PGN8 x2Military Digital Camo.....\$58 ea.
Part: DGN164245-WCM7 x2Faded Glory.....\$58 ea.
Part: DGN164245-WFL2 x2Carbon Fiber.....\$58 ea.
Part: DGN164245-WCB1 x2**8. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery.....\$355 ea.
Part: BAT2004376. HCPCS: E2361NF-22, Gel Cell, MK Battery.....\$432 ea.
Part: BATGEL1008. HCPCS: E2361NF-22, AGM, MK Battery.....\$432 ea.
Part: BAT121346. HCPCS: E2361**9. PATIENT'S WEIGHT****REQUIRED.**Up to 200 lb User Weight
Part: PTOINDV3837201 lb-250 lb User Weight
Part: PTOINDV3861251 lb-300 lb User Weight
Part: PTOINDV3984**10. EXTENDED SEAT TOWERS****Optional. Not available with Base Mounted Foot Platform.**Extended Seat Towers.....No Charge
Part: MEC2105558**11. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick.....Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick.....\$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.**11.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module.....\$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick.....\$840
Part: CTL167724Q-Logic 3 EX Joystick for Full Lights.....\$840
Part: CTL167725Q-Logic 3e Joystick.....\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display.**11.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder.....\$60
Part: ELE2103635

Electronics section continued on following page.

11. ELECTRONICS CONTINUED**11.3. Q-Logic 3 EX Joystick Armor**

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft TouchStandard Part: DGN124307 Also matches Night Train shrouds.	Raspberry Beret..... \$63 Part: DGN124336-PPK4	Lemon Crush..... \$63 Part: DGN124336-PYL4
Red Corvette..... \$63 Part: DGN124336-PRD1	Pretty in Pink \$63 Part: DGN124336-PPK3	Green Machine \$63 Part: DGN124336-PGN11
Back in Black..... \$63 Part: DGN124336-PBK1	Evolution Orange..... \$63 Part: DGN124336-POR7	Military Green \$63 Part: DGN124336-PGN8
Grey Street..... \$63 Part: DGN124336-PSL3	Midnight Blue \$63 Part: DGN124336-PBU10	Military Digital Camo..... \$63 Part: DGN124336-WCM7
Deep Purple..... \$63 Part: DGN124336-PPR5	Electric Blue \$63 Part: DGN124336-PBU9	Faded Glory \$63 Part: DGN124336-WFL2
	Blue Jean \$63 Part: DGN124336-PBU11	Carbon Fiber \$63 Part: DGN124336-WCB1

12. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

☐ Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	☐ Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
☐ Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	☐ Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
☐ T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	☐ Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
☐ Medium Ball Handle (E) \$69 Part: ACCASMB1859. HCPCS: E2323	☐ Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
☐ Soft Ball Handle (F) \$80 Part: ACCASMB1856. HCPCS: E2323	

**13. JOYSTICK MOUNTING BRACKETS**

Optional. Select a mount to be used with a joystick selected in previous section.

13.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617

13.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032

13.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201043. HCPCS: E1032

Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201044. HCPCS: E1032

13.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.



Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.

Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.

13.2. Joystick Mount Location

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226



14. CONTOURED SEATING**REQUIRED.** Select one seat type.**14.1. Occupied Transit, WC-19****Optional.**

WC-19 Occupied Transit Option for Contoured Seat\$368

Part: KIT122609

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**14.2a. Q-Captain Seat**

16"W x 16"D No Charge

Part: SET2005145

18"W x 20"D No Charge

Part: SET2005151

20"W x 22"D No Charge

Part: SET2005909

16"W x 18"D No Charge

Part: SET2005146

18"W x 22"D No Charge

Part: SET2005908

22"W x 18"D No Charge

Part: SET2005155

16"W x 20"D No Charge

Part: SET2005147

20"W x 16"D No Charge

Part: SET2005152

22"W x 20"D No Charge

Part: SET2005156

18"W x 16"D No Charge

Part: SET2005149

20"W x 18"D No Charge

Part: SET2005153

22"W x 22"D No Charge

Part: SET2005157

18"W x 18"D No Charge

Part: SET2005150

20"W x 20"D No Charge

Part: SET2005154

14.2b. Q-Captain Seat with Solid Seat Pan

16"W x 16"D No Charge

Part: SET2005394

18"W x 20"D No Charge

Part: SET2005399

20"W x 22"D No Charge

Part: SET2005911

16"W x 18"D No Charge

Part: SET2005395

18"W x 22"D No Charge

Part: SET2005910

22"W x 18"D No Charge

Part: SET2005403

16"W x 20"D No Charge

Part: SET2005396

20"W x 16"D No Charge

Part: SET2005400

22"W x 20"D No Charge

Part: SET2005404

18"W x 16"D No Charge

Part: SET2005397

20"W x 18"D No Charge

Part: SET2005401

22"W x 22"D No Charge

Part: SET2005405

18"W x 18"D No Charge

Part: SET2005398

20"W x 20"D No Charge

Part: SET2005402

14.3. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559


Highest Seat to Floor

Part: PTOINDV3560

	Cushion	SSP
High	22.25"	18.25"
Mid	21.5"	17.5"
Low	20.5"	16.5"
Q-Captain Seat to Floor Height		

14.4. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



Contoured seating section continued on following page.

14. CONTOURED SEATING CONTINUED**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit.....\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit.....\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

No Top Foam
Reference part: XRF000382

FRG Top Foam
Reference part: XRF000381

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

No Liner
Reference part: XRF000385

Stealth Polyurethane Liner
Reference part: XRF000384

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Stealth Tek Cover
Reference part: XRF000388

Mesh Cover
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam	Liner	Cover	Top Foam	Liner	Cover	Top Foam	Liner	Cover
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	—	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

14.5. Armrests**Optional.** Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Q-Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat arm pads.

Full Length, Left Armrest.....\$90 ea.
Part: MEC2201304. HCPCS: K0020

Full Length, Right Armrest.....\$90 ea.
Part: MEC2201303. HCPCS: K0020

Contoured seating section continued on following pages.

14. CONTOURED SEATING CONTINUED**1b. Drop-In, Height Adjustable Armrests (T-Arm)**

Includes flat armpads.

☐ ☐ Full Length, Drop-In Armrests \$79 ea.
☐ ☐ Left part: XRF000298. Right part: XRF000299. HCPCS: E0973

☐ ☐ Desk Length, Drop-In Armrests \$79 ea.
☐ ☐ Left part: XRF000300. Right part: XRF000301. HCPCS: E0973
1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$79 ea.
☐ ☐ Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests \$79 ea.
☐ ☐ Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973
1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H \$90 ea.
☐ ☐ Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973

☐ ☐ Full Length, Bariatric Armrests, 10"-14"H \$90 ea.
☐ ☐ Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973
1e. Heavy Duty, Height Adjustable, Removable Armrests

Select armrests and armpads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.
☐ ☐ Left part: XRF000253. Right part: XRF000254. HCPCS: E0973

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
☐ ☐ Left part: XRF000257. Right part: XRF000258. HCPCS: E0973
1f. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

☐ ☐ Pediatric, Removable Armrests \$79 ea.
☐ ☐ Left part: ARM142876. Right part: ARM142877. HCPCS: E0973
14.6. Q-Captain Armrest Spacer

Optional.

 Armrest Spacer \$35
 Part: XRF000370
 1" spacer for Q-Captain, Flip-Up, Height Adjustable Armrests only.
14.7. Armrest Height

REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.

1. Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

14.8. Armpads

REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.

1a. Straight Armpads for Single Post Armrests
 Straight, Full Length Left and Right No Charge
 Part: ARM134275
 Straight, Full Length Left, Desk Length Right No Charge
 Part: ARM134281

 Straight, Desk Length Left, Full Length Right No Charge
 Part: ARM134280
 Straight, Desk Length Left and Right No Charge
 Part: ARM134254
1b. Waterfall Armpads for Single Post Armrests
 Waterfall, Full Length Left and Right No Charge
 Part: ARM134285
 Waterfall, Full Length Left, Desk Length Right No Charge
 Part: ARM134289

 Waterfall, Desk Length Left, Full Length Right No Charge
 Part: ARM134288
 Waterfall, Desk Length Left and Right No Charge
 Part: ARM134282
1c. Falcon Gel Armpads for Single Post Armrests
 Falcon Gel, Full Length Left and Right \$311.22
 Part: ARM136438
 Falcon Gel, Full Length Left, Desk Length Right \$287.93
 Part: ARM136442

 Falcon Gel, Desk Length Left, Full Length Right \$287.93
 Part: ARM136439
 Falcon Gel, Desk Length Left and Right \$265.02
 Part: ARM136437
1d. Armpads for Bariatric or Heavy Duty Armrests
 Bariatric Armpads No Charge
 Part: UPSASMB1027
 Only for use with Bariatric Armrests.
 Straight Armpads No Charge
 Part: UPSASMB1021
 Only for use with Heavy Duty Armrests

 Waterfall Armpads No Charge
 Part: UPSASMB1024
 Only for use with Heavy Duty Armrests

Contoured seating section continued on following pages.

14. CONTOURED SEATING CONTINUED**14.9. Quad Release Lever****Optional. Available with Single Post and Pediatric armrests.**Quad Release Lever, Pair..... No Charge
Part: FRMASMB10517**15. LAP BELTS****REQUIRED.** Select one lap belt below.**15.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

15.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$181 Part: ST-IFRB5350	Stealth Positioning, 60" Single Pull, Push Button,\$50 Extra Large Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning, 60" Single Pull, Push Button, Medium.....\$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning, 60" Single Pull, Push Button,\$50 Non-padded, 2" Belt (IFPB5310) Part: ST-IFPB5310-Q
Stealth Positioning, 60" Single Pull, Push Button, Large\$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	

16. LEGRESTS**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**16.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard Part: FRMASMB16077	Omit base mounted foot platform No Charge Part: PTOINDV3674
Bell Shaped, Base Mounted Foot Platform \$247 Part: RIGASMB7120041	

16.1b. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**Basic Center Mount.....No Charge
Part: XRF000323**2. Basic Center Mount Extension and Angle Adjustable Footplates**

Short Extension, 5"x9" Footplates \$180 Part: XRF000312. HCPCS: K0040	Long Extension, 5"x9" Footplates \$180 Part: XRF000314. HCPCS: K0040
Medium Extension, 5"x9" Footplates \$180 Part: XRF000313. HCPCS: K0040	

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount..... \$158 Part: RIG2104369	Large Calf Pads, Basic Center Mount \$158 Part: RIG2202997 Only available with Long Extension.
--	--

4. Basic Center Mount Additional AccessoryBasic Center Mount Footplate Connector \$16
Part: RIG2104240**16.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145 Part: PTOINDV3754. HCPCS: E1012+K0040 Only available with power lift option.	Quick Adjustable Center Mount Foot Platform..... \$263 Part: XRF000311. HCPCS: K0040
---	---

2. Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x8" Footplates Part: RIG2404017
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Medium Extension, 5"x11" Footplates Part: RIG2404018	Long Extension, 5"x11" Footplates Part: RIG2404019
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x6" Footplates Part: RIG151505	

Contoured seating section continued on following pages.

16. LEGRESTS CONTINUED

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Small Calf Pads, Manually Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

Large Calf Pads, Manually Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H \$184
Part: FRMASMB7110007

Large Calf Panel, 24"Wx8"H \$237
Part: RIG123824
Only available on widths of 20" or greater.

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.

Width Extension for 8" \$69
Part: RIG2404919

Width Extension for 11" \$69
Part: RIG2404920

Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.

Front and Rear Cap \$16
Part: RIG2404923
Available with 5"x8" footplates.

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-Q

Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q

Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q

Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q

Footplate Connector \$16
Part: RIGASMB7120023

Footplate Riser \$75
Part: XRF000348

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel \$50
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

16.1d. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
Pair part: XRF000082. Left: XRF000083. Right: XRF000084
Not available with 11"-13" extensions.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
Pair part: XRF000085. Left: XRF000086. Right: XRF000087
For use only with 11"-13" lower extensions.

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Extruded, 13" - 16.5"	!	✓	✓
Extruded, 15.5" - 19.5"	!	!	!
Angle Adjustable, 13" - 16.5"	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	!	!	!
Multi-Axis, 15.5" - 19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

☒ This rigging configuration does not work with this seat-to-floor height.

2. Lower Extension

☐ ☐ ☐ 11"-13" Length
Pair part: XRF000088.
Left: XRF000089.
Right: XRF000090.

☐ ☐ ☐ 13"-16.5" Length
Pair part: XRF000094.
Left: XRF000095.
Right: XRF000096.

☐ ☐ ☐ 15.5"-19.5" Length
Pair part: XRF000097.
Left: XRF000098.
Right: XRF000099.

Contoured seating section continued on following pages.

16. LEGRESTS CONTINUED**3. Angle Adjustable Footplate(s)**

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 Part: INDPART2823. HCPCS: K0108

16.1e. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

P L R Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2a. Lower Extension

P L R 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

P L R 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

P L R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

P L R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

P L R 7.25"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

P L R 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

P L R 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

P L R 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

P L R 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

P L R 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

16.2. Legrest Accessories

Heel Loops (pair) \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper \$37
 Part: XRF000100

Legrests section continued on following pages.

16. LEGRESTS CONTINUED**16.3. Set Legrest Extension Length****REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.****1a. Centermount Extension Length**

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

17. POWER LIFT ELECTRONICS**REQUIRED with 12" Power Lift for Captain Seat with iLevel.**Single Actuator Through Controller and Interface \$2,163
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and Interface.....\$3,040
Part: XRF000336. HCPCS: E2311Single Actuator Through Controller and Interface (AAM).....\$3,040
Part: XRF000068. HCPCS: E2310**18. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.Cup Holder, Q-Captain No Charge
Part: ACC2006355
For use with Flip-Up armrests.Saddle Bag \$40
Part: DWR1020D010Cup Holder, Bariatric Arms or Gel Armpads \$21
Part: ACCASMB2632Large, Black, Weather Cover \$202
Part: ACCCOVR1021Cup Holder, Tube Style Armrests..... \$32
Part: XRF000221
For HD Removable, Single Post, or Pediatric Arms.Quantum Backpack \$168
Part: ACC1905853Oxygen Holder.....\$118
Part: ACC2006353. HCPCS: E2208
Not available with Power Lift optionGrip Lap Board (20005)..... \$38.95
Part: ACC152992Oxygen Holder..... \$180
Part: ACC2400456. HCPCS: E2208
Select with Power Lift option.Grip Drink Holder (20085) \$15.50
Part: ACC152993Cane/Crutch Holder \$110
Part: ACC2006352. HCPCS: E2207
Not available with Power Lift option.Reflector Kit \$53
Part: ACC124670Cane/Crutch Holder \$250
Part: ACC2400455. HCPCS: E2207
Select with Power Lift option.Clothing Guard (Pair)..... \$237
Part: ACC125007Rear Basket \$118
Part: ACC2006351
Not available with Power Lift option.

Available with Flip-Up armrests only. Not available with Elevating Legrests.

Transfer Handles \$105 ea.
Part: MEC2205680. HCPCS: E0705Rear Basket \$124
Part: ACC2400454
Select with Power Lift option.

Available with Flip-Up armrests only. Not available with Elevating Legrests.

Privacy Bag..... \$42
Part: ACC2005723Personal Item Pouch \$84
Part: ACC2101017Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

20. BRANCH PATCHES**Optional.**

Air Force Patch\$21
Part: ACC140240

Army Patch\$21
Part: ACC140241

Coast Guard Patch\$21
Part: ACC140239

Marines Patch\$21
Part: ACC140243

Navy Patch\$21
Part: ACC140244

Purple Heart Patch\$21
Part: ACC140238

**21. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**21.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
Removable Mounting Hardware
Part: ST-UTTP5-Q, HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.
Removable Mounting Hardware
Part: ST-UTTP7-Q, HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q, HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q, HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q, HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
Part: PTOINDV3585

Install on Right
Part: PTOINDV3586

Do Not Install
Part: PTOINDV3719

21.2. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube\$199 ea.
Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube\$199 ea.
Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
Part: POS124067. HCPCS: K0108

21.3. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954

21.4. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q, HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q, HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q, HCPCS: E0951+E0952

22. LAP TRAY

22.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware.....\$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	M	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SML	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

