

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com**Q6 Edge® 3 Order Form with Contoured, Synergy®, and TRU-Balance® Flex Seating**

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____

PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____



TRU-Balance® 3 Power Lift Seating Shown

INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


Thank you for choosing Quantum!**PATIENT INFORMATION**

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____

**1. BASE MODEL**

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 3 3S-SS \$7,895

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 3 3SP-SS \$7,895

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option for power elevating seat.

☐ Quantum Q6 Edge 3 3S-C \$7,895

Part: Q6 EDGE 3 3S-C. HCPCS: K0849

Contoured seat (captain) standard.

☐ Quantum Q6 Edge 3 3SP-C \$7,895

Part: Q6 EDGE 3 3SP-C. HCPCS: K0857

Contoured seat (captain), single power option for 12" lift.

Account Number: _____ Date: _____ Marked for: _____

2. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Patient's Weight: _____

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1. Match Colors

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

3.2a. Standard Shrouds

 Red Corvette Part: DGN1801101-PRD1	 White Satin Part: DGN1801101-PWH2	 Midnight Blue Part: DGN1801101-PBU10	 Military Green Part: DGN1801101-PGN8
 Cherry Bomb Part: DGN1801101-PRD4	 Deep Purple Part: DGN1801101-PPR5	 Electric Blue Part: DGN1801101-PBU9	 Sugar Plum Part: DGN1801101-PPR7
 Back in Black Part: DGN1801101-PBK1	 Raspberry Beret Part: DGN1801101-PPK4	 Blue Jean Part: DGN1801101-PBU11	To Be Determined Follow-up required.
 Night Train Part: DGN1801101-PBK3	 Pretty in Pink Part: DGN1801101-PPK3	 Lemon Crush Part: DGN1801101-PYL4	
 Grey Street Part: DGN1801101-PSL3	 Evolution Orange Part: DGN1801101-POR7	 Green Machine Part: DGN1801101-PGN11	

3.2b. Hydrographic Shrouds

☐  Military Digital Camo \$250 Part: DGN1801101-WCM7	☐  Faded Glory \$250 Part: DGN1801101-WFL2	☐  Carbon Fiber \$250 Part: DGN1801101-WCB1
--	--	---

4. FENDER LIGHTS

Optional. Fender lights operate by pressing included single switch button.

☐ **Switch Fender Lights** No Charge
Part: DGN154117

5. FULL LIGHT PACKAGES

Optional. Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package. Select one light package below.

- ☐ **Base Mounted, LED, Full Light Package** \$995
Part without iLevel: ELE1705999. Part with iLevel: ELE168643
Fender lights will be omitted with Base Mounted, LED, Full Light Package.
- ☐ **Seat Mounted, LED, Full Light Package** \$995
Flex part: ELE131477. Synergy part: ELEASMB5719
Not available with contoured seating.



6. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard
Part: MOT147433

Accu-Trac Motor Technology \$1,295
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7. BASE MOUNTED FOOT PLATFORM

REQUIRED. May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077

Omit base mounted foot platform No Charge
Part: PTOINDV3674

Bell Shaped, Base Mounted Foot Platform \$235
Part: RIGASMB7120041

8. BATTERIES

REQUIRED. The powerchair requires 2 batteries and will be installed if ordered.

☐ **NF-22, Gel Cell, Interceptor Battery** \$310 ea.
Part: BAT1704613. HCPCS: E2361

☐ **NF-22, AGM, MK Battery** \$395 ea.
Part: BAT121346. HCPCS: E2361

☐ **No batteries.** Provider will assume responsibility of batteries.
Part: Configuration selection

☐ **NF-22, Gel Cell, MK Battery** \$395 ea.
Part: BATGEL1008. HCPCS: E2361

9. DRIVE WHEELS**REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**9.1. Drive Wheels**

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires No Charge Part: WHL1713925 x2
4"x14" Knobby, Flat Free Black Tires \$275.50 ea. Part: WHL147601 x2	4"x14" Knobby, Pneumatic Black Tires \$275.50 ea. Part: WHL147599 x2

9.2. Drive Wheel Accent Color

Not available with knobby tires.

Back in Black Standard Part: DGN164245-PBK1 x2	Raspberry Beret No Charge Part: DGN164245-PPK4 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette No Charge Part: DGN164245-PRD1 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Cherry Bomb No Charge Part: DGN164245-PRD4 x2	Evolution Orange No Charge Part: DGN164245-POR7 x2	Sugar Plum No Charge Part: DGN164245-PPR7 x2
Night Train No Charge Part: DGN164245-PBK3 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Military Digital Camo \$50 ea. Part: DGN164245-WCM7 x2
Grey Street No Charge Part: DGN164245-PSL3 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Faded Glory \$50 ea. Part: DGN164245-WFL2 x2
White Satin No Charge Part: DGN164245-PWH2 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$50 ea. Part: DGN164245-WCB1 x2
Deep Purple No Charge Part: DGN164245-PPR5 x2	Lemon Crush No Charge Part: DGN164245-PYL4 x2	To Be Determined Follow-up required.

10. SEAT TOWERS/POWER ADJUSTABLE SEAT HEIGHT**REQUIRED.** Select one option below.**10.1a. Power Base Seat Towers**

Standard Seat Towers Standard Part: Configuration selection, no part number	Extended Seat Towers (increased seat to floor height) ... No Charge Part: MEC120301
--	--

10.1b. Power Adjustable Seat Height

The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

12" Power Lift for Captain Seat with iLevel \$2,600 Part: PTOINDV3857. HCPCS: E2300 For use with K0857 base code. Must select Power Lift Captain Seat in contoured seat section. Includes armrests and lap belt. Available with swing away legrests or manually adjustable center mount foot platform.	Power Elevating Seat with 7" Stroke \$1,895 Contoured seat part: FRMASMB16310. HCPCS: E2300 Comfort contoured seat part: MEC133630. HCPCS: E2300 Synergy seat part: FRMASMB16312. HCPCS: E2300
	Power Elevating Seat with 6" Stroke \$1,895 Contoured seat part: FRMASMB16309. HCPCS: E2300 Comfort contoured seat part: MEC133629. HCPCS: E2300 Synergy seat part: FRMASMB16311. HCPCS: E2300

11. ELECTRONICS**REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

☐ NE 4-Key Joystick Standard Part: CTLDC1554 Not available with power seating.	☐ NE+ 6-Key Joystick (2 Actuators) \$200 Part: CTLDC1560 Control up to 2 actuators through joystick.
☐ NE+ 6-Key Joystick (1 Actuator) No Charge Part: CTLDC1560 Control of one actuator through joystick on single power base model (3SP-SS) only.	☐ NE+ 8-Key Joystick (2 Actuators) \$400 Part: CTLDC1561 Control up to 2 actuators through joystick and full light package.
	☐ No Joystick, Setup for NE+ Joystick No Charge Part: Configuration selection

**11.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

☐ Q-Logic 3 EX Controller Module \$695 Part: See XRef. HCPCS: E2377	☐ Q-Logic 3 Expandable Controller Harness \$425 Part: HAR144005. HCPCS: E2313
---	---

2. Q-Logic 3 Joystick

☐ Q-Logic 3 EX Joystick \$700 Part: CTL167724	☐ Q-Logic 3 EX Joystick for Full Lights \$700 Part: CTL167725
☐ No Joystick, Setup for Q-Logic Part: Configuration selection.	



Electronics section continued on following page.

11. ELECTRONICS CONTINUED**11.2. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft Touch Standard
 Part: DGN124307
 Also matches Night Train shrouds.

Non-matching Armor \$50
 Part: PTOINDV3422
 Select color below

Match Armor to Base Color \$50
 Part: PTOINDV3421
 Part below will be selected with system configurator.

Red Corvette Part: DGN124336-PRD1	White Satin Part: DGN124336-PWH2	Evolution Orange Part: DGN124336-POR7	Lemon Crush Part: DGN124336-PYL4	Military Digital Camo Part: DGN124336-WCM7
Cherry Bomb Part: DGN124336-PRD4	Deep Purple Part: DGN124336-PPR5	Midnight Blue Part: DGN124336-PBU10	Green Machine Part: DGN124336-PGN11	Faded Glory Part: DGN124336-WFL2
Back in Black Part: DGN124336-PBK1	Raspberry Beret Part: DGN124336-PPK4	Electric Blue Part: DGN124336-PBU9	Military Green Part: DGN124336-PGN8	Carbon Fiber Part: DGN124336-WCB1
☒ Night Train* *See Black, Soft Touch	Pretty in Pink Part: DGN124336-PPK3	Blue Jean Part: DGN124336-PBU11	Sugar Plum Part: DGN124336-PPR7	To Be Determined Follow-up required.
Grey Street Part: DGN124336-PSL3				

12. JOYSTICK HANDLES**Optional.** The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.


- | | |
|--|---|
| ☐ Mushroom Handle (B) \$85
Part: ACCASMB1860. HCPCS: E2323 | ☐ Extended I-Handle (G) \$85
Part: ACCASMB1857. HCPCS: E2323 |
| ☐ Chin Cup (C) \$85
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | ☐ Bodypoint 3" Goal Post Handle \$99
Part: PLSKNCP1057. HCPCS: E2323 |
| ☐ T-Bar (D) \$120
Part: ACCASMB1858. HCPCS: E2323 | ☐ Bodypoint 4" Goal Post Handle \$99
Part: PLSKNCP1058. HCPCS: E2323 |
| ☐ Medium Ball Handle (E) \$65
Part: ACCASMB1859. HCPCS: E2323 | ☐ Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| ☐ Soft Ball Handle (F) \$65
Part: ACCASMB1856. HCPCS: E2323 | |

**13. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


13.1a. Fixed, In-line Joystick Bracket

- ☐ **Fixed, In-line, Bracket** Standard
 Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744. SAJ Part: FRMASMB10454
- ☐ **Left Side Mount**
 Part: PTOINDV1225
- ☐ **Right Side Mount**
 Part: PTOINDV1226

**13.1b. Swing-Away (Retractable) Joystick Bracket**


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

- ☐ **Swing-away, Inline, Left Side Mount** \$245
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, 1" Offset, Left Side Mount** \$245
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, Inline, Right Side Mount** \$245
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, 1" Offset, Right Side Mount** \$245
 Part: See XRef. HCPCS: E1028

**13.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

- ☐ **Retract4 Mount** \$350
 NE/NE+ Part: MEC138876. HCPCS: E1028
 Q-Logic Part: MEC138860. HCPCS: E1028
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.
- ☐ **Stealth Swing-Away/Flip Down Bracket** \$294.25
 NE/NE+ Part: ST-TWBM480QL-Q. HCPCS: E1028
 Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028
 SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028
 Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.
- Multi-axis Mount** \$100
 NE/NE+ Part: ACC142528. HCPCS: E1028
 Q-Logic Part: ACC131216. HCPCS: E1028
 SAJ Part: ACC142740. HCPCS: E1028
 Allows height adjustment of joystick.

**Specialty Mount Side**

- ☐ **Left Side Mount**
 Part: PTOINDV1225
- ☐ **Right Side Mount**
 Part: PTOINDV1226

Account Number: _____ Date: _____ Marked for: _____

14. CONTOURED SEATING

Optional. Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

14.1. Occupied Transit, WC-19

WC-19 Occupied Transit Option for Contoured Seat \$350
Part: KIT122609

Only available with Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

14.2a. Captain Style Seating

1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline

Comfort seat, full length armrests must be selected with a depth adjustable high back.

16"Wx16"-18"D w/ 60" Lap Belt \$150	18"Wx18"-20"D w/ 60" Lap Belt \$150	20"Wx18"-20"D w/ 60" Lap Belt \$150
Part: SET131843	Part: SET131844	Part: SET131845

1b. Depth Adjustable High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

18"Wx18"-20"D with 60" Lap Belt No Charge	22"Wx20"-22"D with 70" Lap Belt \$150
Part: SETASMB6895	Part: SETASMB6897
20"Wx18"-20"D with 60" Lap Belt \$150	
Part: SETASMB6896	

1c. Power Lift, Fold Flat Captain Seat

Only available with 12" Power Lift for Captain Seat with iLevel. Includes armrests and 60" lap belt.

16"Wx16"-18"D \$150	18"Wx18"-20"D \$150	20"Wx20"-22"D \$150
Part: SET1701577	Part: SET1701578	Part: SET1701579

1d. High Back with Limited Recline

16"Wx16"D \$150	18"Wx16"D \$150	20"Wx22"D \$150
Part: SETASMB2770	Part: SETASMB2775	Part: SETASMB4889
16"Wx18"D \$150	18"Wx22"D \$150	
Part: SETASMB2771	Part: SETASMB4888	
16"Wx20"D \$150	20"Wx16"D \$150	
Part: SETASMB2772	Part: SETASMB2780	

14.2b. Contoured Seat with Solid Seat Pan

1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline and Solid Seat Pan

Comfort seat, full length armrests must be selected with a depth adjustable high back.

16"Wx16"-18"D with 60" Lap Belt \$150	18"Wx18"-20"D with 60" Lap Belt \$150
Part: SET145814	Part: SET145815

1b. Power Lift, Fold Flat Contoured Seat with Solid Seat Pan

Only available with 12" Power Lift with iLevel. Includes armrests and lap belt.

16"Wx16"-18"D w/ 60" Lap Belt \$150	18"Wx18"-20"D w/ 60" Lap Belt \$150	20"Wx20"-22"D w/ 60" Lap Belt \$150
Part: SET1701580	Part: SET1701581	Part: SET1701582

1c. High Back with Limited Recline and Solid Seat Pan

16"Wx16"D \$150	18"Wx18"D \$150	20"Wx20"D \$150
Part: SETASMB3178	Part: SETASMB2848	Part: SETASMB3185
16"Wx18"D \$150	18"Wx20"D \$150	20"Wx22"D \$150
Part: SETASMB3179	Part: SETASMB3182	Part: SETASMB5008
16"Wx20"D \$150	18"Wx22"D \$150	
Part: SETASMB3180	Part: SETASMB5007	
18"Wx16"D \$150	20"Wx18"D \$150	
Part: SETASMB3181	Part: SETASMB3184	

14.3. Seat to Floor Height

Required with any contoured seating selection.

Lowest Seat to Floor Part: PTOINDV3558	Middle Seat to Floor Part: PTOINDV3559 Not available with Comfort high back seat.	Highest Seat to Floor Part: PTOINDV3560
---	---	--

14.4. Set Depth Selection for Comfort Seat

Required with Comfort Seat selection above.

14" Depth Part: SETDPTH1004	16" Depth Part: SETDPTH1006	18" Depth Part: SETDPTH1008	20" Depth Part: SETDPTH1010
--------------------------------	--------------------------------	--------------------------------	--------------------------------

14.5. Back Height Selection for Comfort Seat

Required with Comfort Seat selection above.

Set Back Height to Lowest Position Part: SETHGHT1100	Set Back Height to Highest Position Part: SETHGHT1101
---	--

14.6. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: _____ Width: _____ Depth: _____

Contoured seating section continued on following page.

14. CONTOURED SEATING CONTINUED**14.7. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests****REQUIRED except with fold flat, power lift captain seat.** Select a left and right armrest or to omit them.

Full Length, Left Armrest.....\$85 Part: FRMASMB8614. HCPCS: E0973	Comfort Seat, Right, Full Length Armrest.....\$85 Part: FRMASMB16315. HCPCS: E0973
Desk Length, Left Armrest.....\$85 Part: FRMASMB8612. HCPCS: E0973	Comfort Seat, Left, Full Length, Armrest\$85 Part: FRMASMB16316. HCPCS: E0973
Full Length, Right Armrest.....\$85 Part: FRMASMB8615. HCPCS: E0973	No Armrest Standard Part: N/A
Desk Length, Right Armrest\$85 Part: FRMASMB8613. HCPCS: E0973	

15. SYNERGY SEATING SYSTEM**Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.**15.1. Synergy Frame Selection**

Synergy Seating System	\$434
Part: PTOINDV2995	

15.2. Synergy Transit Kit**Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

☐ WC-19 Occupied Transit Option for Synergy Seat.....\$275 Part: KITASMB1830
To determine if your chair is WC-19 compliant, please refer to this website: http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs

15.3. Seat Dimension - Width**REQUIRED with Synergy seating.** Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

15.4. Seat Dimension - Width Package**REQUIRED with Synergy seating.** Maximum seat width is 22".

10"-14" Width Package Part: FRMASMB18011	14"-18" Width Package Part: FRMASMB18012	16"-20" Width Package Part: FRMASMB18013	20"-24" Width Package Part: FRMASMB18014 (Maximum width is 22")
---	---	---	--

15.5. Seat Dimension - Depth**REQUIRED with Synergy seating.** Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

15.6. Seat Dimension - Depth Package**REQUIRED with Synergy seating.** Maximum seat depth is 22".

10"-14" Depth Package Part: FRMASMB6981	14"-18" Depth Package Part: FRMASMB17823	15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package Part: FRMASMB7733 (Maximum depth is 22")
--	---	--	---

15.7. Synergy Seat to Floor**REQUIRED with Synergy seating.**

15.875" Part: PTOINDV3652	16.875" Part: PTOINDV3653	17.875" Part: PTOINDV3654	18.875" Part: PTOINDV3761	19.875" Part: PTOINDV3762	20.875" Part: PTOINDV3763
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

15.8. Synergy Back Options**REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Cane Selection**

Synergy Back Cane No Charge Part: FRMASMB15845	2" Endomorph (1" per side) Back Cane\$195 Part: FRMASMB12115
TRU-Comfort Back Canes No Charge Part: N/A, Configuration selection	4" Endomorph (2" per side) Back Cane.....\$195 Part: FRMASMB12139
Flip Forward Back Cane \$150 Part: FRMASMB15846 Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests.	Manual Recline Back \$750 Part: FRMASMB12910. HCPCS: E1226 Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.
2" Mesomorph (1" per side) Back Cane.....\$195 Part: FRMASMB12067	
4" Mesomorph (2" per side) Back Cane.....\$195 Part: FRMASMB12091	

Synergy seating section continued on following pages.

15. SYNERGY SEATING SYSTEM CONTINUED**2. Synergy Back Cane Height**

Omit this selection if a TRU-Comfort back is desired. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes
Part: FRMASMB15832

15" High Canes
Part: FRMASMB15833

16" High Canes
Part: FRMASMB15834

17" High Canes
Part: FRMASMB15835

18" High Canes
Part: FRMASMB15836

19" High Canes
Part: FRMASMB15837

20" High Canes
Part: FRMASMB15838

21" High Canes
Part: FRMASMB15839

22" High Canes
Part: FRMASMB15840

23" High Canes
Part: FRMASMB15841

24" High Canes
Part: FRMASMB15842

25" High Canes
Part: FRMASMB15843

3. Back Selection

Synergy Back..... \$275
Part: See XRef

TRU-Comfort, 20" Width or less..... \$675
Part: See XRef. HCPCS: E2620

TRU-Comfort, 22" Width or greater..... \$775
Part: See XRef. HCPCS: E2621

Sling Back..... \$98
Part: See XRef

Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back..... No Charge
Part: PTOINDV1287

4. Synergy Overall Back Height

Parts: See XRef, configuration selection with back selection.

12" High 15" High 18" High 21" High 24" High

13" High 16" High 19" High 22" High 25" High

14" High 17" High 20" High 23" High

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	-	-	-	-	-	-	-	-
	16"	S	S	-	-	-	-	-	-	-	-
	17"	-	S	S	S	S	S	S	S	S	S
	18"	-	S	S	S	S	S	S	S	S	S
	19"	-	S	S	S	-	S	S	S	S	S
	20"	-	S	S	S TC	-	S TC	S	S TC	S	S
	21"	-	S	S	S TC	S	S TC	S	S TC	S	S
	22"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	23"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S
	24"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC
	25"	-	S	S	S TC	S TC	S TC	S TC	S TC	S	S TC

S Synergy Back S TC Synergy or TRU-Comfort Back

Synergy and TRU-Comfort Back Availability Chart

15.9. Synergy Armrests

REQUIRED with Synergy seating. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

1a. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ Full Length, Flip Back Armrests..... Standard
L R Left Part: FRMASMB7038. Right part: FRMASMB7039

☐ Desk Length, Flip Back Armrests Standard
L R Left part: FRMASMB7040. Right part: FRMASMB7041

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ Full Length, Two-Post Armrests.....\$150 ea.
L R Left Part: FRMASMB7062. Right part: FRMASMB7063. HCPCS: E0973

☐ Desk Length, Two-Post Armrests\$150 ea.
L R Left part: FRMASMB7060. Right part: FRMASMB7061. HCPCS: E0973

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ Full Length, Heavy Duty Armrests, 8"-12" H \$75 ea.
L R Left Part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973

☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$75 ea.
L R Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$75 ea.
L R Left Part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973

☐ Desk Length, Heavy Duty Armrests, 10"-14"H \$75 ea.
L R Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armrests.

☐ Pediatric, Removable Armrests \$75 ea.
L R Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ Full Length, Cantilever Armrests \$97.50 ea.
L R Left Part: FRMASMB16188. Right part: FRMASMB16189

☐ Desk Length, Cantilever Armrests..... \$97.50 ea.
L R Left part: FRMASMB16190. Right part: FRMASMB16191

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ Full Length, Single Post Armrests \$75 ea.
L R Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ Desk Length, Single Post Armrests \$75 ea.
L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

15.10. Synergy Armpads

REQUIRED with armrest selection above. Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard
Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.
Desk/Full part: UPSASMB1023

Waterfall Armpads..... No Charge
Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.
Desk/Full part: UPSASMB1026

No Armpad Selection No Charge
Part: PTOINDV1373

Synergy seating section continued on following page.

15. SYNERGY SEATING SYSTEM CONTINUED**15.11. Joystick Bracket Receiver****REQUIRED.**

- ☐ Standard Fixed Receiver Standard
Part: FRMASMB16062
- ☐ Stealth Height Adjustable Receiver \$214
Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

16. TRU-BALANCE FLEX SEATING SYSTEM

Optional. Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D. **Not available with power adjustable seat height options.**

16.1. TRU-Balance Flex Frame Selection

TRU-Balance Flex Seating System No Charge
Part: PTOINDV3853

16.2. TRU-Balance Flex Seat Dimension - Width

REQUIRED with TRU-Balance Flex seating. Maximum seat width is 22". Part number is SETWDTH plus number following # below.

15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

16.3. TRU-Balance Flex Seat Dimension - Depth

REQUIRED with TRU-Balance Flex seating. Maximum seat depth is 20". Part number is SETDPTH plus number following # below.

15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)
------------------	------------------	------------------	------------------	------------------	------------------

16.4. TRU-Balance Flex Seat to Floor

REQUIRED with TRU-Balance Flex seating.

16.5" at front and rear Part: PTOINDV3459	17.5" at front with 2° posterior angle Part: PTOINDV3460	18.25" at front with 4° posterior angle Part: PTOINDV3805
--	---	--

16.5. TRU-Balance Flex Back Cane Selection

Angle Adjustable Back Cane No Charge
Part: MEC140266

Flip Forward Back Cane \$150
Part: MEC140272

Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$195
Part: MEC1804736

4" Mesomorph (2" per side) Back Cane \$195
Part: MEC1804739

2" Endomorph (1" per side) Back Cane \$195
Part: MEC141455

4" Endomorph (2" per side) Back Cane \$195
Part: MEC141456

Manual Recline Back \$750
Part: FRMASMB12910. HCPCS: E1226

Select Synergy Back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

16.6. TRU-Balance Flex Back Cane Height Selection

16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes
Part: FRMASMB15832

17" High Canes
Part: FRMASMB15835

20" High Canes
Part: FRMASMB15838

23" High Canes
Part: FRMASMB15841

15" High Canes
Part: FRMASMB15833

18" High Canes
Part: FRMASMB15836

21" High Canes
Part: FRMASMB15839

24" High Canes
Part: FRMASMB15842

16" High Canes
Part: FRMASMB15834

19" High Canes
Part: FRMASMB15837

22" High Canes
Part: FRMASMB15840

25" High Canes
Part: FRMASMB15843

16.7. TRU-Balance Flex Back Options

REQUIRED with TRU-Balance Flex seating. See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

1. Back Selection

Synergy Back Cushion \$275
Part: See XRef

Sport Back Cushion, 20" Width or Less \$375
Part: See XRef

Sport Back Cushion, 22" Width \$575
Part: See XRef

TRU-Comfort Cushion, 20" Width or Less \$675
Part: See XRef. HCPCS: E2620

TRU-Comfort Cushion, 22" Width \$775
Part: See XRef. HCPCS: E2621

Sling Back \$98
Part: See XRef

Overall back height must match back cane height selection above.

Omit Back No Charge
Part: PTOINDV1287

		Back Width							
		15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	16"	-	P	T	P	T	P	T	P
	17"	S	S	P	T	S	P	T	S
	18"	S	S	P	T	S	P	T	S
	19"	S	S	P	T	S	P	T	S
	20"	S	S	P	T	S	P	T	S
	21"	S	S	P	T	S	P	T	S
	22"	S	S	P	T	S	P	T	S
	23"	S	S	P	T	S	P	T	S
	24"	S	S	P	T	S	P	T	S
	25"	S	S	P	T	S	P	T	S

S Synergy Back
 P Sport Back
 T TRU-Comfort Back

TRU-Balance Flex Seat Back Availability Chart

TRU-Balance Flex seating section continued on following pages.

16. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED**2. TRU-Balance Flex Overall Back Height**

Refer to the seat back availability chart for stock order back size ranges.

14" Height Part: SETHGHT1027	17" Height Part: SETHGHT1030	20" Height Part: SETHGHT1033	23" Height Part: SETHGHT1036	26" Height Part: SETHGHT1039
15" Height Part: SETHGHT1028	18" Height Part: SETHGHT1031	21" Height Part: SETHGHT1034	24" Height Part: SETHGHT1037	27" Height Part: SETHGHT1040
16" Height Part: SETHGHT1029	19" Height Part: SETHGHT1032	22" Height Part: SETHGHT1035	25" Height Part: SETHGHT1038	28" Height Part: SETHGHT1041

16.8. TRU-Balance Flex Armrests**REQUIRED with TRU-Balance Flex seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU-Balance Flex Seat, Drop-In Armrests**

Includes full length, flat armpads.

☐ ☐ Full Length, Drop-In Armrests No Charge ☐ ☐ Left part: ARM1803253. Right part: ARM156742	☐ ☐ Desk Length, Drop-In Armrests No Charge ☐ ☐ Left part: ARM1803254. Right part: ARM156743
---	---

1b. Flip Back, Cane Mounted Armrests

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests (pair) \$60 ☐ ☐ Left part: FRMASMB7038. Right part: FRMASMB7039	☐ ☐ Desk Length, Flip Back Armrests (pair) \$60 ☐ ☐ Left part: FRMASMB7040. Right part: FRMASMB7041
--	--

1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests \$150 ea. ☐ ☐ Part: See XRef. HCPCS: E0973	☐ ☐ Desk Length, Two-Post Armrests \$150 ea. ☐ ☐ Part: See XRef. HCPCS: E0973
--	--

1d. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests \$75 ea. ☐ ☐ Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973	☐ ☐ Desk Length, Single Post Armrests \$75 ea. ☐ ☐ Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973
---	---

16.9. TRU-Balance Flex Armpads**REQUIRED with armrest selection above except Drop-In Armrests.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134275
Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right No Charge Part: ARM134254

1b. Waterfall Armpads


Waterfall, Full Length Left and Right No Charge Part: ARM134275	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134275
Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134281	Waterfall, Desk Length Left and Right No Charge Part: ARM134254

1c. Falcon Gel Armpads

Falcon Gel, Full Length Left and Right \$296.40 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right \$274.22 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right \$274.22 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$252.40 Part: ARM136437

17. HEADRESTS**Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****17.1. Headrest Pads**

☐ Stealth 8" Comfort Plus (CP450) \$151.94 Part: POS148484. HCPCS: E0955	☐ Stealth Small Contoured (9.75"Wx3"H) \$64.20 (STL-OTSM250) Part: SETHEAD1052. HCPCS: E0955
☐ Stealth 10" Comfort Plus (CP250) (A) \$151.94 Part: SETHEAD1057. HCPCS: E0955	☐ Stealth Large Contoured (11"Wx6"H) (B) \$80.25 (STL-OTLG150) Part: SETHEAD1053. HCPCS: E0955
☐ Stealth 14" Comfort Plus (CP150) \$151.94 Part: INDPART2327. HCPCS: E0955	
☐ Pediatric Comfort Plus (6"Wx4.5"H) \$151.94 (CPN350) Part: SETHEAD1069. HCPCS: E0955	

**17.2. Headrest Brackets**

Stealth's Worlds Best Removable Headrest Bracket (TWB485-1) (C) \$220.50 Part: FRMASMB15999. HCPCS: E1028	Stealth Multi-Axis Removable (D) \$160.50 Headrest Bracket (HMO475P) Part: HDWASMB1084. HCPCS: E1028 Not available with TRU-Balance Flex seating.
Pediatric Offset (HMO475P-17) (E) \$160.50 Part: HDWASMB1166. HCPCS: E1028	

Headrest section continued on following page.

17. HEADRESTS CONTINUED**17.3. Flip-down Headrest Bracket**

Stealth Flip-down Headrest Bracket (FDM350).....\$224.70
 Part: INDPART2219
 Not available with TRU-Balance Flex seating. The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

18. LAP BELTS

REQUIRED except with select contoured seats with included lap belt. Select one lap belt below.

18.1a. Standard Lap Belts

☐ 38" Lap Belt Part: ACCBELT1006 ☐ 50" Lap Belt Part: ACCBELT1000 ☐ 60" Lap Belt Part: ACCBELT1002 ☐ 70" Lap Belt Part: ACCBELT1003 ☐ 80" Lap Belt Part: ACCBELT1004 ☐ 90" Lap Belt Part: ACCBELT1001

18.1b. Optional Lap Belts

Not available with contoured or TRU-Balance Flex seating.

☐ i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q\$86
☐ i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q\$86
☐ i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q\$86
☐ 60" Sure-Lok Lap Belt (Non-retractable).....\$120
 Required with WC-19 Transit Option, Synergy Seat only.
 Part: ACCBELT1016

19. LATERAL EXTENSIONS

REQUIRED with Synergy Back and thoracic lateral selections.

Lateral Extension for Synergy Back (pair).....\$26.25 ea.
 16"OH back part: FRMASMB16274. 17"OH back part: FRMASMB16280. 20"OH back part: FRMASMB16298. 24"OH back part: FRMASMB16281

20. SYNERGY POWER OPTION ELECTRONICS

REQUIRED for Synergy Seat with Power Legrests or Elevating Seat. Select the desired control type for the seating system.

☐ One Actuator Control through NE+ Joystick\$1,895
 Part: See XRef. HCPCS: E2310
☐ One Actuator Control through Q-Logic 3 Electronics\$1,895
 Part: See XRef. HCPCS: E2310
☐ Two Actuator Control through NE+ Joystick\$2,375
 Part: See XRef. HCPCS: E2311
☐ Multiple Actuator Control through Q-Logic 3 Electronics\$2,895
 Part: See XRef. HCPCS: E2311

21. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

21.1. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured seating.**

1. Foot Platform Assembly

Power Articulating Foot Platform\$2,995
 Part: PTOINDV3754. HCPCS: E1012+K0040
 Not available with TRU-Balance Flex seating.

Manually Adjustable Center Mount Foot Platform...\$250
 12"-16" width part: FRMASMB17944. 14"-18" width part: FRMASMB17945. 16"-20" width part: FRMASMB17946.
 HCPCS: K0040

**2. Lower Extension and Angle Adjustable Footplates**

8"-12" Extension, 5"x6" Footplates
 Part: RIG152577
 Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates
 Part: RIG150206

13"-19" Extension, 5"x6" Footplates
 Part: RIG151505

8"-12" Extension, 5"x8" Footplates
 Part: RIG152578
 Not available with power articulating platform.

10"-16" Extension, 5"x8" Footplates
 Part: RIG150203

13"-19" Extension, 5"x8" Footplates
 Part: RIG151506

10"-16" Extension, 5"x11" Footplates
 Part: RIG150212

13"-19" Extension, 5"x11" Footplates
 Part: RIG151508

Center Mount Foot Platform Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
10"-16" Length Extension	!	✓	✓	☒ No interference issues at this seat-to-floor height. ☐ May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
13"-19" Length Extension	!	!	!	

3. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H\$150*
 Part: RIGASMB7120018

Small Calf Panel 18"Wx8"H\$150
 Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H\$150*
 Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H\$200
 Part: RIG123824
 Only available on widths of 20" or greater.

Legrests section continued on following page.

21. LEGRESTS CONTINUED**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates\$50 Part: RIG152517	4" Outside Edge (pair)\$50 Part: RIG152510 Only available with 6"D footplates and 8"D footplate extension.
Width Extension for 11"D Footplates\$50 Part: RIG152518	6" Outside Edge (pair)\$50 Part: RIG152516 Available with 8"D and 11"D footplates and 11"D footplate extension.
2" Bell-Shaped Length Extension\$50 Part: RIG152519 Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.	

5. Center Mount Additional Accessories

Reversed Dartex Calf Pad Cover (pair) \$150 Small part: RIG133717. Large part: RIG133719	Wheel Bumpers (pair)..... \$12.50 ea. Part: WHLASMB1578 x2
Reversed Dartex Panel Cover (pair) \$150 Small part: RIG134673. Large part: RIG134744	Footplate Connector \$15 Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers..... \$205 Part: ST-FPG5008T-Q. HCPCS: K0108	6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$225 Part: ST-FPG6508T-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers..... \$225 Part: ST-FPG5011T-Q. HCPCS: K0108	7"x11" Tapered Gel Footplate Cover for Width Extension \$250 Part: ST-FPG7011T-Q. HCPCS: K0108

21.2. 70° Power Articulating, Swing Away Legrests

Not available with contoured or TRU-Balance Flex seating seating. 4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Articulating Legrest Upper Assembly

☐ ☐ ☐ 70° Articulating, Swing-Away.....\$1,197.50 ea. Legrest with 14.5"-18" Extension Pair Part: FRMASMB9073. HCPCS for pair: E1010+K0040+E0951. Left Part: FRMASMB9074. Right Part: FRMASMB9075	☐ ☐ ☐ 70° Articulating, Swing-Away.....\$1,197.50 ea. Legrest with 11"-13" Extension Pair Part: FRMASMB11444. HCPCS for pair: E1010+K0040+E0951. Left Part: FRMASMB11366. Right Part: FRMASMB11367.
---	---

70° Power Articulating, Swing-Away Legrest Compatibility Matrix	Synergy Seat			
	Low	Medium	High	
Angle Adjustable, 11" - 13"	✓	✓	✓	☒ This rigging configuration causes no interference issues at this seat-to-floor height. ☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
Angle Adjustable, 14.5" - 18"	!	!	!	

2. Angle Adjustable Footplate(s) with Heel Loop and Wheel Bumper

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 4"Wx6"D Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566.	☐ ☐ ☐ 5"Wx8"D Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.	☐ ☐ ☐ 6"Wx11"D ¹ Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570.
☐ ☐ ☐ 5"Wx6"D Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567.	☐ ☐ ☐ 6"Wx8"D ¹ Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569.	☐ ☐ ☐ 8"Wx11"D ² Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571.

3. Power Articulating, Swing-Away Legrest Accessories

Calf Pad (pair) 6"Wx6"H..... No Charge Part: MEC153735	Falcon Gel Padded Knee Adductor Buttons (KBE)(pair)\$125 Part: INDPART2242. HCPCS: E0953^ ^Please verify the appropriate billing code to be used for this item with all third-party payors.
---	---

21.3. Swing-Away Legrests

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest..... \$81 ea. Pair part: RIG120062. Left: FRMASMB17558. Right: FRMASMB17559 Not available with 11"-13" extensions.	☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short \$81 ea. Pair part: RIG136481. Left: RIG132607. Right: RIG132609. For use only with 11"-13" lower extensions.
---	---

1b. 70° Swing-Away Legrest Uppers

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

70° Swing-Away Legrest (pair).....\$162 Pair part: FRMASMB7841 Not available with 11"-13" extensions.	70° Tapered, Swing-Away Legrest, Short (pair).....\$162 Pair part: FRMASMB8052 For use only with 11"-13" lower extensions.
70° Tapered, Swing-Away Legrest (pair)\$162 Pair part: FRMASMB7843 Not available with 11"-13" extensions.	

Legrests section continued on following page.

21. LEGRESTS CONTINUED

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat			
	Low	Medium	High	
Extruded, 13" - 16.5"	!	✓	✓	☒ This rigging configuration causes no interference issues at this seat-to-floor height. ☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. ☒ This rigging configuration does not work with this seat-to-floor height.
Extruded, 15.5" - 19.5"	!	!	!	
Angle Adjustable, 13" - 16.5"	!	✓	✓	
Angle Adjustable, 15.5" - 19.5"	!	!	!	
Multi-Axis, 15.5" - 19.5"	!	!	!	

2. Lower Extension

☐ ☐ ☐ 11"-13" Length Part: See XRef N/A with Multi-axis.	☐ ☐ ☐ 11"-16.5" Length Part: See XRef N/A with Angle Adjustable.	☐ ☐ ☐ 13"-16.5" Length Part: See XRef N/A with Multi-Axis.	☐ ☐ ☐ 15.5"-19.5" Length Part: See XRef
---	---	---	---

3a. Extruded Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Extruded Footplate..... Standard Part pair: FRMASMB7847. Left: FRMASMB7861 Right: FRMASMB7862	☐ ☐ ☐ 8.5"Wx6"D Extruded Footplate Standard Part pair: FRMASMB7848. Left: FRMASMB7863 Right: FRMASMB7864 May cause interference on seat widths less than 18."
---	---

3b. Angle Adjustable Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 4"Wx6"D \$85 ea. Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560 HCPCS: K0040+E0951	☐ ☐ ☐ 6"Wx8"D ¹ \$85 ea. Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563 HCPCS: K0040+E0951
☐ ☐ ☐ 5"Wx6"D..... \$85 ea. Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561 HCPCS: K0040+E0951	☐ ☐ ☐ 6"Wx11"D ¹ \$85 ea. Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564 HCPCS: K0040+E0951
☐ ☐ ☐ 5"Wx8"D \$85 ea. Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562 HCPCS: K0040+E0951	☐ ☐ ☐ 8"Wx11"D ² \$85 ea. Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565 HCPCS: K0040+E0951

3c. Multi-Axis Footplate(s) with Heel Loops

¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 4"Wx6"D \$247.50 ea. Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552 HCPCS: K0040+E0951	☐ ☐ ☐ 6"Wx8"D ¹ \$247.50 ea. Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555 HCPCS: K0040+E0951
☐ ☐ ☐ 5"Wx6"D..... \$247.50 ea. Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553 HCPCS: K0040+E0951	☐ ☐ ☐ 6"Wx11"D ¹ \$247.50 ea. Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558 HCPCS: K0040+E0951
☐ ☐ ☐ 5"Wx8"D \$247.50 ea. Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554 HCPCS: K0040+E0951	☐ ☐ ☐ 8"Wx11"D ² \$247.50 ea. Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559 HCPCS: K0040+E0951

4. Swing-Away Legrest Accessories

Heel Loops (pair)..... \$50 Part: FRMASMB7873. HCPCS: E0951 Standard on angle adjustable and multi-axis footplates.	Wheel Bumper..... \$12.50 ea. Part: See XRef
	Gel Neoprene Footrest Sleeve (pair) \$80 Part: INDPART2243. HCPCS: K0108

5. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea. Part: INDPART2822. HCPCS: E1020 QTY	Therafin 32387 Amputee Support Gel Cover \$155 ea. Part: INDPART2823. HCPCS: K0108 QTY
---	--

21.4. Pediatric, Swing-Away Legrests

Not available with TRU-Balance Flex seating. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Pediatric Legrests Upper

☐ ☐ ☐ 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) \$81 ea. Pair part: FRMASMB8808. Left: FRMASMB8809 Right: FRMASMB8810

2. Angle Adjustable Footplate(s)

☐ ☐ ☐ 4"Wx6"D Angle Adjustable \$85 ea. Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560 HCPCS: K0040	☐ ☐ ☐ 5"Wx6"D Angle Adjustable \$85 ea. Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561 HCPCS: K0040
---	---

3. Accessories

Heel Loops (pair)..... \$50 Part: FRMASMB7873. HCPCS: E0951	☐ ☐ ☐ Wheel Bumper(s) \$12.50 ea. Pair part: WHL169508. Left: WHL169510. Right: WHL169511.
--	--

Legrests section continued on following page.

21. LEGRESTS CONTINUED**21.5. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$162.50 ea.
 Pair part: FRMASMB6483. Left: RIG138012
 Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

☐ ☐ ☐ 11"-14" Length ☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4755 Pair: FRMASMB4756
 Left: RIG138026 Left: RIG138033
 Right: RIG138030 Right: RIG138035
 For reinforced only For reinforced only

☐ ☐ ☐ 12"-14.5" Length ☐ ☐ ☐ 14.5"-18" Length
 Part: See XRef Part: See XRef
 For angle adjustable only For angle adjustable only

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D \$85 ea.
 Pair parts: FRMASMB8812 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$85 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx6"D..... \$85 ea.
 Pair parts: FRMASMB8811 x2.
 HCPCS: K0040

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$85 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$85 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops..... \$85 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Heel Loops (pair)..... \$50
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair)\$125
 Part: INDPART2242. HCPCS: E0953^
 ^Please verify the appropriate billing code to be used for this item with all third-party payors.

Wheel Bumper..... \$12.50 ea.
 Part: See XRef
 Available on angle adjustable footplates.

22. SET LEGREST LENGTH

REQUIRED with matching left and right legrest configuration. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

23. SEAT CUSHION


Optional. Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

23.1. Seat Cushion Model Selection**Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY

Simplicity G, 10"W to 21" W, E2601 \$110

Simplicity G, 22"W, E2602 \$165

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCFSPXXYY

Spectrum Foam SPP, 16"W to 20" W, E2607 \$460

Spectrum Foam SPP, 22"W, E2608 \$530


**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SOLSPXXYY

Solution SPP, 10"W to 21" W, E2607 \$375

Solution SPP, 22"W, E2608 \$430


**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-GLSPXXYY

Glacial SP, 10"W to 21"W, E2603 \$330

Glacial SP, 22"W, E2604 \$440


**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.

Part series: ST-ZENSPXXYY

Zen SP, 10"W to 21" W, E2603 \$315


Zen SP, 22"W, E2604 \$425

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY

TRU-Comfort 2 SPP, 16"W to 20" W, E2607 \$400

TRU-Comfort 2 SPP, 22"W, E2608 \$500

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.

Part series: ST-SPCGSPPXXYY

Spectrum Gel SPP, 14"W to 20" W, E2607 \$460

Spectrum Gel SPP, 22"W, E2608 \$530

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-ESSSPPXXYY

Essence SPP, No Wedge, 16"-20" Wide, E2607 \$400

Essence SPP, No Wedge, 22" Wide, E2608 \$500

Essence SPP, Small Wedge, 16"-20" Wide, E2607 \$435

Essence SPP, Small Wedge, 22" Wide, E2608 \$535

Essence SPP, Medium Wedge, 16"-20" Wide, E2607 . \$465

Essence SPP, Medium Wedge, 22" Wide, E2608 \$565

Essence SPP, Large Wedge, 16"-20" Wide, E2607 \$495

Essence SPP, Large Wedge, 22" Wide, E2608 \$595


**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.

Part series: ST-PREMPXXYY

Premiere P, 16"W to 20" W, E2605 \$325

Premiere P, 22" W, E2606 \$475

**23.2. Seat Cushion Dimensions**

REQUIRED with cushion model selection. See respective [cushion order form](#) for available sizes. Select to match or pick non-matching dimensions.

Match Seat Width and Depth

Cushion Width

10" W	14" W	16" W	18" W	20" W	22" W
12" W	15" W	17" W	19" W	21" W	

Cushion Depth

10" D	14" D	16" D	18" D	20" D	22" D
12" D	15" D	17" D	19" D	21" D	

23.3. Seat Cushion Alternate/Add Fluid Proof Liner

A liner may be added or the standard liner may be substituted with a selection below. *A substitution credit may apply. See the respective cushion order form for more details.

Fluid Proof Urethane Liner \$50-\$70*

Part: URL1709070-XXYY

Not available with Simplicity, Essence, and TRU-Comfort 2.

Fluid Proof Polyurethane Liner \$55-\$75*

Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY

Not available with Simplicity.

23.4. Seat Cushion Alternate Cover

An additional cover may be ordered or the standard cover may be substituted with a selection below. *A substitution credit may apply, see the respective cushion order form for more details. The "prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover No Charge

Part: *prefix-TBLKXXYY

Not available with Simplicity, TRU-Comfort 2, Essence, or Premiere.

Mesh Cover No Charge

Part: *prefix-MBLKXXYY

Not available with Simplicity or TRU-Comfort 2.

24. ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder \$20

Part: ACCASMB2287. Part with 12" lift: ACCASMB2288

Synergy/Flex Seat Cup Holder No Charge

Part: ACCASMB2634

Not available on TRU-Balance Flex seating without two-post armrests.

Oxygen Holder for Contoured Seat \$118

Part: ACCASMB2292. HCPCS: E2208

Oxygen Holder for Synergy Seat \$150

Not available with Synergy Manual Recline.

Part: FRMASMB7030. HCPCS: E2208

Oxygen Holder for Flex Seat \$200

Part: ACC125006. HCPCS: E2208

Not available with Synergy Manual Recline.

Walker Holder \$103

Contoured seat part: ACCASMB1912. Synergy/Flex part: FRMASMB7031

Cane/Crutch Holder \$110

Contoured seat part: ACCASMB2294. HCPCS: E2207

Synergy part: FRMASMB4822. Flex part: ACC152515. HCPCS: E2207

Cell Phone Holder for 12" Lift Captain Seat \$70

Part: ACC1611300

Accessories section continued on following page.

24. ACCESSORIES CONTINUED

Rear Basket\$118 Contoured seat part: ACCASMB1914 Synergy part: FRMASMB4821	Quantum Backpack \$150 Part: ACC1905853
Flag Holder\$25 Contoured seat part: ACCASMB1930 Synergy/Flex part: ACCASMB2415.	Grip Lap Board (20005)\$38.95 Part: ACC152992
Saddle Bag\$35 Part: DWR1020D010	Grip Drink Holder (20085)\$15.50 Part: ACC152993
Stealth Cell Phone/Media Device Holder\$195 Part: ACC7120006	Accessory Bar \$100 12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975. 20"W-22"W part: MEC146018
Stealth Large Tablet Holder \$375 Part: ACC7120005	Large, Black, Weather Cover \$202 Part: ACCCOVR1021

25. BRANCH PATCHES

Optional.

Air Force Patch \$20 Part: ACC140240	Coast Guard Patch \$20 Part: ACC140239	Navy Patch \$20 Part: ACC140244
Army Patch \$20 Part: ACC140241	Marines Patch \$20 Part: ACC140243	Purple Heart Patch \$20 Part: ACC140238

**26. PROGRAMMER**

Optional.

☐ Q-Logic EX/NE Programmer \$545 Part: ELEASMB6976	☐ Q-Logic PC Programming Disk with Cable \$750 Part: ELEASMB5216
---	---

27. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Only available with the Synergy Seating system. Select a type of specialty controls below or to refer to the specialty controls order form.

27.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit.....\$1,720
Part: CTL167726
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft TouchStandard Part: DGN137728	Match Joystick Shroud to Base Color \$50 Part: See color below
Non-matching Joystick Shroud \$50 Select color below	



Red Corvette Part: DGN137729-PRD1	White Satin Part: DGN137729-PWH2	Evolution Orange Part: DGN137729-POR7	Lemon Crush Part: DGN137729-PYL4	Military Digital Camo Part: DGN137729-WCM7
Cherry Bomb Part: DGN137729-PRD4	Deep Purple Part: DGN137729-PPR5	Midnight Blue Part: DGN137729-PBU10	Green Machine Part: DGN137729-PGN11	Faded Glory Part: DGN137729-WFL2
Back in Black Part: DGN137729-PBK1	Raspberry Beret Part: DGN137729-PPK4	Electric Blue Part: DGN137729-PBU9	Military Green Part: DGN137729-PGN8	Carbon Fiber Part: DGN137729-WCB1
Night Train Part: DGN137729-PBK3	Pretty in Pink Part: DGN137729-PPK3	Blue Jean Part: DGN137729-PBU11	Sugar Plum Part: DGN137729-PPR7	To Be Determined Follow-up required.
Grey Street Part: DGN137729-PSL3				

27.2. SCIM Module

Optional.

☐ SCIM Module for Q-Logic 3\$525 Part: CTL168055 Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.
--

**27.3. Q-Logic 3 SAJ**

☐ Q-Logic 3 EX SAJ\$700 Part: CTL143980 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.
--



Alternative controls and electronics section continued on following page.

27. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**27.4. Head Array**

Adult, i-Drive Tri-Array Series, 3-Switch \$4,414

Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array \$4,461

(IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

27.5. Sip and Puff Drive**1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 \$1,895

Part: CTL168056. E2325

Includes module and mounting hardware.

2. Sip and Puff Straw and MountStealth Sip and Puff for Comfort Series, \$294.25
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series, \$294.25
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, \$294.25
Head Array (SP625-U)

Part: ST-SP625-U-Q. E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Sip & Puff Module



Stealth Sip & Puff

27.6. Egg Switches

Red Egg Switch (ESRED) \$80.85

Part: SWTMCRO1077

Blue Egg Switch (ESBLU) \$80.85

Part: SWTMCRO1079

Black Egg Switch (ESBLK) \$80.85

Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$80.85

Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$80.85

Part: SWTMCRO1080



Egg Switch

27.7. Switch Adapter

2 to 1 Switch Adapter \$13.39

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

27.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,400

Part: CTL168057

28. POSITIONING COMPONENTSOptional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).**28.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, Rail Mounts, pair (TWBSTL) \$472.50

Part: INDPART2896. HCPCS: E1028

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5) \$210

Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) \$210

Part: INDPART2664. HCPCS: E0956

Flat Pads, 3"x5", pair (TWBLP3X5) \$210

Part: INDPART2676. HCPCS: E0956

Flat Pads, 4"x6", pair (TWBLP4X6) \$210

Part: INDPART2675. HCPCS: E0956

3. InstallationInstall Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**28.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, \$245 ea.

Removable Mounting Hardware

Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable, \$255 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$95 ea.

Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad \$125 ea.

Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad \$105 ea.

Part: ST-UTP548056100-Q



UniLink Pad & Hardware

Positioning components section continued on following page.

28. POSITIONING COMPONENTS CONTINUED**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link \$75 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link..... \$75 ea. Part: ST-UIL001200A00-Q
---	---

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

28.3. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$236.25 ea. Part: INDPART2804. HCPCS: E1028	Flip Down with LEFT Full Surface Contact \$236.25 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028
Removable Bracket (LP430)..... \$262.50 ea. Part: INDPART2805. HCPCS: E1028	Flip Down with RIGHT Full Surface Contact..... \$236.25 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028
Removable Bracket with 1" Offset \$400 ea. (LP430 & LHW-121) Part: INDPART2854. HCPCS: E1028	

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122).....\$26.25 Part: INDPART2821. HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129)\$53.50 Part: POS140430. HCPCS: K0108 Only available with the flip down full surface bracket.
--	---

3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea. Part: INDPART2801. HCPCS: E0953^/E0956	Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea. Part: INDPART2803. HCPCS: E0953^/E0956
Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea. Part: INDPART2802. HCPCS: E0953^/E0956	

4. Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

28.4. Stealth i-Fit Chest Harness

Upper Classic, Medium (IFPS1103) \$142 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$142 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$142 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$142 ea. Part: ST-IFPS2104-Q. HCPCS: E0960

28.5. Stealth i-Fit Chest Straps

Chest Strap, Medium w/ Belt Anchor (IFCS4113).. \$90 ea. Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large w/ Belt Anchor (IFCS4114)..... \$90 ea. Part: ST-IFCS4114-Q. HCPCS: E0960
---	---

28.6. Stealth Swing-Away Headrest Mounts


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral Hardware

Facial Lateral for Comfort Series.....\$294.25 ea. Right Side (CPS9-R-HW) Part: INDPART2647. HCPCS: E1028	Facial Laterals for Comfort Series \$347.75 ea. Left and Right Side (CPS9-B-HW) Part: INDPART2649. HCPCS: E1028
Facial Lateral for Comfort Series.....\$294.25 ea. Left Side (CPS9-L-HW) Part: INDPART2648. HCPCS: E1028	

2. Facial Lateral Pads and Accessories

1" Ball Kit for Facial Laterals (TWB-1BK) (D) \$80.25 Part: POS150532 Required with facial laterals to mount to comfort headrest brackets.	Facial Pad, 4"x2" (900-Q)..... \$53.50 ea. Part: ST-900-Q. HCPCS: K0108
Egg Switch Mount (ESM) (A) \$52.50 ea. Part: INDPART2815	Facial Pad, 6"x2" (920-Q)..... \$53.50 ea. Part: ST-920-Q. HCPCS: K0108
Figure 8 Bracket with 13" Rod..... \$267.50 ea. and Egg Switch Mount (B) Part: ST-F8675-13-Q	Swing-Away Joystick Mount,..... \$250 ea. 18" Rod (C) Part: ST-SUS18-Q Allows mounting of facial pads and attaches to the headrest



Positioning components section continued on following page.

28. POSITIONING COMPONENTS CONTINUED**28.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket\$183.75 ea.
 Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket\$183.75 ea.
 Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$52.50 ea.
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)\$52.50 ea.
 QTY Part: POS124067. HCPCS: K0108

28.8. Stealth Foot Boxes

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H\$498.75
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H\$498.75
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H\$498.75
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954^



Stealth Foot Boxes

28.9. Stealth Padded Shoe Holder

Stealth Small Padded Shoe Holder,\$266
 3.75"Wx6.25"D (SPLE111) (pair)
 Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder,\$266
 5"Wx10.75"D (SPLE114) (pair)
 Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder,\$266
 4"Wx8.25"D (SPLE112) (pair)
 Part: INDPART2784. HCPCS: E0951+E0952



Stealth Shoe Holders

29. Q-LOGIC 3 MULTIPLIER HARNESS

REQUIRED. Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA) No Charge
 Part: Configuration selection

Secondary Multiplier (SM) \$221.21
 Part: HAR144004

Special Order Configuration (SOF)
 Part: Configuration selection

Arm Mounts	Back Mounts			
		1 Arm	2 Arm	3+ Arm
Joystick (Hand Control) Enhanced Display iAccess USB Charger SAJ	Attendant Control SCIM Sip n Puff ECU	0 Back	NA	SM
		1 Back	NA	SM
		2 Back	SM	SM
		3+ Back	SOF	SOF

Q-Logic 3 Multiplier with Synergy Seat Chart

30. REFERENCE ORDER FORMS

Optional. Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)