



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Order Form with Synergy® Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Named/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3S-SS\$8,870

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge 3 3SP-SS\$8,559

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option for power adjustable seat height.

2. DRIVE MOTORS

REQUIRED.

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7 mph, 4 Pole Drive Motors..... \$1,653
Part: MOT1902455

Requires Q-Logic 3 electronics.

7 mph, 4 Pole, Accu-Trac Motor Technology..... \$1,360
Part: MOT2008501. HCPCS: K0108

Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette
Part: DGN1801101-PRD1

 Deep Purple
Part: DGN1801101-PPR5

 Midnight Blue
Part: DGN1801101-PBU10

 Green Machine
Part: DGN1801101-PGN11

 Back in Black
Part: DGN1801101-PBK1

 Raspberry Beret
Part: DGN1801101-PPK4

 Electric Blue
Part: DGN1801101-PBU9

 Military Green
Part: DGN1801101-PGN8

 Night Train
Part: DGN1801101-PBK3

 Pretty in Pink
Part: DGN1801101-PPK3

 Blue Jean
Part: DGN1801101-PBU11

 Grey Street
Part: DGN1801101-PSL3

 Evolution Orange
Part: DGN1801101-POR7

 Lemon Crush
Part: DGN1801101-PYL4

3.1b. Hydrographic Shrouds

☐  Military Digital Camo..... \$263
Part: DGN1801101-WCM7

☐  Faded Glory \$263
Part: DGN1801101-WFL2

☐  Carbon Fiber..... \$263
Part: DGN1801101-WCB1

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height is 15.875" with 6" stroke and 16.875" with 7" stroke. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

Power Elevating Seat with 7" Stroke \$1,990
Part: XRF000131. HCPCS: E2298/E2300^

Power Elevating Seat with 6" Stroke \$1,990
Part: XRF000132. HCPCS: E2298/E2300^

5. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Fender Lights..... No Charge
Part: DGN154117


Fender lights operate by pressing included single switch button. Fender lighting is not available with Manual Recline.

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129

Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

6. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$21

Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27

Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

7. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

7.1. Drive Wheels

3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair)

3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

4"x14" Knobby, Pneumatic Black Tires..... \$580
Part: XRF000367

Drive wheels section continued on following pages.

7. DRIVE WHEELS CONTINUED**7.2. Drive Wheel Accent Color**

Not available with knobby tires.

Back in Black.....Standard
Part: DGN164245-PBK1 x2Red Corvette.....No Charge
Part: DGN164245-PRD1 x2Night Train.....No Charge
Part: DGN164245-PBK3 x2Grey Street.....No Charge
Part: DGN164245-PSL3 x2Deep Purple.....No Charge
Part: DGN164245-PPR5 x2Raspberry Beret.....No Charge
Part: DGN164245-PPK4 x2Pretty in Pink.....No Charge
Part: DGN164245-PPK3 x2Evolution Orange.....No Charge
Part: DGN164245-POR7 x2Midnight Blue.....No Charge
Part: DGN164245-PBU10 x2Electric Blue.....No Charge
Part: DGN164245-PBU9 x2Blue Jean.....No Charge
Part: DGN164245-PBU11 x2Lemon Crush.....No Charge
Part: DGN164245-PYL4 x2Green Machine.....No Charge
Part: DGN164245-PGN11 x2Military Green.....No Charge
Part: DGN164245-PGN8 x2Military Digital Camo.....\$58 ea.
Part: DGN164245-WCM7 x2Faded Glory.....\$58 ea.
Part: DGN164245-WFL2 x2Carbon Fiber.....\$58 ea.
Part: DGN164245-WCB1 x2**8. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery.....\$355 ea.
Part: BAT2004376. HCPCS: E2361NF-22, Gel Cell, MK Battery.....\$432 ea.
Part: BATGEL1008. HCPCS: E2361NF-22, AGM, MK Battery.....\$432 ea.
Part: BAT121346. HCPCS: E2361**9. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 200 lb User Weight
Part: PTOINDV3837201 lb-250 lb User Weight
Part: PTOINDV3861251 lb-300 lb User Weight
Part: PTOINDV3984**10. SEAT TOWERS****Optional.** Not available with Power Adjustable Seat Height.Extended Seat Towers.....No Charge
Part: MEC120301
Increases seat to floor height by 3".**11. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick.....Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick.....\$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.**11.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module.....\$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick.....\$840
Part: CTL167724Q-Logic 3 EX Joystick for Full Lights.....\$840
Part: CTL167725Q-Logic 3e Joystick.....\$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display.

Electronics section continued on following pages.

11. ELECTRONICS CONTINUED**11.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635**11.3. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick.** Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN124307 Also matches Night Train shrouds.	Pretty in Pink \$63 Part: DGN124336-PPK3	Green Machine \$63 Part: DGN124336-PGN11
Red Corvette \$63 Part: DGN124336-PRD1	Evolution Orange \$63 Part: DGN124336-POR7	Military Green \$63 Part: DGN124336-PGN8
Back in Black \$63 Part: DGN124336-PBK1	Midnight Blue \$63 Part: DGN124336-PBU10	Military Digital Camo \$63 Part: DGN124336-WCM7
Grey Street \$63 Part: DGN124336-PSL3	Electric Blue \$63 Part: DGN124336-PBU9	Faded Glory \$63 Part: DGN124336-WFL2
Deep Purple \$63 Part: DGN124336-PPR5	Blue Jean \$63 Part: DGN124336-PBU11	Carbon Fiber \$63 Part: DGN124336-WCB1
Raspberry Beret \$63 Part: DGN124336-PPK4	Lemon Crush \$63 Part: DGN124336-PYL4	

12. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
Medium Ball Handle (E) \$69 Part: ACCASMB1859. HCPCS: E2323	Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
Soft Ball Handle (F) \$80 Part: ACCASMB1856. HCPCS: E2323	



**13. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**13.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**13.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**13.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**13.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368 Part: XRF000003. HCPCS: E1032 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Multi-Axis Mount \$126 Part: XRF000004 Allows height adjustment of joystick.
--	--

13.2. Joystick Mount LocationLeft Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226

14. SYNERGY SEATING SYSTEM

Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.

14.1. Synergy Seating Frame SelectionSynergy Seating System \$456
Part: PTOINDV2995**14.2. Synergy Seating Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.☐ WC-19 Occupied Transit Option for Synergy Seat \$289
Part: KITASMB1830To determine if your chair is WC-19 compliant, please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**14.3. Seat Width****REQUIRED.** Part number is SETWDTH plus number following # below.

10" W (#1000)	11" W (#1001)	12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

14.4. Seat Width Package**REQUIRED.** Maximum seat width is 22".

10"-14" Width Package Part: FRMASMB18011	14"-18" Width Package Part: FRMASMB18012	16"-20" Width Package Part: FRMASMB18013	20"-24" Width Package (Maximum width is 22") Part: FRMASMB18014
---	---	---	--

14.5. Seat Depth**REQUIRED.** Part number is SETDPTH plus number following # below.

10" D (#1000)	11" D (#1001)	12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

14.6. Seat Depth Package**REQUIRED.** Maximum seat depth is 22".


10"-14" Depth Package Part: FRMASMB6981	14"-18" Depth Package Part: FRMASMB17823	15"-20" Depth Package Part: FRMASMB7732	20"-24" Depth Package (Maximum depth is 22") Part: FRMASMB7733
--	---	--	---

14.7. Seat to Floor Height**REQUIRED.**Lowest Seat to Floor
Part: PTOINDV3558Middle Seat to Floor
Part: PTOINDV3559Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
Synergy Seat		
15.875"	16.875"	17.875"
Extended Towers		
Add 3 inches to above		
Seat to Floor Height		

14.8. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

Synergy seating section continued on following pages.

14. SYNERGY SEATING SYSTEM CONTINUED**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

No Top Foam
Reference part: XRF000382

FRG Top Foam
Reference part: XRF000381

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

No Liner
Reference part: XRF000385

Stealth Polyurethane Liner
Reference part: XRF000384

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Stealth Tek Cover
Reference part: XRF000388

Mesh Cover
Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY

Additional Mesh Cover
Reference Part: XRF000393
QTY

14.9. Synergy Seating Back Cane Selection**REQUIRED.**

Synergy Back Cane No Charge
Part: FRMASMB15845

TRU-Comfort Back Canes No Charge
Part: XRF000133

Flip Forward Back Cane \$184
Part: FRMASMB15846
Not available with oxygen holder, cane mounted armrests, or cantilever armrests.

2" Mesomorph (1" per side) Back Cane \$205
Part: FRMASMB12067

2" Endomorph (1" per side) Back Cane \$205
Part: FRMASMB12115

4" Endomorph (2" per side) Back Cane..... \$205
Part: FRMASMB12139

Manual Recline Back \$788
Part: FRMASMB12910. HCPCS: E1226
Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan.

14.10. Synergy Seating Back Cane Height**REQUIRED except if a TRU-Comfort back is desired.** 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes
Part: FRMASMB15832

17" High Canes
Part: FRMASMB15835

20" High Canes
Part: FRMASMB15838

23" High Canes
Part: FRMASMB15841

15" High Canes
Part: FRMASMB15833

18" High Canes
Part: FRMASMB15836

21" High Canes
Part: FRMASMB15839

24" High Canes
Part: FRMASMB15842

16" High Canes
Part: FRMASMB15834

19" High Canes
Part: FRMASMB15837

22" High Canes
Part: FRMASMB15840

25" High Canes
Part: FRMASMB15843

Synergy seating section continued on following pages.

14. SYNERGY SEATING SYSTEM CONTINUED**14.11. Synergy Seating Back Cushion****REQUIRED.** Select one type of back. See availability chart for sizes.**1a. Synergy Back Selection**

Synergy Back..... \$289
 Part: XRF000363

1b. TRU-Comfort 4 Back Selection

Back size is the product of seat width, back cane, and overall back height selections. See availability chart for sizes.

TRU-Comfort 4 Back, Back Cane Option, 16" to 20" Width \$744
 Part: TC4V10W21WBC. HCPCS: E2620

TRU-Comfort 4 Back, Back Cane Option, 22" Width..... \$855
 Part: TC4V22W32WBC. HCPCS: E2621

1c. ADI Back, AL Series, 4-Point, Active Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
 Part: ALBR16SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622
 Part: ALBR16TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622
 Part: ALBR16XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622
 Part: ALBR18SA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$622
 Part: ALBR18TA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$622
 Part: ALBR18XA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622
 Part: ALBR20SA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622
 Part: ALBR20TA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622
 Part: ALBR20XA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735
 Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735
 Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735
 Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735
 Part: ALBR18SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735
 Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
 Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour**Available with Back Cane Option only.** Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822
 Part: ALBR18SE4. HCPCS: E2615

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822
 Part: ALBR18TE4. HCPCS: E2615

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822
 Part: ALBR18XE4. HCPCS: E2615

1g. Omit Back Selection

Omit Back..... No Charge
 Part: PTOINDV1287

14.12. Synergy Seating Overall Back Height**REQUIRED.**

14" Overall Back Height
 Part: XRF000136

15" Overall Back Height
 Part: XRF000137

16" Overall Back Height
 Part: XRF000138

17" Overall Back Height
 Part: XRF000139

18" Overall Back Height
 Part: XRF000140

19" Overall Back Height
 Part: XRF000141

20" Overall Back Height
 Part: XRF000142

21" Overall Back Height
 Part: XRF000143

22" Overall Back Height
 Part: XRF000144

23" Overall Back Height
 Part: XRF000145

24" Overall Back Height
 Part: XRF000146

25" Overall Back Height
 Part: XRF000147

		Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	15"	S	S	-	-	-	-	-	-	-	-
	16"	S	S	-	-	-	-	-	-	-	TC
	17"	-	S	S	S	S	S	S	S	-	S TC
	18"	-	S	S	S	S	S	S	S	-	S TC
	19"	-	S	S	S	-	S	S	S	S	S TC
	20"	-	S	S	S	-	S TC	S	S	S	S TC
	21"	-	S	S	S	S	S TC	S	S	S	S TC
	22"	-	S	S	S TC	S TC	S TC	S	S TC	S	S TC
	23"	-	S	S	S TC	S TC	S	S	S TC	S	S TC
	24"	-	S	S	S	S	S TC	S	S TC	S	S
	25"	-	S	S	S	S	S TC	S	S TC	S	S

S Synergy TC TRU-Comfort 4

Synergy and TRU-Comfort Back Availability Chart

Synergy seating section continued on following pages.

14. SYNERGY SEATING SYSTEM CONTINUED**14.13. Synergy Armrests****REQUIRED.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

☐ ☐	Full Length, Flip Back Armrests..... Standard	☐ ☐	Desk Length, Flip Back Armrests Standard
L R	Left Part: FRMASMB7038. Right part: FRMASMB7039	L R	Left part: FRMASMB7040. Right part: FRMASMB7041

1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐	Full Length, Two-Post Armrests.....\$158 ea.	☐ ☐	Desk Length, Two-Post Armrests\$158 ea.
L R	Left part: XRF000051. Right part: XRF000052. HCPCS: E0973	L R	Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests

☐ ☐	Full Length, Heavy Duty Armrests, 8"-12" H \$90 ea.	☐ ☐	Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea.
L R	Left Part: XRF000245. Right part: XRF000246. HCPCS: E0973	L R	Left part: XRF000247. Right part: XRF000248. HCPCS: E0973

1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests

☐ ☐	Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.	☐ ☐	Desk Length, Heavy Duty Armrests, 10"-14"H \$90 ea.
L R	Left Part: XRF000249. Right part: XRF000250. HCPCS: E0973	L R	Left part: XRF000251. Right part: XRF000252. HCPCS: E0973

1e. Pediatric, Removable, Height Adjustable Armrests

Only available with desk length armpads.

☐ ☐	Pediatric, Removable Armrests \$79 ea.
L R	Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

1f. Height Adjustable, Cantilever Armrests

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐	Full Length, Cantilever Armrests\$103 ea.	☐ ☐	Desk Length, Cantilever Armrests.....\$103 ea.
L R	Left Part: FRMASMB16188. Right part: FRMASMB16189	L R	Left part: FRMASMB16190. Right part: FRMASMB16191

1g. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐	Full Length, Single Post Armrests \$79 ea.	☐ ☐	Desk Length, Single Post Armrests \$79 ea.
L R	Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973	L R	Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

14.14. Synergy Armpads**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard	Waterfall Armpads..... No Charge
Part: XRF000148	Part: XRF000149
No Armpad Selection No Charge	
Part: PTOINDV1373	

14.15. Joystick Bracket Receiver**REQUIRED.**

☐ Standard Fixed Receiver Standard	☐ Stealth Height Adjustable Receiver..... \$236
Part: FRMASMB16062	Part: XRF000064. HCPCS: K0108

15. HEADRESTS**Optional.** Select a headrest pad and bracket.**15.1. Headrests, Comfort Plus with Cool Core****1. Headrest Pads**

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179	Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955	Part: ST-CP14CCFPM-Q. HCPCS: E0955
Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179	Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955	Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311	Multi-Axis Removable Headrest Bracket\$174
Part: ST-UHB2D33200-Q. HCPCS: E1033	Part: ST-HMO475-17M-Q. HCPCS: E1033

15.2. Flip-down Headrest Bracket

Stealth Flip-down Headrest Bracket (FDM350)..... \$170
Part: INDPART2219

The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

15.3. Q Cam Lock**Optional.**

Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

16. LAP BELTS**REQUIRED.** Select one lap belt below.**16.1a. Standard Lap Belts**

- ☐ **38" Lap Belt** Part: ACCBELT1006
 ☐ **50" Lap Belt** Part: ACCBELT1000
 ☐ **60" Lap Belt** Part: ACCBELT1002
 ☐ **70" Lap Belt** Part: ACCBELT1003
 ☐ **80" Lap Belt** Part: ACCBELT1004
 ☐ **90" Lap Belt** Part: ACCBELT1001

16.1b. Optional Lap Belts

- ☐ **Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50**
 (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q
- ☐ **Stealth Positioning 60" Single Pull, Push Button, \$50**
 Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5315-Q
- ☐ **Stealth Positioning 60" Single Pull, Push Button, Large \$50**
 Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q
- ☐ **60" Sure-Lok Lap Belt (Non-retractable)..... \$120**
 Required with WC-19 Transit Option.
 Part: ACCBELT1016

17. POWER LEGREST CONTROLS**Optional.** Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

- ☐ **Independent Legrests (adds 2 actuator control requirement)**
 Part: XRF000108
- ☐ **Articulating Foot Platform (adds 1 actuator control requirement)**
 Part: XRF000110
- ☐ **Combined Legrests (adds 2 actuator control requirement)**
 Part: XRF000109

18. LEGRESTS**Optional.** Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.**18.1a. Base Mounted Foot Platform**

May interfere with seat mounted legrests. Not available with center mount foot platforms.

- Base Mounted Foot Platform Standard** **Omit base mounted foot platform No Charge**
 Part: FRMASMB16077 Part: PTOINDV3674
- Bell Shaped, Base Mounted Foot Platform \$247**
 Part: RIGASMB7120041

18.1b. Basic Center Mount Foot Platform**1. Basic Center Mount Foot Platform Assembly**

- Basic Center Mount..... No Charge**
 Part: XRF000323

2. Basic Center Mount Extension and Angle Adjustable Footplates

- Short Extension, 5"x9" Footplates \$180** **Long Extension, 5"x9" Footplates \$180**
 Part: XRF000312. HCPCS: K0040 Part: XRF000314. HCPCS: K0040
- Medium Extension, 5"x9" Footplates \$180**
 Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

- Small Calf Pads, Basic Center Mount..... \$158** **Large Calf Pads, Basic Center Mount \$158**
 Part: RIG2104369 Part: RIG2202997
 Only available with Long Extension.

4. Basic Center Mount Additional Accessory

- Basic Center Mount Footplate Connector \$16**
 Part: RIG2104240

18.1c. Center Mount Foot Platform

Select legrest length in following section.

1. Foot Platform Assembly

- Power Articulating Foot Platform \$3,145**
 Part: PTOINDV3754. HCPCS: E1012+K0040

2. Lower Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

- | | | |
|---|--|--|
| Medium Extension, 5"x6" Footplates
Part: RIG150206 | Medium Extension, 5"x11" Footplates
Part: RIG2404018 | Long Extension, 5"x8" Footplates
Part: RIG2404017 |
| Medium Extension, 5"x8" Footplates
Part: RIG2404016 | Long Extension, 5"x6" Footplates
Part: RIG151505 | Long Extension, 5"x11" Footplates
Part: RIG2404019 |

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

- Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard**
 Part: RIGASMB7120018
- Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard**
 Part: RIGASMB7120019
- Small Calf Pads, Manually Adjustable Center Mount Only, \$158**
 6"Wx6"H (pair)
 Part: XRF000156
- Large Calf Pads, Manually Adjustable Center Mount Only, \$158**
 8" taper to 6"Wx8"H (pair)
 Part: XRF000157

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**3b. Center Mount Calf Panel**

Calf pads are standard with a power articulating foot platform.

Small Calf Panel, 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge \$69
Part: RIG2404921
Available with 5"x8" footplates.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923
Available with 5"x8" footplates.Width Extension for 11" \$69
Part: RIG2404920**5. Center Mount Additional Accessories**Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-QLarge Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-QLarge Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-QFootplate Connector \$16
Part: RIGASMB7120023Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108**18.1d. Swing-Away Legrests**

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Uppers

Adjustable in 5° increments.

P	L	R	60°-90° Swing-Away Legrest, Static Seat \$86 ea. Pair part: XRF000082. Left: XRF000083. Right: XRF000084 Not available with 11"-13" extensions.
---	---	---	---

P	L	R	60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea. Pair part: XRF000085. Left: XRF000086. Right: XRF000087 For use only with 11"-13" lower extensions.
---	---	---	---

Swing-Away Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Extruded, 13" - 16.5"	!	✓	✓
Extruded, 15.5" - 19.5"	!	!	!
Angle Adjustable, 13" - 16.5"	!	✓	✓
Angle Adjustable, 15.5" - 19.5"	!	!	!
Multi-Axis, 15.5" - 19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.
☒ This rigging configuration does not work with this seat-to-floor height.

2. Lower Extension

P	L	R	11"-13" Length Pair part: XRF000088. Left: XRF000089. Right: XRF000090.
---	---	---	--

P	L	R	13"-16.5" Length Pair part: XRF000094. Left: XRF000095. Right: XRF000096.
---	---	---	--

P	L	R	15.5"-19.5" Length Pair part: XRF000097. Left: XRF000098. Right: XRF000099.
---	---	---	--

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

P	L	R	5"Wx8"D Footplates with Heel Loops \$90 ea. Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792 HCPCS: K0040+E0951
---	---	---	--

P	L	R	6"Wx11"D1 Footplates with Heel Loops \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797 HCPCS: K0040+E0951
---	---	---	--

P	L	R	6"Wx8"D1 Footplates with Heel Loops \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795 HCPCS: K0040+E0951
---	---	---	---

P	L	R	8"Wx11"D2 Footplates with Heel Loops \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799 HCPCS: K0040+E0951
---	---	---	--

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
Part: INDPART2822. HCPCS: E1020Therafin 32387 Amputee Support Gel Cover \$155 ea.
Part: INDPART2823. HCPCS: K0108

Legrests section continued on following pages.

18. LEGRESTS CONTINUED**18.1e. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Contoured/Synergy Seat		
	Low	Medium	High
Standard, 11"-14"	✓	✓	✓
Standard, 14"-17"	!	!	!
Angle Adjustable, 12.5"-16.5"	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!

☒ This rigging configuration causes no interference issues at this seat-to-floor height.
☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

18.2. Legrest Accessories

Heel Loops (pair)..... \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$19 ea.
 Part: XRF000100

18.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

19. SYNERGY POWER OPTION ELECTRONICS**REQUIRED for Synergy Seat with Power Legrests or Elevating Seat.** Select the desired control type for the seating system.Single Actuator Through Controller Interface \$2,163
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and Interface.....\$3,040
Part: XRF000336. HCPCS: E2311Single Actuator Through Controller and Interface (AAM).....\$3,040
Part: XRF000068. HCPCS: E2310**20. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

21. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.Synergy Seat Cup Holder No Charge
Part: ACCASMB2634Grip Lap Board (20005)..... \$38.95
Part: ACC152992Oxygen Holder for Synergy Seat..... \$158
Not available with Synergy Manual Recline.
Part: FRMASMB7030. HCPCS: E2208Grip Drink Holder (20085) \$15.50
Part: ACC152993Walker Holder \$109
Part: FRMASMB7031Accessory Bar \$132
Part: XRF000128Cane/Crutch Holder \$116
Part: FRMASMB4822. HCPCS: E2207Large, Black, Weather Cover \$213
Part: ACCCOVR1021Rear Basket \$124
Part: FRMASMB4821Privacy Bag..... \$42
Part: ACC2005723Flag Holder \$55
Part: ACCASMB2415Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Saddle Bag \$40
Part: DWR1020D010Q-Logic PC Programming Disk with Cable \$445
Part: ELEASMB5216Stealth Cell Phone/Media Device Holder \$215
Part: ACC7120006Personal Item Pouch \$84
Part: ACC2101017Stealth Large Tablet Holder \$625.28
Part: ST-URPKTAB0000-QTransfer Handle \$105 ea.
Part: MEC2205680 x2. HCPCS: E0705
Available with Flip-Up armrests only.Quantum Backpack \$168
Part: ACC1905853Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496**22. BRANCH PATCHES****Optional.**Air Force Patch \$21
Part: ACC140240Coast Guard Patch \$21
Part: ACC140239Navy Patch \$21
Part: ACC140244Army Patch \$21
Part: ACC140241Marines Patch \$21
Part: ACC140243Purple Heart Patch \$21
Part: ACC140238**23. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**23.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with 2-post, flip up armrests. These laterals are not available with the back cane option.

1. BracketSwing-away, Rail Mounts, pair (TWBSTL) \$562
Part: INDPART2896. HCPCS: E1034Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
Part: ST-ULPK404108B00-Q. HCPCS: E1034
Only available with TRU Comfort backs.

Only available with back cane and omit back configurations.

2. PadsCurved Pads, 3"x5", pair (TWBL-3X5) \$250
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250
Part: INDPART2664. HCPCS: E0956

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**4. Lateral Extensions****REQUIRED with Synergy Back and thoracic lateral selections.**Lateral Extension for Synergy Back (pair)..... \$56
Part: XRF000150**23.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^Flip Down with LEFT Full Surface Contact \$281
Feature (TWBADD-L & LPHW407)
Part: POS140428. HCPCS: E1028/E1034^Flip Down with RIGHT Full Surface Contact..... \$281
Feature (TWBADD-R & LPHW407)
Part: POS140429. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^**2. UniLink Pelvic/Thigh Pads**


^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 12" Pad \$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^UniLink, 4" x 8" Pad..... \$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^**3. UniLink Pelvic/Thigh Additional Links**UniLink, 1.5" Link \$81 ea.
Part: ST-UIL001150A000-QUniLink, 2" Link..... \$81 ea.
Part: ST-UIL001200A000-Q**4. UniLink Pelvic/Thigh Guide Installation**Install on Left
Part: PTOINDV3585Install on Right
Part: PTOINDV3586Do Not Install
Part: PTOINDV3719**23.3. Stealth Positioning Chest Harness**Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960**23.4. Stealth Positioning Chest Straps**Chest Strap, Medium w/ Belt Anchor (IFCS4113) \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor (IFCS4114) \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960**23.5. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. Also select a facial pad(s) with the facial lateral(s).

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right..... \$318 ea.
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower LeftUniLink Facial Lateral, Swing Away, Left \$318 ea.
Part: ST-SAR2A-L-Q. HCPCS: E1028

Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028

UniLink Facial Lateral, Dual Swing Away \$492 ea.
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$434 ea.
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028**2. Facial Lateral Pads and Accessories**Egg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815Swing-Away Joystick Mount,..... \$276 ea.
18" Rod (C)Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrestFacial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea.
Part: ST-920M-Q. HCPCS: K0108

Positioning components section continued on following page.

23. POSITIONING COMPONENTS CONTINUED**23.6. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket\$199 ea.
 Tube Style Arms (ARMS-ES-L)
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket\$199 ea.
 Tube Style Arms (ARMS-ES-R)
 Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
 QTY Part: POS124067. HCPCS: K0108

23.7. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
 (SPLE107) (pair)
 Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
 (SPLE109) (pair)
 Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
 (SPLE108) (pair)
 Part: INDPART2787. HCPCS: E0954

23.8. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211) \$87 ea.
 QTY Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
 QTY Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
 QTY Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

24. LAP TRAY**24.1. Simple Solutions Lap Tray**

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	M	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	M	M	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	M	M	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	M	M	L	L	
Cantilever	S	S	S	S	S	M	M	M	M	L	L	
Q-Captain, Flip-Up						S	M		M	L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	M	L	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

S:Compatible

M:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

25. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below or to refer to the specialty controls order form.**25.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN137728**Deep Purple** \$63
Part: DGN137729-PPR5**Electric Blue** \$63
Part: DGN137729-PBU9**Military Digital Camo** \$63
Part: DGN137729-WCM7**Red Corvette**..... \$63
Part: DGN137729-PRD1**Raspberry Beret**..... \$63
Part: DGN137729-PPK4**Blue Jean** \$63
Part: DGN137729-PBU11**Faded Glory** \$63
Part: DGN137729-WFL2**Back in Black**..... \$63
Part: DGN137729-PBK1**Pretty in Pink**..... \$63
Part: DGN137729-PPK3**Lemon Crush**..... \$63
Part: DGN137729-PYL4**Carbon Fiber** \$63
Part: DGN137729-WCB1**Night Train**..... \$63
Part: DGN137729-PBK3**Evolution Orange**..... \$63
Part: DGN137729-POR7**Green Machine** \$63
Part: DGN137729-PGN11**Grey Street**..... \$63
Part: DGN137729-PSL3**Midnight Blue** \$63
Part: DGN137729-PBU10**Military Green** \$63
Part: DGN137729-PGN8

Enhanced Display

25.2. Q-Logic 3 SAJ☐ **Q-Logic 3 EX SAJ** \$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

25.3. Head Array**i-Drive 4.0 Bluetooth Tri-array Series**,..... \$4,774
3 Switch with Egg Switch
Part: ST-IDHBT200-U-Q. HCPCS: E2330**Adult, Curved, i-Drive Tri-Array with** \$4,774
Mounting Hardware
Part: ST-IDHBT400-U-Q. HCPCS: E2330**i-Drive 4.0 Bluetooth Pro Series**,..... \$4,825
3 Switch with Egg Switch
Part: ST-IDHBT300-U-Q. HCPCS: E2330**Junior, Curved, i-Drive Tri-Array with** \$4,774
Mounting Hardware
Part: ST-IDHBT450-U-Q. HCPCS: E2330

Adult Pro-Series

25.4. Egg Switches**Red Egg Switch (ESRED)**..... \$87
Part: SWTMCRO1077**Green Egg Switch (ESGRN)** \$87
Part: SWTMCRO1078**Blue Egg Switch (ESBLU)** \$87
Part: SWTMCRO1079**Yellow Egg Switch (ESYEL)** \$87
Part: SWTMCRO1080**Black Egg Switch (ESBLK)** \$87
Part: SWTMCRO1081

Egg Switch

25.5. Switch Adapter**2 to 1 Switch Adapter**..... \$15

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

25.6. Environmental Controls**Q-Logic 3 ECU Module Kit** \$2,520
Part: CTL168057