



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU Balance® 3 Positioning

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked for: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____
A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3S-SS\$8,870

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

Quantum Q6 Edge 3 3SP-SS\$8,559

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option.

Quantum Q6 Edge 3 3MP-SS\$8,373

Part: Q6 EDGE 3 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. DRIVE MOTORS**REQUIRED.**

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

7 mph, 4 Pole Drive Motors \$1,653
Part: MOT1902455

Requires Q-Logic 3 electronics.

7 mph, 4 Pole, Accu-Trac Motor Technology \$1,360
Part: MOT2008501. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 **Red Corvette**
Part: DGN1801101-PRD1

 **Deep Purple**
Part: DGN1801101-PPR5

 **Midnight Blue**
Part: DGN1801101-PBU10

 **Green Machine**
Part: DGN1801101-PGN11

 **Back in Black**
Part: DGN1801101-PBK1

 **Raspberry Beret**
Part: DGN1801101-PPK4

 **Electric Blue**
Part: DGN1801101-PBU9

 **Military Green**
Part: DGN1801101-PGN8

 **Night Train**
Part: DGN1801101-PBK3

 **Pretty in Pink**
Part: DGN1801101-PPK3

 **Blue Jean**
Part: DGN1801101-PBU11

 **Grey Street**
Part: DGN1801101-PSL3

 **Evolution Orange**
Part: DGN1801101-POR7

 **Lemon Crush**
Part: DGN1801101-PYL4

3.1b. Hydrographic Shrouds

☐  **Military Digital Camo** \$263
Part: DGN1801101-WCM7

☐  **Faded Glory** \$263
Part: DGN1801101-WFL2

☐  **Carbon Fiber** \$263
Part: DGN1801101-WCB1

4. BASE MOUNTED LIGHTS

Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Accent Lights No Charge
Part: DGN154117


Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129

Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.



Base Mounted Lights



Accent Lights

5. QUANTUM BACKUP CAMERA

Optional. Requires Q-Logic electronics. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21
Part: MEC1910373

Works with .5" x .5" vertical supports only.

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard
Part: PTOINDV3525 (pair)

3"x14" Pneumatic Black Tires No Charge
Part: WHL1713925 x2

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367

Drive wheels section continued on following page.

6. DRIVE WHEELS CONTINUED**6.2. Drive Wheel Accent Color****Optional.** Not available with knobby tires.Back in Black Standard
Part: DGN164245-PBK1 x2Red Corvette No Charge
Part: DGN164245-PRD1 x2Night Train No Charge
Part: DGN164245-PBK3 x2Grey Street No Charge
Part: DGN164245-PSL3 x2Deep Purple No Charge
Part: DGN164245-PPR5 x2Raspberry Beret No Charge
Part: DGN164245-PPK4 x2Pretty in Pink No Charge
Part: DGN164245-PPK3 x2Evolution Orange No Charge
Part: DGN164245-POR7 x2Midnight Blue No Charge
Part: DGN164245-PBU10 x2Electric Blue No Charge
Part: DGN164245-PBU9 x2Blue Jean No Charge
Part: DGN164245-PBU11 x2Lemon Crush No Charge
Part: DGN164245-PYL4 x2Green Machine No Charge
Part: DGN164245-PGN11 x2Military Green No Charge
Part: DGN164245-PGN8 x2Military Digital Camo \$58 ea.
Part: DGN164245-WCM7 x2Faded Glory \$58 ea.
Part: DGN164245-WFL2 x2Carbon Fiber \$58 ea.
Part: DGN164245-WCB1 x2**7. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376. HCPCS: E2361NF-22, AGM, MK Battery \$432 ea.
Part: BAT121346. HCPCS: E2361NF-22, Gel Cell, MK Battery \$432 ea.
Part: BATGEL1008. HCPCS: E2361**8. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 200 lb User Weight
Part: PTOINDV3837201 lb-250 lb User Weight
Part: PTOINDV3861251 lb-300 lb User Weight
Part: PTOINDV3984**9. EXTENDED SEAT TOWERS****Optional.** Only available with TRU Balance Static Seat or Power Recline Only. Not available with base mounted foot platform.Extended Seat Towers No Charge
Part: MEC2107219**10. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**10.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key Joystick Standard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick \$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.**10.1b. Expandable Controls**

Select expandable controls or omit selection if another controller type is desired.

1. Expandable Controls Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected elsewhere, Customer Service will quote the order with these items.

Q-Logic 3 Controller Module \$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness \$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$840
Part: CTL167724Q-Logic 3 EX Joystick for Full Lights \$840
Part: CTL167725Q-Logic 3e Joystick \$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.No Joystick, Setup for Q-Logic
Part: XRF000010
Requires selection of Enhanced Display.**10.2. Q-Logic Magnetic Charger**

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635


Electronics section continued on following pages.

10. ELECTRONICS CONTINUED**10.3. Q-Logic 3 EX Joystick Armor****REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.**

Black, Soft TouchStandard Part: DGN124307 Also matches Night Train shrouds.	Pretty in Pink \$63 Part: DGN124336-PPK3	Green Machine \$63 Part: DGN124336-PGN11
Red Corvette \$63 Part: DGN124336-PRD1	Evolution Orange \$63 Part: DGN124336-POR7	Military Green \$63 Part: DGN124336-PGN8
Back in Black \$63 Part: DGN124336-PBK1	Midnight Blue \$63 Part: DGN124336-PBU10	Military Digital Camo \$63 Part: DGN124336-WCM7
Grey Street \$63 Part: DGN124336-PSL3	Electric Blue \$63 Part: DGN124336-PBU9	Faded Glory \$63 Part: DGN124336-WFL2
Deep Purple \$63 Part: DGN124336-PPR5	Blue Jean \$63 Part: DGN124336-PBU11	Carbon Fiber \$63 Part: DGN124336-WCB1
Raspberry Beret \$63 Part: DGN124336-PPK4	Lemon Crush \$63 Part: DGN124336-PYL4	

11. JOYSTICK HANDLES**Optional.** The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

☐ Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	☐ Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
☐ Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	☐ Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
☐ T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	☐ Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
☐ Medium Ball Handle (E) \$69 Part: ACCASMB1859. HCPCS: E2323	☐ Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
☐ Soft Ball Handle (F) \$80 Part: ACCASMB1856. HCPCS: E2323	

**12. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.**12.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic Standard
Part: FRMASMB10617**12.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.



Swing-away, Q-Logic \$319
Part: ACC1908434. HCPCS: E1032**12.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**12.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount \$368 Part: XRF000003. HCPCS: E1032 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Multi-Axis Mount \$126 Part: XRF000004 Allows height adjustment of joystick.
--	--

12.2. Joystick Mount Location

☐ Left Side Mount Part: PTOINDV1225	☐ Right Side Mount Part: PTOINDV1226
---	--



13. TRU BALANCE® 3 SEATING**REQUIRED.** Maximum seat size is 22"W x 22"D.**13.1. TRU Balance 3 Configuration**

TRU Balance 3 Static Seat.....\$2,410
 Part: PTOINDV3442
 20° of limited recline adjustment (-5/ +15) via gas strut within back cover.

TRU Balance 3 Power Tilt Only\$8,202
 Part: PTOINDV3443. HCPCS: E1002
 50° Center of Gravity (CG) power tilt seating system.

TRU Balance 3 Power Recline Only\$8,234
 Part: PTOINDV3444. HCPCS: E1004
 175° of power recline.

TRU Balance 3 Power Tilt & Recline\$12,615
 Part: PTOINDV3446. HCPCS: E1007
 50° CG power tilt with 175° of power recline. Full drive lockout at 20°.

13.2. 12" Power Lift**Optional.** ^Please verify the appropriate billing code to be used for this item with all third-party payors.

TRU Balance 3, 12" Power Lift with iLevel® \$3,725
 Part: PTOINDV3773. HCPCS: E2298/E2300^
 12" of scissor lift with iLevel System.

**13.3. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

WC-19 Occupied Transit Option..... \$289
 Part: KIT152047

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

13.4. Seat Width**REQUIRED with seat selection.** Maximum seat width is 22". Part number is SETWDTH plus number following # below. 12" Width is not available with 16" and 17" Depth. 13" and 15" widths are only available with non-matching back width.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.5. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package
 Part: XRF000006

14"-16" Width Package
 Part: XRF000007

16"-20" Width Package
 Part: XRF000008

20"-24" Width Package
 Part: XRF000158
 Maximum width is 22".

13.6. Seat Depth**REQUIRED with seat selection.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

13.7. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor
 Part: PTOINDV3558


Middle Seat to Floor
 Part: PTOINDV3559

Highest Seat to Floor
 Part: PTOINDV3560

Low	Med	High
Static Seat, Recline		
16.5"	17.5"	18.5"
Tilt, Tilt & Recline, Lift		
17.5"	18.5"	19.5"
TRU Balance 3 Seat to Floor Heights		

13.8. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.



TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1. Cushion Type**

TRU-Comfort 2 SPP
Reference part: XRF000116. [See cushion order form.](#)
Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPP
Reference part: XRF000118. [See cushion order form.](#)
Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPP
Reference part: XRF000113. [See cushion order form.](#)
Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPP
Reference part: XRF000112. [See cushion order form.](#)
Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPP
Reference part: XRF000117. [See cushion order form.](#)
Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere P
Reference part: XRF000122. [See cushion order form.](#)
Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SP
Reference part: XRF000114. [See cushion order form.](#)
Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SP
Reference part: XRF000115. [See cushion order form.](#)
Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity G
Reference part: XRF000111. [See cushion order form.](#)
Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

2. Essence Cushion Wedge Kit

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

3. Cushion Top Foam

Visco Top Foam
Reference part: XRF000380

FRG Top Foam
Reference part: XRF000381

No Top Foam
Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
Reference part: XRF000383

Stealth Polyurethane Liner
Reference part: XRF000384

No Liner
Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
Reference part: XRF000386

Mesh Cover
Reference part: XRF000387

Stealth Tek Cover
Reference part: XRF000388

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
Reference part: XRF000389
MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
Reference Part: XRF000390
QTY _____

Additional Stealth Polyurethane Liner
Reference Part: XRF000391
QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
Reference Part: XRF000392
QTY _____

Additional Mesh Cover
Reference Part: XRF000393
QTY _____

Additional Stealth Tek Cover
Reference Part: XRF000394
QTY _____

13.9. Back Type Selection**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
Part: XRF000014
12" and 24" back width will use the TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option
Part: XRF000015
Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184
Part: XRF000016
12" and 24" back width will use the TRU Balance 3 Back assembly. Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

	None	Visco/V Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.10. Back Cane Selection (Back Cane Option Only)****REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**Angle Adjustable Back Cane No Charge
Part: MEC140266Flip Forward Back Cane \$53
Part: XRF0002732" Endomorph (1" per side) Back Cane \$74
Part: XRF0002694" Endomorph (2" per side) Back Cane \$74
Part: XRF000271**2. Back Cane Height**16"H Cane
Part: FRMASMB1583420"H Cane
Part: FRMASMB1583824"H Cane
Part: FRMASMB1584217"H Cane
Part: FRMASMB1583521"H Cane
Part: FRMASMB1583925"H Cane
Part: FRMASMB15843**13.11. Back Cushion Configuration****REQUIRED.** Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1a. Sport Back/Omit/After-Market Interface**Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface \$535
Part: ST-SSAPTB4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.**1b. TRU-Comfort Back**TRU-Comfort Back (Static/Flip Forward), 12" Width \$744
Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$744
Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$744
Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855
Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$855
Part: TC4V22W32WBC. HCPCS: E2621

Complete following width package and size.

	Tilt Only Width												Recline Width											
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		12"	14"
Overall Height	14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"
	15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"
	16"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	16"
	17"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	17"
	18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	18"
	19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	19"
	20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	20"
	21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	21"
	22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	22"
	23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	23"
	24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	24"
	25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	25"
	26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	26"
	27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	27"
	28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	28"
	29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	29"
	30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	30"
	31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	31"
	32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	32"

S Sport Back
Q Q-Back
T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1c. ADI Back, AL Series, 4-Point, Active Contour**

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 16" Wide x 19" or 20" Overall Height \$622
Part: ALBR16SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622
Part: ALBR16TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622
Part: ALBR16XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622
Part: ALBR18SA4. HCPCS: E2613

Active Contour, 18" Wide x 21" or 22" Overall Height \$622
Part: ALBR18TA4. HCPCS: E2613

Active Contour, 18" Wide x 23" or 24" Overall Height \$622
Part: ALBR18XA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622
Part: ALBR20SA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622
Part: ALBR20TA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622
Part: ALBR20XA4. HCPCS: E2613

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 16" Wide x 19" or 20" Overall Height \$735
Part: ALBR16SD4. HCPCS: E2615

Deep Contour, 16" Wide x 21" or 22" Overall Height \$735
Part: ALBR16TD4. HCPCS: E2615

Deep Contour, 16" Wide x 23" or 24" Overall Height \$735
Part: ALBR16XD4. HCPCS: E2615

Deep Contour, 18" Wide x 19" or 20" Overall Height \$735
Part: ALBR18SD4. HCPCS: E2615

Deep Contour, 18" Wide x 21" or 22" Overall Height \$735
Part: ALBR18TD4. HCPCS: E2615

Deep Contour, 20" Wide x 19" or 20" Overall Height \$735
Part: ALBR20SD4. HCPCS: E2615

Deep Contour, 20" Wide x 23" or 24" Overall Height \$735
Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 18" Wide x 19" or 20" Overall Height \$822
Part: ALBR18S64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 21" or 22" Overall Height \$822
Part: ALBR18T64. HCPCS: E2615

Extra Deep Contour, 18" Wide x 23" or 24" Overall Height \$822
Part: ALBR18X64. HCPCS: E2615

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20AXX. HCPCS: E2613

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**1g. Q-Back 2.0, AL Series, Deep Contour****Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

1h. Q-Back 2.0, AL Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL14ES. HCPCS: E2620	Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18EXX. HCPCS: E2620
Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703 Part: QB2AL16ES. HCPCS: E2620	Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL20ES. HCPCS: E2620
Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703 Part: QB2AL16ET. HCPCS: E2620	Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL20ET. HCPCS: E2620
Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL16EX. HCPCS: E2620	Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL20EX. HCPCS: E2620
Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703 Part: QB2AL18ES. HCPCS: E2620	Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20EXX. HCPCS: E2620
Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703 Part: QB2AL18ET. HCPCS: E2620	
Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703 Part: QB2AL18EX. HCPCS: E2620	

1i. Q-Back 2.0, AFT Series, Active Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT16AS. HCPCS: E2613	AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AFT18AS. HCPCS: E2613
--	--

1j. Q-Back 2.0, AFT Series, Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643 Part: QB2AFT16DS. HCPCS: E2620	AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643 Part: QB2AFT18DT. HCPCS: E2620
AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT16DT. HCPCS: E2620	AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AFT20DT. HCPCS: E2620
AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AFT18DS. HCPCS: E2620	

1k. Q-Back 2.0, AFT Series, Extra Deep Contour**Not available with Back Cane Option.** Select a back type and the specific overall back height in subsequent section.

AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline.....\$703 Part: QB2AFT16EXX. HCPCS: E2620	
---	--

13.12. TRU Balance 3 Back Shroud**REQUIRED.** The back shroud is a feature of the TRU Balance 3 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette Part: XRF000019	Grey Street Part: XRF000023	Pretty in Pink Part: XRF000027	Electric Blue Part: XRF000030	Green Machine Part: XRF000033
Back in Black Part: XRF000021	Deep Purple Part: XRF000025	Evolution Orange Part: XRF000028	Blue Jean Part: XRF000031	Military Green Part: XRF000034
Night Train Part: XRF000022	Raspberry Beret Part: XRF000026	Midnight Blue Part: XRF000029	Lemon Crush Part: XRF000032	

TRU Balance 3 Seating section continued on following pages.

13. TRU BALANCE 3 SEATING CONTINUED**13.13. Lateral Positioning Wedge****Optional.**

Lateral Positioning Wedge \$16 ea.

Part: SET123465 x2

Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

13.14. Back Width Package**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections

Part: PTOINDV3468

14"W-16"W Package

Part: XRF000037

18"W-22"W Package

Part: XRF000039

10"W-14"W Package

Part: XRF000036

16"W-20"W Package

Part: XRF000038

20"W-24"W Package

Part: XRF000040

13.15. Back Width**REQUIRED** except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width

Part: PTOINDV3562

16" Back Width

Part: PTOINDV3564

19" Back Width

Part: PTOINDV3567

22" Back Width

Part: PTOINDV3570

14" Back Width

Part: PTOINDV3563

17" Back Width

Part: PTOINDV3565

20" Back Width

Part: PTOINDV3568

15" Back Width

Part: PTOINDV3856

18" Back Width

Part: PTOINDV3566

21" Back Width

Part: PTOINDV3569

13.16. Overall Back Height**Static Seat and Power Tilt available height range:** 14"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion. **Power Recline available height range:** 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

14" Height

Part: SETHGHT1027

17" Height

Part: SETHGHT1030

20" Height

Part: SETHGHT1033

23" Height

Part: SETHGHT1036

26" Height

Part: SETHGHT1039

29" Height

Part: SETHGHT1042

15" Height

Part: SETHGHT1028

18" Height

Part: SETHGHT1031

21" Height

Part: SETHGHT1034

24" Height

Part: SETHGHT1037

27" Height

Part: SETHGHT1040

16" Height

Part: SETHGHT1029

19" Height

Part: SETHGHT1032

22" Height

Part: SETHGHT1035

25" Height

Part: SETHGHT1038

28" Height

Part: SETHGHT1041

13.17. Power Tilt Privacy Flap**Optional.** A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap

Part: XRF000351

\$48

13.18. Optional Back Height Extension**Optional.** Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Optional Back Height Extension

Part: XRF000041

\$48

14. TRU BALANCE 3 ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**14.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ ☐	Flip-Back, 8"-16"H Armrests	Standard	☐ ☐	Flip-Back, 8"-16"H Armrests for Static Seat	\$289
☐ ☐	Left part: ARM1806662. Right part: ARM1806663		☐ ☐	Left part: XRF000320. Right part: XRF000321	
Not available with TRU Balance 3 Static Seat.					

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ ☐	Flip-Back Armrest Downpost	\$92 ea.
☐ ☐	Left part: XRF000265. Right part: XRF000266	

14.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests**Only available with TRU Balance 3 Static Seat.** Includes flat armpads.

☐ ☐	Full Length, Drop-In Armrests	No Charge	☐ ☐	Desk Length, Drop-In Armrests	No Charge
☐ ☐	Left part: ARM1803253. Right part: ARM156742		☐ ☐	Left part: ARM1803254. Right part: ARM156743	

14.1c. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ ☐	Full Length, Two-Post Armrest for Power Seating	No Charge	☐ ☐	Desk Length, Two-Post Armrest for Power Seating....	No Charge
☐ ☐	Left part: XRF000047. Right part: XRF000048		☐ ☐	Left part: XRF000049. Right part: XRF000050	

1b. Two-Post Armrest for Static Seating

☐ ☐	Full Length, Two-Post Armrest for Static Seating	\$158 ea.	☐ ☐	Desk Length, Two-Post Armrest for Static Seating	\$158 ea.
☐ ☐	Left part: XRF000051. Right part: XRF000052. HCPCS: E0973		☐ ☐	Left part: XRF000053. Right part: XRF000054. HCPCS: E0973	

TRU Balance 3 armrests section continued on following pages.

14. TRU BALANCE 3 ARMRESTS CONTINUED**14.1d. Height Adjustable, Removable, Heavy Duty Armrests**

Select a type of HD removable armrests and the set height in the following section. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$158
Heavy Duty (Sold in pairs)

Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable, \$158
Heavy Duty (Sold in pairs)

Part: XRF000056. HCPCS w/ Static Seat: E0973

14.2. Armrest Height

REQUIRED. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H

PTOINDV3535

10" H

PTOINDV3537

12" H

PTOINDV3539

14" H

PTOINDV3541

16" H

PTOINDV3543

9" H

PTOINDV3536

11" H

PTOINDV3538

13" H

PTOINDV3540

15" H

PTOINDV3542

2. Right Armrest Height

8" H

PTOINDV3450

10" H

PTOINDV3452

12" H

PTOINDV3454

14" H

PTOINDV3456

16" H

PTOINDV3458

9" H

PTOINDV3451

11" H

PTOINDV3453

13" H

PTOINDV3455

15" H

PTOINDV3457

14.3. Armrest Size

REQUIRED. Select a left (L) and right (R) armrest size. This also affects armpad selection.



10" Length

L part: XRF000057. R part: XRF000123



14" Length

L part: XRF000058. R part: XRF000124



17" Length

L part: XRF000059. R part: XRF000125
Only available with Flip-Back, 8"-16" H
Armrests and straight arm pads.

15. TRU BALANCE 3 ARPADS

REQUIRED. Please note that waterfall and gel armpads may cause interference with the seat back.

15.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge

Part: XRF000060

Not available with Flip-back, 8"-16" H armrests.

Interface for Aftermarket Pad No Charge

Part: XRF000308

Only available with Flip-back, 8"-16" armrests.

15.1b. Foam Armpads

Straight, Armpad Standard

Left part: XRF000061. Right part: XRF000126



Waterfall, Armpad No Charge

Left part: XRF000062. Right part: XRF000127

Not available with Flip-back, 8"-16" armrests.

15.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195

Part: ARM1806807

Falcon Gel, Desk Length Left, Full Length Right.....\$171

Part: ARM1806810

Falcon Gel, Full Length Left, Desk Length Right.....\$171

Part: ARM1806813

Falcon Gel, Desk Length Left and Right\$147

Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22

Part: ARM136438

Falcon Gel, Desk Length Left, Full Length Right.....\$287.93

Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right.....\$287.93

Part: ARM136442

Falcon Gel, Desk Length Left and Right \$265.02

Part: ARM136437

16. TRU BALANCE 3 JOYSTICK BRACKET

REQUIRED.

Standard Fixed Receiver Standard

Part: XRF000063

Stealth Height Adjustable Receiver \$236

Part: XRF000064. HCPCS: K0108

Account Number: _____ Date: _____ Quote Name/Marked for: _____

17. TRU BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

17.1. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1033

Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1033

17.2. Q Cam Lock

Optional.

Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

18. LAP BELTS

REQUIRED. Select a lap belt below.

18.1a. Standard Lap Belts

38" Lap Belt
Part: ACCBELT1006

50" Lap Belt
Part: ACCBELT1000

60" Lap Belt
Part: ACCBELT1002

70" Lap Belt
Part: ACCBELT1003

80" Lap Belt
Part: ACCBELT1004

90" Lap Belt
Part: ACCBELT1001

18.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$181
Part: ST-IFRB5350

Stealth Positioning, 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

Stealth Positioning, 60" Single Pull, Push Button,\$50
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable).....\$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.

19. TRU BALANCE 3 LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (adds 2 actuator control requirement)
Part: XRF000108

Articulating Foot Platform (adds 1 actuator control requirement)
Part: XRF000110

Combined Legrests (adds 2 actuator control requirement)
Part: XRF000109

20. LEGRESTS

Optional. Please note the set lower leg length may vary based on your seat type and seat to floor height selections. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

20.1a. Base Mounted Foot Platform

May interfere with seat mounted legrests. Not available with center mount foot platforms.

Base Mounted Foot Platform Standard
Part: FRMASMB16077

Omit base mounted foot platform No Charge
Part: PTOINDV3674

Bell Shaped, Base Mounted Foot Platform\$247
Part: RIGASMB7120041

20.1b. Basic Center Mount Foot Platform

1. Basic Center Mount Foot Platform Assembly

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates\$180
Part: XRF000312. HCPCS: K0040

Long Extension, 5"x9" Footplates\$180
Part: XRF000314. HCPCS: K0040

Medium Extension, 5"x9" Footplates\$180
Part: XRF000313. HCPCS: K0040

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount.....\$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount\$158
Part: RIG2202997
Only available with Long Extension.

4. Basic Center Mount Additional Accessory

Basic Center Mount Footplate Connector\$16
Part: RIG2104240

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1c. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in later section. Select legrest length in following section.

1. Center Mount Foot Platform TypePower Articulating Foot Platform \$3,145
Part: PTOINDV3754. HCPCS: E1012+K0040Quick Adjustable Center Mount Foot Platform \$263
Part: XRF000311. HCPCS: K0040**2. Extension and Angle Adjustable Footplates**

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates
Part: RIG152577Medium Extension, 5"x11" Footplates
Part: RIG2404018Pediatric Extension, 5"x8" Footplates
Part: RIG2404015Long Extension, 5"x6" Footplates
Part: RIG151505Medium Extension, 5"x6" Footplates
Part: RIG150206Long Extension, 5"x8" Footplates
Part: RIG2404017Medium Extension, 5"x8" Footplates
Part: RIG2404016Long Extension, 5"x11" Footplates
Part: RIG2404019

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

Pediatric	Medium	Long
8.5"-12.5"	11"-16.5"	14"-20"
Quick Adj Center Mount Extension Range		

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019Small Calf Pads, Manually Adjustable Center Mount Only, \$158
6"Wx6"H (pair)
Part: XRF000156Large Calf Pads, Manually Adjustable Center Mount Only, \$158
8" taper to 6"Wx8"H (pair)
Part: XRF000157**3b. Center Mount Calf Panel**Small Calf Panel, 18"Wx8"H \$184
Part: FRMASMB7110012Large Calf Panel, 24"Wx8"H \$237
Part: RIG123825
Only available on widths of 20" or greater.**4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Only available with the 5x6 footplates.Outside Edge for Width Extension \$69
Part: RIG2404922
Available with footplates and width extension.Width Extension for 8" \$69
Part: RIG2404919Front and Rear Cap \$16
Part: RIG2404923Width Extension for 11" \$69
Part: RIG2404920

Available with 5"x8" footplates.

Outside Edge \$69
Part: RIG2404921
Not available with width extensions or 5"x6" footplates.**5. Center Mount Additional Accessories**Small Foam Calf Pad Cover (pair) \$67
Part: ST-CSCS0775-QFootplate Riser \$75
Part: XRF000348Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q

Not available with Long Extension. Select legrest length between 4" to 8" for Power Articulating and 4"-7" with Quick Adjustable, at 1" increments.

Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-QFootplate Riser Calf Panel \$50
Part: FRMASMB7110011Large Foam Panel Cover (pair) \$282
Part: ST-CSCP0927-Q

Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

Footplate Connector \$16
Part: RIGASMB7120023**6. Stealth Gel Footplate Covers**5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K01085"x8" Tapered Gel Footplate Covers for Width Extension \$247
Part: ST-AFPG5008E-Q. HCPCS: K01085"x11" Tapered Gel Footplate Covers \$270
Part: ST-AFPG5011T-Q. HCPCS: K01085"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-AFPG5011E-Q. HCPCS: K0108

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1d. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1a. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

- 60°-90° Swing-Away Legrest, Power Seating..... No Charge**
 Pair part: RIG120062. Left: FRMASMB17558.
 Right: FRMASMB17559
 Available with Power Positioning only. Not available with 11"-11" or 11"-16.5" extensions.
- 60°-90° Swing-Away Legrest, Short, Power Seating.... No Charge**
 Pair part: RIG136481. Left: RIG132607.
 Right: RIG132609.
 Available with Power Positioning only. For use only with 11"-13" or 11"-16.5" lower extensions.
- 60°-90° Swing-Away Legrest, Static Seat \$86 ea.**
 Pair part: XRF000082. Left: XRF000083.
 Right: XRF000084
 Available with Static Seat only. Not available with 11"-13" or 11"-16.5" extensions.
- 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.**
 Pair part: XRF000085. Left: XRF000086.
 Right: XRF000087.
 Available with Static Seat only. For use only with 11"-13" or 11"-16.5" lower extensions.

60°-90° Swing-Away Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Extruded, 11"-13"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Extruded, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Angle Adjustable, 13"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓
Multi-Axis, 11"-16.5"	!	✓	✓	✓	✓	✓	✓	✓	✓
Multi-Axis, 15.5"-19.5"	!	!	!	!	!	!	!	!	✓

- 11"-13" extension causes no interference at all seat-to-floor heights.
- ☒ This rigging configuration causes no interference issues at this seat-to-floor height.
- ☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2. Lower Extension

- 11"-13" Length**
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.
- 13"-16.5" Length**
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.
- 15.5"-19.5" Length**
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20".

- 5"Wx8"D Footplates with Heel Loops \$90 ea.**
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: E0951
- 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.**
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: E0951
- 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.**
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: E0951
- 8"Wx11"D² Footplates with Heel Loops \$90 ea.**
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

- Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020
- Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

20.1e. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

- Power Elevating Legrests, Pair..... \$2,516
 Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High
Power ELR: Static Seat, Recline Only			
8.5"-14.5" Lower Length	8.5"-14"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-14	12.5"-15"	12.5"-16"

Legrest Length Matrix Center Mount	Seat to Floor Option		
	Low	Middle	High
Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
14.5"-21" Lower Length	12.5"-15"	12.5"-16"	12.5"-17"

2. Lower Extension

- 8.5" to 14.5" Lower Length
 Part: XRF000344
- 12.5" to 21" Lower Length
 Part: XRF000347

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

- 4"Wx6"D, Pair
 Part: FRMASMB7448
- 5"Wx8"D, Pair
 Part: FRMASMB7450
- 6"Wx11"D¹, Pair
 Part: FRMASMB7452
- 5"Wx6"D, Pair
 Part: FRMASMB7449
- 6"Wx8"D¹, Pair
 Part: FRMASMB7451
- 8"Wx11"D², Pair
 Part: FRMASMB7453

Legrests section continued on following pages.

20. LEGRESTS CONTINUED**20.1f. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

Elevating Legrest Compatibility Matrix	Power Recline or Static Seat			Power Tilt			Power Lift		
	Low	Med.	High	Low	Med.	High	Low	Med.	High
Reinforced, 11"-14"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reinforced, 14"-17"	!	!	✓	!	✓	✓	✓	✓	✓
Angle Adjustable, 12.5"-16.5"	✓	✓	✓	✓	✓	✓	✓	✓	✓
Angle Adjustable, 15.5"-19.5"	!	!	!	!	!	✓	!	✓	✓

☒ This rigging configuration causes no interference issues at this seat-to-floor height.

☐ This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)

*The 6"x11" footplates may cause interference with seat widths less than 18". *The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D¹ Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² Footplates with Heel Loops \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

20.2. Legrest Accessories

Heel Loops (pair)..... \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$37
 Part: XRF000100

20.3. Set Legrest Extension Length

REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

21. TRU BALANCE 3 POWER POSITIONING ELECTRONICS**REQUIRED** with power seating or legrest options. Select one option from section 21.1 or 21.3.**21.1a. Tilt Only through Toggle**

Single Actuator Control, Power Tilt Only, Through Toggle No Charge

Part: ELE1905907

Available with QLNES. Select toggle side. May require extended lead time when selected with Fender Lights.

1. Toggle Mount Side Selection

Left Side Toggle

Part: PTOINDV1839

Right Side Toggle

Part: PTOINDV1838

21.1b. Actuator Control through Joystick**1. Actuator Control Selection**Single Actuator Through Controller and\$2,163
Interface

Part: XRF000335. HCPCS: E2310

Multiple Actuator Through Controller and \$3,040
Interface

Part: XRF000336. HCPCS: E2311

Single Actuator Through Controller and \$3,040
Interface (AAM)

Part: XRF000068. HCPCS: E2310

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations.

4-Function Push Button

Part: CTL123218

See section 21.2 for switch mounting.

4-Way Toggle

Part: ELEASMB5925

4-Way Button Toggle

Part: ELEASMB5926

3. User Controls Side Selection**Required with 4-Way Toggle and 4-Way Button Toggle.**

Mount User Controls on Left

Part: PTOINDV3996

Mount User Controls on Right

Part: PTOINDV3997

21.1c. Actuator Control with iAccess SwitchSingle Actuator Control Through Flat Button\$2,547
iAccess and Controller and Interface

Part: XRF000337. HCPCS: E2310

Single Actuator Control Through Toggle Button\$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000074. HCPCS: E2310

Single Actuator Control Through Toggle Button\$2,547
iAccess and Controller and Interface

Part: XRF000338. HCPCS: E2310

Multiple Actuator Control Through Flat Button\$3,460
iAccess and Controller and Interface

Part: XRF000339. HCPCS: E2311

Single Actuator Control Through Flat Button\$3,460
iAccess and Controller and Interface (AAM)

Part: XRF000070. HCPCS: E2310

Multiple Actuator Control Through Toggle Button.....\$3,460
iAccess and Controller and Interface

Part: XRF000340. HCPCS: E2311

21.2. 4-Function Push Button/iAccess Mount**REQUIRED with 4-Function Push Button or iAccess options.** Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge

Part: MEC1907500

Attendant Mount \$26

Part: MEC1906564

Under Arm Mount, Flip-Back Armrest..... No Charge

Part: MEC1906531

Under Arm Mount, Flip-Back Armrest, plus iAccess Guard..... \$26

Part: MEC2001392

Under Arm Mount, Tube Style Armrest No Charge

Part: MEC1906530

Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26

Part: MEC2001391

1. Under Arm Mount Side Selection

Left Side Under Arm Mount

Part: PTOINDV3678

Right Side Under Arm Mount

Part: PTOINDV3679

21.3. Power Seating through Switches

Power seating through switches (mounting bracket included) is special order\$495

Part: PTOINDV3932

22. TRU BALANCE 3 ATTENDANT CONTROL**Optional.** *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525

Part: XRF000078. HCPCS: E2331

Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module

Part: XRF000079

23. SEAT MOUNTED LED LIGHTS**Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045


Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.



24. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.**Accessory Bar** \$132
Part: XRF000128

Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder.

Backpack Holder (Pair) (A)..... \$69
Part: ACC147774**Cane/Crutch Holder** \$105
Part: ACC152515. HCPCS: E2207**Cell Phone/ Media Holder, X-Grip (E)**..... \$69
Part: ACC155835
Not available with Heavy Duty armrests.**Clothing Guard (Pair)**..... \$237
Part: ACC125007**Compact Vent Tray** \$1,465
Part: ACC133889. HCPCS: E1029**Cup Holder** No Charge
Part: ACC125003
For use with track-style armrests.**Cup Holder (Tube-style Armrest)**..... No Charge
Part: ACCASMB2634
For use with tube-style armrests.**Flag Kit** \$69
Part: ACC1700443**Flip-Back Armrest, Side Track Accessory Bracket ...** \$11
Part: MEC1809178**Glove Box (C)**..... \$158
Part: ACC125011**Grip Drink Holder (20085)**..... \$15.50
Part: ACC152993**Grip Lap Board (20005)**..... \$38.95
Part: ACC152992**Hydration System (B)**..... \$237
Part: ACC125010**Large, Black, Weather Cover** \$202
Part: ACCCOVR1021**Oxygen Holder**..... \$210
Part: ACC125006. HCPCS: E2208
Includes accessory bar, a credit will be applied PTOIN-DV3769.**Personal Item Hook**..... \$69
Part: ACC125008**Personal Item Pouch** \$84
Part: ACC2101017**Privacy Bag**..... \$42
Part: ACC2005723**Push Handles (pair)**..... \$195
Part: ACC2005027
Not available with oxygen holder or vent tray.**Q-Logic EX/NE Programmer** \$573
Part: ELEASMB6976**Quantum Backpack** \$168
Part: ACC1905853**Rear View Mirror** \$69
Part: ACC155833**Reflector Kit** \$53
Part: ACC124670**Stealth Large Tablet Holder** \$625.28
Part: ACC141818**Transfer Handle (D)**..... \$105 ea.
Part: MEC2205680 x2. HCPCS: E0705**Depth Adjustable Vent Tray** \$595
Part: ACC2302500. HCPCS: E1029
Includes accessory bar, a credit will be applied PTOIN-DV3769.**Remote Stop Switch**..... \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046**Aftermarket Accessory Mounting Bracket, Non-slotted ..** \$53 ea.
Part: ACC2405496**25. BRANCH PATCHES****Optional.****Air Force Patch** \$21
Part: ACC140240**Army Patch** \$21
Part: ACC140241**Coast Guard Patch** \$21
Part: ACC140239**Marines Patch**..... \$21
Part: ACC140243**Navy Patch** \$21
Part: ACC140244**Purple Heart Patch** \$21
Part: ACC140238

26. POSITIONING COMPONENTS**Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**26.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q).... \$664
 Part: ST-ULS4042000-Q. HCPCS: E1034
 Only available with Redesigned back assemblies.

Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664
 Part: ST-ULPK404108B00-Q. HCPCS: E1034
 Only available with Sport or TRU Comfort backs.

Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664
 Part: ST-ULPK404104B00-Q. HCPCS: E1034
 Only available with back cane and omit back configurations.

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332
 Part: ST-ULPK404103L00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332
 Part: ST-ULPK404103R00-Q. HCPCS: E1034
 Only available with Active Contour Q-Backs.

2. Pads

Curved Pads, 3"x5", pair (TWBL3X5)..... \$250
 Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL4X6)..... \$250
 Part: INDPART2664. HCPCS: E0956

Flat Pads, 3"x5", pair (ULP846050320)..... \$250
Coming soon! Part: ST-ULP846050320-Q. HCPCS: E0956

Flat Pads, 4"x6", pair (ULP835050320) \$250
Coming soon! Part: ST-ULP835050320-Q. HCPCS: E0956

3. Installation

Install Laterals
 Part: PTOINDV3499

Do Not Install
 Part: PTOINDV3500

26.2. Stealth Unilink Pelvic/Thigh Guide**1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable, \$265 ea.
 Removable Mounting Hardware
 Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

UniLink, 7" Height Adjustable, \$276 ea.
 Removable Mounting Hardware
 Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Left \$358.87 ea.
 Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^

UniLink Multi-Angle Thigh Guide, Right..... \$358.87 ea.
 Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Left..... \$332 ea.
 Part: ST-UTS25100-L-Q. HCPCS: E1028/E1034^

UniLink Flip Down Hip/Thigh Support, Right \$332 ea.
 Part: ST-UTS25100-R-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.
 Feature (TWBADD-L & LPHW407)
 Part: POS140428. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.
 Feature (TWBADD-R & LPHW407)
 Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea.
 Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 8" Pad..... \$106 ea.
 Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

UniLink, 4" x 12" Pad \$106 ea.
 Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link \$81 ea.
 Part: ST-UIL001150A000-Q

UniLink, 2" Link..... \$81 ea.
 Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left
 Part: PTOINDV3585

Install on Right
 Part: PTOINDV3586

Do Not Install
 Part: PTOINDV3719

26.3. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154
 Part: ST-IFPS1103-Q. HCPCS: E0960

Upper Classic, Large (IFPS1104) \$154
 Part: ST-IFPS1104-Q. HCPCS: E0960

Upper Contour, Medium (IFPS2103)..... \$154
 Part: ST-IFPS2103-Q. HCPCS: E0960

Upper Contour, Large (IFPS2104) \$154
 Part: ST-IFPS2104-Q. HCPCS: E0960

**26.4. Stealth Positioning Chest Straps**

Chest Strap, Medium w/ Belt Anchor, \$97
 (IFCS4113)
 Part: ST-IFCS4113-Q. HCPCS: E0960

Chest Strap, Large w/ Belt Anchor, \$97
 (IFCS4114)
 Part: ST-IFCS4114-Q. HCPCS: E0960



Positioning components section continued on following page.

26. POSITIONING COMPONENTS CONTINUED**26.5. Stealth Swing-Away Headrest Mounts**


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right..... \$318
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left \$318
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$434
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028

Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and AccessoriesEgg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815Facial Pad, 4.5"x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108Facial Pad, 6"x2", Cool Core (920M-Q) \$61 ea.
Part: ST-920M-Q. HCPCS: K0108Swing-Away Joystick Mount, 18" Rod (C) \$276 ea.
Part: ST-SUS18-Q

Allows mounting of facial pads and attaches to the headrest



Headrest Accessories

26.6. Stealth Elbow Support**1. Elbow Support Bracket**Left, Swing-Away, Elbow Support \$199
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Left, UniLink, Fixed, Elbow Support Bracket \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-QRight, Swing-Away, Elbow Support \$199
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Right, UniLink, Fixed, Elbow Support Bracket \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4) \$57 ea.
Part: POS124067. HCPCS: K0108

Elbow Supports

26.7. Stealth Foot BoxesStealth Small Foot Boxes, 5"Wx8"Dx7"H \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954Stealth Large Foot Boxes, 7"Wx12"Dx11"H \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954**26.8. Stealth Positioning Shoe Holder**Small Shoe Holder, 3/4" Webbing (IFSH4211) \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

27. LAP TRAY**27.1. Quantum Lap Tray**

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
 Part: MEC1905328-11B. HCPCS: E0950
 For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
 Part: MEC1905328-13B. HCPCS: E0950
 For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
 Part: MEC1905328-15B. HCPCS: E0950
 For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
 Part: MEC1905328-17B. HCPCS: E0950
 For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
 Part: MEC1905328-19B. HCPCS: E0950
 For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
 Part: MEC1905328-11C. HCPCS: E0950
 For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
 Part: MEC1905328-13C. HCPCS: E0950
 For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
 Part: MEC1905328-15C. HCPCS: E0950
 For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
 Part: MEC1905328-17C. HCPCS: E0950
 For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
 Part: MEC1905328-19C. HCPCS: E0950
 For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
 Part: MEC2104886

27.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	SM	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

:Compatible

:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

28. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**28.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Deep Purple \$63 Part: DGN137729-PPR5	Electric Blue \$63 Part: DGN137729-PBU9	Military Digital Camo \$63 Part: DGN137729-WCM7
Red Corvette \$63 Part: DGN137729-PRD1	Raspberry Beret \$63 Part: DGN137729-PPK4	Blue Jean \$63 Part: DGN137729-PBU11	Faded Glory \$63 Part: DGN137729-WFL2
Back in Black \$63 Part: DGN137729-PBK1	Pretty in Pink \$63 Part: DGN137729-PPK3	Lemon Crush \$63 Part: DGN137729-PYL4	Carbon Fiber \$63 Part: DGN137729-WCB1
Night Train \$63 Part: DGN137729-PBK3	Evolution Orange \$63 Part: DGN137729-POR7	Green Machine \$63 Part: DGN137729-PGN11	
Grey Street \$63 Part: DGN137729-PSL3	Midnight Blue \$63 Part: DGN137729-PBU10	Military Green \$63 Part: DGN137729-PGN8	



Enhanced Display

28.2. SCIM Module**Optional.**☐ **SCIM Module for Q-Logic 3** \$525
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

28.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ** \$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

28.4. Head Array

i-Drive 4.0 Bluetooth Tri-array Series, 3 Switch with Egg Switch \$4,774 Part: ST-IDHBT200-U-Q. HCPCS: E2330	Adult, Curved, i-Drive Tri-Array with Mounting Hardware \$4,774 Part: ST-IDHBT400-U-Q. HCPCS: E2330
i-Drive 4.0 Bluetooth Pro Series, 3 Switch with Egg Switch \$4,825 Part: ST-IDHBT300-U-Q. HCPCS: E2330	Junior, Curved, i-Drive Tri-Array with Mounting Hardware \$4,774 Part: ST-IDHBT450-U-Q. HCPCS: E2330



Adult Pro-Series

28.5. Sip and Puff Drive**1. Sip and Puff Module****Sip and Puff Module for Q-Logic 3** \$1,990
Part: ST-SP625M-U-Q. HCPCS: E2325
Includes module and mounting hardware.

Sip & Puff Module

2. Sip and Puff Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, Left (SP625M-CL) \$318 Part: ST-SP625M-CL-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.	Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U) \$318 Part: ST-SP625-U-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.
--	--

Stealth Sip and Puff for Comfort Series, Right (SP625M-CR) \$318
Part: ST-SP625M-CR-Q. HCPCS: E2326
Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip & Puff

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)** \$35
Part: ST-SP625-S12-Q

Alternative controls and electronics section continued on following page.

28. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED**28.6. Egg Switches**

Red Egg Switch (ESRED)..... \$87
Part: SWTMCRO1077

Blue Egg Switch (ESBLU)..... \$87
Part: SWTMCRO1079

Black Egg Switch (ESBLK)..... \$87
Part: SWTMCRO1081

Green Egg Switch (ESGRN) \$87
Part: SWTMCRO1078

Yellow Egg Switch (ESYEL) \$87
Part: SWTMCRO1080

**28.7. Switch Adapter**

2 to 1 Switch Adapter..... \$15
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

28.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057

28.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
Part: ELE168069

Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

