



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Order Form with TRU Balance® 4 Positioning

Account Number: _____ Date: _____

Provider Name: _____

Contact: _____

Phone: _____ Fax: _____

Email: _____


PO Number: _____

Marked for: _____

Ship to Address: _____

City: _____ State: _____ Zip: _____

EDGE 3



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!


PATIENT INFORMATION

Weight: _____ Height: _____

A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____

B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____

C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and 8 amp off board charger. Dimensions of base without legrests: 24"Wx35.5"L

Quantum Q6 Edge 3 3SP-SS\$8,559

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option.

Quantum Q6 Edge 3 3MP-SS\$8,373

Part: Q6 EDGE 3 3MP-SS. HCPCS: K0861

Multiple actuator power option.

2. DRIVE MOTORS**REQUIRED.**

6 mph Drive Motors Standard
Part: MOT147433

6 mph, Accu-Trac Motor Technology \$1,360
Part: MOT147881. HCPCS: K0108
Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls.

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1801101-PRD1	 Deep Purple Part: DGN1801101-PPR5	 Midnight Blue Part: DGN1801101-PBU10	 Green Machine Part: DGN1801101-PGN11
 Back in Black Part: DGN1801101-PBK1	 Raspberry Beret Part: DGN1801101-PPK4	 Electric Blue Part: DGN1801101-PBU9	 Military Green Part: DGN1801101-PGN8
 Night Train Part: DGN1801101-PBK3	 Pretty in Pink Part: DGN1801101-PPK3	 Blue Jean Part: DGN1801101-PBU11	
 Grey Street Part: DGN1801101-PSL3	 Evolution Orange Part: DGN1801101-POR7	 Lemon Crush Part: DGN1801101-PYL4	

3.1b. Hydrographic Shrouds

☐  Military Digital Camo \$263 Part: DGN1801101-WCM7	☐  Faded Glory \$263 Part: DGN1801101-WFL2	☐  Carbon Fiber \$263 Part: DGN1801101-WCB1
--	--	---

4. BASE MOUNTED LIGHTS


Optional. Select one option. Seat Mounted, LED, Full Light Package option can be found later on this form.

Switch Fender Lights..... No Charge
Part: DGN154117
Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box.

Base Mounted, LED, Full Light Package \$1,045
Part: XRF000129
Requires Q-Logic 3 EX Joystick for Full Lights or Q-Logic 3e Joystick.



Base Mounted Lights



Fender Lights

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Mounting Bracket

Backup Camera Headrest Bracket..... \$21 Backup Camera Track Mount Bracket \$27
Part: MEC1910373 Part: MEC2104227
Works with .5" x .5" vertical supports only.

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

6. DRIVE WHEELS

REQUIRED. Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".

6.1. Drive Wheels

3"x14" Flat Free Black Tires Standard Part: PTOINDV3525 (pair)	3"x14" Pneumatic Black Tires No Charge Part: WHL1713925 x2
4"x14" Knobby, Flat Free Black Tires \$580 Part: XRF000368	4"x14" Knobby, Pneumatic Black Tires \$580 Part: XRF000367

6.2. Drive Wheel Accent Color

Optional. Not available with knobby tires.

Back in Black.....Standard Part: DGN164245-PBK1 x2	Pretty in Pink No Charge Part: DGN164245-PPK3 x2	Green Machine No Charge Part: DGN164245-PGN11 x2
Red Corvette.....No Charge Part: DGN164245-PRD1 x2	Evolution Orange.....No Charge Part: DGN164245-POR7 x2	Military Green No Charge Part: DGN164245-PGN8 x2
Night Train.....No Charge Part: DGN164245-PBK3 x2	Midnight Blue No Charge Part: DGN164245-PBU10 x2	Military Digital Camo..... \$58 ea. Part: DGN164245-WCM7 x2
Grey Street.....No Charge Part: DGN164245-PSL3 x2	Electric Blue No Charge Part: DGN164245-PBU9 x2	Faded Glory \$58 ea. Part: DGN164245-WFL2 x2
Deep Purple No Charge Part: DGN164245-PPR5 x2	Blue Jean No Charge Part: DGN164245-PBU11 x2	Carbon Fiber \$58 ea. Part: DGN164245-WCB1 x2
Raspberry Beret.....No Charge Part: DGN164245-PPK4 x2	Lemon Crush.....No Charge Part: DGN164245-PYL4 x2	

7. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376 x2. HCPCS: E2361

NF-22, Gel Cell, MK Battery \$432 ea.
Part: BATGEL1008 x2. HCPCS: E2361

NF-22, AGM, MK Battery..... \$432 ea.
Part: BAT121346 x2. HCPCS: E2361

8. PATIENT'S WEIGHT

REQUIRED. The patient's weight is required to quote and order the power chair.

Up to 200 lb User Weight
Part: PTOINDV3837

201 lb-250 lb User Weight
Part: PTOINDV3861

251 lb-300 lb User Weight
Part: PTOINDV3984

9. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

9.1. Expandable Controls - Quantum Intelligent Systems™

Select expandable controls or omit selection if another controller type is desired.

1. Expandable Controls Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected elsewhere, Customer Service will quote the order with these items.

Q-Logic 3 Controller Module \$730
Part: XRF000002. HCPCS: E2377

Q-Logic 3 Expandable Controller Harness..... \$458
Part: XRF000264. HCPCS: E2313

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick \$840
Part: CTL167724

Q-Logic 3 EX Joystick for Full Lights \$840
Part: CTL167725

Q-Logic 3e Joystick \$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.

No Joystick, Setup for Q-Logic
Part: XRF000010



9.2. Q-Logic 3 EX Joystick Armor

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2.

Black, Soft Touch Standard
Part: DGN124307
Also matches Night Train shrouds.

Pretty in Pink \$63
Part: DGN124336-PPK3

Green Machine \$63
Part: DGN124336-PGN11

Red Corvette..... \$63
Part: DGN124336-PRD1

Evolution Orange..... \$63
Part: DGN124336-POR7

Military Green \$63
Part: DGN124336-PGN8

Back in Black..... \$63
Part: DGN124336-PBK1

Midnight Blue \$63
Part: DGN124336-PBU10

Military Digital Camo..... \$63
Part: DGN124336-WCM7

Grey Street \$63
Part: DGN124336-PSL3

Electric Blue \$63
Part: DGN124336-PBU9

Faded Glory \$63
Part: DGN124336-WFL2

Deep Purple..... \$63
Part: DGN124336-PPR5

Blue Jean \$63
Part: DGN124336-PBU11

Carbon Fiber..... \$63
Part: DGN124336-WCB1

Raspberry Beret..... \$63
Part: DGN124336-PPK4

Lemon Crush..... \$63
Part: DGN124336-PYL4

10. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323

Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323

Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)

Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323

T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323


Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323

Medium Ball Handle (E)..... \$69
Part: ACCASMB1859. HCPCS: E2323


Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308

Soft Ball Handle (F)..... \$80
Part: ACCASMB1856. HCPCS: E2323



11. JOYSTICK MOUNTING BRACKETS**Optional.** Select a mount to be used with a joystick selected in previous section.**11.1a. Fixed, In-line Joystick Bracket**Fixed, In-line, Bracket for Q-Logic 3 Standard
Part: FRMASMB10617**11.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1032**11.1c. Basic Swing Away Joystick Mount**Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1032Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1032**11.1d. Specialty Joystick Mounts**Retract4 Mount \$368
Part: XRF000003. HCPCS: E1032
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.**11.2. Joystick Mount Location****Required except with basic swing away joystick mount.**Left Side Mount
Part: PTOINDV1225Right Side Mount
Part: PTOINDV1226**12. TRU BALANCE® 4 SEATING****REQUIRED.** Maximum seat size is 22"W x 22"D.**12.1. TRU Balance 4 Configuration****REQUIRED.**TRU-Balance 4 Power Tilt Only \$8,595
Part: PTOINDV4014. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system.TRU-Balance 4 Power Tilt and Recline \$12,790
Part: PTOINDV3954. HCPCS: E1007
50° of posterior tilt and 165° of recline.**12.2. 12" Power Lift****Requires combined selection of TRU Balance 4 Power Tilt and Recline and Articulating Foot Platform. Requires selection in 12.3 Anterior Tilt. ^Please verify the appropriate billing code to be used for this item with all third-party payors.**Power Lift with iLevel \$3,725
Part: PTOINDV3956. HCPCS: E2298/E2300^
12" scissor lift with positional feedback.**12.3. Anterior Tilt****Chest support is required with anterior tilt 20 and 30 degree functions. See chest support section. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.****1. Anterior Tilt Angle****Only available with Power Lift.**10° Anterior Tilt Angle No Charge
Part: MEC200935620° Anterior Tilt Angle \$1,995
Part: MEC2009399. HCPCS: K0108
Requires chest support and knee block.30° Anterior Tilt Angle \$2,595
Part: MEC2009400. HCPCS: K0108
Requires chest support and knee block.0° Anterior Tilt No Charge
Part: PTOINDV3980**12.4. Knee Block Assembly****REQUIRED with 20° or 30° anterior tilt. Optional with 10° and 0° anterior tilt.**Anterior Knee Block Assembly (ST-KBATB4ATC-1-Q) \$314
Part: ST-KBATB4ATC-1-Q. HCPCS: K0108Anterior Knee Block Assembly, Curved (ST-KBATB4ATC-U-Q) .. \$354
Part: ST-KBATB4ATC-U-Q. HCPCS: K0108**12.5. Occupied Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.WC-19 Occupied Transit Option \$289
Part: KIT152047To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**12.6. Seat Width****REQUIRED with seat selection. 12"-15" widths not available with Power Lift.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

12" W (#1002)	13" W (#1003)	14" W (#1004)	15" W (#1005)	16" W (#1006)	17" W (#1007)	18" W (#1008)	19" W (#1009)	20" W (#1010)	21" W (#1011)	22" W (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.7. Seat Width Package**REQUIRED with seat selection.**

12"-16" Width Package

Part: XRF000006

Available with Power Tilt Only.

16"-20" Width Package

Part: XRF000008

20"-24" Width Package

Part: XRF000158

Maximum width is 22".

12.8. Seat Depth**REQUIRED with seat selection. 12"-13" depth not available with Power Lift.** Maximum seat depth is 22". Part number is SETDPTH plus number following # below.

12" D (#1002)	13" D (#1003)	14" D (#1004)	15" D (#1005)	16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------	------------------

12.9. Seat Depth Package**REQUIRED with seat selection.****1a. Power Tilt Depth Package**

12"-16" Depth Package

Part: XRF000011

16"-20" Depth Package

Part: XRF000013

20"-24" Depth Package

Part: XRF000159

Maximum depth is 22".

1b. Power Lift Depth Package

16"-22" Depth Package

Part: XRF000309

12.10. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Low Seat to Floor

Part: PTOINDV3976

Middle Seat to Floor

Part: PTOINDV3559


Highest Seat to Floor

Part: PTOINDV3560

Lowest	Low	Med	High
Q6 Edge 3			
17.50"	18.25"	18.75"	19.25"
TRU Balance 4 Seat to Floor Heights			

12.11. Stealth Seat Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

TRU Balance 4 Seating section continued on following pages.

12. TRU BALANCE 4 SEATING CONTINUED**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
 Part: POP1909006-S
 14"W and smaller

Medium Wedge Positioning Kit.....\$72
 Part: POP1909006-M
 14"W and smaller

Large Wedge Positioning Kit..... \$105
 Part: POP1909006-L
 14"W and smaller

Small Wedge Positioning Kit\$39
 Part: POP1712340-S
 15"W and larger

Medium Wedge Positioning Kit.....\$72
 Part: POP1712340-M
 15"W and larger

Large Wedge Positioning Kit..... \$105
 Part: POP1712340-L
 15"W and larger

3. Cushion Top Foam

Visco Top Foam
 Reference part: XRF000380
 FRG Top Foam
 Reference part: XRF000381

No Top Foam
 Reference part: XRF000382

4. Cushion Liner

Stealth Urethane Liner
 Reference part: XRF000383
 Stealth Polyurethane Liner
 Reference part: XRF000384

No Liner
 Reference part: XRF000385

5. Cushion Cover

CoolCore Cover
 Reference part: XRF000386
 Mesh Cover
 Reference part: XRF000387

Stealth Tek Cover
 Reference part: XRF000388

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	-	✓	-	✓	-	✓	✓	✓	-
Essence	-	✓	-	✓	-	✓	✓	✓	-
Solution	-	✓	-	-	✓	✓	✓	✓	✓
Spectrum Foam	-	✓	✓	-	✓	✓	✓	✓	✓
Spectrum Gel	-	✓	✓	-	✓	✓	✓	✓	✓
Premiere	-	✓	✓	✓	✓	✓	✓	-	-
Glacial	✓	-	-	✓	✓	✓	✓	✓	✓
Zen	-	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	-	-	✓	-	-	-	-	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389
 MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 QTY Reference Part: XRF000390

Additional Stealth Polyurethane Liner
 QTY Reference Part: XRF000391

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 QTY Reference Part: XRF000392
 Additional Mesh Cover
 QTY Reference Part: XRF000393

Additional Stealth Tek Cover
 QTY Reference Part: XRF000394

12.12. Back Type Selection**REQUIRED.** Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back
 Part: XRF000014
 Select a back cushion configuration in the subsequent section.

Back Cane Option
 Part: XRF000015
 Available only with TRU-Balance 3 Power Tilt (1 actuator) selections.
 Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back..... \$184
 Part: XRF000016
 Not available with power recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

12.13. Back Cane Selection (Back Cane Option Only)**REQUIRED when Back Cane Option is selected above.** Select type and height of back canes.**1. Back Cane Type**

Angle Adjustable Back Cane No Charge
 Part: MEC140266
 Flip Forward Back Cane \$53
 Part: XRF000273

2" Endomorph (1" per side) Back Cane \$74
 Part: XRF000269
 4" Endomorph (2" per side) Back Cane..... \$74
 Part: XRF000271

2. Back Cane Height

16"H Cane
 Part: FRMASMB15834
 17"H Cane
 Part: FRMASMB15835

20"H Cane
 Part: FRMASMB15838
 21"H Cane
 Part: FRMASMB15839

24"H Cane
 Part: FRMASMB15842
 25"H Cane
 Part: FRMASMB15843

TRU Balance 4 Seating section continued on following page.

12. TRU BALANCE 4 SEATING CONTINUED**12.14. Back Cushion Configuration**

REQUIRED. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge
Part: XRF000017

Complete following width package and size.

Omit Back Pan/Cushion No Charge
Part: PTOINDV1287

Select a following width package.

After-market Back Interface \$535
Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" to 20" Width \$744

Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855

Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

1c. Q-Back 2.0, AL Series, Active Contour

Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL14AS. HCPCS: E2613

Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL14AT. HCPCS: E2613

Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL16AS. HCPCS: E2613

Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL16AT. HCPCS: E2613

Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL16AX. HCPCS: E2613

Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16AXX. HCPCS: E2613

Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL18AS. HCPCS: E2613

Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL18AT. HCPCS: E2613

Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL18AX. HCPCS: E2613

Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18AXX. HCPCS: E2613

Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AL20AS. HCPCS: E2613

Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$568
Part: QB2AL20AT. HCPCS: E2613

Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$568
Part: QB2AL20AX. HCPCS: E2613

Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20AXX. HCPCS: E2613

1d. Q-Back 2.0, AL Series, Deep Contour

Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL14DS. HCPCS: E2620

Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL14DT. HCPCS: E2620

Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL16DS. HCPCS: E2620

Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL16DT. HCPCS: E2620

Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL16DX. HCPCS: E2620

Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16DXX. HCPCS: E2620

Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL18DS. HCPCS: E2620

Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL18DT. HCPCS: E2620

Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL18DX. HCPCS: E2620

Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18DXX. HCPCS: E2620

Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AL20DS. HCPCS: E2620

Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AL20DT. HCPCS: E2620

Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643
Part: QB2AL20DX. HCPCS: E2620

Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20DXX. HCPCS: E2620

1e. Q-Back 2.0, AL Series, Extra Deep Contour

Select a back type and the specific overall back height in subsequent section.

Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620

Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL16ES. HCPCS: E2620

Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL16ET. HCPCS: E2620

Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620

Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620

Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620

Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620

Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620

Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620

Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620

Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620

Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620

TRU Balance 4 Seating section continued on following page.

12. TRU BALANCE 4 SEATING CONTINUED**12.15. TRU Balance 4 Back Shroud**

REQUIRED. The back shroud is a feature of the TRU Balance 4 seating system and is available in matching shroud colors. The back shroud is not available when Omit Back Pan/Cushion is selected. Reference colors shown on page 2.

Red Corvette
Part: XRF000019

Deep Purple
Part: XRF000025

Midnight Blue
Part: XRF000029

Green Machine
Part: XRF000033

Back in Black
Part: XRF000021

Raspberry Beret
Part: XRF000026

Electric Blue
Part: XRF000030

Military Green
Part: XRF000034

Night Train
Part: XRF000022

Pretty in Pink
Part: XRF000027

Blue Jean
Part: XRF000031

Faded Glory
Part: DGN1714085-WFL2

Grey Street
Part: XRF000023

Evolution Orange
Part: XRF000028

Lemon Crush
Part: XRF000032

Digital Camo
Part: DGN1714085-WCM7

12.16. Lateral Positioning Wedge

Optional. Available with Sport and TRU-Comfort backs. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

Lateral Positioning Wedge \$16
Part: SET123465 x2

12.17. Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

Match Seat Width Selections
Part: PTOINDV3468

14"W-16"W Package
Part: XRF000037

16"W-20"W Package
Part: XRF000038

18"W-22"W Package
Part: XRF000039

12.18. Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

12" Back Width
Part: PTOINDV3562

16" Back Width
Part: PTOINDV3564

18" Back Width
Part: PTOINDV3566

20" Back Width
Part: PTOINDV3568

22" Back Width
Part: PTOINDV3570

14" Back Width
Part: PTOINDV3563

17" Back Width
Part: PTOINDV3565

19" Back Width
Part: PTOINDV3567

21" Back Width
Part: PTOINDV3569

		Tilt Only Width										Recline Width									
		12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"
Overall Height	14"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	14"
	15"	S	T	S	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15"
	16"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	16"
	17"	S	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	17"
	18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	18"
	19"	-	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	S	T	S	S	19"
	20"	-	S	Q	T	S	Q	T	S	Q	T	S	S	S	T	S	S	T	S	S	20"
	21"	-	S	Q	T	S	Q	T	S	Q	T	S	-	S	T	S	S	T	S	S	21"
	22"	-	S	Q	T	S	Q	T	S	Q	T	S	-	S	Q	T	S	S	Q	T	22"
	23"	-	S	T	S	T	S	Q	T	S	Q	T	-	S	Q	T	S	S	Q	T	23"
	24"	-	S	T	S	T	S	Q	T	S	Q	T	-	S	Q	T	S	S	Q	T	24"
	25"	-	S	T	S	Q	T	S	Q	T	S	-	S	Q	T	S	S	Q	T	S	25"
	26"	-	S	T	S	Q	T	S	Q	T	S	-	S	T	S	Q	T	S	S	Q	26"
	27"	-	S	T	S	T	S	T	S	T	S	-	S	T	S	Q	T	S	S	Q	27"
	28"	-	S	T	S	T	S	T	S	T	S	-	S	T	S	Q	T	S	S	Q	28"
	29"	-	S	T	S	T	S	T	S	T	S	-	S	T	S	Q	T	S	S	Q	29"
	30"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	30"
	31"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	31"
	32"	-	-	-	-	-	-	-	-	-	-	-	S	T	S	T	S	T	S	T	32"

S Sport Back
 Q Q-Back
 T TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

12.19. Overall Back Height

Power Recline available height range: 19"-32"H. Odd heights include 6" gap and even heights include 7" gap between seat pan and bottom of back cushion. See the availability chart for stock order back size ranges.

19" Height
Part: SETHGHT1032

22" Height
Part: SETHGHT1035

25" Height
Part: SETHGHT1038

28" Height
Part: SETHGHT1041

31" Height
Part: SETHGHT1102

20" Height
Part: SETHGHT1033

23" Height
Part: SETHGHT1036

26" Height
Part: SETHGHT1039

29" Height
Part: SETHGHT1042

32" Height
Part: SETHGHT1103

21" Height
Part: SETHGHT1034

24" Height
Part: SETHGHT1037

27" Height
Part: SETHGHT1040

30" Height
Part: SETHGHT1043

12.20. Optional Back Height Extension

Optional. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Optional Back Height Extension \$48
Part: XRF000041

13. ARMRESTS**REQUIRED.** Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.**13.1a. Flip-Back, 8"-16" H Armrests**

Choose a left and right armrest and set height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

☐ L	☐ R	Flip-Back, 8"-16"H Armrests Standard
Left part: ARM1806662. Right part: ARM1806663		

Armrest Support**Optional.** Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

☐ L	☐ R	Flip-Back Armrest Downpost \$92 ea.
Left part: XRF000265. Right part: XRF000266		

13.1b. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and set the left and right heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

☐ L	☐ R	Full Length, Two-Post Armrest for Power Seating No Charge	☐ L	☐ R	Desk Length, Two-Post Armrest for Power Seating.... No Charge
		Left part: XRF000047. Right part: XRF000048			Left part: XRF000049. Right part: XRF000050

13.1c. Pediatric, Removable, Height Adjustable ArmrestsHeight adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

☐ L	☐ R	Pediatric, Removable Armrests \$79 ea.
Left Part: ARM142876. Right part: ARM142877		

13.2. Armrest Height**REQUIRED.** Set left and right armrest height. Note the available height range in the option description.**1a. Flip-Back and Two-Post Left Armrest Height**

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

1b. Pediatric, Left Armrest Height

7.75" H PTOINDV3722	8.25" H PTOINDV3724	8.75" H PTOINDV3726	9.25" H PTOINDV3728	9.75" H PTOINDV3730
------------------------	------------------------	------------------------	------------------------	------------------------

2a. Flip-Back and Two-Post Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

2b. Pediatric, Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	8.75" H PTOINDV3725	9.25" H PTOINDV3727	9.75" H PTOINDV3729
------------------------	------------------------	------------------------	------------------------	------------------------

13.3. Armrest Size**REQUIRED.** Select a left (L) and right (R) armrest size. This also affects armpad selection.

☐ L	☐ R	10" Length L part: XRF000057. R part: XRF000123	☐ L	☐ R	14" Length L part: XRF000058. R part: XRF000124	☐ L	☐ R	17" Length L part: XRF000059. R part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
----------------------------	----------------------------	--	----------------------------	----------------------------	--	----------------------------	----------------------------	--

14. ARMPADS**REQUIRED.****14.1a. Omit Armpad**Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests. Not available with chest bar.
---	---

14.1b. Foam Armpads

☐ L	☐ R	Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	☐ L	☐ R	Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
----------------------------	----------------------------	--	----------------------------	----------------------------	--

14.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right.....\$171 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171 Part: ARM1806813	Falcon Gel, Desk Length Left and Right\$147 Part: ARM1806806

Armpads section continued on following page.

Account Number: _____ Date: _____ Marked for: _____

14. ARMPADS CONTINUED

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22
Part: ARM136438

Falcon Gel, Full Length Left, Desk Length Right.....\$287.93
Part: ARM136442

Falcon Gel, Desk Length Left, Full Length Right.....\$287.93
Part: ARM136439

Falcon Gel, Desk Length Left and Right \$265.02
Part: ARM136437

15. ARMREST BRACKETS

REQUIRED.

Standard Fixed Receiver Standard
Part: XRF000063

Stealth Height Adjustable Receiver \$236
Part: XRF000064

16. HEADRESTS

Optional. Select a headrest pad and bracket.

16.1a. Headrests, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink.....\$311
Part: ST-UHB2D33200-Q. HCPCS: E1028

Multi-Axis Removable Headrest Bracket.....\$174
Part: ST-HMO475-17M-Q. HCPCS: E1028

16.2. Q Cam Lock

Q Cam Lock.....\$50
Part: XRF000376

Q Cam Lock Availability (Overall height)	TRU-Comfort/Sport		Q-Back 2	
	Standard	\$50	Standard	\$50
Tilt/Static/Back Canes	14"-21"H	22"-29"H	19"-24"H	25"-26"H
Power Recline Back	19"-21"H	22"-32"H	22"-25"H	26"-29"H

17. LAP BELTS

REQUIRED. Select a lap belt below.

17.1a. Standard Lap Belts

☐ 38" Lap Belt Part: ACCBELT1006 ☐ 50" Lap Belt Part: ACCBELT1000 ☐ 60" Lap Belt Part: ACCBELT1002 ☐ 70" Lap Belt Part: ACCBELT1003 ☐ 80" Lap Belt Part: ACCBELT1004 ☐ 90" Lap Belt Part: ACCBELT1001

17.1b. Optional Lap Belts

Stealth Positioning, 60" Single Pull, Push Button, Medium.....\$50
(7.5") Padded, 1.5" Belt (IFPB5213)
Part: ST-IFPB5213-Q

Stealth Positioning, 60" Single Pull, Push Button, Large\$50
Padded (9.5"), 2" Belt (IFPB5314)
Part: ST-IFPB5314-Q

Stealth Positioning, 60" Single Pull, Push Button,\$50
Extra Large, Padded (11.5"), 2" Belt (IFPB5315)
Part: ST-IFPB5315-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
Part: ACCBELT1016
Required with ISO 7176-19 Compliant Transit Option.

18. LEGRESTS

All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

18.1a. Basic Center Mount Foot Platform

1. Basic Center Mount Foot Platform Assembly

Basic Center Mount..... No Charge
Part: XRF000323

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

2. Basic Center Mount Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates \$180
Part: XRF000312. HCPCS: K0040

Medium Extension, 5"x9" Footplates \$180
Part: XRF000313. HCPCS: K0040

Long Extension, 5"x9" Footplates \$180
Part: XRF000314. HCPCS: K0040

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount..... \$158
Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
Part: RIG2202997
Only available with Long Extension.

4. Basic Center Mount Additional Accessory

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

Legrests section continued on following page.

18. LEGRESTS CONTINUED**18.1b. Center Mount Foot Platform****1. Foot Platform Assembly**

Dual Actuator, Power Articulating Foot Platform.....\$3,595

Part: RIG2009246. HCPCS: E1012+K0040

Positional feedback with independent control of articulation.

Quick Adjustable Center Mount Foot Platform.....\$263

Part: XRF000311. HCPCS: K0040

Not available with power lift.

2. Lower Extension and Angle Adjustable Footplates

Medium length extension range is 11"-17" and long extension range is 14"-20". See compatibility matrix for possible interference.

Pediatric Extension, 5"x8" Footplates

Part: RIG2404015

Medium Extension, 5"x11" Footplates

Part: RIG2404018

Long Extension, 5"x11" Footplates

Part: RIG2404019

Medium Extension, 5"x6" Footplates

Part: RIG150206

Long Extension, 5"x6" Footplates

Part: RIG151505

Medium Extension, 5"x8" Footplates

Part: RIG2404016

Long Extension, 5"x8" Footplates

Part: RIG2404017

Legrest Compatibility Articulating Foot Platform TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Medium Length Extension	11.5"-14.5"	11.5"-15.5"	11.5"-16.5"
Long Length Extension	N/A	14.5"-15.5"	14.5"-16.5"

Legrest Compatibility Quick Adjustable CM TRU Balance 4	Seat to Floor Option		
	Low	Middle	High
Pediatric Length Extension	8.5"-12.5"	8.5"-12.5"	8.5"-12.5"
Medium Length Extension	11"-15"	11"-16"	11"-17"
Long Length Extension	14"-15"	14"-16"	14"-17"

3a. Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard

Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard

Part: RIGASMB7120019

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension\$69

Part: RIG152519

Only available with the 5x6 footplates.

Outside Edge.....\$69

Part: RIG2404921

Available with 5"x8" footplates.

Width Extension for 8".....\$69

Part: RIG2404919

Front and Rear Cap.....\$16

Part: RIG2404923

Available with 5"x8" footplates.

Width Extension for 11"\$69

Part: RIG2404920

5. Center Mount Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair)\$158

Part: RIG133717

Footplate Connector\$16

Part: RIGASMB7120023

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers.....\$222

Part: ST-FPG5008T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension\$247

Part: ST-AFPG5008E-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers.....\$270

Part: ST-AFPG5011T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension\$270

Part: ST-AFPG5011E-Q. HCPCS: K0108

18.2. Legrest Accessories

Wheel Bumper.....\$37

Part: XRF000100

18.3. Legrest Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

11" L (#3478)	11.5" L (#3479)	12" L (#3480)	12.5" L (#3481)	13" L (#3482)	13.5" L (#3483)	14" L (#3484)	14.5" L (#3485)	15" L (#3486)	15.5" L (#3487)	16" L (#3488)	16.5" L (#3489)	17" L (#3490)
------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------	--------------------	------------------

19. POWER POSITIONING ELECTRONICS**REQUIRED.****19.1a. Q-Logic 3 Controls through Joystick**Single Actuator Through Controller and Interface \$2,163
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and Interface.....\$3,040
Part: XRF000336. HCPCS: E2311Single Actuator Through Controller and\$3,040
Interface (AAM)

Part: XRF000068. HCPCS: E2310

19.1b. Q-Logic 3 Controls with iAccess Switch ControlSingle Actuator Control Through Flat Button\$2,547
iAccess and Controller and Interface
Part: XRF000337. HCPCS: E2310Single Actuator Control Through Toggle Button\$2,547
iAccess and Controller and Interface
Part: XRF000338. HCPCS: E2310Single Actuator Control Through Flat Button\$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000070. HCPCS: E2310Single Actuator Control Through Toggle Button\$3,460
iAccess and Controller and Interface (AAM)
Part: XRF000074. HCPCS: E2310Multiple Actuator Control Through Flat Button\$3,460
iAccess and Controller and Interface
Part: XRF000339. HCPCS: E2311Multiple Actuator Control Through Toggle Button.....\$3,460
iAccess and Controller and Interface
Part: XRF000340. HCPCS: E2311**19.2. iAccess Mount****REQUIRED with iAccess options.****1. iAccess Mount Selection**Joystick Mount No Charge
Part: MEC1907500Under Arm Mount, Flip-Back Armrest..... No Charge
Part: MEC1906531Under Arm Mount, Tube Style Armrest No Charge
Part: MEC1906530Attendant Mount \$26
Part: MEC1906564Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26
Part: MEC2001392Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26
Part: MEC2001391**2. Under Arm Mount Side**Left Side Under Arm Mount
Part: PTOINDV3678Right Side Under Arm Mount
Part: PTOINDV3679**20. ATTENDANT CONTROL****Optional.** Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.Attendant Control \$525
Part: XRF000078. HCPCS: E2331


Select this option if no charge option criteria is not met.

Attendant Control with Tilt, Recline, Alternative No Charge
Controls and Q-Logic Display or SCIM Module
Part: XRF000079**21. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.Seat Mounted, LED, Full Light Package \$1,045
Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

22. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur.Accessory Bar \$132
Part: XRF000128Horizontal rail attached to rear of seat. Standard with
the Trilogy vent tray and oxygen holder.Backpack Holder (Pair) (A)..... \$69
Part: ACC147774Cane/Crutch Holder \$105
Part: ACC152515. HCPCS: E2207Clothing Guard (Pair)\$237
Part: ACC125007Compact Vent Tray\$1,465
Part: ACC133889. HCPCS: E1029Cup Holder No Charge
Part: ACC125003Flag Kit \$69
Part: ACC1700443Flip-Back Armrest, Side Track Accessory Bracket ...\$11
Part: MEC1809178Front Mounted, Cell Phone/ Media Holder, \$69
X-Grip (E)

Part: ACC155835

Glove Box (C)..... \$158
Part: ACC125011Grip Drink Holder (20085) \$15.50
Part: ACC152993Grip Lap Board (20005)..... \$38.95
Part: ACC152992Hydration System (B).....\$237
Part: ACC125010Large, Black, Weather Cover \$202
Part: ACCCOVR1021Oxygen Holder..... \$210
Part: ACC125006. HCPCS: E2208Includes accessory bar, a credit will be applied PTOIN-
DV3769.

Accessories section continued on following page.

22. ACCESSORIES CONTINUED

Personal Item Hook.....	\$69
Part: ACC125008	
Personal Item Pouch	\$84
Part: ACC2101017	
Privacy Bag.....	\$42
Part: ACC2005723	
Push Handles (pair).....	\$195
Part: ACC2005027	
Not available with oxygen holder or vent tray.	
Q-Logic EX/NE Programmer	\$573
Part: ELEASMB6976	
Quantum Backpack	\$168
Part: ACC1905853	

Rear View Mirror	\$69
Part: ACC155833	
Reflector Kit	\$53
Part: ACC124670	
Stealth Large Tablet Holder	\$625.28
Part: ACC141818	
Transfer Handle (D).....	\$105 ea.
Part: MEC1809930 x2. HCPCS: E0705	
Depth Adjustable Vent Tray	\$595
Part: ACC2302500. HCPCS: E1029	
Includes accessory bar, a credit will be applied PTOINDV3769.	
Remote Stop Switch	\$1,295
Part: XRF000379	

Aftermarket Accessory Mounting Bracket

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted	\$53 ea.
Part: ACC2407046	

Aftermarket Accessory Mounting Bracket, Non-slotted ..	\$53 ea.
Part: ACC2405496	

23. BRANCH PATCHES

Optional.

Air Force Patch	\$21
Part: ACC140240	
Army Patch	\$21
Part: ACC140241	

Coast Guard Patch	\$21
Part: ACC140239	
Marines Patch	\$21
Part: ACC140243	

Navy Patch	\$21
Part: ACC140244	
Purple Heart Patch	\$21
Part: ACC140238	

**24. POSITIONING COMPONENTS**Optional. More options are available on the [Positioning Components Order Form](#).**24.1. Stealth Unilink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

^ Please verify the appropriate billing code depending upon the placement of the support pad.

UniLink, 5" Height Adjustable & Removable Hardware	\$265 ea.
Part: ST-UTTP5-Q. HCPCS: E1028/E1034^	
Flip Down with LEFT Full Surface Contact	\$281 ea.
Feature (TWBADD-L & LPHW407)	
Part: POS140428. HCPCS: E1028/E1034^	

UniLink Multi-Angle Thigh Guide, Left	\$358.87 ea.
Part: ST-UTT8A11602-L-Q. HCPCS: E1028/E1034^	

UniLink, 7" Height Adjustable & Removable Hardware	\$276 ea.
Part: ST-UTTP7-Q. HCPCS: E1028/E1034^	
Flip Down with RIGHT Full Surface Contact.....	\$281 ea.
Feature (TWBADD-R & LPHW407)	
Part: POS140429. HCPCS: E1028/E1034^	

UniLink Multi-Angle Thigh Guide, Right.....	\$358.87 ea.
Part: ST-UTT8A11602-R-Q. HCPCS: E1028/E1034^	

2. UniLink Pelvic/Thigh Pads

^ Please verify the appropriate billing code to be used with all third-party payors.

UniLink, 4" x 4" Pad	\$106 ea.
Part: ST-UTP544056100-Q. HCPCS: E0956/E0953^	
UniLink, 4" x 8" Pad.....	\$106 ea.
Part: ST-UTP548056100-Q. HCPCS: E0956/E0953^	

UniLink, 4" x 12" Pad	\$106 ea.
Part: ST-UTP54D056100-Q. HCPCS: E0956/E0953^	

3. UniLink Pelvic/Thigh Guide Installation

Install on Left
Part: PTOINDV3585

Install on Right
Part: PTOINDV3586

Do Not Install
Part: PTOINDV3719

24.2. Stealth Thoracic Laterals

Swing-away laterals may interfere with flip back armrests. These laterals are not available with the back cane option.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000-Q)....	\$664
Part: ST-ULS4042000-Q. HCPCS: E1034	
Only available with Redesigned back assemblies.	
Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ...	\$664
Part: ST-ULPK404108B00-Q. HCPCS: E1034	
Only available with Sport or TRU Comfort backs.	
Swing-away, UniLink, Back Cane, pair (ULPK404104B00)	\$664
Part: ST-ULPK404104B00-Q. HCPCS: E1034	
Only available with back cane and omit back configurations.	

Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00)	\$332
Part: ST-ULPK404103L00-Q. HCPCS: E1034	
Only available with Active Contour Q-Backs.	
Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00)	\$332
Part: ST-ULPK404103R00-Q. HCPCS: E1034	
Only available with Active Contour Q-Backs.	
UniLink Flip Down Hip/Thigh Support, Left.....	\$332 ea.
Part: ST-UTS25100-L-Q. HCPCS: E1034	
UniLink Flip Down Hip/Thigh Support, Right	\$332 ea.
Part: ST-UTS25100-R-Q. HCPCS: E1034	

Positioning components section continued on following page.

24. POSITIONING COMPONENTS CONTINUED**2. Pads**Curved Pads, 3"x5", pair (TWBL3X5)..... \$250
Part: INDPART2663. HCPCS: E0956Curved Pads, 4"x6", pair (TWBL4X6)..... \$250
Part: INDPART2664. HCPCS: E0956**3. Installation**Install Laterals
Part: PTOINDV3499Do Not Install
Part: PTOINDV3500**24.3. Chest Bar Assembly****Chest support is required with anterior tilt 20 and 30 degree functions. Omitting a chest support selection requires the provider to facilitate aftermarket chest supports that best fits consumer needs.**Small Chest Bar, 14" Foam Pad (ST-CBATB4ATT-S-Q)..... \$368
Part: ST-CBA-TB4-S. HCPCS: K0108
Available with 14" to 20" back widths.Medium Chest Bar, 16" Foam Pad (ST-CBATB4ATT-M-Q)..... \$368
Part: ST-CBA-TB4-M. HCPCS: K0108
Available with 17" to 22" back widths.**24.4. Stealth Positioning Chest Harness**Upper Classic, Medium (IFPS1103) \$154 ea.
Part: ST-IFPS1103-Q. HCPCS: E0960Upper Contour, Medium (IFPS2103)..... \$154 ea.
Part: ST-IFPS2103-Q. HCPCS: E0960Upper Classic, Large (IFPS1104) \$154 ea.
Part: ST-IFPS1104-Q. HCPCS: E0960Upper Contour, Large (IFPS2104) \$154 ea.
Part: ST-IFPS2104-Q. HCPCS: E0960**24.5. Stealth Positioning Chest Straps**Chest Strap, Medium w/ Belt Anchor, (IFCS4113)..... \$97 ea.
Part: ST-IFCS4113-Q. HCPCS: E0960Chest Strap, Large w/ Belt Anchor, (IFCS4114)..... \$97 ea.
Part: ST-IFCS4114-Q. HCPCS: E0960**24.6. Stealth Swing-Away Headrest Mounts**

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral HardwareUniLink Facial Lateral, Swing Away, Right..... \$318
Part: ST-SAR2A-R-Q. HCPCS: E1028UniLink Facial Lateral, Swing Away, Left \$318
Part: ST-SAR2A-L-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Right
Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$492
2 Upper & 1 Lower Left
Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028UniLink Facial Lateral, Dual Swing Away \$434
2 Upper
Part: ST-SAR4A-DBU-Q. HCPCS: E1028**2. Facial Lateral Pads and Accessories**Egg Switch Mount (ESM) (A) \$57 ea.
Part: INDPART2815Facial Pad, 4".5x2", Cool Core (900M-Q) \$61 ea.
Part: ST-900M-Q. HCPCS: K0108Facial Pad, 6"x2", Cool Core (920M-Q)..... \$61 ea.
Part: ST-920M-Q. HCPCS: K0108Swing-Away Joystick Mount, 18" Rod (C)..... \$276 ea.
Part: ST-SUS18-Q
Allows mounting of facial pads and attaches to the headrest


Contoured and Classic Chest Harness




Chest Strap



Left and Right Facial Lateral with Pads



Headrest Accessories

24.7. Stealth Elbow Support**1. Elbow Support Bracket**Left, Swing-Away, Elbow Support \$199 ea.
Bracket Tube Style Arms (ARMS-ES-L)
Part: POS124064. HCPCS: E1028Left, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q)
Part: ST-UETB4021FL-QRight, Swing-Away, Elbow Support..... \$199 ea.
Bracket Tube Style Arms (ARMS-ES-R)
Part: POS124066. HCPCS: E1028Right, UniLink, Fixed, Elbow Support Bracket..... \$211
for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q)
Part: ST-UETB4021FR-Q**2. Elbow Support Pad**3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.
Part: POS124065. HCPCS: K01084"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.
Part: POS124067. HCPCS: K0108

Elbow Supports

24.8. Stealth Foot BoxesStealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539
(SPLE107) (pair)
Part: INDPART2786. HCPCS: E0954Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539
(SPLE109) (pair)
Part: INDPART2788. HCPCS: E0954Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539
(SPLE108) (pair)
Part: INDPART2787. HCPCS: E0954**24.9. Stealth Positioning Shoe Holder**Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.
Part: ST-IFSH4211-Q. HCPCS: E0951+E0952Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.
Part: ST-IFSH4213-Q. HCPCS: E0951+E0952Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.
Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

Positioning components section continued on following page.

25. LAP TRAY

25.1. Quantum Lap Tray

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

Black Lap Tray, 11" W \$200
Part: MEC1905328-11B. HCPCS: E0950
For use with 12"-13" back width.

Black Lap Tray, 13" W \$200
Part: MEC1905328-13B. HCPCS: E0950
For use with 14"-15" back width.

Black Lap Tray, 15" W \$200
Part: MEC1905328-15B. HCPCS: E0950
For use with 16"-17" back width.

Black Lap Tray, 17" W \$200
Part: MEC1905328-17B. HCPCS: E0950
For use with 18"-19" back width.

Black Lap Tray, 19" W \$200
Part: MEC1905328-19B. HCPCS: E0950
For use with 20"-22" back width.

1b. Clear Lap Tray

Clear Lap Tray, 11" W \$200
Part: MEC1905328-11C. HCPCS: E0950
For use with 12"-13" back width.

Clear Lap Tray, 13" W \$200
Part: MEC1905328-13C. HCPCS: E0950
For use with 14"-15" back width.

Clear Lap Tray, 15" W \$200
Part: MEC1905328-15C. HCPCS: E0950
For use with 16"-17" back width.

Clear Lap Tray, 17" W \$200
Part: MEC1905328-17C. HCPCS: E0950
For use with 18"-19" back width.

Clear Lap Tray, 19" W \$200
Part: MEC1905328-19C. HCPCS: E0950
For use with 20"-22" back width.

2. Mounting Hardware

Mounting Bracket, Chest Bar/Lap Tray \$75
Part: MEC2104886

25.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	M	M	M	M	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	M	M	M	M	L	L	
Drop-In, T-Arm	S	S	S	S	S	M	M	M	M	L	L	
Two-Post, Flip-Back	S	S	S	S	S	M	M	M	M	L	L	
Cantilever	S	S	S	S	S	M	M	M	M	L	L	
Q-Captain, Flip-Up						M		M		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	M	M	M	M	L	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.

26. ALTERNATIVE CONTROLS AND ELECTRONICS**Optional.** Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).**26.1. Enhanced Display****Q-Logic 3 Enhanced Display Kit**.....\$1,980
Part: CTL167726

Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

Q-Logic 3 Enhanced Display Shroud

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

Black, Soft Touch Standard Part: DGN137728	Deep Purple \$63 Part: DGN137729-PPR5	Electric Blue \$63 Part: DGN137729-PBU9	Military Digital Camo \$63 Part: DGN137729-WCM7
Red Corvette \$63 Part: DGN137729-PRD1	Raspberry Beret \$63 Part: DGN137729-PPK4	Blue Jean \$63 Part: DGN137729-PBU11	Faded Glory \$63 Part: DGN137729-WFL2
Back in Black \$63 Part: DGN137729-PBK1	Pretty in Pink \$63 Part: DGN137729-PPK3	Lemon Crush \$63 Part: DGN137729-PYL4	Carbon Fiber \$63 Part: DGN137729-WCB1
Night Train \$63 Part: DGN137729-PBK3	Evolution Orange \$63 Part: DGN137729-POR7	Green Machine \$63 Part: DGN137729-PGN11	
Grey Street \$63 Part: DGN137729-PSL3	Midnight Blue \$63 Part: DGN137729-PBU10	Military Green \$63 Part: DGN137729-PGN8	

26.2. SCIM Module**Optional.**☐ **SCIM Module for Q-Logic 3** \$525
Part: CTL2001963

Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

26.3. Q-Logic 3 SAJ**Q-Logic 3 EX SAJ** \$840
Part: CTL143980

Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick.



EX SAJ

26.4. Head Array**i-Drive 4.0 Bluetooth Pro Series, 3 Switch with Egg Switch** \$4,825
Part: ST-IDHBT300-U-Q. HCPCS: E2330

Adult Pro-Series

26.5. Sip and Puff Drive**1. Sip and Puff Module****Sip and Puff Module for Q-Logic 3** \$1,990
Part: CTL168056. HCPCS: E2325
Includes module and mounting hardware.

Sip & Puff Module

2. Sip and Puff Straw and Mount**The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.****Stealth Sip and Puff for Comfort Series, Left Mount (SP625-CL)** \$318
Part: ST-SP625-CL-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U) \$318
Part: ST-SP625-U-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.

Stealth Sip and Puff for Comfort Series, Right Mount (SP622-CR) \$318
Part: ST-SP625-CR-Q. HCPCS: E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.



Stealth Sip & Puff

3. Sip & Puff Accessories**Sip & Puff Straws, 12 pack (SP625-S12)** \$35
Part: ST-SP625-S12-Q**26.6. Egg Switches****Red Egg Switch (ESRED)** \$87
Part: SWTMCRO1077**Blue Egg Switch (ESBLU)** \$87
Part: SWTMCRO1079**Black Egg Switch (ESBLK)** \$87
Part: SWTMCRO1081**Green Egg Switch (ESGRN)** \$87
Part: SWTMCRO1078**Yellow Egg Switch (ESYEL)** \$87
Part: SWTMCRO1080

Egg Switch

Alternative controls and electronics section continued on following page.

26. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

26.7. Switch Adapter

2 to 1 Switch Adapter \$15
 Part: DWR1361H037
 Connects 2 single switches to a stereo jack to act as a dual switch.

26.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
 Part: CTL168057

26.9. USB Mobile Device Charger

USB Mobile Device Charger \$184
 Part: ELE168069
 Can charge 1 device. A multiplier harness is required and must be added when a Static Seat or Power Tilt Only seating system is chosen. The selection of more than 3 modules requires an additional multiplier.

