



Quantum Rehab®

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: quantumorders@pridemobility.com

Q6 Edge® 3 Stretto Order Form with Contoured Seating

Account Number: _____ Date: _____
Provider Name: _____
Contact: _____
Phone: _____ Fax: _____
Email: _____
PO Number: _____
Quote Name/Marked For: _____
Ship to Address: _____
City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. Select a group 3 power base. 300 lb. weight capacity with standard 6mph programmable motors and off board charger. Dimensions of base without legrests: 20.5"Wx35.5"L

Quantum Q6 Edge 3 Stretto 3S-C.....\$8,870
Part: Q6 EDGE 3 STRETTO 3S-C. HCPCS: K0849
Contoured seat (captain) standard.

Quantum Q6 Edge 3 Stretto 3S-SS.....\$8,870
Part: Q6 EDGE 3 STRETTO 3S-SS. HCPCS: K0848
Solid seat pan standard with no power option.

Quantum Q6 Edge 3 Stretto 3SP-C.....\$8,559
Part: Q6 EDGE 3 STRETTO 3SP-C. HCPCS: K0857
Contoured seat (captain), single power option for 12" lift.

Quantum Q6 Edge 3 Stretto 3SP-SS.....\$8,559
Part: Q6 EDGE 3 STRETTO 3SP-SS. HCPCS: K0856
Solid seat pan standard, single power option for 12" lift.

2. DRIVE MOTORS

Optional. Requires Q-Logic 3 electronics. The standard 6 mph motor can be upgraded to Accu-Trac motors that delivers enhanced traction and tracking, especially with switched drive controls.

Accu-Trac Motor Technology \$1,360
Part: XRF000268. HCPCS: K0108

3. BASE COLOR

REQUIRED. Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

3.1a. Standard Shrouds

 Red Corvette Part: DGN1900199-PRD1	 Deep Purple Part: DGN1900199-PPR5	 Midnight Blue Part: DGN1900199-PBU10	 Green Machine Part: DGN1900199-PGN11
 Back in Black Part: DGN1900199-PBK1	 Raspberry Beret Part: DGN1900199-PPK4	 Electric Blue Part: DGN1900199-PBU9	 Military Green Part: DGN1900199-PGN8
 Night Train Part: DGN1900199-PBK3	 Pretty in Pink Part: DGN1900199-PPK3	 Blue Jean Part: DGN1900199-PBU11	
 Grey Street Part: DGN1900199-PSL3	 Evolution Orange Part: DGN1900199-POR7	 Lemon Crush Part: DGN1900199-PYL4	

3.1b. Hydrographic Shrouds

Military Digital Camo.....\$263
Part: DGN1801101-WCM7



Faded Glory.....\$263
Part: DGN1801101-WFL2

4. POWER ADJUSTABLE SEAT HEIGHT

Optional. The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit. ^Please verify the appropriate billing code to be used for this item with all third-party payors.

12" Power Lift for Captain Seat with iLevel..... \$3,725

Part: PTOINDV3857. HCPCS: E2298/E2300^

Available with swing away legrests or manually adjustable center mount foot platform.


5. BASE MOUNTED LIGHTS

Optional. Select one option.

Switch Fender Lights..... No Charge

Part: DGN1810041

Fender lights operate by pressing included single switch button. Fender lighting is not available with Single Actuator Control through Toggle Box. Fender lights come standard with track mount hardware.



Fender Lights

6. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368
Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Captain's Headrest Bracket \$21
Part: MEC1901231

Backup Camera Track Mount Bracket \$27
Part: MEC2104227

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105
Part: MEC2106123



Backup Camera

7. DRIVE WHEELS

REQUIRED. Select a type of drive wheel. Accent colors are available for 14" tires.

7.1. Drive Wheels

3"x14" Pneumatic Tires No Charge
Part: WHL1713925

3"x14" Flat Free Tires No Charge
Part: WHL1704058

4"x14" Knobby, Flat Free Black Tires \$580
Part: XRF000368

2"x12.5" Pneumatic Tires No Charge
Part: WHL1810666

2"x12.5" Flat Free Tires No Charge
Part: WHL1810665

4"x14" Knobby, Pneumatic Black Tires \$580
Part: XRF000367

Drive wheels section continued on following page.

7. DRIVE WHEELS CONTINUED**7.2. Drive Wheel Accent Color for 14" Tires**

Not available with 12.5" tires.

Back in Black..... No Charge
Part: DGN164245-PBK1 x2Red Corvette..... No Charge
Part: DGN164245-PRD1 x2Night Train..... No Charge
Part: DGN164245-PBK3 x2Grey Street..... No Charge
Part: DGN164245-PSL3 x2Deep Purple..... No Charge
Part: DGN164245-PPR5 x2Raspberry Beret..... No Charge
Part: DGN164245-PPK4 x2Pretty in Pink..... No Charge
Part: DGN164245-PPK3 x2Evolution Orange.... No Charge
Part: DGN164245-POR7 x2Midnight Blue No Charge
Part: DGN164245-PBU10 x2Electric Blue No Charge
Part: DGN164245-PBU9 x2Blue Jean No Charge
Part: DGN164245-PBU11 x2Lemon Crush No Charge
Part: DGN164245-PYL4 x2Green Machine No Charge
Part: DGN164245-PGN11 x2Military Green No Charge
Part: DGN164245-PGN8 x2Military Digital Camo... \$58 ea.
Part: DGN164245-WCM7 x2Faded Glory \$58 ea.
Part: DGN164245-WFL2 x2**8. BATTERIES****Optional.** The powerchair requires 2 batteries and will be installed if ordered.NF-22, Gel Cell, Interceptor 2 Battery \$355 ea.
Part: BAT2004376. HCPCS: E2361
Uses 8 AMP chargerGel Cell, 40AH, Interceptor 2 Battery.....\$191 ea.
Part: BAT2100316. HCPCS: K0108
Uses 5 AMP chargerNF-22, AGM, MK Flat Top Battery \$432 ea.
Part: BAT2003982. HCPCS: E2361
Uses 8 AMP charger**9. PATIENT'S WEIGHT****REQUIRED.** The patient's weight is required to quote and order the power chair.Up to 150 lb User Weight
Part: PTOINDV3977151-250 lb User Weight
Part: PTOINDV3985250-300 lb User Weight
Part: PTOINDV3984**10. EXTENDED SEAT TOWERS****Optional.**Extended Seat TowersNo Charge
Part: MEC2105559**11. ELECTRONICS****REQUIRED.** Select one type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

QLNES 6 Key JoystickStandard
Part: CTL2205379
Controls up to 2 actuators.QLNES 10 Key Joystick\$420
Part: CTL2205380
Controls up to 2 actuators through joystick and seat mounted, LED, full light package.

QLNES Joystick

11.1b. Expandable Controls

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

1. Q-Logic 3 Required Components

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module\$730
Part: XRF000002. HCPCS: E2377Q-Logic 3 Expandable Controller Harness.....\$458
Part: XRF000264. HCPCS: E2313**2. Q-Logic 3 Joystick**

Q-Logic 3 EX Joystick includes XLR USB charger.

Q-Logic 3 EX Joystick\$840
Part: CTL167724Q-Logic 3 EX Joystick for Full Lights\$840
Part: CTL167725Q-Logic 3e Joystick \$735
Part: CTL1910600
Can operate full light package. Not available with alternative or specialty controls.No Joystick, Setup for Q-Logic 3
Part: XRF000010
Requires selection of Enhanced Display

Q-Logic 3 EX Joystick

11.2. Q-Logic Magnetic Charger

Available with QLNES. Included with Q-Logic 3e Joystick.

Magnetic to XLR Charger Adapter with Wall Holder \$60
Part: ELE2103635

Electronics section continued on following page.

11. ELECTRONICS CONTINUED**11.3. Q-Logic 3 Joystick Armor**

REQUIRED with Q-Logic 3 EX Joystick selections above. Not available with Q-Logic 3e Joystick. Reference colors shown on page 2. *Night train joystick armor will be the black, soft touch.

Black, Soft TouchStandard Part: DGN124307 Also matches Night Train shrouds.	Raspberry Beret..... \$63 Part: DGN124336-PPK4	Lemon Crush..... \$63 Part: DGN124336-PYL4
Red Corvette..... \$63 Part: DGN124336-PRD1	Pretty in Pink \$63 Part: DGN124336-PPK3	Green Machine \$63 Part: DGN124336-PGN11
Back in Black..... \$63 Part: DGN124336-PBK1	Evolution Orange..... \$63 Part: DGN124336-POR7	Military Green \$63 Part: DGN124336-PGN8
Grey Street \$63 Part: DGN124336-PSL3	Midnight Blue \$63 Part: DGN124336-PBU10	Military Digital Camo..... \$63 Part: DGN124336-WCM7
Deep Purple \$63 Part: DGN124336-PPR5	Electric Blue \$63 Part: DGN124336-PBU9	Faded Glory \$63 Part: DGN124336-WFL2
	Blue Jean \$63 Part: DGN124336-PBU11	

12. JOYSTICK HANDLES

Optional. The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.

☐ Mushroom Handle (B) \$90 Part: ACCASMB1860. HCPCS: E2323	☐ Extended I-Handle (G) \$90 Part: ACCASMB1857. HCPCS: E2323
☐ Chin Cup (C) \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	☐ Bodypoint 3" Goal Post Handle \$104 Part: PLSKNCP1057. HCPCS: E2323
☐ T-Bar (D) \$126 Part: ACCASMB1858. HCPCS: E2323	☐ Bodypoint 4" Goal Post Handle \$104 Part: PLSKNCP1058. HCPCS: E2323
☐ Medium Ball Handle (E)..... \$69 Part: ACCASMB1859. HCPCS: E2323	☐ Q-Logic 3 Carrot Handle (A) No Charge Part: DGN124308
☐ Soft Ball Handle (F)..... \$80 Part: ACCASMB1856. HCPCS: E2323	



Joystick Handles

13. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

13.1a. Fixed, In-line Joystick Bracket

Fixed, In-line, Bracket for Q-Logic 3Standard
Part: FRMASMB10617

13.1b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1028

13.1c. Basic Swing Away Joystick Mount

Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1028

Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1028

13.1d. Specialty Joystick Mounts

Select a bracket and side to mount it.

Retract4 Mount..... \$368
Part: XRF000003. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.


Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick.

Stealth Swing-Away/Flip Down Bracket \$263
Part: XRF000005. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.


13.2. Joystick Mount Location

Left Side Mount
Part: PTOINDV1225

Right Side Mount
Part: PTOINDV1226



Swing Away Joystick



Retract4



Multi-Axis Joystick Mount

14. CONTOURED SEATING**Optional.** Select only one seating type.**14.1. Occupied Transit, WC-19****Optional.**

WC-19 Occupied Transit Option for Contoured Seat \$368

Part: KIT122609

To determine if your chair is WC-19 compliant (occupied), please refer to this website: <https://wc-transportation-safety.umtri.umich.edu/>**14.2a. Q-Captain Seat**

16"W x 16"D No Charge

Part: SET2005145

18"W x 16"D No Charge

Part: SET2005149

20"W x 16"D No Charge

Part: SET2005152

16"W x 18"D No Charge

Part: SET2005146

18"W x 18"D No Charge

Part: SET2005150

20"W x 18"D No Charge

Part: SET2005153

16"W x 20"D No Charge

Part: SET2005147

18"W x 20"D No Charge

Part: SET2005151

20"W x 20"D No Charge

Part: SET2005154

14.2b. Q-Captain Seat with Solid Seat Pan

16"W x 16"D No Charge

Part: SET2005394

18"W x 16"D No Charge

Part: SET2005397

20"W x 16"D No Charge

Part: SET2005400

16"W x 18"D No Charge

Part: SET2005395

18"W x 18"D No Charge

Part: SET2005398

20"W x 18"D No Charge

Part: SET2005401

16"W x 20"D No Charge

Part: SET2005396

18"W x 20"D No Charge

Part: SET2005399

20"W x 20"D No Charge

Part: SET2005402

14.3. Seat to Floor Height**REQUIRED.**

Lowest Seat to Floor

Part: PTOINDV3558

Middle Seat to Floor

Part: PTOINDV3559


Highest Seat to Floor

Part: PTOINDV3560

	Cushion	SSP
High	22.5"	18.5"
Mid	21.5"	17.5"
Low	20.5"	16.5"
Q-Captain Seat to Floor Height		

14.4. Stealth Cushion

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions will be found on the respective complete [Stealth Cushion Order Forms](#). Part numbers in this section are for reference only. This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. The complete cushion part number and MSRP is a product of the cushion model, top foam, liner, cover, width, and depth as shown below. Reference the subsequent chart for available top foam, liner, and cover options for each cushion type. See specific cushion order form for all sizes, top foam, liners, covers, parts numbers, and MSRP.

**1. Cushion Type****TRU-Comfort 2 SPP**Reference part: XRF000116. [See cushion order form.](#)

Part series: TC2-###XXYY. MSRP: \$441-\$651. HCPCS: E2607/E2608

Essence SPPReference part: XRF000118. [See cushion order form.](#)

Part series: ESS-###XXYY. MSRP: \$441-\$630. HCPCS: E2607/E2608

Solution SPPReference part: XRF000113. [See cushion order form.](#)

Part series: SOL-###XXYY. MSRP: \$414-\$533. HCPCS: E2607/E2608

Spectrum Foam SPPReference part: XRF000112. [See cushion order form.](#)

Part series: SPCF-###XXYY. MSRP: \$507-\$709. HCPCS: E2607/E2608

Spectrum Gel SPPReference part: XRF000117. [See cushion order form.](#)

Part series: SPCG-###XXYY. MSRP: \$507-\$602. HCPCS: E2607/E2608

Premiere PReference part: XRF000122. [See cushion order form.](#)

Part series: PREM-###XXYY. MSRP: \$358-\$603. HCPCS: E2605/E2606

Glacial SPReference part: XRF000114. [See cushion order form.](#)

Part series: GLA-###XXYY. MSRP: \$364-\$564. HCPCS: E2603/E2604

Zen SPReference part: XRF000115. [See cushion order form.](#)

Part series: ZEN-###XXYY. MSRP: \$348-\$547. HCPCS: E2603/E2604

Simplicity GReference part: XRF000111. [See cushion order form.](#)

Part series: SIM-###XXYY. MSRP: \$122-\$182. HCPCS: E2601/E2602

Contoured seating section continued on following pages.

14. CONTOURED SEATING CONTINUED**2. Essence Cushion Wedge Kit**

Only available with Essence cushion selection above.

Small Wedge Positioning Kit\$39
 Part: POP1909006-S
 14"W and smaller

Medium Wedge Positioning Kit.....\$72
 Part: POP1909006-M
 14"W and smaller

Large Wedge Positioning Kit..... \$105
 Part: POP1909006-L
 14"W and smaller

Small Wedge Positioning Kit\$39
 Part: POP1712340-S
 15"W and larger

Medium Wedge Positioning Kit.....\$72
 Part: POP1712340-M
 15"W and larger

Large Wedge Positioning Kit..... \$105
 Part: POP1712340-L
 15"W and larger

3. Cushion Top Foam

Visco Top Foam
 Reference part: XRF000380

No Top Foam
 Reference part: XRF000382

FRG Top Foam
 Reference part: XRF000381

4. Cushion Liner

Stealth Urethane Liner
 Reference part: XRF000383

No Liner
 Reference part: XRF000385

Stealth Polyurethane Liner
 Reference part: XRF000384

5. Cushion Cover

CoolCore Cover
 Reference part: XRF000386

Stealth Tek Cover
 Reference part: XRF000388

Mesh Cover
 Reference part: XRF000387

	None	Visco/V.Elastic	FRG	None	Urethane	Polyurethane	CoolCore	Mesh	Stealth Tek
	Top Foam			Liner			Cover		
TRU-Comfort 2	—	✓	—	✓	—	✓	✓	—	—
Essence	—	✓	—	✓	—	✓	✓	✓	—
Solution	—	✓	—	—	✓	✓	✓	✓	✓
Spectrum Foam	—	✓	✓	—	✓	✓	✓	✓	✓
Spectrum Gel	—	✓	✓	—	✓	✓	✓	✓	✓
Premiere	—	✓	✓	✓	✓	✓	✓	—	—
Glacial	✓	—	—	✓	✓	✓	✓	✓	✓
Zen	—	✓	✓	✓	✓	✓	✓	✓	✓
Simplicity	✓	—	—	✓	—	—	—	—	✓

6. Optional Base Insert

Option selection available with Glacial, Premiere, Simplicity, Solution, Spectrum Foam, Spectrum Gel, and Zen cushions.

ABS Flat Rigidizer
 Reference part: XRF000389
 MSRP: \$48-\$69

7. Additional Cushion Liner

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional Stealth Urethane Liner
 Reference Part: XRF000390
 QTY _____

Additional Stealth Polyurethane Liner
 Reference Part: XRF000391
 QTY _____

8. Additional Cushion Cover

Up to 5 may be selected below. See respective cushion order form for MSRP and part numbers.

Additional CoolCore Cover
 Reference Part: XRF000392
 QTY _____

Additional Stealth Tek Cover
 Reference Part: XRF000394
 QTY _____

Additional Mesh Cover
 Reference Part: XRF000393
 QTY _____

14.5. Armrests**Optional.** Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. Q-Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat arm pads.

Full Length, Left Armrest\$90 ea.
 Part: MEC2201304. HCPCS: K0020

Full Length, Right Armrest.....\$90 ea.
 Part: MEC2201303. HCPCS: K0020

1b. Drop-In, Height Adjustable, Armrests (T-Arm)

Includes flat arm pads.

☐ ☐ Full Length, Drop-In Armrests.....\$79 ea.
 L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973

☐ ☐ Desk Length, Drop-In Armrests\$79 ea.
 L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973

1c. Single Post, Height Adjustable, Removable Armrests

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests\$79 ea.
 L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests.....\$79 ea.
 L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973

1d. Bariatric, Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Bariatric Armrests, 8"-12" H\$90 ea.
 L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973

☐ ☐ Full Length, Bariatric Armrests, 10"-14"H\$90 ea.
 L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973

1e. Heavy Duty, Height Adjustable, Removable Armrests

Select armrests and arm pads. Armrests are sold in same-size pairs.

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H\$90 ea.
 L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H.....\$90 ea.
 L R Left Part: XRF000257. Right part: XRF000258. HCPCS: E0973

Contoured seating section continued on following page.

14. CONTOURED SEATING CONTINUED**1f. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

L R	Pediatric, Removable Armrests \$79 ea. Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973
---	--

14.6. Q-Captain Armrest Spacer

Optional.

Armrest Spacer	\$35
-----------------------------	-------------

Part: XRF000370

1" spacer for Q-Captain, Flip-Up, Height Adjustable Armrests only.

14.7. Armrest Height**REQUIRED with Bariatric, Heavy Duty, and Pediatric Height Adjustable armrests.****1. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

14.8. Armpads**REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.****1a. Straight Armpads for Single Post Armrests**

Straight, Full Length Left and Right No Charge Part: ARM134275	Straight, Desk Length Left, Full Length Right No Charge Part: ARM134280
--	---

Straight, Full Length Left, Desk Length Right No Charge Part: ARM134281	Straight, Desk Length Left and Right No Charge Part: ARM134254
---	--

1b. Waterfall Armpads for Single Post Armrests

Waterfall, Full Length Left and Right No Charge Part: ARM134285	Waterfall, Desk Length Left, Full Length Right No Charge Part: ARM134288
---	--

Waterfall, Full Length Left, Desk Length Right No Charge Part: ARM134289	Waterfall, Desk Length Left and Right No Charge Part: ARM134282
--	---

1c. Falcon Gel Armpads for Single Post Armrests

Falcon Gel, Full Length Left and Right\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right\$287.93 Part: ARM136439
--	---

Falcon Gel, Full Length Left, Desk Length Right\$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437
---	---

1d. Armpads for Bariatric or Heavy Duty Armrests

Bariatric Armpads No Charge Part: UPSASMB1027 Only for use with Bariatric Armrests.	Waterfall Armpads No Charge Part: UPSASMB1024 Only for use with Heavy Duty Armrests
--	--

Straight Armpads No Charge Part: UPSASMB1021 Only for use with Heavy Duty Armrests

14.9. Quad Release Lever**Optional. Available with Single Post and Pediatric armrests.**

Quad Release Lever, Pair	No Charge
---------------------------------------	------------------

Part: FRMASMB10517

15. LAP BELTS**REQUIRED. Select one lap belt below.****15.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

15.1b. Optional Lap Belts

Stealth Positioning, 60" Single Pull, Push Button, Medium \$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning, 60" Single Pull, Push Button, Extra Large \$50 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
---	--

Stealth Positioning, 60" Single Pull, Push Button, Large \$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	Stealth Positioning, 60" Single Pull, Push Button, Non-padded ... \$50 2" Belt (IFPB5310) Part: ST-IFPB5310-Q
--	--

16. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

16.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Type

Basic Center Mount No Charge
Part: XRF000323

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

Short Extension, 5"x9" Footplates No Charge Medium Extension, 5"x9" Footplates No Charge
Part: MEC1906675 Part: MEC1906680

3. Basic Center Mount Calf Pads

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount \$158
Part: RIG2104369

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
Part: RIG2104240

16.1b. Center Mount Foot Platform

Select legrest length in following section.

1. Center Mount Foot Platform Assembly

Power Articulating Foot Platform \$3,145 Quick Adjustable Center Mount Foot Platform \$263
Part: PTOINDV3754. HCPCS: E1012+K0040 Part: XRF000311. HCPCS: K0040
Only available with power lift option.

2. Center Mount Lower Extension and Angle Adjustable Footplates

Pediatric extensions are not compatible with Power Articulating Foot Platform.

Pediatric Extension, 5"x6" Footplates Part: RIG152577	Medium Extension, 5"x8" Footplates Part: RIG2404016	Long Extension, 5"x8" Footplates Part: RIG2404017
Pediatric Extension, 5"x8" Footplates Part: RIG2404015	Medium Extension, 5"x11" Footplates Part: RIG2404018	Long Extension, 5"x11" Footplates Part: RIG2404019
Medium Extension, 5"x6" Footplates Part: RIG150206	Long Extension, 5"x6" Footplates Part: RIG151505	

3a. Center Mount Calf Pads

Calf pads are standard and required with a power articulating foot platform.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard Part: RIGASMB7120018	Small Calf Pads, Quick Adjustable Center Mount Only, \$158 6"Wx6"H (pair) Part: XRF000156
--	---

3b. Center Mount Calf Panel

Small Calf Panel, 18"Wx8"H \$184
Part: FRMASMB7110012

4. Center Mount Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

2" Bell-Shaped Length Extension \$69 Part: RIG152519 Only available with the 5x6 footplates.	Outside Edge \$69 Part: RIG2404921 Available with 5"x8" footplates.
Width Extension for 8" \$69 Part: RIG2404919	Front and Rear Cap \$16 Part: RIG2404923 Available with 5"x8" footplates.
Width Extension for 11" \$69 Part: RIG2404920	

5. Center Mount Additional Accessories

Small Foam Calf Pad Cover (pair) \$67 Part: ST-CSCS0775-Q	Footplate Riser \$75 Part: XRF000348 Select legrest length between 2" to 9" at 1" increments.
Small Foam Panel Cover (pair) \$176 Part: ST-CSCP0920-Q	Footplate Riser Calf Panel \$50 Part: FRMASMB7110011 Only available and required with Footplate Riser selection and extension lengths of 6" to 9".
Footplate Connector \$16 Part: RIGASMB7120023	

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222 Part: ST-FPG5008T-Q. HCPCS: K0108	5"x8" Tapered Gel Footplate Covers for Width Extension \$247 Part: ST-AFPG5008E-Q. HCPCS: K0108
5"x11" Tapered Gel Footplate Covers \$270 Part: ST-AFPG5011T-Q. HCPCS: K0108	5"x11" Tapered Gel Footplate Cover for Width Extension \$270 Part: ST-AFPG5011E-Q. HCPCS: K0108

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED**16.1c. Swing-Away Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. 60°-90° Swing-Away Legrest Upper

Adjustable in 5° increments.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Static Seat \$86 ea.
 Pair part: XRF000082. Left: XRF000083. Right: XRF000084
 Not available with 11"-13" extensions.

☐ ☐ ☐ 60°-90° Swing-Away Legrest, Short, Static Seat \$86 ea.
 Pair part: XRF000085. Left: XRF000086. Right: XRF000087
 For use only with 11"-13" lower extensions.

2. Lower Extension

☐ ☐ ☐ 11"-13" Length
 Pair part: XRF000088.
 Left: XRF000089.
 Right: XRF000090.

☐ ☐ ☐ 13"-16.5" Length
 Pair part: XRF000094.
 Left: XRF000095.
 Right: XRF000096.

☐ ☐ ☐ 15.5"-19.5" Length
 Pair part: XRF000097.
 Left: XRF000098.
 Right: XRF000099.

3b. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D¹ with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Amputee Supports with Swing-Away Upper

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

Therafin 31515 Multi-Axis Amputee Support \$230 ea.
 QTY Part: INDPART2822. HCPCS: E1020

Therafin 32387 Amputee Support Gel Cover \$155 ea.
 QTY Part: INDPART2823. HCPCS: K0108

16.1d. Elevating Legrests

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Elevating Legrests Upper

☐ ☐ ☐ Elevating Legrest Assembly \$171 ea.
 Pair part: FRMASMB6483. Left: RIG138012. Right: RIG138013. HCPCS: E0990
 Select lower extension and footplates

2a. Lower Extension

☐ ☐ ☐ 11"-14" Length
 Pair: FRMASMB4755. Left: RIG138026. Right: RIG138030
 For reinforced only

☐ ☐ ☐ 14"-17" Length
 Pair: FRMASMB4756. Left: RIG138033. Right: RIG138035
 For reinforced only

☐ ☐ ☐ 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length
 Pair part: XRF000101. Left: XRF000102. Right: XRF000103
 For angle adjustable only

☐ ☐ ☐ 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length
 Pair part: XRF000104. Left: XRF000105. Right: XRF000106
 For angle adjustable only

2b. Residual Limb Support

Do not select a footplate or accessory.

Left Residual Limb Support with Calf Pad..... \$342
 Part: FRMASMB5301. HCPCS: E1020

Right Residual Limb Support with Calf Pad \$342
 Part: FRMASMB5300. HCPCS: E1020

3a. Reinforced Cast Footplate(s)

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056

☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate Standard
 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062
 May cause interference on seat widths less than 18."

3b. Angle Adjustable Footplate(s)¹The 6"x11" footplates may cause interference with seat widths less than 18". ²The 8"x11" footplates may cause interference with seat widths less than 20".

☐ ☐ ☐ 5"Wx8"D with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D¹ with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D¹ with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D² with Heel Loops..... \$90 ea.
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799
 HCPCS: K0040+E0951

4. Elevating Legrest Accessories

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) \$131.25
 Part: INDPART2242. HCPCS: E0953

16.2. Legrest Accessories

Heel Loops (pair)..... \$53
 Part: FRMASMB7873. HCPCS: E0951
 Standard on angle adjustable and multi axis footplates when noted.

Wheel Bumper..... \$19 ea.
 Part: XRF000100

Legrests seating section continued on following page.

16. LEGRESTS CONTINUED**16.3. Set Legrest Extension Length****REQUIRED with center mounts and matching left and right rigging configurations. Part number is PTOINDV plus number following # below.****1a. Centermount Extension Length**

2" L (#3989)	4" L (#3991)	6" L (#3993)	8" L (#3680)	9" L (#3682)	10" L (#3684)	11" L (#3686)	12" L (#3688)	13" L (#3690)	14" L (#3692)	15" L (#3694)
3" L (#3990)	5" L (#3992)	7" L (#3994)	8.5" L (#3681)	9.5" L (#3683)	10.5" L (#3685)	11.5" L (#3687)	12.5" L (#3689)	13.5" L (#3691)	14.5" L (#3693)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)
7.5" L (#3471)	8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	

17. POWER ELECTRONICS**REQUIRED with actuator seating.**Single Actuator Through Controller and Interface \$2,163
Part: XRF000335. HCPCS: E2310Multiple Actuator Through Controller and Interface.....\$3,040
Part: XRF000336. HCPCS: E2311Single Actuator Through Controller and Interface (AAM).....\$3,040
Part: XRF000068. HCPCS: E2310**18. SEAT MOUNTED LED LIGHTS****Optional.** Requires Q-Logic 3 Joystick for Full Light Package or Q-Logic 3e Joystick.

Seat Mounted, LED, Full Light Package \$1,045

Part: XRF000130

Not available with Base Mounted, LED, Full Light Package.

19. ACCESSORIES**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.Cup Holder, Q-Captain No Charge
Part: ACC2006355Grip Drink Holder (20085) \$15.50
Part: ACC152993Front Mounted Cell Phone/ Media Holder, X-Grip, Q-Captain..... \$74
Part: ACC2103691Large, Black, Weather Cover \$202
Part: ACCCOVR1021Rear View Mirror, Q-Captain \$69
Part: ACC2103692Reflector Kit \$53
Part: ACC124670Oxygen Holder.....\$118
Part: ACC2006353. HCPCS: E2208
Not available with Power Lift option.Clothing Guard (Pair)..... \$237
Part: ACC125007
Not available with Heavy Duty armrests.Oxygen Holder..... \$180
Part: ACC2400456. HCPCS: E2208
Select with Power Lift option.Transfer Handles\$105 ea.
Part: MEC2205680. HCPCS: E0705
Not available with Heavy Duty armrests.Cane/Crutch Holder\$110
Part: ACC2006352. HCPCS: E2207
Not available with Power Lift option.Q-Logic EX/NE Programmer \$573
Part: ELEASMB6976Cane/Crutch Holder \$250
Part: ACC2400455. HCPCS: E2207
Select with Power Lift option.Q-Logic PC Programming Disk with Cable \$445
Part: ELEASMB5216Rear Basket\$118
Part: ACC2006351
Not available with Power Lift option.Privacy Bag\$84
Part: ACC2005723Rear Basket\$124
Part: ACC2400454
Select with Power Lift option.Personal Item Pouch, Rail Mounted.....\$84
Part: ACC2101017Saddle Bag\$40
Part: DWR1020D010Personal Item Pouch, Arm Mounted\$84
Part: ACC2103694Quantum Backpack \$168
Part: ACC1905853Glove Box, Rail Mounted \$158
Part: ACC125011Grip Lap Board (20005).....\$38.95
Part: ACC152992Glove Box, Arm Mounted..... \$158
Part: ACC2103693Remote Stop Switch \$1,295
Part: XRF000379**Aftermarket Accessory Mounting Bracket**

Select a quantity up to 5.

Aftermarket Accessory Mounting Bracket, Slotted \$53 ea.
Part: ACC2407046Aftermarket Accessory Mounting Bracket, Non-slotted .. \$53 ea.
Part: ACC2405496

20. BRANCH PATCHES**Optional.**

Air Force Patch\$21

Part: ACC140240

Coast Guard Patch\$21

Part: ACC140239

Navy Patch\$21

Part: ACC140244

Army Patch\$21

Part: ACC140241

Marines Patch\$21

Part: ACC140243

Purple Heart Patch\$21

Part: ACC140238

**21. POSITIONING COMPONENTS****Optional.** More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.**21.1. Stealth Unilink Pelvic/Thigh Guide**

Available with Flip-Up armrests only. Not available with Elevating Legrests.

1. UniLink Pelvic/Thigh Hardware

^ Please verify the appropriate billing code depending upon the placement of the support pad.

Unilink, 5" Height Adjustable, \$265 ea.

Removable Mounting Hardware

Part: ST-UTTP5-Q. HCPCS: E1028/E1034^

Flip Down with LEFT Full Surface Contact \$281 ea.

Feature (TWBADD-L & LPHW407)

Part: POS140428. HCPCS: E1028/E1034^

Unilink, 7" Height Adjustable, \$276 ea.

Removable Mounting Hardware

Part: ST-UTTP7-Q. HCPCS: E1028/E1034^

Flip Down with RIGHT Full Surface Contact..... \$281 ea.

Feature (TWBADD-R & LPHW407)

Part: POS140429. HCPCS: E1028/E1034^

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Unilink, 4" x 4" Pad \$106 ea.

Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^

Unilink, 4" x 12" Pad \$106 ea.

Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^

Unilink, 4" x 8" Pad..... \$106 ea.

Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^

3. UniLink Pelvic/Thigh Additional Links

Unilink, 1.5" Link \$81 ea.

Part: ST-UIL001150A000-Q

Unilink, 2" Link..... \$81 ea.

Part: ST-UIL001200A000-Q

4. UniLink Pelvic/Thigh Guide Installation

Install on Left

Part: PTOINDV3585

Install on Right

Part: PTOINDV3586

Do Not Install

Part: PTOINDV3719

21.2. Stealth Elbow Support**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube\$199 ea.

Style Arms (ARMS-ES-L)

Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket, Tube\$199 ea.

Style Arms (ARMS-ES-R)

Part: POS124066. HCPCS: E1028

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3) \$57 ea.

Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea.

Part: POS124067. HCPCS: K0108

21.3. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539

(SPLE107) (pair)

Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539

(SPLE109) (pair)

Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539

(SPLE108) (pair)

Part: INDPART2787. HCPCS: E0954

21.4. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.

Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) \$87 ea.

Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) \$87 ea.

Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

22. LAP TRAY

22.1. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q, HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309
Part: ST-SSLXL010001171-Q, HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309
Part: ST-SSLXM010001171-Q, HCPCS: E0950

2. Mounting Hardware

Swing Away Mounting Hardware..... \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	SM	ML	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	SM	SM	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	SM	SM	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	SM	SM	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						SM		ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	SM	SM	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q

■:Compatible

■:Armrest configuration may not hold tray position

*Some armrests listed may not be available on this order form.